

Instructions for use

EXPERTmatic contra-angle E31 C - REF 1.008.4998



KaVo. Dental Excellence.

Distributed by:

KaVo Dental GmbH
Bismarckring 39
D-88400 Biberach
Phone +49 7351 56-0
Fax +49 7351 56-1488

Manufacturer:

Kaltenbach & Voigt GmbH
Bismarckring 39
D-88400 Biberach
www.kavo.com



Table of contents

1	User instructions	4
2	Safety	6
2.1	Description of safety instructions	6
2.2	Safety instructions	8
3	Product description	12
3.1	Purpose – Intended use	13
3.2	Technical Specifications	15
3.3	Transportation and storage conditions	16
4	Start up and shut down	18
5	Operation	19
5.1	Attach the medical device	19
5.2	Remove the medical device	21

5.3	Insert the rubber cups or polishing brushes	21
5.4	Removing the rubber cup or polishing brushes	25
6	Preparation methods according to ISO 17664	26
6.1	Preparations at the site of use	26
6.2	Cleaning	27
6.2.1	Cleaning: Manual cleaning - external	27
6.2.2	Cleaning: Automated external cleaning	28
6.2.3	Cleaning: Manual cleaning - internal	29
6.2.4	Cleaning: Automated internal cleaning	31
6.3	Disinfection	32
6.3.1	Disinfection: Manual disinfection - external	33
6.3.2	Disinfection: Manual disinfection - internal	34
6.3.3	Disinfection: Machine disinfection - external and internal	35
6.4	Drying	36

6.5	Care products and systems - Servicing	37
6.5.1	Care products and systems - Servicing: Care with KaVo Spray	38
6.5.2	Care products and systems - Servicing: Care with the KaVo SPRAYrotor	39
6.5.3	Care products and systems - Servicing: Servicing with KaVo QUATTROcare 2104 / 2104A	40
6.5.4	Care products and systems - Servicing: Servicing with KaVo QUATTROcare PLUS	41
6.6	Packaging	42
6.7	Sterilisation	43
6.8	Storage	46
7	Tools	47
8	Warranty terms and conditions	49

1 User instructions

Dear User

Congratulations on purchasing this KaVo quality product. By following the instructions below you will be able to work smoothly, economically and safely.

© Copyright by KaVo Dental GmbH

Symbols

	Refer to the chapter on Safety/Warning symbol
	Important information for users and service technicians

	Action request
	CE mark (European Community). A product bearing this mark meets the requirements of the applicable EC directive.
	Can be steam-sterilised at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4 °F)
	Thermodisinfected

Target group

This document is intended for dentists and their assistants. The section on starting up is also intended for service technicians.

2 Safety

2.1 Description of safety instructions



Warning symbol

Structure



 **DANGER**

The introduction describes the type and source of the hazard.

This section describes the potential consequences of non-observance.

- ▶ The optional step includes necessary measures for hazard prevention.

Description of hazard levels

The safety instructions listed here, together with the three levels of danger will help avert property damage and injury.



CAUTION

CAUTION

indicates a hazardous situation that can cause damage to property or mild to moderate injuries.



WARNING

WARNING

indicates a hazardous situation that can cause death or serious injury.

**⚠ DANGER****DANGER**

indicates a hazardous situation that can directly cause death or serious injury.

2.2 Safety instructions

**⚠ WARNING****Hazard to the care provider and patient**

Damage, irregular noise during operation, excessive vibration, unusual build-up of heat or if the insert cannot be firmly held.

- ▶ Stop work and seek service support.

**⚠ CAUTION**

Burning hazard from hot instrument head and instruments cover.
If the instrument overheats, burns may arise in the oral area.

- ▶ Never contact soft tissue with the instrument head.

**⚠ CAUTION**

Premature wear and malfunctioning from improper storage during long periods of nonuse.
Reduced product life.

- ▶ The medical device should be cleaned, serviced and stored in a dry location, according to instructions, before long periods of non-use.

**Note**

The head housing is protected by special seals against the penetration of polishing paste.

**Note**

Before each treatment, insert a dental dam for safety reasons.

The following individuals are authorized to repair and service KaVo products:

- Technicians at KaVo branches throughout the world
- Technicians specially trained by KaVo

To ensure proper function, the medical device must be set up according to the reprocessing methods described in the KaVo Instructions for Use, and the care products and care systems described therein must be used.

KaVo recommends specifying a service interval at the dental office for a licensed shop to clean, service and check the functioning of the medical device. This service interval depends on the frequency of use and should be adjusted accordingly.

Service may only be carried out by KaVo-trained repair shops using original KaVo replacement parts.

3 Product description



EXPERTmatic contra angle E31 C (**Mat. no.1.008.4998**)

3.1 Purpose – Intended use

Purpose:

This medical device is

- intended for dental treatment only. Any other type of use or alteration to the product is impermissible and can be hazardous.
In connection with a polishing tool, the EXPERTmatic contra angle can be used to support the periodontal treatment by means of rotating polishing.
- A medical device according to relevant national statutory regulations.

Proper use:

According to these regulations, this medical device may only be used for the described application by a knowledgeable user. The following must be observed:

- the applicable health and safety regulations
- the applicable accident prevention regulations
- these instructions for use

According to these regulations, it is the responsibility of the user to:

- only use equipment that is operating correctly,
- use the equipment for the proper purpose,
- protect him or herself, the patient and third parties from danger, and
- avoid contamination from the product.

3.2 Technical Specifications

Drive speed	max. 10,000 rpm
identification	1 green ring
Transmission ratio	5.4 : 1



Note

Operate the product with clockwise rotation.

The rubber caps and polishing brushes or mandrels with snap-on function in accordance with DIN EN ISO 13295, type 5 mandrels available in the dentistry field can be applied.

The contra-angle handpiece can be mounted on all INTRAmatic motors, and motors with a connection in accordance with ISO 3964 / DIN 13940.

**Note**

Do not operate at maximum speed.
Adjust the speed via the motor control.

3.3 Transportation and storage conditions**⚠ CAUTION**

It is hazardous to start up the medical device after it has been stored strongly refrigerated.

This can cause the medical device to malfunction.

- ▶ Prior to start-up, very cold products must be heated to a temperature of 20°C to 25°C (68°F to 77°F).

	Temperature: -20°C to +70°C (-4°F to +158°F)
	Relative humidity: 5% RH to 95% RH absence of condensation
	Air pressure: 700 hPa to 1060 hPa (10 psi to 15 psi)
	Protect from moisture

4 Start up and shut down



WARNING

Hazard from nonsterile products.

Infection danger to the care provider and patient.

- ▶ Before first use and after each use, prepare and sterilise the medical device if needed.



WARNING

Disposal of the product in the appropriate manner.

Prior to disposal, the product must be appropriately prepared or sterilised if this is necessary.

5 Operation

5.1 Attach the medical device



WARNING

Release of the medical device during treatment.

A medical device that is not properly locked in place can release from the motor coupling and fall off.


- ▶ Carefully pull on it before each treatment to ensure that the medical device is securely locked on the motor coupling.

**⚠ CAUTION**

Removing and attaching the contra-angle handpiece while the drive motor is rotating.

Damage to the drive.

- ▶ Never attach or remove the contra-angle handpiece while the drive motor is rotating!



- ▶ Place the medical device on the (LUX) motor coupling and lock it into place.

- ▶ Before each treatment, pull on the medical device to see if it is securely seated in the motor coupling.

5.2 Remove the medical device

- ▶ Unlock the medical device from the motor coupling by twisting it slightly and then pulling it along its axis.

5.3 Insert the rubber cups or polishing brushes



Note

Only use rubber calyces and polishing brushes that are conformant with DIN EN ISO 13295.

**⚠ WARNING****Use of impermissible inserts.**

Injury to the patient or damage to the medical device.

- ▶ Observe the instructions for use, and use the insert properly.

**⚠ CAUTION****Danger of injury from insert.**


Infections or cuts.

- ▶ Wear gloves or fingerstalls.

**⚠ CAUTION****Hazard to patients**

A loosened insert can lead to injury.

- ▶ Check that the insert is firmly attached to the head by pulling on it.
- ▶ Block drive with holder..





- ▶ Push or screw the insert onto/into the head.
- ▶ Remove the holder.

5.4 Removing the rubber cup or polishing brushes



WARNING

Hazard from rotating insert.

Lacerations.

- ▶ Do not touch rotating insert!
 - ▶ Remove the insert from the contra-angle handpiece after treatment to avoid injury and infection during storage.
- ▶ After movement has stopped, unscrew or remove the insert.

6 Preparation methods according to ISO 17664

6.1 Preparations at the site of use



WARNING

Hazard from nonsterile products.

There is a risk of infection from contaminated medical devices.

- ▶ Take suitable personal protective measures.

- ▶ Remove all residual cement, composite or blood without delay.
- ▶ Recondition the medical device as soon as possible after treatment.
- ▶ Remove the insert from the medical device.
- ▶ The medical device must be dry when transported for reconditioning.
- ▶ Do not place it in a solution or similar.

6.2 Cleaning



⚠ CAUTION

Malfunctions from cleaning in the ultrasonic unit.


Defects in the product.

- ▶ Only clean manually or in a thermodisinfectator.

6.2.1 Cleaning: Manual cleaning - external

Accessories required:

- Tap water $30\text{ °C} \pm 5\text{ °C}$ ($86\text{ °F} \pm 10\text{ °F}$)
- Brush, e.g. medium-hard toothbrush



- ▶ Brush off under flowing tap water.

6.2.2 Cleaning: Automated external cleaning

KaVo recommends thermodisinfectors in accordance with EN ISO 15883-1, which are operated with alkaline cleaning agents with a pH value of max. 10 (e.g. Miele G 7781 / G 7881 – Validation was carried out with Programme "VARIO-TD", cleaning agent "neodisher® mediclean", neutralisation agent "neodisher® Z" and rinsing agent "neodisher® miel-clear" and only applies to the material compatibility with KaVo products).



- ▶ For program settings as well as cleansers and disinfectants to be used, please refer to the Instructions for Use of the thermodisinfectant (complying with max. pH value of 10).
- ▶ In order to prevent negative effects on the medical device, make sure that the interior and the exterior of the medical device are dry, and then lubricate immediately with care agents from the KaVo care system.

6.2.3 Cleaning: Manual cleaning - internal

Can only be done with KaVo CLEANspray or KaVo DRYspray.

- ▶ Cover the medical device with the KaVo CLEANpac bag, and place it on the corresponding care adapter. Press the spray button three times for 2 seconds each time. Remove the medical device from the spray attachment and let the cleaner work for one minute.

- ▶ Afterwards, rinse for 3-5 seconds with KaVo DRYspray.

See also: KaVo CLEANspray / KaVo DRYspray Instructions for Use



Note

KaVo CLEANspray and KaVo DRYspray for manual interior cleaning are only available in the following countries:

Germany, Austria, Switzerland, Italy, Spain, Portugal, France, Luxembourg, Belgium, Netherlands, United Kingdom, Denmark, Sweden, Finland and Norway.

In other countries interior cleaning can only be carried out with thermosinfectors in accordance with EN ISO 15883-1.



6.2.4 Cleaning: Automated internal cleaning

KaVo recommends thermodisinfectors in accordance with EN ISO 15883-1, which are operated with alkaline cleaning agents with a pH value of max. 10 (e.g. Miele G 7781 / G 7881 – Validation was carried out with Programme "VARIO-TD", cleaning agent "neodisher® mediclean", neutralisation agent "neodisher® Z" and rinsing agent "neodisher® miel-clear" and only applies to the material compatibility with KaVo products).

- ▶ For program settings as well as cleansers and disinfectants to be used, please refer to the Instructions for Use of the thermodisinfectant (complying with max. pH value of 10).

- ▶ In order to prevent negative effects on the medical device, make sure that the interior and the exterior of the medical device are dry, and then lubricate immediately with care agents from the KaVo care system.

6.3 Disinfection




CAUTION

Malfunctioning from using a disinfectant bath or disinfectant containing chlorine.

Defects in the product.

- ▶ Only disinfect in a thermodisinfectant or manually.

6.3.1 Disinfection: Manual disinfection - external



KaVo recommends the following products based on material compatibility. The microbiological efficacy must be ensured by the disinfectant manufacturer.

- Mikrozyd AF Liquid made by Schülke & Mayr
- FD 322 made by Dürr

- CaviCide made by Metrex

Consumables required:

- Cloths for wiping off the medical device.

- ▶ Spray the disinfectant on a cloth, then wipe down the medical device and allow the disinfectant to act according to the instructions of the disinfectant manufacturer.
- ▶ Follow the instructions for use of the disinfectant.

6.3.2 Disinfection: Manual disinfection - internal

The efficacy of manual internal disinfection must be demonstrated by the manufacturer of the disinfection agent. With KaVo products, use only disinfection agents that have been released by KaVo with respect to the compatibility of materials (e.g. WL-cid / made by ALPRO).

- ▶ Follow the instructions for use of the disinfectant.

- ▶ Immediately after internal disinfection, lubricate the KaVo medical device immediately with care agents from the KaVo care system.

6.3.3 Disinfection: Machine disinfection - external and internal



KaVo recommends thermodisinfectors in accordance with EN ISO 15883-1, which are operated with alkaline cleaning agents with a pH value of max. 10 (e.g. Miele G 7781 / G 7881 – Validation was carried out with Programme "VARIO-TD", cleaning agent "neodisher® mediclean", neutralisation agent "neodisher® Z" and rinsing agent "neodisher® miel-clear" and only applies to the material compatibility with KaVo products).

- ▶ For program settings as well as cleansers and disinfectants to be used, please refer to the Instructions for Use of the thermodisinfectant (complying with max. pH value of 10).

- ▶ In order to prevent negative effects on the medical device, make sure that the interior and the exterior of the medical device are dry, and then lubricate immediately with care agents from the KaVo care system.

6.4 Drying

Manual Drying

- ▶ Blow off the outside and inside with compressed air until water drops are no longer visible.

Automatic Drying

The drying procedure is normally part of the cleaning program of the thermodisinfectant.

- ▶ Follow the instructions for use of the thermodisinfectors.

6.5 Care products and systems - Servicing



CAUTION

Premature wear and malfunctions from improper servicing and care.
Reduced product life.

- ▶ Perform proper care regularly!



Note


Remove the insert for care.

**Note**

KaVo only guarantees that its products will function properly when the care products used are those listed as accessories, as they were tested for proper use on our products.


6.5.1 Care products and systems - Servicing: Care with KaVo Spray

KaVo recommends servicing the product after each time it is used, i.e. after each automatic cleaning and before each sterilisation.



- ▶ Remove the insert.
- ▶ Cover the product with the Cleanpac bag.
- ▶ Place the product on the cannula, and press the spray button for one second.

6.5.2 Care products and systems - Servicing: Care with the KaVo SPRAYrotor




KaVo recommends servicing the product after each time it is used, i.e. after each automatic cleaning and before each sterilisation.

- ▶ Place the product on the appropriate coupling of the **KaVo SPRAYrotor** and cover it with a CLEANpac bag.
- ▶ Service the product.

See also: Instructions for use **KaVo SPRAYrotor**

6.5.3 Care products and systems - Servicing: Servicing with KaVo QUATTROcare 2104 / 2104A

Cleaning and care unit with expansion pressure for effective cleaning and care.




KaVo recommends servicing the product after each time it is used, i.e. after each automatic cleaning and before each sterilisation.

- ▶ Removing the insert.
- ▶ Servicing the product.

See also: Instructions for use **KaVo QUATTROcare 2104 / 2104A**

6.5.4 Care products and systems - Servicing: Servicing with KaVo QUATTROcare PLUS



KaVo recommends servicing the product after each time it is used, i.e. after each automatic cleaning and before each sterilisation.

- ▶ Removing the insert.
- ▶ Servicing the product in QUATTROcare PLUS.

See also: Instructions for use **KaVo QUATTROcare PLUS**

6.6 Packaging



Note

The sterilisation bag must be large enough for the handpiece so that the bag is not stretched.

The quality and use of the sterilisation packaging must satisfy applicable standards and be suitable for the sterilisation procedure!

- ▶ Individually seal the medical device in the sterilised item packaging.

6.7 Sterilisation

Sterilisation in a steam steriliser (autoclave) in accordance with EN 13060 / ISO 17665-1 (e.g. KaVo STERlclave B 2200 / 2200 P)



CAUTION

**Premature wear and malfunctions from improper servicing and care.
Reduced product life.**

- ▶ Before each sterilisation cycle, service the medical device with KaVo care products.

**⚠ CAUTION**

Contact corrosion due to moisture.

Damage to product.

- ▶ Immediately remove the product from the steam steriliser after the sterilisation cycle!

135 °C



The KaVo medical device has a maximum temperature resistance up to 138 °C (280.4 °F).

(Depending on the available autoclave,) select a suitable procedure from the following sterilisation processes:

- Autoclave with three times initial vacuum:
 - at least 3 minutes at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4 °F)
 - Autoclave using the gravitation method:
 - at least 10 minutes at 134 °C -1 °C / +4 °C (273 °F -1.6 °F / +7.4 °F) or alternatively
 - at least 60 minutes at 121 °C -1 °C / +4 °C (250 °F -1.6 °F / +7.4 °F)
- ▶ Use according to the manufacturer's Instructions for Use.

6.8 Storage

- ▶ Reprocessed products should be stored protected from dust with minimum exposure to germs in a dry, dark and cool space.
- ▶ Comply with the expiry date of the sterilised items.

7 Tools

Available from dental suppliers.

Material summary	Mat. no.
Holder	1.004.0596
Instrument stand 2151	0.411.9501
CLEANpac 10 units	0.411.9691
Cellulose pad 100 units	0.411.9862

Material summary	Mat. no.
Adaptor INTRAmatic (CLEANspray and DRYspray)	1.007.1776
KaVo CLEANspray 2110 P	1.007.0579
KaVo DRYspray 2117 P	1.007.0580
KaVo Spray 2112 A	0.411.9640

Material summary	Mat. no.
ROTAspray 2 2142 A	0.411.7520
QUATTROcare plus Spray 2140 P	1.005.4525

8 Warranty terms and conditions

The following warranty conditions apply to this KaVo medical device:

KaVo provides the end customer with a warranty of proper function and guarantees zero defects in respect of material and processing for a period of 18 months from the date of the invoice, subject to the following conditions:

In case of justified complaints, KaVo will honour its warranty with a free replacement or repair. Other claims of any nature whatsoever, in particular with respect to compensation, are excluded. In the event of default, gross negligence or intent, this shall only apply in the absence of mandatory legal regulations to the contrary.

KaVo shall not be liable for defects and their consequences that have arisen or may arise from natural wear, improper handling, cleaning or maintenance, non-compliance with operating, maintenance or connection instructions, calcination or corrosion, contaminated air or water supplies

or chemical or electrical factors deemed abnormal or impermissible in accordance with KaVo's instructions for use or other manufacturer's instructions. The warranty granted does not usually extend to lamps, light conductors made of glass and glass fibres, glassware, rubber parts, and the colourfastness of plastic parts.

All liability is excluded if defects or their consequences originate from manipulations or changes to the product made by the customer or a third party that is not authorised by KaVo.

Warranty claims will only be accepted if the product is submitted along with proof of purchase in the form of a copy of the invoice or note of delivery. The dealer, purchase date, type, and serial number must be clearly evident from this document.

1.008.5577 · kb · 20120622 - 02 · en



KaVo. Dental Excellence.