



Transposagen

Better Models. Better Results.™

piggyBac™ Transposase

User Manual, Version 3.04.16

Transposagen Biopharmaceuticals, Inc.
535 W. Second Street, Suite L102. Lexington, KY 40508
info@transposagenbio.com

Contents

1.	Introduction and Background.....	3
A.	The <i>piggyBac</i> [™] Transposon/Tranposase Genetic Modification System	3
2.	Design and Construction of <i>piggyBac</i> [™] Transposon and Transposase Vectors	5
A.	The <i>piggyBac</i> [™] transposon and transposase vectors	5
B.	The <i>piggyBac</i> [™] selection cassettes	6
3.	Integration of the <i>piggyBac</i> [™] Transposon	6
A.	Basic Protocol	6
4.	Excision of <i>piggyBac</i> [™] Transposon Cargo	7
A.	Basic Protocol	7
5.	Confirmation of Integrated <i>piggyBac</i> [™] Transposon Cargo by Splinkerette PCR.....	8
6.	Frequently Asked Questions.....	8
7.	Troubleshooting.....	9
8.	References	9
9.	Limited Use	9



1. Introduction and Background

A. The piggyBac™ Transposon/Transposase Genetic Modification System

Transposagen's *piggyBac*™ (PB) transposon genetic modification system is a highly efficient, non-viral means of DNA integration into target genomes. To date, the *piggyBac*™ transposon has been utilized in gene therapy, regenerative medicine, cell line engineering, and animal model creation. The *piggyBac*™ genetic modification system enables you to:

- » Use simple transfection to alter a variety of animal genomes
- » Engineer cell lines for high level protein production
- » Complete non-viral gene delivery efficiently and economically
- » Revert modifications to the genome in a scarless, Footprint-Free™ manner

The *piggyBac*™ transposon is a mobile genetic element that efficiently transposes between vectors and chromosomes via a “cut-and-paste” mechanism (Figure 1). During transposition, the *piggyBac*™ transposase recognizes transposon-specific inverted terminal repeat sequences (ITRs) located on each end of the transposon. Subsequent to ITR recognition, the *piggyBac*™ transposase excises (“Cut”) the transposon from the transposon vector which is then efficiently integrated (“Paste”) into random TTAA genetic locations. The *piggyBac*™ transposons have an enormous cargo limit with over 200 kb demonstrated.




Figure 1: Mechanism of the *piggyBac*™ transposase/transposon gene modification system.

Genomes containing an inserted *piggyBac*™ cargo can be transiently re-transfected with the excision-only *piggyBac*™ (PBx) vector to remove the cargo in a scarless, or Footprint-Free™ manner (Figure 2). The powerful activity of the *piggyBac*™ transposon system enables any DNA cargo to be easily integrated into target genomes.




Figure 2: Scarless, Footprint-Free[™] removal of cargo with the excision-only *piggyBac*[™] vector (PBx)

Using Transposagen’s *piggyBac*[™] transposon selection technology in conjunction with site-specific nucleases, including CRISPR/Cas9 and TALEN, provides a clean and efficient method to select for gene editing down to a single base pair (Figure 3).




Figure 3: Editing DNA using Transposagen’s *piggyBac*[™] transposon selection technology using a site-specific nuclease to stimulate homologous recombination, creation of donor DNA with the edited gene and Multivector selection cassette, and Footprint-Free[™] removal of the cassette by excision-only *piggyBac*[™] vector (PBx)

▶▶ Learn more about Footprint-Free[™] Gene Editing at <http://www.transposagenbio.com/gek-footprint-free>





2. Design and Construction of *piggyBac*[™] Transposon and Transposase Vectors

Transposagen sells many *piggyBac*[™] transposon vectors that can be modified by customers to suit your research needs (shown below). In addition, we can provide fee-for-service cloning to create custom vectors with your gene cargo of interest (Cat # CPB-001). Contact us for instructions on how to order your custom *piggyBac*[™] vector. After submission of the DNA sequence, our scientists will design and synthesize your custom vector.

A. The *piggyBac*[™] transposon and transposase vectors

Transposagen has many validated, high purity *piggyBac*[™] transposon and transposase vectors available for research. These include, but are not limited to:

Catalog #	Description
<i>Transposon Vectors</i>	
SPB-007	This customizable <i>piggyBac</i> [™] transposon vector contains an insulator-flanked large multiple cloning site (MCS), allowing the flexibility to clone in your promoter and DNA cargo of choice. The core insulators protect cargo from genomic position effects and gene silencing. 10 µg size. [PB007]
SPB-008	Similar to SPB-007, but SPB-008 contains a GFP gene preceded by cloning sites that allow cloning of your promoter and gene of interest in frame to monitor expression. 10 µg size. [PB007G]
SPB-006	This customizable <i>piggyBac</i> [™] transposon vector contains a Polr2A promoter for ubiquitous gene expression. this promoter contains an intron and has demonstrated excellent expression levels <i>in vivo</i> and in transgenic animal models. 10 µg size. [PB001]
SPB-009	Similar to SPB-006, but SPB-009 provides the flexibility to clone any gene or cDNA with an in-frame GFP-fusion marker. This vector can also serve as a positive control for transposition. 10 µg size. [PB001G]
SPB-010	This customizable <i>piggyBac</i> [™] transposon vector contains the strong CAG promoter driving GFP expression. This vector also features an MCS to clone an expression cassette of your choice. The GFP marker will be driven independently from your cargo. 10 µg size. [PB002G]
SPB-011	Similar to SPB-010, but SPB-011 replaces the CAG promoter with a CMV promoter to express GFP. 10 µg size. [PB003G]
SPB-005	This customizable <i>piggyBac</i> [™] transposon vector contains an SV40 promoter-driven neomycin resistance gene. An insulator-flanked MCS allows for cloning in elements of interest, including user chosen genes driven by any selected promoter. The core insulators will protect cargo from genomic position effects and gene silencing. 10 µg size. [PBNeo]
<i>Transposase Vectors</i>	
SPB-001-10	Super <i>piggyBac</i> [™] transposase (sPBo) expression vector is a high purity, transfection ready, hyperactive <i>piggyBac</i> [™] transposase that facilitates stable integration of small and large genes into target genomes. This vector is non-renewable, and cannot be propagated in bacteria. Available in 10, 25, and 50 µg quantities. Also available as mRNA (Cat # SPB-003). [sPBo]
SPB-001-25	
SPB-001-50	
SPB-002-10	Excision-only <i>piggyBac</i> [™] transposase (PBx) expression vector is a high purity, transfection ready, transposase for excision of <i>piggyBac</i> [™] selection cassettes and transgenes for Footprint-Free [™] gene editing or phenotype reversions. This vector is non-renewable, and cannot be propagated in bacteria. Available in 10, 25, and 50 µg quantities. [PBx]
SPB-002-25	
SPB-002-50	

B. The *piggyBac*[™] selection cassettes

Transposagen has many validated *piggyBac*[™] Multivector plasmids that you can use to make your own targeting vectors. They include a choice of the most widely used antibiotic selection genes and/or GFP, and flexible cloning options for your arms of homology to facilitate clonal isolation after gene editing using a targeted donor vector insertion strategy. Flanking *piggyBac*[™] ITR excision sequences enable seamless cassette removal with Excision-only *piggyBac*[™] transposase (Cat # SPB-002). These selection cassettes include, but are not limited to:

<i>Available Multivector Options</i>		
Catalog #	Selection Cassette	Promoter
SGK-001 [MV-CMV-GFP-Puro]	GFP/Puromycin	CMV
SGK-002 [MV-PGK-GFP-Puro]	GFP/Puromycin	PGK
SGK-003 [MV-PGK-Neo]	Neomycin	PGK
SGK-004 [MV-PGK-Neo-TK]	Neomycin/Thymidine Kinase (TK)	PGK
SGK-005 [MV-PGK-Puro-TK]	Puromycin/TK	PGK
SGK-006 [MV-Hygro-TK]	Hygromycin/TK	PGK
SGK-007 [MV-EF1 α -Puro-TK]	Puromycin/TK	EF1 α

▶ Learn more about pricing and ordering at <http://www.transposagenbio.com>

3. Integration of the *piggyBac*[™] Transposon

Transposagen suggests following manufacturer's protocols for the introduction of DNA/RNA into your cell type of interest. The *piggyBac*[™] system is titratable. More or fewer integrations can be achieved by varying the amounts of transposase and transposon, as appropriate.

A. Basic Protocol

Note: Below is an example integration protocol using Transposagen's *piggyBac*[™] transposase vector and any *piggyBac*[™] transposon vector you choose. Every cell line will respond differently to transfection and the introduction of foreign DNA. Be sure to determine the optimal transfection efficiency in your cell type prior to beginning the integration protocol.

Co-transfect the Super *piggyBac*[™] transposase with the *piggyBac*[™] transposon vector:

1. Clone the desired cargo into the appropriate *piggyBac*[™] transposon vector.
2. Sequence verify the clones.
3. Grow target cells to 60–80% confluency.
4. Prepare the transfection mixture. For one well of a 6-well dish combine:





Amount	Item	Recommendation
1.0 µg	PB Transposon vector	Start with a 1:1 molar ratio of transposase to transposon to determine baseline integration efficiency.
1.0 µg	PB Transposase vector (or PB mRNA*)	Calculate the molar ratio using your final transposon size.
x µL	Transfection reagent of choice	Follow manufacturer's protocol for reagent amounts
50 µL	serum-free or reduced-serum DMEM	

* Note: If using pB mRNA, alter protocol to manufacturer's recommendations for RNA and/or consider electroporation.

- Mix by pipetting.
- Allow complex formation by incubating the mixture for 15–30 minutes at room temperature.
- Transfer the transfection mixture drop-wise to cells in culture wells and swirl to disperse.
- Incubate the cells and the transfection mixture at 37°C in a CO₂ incubator. If lipofection was used, change media after 24 hours.
- Check for positive integrations after 72 hours.
 - If you can, apply antibiotic selection or isolate cells by fluorescence-activated cell sorting (FACS), if applicable, to select for positive clones and measure expression of integrated cargo by standard RT-PCR methods.
 - If you did not use antibiotic or fluorescent markers, allow appropriate time for cells to remove trace amounts of remaining episomal transposon cargo before measuring expression of cargo by standard RT-PCR methods. Use non-transfected cells as negative control in RT-PCR.

4. Excision of *piggyBac*[™] Transposon Cargo

We recommend using excision-only *piggyBac*[™] (Cat # SPB-002) to remove inserted cargo (see Figure 2). Increasing amounts of PBx will increase excision efficiencies. Negative selection, such as thymidine kinase/ganciclovir or loss of a fluorescent marker, can also be utilized to enhance excision efficiencies.

A. Basic Protocol

- Note:** Below is an example excision protocol with Transposagen's Puromycin/Thymidine Kinase (Puromycin/TK) selection cassette. The protocol assumes that you have identified puromycin-resistant clones. Before beginning any selection protocol, you need to perform kill curves using your cell line and selection method of choice. Every cell line will respond differently to drug selection.
- Mix approximately 2 x 10⁶ cells with 10 µg excision-only *piggyBac*[™] plasmid (PBx) and transfect according to the manufacturer's protocol using either a nucleofection device or lipofection reagents.

Note: 5 µg is an appropriate starting point for induced pluripotent stem cells (iPS cells) and other difficult-to-transfect cells; for transformed cell lines use approximately 250 ng for 1 x 10⁶ cells.
 - Plate the cells onto a 6-well plate in 1:2, 1:4, and 1:6 dilutions in media without puromycin. If lipofection was used, change media after 24 hours.

3. On day 4, count cells and plate 1×10^4 cells onto a 10-cm dish in media containing ganciclovir (Sigma Cat # G2536) at a concentration determined from your previously performed kill curve. Ganciclovir should not affect cells without the thymidine kinase (TK) gene.
4. Change media every other day. Ganciclovir killing takes approximately 2–4 days for fast growing cells, but may take up to two weeks in some cell lines.
5. Check surviving clones for puromycin sensitivity and confirm appropriate cassette removal by standard molecular analysis methods.

Note: We recommend that you perform a literature search to find the appropriate conditions for your cell type prior to undertaking any gene editing project. Below are some published examples of the Footprint-Free[™] Gene Editing Kit system:

- ▶ Seamless genome editing in human pluripotent stem cells using custom endonuclease-based gene targeting and the *piggyBac*[™] transposon: <http://www.nature.com/nprot/journal/v8/n10/abs/nprot.2013.126.html> and <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC4559351/>
- ▶ Seamless gene correction of β -thalassemia mutations in patient-specific iPSCs using CRISPR/Cas9 and *piggyBac*: <http://genome.cshlp.org/content/early/2014/07/30/gr.173427.114>

5. Confirmation of Integrated *piggyBac*[™] Transposon Cargo by Splinkerette PCR

The Splinkerette PCR assay locates and maps the exact location of each transposition event in the genome, and provides a copy number for the number of transposon integrations. Transposagen offers a fee-for-service Splinkerette Assay, or an off-the-shelf Splinkerette kit (Cat # KPB-001) that includes the full protocol.

The full reference for the Splinkerette PCR procedure is published:

- ▶ A high-throughput splinkerette-PCR method for the isolation and sequencing for retroviral insertion sites, *Nature Protocols*, Col.4 No5. 2009, Page 789–798. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3627465/>

6. Frequently Asked Questions

Q. Is excision-only *piggyBac*[™] (PBx) mutagenic?

A. No, the excision-only *piggyBac*[™] is not mutagenic. It will only remove cargo flanked by the specific ITR sequences inserted previously. Removal of inserts by PBx is a scarless mechanism which restores the original wildtype sequence.

Q. How many copies of my integrated gene will be present?

A. The number of integration events will largely depend on the amounts and ratios of *piggyBac*[™] transposase and transposon. Titrating these amounts can yield of a majority of single integration events, or multiple events depending on your preference. Confirmation of integration events is also possible to verify the exact copy number present.

Q. How efficient are *piggyBac*[™] expression vectors?

A. While integration using Super *piggyBac*[™] transposase (sPBo) is highly efficient, excision with PBx is somewhat less so. To monitor excision, Transposagen recommends inclusion of the TK gene for counterselection of unexcised cells or use of a fluorescent marker.





7. Troubleshooting

If you see...	Then try this...
A lower-than-expected integration efficiency	Increase the amounts of transposase
A higher-than-expected integration efficiency	Decrease the amounts of transposase
Incomplete transposon excision across all cells	Multiple rounds of transfection may be required to completely remove all integrations.
Cell toxicity after transfection	Decrease DNA and transfection reagent amounts. Lower the amount of transposase used, as integration may have occurred in essential genomic locations.

8. References

1. Seamless genome editing in human pluripotent stem cells using custom endonuclease-based gene targeting and the piggyBac transposon. *Nat Protoc.* 2013; 8: 2061–2078. <http://www.nature.com/nprot/journal/v8/n10/abs/nprot.2013.126.html>
2. Seamless gene correction of β -thalassemia mutations in patient-specific iPSCs using CRISPR/Cas9 and piggyBac. *Genome Res.* 2014; 24: 1526–1533. <http://genome.cshlp.org/content/early/2014/07/30/gr.173427.114>
3. A high-throughput splinkerette-PCR method for the isolation and sequencing for retroviral insertion sites, *Nat Protoc.* 2009; 4(5): 789–798. <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3627465/>

9. Limited Use

Information in this document is subject to change without notice.

DISCLAIMER: TRANSPOSAGEN AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS PRODUCT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL TRANSPOSAGEN AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THE USE OF THIS PRODUCT.

LIMITED USE LABEL LICENSE FOR PIGGYBAC™ VECTORS AND TRANSPOSASE: The purchase of this product conveys to the buyer a limited, non-exclusive, non-transferable license under Transposagen Intellectual Property to use this product and components of this product only to perform internal research for the sole benefit of the buyer for a period of 12 (twelve) months from date of purchase. Following the 12 month trial period, purchaser will be required to obtain a nominal technology license in order to continue to use the piggyBac™ technology. If applicable, the buyer may also use standard molecular biology techniques to make additional copies of this product for purposes of internal research for the sole benefit of the buyer. Any buyer that is a for-profit entity may not use this product, its components or materials or cells made using this product or its components to generate any animal without first obtaining additional rights from Transposagen. The buyer cannot sell or otherwise transfer (a) this product, (b) its components, or (c) materials, cells, or organisms made using this product to a third party, or otherwise use this product, its components, or materials, cells, or organisms made using this product or its components for any Commercial Purpose or Development Purpose with the sole exception that buyer may transfer this product, its components, and/or materials, cells, or organisms made using this product or its components to (i) the buyer's legal affiliates, (ii) contract service provider and/or (iii) a scientific collaborator, provided that each such legal affiliate, contract service provider and/or scientific collaborator agrees in writing (1) not to transfer such product, its components, or materials, cells, or organisms made using such product or its components to any third party, and (2) to use such product, its components, and materials, cells, and organisms made using such product and/or its components solely for research as set forth in this limited use label license and not for Commercial Purposes or Development Purposes. If this product is subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control.

“Commercial Purpose” means any activity for consideration, including but not limited to: (a) any use, directly or indirectly, in manufacturing, production, production of biomolecules or quality control; (b) any use to provide a service, information, or data for consideration; (c) any use for therapeutic, diagnostic or prophylactic purposes; or (d) any sale, resale, leasing, or licensing, whether or not for research purposes.

“Development Purpose” means any: (a) clinical activity following IND-enabling preclinical toxicological studies or equivalents thereof; (b) activity in any agricultural field trial, including any such field trial that would be subject to regulation by the United States Department of Agriculture if the plants were to constitute genetically engineered organisms under 7 C.F.R. § 340 (or any successor regulation); (c) activity directed towards the submission of data to the United States Department of Agriculture, or any equivalent regulatory agency outside of the United States, in support of an application for clearance, approval or deregulation by such agency; (d) scale-up activities, the primary focus of which is to increase from small-scale to production-scale.

“Transposagen Intellectual Property” means: (a) certain patents and/or patents pending and (b) technical information, trade secrets, protocols, practices, methods, DNA sequences and other knowledge, information and skills controlled by Transposagen (“Transposagen Know-How”), that are necessary and/or useful in buyer’s use of this product.

For information on obtaining additional rights to the product for any use not permitted herein, please contact us at info@transposagenbio.com.

