

Thawing and Quality Control Protocol for Human Peripheral Blood Mononuclear Cells (PBMCs), Bone Marrow Mononuclear Cells (BMMCs), & Ascites Cells

PURPOSE:

This document establishes a standard protocol for quality control testing of our cryopreserved cell samples.

THAWING AND QUALITY CONTROL PROCEDURE:

All steps involving pipetting are to be executed in a Class II biosafety cabinet. Cells are counted using a Nexcelom Bioscience Cellometer Vision for quality control purposes.

A. Preparation of media, reagents and equipment

1. Using a micropipette, aliquot 20 μ L of AOPI (acridine orange/propidium iodide) solution into one 0.5 mL black, light-protected microcentrifuge tube for each sample to be counted. Store at 4°C until ready to use. (AOPI Staining Solution catalog number CS2-0106-5mL, Nexcelom Bioscience)
2. For post-thaw QC, prepare Dulbecco's Modified Eagle Medium (DMEM) + 10% Fetal Bovine Serum (FBS) solution as follows:
 - a. Mix the following:
 - 500mL DMEM (VWR catalog number 10-090-CV)
 - 50mL heat-inactivated FBS (VWR catalog number 45001-108)
 - 5.5mL pen/strep (pen: 10,000U/mL, strep: 10,000U/mL) (VWR catalog number 16777_164)
 - 5.5mL Fungizone (250 μ g/mL) (VWR catalog number 82026-728)
 - b. Aliquot DMEM + 10% FBS into 1.8mL aliquots in 15mL conical tubes to be used to count post-thaw cells. Store tubes at 4°C until ready for use.
3. 37°C water bath.

4. Dry ice in a plastic storage cooler.
- B. Viability and cell number measurement for post-thaw QC
1. Turn on water bath and set at 37°C.
 2. Pull QC vials from liquid nitrogen storage and keep them on dry ice until ready for use.
 - a. Samples must have been stored in the vapor phase of liquid nitrogen for a minimum of 24 hours before being pulled for QC purposes.
 - b. QC vials will be 200µL aliquots.
 3. Take out the number of 1.8mL aliquots of DMEM + 10% FBS needed to perform post-thaw QCs from the 4°C storage.
 - a. Place tubes in a tube rack and warm to 37°C.
 - b. Once warmed, remove from heat source and spray rack and tubes with 70% ethanol.
 - c. Place tubes and rack in hood.
 4. From the dry ice, remove a few of the QC vials and thaw the samples using the 37°C water bath. Work in small batches of no more than 4 or 5 aliquots at a time.
 - a. Quickly move the sample vials through the water to speed thawing.
 - b. Thaw until only a 2mm frozen crystal remains in the sample.
 5. For each 200µL QC aliquot, do the following:
 - a. Get one of the 1.8mL conical tubes and take ~300µL of DMEM + 10% FBS from the tube.
 - b. Add the ~300µL DMEM + 10% FBS to the cryovial of containing 200µL of cells. Constantly stir while adding the media.
 - c. Collect cells and media from the cryovial and slowly place back into the 15mL conical with constant stirring.
 - d. Rinse the QC vial with ~300µL from the 15mL conical and put back into the 15mL conical to be sure to get all of the cells from the QC cryovial.

- e. 200 μ L of cells + 1.8mL DMEM/FBS gives a 1:10 dilution of cells to media.
- 6. Add 20 μ L of thawed cell suspension to the AOPI solution to stain the cells. Mix well by gently pipetting up and down at least 3 times.
 - a. The final concentration of acridine orange will be 0.5 μ g/mL when mixed with the cells.
 - b. The final concentration of propidium iodide will be 10 μ g/mL when mixed with the cells.
- 7. Calculate the dilution factor.
 - a. The dilution factor is 1:2 of cells to fluorescent stain.
 - b. The dilution factor of the cell suspension is 1:10.
 - c. By multiplying the two dilutions, the final dilution of the stained cells is 1:20.
- 8. Peel the plastic film off of both sides of a Cellometer Counting Chamber (Nexcelom Bioscience catalog number SD100).
- 9. Load 20 μ L into one port of the counting chamber, and insert the loaded slide into the instrument and count.