

## Thawing and Quality Control Protocol for Dissociated Tumor Cells (DTCs)

### PURPOSE:

This document establishes a standard protocol for quality control testing of our cryopreserved cell samples.

### THAWING AND QUALITY CONTROL PROCEDURE:

All steps involving pipetting are to be executed in a Class II biosafety cabinet. Cells are counted using a Nexcelom Bioscience Cellometer Vision for quality control purposes.

#### A. Preparation of media, reagents and equipment

1. Using a micropipette, aliquot 20 $\mu$ L of AOPI (acridine orange/propidium iodide) solution into one 0.5 mL black, light-protected microcentrifuge tube for each sample to be counted. Store at 4°C until ready to use (AOPI Staining Solution catalog number CS2-0106-5mL, Nexcelom Bioscience).
2. For post-thaw QC, prepare Dulbecco's Modified Eagle Medium (DMEM) + 10% Fetal Bovine Serum (FBS) solution as follows:
  - a. Mix the following:
    - 500mL DMEM (VWR catalog number 10-090-CV)
    - 50mL heat-inactivated FBS (VWR catalog number 45001-108)
    - 5.5mL pen/strep (pen: 10,000U/mL, strep: 10,000U/mL) (VWR catalog number 16777\_164)
    - 5.5mL Fungizone (250 $\mu$ g/mL) (VWR catalog number 82026-728)
  - b. Aliquot DMEM + 10% FBS into 1.0mL aliquots in 15mL conical tubes to be used to count post-thaw cells. Store tubes at 4°C until ready for use.
3. 37°C water bath.

4. Dry ice in a plastic storage cooler.
- B. Viability and cell number measurement for post-thaw QC
1. Turn on water bath and set at 37°C.
  2. Pull QC vials from liquid nitrogen storage and keep them on dry ice until ready for use.
    - a. Samples must have been stored in the vapor phase of liquid nitrogen for a minimum of 24 hours before being pulled for QC purposes.
    - b. QC vials will be 1.0mL aliquots.
  3. Take out the number of 1.0mL aliquots of DMEM + 10% FBS needed to perform post-thaw QCs from the 4°C storage.
    - a. Place tubes in a tube rack and warm to 37°C.
    - b. Once warmed, remove from heat source and spray rack and tubes with 70% ethanol.
    - c. Place tubes and rack in hood.
  4. From the dry ice, remove a few of the QC vials and thaw the samples using the 37°C water bath. Work in small batches of no more than 4 or 5 aliquots at a time.
    - a. Quickly move the sample vials through the water to speed thawing.
    - b. Thaw until only a 2mm frozen crystal remains in the sample.
  5. For each 1.0mL QC aliquot, do the following:
    - a. Slowly add the 1.0mL DTC sample to the 1.0mL of pre-warmed media in the 15mL conical while constantly stirring.
    - b. Rinse the QC vial with ~300uL from the 15mL conical to be sure to get all of the cells from the QC cryovial.
    - c. 1.0mL of cells + 1.0mL DMEM/FBS gives a 1:2 dilution of cells to media.
  6. Add 20µL of thawed cell suspension to the AOPI solution to stain the cells. Mix well by gently pipetting up and down at least 3 times.

- a. The final concentration of acridine orange will be 0.5  $\mu\text{g/mL}$  when mixed with the cells.
  - b. The final concentration of propidium iodide will be 10  $\mu\text{g/mL}$  when mixed with the cells.
7. Calculate the dilution factor.
- a. The dilution factor is 1:2 of cells to fluorescent stain.
  - b. The dilution factor of the cell suspension is 1:2, since the 1.0mL cells were added to 1.0mL of warmed media.
  - c. By multiplying the two dilutions, the final dilution of the stained cells is 1:4.
8. Peel the plastic film off of both sides of a Cellometer Counting Chamber (Nexcelom Bioscience catalog number SD100).
9. Load 20 $\mu\text{L}$  into one port of the counting chamber, and insert the loaded slide into the instrument and count.