



Shear Bond Strength of Ceramic Primers with Lithium Disilicate

A. Heleba^{1*}, T. Hill², S. Singhal², P. McCabe² and G. W. Tysowsky²



¹Cornell University, Ithaca, NY, USA; ²Ivoclar Vivadent, Inc., Amherst, NY, USA.

Poster #1083



INTRODUCTION

Adhesive cementation plays an important part in the clinical success of esthetic restorations. The etching and priming step is a critical step during adhesive cementation which performs the following processes - creation of a micro-retentive surface and the activation of the surface for the subsequent chemical bonding procedure using a bonding agent. This bonding agent, a silane methacrylate, ensures a durable chemical bond between the glass-ceramic and the methacrylate-based luting composite. The longevity of the restoration can be dependent on the stability of this bond in an aqueous environment during stress. One method for testing the longevity of the bond is thermocycling.

OBJECTIVE

To compare the shear bond strength of resin cements with lithium-disilicate as a function of ceramic primers and artificial aging.

MATERIALS AND METHODS

Materials

	Monobond® Etch & Prime Ivoclar Vivadent Inc. Lot # U12508		IPS e.max® CAD Ivoclar Vivadent, Inc. Lot# 522131
	Scotchbond™ Universal / RelyX™ Ultimate 3M ESPE Lot # 577566 / 602839		Adhese® Universal / Variolink® Esthetic DC Ivoclar Vivadent Inc. Lot # SM0036 / T27196
	Optibond™ XTR / NX3 Nexus™ Kerr Lot # 602839 / 5314469		Tetric EvoCeram® Bulk Fill (IVA) Ivoclar Vivadent, Inc. Lot # U16982
	bluephase® G2 Ivoclar Vivadent Inc. Serial # 110003235		

Experimental Method:

A) Specimen Preparation




Eighty specimens of lithium-disilicate (IPS e.max CAD) were sectioned, crystallized and mounted. Specimens were polished through 400 grit SiC paper. Specimens (n=10) were distributed in eight groups outlined in the table below [Group 1A-4A (24h) and Group 1T-4T (Thermocycled)].

Experimental Groups: Stored 24h (A)/Thermocycled (T)

	Surface Treatment		Resin Cement
	Lithium-disilicate	Pre-cured Composite Rod	
Groups-1A/1T	Monobond Etch & Prime	Adhese Universal	Variolink Esthetic DC
Groups-2A/2T	Monobond Plus	Adhese Universal	Variolink Esthetic DC
Groups-3A/3T	HF (5%) 20 s	Scotchbond Universal	RelyX Ultimate
Groups-4A/4T	HF (5%) 20 s	Optibond XTR	Nexus 3


B) Resin Composite Rod (RCR) Preparation



Pre-cured resin composite rods (Diameter - 2.38 mm) were air-abraded (Silica/50µm/15psi), adhesive was applied and light cured (10 seconds).

MATERIALS AND METHODS, Cont.

C) Cementation and shear bond measurement



The RCR were cemented using adhesive resin cement under constant load (113.4 grams). Specimens were light-cured per manufacturer's instruction followed by storing for 24 hours at 37°C/100% humidity before shear-bond testing using universal testing device (Instron / crosshead speed-1.0 mm/min).

D) Thermocycling and shear bond measurement




The other half of the specimens for each treatment group were thermocycled in a SD Mechatronik Chewing Simulator from 5-55°C for 10,000 cycles using 40s dwell time before loading to failure (Instron / crosshead speed-1.0 mm/min).

RESULTS

Data was analyzed using one-way analysis of variance (ANOVA) and Tukey's post hoc analysis to determine statistical difference ($\alpha=0.05$) between or within each group.

Shear Bond Strength (24 hours/10K Thermocycling)



*Means with different letters are statistically different

- ❖ There was no significant difference in shear bond strength when using a new self-etch ceramic primer compared to Group 2.
- ❖ No significant differences in bond strength were found after thermocycling for Groups 1 and 2.
- ❖ Mean bond strength of Groups 3 and 4 decreased statistically after thermocycling.
- ❖ Group 4 had lower shear bond strengths before and after thermocycling than all other groups.

DISCUSSION

Within the bounds of this study that the single step etching and priming surface treatment in Group 1 produced similar results to a conventional technique of 5% HF and silane containing primer in Group 2. The use of universal types of adhesives as ceramic priming agents, Groups 3 and 4, appeared to have similar degradation of the bond to the ceramic surface over time to varying degrees in this study.

CONCLUSION

Within the limitation of this study, for 24 hour storage Groups 1A-3A produced statistically higher bond strengths than 4A. After thermocycling, Groups 3T and 4T produced statistically lower bond strengths than for 24h storage (3A and 4A).