HOW to STORE and MAINTAIN NOVASINA humidity standards

Novasina offers saturated analytical grade salt solutions, which can be used over years for verification and calibration of Novasina measurement instruments provided they are stored and maintained correctly. A humidity standard is fully saturated, when a visible solid remaining (accumulation of salt crystals) as well as a liquid phase e.g. as a salt solution film over the solid remaining is present. The only exceptions are SAL-T/SAL-SC 4%, 6% as well as SAL-T 11%, 33% and 58%. They might have a jellylike consistency (solid mixed with solution). Are these prerequisites met, a correctly generated humidity value can be assumed.

Storage :

it is recommended to store the salts in a cabinet as used for non-hazardous chemicals. The salts of different humidity values should not be kept in a common container, otherwise there will be a humidity exchange in order to reach equilibrium. Salts with values below the ambient air humidity, mainly 4%, 6%, 11% and 33% tend to absorb water. Salts above the ambient air humidity though, such as 84%, 90%, 97% and 98%, tend to lose water molecules and dry out. All salts are supplied in a well sealable container, thus keeping a humidity exchange with the ambient air as low as possible.

Note :

Shake the salts regularly (once a month) even if not used. Before using in a measurement instrument, shake the salt vigorously. After use, the salts should immediately be placed back into the container and closed and never let left on a table. This way the optimal lifetime of typically 2 to 5 years can be reached.

Tropical climate :

it is recommend to store low humidity salts in tightly screwed container and additionally in a well closed plastic bag.

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Regeneration of dry salts :

- in the case of the cylindric Sensor Checks SC a screw on the top can be opened and distilled water injected using a syringe until the liquid is visible over the solid remaining.
- for SAL-T no regeneration is intended (no opening available)

Both humidity standards type (SAL-SC & SAL-T) are easy to use, as there is no direct contact with the salt. Only the water molecules diffuse through the special membrane, but droplets are stopped. Low value salts, mainly 4%, 6%, 11% and 33% absorb more and more water which eventually dilutes the salt, until the solid remaining is dissolved. In this case, the humidity standard is useless and must be replaced before all crystals are dissolved.

For accurate results, the salt must appear like explained above. Furthermore it is recommended to shake the salts regularly (see chapter notes). This way an optimal life time and reproducible humidity values can be achieved.

Further information as well as the related matters safety data sheets (MSDS) can be downloaded from www.novasina.com

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