

Aging Concerns: New Solutions with Today's Technology

Ginger Walls

Brenlee Mogul-Rotman

permobil





1

Ginger Walls, PT, MS,
NCS, ATP/SMS

Brenlee Mogul-Rotman
OT Reg. (Ont.), ATP/SMS

permobil



2

Aging Concerns: New Solutions with Today's Technology

1

*Just because your patients
are aging, doesn't mean you
have to use old technology!*

—
People who are aging with a
disability have the right to
maximize their health and
independence using today's
technology!

permobil



3




The Aging Process

©Permobi

5

Aging

- Aging INTO disability
 - acquire conditions and limitations as the person grows old.
- Aging WITH a disability
 - Face unique healthcare needs and challenge conventional notions of 'healthy aging'



©Permobil

6

Healthy Aging and Age-Related Changes



©Permobil

- Movement Functions
- Sensory Function
- Skin and Tissues
- Cognitive Function
- Immune Function
- Sexuality


7

Multimorbidity: The presence of multiple chronic conditions at the same time


- Movement/musculoskeletal
- Sensory
- Cognitive
- Immune
- Skin
- Cardio-pulmonary

- CVA
- Falls/orthopedic issues
- COPD
- Parkinson's disease
- RA
- Diabetes
- Alzheimer's/dementia

Multiple diseases lead to disease interaction and polypharmacy, which leads to further disability




World Report of Ageing and Health, World Health Organization, 2015
https://apps.who.int/iris/bitstream/handle/10665/186463/9789240694811_eng.pdf?sequence=1




8


Healthy Ageing

- **Process of developing and maintaining the functional ability that enables well-being in older age**
- Functional ability-health related attributes and interaction of capacity and environment
- Intrinsic capacity-physical and mental capacities
- Environments-all factors in the extrinsic world
- Well Being- happiness, satisfaction, fulfillment





World Report of Ageing and Health, World Health Organization, 2015
https://apps.who.int/iris/bitstream/handle/10665/186463/9789240694811_eng.pdf?sequence=1



9

Aging Concerns: New Solutions with Today's Technology

4

Accelerated Aging

- Altered aging trajectory where the rate and the effects of aging are accelerated
- Health conditions occur *earlier and/or more frequently* than would otherwise be observed, leading to a narrow margin of health.
 - Due to physiologic changes due to SCI and impairments that lead to immediate and long term effects on the body
 - Susceptibility of those with SCI to numerous medical conditions that impart a health hazard
 - Disproportionate % of deaths as a result of preventable causes, including septicemia

Groah, et al. *Am. J. Phys. Med. Rehabil.* Vol 91, No 1, Jan 2012

permobil

©Permobil

10

Plan Ahead!

The diagram illustrates the interconnectedness of various factors in health and disability. At the top is 'Health Condition (disorder / disease)'. Below it are three boxes: 'Body Functions & Structure', 'Activity', and 'Participation', which are connected by double-headed horizontal arrows. Below these are 'Environmental factors' and 'Personal factors', also connected by a double-headed horizontal arrow. Vertical double-headed arrows connect 'Health Condition' to each of the three middle boxes, and each of the three middle boxes to the bottom horizontal arrow. Additionally, a vertical double-headed arrow connects 'Activity' to the bottom horizontal arrow.

Predictable issues with Body Function/Structure, MRADLs, & Participation are common to Progressive Disorders, Aging with a Disability AND Healthy Aging. PLAN AHEAD!

permobil

©Permobil

11

- “The consequences of poorly fitted equipment can cause an accelerated decline in function, pain and fatigue from poor posture, and increase in need for caregiver assistance or alternate mobility devices...”
- Requejo, P.S. et al. *Top Geriatric Rehabil* 2015;31(1):26-41

permobil

- Seating Considerations
- Manual Mobility Considerations
- Power Assist Considerations
- Power Mobility Considerations
- Training Considerations

©Permobil

12


Seating Considerations:
Posture, Skin and Tissue

permobil



©Permobil

13



14

Skin layer changes and related risks

Thinning of layers	Decreased vascularization	Loss of viscoelasticity	Nerve Ending changes	Decreased glands, fatty layer, connective tissue
<ul style="list-style-type: none">Slower healing time (30-50%)	<ul style="list-style-type: none">Decreased blood flow, O2 and nutrients to tissues	<ul style="list-style-type: none">Weaker tissue. Deep tissue at more risk for shear and deformation	<ul style="list-style-type: none">Less effectiveness to perceive input from pressure, heat and cold	<ul style="list-style-type: none">Brittle skin, looser skin, decreased thermoregulation

permobil

©Permobil

15

Incontinence and mobility

- “Impaired mobility can promote urge and mixed urinary incontinence.”
- Women with impaired mobility take longer to reach the toilet, increasing the risk of leakage when the urge to urinate is strong
- Women with urge incontinence may be more likely to limit their activities so that they are always near a toilet

Fritel X, Lachal L, Cassou B, Fauconnier A, Dargent-Molina P. Mobility impairment is associated with urge but not stress urinary incontinence in community-dwelling older women: results from the Ossébo study [published online ahead of print June 10, 2013]. BJOG. doi:10.1111/1471-0528.12316.

permobil

©Permobil

16

The diagram illustrates common areas of skin/tissue breakdown for two scenarios. The top part shows a person in a wheelchair with labels for OCCIPUT, THORACIC SPINE, HIPS/GREATER TROCHANTERS, SACRUM/COCCYX, and HEEL. The bottom part shows a person lying down with labels for THORACIC SPINE, HIPS/GREATER TROCHANTERS, SACRUM/COCCYX, SCAPULA, RIBS, and ISCHIAL TUBEROSITY.

Common Areas of Skin/Tissue Breakdown

- Wound development at any of these areas can be classified in one of the stages of pressure injuries as defined by the NPIAP

permobil

©Permobil

17

Shear and Friction in Bed

Friction = Force caused when two touching surfaces move in opposite directions

Shear effect = Skin stays in one place as the deep fascia & skeletal muscle slide down with gravity

permobil

©Permobil

18

Seating System impact on posture and function

Cushion+

Back Support +

Mobility Base

= Optimal Seating System

permobil

- Together they allow for maximum contact with the seat and back surface to:
 - Achieve postural **stability**
 - Allow for **pressure redistribution**
 - Attain postural **alignment**
 - Allow for **optimal mobility**

19

Seating Components

‘Clients area able to participate in more numerous and challenging activities if they have hands-free balance, feel stable in their seating and mobility base, and can reach beyond the limits of their arm length’

Lange, Minkel, 2018

permobil

Postural Stability and Alignment



20



Wheelchair Seating is more than just the cushion

permobil

©Permobil

21

Cushion Design, Material, Cover

Linear

Anti-Thrust

Wedge

Front of

22

Cushion features → function

Stability

Decreased shear, friction

Ability to reach

Ability to lean L/R, A/P


Less need to correct position during function

Ability to manage w/c skills

Design, medium, cover

©Permobil

23



24



25

Back Supports

Diagram illustrating the components of Back Support:

- Skin Integrity
- Aesthetics
- Function
- Stability
- Positioning
- Comfort, increased sitting tolerance

Back Support Features


- Angle
- Depth
- Height/Length
- Contour
- Medium

©Permobil

26



27



28

Secondary Supports

Four photographs showing different users in wheelchairs with secondary supports: an elderly man, a woman, a young boy, and a young man.

- Stability
- Shear/friction
- Ability to move to propel
- Ability to reach and lean

permobil

©Permobil

29

Manual Mobility

“independence in mobility is a key factor in life satisfaction ad contributes to maintaining the quality of life for the elderly and aging MWC user”

Requejo, R.S. et al, 2015




30



32

Benefits of optimally configured MWC



- Reduce pain level
- Increase activity level
- Extend number of active years in a MWC

If it's not equally configured, it's not equally effective!


permobil

©Permobil

33

Manual Wheelchair Frame Materials

STEEL ALUMINUM TITANIUM



DECREASES IN

- Weight
- Fatigue
- Corrosion
- Vibration

permobil

©Permobil

34


Clinical/Function Implications of Frame Materials

Chenier & Aissaoui, 2014

- Most WCs users experienced more than the acceptable standard of vibration

Wolf et al. 2007

- Whole body vibrations can cause low back pain, neck pain, fatigue, muscle ache



permobil

©Permobil

35

Rotational inertia is correlated with mass and w/c configuration

- Mass distribution changes
 - Heavier mag wheels, which focus mass toward outer edges of chair, increase the inertia and have an even greater impact than the accessories on rotational inertia
- Removing the accessories
 - Decreases mass of the chair
- Changing the axle position from rearward to forward
 - Caspall 2013





permobil

©Permobil

36

Struggling to move



permobil

©Permobil

37

...vs. ease of movement



permobil

©Permobil

38

What Impacts Rolling Resistance

- **Surface interaction**
 - Hard Smooth surface vs. Carpet
- **Material of tire**
 - Durometer (hardness of materials)
 - Solid
 - Pneumatic with airless inserts (foam)
 - Pneumatic – regular vs. high pressure
- **Tread Profile**
 - Lower tread profile = lower rolling resistance.
- **Camber** – increased camber decreases rolling resistance
- **Weight of Tire**
 - Increased material weight leads to increased deformation/strain
- **Deformation**
 - Increased deformation = higher rolling resistance
 - Deformation includes the weight of the client.
 - More mass...harder to push
 - Solid tires wear faster with load increases

permobil

©Permobil

39

Manual Wheelchair Configuration

- Seat depth/seat width
 - Frame depth
 - Center of Gravity- horizontal and vertical
 - Seat to floor height- front and rear (seat slope)
-
- “a wheelchair that is poorly matched to the individual, adversely affects potential activities and participation, lifestyle goals, health status and are costly.”
 - Requejo R.S. et al, 2015

permobil

©Permobil

40

“Participants who received individualized MWC found they had less difficulty independently propelling their wheelchairs and improved postural stability, which increase their ability to lean forward and reach”
Requejo, P.S. et al, 2015



permobil

©Permobil

41

Navigating the Barriers of the HealthCare System

—

permobil



42

So what do I do?

Plan of Care

- If the client will be d/c'd prior to receiving the ultra-lightweight w/c, the clinician must include in the Plan of Care for the client to follow up with the next clinician in the continuum and the supplier to obtain the ultra-lightweight MWC that is recommended

Educate your client

- Empower your client by explaining that they are going home in a rental MWC that will turn into a purchase in the next 12 months
- Encourage client to follow up with plan and the next therapist in the continuum about getting a better ultra-lightweight MWC


Document

- Rule out lesser wheelchairs as not effective. If it's not equally configured it's not equally effective!

45

Own the continuum!

Although shorter stays make it challenging to provide clients with the most appropriate wheelchairs, IT IS POSSIBLE when we learn the process and communicate across the continuum of care and empower our clients!

A photograph of a middle-aged man with grey hair, wearing a red long-sleeved shirt and dark pants, sitting in a black Permobil wheelchair. He is in a home environment with a stone fireplace wall behind him and a wooden floor.

permobil

46



47


6/30/98

Dear Ginger:

All day long
Day after day
I think of you
as I ride on my way
In my powerful little
one horse shay


And I thank you again
+ again for making
it possible.

Blessings + love
Shandra




permobil

Power Assist



The mobility challenges are the same.
The technology has come a long way!



©Permobil

48

Better solution with today's technology

A photograph of Alvin, a 90-year-old man with a white beard, sitting in a standard manual wheelchair. He is wearing a dark grey hoodie and dark pants. He is in a hallway with a wooden floor and a white wall with a handrail.


permobil

Meet Alvin

- 90-years-old
- Right AKA
- Limited upper extremity range of motion – so much that he was unable to independently propel MWC
- Anything more then a few feet, and he was dependent on caregiver
- He didn't want a powerchair because he had no way of transporting it

©Permobil

49

A photograph of Alvin, the same 90-year-old man, now sitting in a TiLite AeroX wheelchair with SmartDrive. He is in the same hallway, but the wheelchair is a different model, designed for independent use with a switch control system.

permobil

Meet Alvin

- TiLite AeroX with SmartDrive
- The switch control option has given him independent mobility!
- He was not able to use the SD with just the PushTracker.
- When asked what he liked most about the SD, he said, "I don't have to push the damn wheels anymore!"

©Permobil

50




51

When to consider power assist					
	Body Function & Structures	Activity	Participation	Environmental Factors	Personal Factors
	<div><div>Pain</div><div>Neuromusculoskeletal<ul style="list-style-type: none">• Range of motion• Strength• Muscle tone• Coordination</div><div>Cognitive Functions</div><div>Fatigue</div><div>Cardiovascular</div><div>Respiration</div><div>Postural Impairment</div><div>Shoulder joint integrity</div><div>permobil</div></div>	<div><div>Mobility<ul style="list-style-type: none">• Carrying or pushing items while propelling• Using transportation• Driving</div><div>Self-care</div><div>Domestic Life<ul style="list-style-type: none">• Shopping• Housework</div></div>	<div><div>Work</div><div>School</div><div>Recreation</div><div>Life roles</div><div>Relationships</div></div>	<div><div>Ramps</div><div>Hills</div><div>Distances</div><div>Varied terrain</div><div>Transportation</div><div>Support</div><div>Attitudes</div></div>	<div><div>Gender</div><div>Age</div><div>Fitness</div><div>Lifestyle</div><div>Life experiences....</div></div>

52

Questions to consider

	No	Yes
1. Will the client be a long-term manual wheelchair user?		
2. Does the client have a history of upper limb pain or dysfunction?		
3. Does the client experience fatigue that limits participation in daily life?		
4. Does the client experience changes in oxygen saturation with increased physical activity?		
5. Does the client require more than manual mobility for independence in all environments, but does not wish to pursue a power wheelchair?		
6. Does inefficiency in manual wheelchair propulsion inhibit participation in mobility related activities of daily living? Consider varying propulsion techniques including hemi-propulsion.		
7. Is the client dependent for mobility, and the caregiver is having a difficult time pushing the wheelchair in necessary environments/terrains?		
8. Is the client at risk for development of postural asymmetry over time, and does posture during propulsion increase this risk?		




permobil

©Permobil

53

Challenges inherent to hemi propulsion

- Many persons post-stroke utilize unilateral upper & lower extremity propulsion
- Kirby, et al (2005) – whether an individual had a neurologic impairment or not, for skills that had a higher rolling resistance or a high balance demand, the hemi propulsion style was difficult



permobil

©Permobil

54

Switch control for power assist

- Momentary switch control for foot propulsion or hemi propulsion
- Positive Impact on posture with effort of foot or hemi-propulsion




permobil

©Permobil

55

Manual tilt-in-space

- Caregivers can also benefit



permobil

©Permobil

56

Justification for Power Assist Use Objective & MRADL examples

[CLIENT] requires the use of SmartDrive power assist to accomplish the following Mobility Related Activities of Daily Living (MRADLs):

[SELECT THOSE THAT APPLY]: taking out the garbage, doing laundry, navigating small thresholds to enter the bathroom/garage/front door, ramp negotiation, childcare, housecleaning, meal preparation, propelling on carpet, etc.]

permobil

The image shows a screenshot of the SmartDrive mobile application. At the top, the status bar shows the time 2:58, LTE signal, and battery level. The app header is 'SmartDrive'. Below the header, there are two main sections: 'Clinical Tools' and 'Demos'. The 'Clinical Tools' section has a sub-header 'Training how to use SmartDrive' and two icons: 'Begin SmartDrive Evaluation' (a clipboard with a question mark) and 'SmartDrive Training' (a person at a computer). The 'Demos' section has a sub-header 'Demo location and management' and two images of the SmartDrive device. Below each image is a location: '216375 Orlando, Florida 32819, United States' and '220526 San Antonio, Texas 78217, United States'. The Permobil logo is in the bottom right corner.

57

Example – objective justification of power assist

With the use of the SmartDrive power assist, client is also able to complete MRADLs in a timely and efficient manner, as demonstrated by:

- Ex: to take laundry to the laundry room, it took 7 minutes and 25 pushes; however with the SmartDrive, it took 2 minutes and 4 pushes
- Ex: In addition, without the SmartDrive, client had to reposition the laundry basket 5 times and with the SmartDrive, there was no need to stop and reposition the basket.

permobil

58

BM1

Justification Tips

- Connect impairment in body functions and structures to inability to fully participate in MRADLs and life participation
- Include objective measures in justification:
 - Wheelchair Skills Test
 - Wheelchair Propulsion Test
 - Oxygen saturation
 - Fatigue Severity Scale
 - Pain scale



permobil

©Permobil

59

Power Wheelchair Mobility and Power Seat Functions

permobil



60

1968

Plumber, Al Theime invented the Amigo scooter in his garage in Michigan for a family member with MS



permobil


Nothing has changed much since... at least with POVs

- Limited seating options
- Stability issues
- Transfer challenges
- Table access issues
- Turning radius for in-home use?
- Ability to meet long term needs considering diagnosis/prognosis?

62

What about captain's seat power wheelchairs?

- Poor posture
- Decreased sitting tolerance
- Increased risk for Pain & PI
- Poor position for function
- Decreased independence
- Decreased participation
- Less time up in the chair



permobil

©Permobil

63



64

Actuators

SMART	NON-SMART
Position Feedback Sensors	Housing limit/Soft pontentiometer
Infinite position sensing	Binary position sensing
Allows for position limiting	Allows for position limiting
Allows for memory positioning	Cannot do memory positioning


©Permobi

65

Smart Actuator Functional Implications

- Memory seating
- Anterior sequence
- Standing sequence
- Independent repositioning sequence (IRM)
- Latched or Momentary
- Assign short cut key
- Ability to read & coach actuator position in real time

permobil



©Permobil

66



Sequential
35° of tilt
135° legrest extension
120° recline
(Kreutz,1997)

Muscle versus skin
reperfusion (Jan et al., 2013)

Significant muscle
reperfusion achieved by:
≥ 25° tilt and 120° recline;
*Unable to be achieved with
35° tilt alone*

©Permobi

67



Memory Positions

Save a frequently accessed position

Function

- MRADLs
- Functional Reach
- Transfers

Positioning

- Tone Management

Communication

Participation

Short Cut Keys!

©Permobil

68

Falls and Wheelchair Users

- Elderly w/c users in SNF fall frequently due to:
 - Their general impairment
 - Inability to compensate for environmental barriers
 - Unavailability of caregiver at the time of fall


Kiely, DK, Kiel, DP, Burrows, AB, Lipsitz, LA. (1998)
- Fall circumstances are multifactorial – action related (transferring); location related (bathroom); or fall attributions (surface conditions)
- Fall prevention strategies must be multifactorial – includes interdisciplinary assessment, transfer and driver skill education, seated postural control, appropriate match of technology solution

Rice, L, et al. (2017)

permobil

©Permobil

69

A woman with short dark hair, wearing a grey patterned top and a tan skirt, is seated in a wheelchair. She is leaning forward over a kitchen counter, focused on a task. The counter has various items on it, including a coffee maker and some containers. The background shows a typical kitchen setting with wooden cabinets.

Practice activities of daily living


PWC skills training can be provided to address some of the challenges new users face and to mitigate the potential risk to the user or others in his or her environment.

(Smith 2019, Canadian Journal of OT; Kirby, 1995)

©Permobil

70

Self-care – getting clothes in/out of the closet

A series of four small images showing a woman in a wheelchair performing self-care tasks. 1. She is sitting on the floor, leaning over to pick up a dark bag. 2. She is sitting in her wheelchair, looking towards a doorway. 3. She is in her wheelchair, reaching into a closet to handle clothes. 4. She is in her wheelchair, standing at the open doorway of a closet, looking inside. In all images, she is wearing the same grey patterned top and tan skirt. A green light is visible on her wheelchair in the third and fourth images.

©Permobi


71

Meet Anne

Self Care

- Sink access
- Brushing Teeth
- Washing face

permobil


A side-view photograph of Anne, a woman with short grey hair and glasses, wearing a blue sweater and a patterned scarf. She is seated in a black Permobil wheelchair, positioned in a kitchen. Her hands are resting on her lap, and she is looking towards the left side of the frame.

72



73

Smart Actuators
Program transfer position & assign short cut switch



permobil

74

Wheelchair Skills Training
for the Older Adult

permobil



75

Wheelchair Self Efficacy and Life Space Mobility

- Wheelchair skills accounted for 91.3% of the direct effect of self-efficacy on life-space mobility
- Self-efficacy is the belief that people have in their ability to use wheelchairs in various situations.

Sakakibara, et al 2015

```
graph TD; A[Lack of confidence] --> B[Decreased Function]; B --> C[Isolation]; C --> D[Less time using the chair]; D --> A;
```

permobil

©Permobil

76

Power wheelchair driving challenges – wheelchair user’s perspective

Difficulties with indoor mobility

- Moving in a constrained space
- Doorways
- Maneuvering around obstacles/people
- Carrying things while maneuvering

Maneuvering on streets, sidewalks

- Curb cuts, gaps, uneven surfaces
- Ramps
- Crowds
- Accessing transportation

Maneuvering for specific tasks

- Accessing sinks
- Activities of daily living

Managing unforeseen barriers

- Most difficulties reported were specific to the execution of a daily task and/or to the environmental context.

Torkia, C, et al 2015, Disabil Rehabil Assist Technol.


permobil

©Permobil

77

Training Considerations


CONTENT	TECHNIQUES	DOSE	MODIFICATIONS
WHAT TO TEACH Chair operation Basic driving skills Community driving skills Cognitive and perceptual skills related to driving Social skills related to driving Problem solving	HOW TO TEACH Demonstration Trial-and-Error Games Verbal and visual cueing Grading Chaining Limiting Chair Functions Feedback	HOW MUCH & HOW OFTEN Client-centred approach Depends on learning needs Several times per week Individual sessions Consider fatigue Time of day Sitting tolerance Length of total training Consider learning needs	ADDRESS TRAINING CHALLENGES Memory, cognition needs? Progressive condition? Correct drive configuration? Access method? Optimal positioning? Do they have anxiety? Programming needs?

(Smith, 2019)


©Permobil

78

Training for people with cognitive or memory impairment



Smith, 2019. Feasibility RCT protocol evaluating a powered-wheelchair training program for older adults. Canadian Journal of Occupational Therapy. Vol 86 (3) 232-242.



©Permobil

- PWC skills training based in trial-and-error methods may not provide the best learning opportunities for individuals with cognitive and memory impairments
- *Errorless learning and shared control strategies* may provide opportunities for safe and effective powered wheelchair skills training

79

permobil

PIDA ©2006

POWER-MOBILITY INDOOR DRIVING ASSESSMENT

MOBILITY DEVICE AND DRIVER EXPERIENCE CHECK LIST

Driver Name: _____

Assessment Date: _____

Assessor's Name: _____

TYPE OF MOBILITY DEVICE:

- Trial device: ☐ Client owned: ☐
- Make and model: _____
- Seating system: _____
- Type of controls: _____
- Special adaptations: _____

USE OF DEVICE SAFETY ACCESSORIES:

Please check all accessories that are currently used. Indicate those that are needed with a star.

_____ flag	_____ lights	_____ anti-tippers
_____ reflectors	_____ seat belt	_____ horn
_____ rear view mirror	_____ other (specify) _____	

Is the driver able to explain or demonstrate how each accessory is used: YES ☐ NO ☐

DRIVING EXPERIENCE:

- Years _____
- Devices used currently in the past _____
- Environment facility (briefly describe) _____

CAN THE CLIENT...

• Turn device on/off	YES	NO	N / A
• Utilize braking system	YES	NO	N / A
• Disengage braking system	YES	NO	N / A
• Use speed control switch	YES	NO	N / A
• Use special features of device	YES	NO	N / A
• Request assistance if necessary	YES	NO	N / A
• Independent transfer on/off	YES	NO	N / A

©Davison, Kaiserman, Chan & Gleason, 1995, 2006

80

permobil

©Permobil

81


Elderly and aging wheelchair users have increased risk for accelerated loss of function and mobility that greatly limits independence and affects quality of life.

Requejo, et al, 2015

permobil





82







Thank you for attending today's webinar!

Questions?

Please contact **education@permobil.com** with webinar questions and inquiries

permobil.com |    

83

References

Boland, P., Levack, W., Perry, M., & Graham, F. (2017). Equipment provision after stroke: A scoping review of the use of personal care and mobility aids in rehabilitation. *British Journal of Occupational Therapy*, 80(2), 73-88.

Caspall, J. J., Seligsohn, E., Dao, P. V., & Sprigle, S. (2013). Changes in inertia and effect on turning effort across different wheelchair configurations. *Journal of Rehabilitation Research and Development*, 50(10), 1353-1362. doi: 10.1682/jrrd.2012.12.0219.

Centers for Disease Control and Prevention (CDC). (2017). Stroke facts. Retrieved from <https://www.cdc.gov/stroke/facts.htm>

Charbonneau, R., Kirby, R. L., & Thompson, K. (2013). Manual wheelchair propulsion by people with hemiplegia: within-participant comparisons of forward versus backward techniques. *Archives of Physical Medicine and Rehabilitation*, 94(9), 1707-1713.

Dicanno, B., Liberman, J., Schmeler, M., Schuler, A., Cooper, R., Lagne, M., Liu, H., Jan, Y.(2015). *RESNA Position on the Application of Tilt, Recline, and Elevating Legrests for Wheelchairs: Literature Update*, [position paper]. RESNA. Retrieved September 11, 2015 from http://www.resna.org/sites/default/files/legacy/resources/position-papers/RESNA%20PP%20on%20Tilt%20Recline_2015.pdf



©Permobil

85

References


Kirby, R. L., Adams, C. D., MacPhee, A. H., Coolen, A. L., Harrison, E. R., Eskes, G. A., ... & Dupuis, D. J. (2005). Wheelchair-skill performance: controlled comparison between people with hemiplegia and able-bodied people simulating hemiplegia. *Archives of Physical Medicine and Rehabilitation*, 86(3), 387-393.

Kwakkel, G., & Kollen, B. J. (2013). Predicting activities after stroke: what is clinically relevant?. *International Journal of Stroke*, 8(1), 25-32.

Lange, M., Minkel, J. (2018) *Seating and Wheeled Mobility – A Clinical Resource*. Chapter 18, Considerations when working with the geriatric population, 297-316.

Lord, S. E., McPherson, K., McNaughton, H. K., Rochester, L., & Weatherall, M. (2004). Community ambulation after stroke: how important and obtainable is it and what measures appear predictive? *Archives of Physical Medicine and Rehabilitation*, 85(2), 234-239.

Mountain, A. D., Kirby, R. L., MacLeod, D. A., & Thompson, K. (2010). Rates and predictors of manual and powered wheelchair use for persons with stroke: a retrospective study in a Canadian rehabilitation center. *Archives of Physical Medicine and Rehabilitation*, 91(4), 639-643.



©Permobil

86

References

- Nilsson, L., Eklund, M., Nyberg, P., & Thulesius, H. (2011). Driving to learn in a powered wheelchair: the process of learning joystick use in people with profound cognitive disabilities. *American Journal of Occupational Therapy*, 65(6), 652-660.
- Punt, T. D., Kitadono, K., Hulleman, J., Humphreys, G. W., & Riddoch, M. J. (2011). Modulating wheelchair navigation in patients with spatial neglect. *Neuropsychological Rehabilitation*, 21(3), 367-382.
- Regier, A. D., Berryman, A., Hays, K., Smith, C., Staniszewski, K., & Gerber, D. (2014). Two approaches to manual wheelchair configuration and effects on function for individuals with acquired brain injury. *NeuroRehabilitation*, 35(3), 467-473.
- Rice, L. A., Yarnot, R., Mills, S., & Sonsoff, J. (2019). A pilot investigation of anterior tilt use among power wheelchair users. *Disability and Rehabilitation: Assistive Technology*, 1-8. doi: 10.1080/17483107.2019.1644676.
- Sakakibara, B., Miller, W., Eng, J., Backman, C., Routhier, F. (2014). Influences of Wheelchair-Related Efficacy on Life-Space Mobility in Adults Who Use a Wheelchair and Live in the Community. *Physical Therapy*. 94(11), 1604-1613.



©Permobil

87

References

- Serra, M. C., Treuth, M. S., Hafer-Macko, C. E., & Ryan, A. S. (2016). Increased Energy Cost of Mobility in Chronic Stroke. *Journal of Gerontology & Geriatric Research*, 5(6), 356.

Source: Online survey conducted in 2017 by PatientsLikeMe on 121 Power Wheelchair Users diagnosed with ALS, MS or SCI.

- Smith, E. M. (2019). *Powered wheelchair skills training for older adults with cognitive impairment: using shared control to facilitate independence* (T). University of British Columbia. Retrieved from <https://open.library.ubc.ca/collections/ubctheses/24/items/1.0378852>

- Smith EM, Miller WC, Mortenson WB, Mihailidis A. (2018) [*Feasibility Randomized Controlled Trial Protocol Evaluating an Errorless Learning Powered Wheelchair Skills Program for Older Adults with Cognitive Impairment*](#). Canadian Journal of Occupational Therapy. 86(3): 232-242.



©Permobil

88

References

Smith EM, Rismani S, Mortenson WB, Mihailidis A, Miller WC. (2019) *"A Chance to Try": Shared Control for Powered Wheelchair Assessment and Training*. American Journal of Occupational Therapy. [Accepted]

Van de Port, I. G., Kwakkel, G., Schepers, V. P., & Lindeman, E. (2006). Predicting mobility outcome one year after stroke: a prospective cohort study. *Journal of Rehabilitation Medicine*, 38(4), 218-223.

Woodman, P., Riazi, A., Pereira, C., & Jones, F. (2014). Social participation post stroke: a meta-ethnographic review of the experiences and views of community-dwelling stroke survivors. *Disability and Rehabilitation*, 36(24), 2031-2043.

permobil

©Permobil