SILVERLON

safe, strong, simple

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
Silverlon® Lifesaver® IV/Catheter Dressings

**EC REP**

**Emergo Europe,
Prinsessegracht 20,
2514 AP The Hague,
Netherlands
Tel.: +(31) (0) 70.345.8570**

LC-CE-IFU-IVCD ; Rev 06

Argentum Medical
2571 Kaneville Court
Geneva, IL 60134 USA
Tel. (+1) 888.551.0188
Fax. (+1) 888.551.0188
www.silverlon.com

SILVERLON

safe, strong, simple

Silverlon® Lifesaver® IV/Catheter Dressings

English

Device Description

Silverlon® Lifesaver® IV/ Catheter Dressings are sterile, multi-layer, non-adherent, absorbent, antimicrobial barrier dressings, intended for topical use in wound management, to reduce risk of wound infection.

Available as:

Code	Pad Size
Lifesaver™ Ag IV/CD-15	2,5 cm circle with 1.5 mm center hole
Lifesaver™ Ag IV/CD-40	2,5 cm circle with 4.0 mm center hole
Lifesaver™ Ag IV/CD-70	2,5 cm circle with 7.0 mm center hole

Silverlon® Lifesaver® IV/ Catheter Dressings contain silver to provide an ancillary, local antimicrobial effect to reduce the risk of wound infection by organisms sensitive to silver. Ionic silver is a broad spectrum antimicrobial, which has been shown to be effective against a range of pathogens associated with wound infection. Silverlon® dressings have been tested in vitro and found effective against microorganisms such as: Enterococcus faecalis, Staphylococcus aureus, Pseudomonas aeruginosa, Acinetobacter Baumannii, Enterobacter cloacae, Staphylococcus epidermidis, and Klebsiella pneumoniae.

In vitro testing has shown Silverlon® dressings to be MRI compatible.

Indications

Silverlon® Lifesaver® IV/ Catheter Dressings are professional use wound dressings indicated for local management of:

- Vascular access sites (including around IV catheters, central venous lines, arterial catheters, dialysis catheters);
- Drainage tube access sites;
- Epidural catheter access sites;
- Infected Wounds*

*Silverlon® Lifesaver® IV/Catheter Dressings may be used on those infected wound types listed above. Where used on infected wounds, the infection should be treated as per local clinical protocol.

Silverlon® Lifesaver® IV/ Catheter Dressings are indicated for a total duration of contact (i.e. where treatment involves consecutive application of individual dressings) of up to 30 days.

Warnings

- Do not use past expiration date on the product packaging.
- Do not use if pouch is damaged or open.
- Clinicians / Healthcare Professionals should be aware that there are very limited data on the prolonged and repeated use of silver-containing dressings, particularly in children and neonates.
- Silverlon® dressings should not be used during pregnancy or lactation because of the lack of data concerning the effects of silver on reproduction.

- Do not use petroleum-based ointments or creams under Silverlon® Lifesaver® IV/ Catheter Dressing.
- Do not soak Silverlon® Lifesaver® IV/Catheter Dressing in hydrogen peroxide, or povidone iodine.
- Silverlon® Lifesaver® IV/Catheter Dressings are intended for single use only and should not be re-used. Re-use may adversely affect the dressing's performance characteristics and also presents a risk of infection to the patient.
- Total duration of use of Silverlon® Lifesaver® IV/Catheter Dressing (i.e. where treatment involves consecutive application of individual dressings) should not exceed 30 days.

Contraindications

Avoid using Silverlon® Lifesaver® IV/Catheter Dressings on patients with known sensitivity to silver or nylon.

Adverse Reactions

Silverlon® Lifesaver® IV/ Catheter Dressings have been subjected to independent standard in vitro and in vivo biocompatibility tests, including cytotoxicity, sensitization and intracutaneous reactivity. All tests were performed in accordance with the International Standard Organization (ISO) 10993 Standard Series for Biological Evaluation of Medical Devices. The results of these studies indicated that Silverlon® Lifesaver® IV/ Catheter Dressings are safe for their intended use.

Instructions for Use

- Prepare the skin around wound site according to hospital protocol.
- Select the Silverlon® Lifesaver® IV/ Catheter Dressing with a hole size that is slightly larger in diameter than the catheter (or other device) protruding from the wound.
- Activate Silverlon® Lifesaver® IV/ Catheter Dressing by wetting it with sterile water, distilled water, or normal saline.
- Place the Silverlon Lifesaver® IV/ Catheter Dressing around the device protruding from the wound, so that the silver side of the dressing is in direct contact with the skin.
- Secure Silverlon® Lifesaver® IV/ Catheter Dressing in position with a transparent film dressing.
- Silverlon® Lifesaver® IV/ Catheter Dressing may be used for up to 7 days, but may require more frequent changing depending on wound condition and exudate buildup.
- To remove Silverlon® Lifesaver® IV/ Catheter Dressing, first remove the outer dressing per local protocol, then gently lift the edges of the Silverlon® Lifesaver® IV/ Catheter Dressing.
 - If the dressing sticks to the skin, moisten the dressing as needed with sterile water, distilled water, or normal saline, until it can be easily removed by gently lifting the corners.

SILVERLON [®]	Compressas Silverlon® Lifesaver® para IV/Cateter	Português	SILVERLON [®]	Silverlon® Lifesaver® IV-/katheterverbanden	Nederlands
safe, strong, simple			safe, strong, simple		

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