

# RET Fusion Factsheet



**CrownBio**  
CONNECTING SCIENCE TO PATIENTS

## Accelerate RET inhibitor research with novel RET fusion PDX models and derived early drug development resources

Constitutively active RET fusion oncogenes are important targets for oncology research. RET is a proto-oncogene that encodes a receptor tyrosine-protein kinase involved in numerous cellular mechanisms including cell proliferation, migration, and differentiation. Fusions of RET with different nearby genes on chromosome 10 (at 10q11.2) result in a constitutively active gene – permanently switched on to drive cell growth and tumor development. RET fusion oncogenes are therefore highly ‘druggable’ and important targets for cancer research.

Historically, RET fusions have been found in 10% to 20% of sporadic papillary thyroid cancers<sup>(1)</sup>; however, more recently a range of RET fusions were found in approximately 1% of lung adenocarcinomas<sup>(2,3)</sup>. This incidence can increase to as much as 15% in never-smokers with lung adenocarcinomas that lack other known driver oncogenes<sup>(4)</sup>. RET fusions (specifically CCDC6-RET and NCOA4-RET) have also been identified in patients with advanced colorectal cancer (CRC). Within CRC, RET fusions represent a novel class of oncogenic drivers, which occur at 0.2% frequency and without other known driver mutations e.g. KRAS, NRAS, BRAF, PIK3CA, or other fusion tyrosine kinases<sup>(5)</sup>.

Currently, a specific RET-only inhibitor is not available. Multi-targeted tyrosine kinase inhibitors which are known to have anti-RET activity are already approved for use in medullary thyroid cancer (e.g. cabozantinib and vandetanib), certain leukemias (e.g. ponatinib), and renal cell carcinoma (e.g. sunitinib and sorafenib). Many of these agents are now in clinical trials for cancer types including RET-rearranged NSCLC<sup>(6)</sup>. Targeted therapies against RET fusions may provide novel treatments for a patient subset lacking other treatment options, particularly in advanced CRC where there is an unmet need for molecularly directed therapies<sup>(5)</sup>.

### CrownBio RET Fusion Resources

CrownBio has a range of xenograft and cell line RET fusion platforms, to provide a critical link between *in vitro* testing and clinical efficacy in drug development.

CrownBio has the largest commercially available collection of over 2,500 patient-derived xenograft (PDX) models. The **HuPrime**® PDX collection closely reflects patient tumor histopathological profiles and is supported by a suite of curated patient and sample data including clinical diagnosis, patient clinical history, histopathology, and mutational status. Available gene expression data (measured by RNA sequencing or using Affymetrix® Human Genome U219 Array Plate), copy number (measured using Affymetrix Genome-Wide Human SNP Array 6.0), and hotspot mutation analysis for 10 genes for our **HuPrime** PDX models can be found in **HuBase**™, our free to access, curated online PDX database. **HuBase** can be accessed directly from the homepage of our website at [www.crownbio.com](http://www.crownbio.com). Within the **HuPrime** collection, we have identified 3 unique colorectal adenocarcinoma PDX models containing RET fusions: CR1520 with NCOA4-RET fusion, CR2518 with CCDC6-RET fusion (which within our models were

the first reported oncogenic RET fusions in CRC patient samples)<sup>(7)</sup>, and CR2545 which is a RET inhibitor resistant derivative of CR2518.

A challenge encountered by performing PDX *in vivo* studies is the lack of a corresponding *in vitro* cell culture system for cost-effective and high-throughput drug screening and model selection. CrownBio has therefore established our *in vitro* **PrimePanel**™ cell line collection<sup>(8)</sup>, derived from mouse stromal cell-depleted primary cancer cell cultures from our PDX tumors. The cells are all early passage (<10) and maintain essential histopathological features and genetic profiles of the original patient tumors including genomic mutational status, biochemical signaling, and response to tumor cell autonomously targeted therapeutics. Therefore, **PrimePanel** provides a higher throughput, faster turnaround, more versatile, and lower cost platform to assess drug efficacy prior to *in vivo* study, to fulfill critical research needs in early stage drug discovery.

CrownBio has derived a **PrimePanel** cell line from the CR1520 NCOA4-RET fusion PDX model, with both 2D and 3D *in vitro* assays validated for the cell line. An associated **PrimeXeno**™ model has also been developed. **PrimeXeno** *in vivo* models are established from our **PrimePanel** cell lines and are ideal for early stage drug discovery, where a robust system is needed to screen large numbers of compounds in an assay such as PK/PD analysis.

### HuPrime RET Fusion PDX Models: CR1520 and CR2518 Model Background

The available patient background information, tumor pathology diagnosis, and completed sequencing data for CR1520 and CR2518 are shown in **Table 1**. The mutation status of 10 common oncogenes and tumor suppressor genes has also been verified, and the genes that have been evaluated for each specific model are included in **Table 1**. Both PDX are negative for hotspot mutations in KRAS, BRAF, and PI3K (CR2518 has a mutation in CTNNB1 identified by RNAseq), suggesting that RET fusion may be the primary oncogenic driver in these tumors<sup>(7)</sup>. The pathology diagnosis for each tumor was QC'd by CrownBio, and also established for the PDX model at a range of low passage numbers. QC data are included in **Table 1** and representative H&E stain images (400x) are included within **Table 2**; further images from other PDX passage numbers are available within **HuBase**.







# RET Fusion Factsheet

**Table 1: Summary of HuPrime RET Fusion PDX Models Background**

HuPrime Identifier	Patient Background	Tumor Pathology Diagnosis	PDX Tumor Pathology QC	Genomic Profiling	IHC	Treatment History	Examples of Oncogene Mutation Status
CR1520 (NCOA4-RET)	Asian male, aged 82 years	Moderately-poorly differentiated adenocarcinoma	Pa, P2: Moderately-poorly differentiated adenocarcinoma	Affy U219 and SNP 6.0 data (P2) RNAseq (P3)	NA	Naive	WT: AKT, BRAF, c-MET, CTNNB1, EGFR, ERK, KRAS, PIK3CA, PTEN, TP53
CR2518 (CCDC6-RET)	Asian female, aged 82 years	Adenocarcinoma of ascending colon, ulcerative type, poorly differentiated, invaded through intestinal wall to peripheric adipose tissue. No malignant cells adjacent to superior and inferior stump. Regional lymph nodes (LN): paraintestinal LN (0/17)	Pa, P3: Poorly differentiated adenocarcinoma	Affy U219 and SNP 6.0 data (P3) RNAseq (P2)	CK7(+/-) CK20(++) Syn(+++) CD56(-) Ki-67(60%+) CerbB-2(-) VEGF(-) D2-40 (lymphatic vessel +)	Naive	WT: AKT, BRAF, EGFR, ERK, KRAS, PIK3CA, PTEN, TP53  Mutation: CTNNB1 Arg565His

**Table 2: HuPrime RET Fusion PDX Models Pathology**

H&E stained pathology images (400x) from the original patient tumor, and from a low passage number PDX model.

HuPrime Identifier	CR1520	CR2518
Patient Tumor		
PDX		

## HuPrime RET Fusion PDX Model Expression Levels

The expression levels of RET in PDX models CR1520 and CR2518 are shown in **Figure 1**, in comparison with all HuPrime CRC models, and just showing CRC models with expression levels >3. **Figure 1** shows Affymetrix Human Genome U219 Array Plate expression data, RNAseq expression data is also available in [HuBase](#). Both models show high expression of RET, as does model CR2545 which is derived from model CR2518.


## HuPrime RET Fusion PDX Model Genomic Profiling

RET fusions within the PDX models were identified by RNA sequencing (full genomic and sequencing data for the models and the fusions are available within [HuBase](#)). The CR1520 model has an in-frame NCOA4-RET fusion at exon 9 of the NCOA4 gene and exon 12 of the RET gene, detailed in **Figure 2**. RT-PCR analysis of the patient tumor and passage 3 of CR1520, and direct sequencing of the model confirmed the fusion protein was present (**Figure 3**). The CR2518 model has an in-frame CCDC6-RET fusion at exon 8 of the CCDC6 gene and exon 12 of the RET gene, detailed in **Figure 4**. RT-PCR analysis on passage 2 and 5 of CR2518, and direct sequencing of the model confirmed the fusion protein was present (**Figure 5**).

# RET Fusion Factsheet


**Figure 1: HuPrime RET Fusion PDX Model Expression Data**

Gene expression levels of RET (using Affymetrix Human Genome U219 Array Plate, probeset 5979\_at) in all HuPrime CRC PDX models, and CRC models with expression levels >3 (overlaid panel).



**Figure 2: HuPrime RET Fusion PDX Model Genomic Profiling: CR1520 Gene Fusion**


The chromosome 10 position of the NCOA4-RET fusion and associated fusion protein sequence.



PDX model	Up-stream Fusion Gene	Up-stream Chromosome	Up-stream Strand	Upstream Genome Junction Position	Down-stream Fusion Gene	Down-stream Chromosome	Down-stream Strand	Downstream Genome Junction Position
CR1520	NCOA4	10	+	51586411	RET	10	+	43612032


**Figure 3: HuPrime RET Fusion PDX Model Genomic Profiling: CR1520 RT-PCR Analysis and Sequencing**

RT-PCR and direct sequencing of CR1520 confirms the presence of the NCOA4-RET fusion protein.



**Figure 4: HuPrime RET Fusion PDX Model Genomic Profiling: CR2518 Gene Fusion**

The chromosome 10 position of the CCDC6-RET fusion and associated fusion protein sequence.




PDX model	Up-stream Fusion Gene	Up-stream Chromosome	Up-stream Strand	Upstream Genome Junction Position	Down-stream Fusion Gene	Down-stream Chromosome	Down-stream Strand	Downstream Genome Junction Position
CR2518	CCDC6	10	-	61554231	RET	10	+	43612032

# RET Fusion Factsheet

**Figure 5: HuPrime RET Fusion PDX Model Genomic Profiling: CR2518 RT-PCR Analysis and Sequencing**


RT-PCR and direct sequencing of CR2518 confirms the presence of the CCDC6-RET fusion protein.




## HuPrime RET Fusion PDX Model Treatment Data

Treatment data with 5-FU for CR1520, and 5-FU and irinotecan for CR2518 is stored within HuBase. The CR1520 model is sensitive to 5-FU, whereas CR2518 is only partially responsive to 5-FU and sensitive to irinotecan. The PDX models have also been treated with a RET inhibitor, with both CR1520 and CR2518 shown to be sensitive to the agent (Figure 6).

**Figure 6: HuPrime RET Fusion PDX Models CR1520 and CR2518 are Sensitive to RET Inhibitor Treatment** A: CR1520 B: CR2518.




**Figure 6: Continuation.**



## Development of a RET Fusion PDX Model Resistant to RET Inhibition

As shown above, model CR2518 is responsive to RET inhibition. However, following consecutive cycles of treatment with the RET inhibitor AP24534, resistance to the agent is developed and the RET inhibitor resistant PDX variant CR2545 is produced (Figure 7). We have validated the model to confirm that the CCDC6-RET fusion is still present. RNAseq data revealed that the acquired resistance is due to a new mutation (V804M) at the RET kinase domain, that is known to cause significant changes in the binding domain for the inhibitor.

**Figure 7: HuPrime RET Fusion PDX Model CR2545 is Resistant to RET Inhibitor Treatment**




## PrimePanel RET Fusion Cell Line: CR1520

The PrimePanel CR1520 primary cancer cell line was generated from the CR1520 PDX model, and validated to confirm that the NCOA4-RET fusion is still present (Figure 8). The CR1520 cell line has been shown to be sensitive to a range of RET inhibitors (Figure 9).


# RET Fusion Factsheet

**Figure 8: PrimePanel RET Fusion Cell Line Genomic Profiling: CR1520 RT-PCR Analysis**

NCOA4-RET fusion product (271bp) confirmed by RT-PCR.



**Figure 9: PrimePanel RET Fusion Cell Line CR1520 is Sensitive to RET Inhibitor Treatment**




RET Inhibitor	Sorafenib	Sunitinib	XL-184	ZD6474
IC <sub>50</sub> (μM)	3.87	3.45	4.58	1.46

## PrimeXeno RET Fusion Xenograft: CR1520

A PrimeXeno *in vivo* model has been developed from the CR1520 PrimePanel cell line. The subcutaneous xenograft model shows robust growth following inoculation (Figure 10).

**Figure 10: PrimeXeno RET Fusion Xenograft Model CR1520 Growth Curve**



## Conclusions

The recent identification of RET fusions in CRC and lung cancer has led to a need for RET fusion models across all stages of drug discovery. CrownBio provides a range of xenograft and cell line RET fusion platforms, to provide a critical link between *in vitro* testing and clinical efficacy in drug development.

Our HuPrime PDX models provided the first description of oncogenic RET fusions in CRC patient samples<sup>(7)</sup>, NCOA4-RET fusion and CCDC6-RET fusion, which have been further observed in patients with advanced CRC<sup>(5)</sup>. Our models provide clinically relevant xenografts for evaluation of therapeutics targeting RET fusions, with models responsive and resistant to RET inhibitors available.

Our PDX models have also been used to derive a PrimePanel cell line to provide a corresponding *in vitro* cell culture system for cost-effective and high-throughput drug screening and model selection. The related PrimeXeno CDX model completes our full resource platform of cell line, CDX, and PDX for all research requirements around RET fusions.

CrownBio can be contacted at [busdev@crownbio.com](mailto:busdev@crownbio.com) for any further questions or information required on our RET fusion resources, or for information on other CrownBio products and services.

## References

- My Cancer Genome® RET Fusions in Thyroid Cancer <http://www.mycancergenome.org/content/disease/thyroid-cancer/ret/127> Accessed 08 April 2016.
- Kohno T, Ichikawa H, Totoki Y *et al.* KIF5B-RET fusions in lung adenocarcinoma. *Nature Medicine* 2012; 18(3): 375-377.
- Takeuchi K, Soda M, Togashi Y *et al.* RET, ROS1 and ALK fusions in lung cancer. *Nature Medicine* 2012; 18: 378-381.
- Drillon A, Wang L, Hasanovic A *et al.* Screening for RET and ROS1 fusions in an enriched cohort of pan-negative never-smokers with advanced lung adenocarcinomas to identify patients for treatment in targeted therapy trials. *Journal of Clinical Oncology* 2013; 31(suppl abstr 8067).
- Le Rolle A-F, Klempner SJ, Garrett CR *et al.* Identification and characterization of RET fusions in advanced colorectal cancer. *Oncotarget* 2015; 6(30):28929-28937.
- Shaw AT, Hsu PP, Awad MM *et al.* Tyrosine kinase gene rearrangements in epithelial malignancies. *Nature Reviews Cancer* 2013; 13(11):772-787.
- Gozgit JM, Chen T-H, Clackson T *et al.* RET fusions identified in colorectal cancer PDX models are sensitive to the potent RET inhibitor ponatinib [abstract]. In: Proceedings of the 105th Annual Meeting of the American Association for Cancer Research; 2014 Apr 5-9; San Diego, CA. Philadelphia (PA): AACR; *Cancer Research* 2014; 74(19 Suppl):Abstract nr 2726.
- Zhang Y, Ge Y, Liu Y *et al.* PrimePanel provides a high throughput *in vitro* drug screening platform that intimately links to *in vivo* pharmacological analysis in PDX models [abstract]. In: Proceedings of the 104th Annual Meeting of the American Association for Cancer Research; 2013 Apr 6-10; Washington, DC. Philadelphia (PA): AACR; *Cancer Research* 2013; 73(8 Suppl):Abstract nr 2787.