Safety and
Environmental Protection

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Storage of hazardous materials in 21st Century laboratories

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»» Why settle for less? The storage of hazardous materials according to European Standards (EN 14470 Part 1 and Part 2) is safe and flexible!

Sascha Kunkel, International Sales Manager,
expert in safe storage and handling of hazardous materials

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Maximum safety for you

Don't compromise when storing hazardous materials.
Trust in state of the art technology.

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In 21st century laboratories flexible layouts are of the utmost importance. The operational purposes of a laboratory change quite frequently and it is important to consider a sustainable storage concept for hazardous materials, such as flammable liquids, during the initial planning stage. A central fire resistant storage room is a thing of the past. Today fire resistant storage cabinets play a major role in realizing flexible concepts in modern laboratories, while maintaining a high degree of safety and preventing higher costs in future lab refurbishments.

An efficient and appealing room concept and technical equipment that adapts itself to changing laboratory procedures – is the aim of every laboratory. The use of suitable solutions when furnishing means that no compromises need to be made with regard to functionality and above all safety. The following pages provide detailed information.

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It can happen to anyone

Many different dangers threaten when handling hazardous materials

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The use of hazardous and combustible materials is part of the daily routine in most laboratories and is unavoidable. What is often quickly forgotten in daily handling: the improper storage of such materials puts people, the environment and property in danger.

- **Fires and explosions** caused, for example, by the improper storage of combustible liquids
- **Pollution of soil, groundwater and surface water** e.g. by contaminated fire fighting water or leakages
- **Downtimes**, for example due to laboratories destroyed by fires
- **Endangerment to human life**

Hazardous materials are therefore always very much a current topic for every laboratory. Apart from the proper handling of these health-endangering substances, particularly strict storage regulations must be observed.

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The legal basis concerning the handling of hazardous materials is defined

- in global regulations (e.g. GHS – globally harmonized system)
- in directives and guidelines of state confederations (e.g. the Europe Union)
- and finally in rules and laws of individual countries

- Do you know the legislation for handling hazardous materials in your country?
- Do you apply these rules?
- Do you fulfill your duty of care for humanity and environment?
- Do you already use safety storage cabinets with the highest protection level?

Only if you have answered all these questions with „Yes“ you are on the safe side!

But nothing has ever happened so far...

„Life long experience has shown that one has to face the possibility of an outbreak of a fire at any time. The fact that in most buildings there was no outbreak of a fire for decades is no proof that this danger does not exist. It is rather long-lasting luck which could end at anytime.“

Quote of a decision of the higher administrative court in Münster, Germany, Az.: 10 A 363/86, 11.12.87

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Hazardous materials – terms, definitions and labelling

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ALCOHOLS
(PROPAN-1-OL)

ALCOHOLS, N.O.S.
(PROPAN-1-OL, ETHANOL)

Classification and labelling of hazardous materials

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All chemicals are subject to compulsory classification and labelling prior to being brought onto the market. The goal of this obligation is to identify the properties of the substances and mark them with the relevant hazard symbols. People and the environment are therefore protected against harmful effects when handling chemicals.

In the past, worldwide systems for the classification and labelling of chemicals were very different. It was possible for a substance or a mixture of substances to be classified as hazardous and handled as such in one country and not in another. This led to problems in terms of health and safety at work.

With the development of the GHS (Global Harmonised System) in September 2008, a worldwide uniform system was created for the classification and labelling of chemicals. Wherever this global harmonised system is introduced, be it in China, India, the USA or Europe, chemicals are classified and labelled according to uniform criteria.

The goals of the GHS are to enhance protection of mankind and the environment, by:

- Providing an internationally comprehensible system for hazard communication.
- Providing a recognized framework for those countries without an existing system.
- Facilitating international trade in chemicals whose hazards have been properly assessed and identified on an international basis.
- Reducing the need for testing and evaluation of chemicals.

The GHS classification system is very complex. It consists of multiple physical, health and environmental hazards. The most important standardized label elements included in the GHS are:

- **Symbols (hazard pictograms):**
- **Signal Words** ("Danger" or "Warning")
- **Hazard Statements**

The diagram provides an overview of the highly visible pictograms with a red border:

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GHS/CLP New

	Explosive
	Compressed gases
	Flammable Category 1 Flammable Category 2 Flammable Category 3
	Oxidising Category 1, 2, 3
	Acute toxicity Category 1 Acute toxicity Category 2 Acute toxicity Category 3
	Corrosive Category 1
	Acute toxicity Category 4 Corrosive Category 2 Skin irritant Specific target organ toxicity Category 3
	CMR* Category 1A, 1B, 2 Respiratory sensitization, specific target organ toxicity Category 1, 2
	Environmental pollutant Hazardous to the aquatic environment Damaging to the ozone layer

CMR*: C = carcinogenic, M = mutagenic, R = toxic to reproduction

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Storage:

Storage is the keeping of hazardous materials for later use and for dispensing to others. Hazardous materials are to be kept or stored in such a way that they do not endanger human health and the environment. Abuse or incorrect use must be prevented.

Flash point:

The temperature at which enough vapour escapes from a liquid that an ignitable mixture is created with the ambient air (EN 1127-1).

Definition auto-ignition point

The auto-ignition point is defined as the lowest temperature at which vapours, from a flammable liquid in contact with air or a hot object, ignite themselves. Many widely-used chemicals (e.g. Diethylether) have an auto-ignition point of 200 °C (392 °F) or greater.

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
Definition of explosion limits

The Lower Explosion Limit (LEL) and Upper Explosion Limit (UEL) are respectively the lower and upper value limits for the concentration of a combustible substance in a mixture of gases, vapours, mists and/or dusts. Following the ignition, a flame that is independent of the source of ignition becomes unable to independently propagate itself.

Example of explosion limits

Substance	LEL in Vol.-%	UEL in Vol.-%
Hydrogen	4	75
Propane	2.1	9.5 - 10.1
Methane	4.4 - 5	15 - 17
Ethanol	3 - 3.3	19
Silane	1.5	98

Flammable liquids



Liquids with a flash point equal to or greater than 23 °C (73.4 °F) and less than or equal to 60 °C (140 °F) and which support combustion when tested at 60 °C (140 °F).

Highly flammable liquids

Liquids with a flash point lower than 23 °C (73.4 °F) and with a boiling point higher than 35 °C (95 °F).

Extremely flammable liquids

Liquids with a flash point lower than 0 °C (32 °F) and with a boiling point lower than or equal to 35 °C (95 °F).

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Fire triangle – 3 requirements for an explosion

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The basis for the handling of hazardous materials is knowledge of the fire triangle. This illustrates in an unforgettable way that a fire or an explosion is only possible if three fundamental requirements are satisfied or present:

Source of ignition (e.g. mechanical or electrical spark, chemical reaction, static charge, hot surface, open flame)

Combustible substance (e.g. gases, liquids, solids)

Oxygen (the correct mixing ratio of combustible material and oxygen)

To prevent fires it must be made sure that at least one of the three requirements can be excluded.

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
Storage of hazardous materials

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Storage of hazardous materials – past and present

A comparison of methods – central and local storage

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Looking back over recent years, developments in laboratory safety have greatly advanced. A few years ago it was still usual to centrally store hazardous materials in a separate hazardous materials store.

The current state of the art storage solutions, on the other hand, provide local storage, close to the workplace, in modern safety storage cabinets. In addition to maximum safety and flexibility, the requirements for a modern laboratory include the most efficient work procedure with the shortest routes possible.

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The advantages of a safety storage cabinet compared to a central hazardous material store:

Transportation/Working hours

- Avoiding internal transport of hazardous material from the storage room to the workplace.
- All hazardous materials for daily use can be stored safely, directly at the workplace.
- Efficient use of working time, with no need to travel far every day just to bring flammables in and out from the central storage room.
- Quick access to bottles and containers.

Maximum safety for the user

- Only minimum quantities needed for a specific job are in use at the workplace. All other hazardous materials are stored safely in the fire resistant cabinet.
- Maximum fire protection due to the centralized storage of all flammables at the workplace in a fire resistant safety storage cabinet (Type 90 – EN 14470-1).
- Explosive atmospheres are avoided when hazardous materials are stored in technically ventilated safety storage cabinets.

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
Safety storage cabinets – overview

What distinguishes modern safety storage cabinets, what can they do?

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Just like a central hazardous material store, the purpose of a safety storage cabinet is the safe keeping of combustible chemicals or compressed gas cylinders. The areas of use are varied: laboratory, industry, handicraft, universities, hospitals, public organisations, etc.

The primary goal of a modern safety storage cabinet is the minimisation of the fire risk when storing combustible liquids or compressed gases. In case of fire the cabinets prevent the substances stored within from assisting the spread of the fire. A special fire-retardant construction makes this possible.

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Looking back – strict standards for the storage of hazardous materials were only known in Germany. As far back as the 1980s there was a standard – DIN 12925 – for the storage of combustible liquids (part 1) and compressed gas cylinders (part 2). These high standards, unsurpassed in international comparison, were the basis for the present-day European standard EN 14470, which has defined the requirements with regard to the design and fire resistance of safety storage cabinets since 2004.

Altogether the EN distinguishes between four cabinet type classes: 15, 30, 60 and 90. The respective type class indicates how long the cabinet protects the hazardous materials stored inside against fire. A type 90 cabinet, as used in Germany, for example, protects the chemicals stored inside it for one and a half hours in case of fire: additional time to allow employees to reach safety or for the rescue services to initiate fire-fighting operations.

Safety storage cabinets are distinguished by the type of substances to be stored. There are cabinets for the storage of

- Combustible liquids
- Compressed gas cylinders
- Acids and alkalis
- Non-combustible Chemicals

Cabinet models are available on the market in various widths and heights, with different door and door arrest systems and innumerable equipment variations. There are also very big differences in the cabinet designs with single-wall steel cabinets that are sold as safety cabinets up to type 90 EN cabinets.

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What must a safety storage cabinet be able to do?

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The most important protection aim of a safety storage cabinet is to reliably shield flammables against a critical temperature increase in the event of a fire over a defined period of time.

- avoid explosions accelerating the fire
- adequate time for personnel to escape
- adequate time for rescue/fire services for fire-fighting and rescue missions

What happens if 200 °C (392 °F) is reached? Explosions can occur if the temperature rises by more than 200 °C (392 °F), because many flammable chemicals start to reach their auto ignition point.


The temperature curve shows that after only 5 minutes a fire has reached a temperature of 576 °C (1068.8 °F).

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Comparison of fire resistant safety times for different cabinet constructions

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1. Double wall steel cabinet

How long will the stored flammables remain safe? The illustration makes the difference very clear: the double wall steel cabinet offers no more than 3 minutes of safety before the temperature increase exceeds 180 °C (356 °F) inside the cabinet, when subjected to a fire test according to the European standard.

2. Type 30 EN safety storage cabinet

A Type 30 EN cabinet from asecos guarantees 30 minutes of safety before the temperature increase on the inside becomes critical. 10 times the amount of safety time provided by a double wall steel cabinet!

3. Type 90 EN safety storage cabinet

asecos Type 90 EN cabinets offer maximum safety. It takes more than 90 minutes of exposure to a fire before the temperature inside rises by more than 180 °C (356 °F). 30 times the amount of safety time! This gives plenty of time for the rescue services and fire brigade to respond.

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


- Single wall or double wall steel cabinets without any fire protection measures do not fulfil the requirements of EN 14470-1.
- 30 minutes of fire resistance provides the basic protection for rescue services and personnel
- 90 minutes of fire resistance means optimum quality and therefore the highest possible safety
- Cabinets meeting EN 14470-1 give your investment maximum security

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The current state of the art safety storage cabinets

Safety storage cabinets designed and tested according to European standard 14470 offers maximum protection

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Although there have always been strict regulations in Germany for the storage of hazardous materials, these were substantially tightened - with regard to the design and testing of safety storage cabinets - in 2004 with the European standard 14470. The following points describe the basic and main safety, test and construction requirements.

The range of applications

- The standard specifies the performance required of the design and fire resistance of safety storage cabinets used for storing flammable liquids in working areas.

The principal safety requirements

- Minimisation of the fire risk associated with the storage of flammable materials.
- Protection of the cabinet contents in the event of fire for a known (and tested) period of time.
- Minimisation of the fumes released to the working environment.
- Retention of any possible leakage within the cabinet.
- Provision of enough time, in the event of fire, for personnel to leave the room, and sufficient time for fire service personnel to enter the building before the stored materials turn a small fire into an uncontrollable blaze.

Fire protection

- In the event of a fire the cabinet must ensure that, over a period of time defined by the manufacturer (but at least 15 minutes), its contents do not present an additional risk that the fire will spread.

Doors

- The cabinet doors must close entirely, starting from any position (closing time max. 20 seconds).
- Immobilising equipment fitted must release the locked doors at a temperature of 50 (-10) °C (122 (-50) °F).
- Avoiding risk of injury: the closing force of the doors must not exceed 100 N.
- One-handed operation must be possible and the doors must close entirely even if open and locked.

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Side and rear walls

- The side and rear walls of the cabinet must have the same thickness and comparable structures.

Air inlet and outlet openings

- The cabinets must have openings for air inlet and outlet (to connect the cabinet to an exhaust system).
- The ventilation openings must close automatically at a temperature of 70 °C (158 °F).

Storage locations (shelves or drawers)

- Storage surfaces must be able to support the loading specified by the manufacturer over the period of the test in the furnace.

Spill containment sump

- The spill containment sump must retain its ability to function after the fire resistance test. This is to be checked visually by filling the spill containment sump with water.

Fire resistance

- Must be investigated by tests on a design sample.
- The fire resistant cabinet is exposed to flames in a suitable furnace.
- The doors, walls and ceiling of the cabinet being tested must be exposed to the same heating conditions.
- Cabinets must be tested as freestanding single cabinets. The example being tested must be positioned with its rear wall at least 100 mm from the furnace wall.
- The flame exposure is carried out in accordance with the standard temperature curve of EN 1363-1.
- The temperature rise is measured inside the cabinet.
- The cabinet must then be classed as type 15, 30, 60 or 90, according to the time that has elapsed before the temperature rose by 180 °C.

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The Advantages of fire rated safety storage cabinets

The comparison of a double-wall steel cabinet with an EN safety storage cabinet demonstrates impressive differences

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1

Cabinet structure

EN 14470-1 CABINETS: Multiple layers of fire-proof mineral fibre insulation material (Calcium sulphate), encased in sheet steel.

STEEL CABINETS: Only double wall sheet steel.

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2

Fire test

EN 14470-1 CABINETS: The cabinet is fire tested for 15, 30, 60 or 90 minutes, depending on the fire resistance class. Longer resistance against a temperature increase means longer protection for personnel and more time for the fire fighters/rescue services.


STEEL CABINETS: Cabinets are fire tested for only 10 minutes.

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3

Maximum temperature increase allowed

EN 14470-1 CABINETS: The temperature increase must not exceed 180 K (relative temperature increase in Celsius) at any measuring point during the entire duration of the fire test. Many flammable liquids reach their auto ignition point at approximately 392 °F (approx. 200 °C).

STEEL CABINETS: The temperature increase must not exceed 325 °F (approx. 163 °C) in the first 10 minutes of exposure to a fire.

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4

Self-closing doors

EN 14470-1 CABINETS: Mandatory. If no door open arrest system is used the doors must close automatically from every position when released.

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5


Door open arrest system

EN 14470-1 CABINETS: Possible, doors must close at 122 °F maximum (approximately 50 °C), the closing sequence must be completed within 20 seconds, starting from any position.

STEEL CABINETS: Possible, doors must close at 212 °F (100 °C) maximum.

6

Perforated metal plate insert

EN 14470-1 CABINETS: Available for sumps and drawers. If sumps/drawers are used as a storage level, the stored containers reduce the retention volume of the sump, which then could lead to an overflow if liquids are spilled. For this reason perforated inserts are put into the sumps to hold the entire retention volume.

STEEL CABINETS: Not specified.

7

Bottom collecting sump

EN 14470-1 CABINETS: Mandatory, capacity has to be at least 10 % of the volume of all containers stored or a minimum of 110 % of the largest container stored, whichever volume is greater.

STEEL CABINETS: Mandatory, no capacity mentioned, only minimum height specified, regardless of collection capacity.

8

Fusible Link

EN 14470-1 CABINETS: Fusible links at door open arrest system, ventilation spigots and drawers. A mechanical fusible link is a device consisting of two strips of metal soldered together with a fusible alloy.

STEEL CABINETS: Fusible links only at door open arrest system.

9

Insulation

EN 14470-1 CABINETS: Multiple layers of fire-proof mineral fibre insulation material (Calcium sulphate). A steel corpus protects against damage from outside but has no heat retardant properties; in fact it is a strong heat conductor.

STEEL CABINETS: Plain air insulation.

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10

Seals

EN 14470-1 CABINETS: All gaps become totally sealed between the door and frame with special intumescent seal protection strips which, when exposed to a fire, expand and seal hermetically to avoid heat entering the cabinet.

STEEL CABINETS: This is not a specified requirement.

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11

Exhaust air outlets /spigots

EN 14470-1 CABINETS: Exhaust air outlets/spigots are mandatory. The recommendation is for cabinets to be connected to a 10:1 air change ventilation system.

STEEL CABINETS: mandatory, but no recommendation regarding the technical ventilation of the cabinet.

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12

Ground connection

EN 14470-1 CABINETS: Possible. Limits the built-up of static electricity (possible ignition source!) when connected to the Earth.

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13

Storage volume

EN 14470-1 CABINETS: Actual usable storage space for each shelf/drawer is stated. The total interior volume of any cabinet is always much higher than the actual usable volume. Indicating the volume per shelf/drawer is a much more practical and realistic approach, which enables the user to exactly calculate whether a particular cabinet model/size is sufficient for their requirements.

STEEL CABINETS: The total internal volume is stated in most cases e.g. 120 gallons (454 litres).

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14

Third party testing

EN 14470-1 CABINETS: Yes. Only an independent third party test can ensure that all safety features comply to the standards.

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DOPPELWANDIGE STAHL-SCHRÄNKE: Ja.

Storage of hazardous materials in laboratories / European EN standard vs. double wall steel cabinets

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Statement extracts of Klaus Söhngen, CEO of the international active HVAC and laboratory planning office ERETEC, regarding the topic "storage of hazardous materials in laboratories / European EN standard vs. double wall steel cabinets"

... During our international planning activities, we observe that the opinion concerning safe storage of flammable liquids and pressurised gases is often different to our opinion, particularly outside European countries ...

... due to our European understanding of safety, it is incomprehensible for us that flammable liquids, hazardous materials and pressurised gases are stored in presumably safe double wall steel cabinets. On closer examination, these cabinets, guarantee only insufficient to no fire protection in the event of an emergency. Laboratory users and operators are led to believe that safety is provided, but in case of emergency it is not...

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... in the following illustration (ground plan of a laboratory) one can easily see, that nowadays laboratories are planned with open workplaces and service areas. Especially in the area of core facilities, continuous areas of several hundred square metres are common...

...these types of ground plans are not limited to Europe. Especially in American-influenced countries, these types of laboratories are also planned and operated. From our point of view, it is highly questionable to use steel cabinets in laboratories where many people work. These single-walled resp. double-walled cabinets including an air insulation system do not provide a fire resistance classification and offer only a few minutes of fire protection applying European-Standard criteria...

...apart from automatic fire extinguishing systems and other efforts to raise fire protection, an intelligent and reliable safety concept for laboratories should always include equipment with fire-resistant cabinets (type 90)...

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From the point of view of the HVAC planner and architect

...as hazardous materials can be stored safely in fire-resistant cabinets (type 90) directly at the workplace, the planning of a central storage room is no longer necessary in all cases.

Therefore, the investment cost for such rooms and their equipment can be saved.

If the architectural concept includes these rooms anyway, additional space can be gained, as the rooms can be used for other purposes...

...as safety storage cabinets can be planned and flexibly installed within the ground plan, a higher flexibility of availability is guaranteed for the workflow...

...the use of fire-resistant (type 90) cabinets can also contribute to a leaner and less extensive constructional fire protection, as the required fire protection is already integrated in the cabinet. The costs for more expensive constructional fire protection components, such as sprinkler systems, for instance, can be completely saved...

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From the point of view of the user, fire protection and emergency staff

...Due to the use of fire-resistant cabinets (type 90) the employees of a laboratory have on the one hand much more time to leave the danger zone in case of fire. On the other hand, they can be sure that the protection mechanism of the cabinets will avoid an acute or short-term deterioration of the situation. The fire brigades and paramedics also profit from this fact, as the fire-resistant cabinets (type 90) avoid an additional risk regarding an uncontrolled or accelerated fire spread during the fire resistance time of the cabinet. Therefore, the emergency-services are faced with an easier situation on-site, as the circumstances and risks can be assessed more reliably... thus the installation of fire-resistant cabinets (type 90) causes a minimized economic risk (risk of total damage). Especially, in countries where the cabinets according to European standards are not obligatory, this fact can lead to reduced insurance costs...in addition, the use of fire-resistant cabinets directly at the workplace (type 90) drastically reduces the transportation of hazardous materials and gases and transport-accidents are almost non-existent...

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Safety storage cabinets according to EN 14470 – the perfect choice

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Statement extracts of Jun.- Prof. Dr. Ing. Marcus Marx, state-certified safety engineer and junior professor for fire and explosion protection at the Otto-von-Guericke-University Magdeburg.

...safety storage cabinets with a proven fire resistance class in accordance with EN 14470 enable modern laboratories – depending on the tested fire resistance duration of the cabinet – a less risky, local storage of

- toxic,
- flammable or
- potentially explosive substances
- and/or of substances with combinations of the stated dangerous characteristics.

Apart from the purely architectonic advantages, such as additional space, and the laboratory-specific benefits like time saving in the supply of materials, the use of safety storage cabinets according to EN 14470 also signifies an immense safety advantage regarding preventive and defensive fire protection.

The increasingly flexible design of laboratories and the fast change of applications makes a decentralised fire protection concept desirable.

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
This concept should be quickly adaptable without any complex structural changes and without downgrading safety...

...in case of fire, the risk which could arise from the substances stored in the cabinet is almost non-existent during the fire-resistant time of the cabinet. If the substances are stored properly in the safety storage cabinet, there is no need for preferential consideration of the fire accelerating or even potentially explosive substances during the fire-fighting and rescue operations. Thus, the emergency services can primarily concentrate efforts to control the fire and to rescue injured persons...

...also from the point of view of insurance, a certified safety storage cabinet has positive aspects. The anticipated risk in case of fire is minimized by less severe consequences. Therefore the classification according to EN-standards makes a more precise risk assessment possible...

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...these safety-related statutory specifications do not exist in all countries. Laboratory staff using double wall steel cabinets are exposed to a higher risk, as these cabinets do not fulfil strict European specifications...


...the illustration shown below demonstrates the comparison between a safety storage cabinet (according to EN 14470) and a storage room for hazardous materials and flammable liquids designed in accordance to European standards.

This illustration clearly shows that the safety storage cabinet has comparable safety-related characteristics to the storage room. According to European standards, the safety storage cabinet thus becomes the "storage room within the laboratory". On the one hand the substances are protected according to the fire resistance duration and on the other hand the laboratory is protected against the substances stored in the cabinet...

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...the insulated safety storage cabinet according to EN helps to restrict local damage. Safety storage cabinets with other construction standards often have no insulation, but there are sprinkler systems in the rooms where they are installed to limit the temperatures in case of fire. In this case the rapid control of the fire is contradicted not only by typical water damage, but also by the potential spread of hazardous materials mixing with the fire fighting water. Furthermore, the sprinkler installation is a so-called active safety system, which only reacts in the event of fire. This means that you only know in the event of fire, whether the system is working or not...

...the European safety storage cabinet excels with its more significant passive safety provided by its fire resistance...this system does not require activation in the event of a fire, as it is immanent...

...The equipment of laboratories worldwide with safety storage cabinets in accordance with European safety standards is desirable, not only with regard to personal safety but in the event of a fire, far less property damage can be expected. This emphasises the practical added value of the cabinets compared to less expensive cabinets which do not comply with the requirements of EN 14470...

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Comparison Safety Storage Cabinet Fire resistance Type 90 / Storage room Fire resistance F90/T30

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Cabinet walls with a fire resistance 90 minutes
Automatic closing of doors in case of fire - fire resistance 90 minutes
Technical ventilation
spill containment sump system
stored hazardous materials

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Walls with a fire resistance 90 minutes
Technical ventilation
stored hazardous materials
spill containment sump system
Automatic closing of doors in case of fire - fire resistance 30 minutes

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Storage of compressed gas cylinders – central or local?

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When storing compressed gases the question also arises: Should you store them in a gas cylinder cabinet directly at the place of work or in central storage rooms? A comparison of the two models can help to make this decision.

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With their extensive requirements and fields of application, gas supply systems (particularly for special gases) also place special demands on the components for the storage, distribution, control and monitoring of the gases. Before installing new gas supply systems or extending existing ones, the supply concept that offers the most advantages should always be determined. In principle, a distinction must be made between a central and a local gas supply.

Central gas supplies are systems in which the required gas containers are installed in a single place and the gases are fed to the points of consumption via a permanently installed piping system. With central gas supplies a further distinction can be made between outdoor installation and installation in a storeroom.

Local gas supplies are individual supply units that are installed (usually inside the building) directly at, or in direct proximity to, the point of consumption.

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When faced with the choice of one of these types of gas supply, current technical developments should always be taken into account. Today, gas cylinder cabinets tested in accordance with EN 14470-2 offer the same level of safety as safety storage cabinets for the storage of flammable liquids – with up to 90 minutes of fire resistance.

The local provision and withdrawal of compressed gases is often the more economical solution, because, particularly in the case of corrosive gases, the price for the long pipelines required is very high.

With local storage the necessary pipelines are reduced to a minimum. This leads, for example, to a considerable reduction in the consumption of purging gases.

If the flexible local use of compressed gases is necessary, then the more appropriate choice would be the type G90 safety storage cabinets.

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The pros and cons of central and local storage

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Central storage (external store and storeroom):

- ⊖ Protective zones in accordance with TRG 280 must be maintained around the cylinders (external store/storeroom)
- ⊖ Expensive piping system up to the points of use (external store/storeroom)
- ⊖ Partitions must be at least fire-retardant (storeroom)
- ⊖ Doors must be self-closing (storeroom)
- ⊖ Floor covering must be at least fire-resistant (storeroom)
- ⊖ Adequate ventilation openings leading directly to the open air are required (storeroom)
- ⊖ No additional storage of flammable liquids (storeroom)
- ⊕ Saving space at the points of consumption
- ⊕ Central monitoring of the entire gas stock
- ⊕ Low in-plant transport costs
- ⊕ Better monitoring of the cylinder stock levels
- ⊕ Smaller cylinder pool thanks to fewer reserve cylinders

Gas cylinder cabinets in accordance with EN 14470-2

- ⊕ Protective zones around the cylinders are unnecessary
- ⊕ Short routes for necessary pipelines to the points of use
- ⊕ Hence low installation and usage costs
- ⊕ Monitoring of gas stock directly at the place of work by integrated low gas indicators
- ⊕ Increased fire protection for up to 90 minutes in accordance with DIN EN 14470-2
- ⊕ Offers the same safety level as safety storage cabinets for the storage of flammable liquids
- ⊕ Flexible and uncomplicated adjustments for changes of use

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Safety with a seal of approval

EN-tested safety storage cabinets from asecos live up to their promise – independent testing institutes also confirm this

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When purchasing safety storage cabinets you should make sure that the approval documentation is complete. Only then can you be sure that the cabinet actually does what it promises to do.

Safety in case of fire

All asecos cabinet models, with a fire-resistant classification, have been successfully tested by an independent material testing body in a fire chamber (type tested). These cabinets possess a type-tested fire resistance according to EN 14470-1 or 2.

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DIN EN

Safety in daily use

All cabinets marked with the GS logo have been successfully tested according to the requirements of the German Equipment Safety Act – documented by the test certificates according to the Equipment Safety Act by independent testing organisations. The GS certificate guarantees that, in addition to fire protection, the cabinets also function perfectly in everyday use. During the GS certification the safety storage cabinet passes through numerous tests. The cabinet is examined to determine whether there are squeezing and trapping points or risks of injury. In addition the cabinet and the door mechanism are subjected to a tough fatigue test: the GS certificate is only issued if 50,000 door closing cycles are completed without damage, which is a real test of endurance and safety!

Quality control

The manufacture of asecos safety storage cabinets is subject to voluntary, permanent external supervision by a recognised independent material testing body.

Cabinets are regularly submitted to a repeat test in the fire chamber.

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EN-tested safety storage cabinet in hardness test

Type 90 safety storage cabinet proved in a real-life situation, what the test documents certify

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An industrial enterprise from Bamberg in Bavaria suffered what all companies fear: following a technical defect that was no-one's fault, fire broke out in its production halls. Fortunately, people and the neighbouring buildings came to no harm. Despite extensive attempts by the local fire brigade, however, it was not quite possible to extinguish the fire and the hall had to be allowed to burn down overnight in a controlled manner.

During the inspection of the burnt-down hall, an asecos type 90 safety storage cabinet for the storage of combustible liquids was examined. Although its outer body looked very much the worse for wear following the fire, which lasted several hours, the cabinet had retained its safety function despite the temperatures of one to two thousand degrees Celsius that had developed. The fire protection insulation and seals had performed their tasks well.

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After opening the cabinet, it was possible to salvage all of the hazardous materials inside. Both the plastic containers and the flammable liquids contained in them had been shielded against the heat by the safety cabinet.

What would have happened if the flammable hazardous materials had also been caught up in the fire is anyone's guess. One thing is sure however: the current type 90 technology, with at least 90 minutes of fire resistance, more than lives up to its promise, because the duration of 90 minutes of heat acting on the cabinet was exceeded by far in the fire in Bamberg. Type 90 safety cabinets are therefore a mature and proven preventative measure for protection against the fast propagation of fire.

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Technical ventilation of safety storage cabinets

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In Germany, the law allows the operation of safety storage cabinets without permanent effective technical ventilation. In this case, however, other measures must be taken to guarantee comparable protection against explosion. If corrosive or toxic materials are to be stored in the safety storage cabinet, for example, or in the case of the employees being disturbed by unpleasant odours, then technical ventilation should always be provided.

There are various possibilities to implement the necessary technical extraction:

Plug-in extraction units

These need only be placed on the roof of the safety storage cabinet and connected to the existing extraction system. A wide range of extraction units are available with and without exhaust air monitoring, with potential-free alarm contact, for wall mounting, etc.

Radial ventilators

Powerful radial fans can be integrated into exhaust air systems and provide reliable extraction from the safety storage cabinets.

Recirculating air filter systems

Recirculating air filter systems are simply mounted on top of the safety storage cabinet and connected to its exhaust ducts. They are equipped with an active carbon filter, which removes the hydrocarbons contained in the exhaust air within the cabinet. The cleaned air is returned to the installation room. No customer exhaust air system is required.

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Pro and cons...

... without technical extraction from the safety cabinet?

- + no cost-intensive exhaust air ducts necessary
- definition of explosion zones in and around the safety storage cabinet (illustration/description)
- possibly dangerous to the health of the employees

... with technical extraction (min. 10 times/h) and exhaust air fed outside?

- + no danger to the health of the employees
- + all hazardous substances (including acid and alkali vapours) are eliminated
- cost-intensive exhaust air ducts with possible wall/floor breakthroughs are necessary
- installation of the safety storage cabinet is no longer flexible
- high energy consumption due to heat loss in the exhaust air

... with technical extraction (min. 10 times/h) with a recirculating air filter system without feeding the exhaust air to the outside?

- + no danger to the health of the employees
- + no cost-intensive exhaust air ducts necessary
- + flexible installation of the safety storage cabinet
- + no loss of energy due to feeding the exhaust air to the outside
- + permanent electronic monitoring of the exhaust air quantity
- + tested and certified system
- only solvent vapours (hydrocarbons) can be filtered

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Recirculating air filter system
Umluftfilteraufsatz

Radial ventilator
Radialventilator

Storage of non-flammable chemicals

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In addition to combustible liquids and compressed gases, a large number of other chemicals are required in laboratories and these also have to be safely stored.

The chemical storage cabinets offer the perfect solution for all requirements. Cabinet models in three different sizes are available, depending on the number of containers to be stored. The doors can also be individually selected - wing doors, space-saving folding doors or doors with glass inlays. You even have a choice of colour: there are 7 door colours to choose from!

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Storage of acids and alkalis

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Acids and alkalis can give off aggressive vapours, mists or smoke. Therefore, these substances must be stored in such a way that neither people nor the environment are endangered.

Acids and alkalis can be stored directly at the workplace – separately if possible – in special so-called “acid/alkali cabinets” made mostly of corrosion-resistant material. In Germany, however, there is no DIN or EN standard for these storage cabinets. The only legal requirements for these storage cabinets arise from various laws.

- Equipped with sumps to prevent leakages. They must be able to retain 10% of the stored amount or at least the contents of the largest container.
- The cabinet must be resistant to corrosion and compatible with the stored acids and alkalis.
- The cabinet must be connected to an exhaust air blower, which constantly feeds the escaping vapours to the open air.
- Strong acids and strong alkalis must be stored and extracted in separate areas of the cabinet (separate from each other) because their vapours/mists can react and create heat.

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Laboratory planning using EN safety storage cabinets

Different cabinet sizes allow for tailored laboratory layouts

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The requirements of laboratories are very individual. Planning the laboratory must be just as individual. Cabinets for the storage of flammable liquids and gas cylinders are available in different heights and widths, making the integration into any new laboratory very easy.

EN approved **tall cabinets** are available in:
3 widths: 120 cm, 90 cm and 60 cm (47.24 in, 35.43 in, 23.26 in)
2 heights: 129 cm, 196 cm (50.79 in, 77.17 in)

EN approved **underbench cabinets** are available in:
4 widths: 140 cm, 110 cm, 89 cm and 59 cm (55.12 in, 43.31 in, 35.04 in, 23.23 in)


EN approved **gas cylinder cabinets** are available in:
4 widths: 140 cm, 120 cm, 90 cm, 60 cm (55.12 in, 47.24 in, 35.43 in, 23.26 in)
2 heights: 205 cm, 145 cm (80.71 in, 57.09 in)

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

Wing Doors vs. Folding Doors

Flügeltüren vs. Falttüren

It is also true for a Type 90 safety storage cabinet: Keep the doors closed when not in use. Only that way maximum safety can be guaranteed.

Compared to double-walled steel cabinets the cabinets according to the European Standard offer considerably more protection – up to 90 minutes fire resistance – and ease of use. It begins at the doors: Wing doors or space saving folding doors can be chosen.

Folding doors are the best choice whenever space is limited, like in the laboratory. They make it possible to pass the cabinet even when the doors are fully open. Furthermore there are cabinets available where you can open both doors with just one hand and even fully-automatic cabinet doors with sensor detection are available. A real benefit with regard to safety when the containers can be held with both hands at all times.

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Frequently changing requirements on the laboratory are no problem for modern safety storage cabinets. Relocating the cabinets e.g. when moving or refurbishing the lab is very easy.

asecos safety storage cabinets can be equipped with an optional base. Installation or relocation of these cabinets can be easily done with a pallet truck. The latest cabinet models can even be moved through standard room doors using a special transport device (Q-mover). All underbench cabinets from asecos that come with a base are equipped with sturdy castors making it very easy to move the cabinets e.g. for service works at the fume cupboards. If the cabinets are also equipped with a recirculating air filter system even the connection to the ventilation system becomes obsolete.

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Range of versions in all areas

The equipment variations installed in modern safety storage cabinets are as diverse as the demands on a modern laboratory

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Safety storage cabinets are not only available in different widths, heights and depths, with different door variants, they can also be equipped with a multitude of different interior fittings in an individual way so that they meet all laboratory requirements. Some manufacturers of single-wall cabinets specify the maximum storage capacity in litres or gallons. Although this corresponds to the internal volume of the cabinet, it does not say anything about the actual possible storage capacity. It is a theoretical value that is never achieved in practice. In the case of EN cabinets it is possible, using fully extractable drawers and taking into account customer-specific height settings, to optimally use the interior space without any waste.

When planning cabinets for the storage of combustible liquids, the following fittings are among those available:

- Shelves
- Fully extractable drawers
- Heavy-duty drawers
- Storage boxes
- Spill containment sumps
- Perforated metal plate inserts
- Insertable sumps
- Rubber mats
- Bottle racks
- Different door closing alternatives

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Natural gas cylinder cabinets can also be individually planned. For the storage of smaller compressed gas cylinders (for example 2 or 3-litre cylinders), fittings with several storage levels adapted to the cylinder size are possible.

The following fittings are available for gas cylinder cabinets:

- Rolling ramp, standard
- Rolling ramp, comfort
- Side cylinder holder
- Shelves
- Heavy-duty shelves
- Retaining chains
- Ratchet straps

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
Safety storage cabinets with shelves are the ideal choice when storing bulk containers. The shelves can easily be height-adjusted in standardized increments within the cabinet. Oftentimes only small quantities of a certain hazardous substance are required. In order to accommodate small bottles and containers in a most efficient and space saving way EN safety storage cabinets can be equipped with fully-extractable drawers.

Drawers can be installed with much less clearance than shelves because the containers can be accessed from the top as the drawers can be pulled out entirely. This way no space is wasted. In addition every drawer is a collecting sump thus making a bottom collecting sump obsolete.

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

Shelf available in powder-coated sheet steel or stainless steel

W x D x H (mm) 1043 x 435 x 30
(41.06 x 17.13 x 1.18 in)
Maximum load 75 kg
Increments for height adjustment 32 mm (1.26 in)

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

Drawer available in powder-coated sheet steel or stainless steel

W x D x H (mm) 865 x 510 x 90
(34.06 x 20.08 x 3.54 in)
Maximum load 60 kg
Retention volume 30 l
Max. size of stored container 27 l

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
Heavy-duty drawer available in powder-coated sheet steel or stainless steel

W x D x H (mm) 865 x 510 x 162
(34.06 x 20.08 x 6.38 in)
Maximum load 100 kg
Retention volume 66 l
Max. size of stored container 60 l

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asecos

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asecos

asecos – your partner for safety and environmental protection

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The name asecos stands for 'security' and 'ecology' – for safety and environmental protection in the handling and storage of hazardous materials. Since its establishment in 1994 the company, which is based in Gründau (in the German state of Hessen), has dedicated itself to the development and production of innovative products in this field.

Today asecos is Europe's leading manufacturer of type-tested safety storage cabinets and is well known for its trend-setting developments in the field of the storage and handling of hazardous substances.

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asecos has more than 1,700 technical options for the storage, safe transport, filling and refilling and handling of hazardous materials, both inside and outside buildings and in conformity with the law: The product range reaches from safety cabinets, cabinets for gas cylinders and hazardous material workplaces to extraction units, absorbents and sump systems through to hazardous material containers. Current developments encompass product solutions in the fields of non-smokers protection and cleaning room air.


Region-wide distribution network

Production is carried out at the headquarters in Gründau. Specialised trade dealers and partners handle marketing in Germany. Outside Germany asecos is represented by its own subsidiaries in the Netherlands, France, Spain and Great Britain; selected partner companies sell asecos products worldwide.

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The asecos think tank – non-stop trends

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asecos products are characterised by trend-setting, consumer-orientated solutions. We are of the opinion that safety has to be the overriding feature of our safety storage cabinets followed closely by convenience. Based on this conviction we develop intelligent, unique product solutions.

asecos holds more than 40 brand, industrial and design rights and makes considerable annual investments into research and development – that's how trends are created! The world's first type 90 safety cabinet was only one invention by asecos and was followed by many more. It caused a wave of innovations, setting the trend for future development of the safety levels for hazardous material storage.

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Safety and Environmental Protection – is our profession

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For asecos, the careful use of natural resources is as important as social obligation. State-of-the-art production processes sustainably preserve the environment. Local raw material suppliers guarantee not only fast delivery but also short delivery distances.

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asecos stands for quality 'made in Germany'

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asecos is a German manufacturer that always strives for the highest quality. That's why all our safety storage cabinets and smoker's cabins are produced at our headquarters in Germany. At the Gründau site motivated people invent and manufacture excellent products, which are widely distributed and used throughout Europe and beyond.

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Safe handling, approved storage and environmental protection:
asecos offers a wide range of safety products.

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Extraction of hazardous materials

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Handling of hazardous materials, cleaning, disposal & prevention

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Air purification concepts

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The asecos safety concept

You're always within safety limits with asecos' services

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Our expertise as developers and manufacturers of safety installations is also reflected in our range of services: We offer you a complete and professional range of services that leaves nothing to be desired.

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Why is the safety inspection so important?

The law specifies the regular inspection of safety-relevant installations. If safety installations are not inspected properly there is a risk of the loss of insurance cover in the case of a claim.

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The asecos academy

Die asecos academy

In our experimental lectures we sensitise you and your customers to the dangers posed by the improper storage of hazardous materials. Up close, but at a safe distance, we demonstrate deflagration, explosion and other chemical reactions with exciting experiments and our expert provides a gripping and informative introduction to the basic background knowledge. You experience the behaviour of flammable gases and liquids and learn about the various sources of ignition. Beyond that you learn how you can take protective and preventative action.

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47798

Promic

halogenated
Diethyl ether, Distilled

DANGER

University of Cambridge, United Kingdom

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Area of use: Whiffen Laboratory
Requirement: Replacement of the single-wall cabinets and optimal use of space
Solution: The dimensions of the cabinets were adapted to the available space; state-of-the-art type 90 cabinets

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Justus Liebig University, Giessen

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Area of use: Chemical store, Institute of Organic Chemistry
Requirement: Flexible organisation of the respective rooms and clearly arranged user-friendliness
Solution: 50 safety storage cabinets connected to the exhaust air system enable the uncomplicated storage and withdrawal of chemicals

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Johannes Kepler University, Linz

Johannes Kepler Universität, Linz

Area of use: Chemical store, laboratory, Institute of Chemistry
Requirement: State of the art; extraction from each cabinet must be ensured
Solution: Type 90 safety storage cabinets were installed with connecting lines and withdrawal stations

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Office of Chemical and Veterinarian Investigation, Krefeld

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Area of use: Preparation and analysis of food samples
Requirement: Modernisation of the analysis laboratory, protection of health and ergonomic working
Solution: Hazardous substance workplace with integrated supply

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GlaxoSmithKline Biologicals, Dresden

GlaxoSmithKline Biologicals, Dresden

Area of use: Laboratory area
Requirement: Modernisation of the existing laboratory equipment
Solution: Equipment with hazardous substance workplace, underbench cabinet, refrigerated underbench cabinet and combination cabinet with type 90 compartment and acid/alkali compartment

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Area of use: Outdoor area
Solution: Use of LG cabinets for the safe storage of compressed gases in the outdoor area

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Universities in Germany (e.g. Munich, Frankfurt, Jena, Chemnitz, Dresden, Stuttgart, Potsdam, Regensburg, Hohenheim, Ulm ...)

Universities in Spain (e.g. Barcelona, Almería, Madrid, Alicante, la Coruña, Cádiz, León, Murcia, Sevilla ...)

Universities in UK (e.g. Newcastle, Bradford, Glasgow, Manchester ...)

Summary and conclusion

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Summary and conclusion

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A modern laboratory without hazardous materials is inconceivable – just as these chemicals, compressed gases, acids or alkalis are inconceivable without the dangers they pose. However, you can reduce the risks associated with the storage of these hazardous materials to a minimum with safety storage cabinets designed and fire-chamber-tested in accordance with European standards EN 14470 Part 1 and Part 2!

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We hope this brochure has given you a first insight into the European concept of storing hazardous materials in modern laboratories. The technologies in the lab continue to advance very quickly and so we constantly strive for innovations to meet these new challenges. It would be a pleasure for us to provide further information on the subject and to discuss about the solutions asecos can provide. Please do not hesitate to get in touch with us!

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asecos – Knowledge Base

asecos – Wissen kompakt!

Oops, the DVD 'asecos – Knowledge Base' is missing?
Then simply request your own personal copy:

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Phone +49 60 51 92 20-0
Fax +49 60 51 92 20-10
info@asecos.com

asecos Ltd.

Safety and Environmental Protection
c/o Park Associates
Gretton House, Waterside Court
Third Avenue, Centrum 100
Burton-upon-Trent, Staffordshire
GB-DE14 2WQ

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Veiligheid en milieubescherming
Tuinderij 15
NL-2451 GG Leimuïden

Phone +31 172 50 64 76
Fax +31 172 50 65 41
info@asecos.nl

For all other countries please contact:

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asecos GmbH