

Protocol - Flow Cytometric Analysis of Dissociated Tumor Cells

Reagents:

- **FACS Buffer** (PBS + 2% FBS, or similar)
- **Fc Blocking Solution** (BioLegend Human TruStain FcX, or similar)
- **Fluorescently-Tagged Antibodies**
- **Viability Dye** (Propidium Iodide, DAPI, 7AAD, or similar)
- **Conversant Bio Dissociated Tumor Cells** (DTCs)
 - All Conversant Bio Viable Cell Products are shipped on dry ice and should be used immediately or stored in liquid nitrogen. For more information, please refer to our document [Recommended Handling and Storage for Biological Specimens](#).

Equipment:

- **Flow Cytometer**
- **Centrifuge**
- **Micropipettes**
- **Pipettors**

Protocol:

1. Thaw DTC samples as described in [Thawing Viable Cell Products](#).
2. Resuspend cells at $0.5-1 \times 10^6$ cells/mL in FACS Buffer.
3. For each stain, aliquot 0.5-1 mL cell suspension into 5 mL round bottom tubes.
 - a. For rare cell populations, staining greater than 1×10^6 cells may be required.
4. Centrifuge at 300xg for 5 minutes at room temperature.
5. Decant the supernatant and gently wick away any residual buffer.
6. Resuspend cells in Fc Blocking Solution according to the manufacturer's protocol.
 - a. For BioLegend Human TruStain FcX, prepare a stock solution of 25 μ l FACS Buffer + 2.5 μ l Human TruStain FcX per sample to be stained.
 - b. Resuspend cells in 25 μ l Human TruStain FcX Solution.
 - c. Incubate 15 minutes at room temperature.
7. Stain cells with fluorescently-tagged antibodies.
 - a. If multiple samples are to be stained, generate a master mix with 25 μ l FACS Buffer + the antibodies of interest per sample.
 - i. All antibodies should be titrated to ensure optimal staining concentration.
 - b. Add 25 μ l staining master mix to each sample. Mix gently.
 - c. Incubate 30-45 minutes at 4°C.
8. Add 1 mL FACS Buffer to each tube.
9. Centrifuge at 300xg for 5 minutes at room temperature.
10. Decant the supernatant and gently wick away any residual buffer.
11. Resuspend cells in 300 μ l FACS Buffer containing a viability dye at the recommended concentration.
12. Proceed to flow cytometry analysis.