



Ghost Peaks: **How to Fix a Haunting Problem**

Jacob A. Rebholz
Teledyne Tekmar
VOC Product Line Manager




TELEDYNE TEKMAR
Everywhere you look™

Outline of Topics

- What is a “ghost” peak?
- How to identify ghost peaks and their source
- How to correct various sources of contamination
- How to prevent future “hauntings”!


What is a “Ghost” Peak?

- Simply put: any peak that does not belong in your chromatogram.



What is a “Ghost” Peak?

- Simply put: any peak that does not belong in your chromatogram.




What is a “Ghost” Peak?

- Less simply put: contamination can present in several ways
 - Elevated baseline
 - Contribution to target peaks
 - Tailing/ “splitting”

What is a “Ghost” Peak?

- Less simply put: contamination can present in several ways
 - Elevated baseline
 - Contribution to target peaks
 - Tailing/ “splitting”




TELEDYNE TEKMAR
Everywhere you look™

What is a “Ghost” Peak?


- Less simply put: contamination can present in several ways

- Elevated baseline
- Contribution to target peaks
- Tailing/ “splitting”



What is a “Ghost” Peak?

- Less simply put: contamination can present in several ways
 - Elevated baseline
 - Contribution to target peaks
 - Tailing/ “splitting”



TELEDYNE TEKMAR
Everywhere you look™

How to Identify a “Ghost” Peak

- Visually detect stray peaks in chromatogram
 - Helps to have an example of a clean chromatogram, or old chromatograms from when the system was clean
 - Helps identify extra peaks, elevated baselines and potentially tailing/split peaks
- Poor calibration results
 - Can help ID contribution to peaks and split/tailing peaks

Troubleshooting Tips

- Take a moment to review manuals
- Examine the facts and use valid reasoning
- Identify the “root cause” of the problem
- Avoid quick fixes and shotgun approach


BE PATIENT!



TELEDYNE TEKMAR
Everywhere you look™

Exorcising the Demon

- Step 1: locate the source of the contamination
 - Break system down into major components:
 - Autosampler components, Detector, GC inlet, GC column



A Word on Method Parameters

■ GC cycle time

- Heavy compounds missing from the end of chromatogram
- Show up in following blank
- Fix: increase final temperature and/or hold time

A Word on Method Parameters

- Inlet temperature
 - Typically heavy compound carryover, potentially others
 - Check temperature of inlet and/or any transfer lines
 - Fix: ensure temperature is adequate to prevent condensation of compounds, increase if necessary

Exorcising the Demon

- System blanks
 - Run a full system blank
 - Run GC only blank (no autosampler)
 - Run Detector only Blank (no column connected if possible)

Solvents/Standards

■ Solvent

- Methanol, acetone, etc.
- Direct inject of solvent used to prepare standards
- Compare results to a second source of solvent
- Don't forget that syringe could contribute
- Second GC system to confirm

Solvents/Standards

- Solvent ruled out
- Standards
 - Crack new ampules of internal and surrogate standards
 - Direct inject
 - Crack new ampules of calibration standard
 - Direct inject
 - Do not forget syringes as source of contamination

Carrier Gas

- Gas filters
- Gas supply tubing
- Gas regulators
- Carrier gas supply



<http://www.restek.com/catalog/view/3445/22019>



Carrier Gas

- Allow GC oven to sit at or below starting temperature for up to several hours
 - Run GC blank (no injection)
 - If carrier gas contamination is focusing on cold column, peaks will be larger
 - If no change in peak size, continue with trouble shooting

Carrier Gas

- Verify filter is not expired (color change indicator)
 - Check for contamination
 - Replace if uncertain
- Gas supply
 - Tank tubing and regulator
 - Replace all, verify issue resolved, then isolate issue to source
 - If gas contaminated, other components may be contaminated as well

Detectors

■ Mass Spectrometers

- Check the tune report, and/or air-water report
- Is the detector working properly?
- Is there sufficient vacuum?
- Is the source clean?



Detectors

- Other detectors (ECD, TCD)
 - Remove column and cap detector
 - Run a blank and compare detector response to that of blank with column installed
 - Improved: contamination from GC or autosampler
 - Same: contamination from make up gas, detector itself

Inlet Maintenance

- Detector has been ruled out
- Replace inlet liner, septum
- Inspect seals for damage/discoloration and replace
- Clip the column
- Verify correct column depth and re-install
- Re-run system blank



Other Inlet Issues

■ Split ratio

- Verify that split ratio valve is functioning properly
- Sticking or non-opening valve can cause overloading of inlet, column and detector

■ Septum purge

- Faulty septum purge can allow septum off-gassing to reach detector.
- Tubing can collect contamination

Column Issues

- Detector and inlet issues ruled out
- Column clipped during inlet maintenance
- Swap column with a new one or other column of similar phase
 - Issue resolved, replace column

Autosamplers

- Liquid Injection
- Headspace
- Purge and Trap



Liquid Injection Autosamplers

- Syringe
 - Replace/clean- verify no change
- Cleaning solvent
 - Verify adequate number of cleans between samples with uncontaminated solvent
- Vials
 - Verify contamination-free sample vials



Headspace Autosamplers

- Needle
- Transfer line
- Sample loop
- Vials
- 6-port valve



Headspace Autosamplers

- Vials
 - Verify clean, compare to a different lot/manufacture
- Needle
 - Often the source of contamination
 - Clean/replace and re-run
- Trap
 - In dynamic systems, replace trap and re-run to verify

Headspace Autosamplers

- Transfer line, valve oven, sample loop
 - Verify temperatures sufficient to prevent condensation of compounds
 - If condition persists, active site or more persistent contamination

Purge and Trap Systems

- Two main components

- Concentrator

- Sparger, sample mount, analytical trap, condensate trap, 6-port valve, transfer line

- Autosampler

- Sample needle, syringe/sample loop, transfer tubing, standards injector



TELEDYNE TEKMAR
Everywhere you look™

Concentrator

- Typically carryover is just late eluters, but contamination can be anything
- Check temps during Desorb and Bake
 - Ensure desorb of trap and cleanup between samples

Concentrator

- Desorb only
 - Rules out gas pathway
 - If contaminated, change trap and re-run
 - Still have issues, likely an active site in gas pathway, including tubing, 6 port, transfer line, and/or condensate trap

Concentrator

- Disconnect autosampler/bypass autosampler functions
- Liquid pathway
 - Sparger, sample introduction valve, drain valve
- Sparger
 - Check for discoloration. Clean or replace, then run blank
 - Check drain valve for proper function

Autosampler

- Check hot water and/or MeOH Rinse
 - Ensure cleanup of liquid path between samples
- Possible sources of contamination
 - Sample vials
 - Multi-stage needle
 - Syringe/sample loop
 - Standard addition system



Autosampler

■ Soil autosampler

- Contamination in water mode or soil mode? Run each to verify
- Soil mode only: needle, soil valve, transfer line
- Water mode only: liquid handling pathway
- Both: rinse water, syringe, needle



Autosampler

- Syringe based system
 - Bypass needle by placing transfer tubing directly into vial of water
 - Clean/replace syringe
 - Bypass standard addition system
 - Still have issues, multi-port valve

Autosampler

- Loop and valve system
 - Verify loop rinse times, ensure adequate to rinse loop between samples
 - Clean/replace loop
 - Bypass standard addition system
 - Clean/replace transfer tubing



Preventative Actions

- Maintain proper methods
 - Before making method changes, ask self “Why?”
 - Retain records of changes made and why
 - Keep a copy of last known “good” method
 - If its not broke, don’t fix it!



Preventative Actions

- Regular preventative maintenance
 - GC maintenance
 - Inlet maintenance, clip column, regular detector cleaning/tuning
 - Column life
 - Condition columns properly and replace on a regular basis as needed
 - Don't wait until problems arise, waste valuable time troubleshooting!

Preventative Actions

- Autosampler system maintenance
 - Clean/replace syringes, rinse solvents/water and other consumables on a regular basis
 - Analytical traps
 - Same as columns, condition properly and replace regularly as needed
 - Monitor purge and bake pressures, spike in pressure may signal time to replace

Preventative Actions

■ Consumables

- Purchase standards, solvents and vials from reputable manufacturers
- Refresh standards and solvents regularly to prevent contamination
- Use good practice when cleaning syringes, flasks and other glassware

Preventative Actions

- Know your samples
 - Consider pre-screening for P&T samples
 - Headspace, prevents contamination from high conc. samples
 - Consider in-vial purge for difficult matrices
 - Consider headspace to prevent column fouling
 - Increase maintenance frequency if/when necessary

As Always, We're Here to Help!

- Don't let troubleshooting scare you, we're here to help!

- Tech Support 1-800-874-2004
- Tekmar_IntlTech@Teledyne.com
- tekmarsupport@teledyne.com
- www.teledynetekmar.com
- Jaocb.Rebholz@Teledyne.com



TELEDYNE TEKMAR
Everywhere you look™