

JC

HIGH NOBLE YELLOW TYPE 3 CROWN & BRIDGE ALLOY

JC was one of the first alloys developed by Jensen Dental more than 35 years ago. A high noble, yellow crown and bridge alloy, JC is perfect for most cast-gold restorative work including single and partial crowns, onlays, and short-span bridges. JC offers a stunning 18 carat yellow color and is well known by industry professionals to be a strong alloy that melts, casts, finishes and polishes with ease and accuracy.

PROPERTIES		CHEMISTRY	
Melting Range	1670 to 1760°F (910° to 960°C)	Gold	75%
Density	15.2 g/cm ³	Silver	11%
Grain Size	24 microns	Copper	10%
Hardness	125 HV	Palladium	3 %
Tensile Elongation	50%	Contains less than 1% Zinc, Indium, Iridium	
Tensile Yield Strength	39,500 psi (270 MPa)	Au & Pt group - 78%	
Ultimate Tensile Strength	62,000 psi (430 MPa)	Classification - High Noble	

PROCESSING TECHNIQUE

SPRUIING

The indirect method is recommended for multi-units. Use an 8 gauge runner bar with 10 gauge connectors. If preferred, the direct method may be used on both single units and small bridges. Use a 10 gauge sprue 1/4" (6mm) to 3/8" (9mm) long. Sprues longer than 3/8" (9mm) should have a reservoir 1/16" (1.5mm) from pattern. Patterns should be a maximum of 1/4" (6mm) from top of investment.

INVESTMENT AND BURNOUT

Either gypsum or phosphate bonded investment may be used following the manufacturer's instructions. The burnout temperature should be at least 900°F (480°C) and should not exceed 1200°F (650°C).

MELTING AND CASTING

Extra winds of the casting arm are not required. A gas/compressed air or gas/oxygen flame with 5 psi gas and 10 psi oxygen is recommended. The alloy will fully puddle and form a ball before it is ready to cast. The casting temperature is 1850°F (1010°C). Bench cool to a dull red heat before quenching.

DEVESTING AND FINISHING

Blast with aluminum oxide to remove investment particles and oxidation. Finish and polish using standard techniques.

SOLDER AND FLUX

Solder: 650 Fine Solder
Flux: Brown Fluoride Flux

5003Y r1

