

JP-II

HIGH NOBLE YELLOW PORCELAIN ALLOY

JP-II is a high noble, yellow porcelain dental alloy that is perfect for just about any type of restoration, from single units to small-span bridges. Silver-free, JP-II offers a pleasant-yellow color with a very light oxide for superior shade control. The opaque layer will cover the coping easily and completely. JP-II is easier to polish than most yellow gold porcelain alloys and finishes to a fine luster.

PROPERTIES	
Melting Range	1990° to 2190°F (1085° to 1200°C)
Coefficient of Thermal Expansion	
from 25°C to 500°C:	14.2x10-6C ⁻¹
from 25°C to 600°C:	14.4x10-6C ⁻¹
Density	18.5 g/cm ³
Grain Size	15 microns
Hardness	155 HV
Tensile Elongation	17%
Tensile Yield Strength	50,700 psi (350 MPa)
Ultimate Tensile Strength	64,200 psi (440 MPa)

CHEMISTRY	
Gold	87%
Platinum	7%
Palladium	4.5%
Contains less than 1% Tin, Indium, Iron, Rhenium	
Au & Pt group – 98.5%	
Classification - High Noble	

PROCESSING TECHNIQUE

- WAXING** Wax to a minimum thickness of .3mm for single units and .5mm for bridge work. Avoid sharp angles and corners.
- SPRUIING** The indirect method is recommended for multi-units. Use an 8 gauge runner bar with 10 gauge connectors. If preferred, the direct method may be used on both single units and small bridges. Use a 10 gauge sprue 1/4" to 3/8" long. Sprues longer than 3/8" should have a reservoir 1/16" from pattern. Patterns should be a maximum of 1/4" (6mm) from top of investment.
- INVESTMENT** A phosphate-bonded, high heat investment with or without carbon content is recommended.
- BURNOUT** 1350°F (735°C)
- MELTING AND CASTING** Extra winds of the casting arm are not required. Use a multi-orifice torch with 10 psi fuel and 20 psi oxygen. The alloy will fully puddle and form a ball before it is ready to cast. **DO NOT OVERHEAT. DO NOT USE CASTING FLUX.** The casting temperature is 2290°F (1255°C).
- DEVESTING AND FINISHING** Blast with aluminum oxide to remove investment particles. Finish with carbides or aluminum oxide stones (Refer to our *Finishing and Polishing Guide – Yellow Ceramic Alloys* for more information). Reblast porcelain receiving surface with nonrecycled aluminum oxide. Clean in ultrasonic for 10 minutes in distilled water.
- OXIDATION** Oxidize from 1200°F (650°C) to 1850°F (1010°C) at 145°F/min (80°C/min) in air with 5 minutes hold time. Bench cool. Proceed with normal opaque technique.
- SOLDERS AND FLUX** Pre-Solder: Spirit Solder or LX Solder
Post-Solder: 615 Solder
Flux: Brown Fluoride Flux for both pre and post soldering

5015Y r2

