

JP-I

HIGH NOBLE WHITE SILVER FREE PORCELAIN ALLOY

Our interpretation of the classic 52% gold, silver free porcelain dental alloy, **JP-I** has been one of Jensen's most popular PFM alloys for over 30 years. Suitable for all types of restorations, our customers continually tell us that JP-I offers sharper, more stable margins than competitive alloys of this type. JP-I gives the technician a moderate gray, uniform and consistent oxide layer, reliable bonding, and excellent cast ability. JP-I is used with virtually all conventional dental PFM ceramics, but prefers low to moderate expansion porcelains that tend to be silver sensitive. A great alloy for Pulse press-to-metal.

PROPERTIES	
Melting Range	2075° to 2220°F (1135° to 1215°C)
Coefficient of Thermal Expansion	
from 25°C to 500°C:	13.8x10-6C ⁻¹
from 25°C to 600°C:	14.1x10-6C ⁻¹
Density	14 g/cm ³
Grain Size	38 microns
Hardness	240 HV
Tensile Elongation	25%
Tensile Yield Strength	73,000 psi (505 MPa)
Ultimate Tensile Strength	110,000 psi (760 MPa)

CHEMISTRY	
Gold	51.5%
Palladium	38.5%
Indium	8.5%
Gallium	1.5%
Contains less than 1% Ruthenium	
Au & Pt group - 90%	
Classification - High Noble	

PROCESSING TECHNIQUE

- WAXING** Wax to a minimum thickness of 0.3mm for single units and 0.5mm for bridge work. Avoid sharp angles and wax to provide for an even thickness of porcelain.
- SPRUIING** The indirect method is recommended for multi-units. Use an 8 gauge runner bar with 10 gauge connectors. If preferred, the direct method may be used on both single units and small bridges. Use a 10 gauge sprue 1/4" (6mm) to 3/8" (9mm) long. Sprues longer than 3/8" (9mm) should have a reservoir 1/16" (1.5mm) from pattern. Patterns should be a maximum of 1/4" (6mm) from top of investment.
- INVESTMENT** A phosphate-bonded, high heat investment with or without carbon content is recommended.
- BURNOUT** 1500°F (815°C)
- MELTING AND CASTING** Wind casting arm one turn more than used for casting gold. Use a multi-orifice torch with 10 lbs. gas and 20 lbs. oxygen. As JP-I melts, a cloudy surface will appear. Continue heating until the cloudy surface clears, before releasing the casting arm. **DO NOT OVERHEAT.** The casting temperature is 2300°F (1260°C). **DO NOT USE CASTING FLUX.**
- DEVESTING AND FINISHING** Blast with aluminum oxide to remove investment particles. Shape and finish down metal with aluminum oxide stones. Blast outer surface with non-recycled aluminum oxide (50 micron-white preferred). Clean in ultrasonic for 10 minutes in distilled water.
- OXIDATION** Oxidize from 1200° to 1850°F (650° to 1010°C) at 145°F/min (80°C/min) in air. Hold for 5 minutes at 1850°F (1010°C) in air. Bench cool. Proceed with opaque following porcelain manufacturer's instructions.
- SOLDERS AND FLUX** Pre-Solder: Spirit Solder or LX Solder
Post-Solder: 1400 Solder
Flux: Brown Fluoride Flux for both pre and post soldering

5014Y r2

