


WHEELCHAIR SEATING & POSITIONING GUIDE




MEET THE COMFORT COMPANY CLINICAL EDUCATION TEAM



Stacey Mullis OTR/L, ATP
Director of Clinical Education



Ana Endsjo MOTR/L, CLT
Clinical Education Manager
LTC Division



Lindsey Sharpe PT, DPT
Clinical Education Manager
Rehab Division

COMFORT COMPANY MISSION

Comfort Company is committed to developing high quality resources and education related to wheelchair seating and positioning. We understand the challenges that healthcare professionals face in the provision of wheeled mobility and positioning. Our goal is to simplify the process and provide the tools needed to help you meet the needs of end users.



PURPOSE OF THIS GUIDE


This guide is intended to demystify the provision of wheelchair seating and positioning for healthcare professionals. The aim is to assist healthcare professionals in achieving optimal individualized wheelchair seating to meet the goals for each client by:

- Understanding the qualifications for both manual and power wheelchairs
- Identifying the abnormal postures and understanding how ill-fitting wheelchair components lead to abnormal postures and pressure injuries
- Understanding how to select the proper seating components
- Identifying the key components to completing a successful wheelchair evaluation

*** Look for these blue boxes throughout the guide. They're our quick tips or takeaways for that section.**



This guide uses common therapy terms to keep the concepts as relatable as possible. For the most up-to-date, standardized wheelchair terms and definitions, see reference #9, *Glossary of Wheelchair Terms and Definitions*, as listed on page 88.



* Despite wanting what is best for our patients, the above factors may prevent us from doing it.

TABLE OF CONTENTS

GETTING STARTED

The Process of Getting Seating Equipment	3
DME vs CRT	5
The Team	6
Identifying the Need	8

MANUAL WHEELCHAIRS

Manual Wheelchair Universal Terms	9
Manual Wheelchair Justification	10
Wheelchair Options: Features & Limitations	12
Standard Wheelchair Comparison Chart	15

POWER WHEELCHAIRS

Power Wheelchair Justification	17
Reasons to Choose Group 3 Power	20

SEATING AND POSITIONING

Looking at Abnormal Postures	21
Fixed vs Flexible Postural Abnormalities	30
Issues with the Current Wheelchair System Can Cause Abnormal Postures	31
The Therapy Evaluation: The Pieces Required	43
Measuring for Properly Fitting Wheelchair Components	45
Choosing a Wheelchair Seat Cushion	56
When to Choose a Separate Off-The-Shelf Back Support	70
Choosing a Back Support	72
Accessories	78

COMMON MISCONCEPTIONS

The Elevating Leg Rest Myth	81
The Wedge Misconception	84
The Problem of Sliding Out of the Wheelchair	85

CONCLUSION

87

GLOSSARY

88


THE PROCESS OF GETTING SEATING EQUIPMENT




WHERE DO WE START?



First, let's look at the big picture of how to get wheelchair seating and positioning equipment:





*** This guide will give you tools to be able to confidently walk through this process.**

DME VS CRT



**LET'S START WITH THE BASICS:
IS ALL WHEELED MOBILITY AND
SEATING EQUIPMENT CLASSIFIED
THE SAME?**



NO! When it comes to seating and wheeled mobility, products are divided into two groups: Durable Medical Equipment (DME) and Complex Rehabilitation Technology (CRT).

Patients will qualify for certain equipment based on the severity or complexity of their condition. Most products are covered based on diagnosis, and some are covered because lesser products have been tried and were found to be inappropriate.

"Complex Rehab Technology (CRT) products and services are significantly different than standard Durable Medical Equipment (DME). The DME benefit was created over forty years ago to address the medical equipment needs of elderly individuals. Over the years technology has advanced and now includes complex rehab power wheelchairs, highly configurable manual wheelchairs, adaptive seating and positioning systems, and other specialized equipment. These products are classified as Complex Rehab Technology. Suppliers who furnish CRT provide highly specialized products and services which are much different than standard DME."

NCARTblogger. "Why a Separate Benefit Category Is Necessary for CRT." NCART, NCART, 25 June 2015, blog.ncart.us/why-a-separate-benefit-category-is-necessary-for-crt/.

THE TEAM




**IS IT ONLY UP TO THE PHYSICIAN
AND/OR THERAPIST TO DECIDE
ON A PATIENT'S EQUIPMENT?**



NO! The most important part of successfully providing appropriate equipment is understanding that it's a team effort.


MEET THE TEAM: Each member plays a role in the selection and attainment of proper equipment.




Client/Caregiver



Physician




Nurse



Therapist



ATP/Supplier/Dealer



Manufacturer



What does each team member bring to the table?

CLIENT/CAREGIVER

Has knowledge of their own body, functional needs, what works best for them in their day-to-day lives

PHYSICIAN

Determines the medical need for seating and wheeled mobility

THERAPIST

Provides clinical knowledge of body, tone, debility, stability vs mobility, function vs exercise, and patient advocacy

SUPPLIER/DEALER

Has knowledge of currently available equipment, billing/insurance issues, qualification requirements

NURSE

Adds wound care expertise

ATP - ASSISTIVE TECHNOLOGY PROFESSIONAL

Must go through an in-depth certification process

Has specialized knowledge of complex rehab products/equipment, qualification requirements, is competent in analyzing the needs of consumers with disabilities, aids in the selection of appropriate assistive technology for the consumer's needs, and provides training in the use of the devices

MANUFACTURER

Offers knowledge of their products, clinical applications, integration of their products with other technology, and all the pros/cons of products for various patient presentations



It is critical for health care professionals to develop a relationship with a reputable dealer who employs an ATP.

IDENTIFYING THE NEED



HOW DO I KNOW THAT MY PATIENT NEEDS WHEELED MOBILITY?




Consider your patient's quality of life. The list below provides good indicators that your patient would benefit from wheeled mobility.

1. Patient is **non-ambulatory**
2. **Demonstrates decreased safety** with ambulation or is **at risk for falls** within the home. Ask about history of falls; perform an objective balance assessment i.e. BERG, DGI
3. Requires **assistance for ambulation** within the home and wheeled mobility would allow independence
4. Requires **increased time for ambulation** within the home. Perform a gait speed test; think about performing ADLs in a reasonable amount of time
5. **Unable to consistently ambulate throughout the day** in the home which affects their ADLs. Look at a 24 hour period

* Here are scenarios where wheeled mobility could significantly increase a person's quality of life:

- They can ambulate but are at high risk of falls
- They have frequent urge incontinence because they are unable to get to the restroom on time
- Their O₂ saturations drop below or heart rate increases above a safe range with ambulation
- Their day consists of sitting in a recliner and transferring to a bedside commode as needed

MANUAL WHEELCHAIR UNIVERSAL TERMS



MANUAL WHEELCHAIR JUSTIFICATION




NOW THAT I KNOW MY PATIENT WILL BENEFIT FROM A WHEELCHAIR, WHAT DO I DO?



Wheelchair selection requires evidence of medical necessity. Step one is a physician's visit with notes that state:

- Mobility related diagnosis
- Symptoms that affect mobility
- MRADLs affected by the mobility limitation
- Current ambulation status

Then, a referral is made to PT/OT and the fun begins! It becomes our task to evaluate the patient and determine what level of wheeled mobility they need to lead safe, functional lives.





HOW DO I JUSTIFY MY CLIENT'S NEED FOR A WHEELCHAIR?



Prior to choosing the level of wheelchair, the PT/OT needs to justify the need for a manual wheelchair. Ask yourself the following questions, and the answers will begin to guide you towards the right wheelchair:

1. Does your patient have a mobility limitation that significantly impairs his/her ability to participate in **one or more MRADLs** in the home?

*Does it **prevent** them from doing MRADLs?*

*Are they **unsafe** to perform MRADLs?*

*Can they perform MRADLs in a **reasonable time frame**?*

2. Can the mobility limitation be resolved by a cane or walker?

3. Do they have the desire or capability to propel a wheelchair?

If they can't propel, do they have a willing caregiver?

4. Does the patient's home have the space/layout for functional wheelchair use?

Measure doorways and ask your ATP for required measurements to get through doorways based on the wheelchair model selected

*** Always document how the right equipment allows them to perform routine tasks more independently.**

WHEELCHAIR OPTIONS: FEATURES & LIMITATIONS

? NOW THAT I'VE IDENTIFIED THE NEED, AND KNOW HOW TO JUSTIFY A MANUAL WHEELCHAIR, WHICH DO I CHOOSE?

STANDARD WHEELCHAIR OPTIONS			
Very Minimal Adjustability	Minimal Adjustability	Most Basic Adjustability	Most Customizable
K0001/K0002	K0003	K0004	K0005

WHEELCHAIR KEY

DME

K0001: Standard Wheelchair

K0002: Standard Hemi Height Wheelchair

K0003: Lightweight Wheelchair

K0004: High Strength, Lightweight Wheelchair

CRT

K0005: Ultra Lightweight Wheelchair

Tilt-in-Space: Dependant Manual Wheelchair



The first question to ask is: Will my patient need this **short term** or **long term**?



SHORT TERM USE: K0001 - K0003

- These chairs are heavier, basic manual wheelchairs with minimal adjustability and are most appropriate for short term use.

LONG TERM USE: K0004 - K0005 - TILT K0004

- Patient qualifies for a basic manual chair, but requires a seat width, depth, or height that can't be accommodated by a basic MWC
- They are unable to perform all ADLs efficiently in a basic MWC
- They use a wheelchair more than two hours per day



K0005

- Patient is a full time wheelchair user
- Patient requires customization such as axle configuration, wheel camber angle, or seat-to-back angle that can't be accommodated by a K0001-K0004
- This requires an evaluation by a PT/OT, a letter of medical necessity, and the involvement of an ATP in the equipment selection process



MANUAL TILT-IN-SPACE



WHEN IS TILT-IN-SPACE APPROPRIATE?

- Patient is dependent in mobility
- Patient is unable to perform independent pressure relief
- Patient requires gravity-assisted positioning/repositioning
- Patient requires postural support, head and trunk control, and accommodation of postural asymmetries
- The goal is to increase sitting tolerance/endurance
- Patient needs improved line of sight due to forward head posture
- Patient will benefit from trunk support and open thoracic posture for increased respiratory function
- Patient requires safe positioning for feeding/gravity-assisted swallowing

*** A manual tilt-in-space chair requires an PT/OT evaluation, justification that other manual wheelchairs are not appropriate, and have an ATP involved in the process.**

STANDARD WHEELCHAIR COMPARISON CHART

	K0001: Standard	K0002: Standard Hemi Height	
Dimensions			
Chair Weight without Leg Rests	>35 lbs	>35 lbs	
Seat Width: Standard	16", 18", 20"	16", 18", 20"	
Seat Depth: Standard	16"	16"	
Weight Capacity	300 lbs	300 lbs	
Back Height	18"	18"	
Lowest Achievable Seat-to-Floor Height	21"	19"	
Adjustability to accommodate for postural abnormality:			
Arm Rest Height	No	No	
Back Height	No	No	
Seat-to-Back Angle	No	No	
Ability to create a Fixed Tilt	No	No	
Seat-to-Floor Height (STFH)	No	Hemi Height Only	
Arm Rest Options			
Desk Length	Yes	Yes	
Full Length	Yes	Yes	
Leg Rest Options			
Standard	Yes	Yes	
Elevating (ELR)	Yes	Yes	
Swing-Away	Yes	Yes	
Meant for Long Term Sitting	No	No	

**Features vary according to model.*

	K0003: Lightweight	K0004: High Strength, Lightweight	K0005: Ultra Lightweight	Tilt-in- Space
	33 - 35 lbs	30 - 34 lbs	<30 lbs	>45 lbs
	16", 18", 20"	16", 18", 20", 22"	Customizable	16", 18", 20"
	16", 18"	16", 18", 20"	Customizable	16", 18", 20"
	300 lbs	300 lbs	Customizable	300 lbs
	18"	16" to 20"	Customizable	24"
	17"	13.5" (most 14.5")	Customizable	17"
	No	Yes	Yes	Yes
	No	Yes	Yes	Yes
	No	Yes	Yes	Yes
	No	Yes	Yes	Yes
	Yes	Yes	Yes	17" - 21"
	Yes	Yes	Yes	Yes
	Yes	Yes	Yes	Yes
	Yes	Yes	Yes	Yes
	Yes	Yes	Yes	Yes
	Yes	Yes	Yes	Yes
	No	Yes	Yes	Yes

POWER WHEELCHAIR JUSTIFICATION



WHY CHOOSE A POWER WHEELCHAIR?



All of the basic criteria for a manual wheelchair apply **AND** they must meet the basic power chair criteria below:

1. They are unable to propel a manual wheelchair due to upper extremity limitations:
 - Strength
 - Coordination
 - Pain
 - Range of motion
2. The home has adequate access for maneuvering of the power chair
3. Use of the power wheelchair will significantly improve their ability to perform MRADLs

*** These all require a physician face-to-face appointment, an OT/PT evaluation, and an ATP directly involved in the equipment choice.**

SCOOTER

Meets all of the basic power wheelchair criteria

SCOOTERS VS HIGHER LEVEL BASES:

A scooter will be often ruled out due to:

- Poor trunk strength
- Difficulty with transfers
- Poor endurance
- Existing postural abnormalities
- Large turning radius

*** Very rarely does a scooter fit in the home so a Group 2 power wheelchair is the next base to consider!**

GROUP 2 POWER

Meets all of the basic power wheelchair criteria and:

- Group 2 models are typically characterized by "captain's seating"
- The patient is unable to safely transfer, operate, and maintain postural stability in scooter
- The home does not provide adequate access for operating scooter
- The patient is able to safely operate a power wheelchair OR has a caregiver willing to operate, who cannot push a manual wheelchair
- Group 2 will improve the ability to participate in MRADLs in the home
- Some Group 2 models have multi power options available if they qualify for a tilt/recline system, or if they use a ventilator that is mounted on the chair



GROUP 3 POWER



Meets all of the criteria for a basic power wheelchair and Group 2 power.

Additional Criteria:

- The patient has a neurological condition, myopathy, or congenital skeletal deformity

What's different than a Group 2?

- Group 3 allows integration of rehab seating, and typically has better drive wheel suspension
- Group 3 models offer multiple alternative drive control options
- Group 3 models offer the option of power tilt, recline, and/or power legs, and power elevate
- Tighter turning radius and increased speeds

REASONS TO CHOOSE GROUP 3 POWER



WHY WOULD MY PATIENT NEED A GROUP 3 RATHER THAN A GROUP 2?



Neurological conditions are the main reason a patient would need a Group 3 power chair:

- Neurological conditions are often associated with abnormal tone and/or reflexes. Drive wheel suspension becomes important to minimize jarring forces that can trigger tone or make the user unstable
- When using a power chair all day, the distance per charge is important
- Group 3 suspension will provide the safest navigation over uneven terrain, and minimize jarring forces
- If the patient has impaired sensory-motor function, they may not sense a Group 2 base tipping during an incline or decline, which puts them at risk

LOOKING AT ABNORMAL POSTURES



HOW DO I KNOW IF MY PATIENT NEEDS A CUSHION OR BACK SUPPORT?



Look at the patient's posture and ask: what is the body doing that it shouldn't be doing?

With prolonged sitting, patients begin to shift their bodies into what is known as abnormal postures to seek stability and/or alleviate pain and pressure.


Abnormal Postures are grouped into 5 categories and are as follows:

- **Anterior Pelvic Tilt** (page 24) with lumbar lordosis with or without neck hyperextension.
- **Posterior Pelvic Tilt** (page 25-26) with thoracic kyphosis with or without forward neck flexion. Referred to as Sacral Sitting.
- **Pelvic Obliquity** (page 27) with scoliosis with or without lateral neck flexion.
- **Pelvic Rotation** (page 28) with rotation of the spine with or without lateral neck flexion.
- **Windswept Posture** (page 29)


*** Understand what you are looking at and decide whether your goal is to *correct* or *accommodate* for the abnormal posture. (Page 30)**

ANATOMY REFRESHER


RIGHT LATERAL VIEW OF SPINE



RIGHT LATERAL VIEW OF PELVIC BONE



POSTERIOR VIEW OF PELVIS



NEUTRAL PELVIC POSTURE



NEUTRAL PELVIC POSTURE

What is the pelvis doing?

Pelvis in midline.
ASIS & PSIS at equal height: no pelvic tilt
L ASIS & R ASIS at equal height: no obliquity
L ASIS & R ASIS at equal depth: no rotation

What is the spine doing?

Spine is balanced and upright, no rotation, no lateral curvature.
Normal lordosis in cervical and lumbar spine and normal kyphosis in thoracic spine creating the desired “S” shape.

What is the head doing?

Head is functionally upright with only mild forward/lateral flexion or rotation.

Goals

Maintain proper alignment with an appropriate cushion and back support.

ANTERIOR PELVIC TILT



ANTERIOR PELVIC TILT

What is the pelvis doing?

Pelvis sits with PSIS higher than the ASIS.

What is the spine doing?

Excessive lordosis of the lumbar and cervical spine: the patient hyperextends his or her back over the back of the chair, placing him or her at risk to tip the chair backwards.

What is the head doing?

Excessive lordosis of cervical spine causes hyperextension of the neck and upward eye gaze.

Goals

Utilize a cushion and back support that maximizes contact with the seat surface for optimal pelvic and spinal stability and pressure redistribution.

Stability is the goal, so provide a back support that is tall enough for the patient. Measure from seat surface to the top of shoulder.


Look for a moldable back support to conform to the curvature of the spine.

POSTERIOR PELVIC TILT




RIGHT LATERAL VIEW OF PELVIS

Neutral Pelvic Position



Posterior Pelvic Tilt



POSTERIOR PELVIC TILT

What is the pelvis doing?

Pelvis sits with ASIS higher than the PSIS, resulting in the posterior pelvic tilt, aka sacral sitting posture.

What is the spine doing?

Excessive thoracic kyphosis, which produces a “C” shaped spine.

“Flattening out” of the lordosis of the cervical spine.

“Flattening out” of the lordosis of the lumbar spine.

What is the head doing?

Decreased lordosis in cervical spine, causing forward neck flexion & downward eye gaze to the floor or lap.

Goals

Use a cushion with medial and lateral contour to promote LE alignment and pelvic stability.

Ensure appropriate cushion depth to prevent patient from sliding forward, seeking reduced pressure behind the knees.

Add a rigid insert to prevent hammocking of the seat and cushion and keep the pelvis from collapsing into a posterior pelvic tilt.

If FLEXIBLE: Try a cushion with tapered adductors to load the trochanters, stabilizing the pelvis in the patient’s most neutral alignment.

If FLEXIBLE: Use a cushion with an anti-thrust component to reduce forward sliding of the pelvis into posterior pelvic tilt.

If FIXED: Use an immersion style cushion that contours to the shape of the patient to maximize pressure redistribution and minimize peak pressures on the ITs, sacrum, and coccyx.

If FIXED: Consider opening seat-to-back angle in conjunction with a fixed tilt in the wheelchair, to match the patient’s ROM limitations and minimize forward sliding.

Refer to Page 30 for more information on Fixed and Flexible postures

PELVIC OBLIQUITY



PELVIC OBLIQUITY

What is the pelvis doing?

Pelvis sits with the L or R ASIS higher than the other, causing one hip to raise.

What is the spine doing?

When one side of the pelvis is raised higher than the other, the thoracic spine curves away from the higher side, creating a scoliosis over time.

What is the head doing?

The neck will go into lateral flexion as if the person is dropping the ear to his or her shoulder. The lateral flexion will usually be towards the side where the hip is higher.


Goals

Pressure redistribution is the goal whether the deformity is FIXED or FLEXIBLE.

If FLEXIBLE: level the pelvis by building up the lower side.

If FIXED: accommodate for the deformity. Protect the bony prominences from pressure by "filling in" the higher side and immersing the lower side IT.

PELVIC ROTATION



PELVIC ROTATION

What is the pelvis doing?

Pelvis sits with L or R ASIS more forward than the other, producing the rotation in the hips.

What is the spine doing?

The thoracic spine follows and rotates in the same direction as the pelvis. Therefore, if the right side of the pelvis is rotated more forward, the right side of the spine is rotated more forward as well.

What is the head doing?

The neck will go into lateral flexion as if the person is dropping the ear to his or her shoulder.

Goals


Stabilize the pelvis in the patient's most neutral position—adjusting for flexible or fixed postures—to prevent further pelvic rotation.

Consider an anti-thrust cushion to reduce forward migration of the pelvis.

If FLEXIBLE: use tapered adductors and a medial abductor to create leg troughs for midline LE alignment and to maintain a pelvic neutral position.

If FIXED: look for a cushion with less aggressive contouring. Protect the bony prominences with an immersion style cushion.

WINDSWEPT POSTURE



WINDSWEPT POSTURE

What are the pelvis and LEs doing?	Usually goes hand in hand with a rotation or an obliquity and causes one leg to adduct and the other leg to abduct, presenting as if the legs were “swept” away together to one side of the chair.
What is the spine doing?	Will present usually with a rotation or scoliosis, depending on the position of the pelvis.
What is the head doing?	The neck will go into lateral flexion and may be accompanied by forward flexion.
Goals	Look for a cushion that can be adjusted (built-up or minimized) to accommodate for internal rotation and adduction of one LE and external rotation and abduction of the opposite LE.
	Avoid use of ELRs that decrease femoral contact with seat surface and promote windswept posture.
	Windswept posture often accompanies either a rotation or an obliquity, so follow the goals for the underlying causes of this posture.

FIXED VS FLEXIBLE POSTURAL ABNORMALITIES



HOW DO I KNOW WHEN TO CORRECT THE POSTURAL ABNORMALITY AND WHEN TO ACCOMMODATE FOR IT?



Each abnormality is either **FIXED** or **FLEXIBLE**. Understanding whether the postural abnormality is fixed or flexible helps decide if we need to **ACCOMMODATE** for or to **CORRECT** the postural abnormality.

FIXED = ACCOMMODATE

When a posture is **FIXED**, the skeleton does not move out of that posture. The patient needs equipment that will **ACCOMMODATE** the fixed posture, providing optimal support and pressure redistribution.

GOAL = Preventing further progression

FLEXIBLE = CORRECT

When a posture is **FLEXIBLE**, the skeleton still moves and equipment should be chosen with the goal to **CORRECT** the abnormal posture.


GOAL = Prevention from becoming a fixed posture

* FYI- The standardized terms in the wheelchair seating world are Non-reducible for *Fixed* and Reducible for *Flexible*.

- Why? The term *non-reducible* (vs *fixed*) accounts for the potential worsening of a postural abnormality.

ISSUES WITH THE CURRENT WHEELCHAIR SYSTEM CAN CAUSE ABNORMAL POSTURE

HOW DOES AN IMPROPERLY FITTING WHEELCHAIR IMPACT MY PATIENT’S POSTURE?

 The wheelchair can definitely be at fault when the patient is sitting in one of the five abnormal postures. Here is a closer look into seating dimensions and how an ill-fitting wheelchair component can be the culprit behind the body moving into unwanted, harmful postures.

- **Seat Dimension Issues:** Pages 32-34
- **Seat-to-Floor Height Issues:** Pages 35-36
- **Back Support Issues:** Pages 37-38
- **Leg Rest Issues:** Pages 39-40
- **Arm Rest Issues:** Page 41
- **Head Support Issues:** Page 42

What is going on?	What can you do?
Use this column as your problem list in your evals.	This column gives you ways to address the issues present. Letters referenced correspond to measuring guide on pages 46-54. Use those measurements to find the appropriate wheelchair dimension to fit your patient’s body size.

SEAT DIMENSION

Wheelchair seat width and depth are crucial to proper pelvic and LE alignment

WIDTH: TOO WIDE

What is going on?	Negative Result	What can you do?
Patient leans to one side to increase stability. Pelvis will be lower on that side	Pelvic obliquity	Measure hip width (C)
Pelvis collapses	Posterior pelvic tilt, aka sacral sitting	Measure hip width (C)
Pelvic collapse causes hips to internally rotate and LEs to excessively adduct	Strain & contracture risk at hip joints Pressure injury risk at medial knees where knees rub together	Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment
LEs "sweep" to one side when LE weakness is present	Windswept positioning of LEs	Measure hip width (C) Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment

WIDTH: TOO NARROW

What is going on?	Negative Result	What can you do?
Excess pressure at trochanter from contact with the chair	Pressure injury risk at the trochanter	Measure hip width (C) Look for a skin protection cushion that immerses the trochanters
Patient rotates hips to "fit" into the chair	Pelvic rotation	Measure hip width (C)
Patient "sweeps" LEs to one side trying to avoid a buildup of pressure on the trochanters	Windswept posture of LEs	Measure hip width (C) Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment

SEAT DIMENSION

HAMMOCKING SEAT SLING

What is going on?	Negative Result	What can you do?
Pelvis collapses	Posterior pelvic tilt, aka sacral sitting	Add Rigid Insert
Patient seeks out one side of wheelchair for stability	Pelvic obliquity	Add Rigid Insert
Patient compensates by rotating pelvis for stability	Pelvic rotation	Add Rigid Insert
LEs "sweep" to one side	Windswept posture of LEs	Add Rigid Insert

DEPTH: TOO DEEP


What is going on?	Negative Result	What can you do?
Seat sling digs into the back of legs, causing pain Patient will slide forward to alleviate the pain	Posterior pelvic tilt, aka sacral sitting	Measure upper leg length (K) minus 2"
Seat sling digs into the back of legs, decreasing circulation, increasing LE edema Patient slides forward to alleviate numbness	Posterior pelvic tilt, aka sacral sitting	Measure upper leg length (K) minus 2"
Foot propulsion more difficult, patient slides forward for better heel strike	Posterior pelvic tilt, aka sacral sitting	Measure upper leg length (K) minus 2"
Patient slides forward immediately after repositioning	Posterior pelvic tilt, aka sacral sitting	Measure upper leg length (K) minus 2"

*** At a minimum, use a MWC that has seat depth adjustability: K0004**

DEPTH: TOO SHALLOW

What is going on?	Negative Result	What can you do?
Decreased femoral contact and LE support	Windswept posture, abduction, or adduction of LEs	Measure upper leg length (K) minus 2" Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment
Reduced area for pressure redistribution results in increased pressure at the ITs, sacrum, and coccyx	Pressure injury risk at ITs, sacrum, and coccyx	Measure upper leg length (K) minus 2" Look for a skin protection cushion that immerses and/or offloads the bony prominences
Pelvis collapses inward due to lack of LE support to ensure pelvic alignment Hips internally rotate and LEs excessively adduct	Strain & contracture risk at hip joints Pressure injury risk at medial knees where knees rub together	Measure upper leg length (K) minus 2" Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment

*** When LE alignment is an issue: think CONTOURS when choosing a cushion!**



SEAT-TO-FLOOR HEIGHT (STFH)

STFH is crucial for proper pelvic alignment, LE alignment, and heel strike for self-propulsion

Below is a quick reference to match lower leg length to wheelchair model for proper STFH

Wheelchair Model	Achievable STFH	Lower Leg Length Required to Fit MWC
K0001	21"	23" or longer
K0002	19"	21" or longer
K0003	17" - 19"	19" or longer
K0004	As low as 13.5" (most 14.5")	15.5" or longer (16.5")
K0005	Completely Customizable	Any lower leg length

*** At minimum, use a MWC with STFH adjustability: K0004**

STFH: TOO LOW

What is going on?	Negative Result	What can you do?
Leg rests are shortened to compensate for lack of threshold clearance Knees are higher than the hips, causing decreased femoral contact with the seat surface and LEs "sweep" to one side	Windswept positioning of LEs	Measure lower leg length (L) Use a higher profile cushion Use a wedge cushion if they can tolerate the hip flexion Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment
Knees are higher than hips, increasing peak pressure at ITs, sacrum, and coccyx	Pressure injury risk at ITs, sacrum, and coccyx	Measure lower leg length (L) Use a higher profile cushion Look for a skin protection cushion that immerses and/or offloads the bony prominences
On a chair without leg rests, feet can drag and get caught under chair during propulsion or transport	Patient can be thrown out of chair or injury to LEs can occur	Measure lower leg length (L) Use a higher profile cushion

STFH: TOO HIGH

What is going on?	Negative Result	What can you do?
Foot propulsion is more difficult so patient slides forward for better heel strike	Posterior pelvic tilt, aka sacral sitting	Measure lower leg length (L) Use a lower profile cushion
Increased pressure at distal thigh, increasing risk of LE edema, patient slides forward to reduce numbness and pain	Posterior pelvic tilt, aka sacral sitting	Measure lower leg length (L) Use a lower profile cushion
Patient rotates forward on one side for better heel strike of one foot	Pelvic rotation	Measure lower leg length (L) Use a lower profile cushion
Feet “dangle”	Increased internal rotation and adduction of the hip	Measure lower leg length (L) Use a lower profile cushion Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment
Patient slides forward to alleviate strain on pelvis and knees	Posterior pelvic tilt, aka sacral sitting	Measure lower leg length (L) Use a lower profile cushion

*** Choose a higher or lower profile cushion to compensate for incorrect STFH if switching out the MWC just isn't an option!**


TOO HIGH



TOO LOW



GOOD STFH



Refer to Page 53 for more information on measuring for proper STFH

BACK SUPPORT

An appropriate back support can: align the spine, provide pelvic stability, decrease risk of pressure injury development and pain, and improve function.

When choosing appropriate back support height, define your **GOAL** for the patient:

- If **positioning** is the objective, then measure seat to top of shoulder height (G) to choose a back support that can provide optimal trunk stabilization.
- If **self-propulsion** for a patient with good trunk strength is the goal, measure seat to inferior angle of scapula (H) for increased freedom of movement and access to the wheels.

HEIGHT: TOO LOW


This causes insufficient stability for a patient with poor trunk strength and balance

What is going on?	Negative Result	What can you do?
Patient slides down in the chair, seeking more surface area for increased stability	Posterior pelvic tilt with kyphosis of the thoracic spine	Measure seat to top of shoulder (G) or Seat to inferior angle of scapula (H) Provide an appropriately sized contoured back support for added stability that allows for immersion and envelopment of patient's curvature
Patient seeks out one arm rest to gain more stability	Pelvic obliquity with scoliosis of the spine	Measure seat to top of shoulder (G) or Seat to inferior angle of scapula (H) Provide an appropriately sized contoured back support for added stability that allows for immersion and envelopment of patient's curvature
Patient rotates spine and pelvis to seek out more stability on one side of the body	Pelvic rotation with rotation of the spine	Measure seat to top of shoulder (G) or Seat to inferior angle of scapula (H) Provide an appropriately sized contoured back support for added stability that allows for immersion and envelopment of patient's curvature

HEIGHT: TOO HIGH/TOO UPRIGHT

What is going on?	Negative Result	What can you do?
Patient with poor core strength will slide down to alleviate fatigue	Posterior pelvic tilt with kyphosis of the thoracic spine	<p>Measure seat to top of shoulder (G) or Seat to inferior angle of scapula (H)</p> <p>Open seat-to-back angle either through back support hardware or by adjusting the back canes in a MWC with that option</p> <p>Provide an appropriately sized contoured back support for added stability that allows for immersion and envelopment of patient's curvature</p>
Sitting too upright pushes the trunk forward resulting in instability	Anterior pelvic tilt with lumbar lordosis initially. Eventually slides into a posterior pelvic tilt	<p>Measure seat to top of shoulder (G) or Seat to inferior angle of scapula (H)</p> <p>Open seat-to-back angle either through back support hardware or by adjusting the back canes in a MWC with that option</p> <p>Provide an appropriately sized contoured back support for added stability that allows for immersion and envelopment of patient's curvature</p>

*** Another reason to choose a K0004 MWC at a minimum, is the ability to adjust seat-to-back angle. This helps to alleviate fatigue of the trunk muscles and reduce the sensation of being pushed out of the chair.**



Refer to Pages 51-52 for more on measuring for a proper back support

LEG REST

Properly fitting leg rests not only provide a place to rest the feet, they:

- Stabilize the LEs for optimal pelvic and spine alignment
- Promote femoral contact for pressure redistribution away from the bony prominences to reduce risk of pressure injuries

USE OF ONLY ONE LEG REST

What is going on?	Negative Result	What can you do?
Pelvis is higher on the side with the leg rest	Pelvic obliquity	Add second leg rest if positioning is your goal
One-legged heel strike for propulsion	Pelvic rotation	Add second leg rest if positioning is your goal
LEs “sweep” toward the side with the leg rest	Windswept positioning of LEs	Add second leg rest if positioning is your goal

*** Sometimes only one leg rest is necessary, but having only one can cause the above issues. We can’t have our cake and eat it too. When considering the addition of a second leg rest, decide what is more important: one-legged propulsion OR the prevention of a postural abnormality.**

UNEQUAL FOOT PLATE HEIGHTS

What is going on?	Negative Result	What can you do?
Unequal foot plate heights cause one side of the pelvis to be higher	Pelvic obliquity	Adjust foot plate heights to be even

LEG REST TOO SHORT/FOOT PLATE TOO HIGH

What is going on?	Negative Result	What can you do?
Decreased femoral contact, which increases peak pressure at the ITs, sacrum, and coccyx	Pressure injury risk at ITs, sacrum, and coccyx	Measure lower leg length (L) Lengthen leg rest/lower foot plate Look for a skin protection cushion that immerses and/or offloads the bony prominences
Decreased femoral contact causes LEs to “sweep” to one side	Windswept positioning of LEs	Measure lower leg length (L) Lengthen leg rest/lower foot plate Find a cushion with medial abduction and lateral adduction contours to maintain LE alignment
Decreased femoral contact increases hip flexion, pulling downward on the pelvis	Posterior pelvic tilt, aka sacral sitting Shortened hamstrings	Measure lower leg length Lengthen leg rest/lower foot plate

LEG REST TOO LONG/FOOT PLATE TOO LOW

What is going on?	Negative Result	What can you do?
Patient slides forward and stretches legs to reach foot plates	Posterior pelvic tilt, aka sacral sitting	Measure lower leg length (L) Shorten leg rest/raise foot plate
Promotes sacral sitting, adding pressure directly onto the ITs, sacrum, and coccyx	Pressure injury risk at ITs, sacrum, and coccyx	Measure lower leg length (L) Shorten leg rest/raise foot plate Look for a skin protection cushion that immerses and/or offloads the bony prominences
Patient stretches foot & ankle downward to reach a foot plate that is too low, promoting ankle plantar flexion & inversion	Contracture risk of ankle joint	Measure lower leg length (L) Shorten leg rest/raise foot plate Use a single or double foot support
Patient stretches to reach the foot plate, only the ball of the foot; makes contact with the foot plate	Risk of eliciting abnormal reflexes and tone	Measure lower leg length (L) Shorten leg rest/raise foot plate Use a single or double foot support

ARM REST

Properly adjusted arm rests serve many functions:

- A place to rest UEs
- Contribute to overall trunk stability and upright posture when muscle weakness is present
- Facilitates ease of transfers

HEIGHT: TOO LOW

What is going on?	Negative Result	What can you do?
Patient slides down in chair to make contact with arm rests	Posterior pelvic tilt with kyphosis of the thoracic spine	Measure seat to elbow (I) and adjust arm rest height accordingly
Patient leans to one side, seeking more support: pelvis on that side will be lower	Pelvic obliquity with scoliosis of the spine	Measure seat to elbow (I) and adjust arm rest height accordingly
Excessive shoulder depression to make contact with arm rest allows gravity to pull at shoulder joint	UEs are placed in in lap	Measure seat to elbow (I) and adjust arm rest height accordingly

HEIGHT: TOO HIGH

What is going on?	Negative Result	What can you do?
Excessive elevation of shoulder to place arms on arm rests	Shoulder internal rotation and abduction causing pain, numbness, and contracture risk	Measure seat to elbow (I) and adjust arm rest height accordingly
Patient leans against arm rest seeking stability versus placing arm on arm rest	Lateral leaning of trunk	Measure seat to elbow (I) and adjust arm rest height accordingly
Patient places UEs in lap to avoid discomfort of shoulder elevation	Posterior pelvic tilt with kyphosis of the thoracic spine	Measure seat to elbow (I) and adjust arm rest height accordingly


*** Use a MWC with arm rest adjustability.
Choose a K0004 at a minimum.**

HEAD SUPPORT

Proper head support is important for:

- Socialization and Communication
- Safe Swallowing
- Respiration
- Attention to Task
- Mobility

A standard posterior head support pad should be positioned in the suboccipital area with the goal to prevent excessive extension, rotation, or lateral flexion.



TOO HIGH		
What is going on?	Negative Result	What can you do?
The pad will rest on the occipital area or above, causing strain on the neck	Patient will try to adjust by moving the head away from the pad	Place pad in suboccipital area

TOO LOW		
What is going on?	Negative Result	What can you do?
The pad will rest on the cervical spine, resulting in poor posterior support	Head extension	Place pad in suboccipital area

*** Look for a head support that will simultaneously support the head posteriorly and laterally and is height and angle adjustable.**

THE THERAPY EVALUATION: THE PIECES REQUIRED TO JUSTIFY WHEELED MOBILITY AND SEATING EQUIPMENT



**HOW IS A SEATING EVALUATION
DIFFERENT THAN A TYPICAL
OT/PT EVAL?**



**Don't let the term "Seating Evaluation"
scare you!!**

In addition to your typical evaluation, you **MUST**
assess and document the following:

HOME ENVIRONMENT:

Measurements of doorways,
table height, bed height,
stairs, etc

TRANSPORTATION:

Document what kind of car
they drive and/or ride in

MOBILITY SKILLS:

How do they currently get
around to perform their
daily tasks? Be descriptive
(efficiency, ability)

CURRENT SEATING/ MOBILITY:

Describe ALL current
equipment and the
condition of each

MEASUREMENTS:

Complete the 8 key patient
measurements to ensure
proper fitting equipment

EQUIPMENT TRIALS:

Explain why a lesser product
will not work and how the
selected product increases
independence in the home

*** When documenting your evaluation, be as descriptive as possible. Describe their posture while sitting, what posture is required for optimal independence with ADLs, and how the proper equipment can maximize function and minimize pain.**

- Understand that you are doing your typical OT/PT evaluation and **adding a few components** to your documentation
- Assess limitations in ROM, strength, balance, and endurance and how they affect posture in the wheelchair
- Choose an appropriate wheelchair model that allows for adjustability, seat cushion, back support, head support, and accessories to address the patient's physical limitations
- Schedule a visit with your trusted ATP/Dealer to trial equipment with your patient

*** Is there a postural abnormality? Decide if that postural abnormality is flexible or fixed. Should your goal be to correct or accommodate it?**



Goals of a therapist when fitting a patient for a wheelchair:

- Provide pelvic and trunk stability
- Optimize function for ADLs
- Protect the skin to prevent and/or heal existing pressure injuries
- Maximum comfort for patient
- Minimize unwanted movement
- Correct or accommodate for postural abnormalities
- Prevent the progression of postural abnormalities

HOW TO MEASURE FOR PROPERLY FITTING WHEELCHAIR COMPONENTS:



WHAT SHOULD I KNOW TO TAKE THE MEASUREMENTS NEEDED?



Measuring Tips

- All measurements should be taken while sitting on a **firm, flat surface** in the **posture the therapist is trying to accomplish** for the patient.
- Measure using a **hard measuring tape**, not one that can bend and wrap around the patient's body. This will lead to measurement errors, adding circumference or length.
- All measurements should be taken straight across the body.
- Use a measurement form to remind you of the measurements needed and to record them properly.

*** Don't forget to include excess tissue in all measurements and measure the patient in their ideal posture**

*** Don't forget to look for hip contractures and tight hamstrings during the evaluation. They are often the culprits of poor pelvic positioning!**



WHAT ARE THE MEASUREMENTS THAT I NEED TO TAKE?





You will need to perform the following required measurements when choosing the best fitting equipment. Copy the easy to use measurement form on pages 47-48 to document your findings in the field.



Measuring Guide

Letters refer to the measurement form on 47-48


- **SHOULDER WIDTH (A)**
- **CHEST WIDTH (B)**
- **HIP WIDTH (C)**
- **BETWEEN THE KNEES (D)**
- **TOP OF HEAD (E)**
- **OCCIPUT (F)**
- **SEAT TO TOP OF SHOULDER (G)**
- **INFERIOR ANGLE OF SCAPULA (H)**
- **SEAT TO ELBOW (I)**
- **ELBOW TO TIP OF FINGERS (J)**
- **UPPER LEG LENGTH (K)**
- **LOWER LEG LENGTH (L)**
- **FOOT LENGTH (M)**



Don't worry!!! Keep reading, we will walk you through the list. And if you're short on time, we offer a solution.




Measuring Guide



A - Shoulder Width *	
B - Chest Width *	
C - Hip Width *	
D - Between Knees	

*There are 8 MUST HAVE Measurements



E - Top of Head	
F - Occiput	
G - Seat to Top of Shoulder *	
H - Inferior Angle of Scapula *	
I - Seat to Elbow *	
J - Elbow to Tip of Fingers	
K - Upper Leg Length *	
L - Lower Leg Length *	
M - Foot Length	




WHAT IF I DON'T HAVE ENOUGH TIME TO TAKE ALL THE MEASUREMENTS?



Since time is a constraint, the top 8 measurements you **MUST** do are:

- **Shoulder Width (A)**
- **Chest Width (B)**
- **Hip Width (C)**
- **Seat to Top of Shoulder (G)**
- **Inferior Angle of Scapula (H)**
- **Seat to Elbow (I)**
- **Upper Leg Length (K)**
- **Lower Leg Length (L)**


These will allow you to match a wheelchair, including a cushion and back support, to your patient.




REMEMBER TO HAVE PATIENT SITTING ON A FLAT, FIRM SURFACE FOR ACCURATE MEASUREMENTS



SEAT DIMENSION



WIDTH



DEPTH: K - 2"

HOW TO MEASURE **WIDTH**:


- Measure from widest point to widest point of the hips, including all residual tissue

HOW TO MEASURE **DEPTH**:

- Measure from the buttocks, including any excess tissue, across the femur to the popliteal fossa
- Then subtract 2"



BACK SUPPORT WIDTH



HOW TO MEASURE **SHOULDER WIDTH (A)**:


- Measure from humeral head to humeral head, incorporating excess tissue

HOW TO MEASURE **CHEST WIDTH (B)**:

- Measure from axilla to axilla



BACK SUPPORT HEIGHT



HOW TO MEASURE HEIGHT:


- Measure the person in the position you are striving to achieve

* Measure both sides L & R in case of scoliosis and/or a pelvic obliquity

* An average back support height for a self propeller is 16" tall, unless you have an active client who requires more freedom of movement and is more independent.



SEAT-TO-FLOOR HEIGHT (STFH, AKA LOWER LEG LENGTH)



HOW TO MEASURE:


- Measure from the top of the patient's seat (top of cushion if applicable) to the floor
- Measure with feet on floor (with shoes that they normally wear or barefoot if they do not normally wear shoes) and measure bottom of foot/shoe to popliteal fossa

*** There must be a 2" clearance between foot plate and floor to clear thresholds**

*** If providing the patient with a new cushion, remember to consider the thickness of the desired cushion and adjust STFH accordingly.**



ARM REST



HOW TO MEASURE **HEIGHT (I):**

- Seat your patient with shoulder in neutral and elbow bent at his or her side to 90°
- Measure from top of seat to under forearm/ elbow

*** Measure both side in case of pelvic obliquity/ scoliosis. This will justify the need for a model of wheelchair with adjustable height arm rests.**

- If patient needs a specialized arm length (J), measure from end of elbow to where the therapist wants the desired replacement arm support to end




HEAD SUPPORT

HOW TO MEASURE HEIGHT:

There is no miracle formula for fitting a patient for a proper head support

- The rest of the body must be in the best alignment possible before trying to fit for a head support
- Remember that head position changes constantly with the slightest body movement
- Aim to position the pad in the suboccipital area
- Choose a head support with lateral, anterior/posterior, height, and angle adjustability



*** Head support hardware that has offset capabilities accommodates patients with head position out of midline.**

CHOOSING A WHEELCHAIR SEAT CUSHION



HOW DO I KNOW WHAT CUSHION TO CHOOSE FROM THE HUNDREDS OUT THERE?



The first step is knowing **what your client qualifies for**. This is determined by diagnosis and documented postural presentation. There are resources that can help you find the right code but don't be afraid to ask for help from your trusted ATP/dealer!

Once you know the code, you still have so many options. What then?

There are three areas to understand when choosing a wheelchair seat cushion:



- **Methods of Pressure Redistribution**
- **Cushion Geometry**
- **Cushion Medium**

* When choosing a cushion, always keep in mind your goals for the pelvis. For example, if the pelvis is in a *FIXED* posture, your goal is to ACCOMMODATE that posture through cushion geometry and medium. If the pelvis is *FLEXIBLE*, you will look for contours and a medium that can CORRECT the pelvic position.

METHODS OF PRESSURE REDISTRIBUTION



Offloading

The principle of **taking pressure off of a small surface area and loading it onto a greater surface area** that can withstand more pressure and prevent unwanted skin breakdown



Immersion

The principle of **conforming** to the person's curvature by "sinking the body in." We allow the cushion and/or back support to take the body's shape, **alleviating the bony prominences** from **unwanted peak pressure** to maximize **pressure redistribution**



BENEFITS AND CONSIDERATIONS OF USING AN OFFLOADING CUSHION

Benefits	Considerations
Firmer surface, creating a “safer” end feel for the patient	May not be suitable for someone with hypersensitivity due to the firm end feel; depends on personal preference
Partially or completely suspends the ITs, sacrum, and coccyx for pressure relief	Since the nature of the contours adds more stability, it may be more difficult for independent transfers
Redistributes pressure to the trochanters and femurs, areas with greater surface area that can withstand more pressure	Someone with significant contractures may not “fit” the pre-contoured shape
Provides more stability to someone lacking trunk or pelvic strength to prevent unwanted movement	Not good for patients with trochanter wounds since pressure is redistributed to that area

BENEFITS AND CONSIDERATIONS OF USING AN IMMERSION CUSHION

Benefits	Considerations
Less firm surface, creating a “softer” end feel for the patient, may be an ideal solution for pain management	Less stable surface that may not be suitable for someone with trunk and pelvic weakness, causing fatigue with prolonged, upright posture
Envelops or captures the exact curvature of the pelvis to increase pressure redistribution over the largest area possible	May not be suitable for someone who desires proprioceptive input from a firmer surface to maintain an upright, midline posture
Good for any type of pressure injury at any stage, depending on the medium	
May be easier for independent transfers	Some immersion style cushions require more maintenance, and if not properly maintained, can increase risk of a pressure injury
Reduces peak pressures at the bony prominences (ITs, sacrum, and coccyx)	

CUSHION GEOMETRY

Linear

- Refers to a flat seat cushion
- Can be foam layers, foam with gel, or air

LINEAR



*** When choosing a linear cushion, the medium is CRITICAL! Choose a high quality medium that will allow for immersion and envelopment of the bony prominences.**

Contoured

Contoured cushions are made of up various geometric components:


1 Adductors: Tapered or Straight

2 Medial Abductor

3 Posterior Pelvic Well

4 Anti-Thrust


5 Wedge (Not Pictured)




Refer to page 65 for the benefits and considerations of each style of cushion geometry

1 ADDUCTORS

Front



Rear



Front View

Straight

- Can be rear, front, or entire length of cushion
- Act as a boundary to assist with LE alignment
- Minimize abduction of LEs

Benefits of using Straight

Adductors:

- Provides proprioceptive input and a boundary to the lateral surface of the leg to help with LE alignment

Front



Rear



Rear View

Tapered


- Wider in the rear to form a shelf where the trochanters sit
- Allows for loading of trochanters and femurs for pressure redistribution
- Locks the head of the femur into the acetabulum, stabilizing the pelvis
- Promotes offloading of the ITs, sacrum, and coccyx

Trochanter

Benefits of using Tapered Adductors:


- Optimal offloading
- Pelvic stability for patients with trunk and pelvic weakness
- Prevents unwanted shifting in chair PLUS pressure redistribution

② MEDIAL ABDUCTOR




Abductor

- Minimizes adduction of LEs and promotes LE alignment
- Helps create leg troughs when combined with adductors



* Always consider the anatomical shape of your patient and whether he or she will fit in the contours.

3 POSTERIOR PELVIC WELL



What does a Posterior Pelvic Well do?

- Delineates where the pelvis is supposed to be on the cushion for stabilization and pressure redistribution
- Works in combination with lateral tapered adductors to create a deeper pelvic well, which will offload and suspend the ITs, sacrum, and coccyx


4 ANTI-THRUST

LINEAR VS ANTI-THRUST

LINEAR



ANTI-THRUST



What does an Anti-Thrust do?

- Provides a boundary to minimize forward migration of the ITs into a posterior pelvic tilt
- Often used when extensor tone is present
- Standard placement is half the depth of the cushion to allow for 1" – 1½" from ITs to the front edge of anti-thrust

5 WEDGE

LINEAR VS WEDGE

LINEAR



WEDGE




What does a Wedge do?

- Used to close the seat-to-back angle and reduce extensor tone
- Used to maintain pelvic positioning for certain pelvic abnormalities
- Can be used to increase STFH for someone with long LEs

Now take a look at what cushion contouring does for pressure redistribution

Contours increase **surface area contact** with the patient minimizing peak pressures.

WITHOUT CONTOUR




Front View




Rear View

WITH CONTOUR



Front View



Rear View

BENEFITS AND CONSIDERATIONS OF USING A LINEAR CUSHION

Benefits	Considerations
Allows freedom of movement	If using traditional gel, it disperses over time, exposing the bony prominences to peak pressures
Fits any body type	Over time, low quality foam compresses and loses its shape resulting in increased peak pressures
Can accommodate for contractures if using an immersion style	If using air, look for a cushion that provides stability through its construction and does not require frequent maintenance
Immersion style provides envelopment of bony prominences	Decreased stability resulting in fatigue over time
	Often needs to be used with a rigid insert in a sling seat

BENEFITS AND CONSIDERATIONS OF USING A CONTOURED CUSHION

Benefits	Considerations
Promotes pressure redistribution to greater surface areas, reducing peak pressures	Various body shapes & sizes may not always fit in the contours
	May not fit someone with significant contractures
Increased stability, securing pelvis in optimal position	May restrict freedom of movement during ADLs
Can accommodate or correct abnormal postures	Aggressive contours may make independent transfers more difficult

CUSHION MEDIUM

Understand the available mediums when choosing a cushion. No matter the contours, the quality of the medium will affect the application of the cushion. Each has benefits and considerations, and we will look at each one.

Cushion Mediums:

- **Foam**
- **Air**
- **Honeycomb**
- **Gel/Fluid Inserts**
- **Combination**

Considerations



Weight of the cushion- Each medium has different properties, making them lighter or heavier.



Support/Stability Needs- Each medium has different qualities that make them more or less supportive, which impacts how long a patient can sit in their chair, ease of transfers, and pelvic posture when weakness is present. Remember to think of the functional and postural goals of your patient!

*** Foam is not foam is not foam! The quality and properties of foam affect the application and effectiveness of a cushion. Be an advocate for your patients and CHOOSE QUALITY PRODUCTS!**

FOAMS TYPICALLY USED IN CUSHIONS

HIGH RESILIENCY (HR) FOAM

- Instantly "springs" back to its original shape post compression
- Durable foam, can withstand repeated compression and maintains resiliency over time
- Provides structure and stability
- Great for use as the **base** of a cushion



VISCO ELASTIC FOAM

- Memory foam, meaning it takes time to resume the original shape post compression
- Allows for the most immersion and envelopment, contouring to the unique shape of the user
- Great to use as top layer of a cushion



*** Density is the measure of quality for foam. A density of 1-3 pounds is considered good for heavy use. What you need to know is, you get what you pay for. WHEN THE COST OF A CUSHION SEEMS TOO GOOD TO BE TRUE, IT IS!!!**

FOAM

Foam Benefits	Foam Considerations
Can be designed to be as supportive/contoured as needed	Foam can be heavy. Consider the weight of the cushion
Can allow for offloading or immersion	Need to protect the foam from incontinence
Low maintenance	Most foams inherently retain heat and moisture

AIR

There are different types of air cushions:
Individual air cells and air bladders

Air Benefits	Air Considerations
Offer envelopment and immersion of bony prominences, increasing pressure redistribution	Certain styles require frequent maintenance
Perceived as “soft” & comfortable	Certain styles are not stable, resulting in fatigue, decreased pelvic stability, and difficulty with transfers
Can be adjustable	Easily affected by altitude: higher altitude = firmer

HONEYCOMB

These cushions are made of a thermoplastic material, using a hexagonal honeycomb shape

Honeycomb Benefits	Honeycomb Considerations
Inherently breathable, addresses microclimate	Not adjustable and may not accommodate orthopedic deformities
Lightweight	
Allows freedom of movement	Balance and endurance of patient
Low maintenance	Can be perceived as “firm” and this is not adjustable

GEL/FLUID INSERTS

There is a perception that gel is the best medium for skin protection: NOT TRUE!

Gel/Fluid Benefits	Gel/Fluid Considerations
Can allow for envelopment and immersion of bony prominences for increased pressure redistribution	Their effectiveness is dependent on the geometry and stability of the cushion
	Not a stable surface, so provides minimal pelvic stability; transfers may be difficult
	Gel can migrate and result in high peak pressures
Some gels state that they are “cooling” to manage microclimate	Certain styles require daily maintenance or with every reposition, need to be kneaded and readjusted to provide pressure relief
New gel technology offers solid gel options, eliminating the problem of gel dispersion	With sun exposure, can retain heat
	More susceptible to damage than other mediums, can leak

COMBINATION

These are a combination of foam and air or foam and gel; style is dependent on the manufacturer

Combination Benefits	Combination Considerations
Foam provides a stable base for transfers with the skin protection of air/gel	Can be heavy depending on medium/manufacturer
Contours are supportive of pelvis in optimal position and provide LE alignment	With well cut outs, a migration of gel or air loss can result in ITs resting on edge of hard foam, causing peak pressure
	May need frequent maintenance to ensure gel packs are in position or air cells are properly inflated

WHEN TO CHOOSE A SEPARATE OFF-THE-SHELF BACK SUPPORT



HOW DO I KNOW IF MY PATIENT NEEDS IT?



There are key questions to ask yourself when deciding if the manufacturers back support should be replaced with an off-the-shelf back support:

- Is your patient immobile?
- Does your patient have a postural abnormality that affects pelvic and trunk stability and alignment?
- Does your patient require assistance to sit unsupported?
- Does your patient complain of pain while sitting in the wheelchair?

If you answered yes to any of these questions, then, without a doubt, replace the sling back with a separate off-the-shelf back support.

*** If your patient qualifies for any model wheelchair, they automatically qualify for at least a general use back support. Always consider replacing the sling back with an off-the-shelf option to provide better pelvic and trunk positioning.**



The right back support is an essential part of an **optimal seating system**. Here are some of the important **ADVANTAGES** of an off-the-shelf back support:

- Provides depth and angle adjustability through **hardware**
- Provides **trunk stability** to promote functional sitting
- Works with the cushion to provide **pelvic stability**
- Provides increased **pressure redistribution** at the trunk and pelvis
- It can **minimize the progression of abnormal postures**
- **Decreases pain** by increasing stability and pressure redistribution

* Even when dealing with a more "typical" spinal posture, a solid back support can provide stability to conserve energy for self-propulsion.



Optimal Seating System

CHOOSING A BACK SUPPORT



WHAT DO WE NEED TO CONSIDER WHEN CHOOSING A BACK SUPPORT?



There are three areas to understand when choosing a wheelchair back support:

- **Shape of the Shell**
- **Medium**
- **Adjustability**




Back Support Goals:

- The patient will function at his/her optimal capacity
- Stabilize the pelvis to maximize function of the cushion
- Provide posterior thoracic stability
- Facilitate maximum pressure redistribution to minimize high peak pressure areas
- Maintain and support natural curves of the spine

SHAPE OF THE SHELL

Back support shells can vary in their geometry:

- **Lateral Contoured Shell**
- **Posterior Open Shell**
- **Posterior Planar Shell**



LATERAL CONTOURED SHELL

Benefits	Considerations
Designed to fit the anatomical shape of the user	Contour is lateral, not on the posterior shell
Various options depending on patient goals. Deep, upper thoracic or lower thoracic contours	
Look at the contour of the upper shell for scapular cut outs if your patient is a self-propeller	Look for a contoured shell with the ability to conform posteriorly

POSTERIOR OPEN SHELL

Benefits	Considerations
New concept	
Allows for immersion posteriorly with technology that conforms to the shape of the end user	May not be able to adjust for a significant lordosis
Lightweight	
Can accommodate and/or correct for multiple postural abnormalities	

POSTERIOR PLANAR SHELL

Benefits	Considerations
Certain planar styles can be adjustable for conforming to shape	Most have no ability to conform to the shape of the user, resulting in the progression of postural abnormalities
Flat shell	
Room to attach harnesses, pelvic belts, lateral trunk supports	

MEDIUM

Understand the available mediums when choosing a back support. No matter the contours, the quality of the mediums will affect the application of the back support. Each has benefits and considerations, and we will look at each one.

BACK SUPPORT MEDIUMS:

- **Foam**
- **Air**
- **Combination**

FOAM

Benefits	Considerations
Most common	Can be hot
A high quality visco elastic foam allows immersion and envelopment of a person's shape	Consider quality and type of foam used as it affects immersion and envelopment and pressure redistribution

AIR

Benefits	Considerations
Air cells that allow the transfer of air, mimicking the pressure-relieving properties of water	May require maintenance
Individual air cells that shift according to the shape of the user	
Allows for optimal envelopment and pressure redistribution	
Adjustable to configure the back support for optimal pressure redistribution.	

COMBINATION

Benefits	Considerations
Stability of foam with the pressure relieving properties of air	May require maintenance

ADJUSTABILITY



WHY IS ADJUSTABILITY IMPORTANT IN BACK SUPPORTS?



Adjustability is important to wheelchair back supports for the following reasons:

- The human skeleton has natural curves that allow for optimal function.
- If these curves are not supported well, they will either flatten or change into an abnormal posture.
- This is especially important when providing a back support for a child to promote natural curves of the spine as he or she grows.




What are ways that we can adjust a back support?

Angle and Depth

- Using the mounting hardware

Posterior Contours


- With moldable stays
- Using The Boa® System
- Through mounted wings or foam wedges



* 90/90/90 will NOT allow us to match the patient's joint angles.



WHY IS ADJUSTING SEAT-TO-BACK ANGLE IMPORTANT?



- Back Support Surface
- Surface Contact
- Lack of Surface Contact

Adjusting seat-to-back angle allows for:

- Maximum contact between the back and seat surface
- The pelvis to reside in the intended area of the cushion
- Pressure redistribution off the bony prominences
- Stabilization at the pelvis and spine
- Accommodation of limited hip flexion

*** Hardware is key when choosing a back support. Always consider what kind of adjustability it has so you can achieve optimal pressure redistribution.**




WHAT SHOULD WE AIM FOR WITH ADJUSTABILITY?




Look at what the patient's body is doing and consider that seat angles depend on the individual's ROM, muscle integrity, and joint integrity.


- Open or close the seat-to-back angle to accommodate for your patient's degree of hip ROM.
- Change seat-to-back angle to allow for more contact with the back support and seat surface.
- Use methods such as moldable stays, wings, or wedges to conform to the patient's individual curvature.



Fixed 90°
Peak Pressure
at the Apex



Open STBA
More Surface Contact
Match Patient ROM



Open STBA + Contour
Optimal Pressure
Redistribution

*** When more surface area makes contact with the patient, stability and comfort increase.**

ACCESSORIES



WHEN WOULD I NEED TO ADD ACCESSORIES?



As in most things, ***LESS IS MORE*** and this is certainly true in wheelchair positioning. If you take your time in selecting the right model and size wheelchair and then select an appropriate cushion and back support, accessories may not be needed. Too often accessories are used as “band-aids” for a less than optimal wheelchair and seating system. When used the proper way, they can be the finishing touch to achieve the best posture possible for your patient!

UPPER EXTREMITY SUPPORTS

Can be added onto the existing arm rest to:

- Add length
- Add additional postural support
- Decrease contracture risk
- Decrease dependent edema in UE
- Increase comfort while sitting in wheelchair
- Stabilize a weak UE in an optimal position

*** Most of us don't sit with our elbows bent to 90° with the forearm straight out. So why do we position our patients this way, especially when dealing with abnormal tone or decreased range of motion? Choose an upper extremity support that has angle adjustability!!**

TYPES OF UPPER EXTREMITY SUPPORTS:

You can choose a basic non-adjustable support or you can choose a pad and hardware style to suit the needs of your patient.



HARDWARE TYPES

- **SLIDE ON:** Pad will be aligned with existing arm rest
- **BOLT ON:** Pad will be aligned with existing arm rest
- **SWING-AWAY:** Pad can move towards or away from body as needed
- **ARTICULATING:** Pad can be positioned in any plane; most adjustable option

PAD TYPES

- Standard trough that keeps arm from sliding off the pad in any direction
- Minimal contour for low profile support
- Moldable option at the wrist for flexion or extension of the wrist and finger abduction

LOWER EXTREMITY SUPPORTS

Can be added to a chair to:

- Support the lower leg, foot, and ankle in optimal position
- Decrease risk of contracture
- Accommodate for contractures and protect from injury
- Manage abnormal tone and reflexes
- Support an amputated limb
- Support a post surgical limb on leg rest

* An external fixator support provides great stability and support for a LE post hip/knee surgery, or when a LE immobilizer is in use.

TYPES OF LE SUPPORTS:

- Single Foot Support
- Double Foot Support
- Double Foot Support with Separator
- Calf Pads
- External Fixator Support
- Amputee Support

*** Don't forget that quality leg supports have a range of sizes or can be customized. Consider where you want the support and measure to ensure you get the right fit.**

LATERAL TRUNK SUPPORTS

- Prevent lateral leaning caused by decreased trunk strength and balance
- Correct or accommodate for a scoliosis or trunk rotation and prevent further progression of the postural abnormality
- Decrease pressure points on the apex of the curvature that could possibly lead to skin breakdown and cause pain to the patient

TYPES OF LATERAL SUPPORTS:

- Mounted swing-away or fixed lateral supports
- Traditional foam lateral supports

*** If your patient has tone and/or is a heavy leaner, choose mounted lateral trunk supports. Traditional foam supports cannot withstand the sustained pressure, and the foam will collapse.**

THE ELEVATING LEG REST MYTH!

ELEVATING LEG RESTS (ELRs) ARE A GOOD OPTION, RIGHT?

NO!

Here are some common myths about ELRs that you **NEED ABSOLUTE CLARIFICATION** on to help you understand why standard leg rests are usually the better option:

See Page 83 for visual

1 ELRs keep hips back in chair: Wrong!

- ELRs do just the opposite
- When we elevate the legs, we pull on the already tightened hamstrings of the patient
- The shortened hamstrings cause a greater pull on the pelvis, bringing the pelvis into a posterior pelvic tilt
- The posterior pelvic tilt makes the pelvis slide forward, resulting in sacral sitting
- The patient begins to slide forward out of the chair, the exact opposite of keeping the hips back in the chair

2 ELRs decrease edema: Wrong!

- Contrary to popular belief ELRs on K0001-K0004 MWCs unfortunately **CANNOT** reduce edema

- In order to decrease edema, the legs must be 30 cm above heart level. The only way to achieve that degree of elevation with an ELR is when used in combination with tilt and recline, often seen in manual tilt-in-space or power chairs
- ELRs on a manual chair can actually decrease optimal circulation by cutting off blood flow at the groin area and inhibiting flow to the lower extremities


3 ELRs decrease pressure on the pelvis: Wrong!

- ELRs actually position the pelvis in a posterior pelvic tilt
- The forced posterior pelvic tilt increases pressure on the bony prominences of the ITs, sacrum, and coccyx
- Pressure injury development risk is increased in those areas
- Elevating the leg rest promotes knee flexion, leading to decreased femoral contact, shifting the pressure back onto the ITs and sacrum/coccyx

4 ELRs help with LE alignment: Wrong!

- ELRs prevent full femoral contact with the seat surface
- Since less of the leg is making contact with the seat surface, it is easier for the leg to internally/externally rotate or abduct/adduct
- Legs then fall off leg rests more easily
- Windswept posture is more prevalent
- ELRs promote flexion of knees, hips, and ankles, increasing risk of contractures at those joints

RIGHT LATERAL VIEW - WHAT IS REALLY HAPPENING WHEN USING AN ELR



*** AVOID the use of ELRs as much as possible for optimal positioning in a wheelchair.**

THE WEDGE MISCONCEPTION



Often a wedge is used to prevent someone sliding out of a chair

However, examine what a wedge actually does:

- Closes seat-to-back angle
- Pulls shortened hamstring muscles, resulting in a posterior pelvic tilt and sliding forward
- Increases the risk of pressure injury development due to peak pressures on the ITs, sacrum, and coccyx that now rest directly on the incline of the wedge



Using a Wedge

So who should we AVOID using a wedge with?

- Someone with tightened hamstrings
- Someone who cannot tolerate a 90° or less (85°, etc) seat-to-back angle

Better for use with:

- Someone with the goal of reducing extensor tone
- Someone with full ROM at the hips and knees
- Someone with long LEs to increase STFH

THE PROBLEM OF SLIDING OUT OF THE WHEELCHAIR

? WHAT SHOULD I DO WHEN MY PATIENT KEEPS SLIDING OUT OF THE CHAIR?



First, we need to find out WHY...

- Check for **tightened hamstrings** as they can pull the pelvis forward in the seated posture
- Check the **hip ROM** and see if they are trying to increase the seat-to-back angle by sliding forward
- Check the **seat depth** and see if it is too deep and digging in to the back of their legs
- Check the **seat-to-floor height**, and see if they are sliding to try and reach the floor for propulsion
- Check **trunk stability and strength**. If weakness is present, the patient may slide forward for stability



Normal Hamstring



Shortened Hamstring



When we know why, we can look at possible SOLUTIONS:

- Accommodate tightened hamstrings by adapting the seat depth and back angle to match the patient's posture. Use an immersion style cushion to protect the sacrum and coccyx from pressure injury
- Adjust the seat-to-back angle either through the back canes or the hardware of the back support. Find their optimal seat-to-back angle based on your mat evaluation results
- Measure upper leg length (Page 50) and get a chair with the appropriate seat depth
- Lower STFH at the wheel axle, or use a drop seat
- Try a cushion with tapered adductors that will stabilize the pelvis and provide increased support for the trunk

*** A wedge cushion is not the solution to sliding forward. Be a detective and find out WHY the patient is sliding forward, and this will lead you to the appropriate solution!**

CONCLUSION

*** Remember, there is NO cookbook solution to wheelchair seating and positioning. Use your clinical reasoning skills to figure out what's best within the given parameters.**



Keep the following goals in mind, and do the best you can!

- Maximize patient comfort
- Minimize pain
- Correct deformity when possible
- Prevent further deformity
- Prevent pressure injuries
- Increase functional safety and independence

REFERENCES

1. "Adult Wheelchairs Categorized and Explained - Northwest Pompano." Active Mobility Center, 10 May 2015. Web. 19 July 2016. <http://nearsay.com/c/44448/33117/adult-wheelchairs-categorized-and-explained>
2. Dicianno, Brad E., Juliana Arva, Jenny M. Lieberman, Mark R. Schmeler, Ana Souza, Kevin Phillips, Michelle Lange, Rosemarie Cooper, Kim Davis, and Kendra L. Betz. "RESNA Position on the Application of Tilt, Recline, and Elevating Legrests for Wheelchairs." *Assistive Technology* 21.1 (2009): 13-22. RESNA.org. RESNA, 25 Feb. 2015. Web. 21 July 2016. http://www.resna.org/sites/default/files/legacy/resources/position-papers/RESNA%20PP%20on%20Tilt%20Recline_2015.pdf
3. Cole, Elizabeth. "Causes and Prevention Skin Breakdown." US Rehab Seating and Mobility Master Program. 2016 Lecture.
4. Cole, Elizabeth. "The Biomechanics of Seating." US Rehab Seating and Mobility Master Program. 2016 Lecture.
5. Cole, Elizabeth. "Matching the Evaluation Results to Seating Interventions Part 1, 2, 3." US Rehab Seating and Mobility Master Program. 2016 Lecture.
6. Cole, Elizabeth. "How to Measure for Seating and Positioning Products." US Rehab Seating and Mobility Master Program. 2016 Lecture.

GLOSSARY



ATP: Assistive Technology Professional

NP: Nurse Practitioner

PA: Physician Assistant

LMN: Letter of Medical Necessity

DME: Durable Medical Equipment

CRT: Complex Rehab Technology

ROM: Range of Motion

ADLs: Activities of Daily Living

MRADLs: Mobility Related ADLs

MWC: Manual Wheelchair

PWC: Power Wheelchair

PPT: Posterior Pelvic Tilt

PSIS: Posterior Superior Iliac Spine

ASIS: Anterior Superior Iliac Spine

IT: Ischial Tuberosity

ELR: Elevating Leg Rest


STFH: Seat-to-Floor Height

LE: Lower Extremity

UE: Upper Extremity

REFERENCES CONTINUED

7. Mounic, Gerard, OTD,OTR/L,C/NDT,COS-C. "Evaluation and Seating Strategies with Regulatory Guidelines, Wheelchair Seating and Positioning Improving Functional Mobility and Patient ADLs." Summit Professional Education Continuing Education Course. North Carolina, Asheville. 2014 Lecture.
8. VGM Course Anatomy and Physiology. US Rehab Seating and Mobility Master Program. 2016 Lecture. https://univ.vgmu.com/opcv/pages/course/run_course.jsp?id=282740
9. Waugh, K; Crane, B; Taylor, S; Davis, K; Cwertnia, S; Brown, L; Saftler, F; Christie, S. Glossary of Wheelchair Terms and Definitions. Published by the University of Colorado, Assistive Technology Partners, through a grant from the Paralyzed Veterans Association. December, 2013.
10. Waugh, K; Crane, B; Taylor, S; Davis, K; Cwertnia, S. A Clinical Application Guide to Standardized Wheelchair Seating Measures of the Body and Seating Support Surface. Published by the University of Colorado, Assistive Technology Partners, through a grant from the Paralyzed Veterans Association. 2012, revised July, 2013.
11. Fedor, Dan. "No More Addendums...Getting it right the first time!". VGM/US Rehab. 2016 Lecture.
12. Pedersen, Jessica Presperin, and Michelle L Lange. "Positioning Chart." POSITIONING CHART, 29 June 2017.



www.comfortcompany.com

education@comfortcompany.com

1.800.564.9548