

Sampling, measurement and analysis of VOCs & SVOCs: What are the best tools for the job?

Nicola Watson

Sale Support Manager - Americas

Markes International

nwatson@markes.com

Agenda

- Sampling VOCs in air – The tools of the trade
Focusing on:
 - Canisters and online sampling
 - Sorbent tubes
 - Active sampling
 - Passive sampling
 - When should each technique be used?

How on-line/canister analysis works

1. Air/gas focusing




Quantitative retention of ultra-volatiles from up to 1.5 L volumes and efficient low-flow, splitless desorption ensure low detection limits.

Peltier-cooled focusing trap eliminates ice-plug formation, while fast trap cooling minimises cycle times.

Samples are introduced directly onto the electrically-cooled, sorbent-packed focusing trap of the UNITY-xr thermal desorber, typically held between ambient and -30°C .

How on-line/canister analysis works


2. Trap desorption and outlet split



Focusing trap rapidly heated (up to 100° C/s) in a reverse flow of carrier gas ('backflush' operation), to transfer the analytes to the GC column.

- Fully automated sequences of air/gas (from a sample stream, calibration gas or zero air/gas) can be set at user-defined frequencies.
- During trap desorption, the flow of analytes can be split and re-collected onto a clean sorbent tube.
- Tubes and traps can contain multiple sorbents, for analysis of an extended range of analytes.


Analysis of Air Toxic compounds – Canisters




Application: Monitoring of C₂–C₁₀ hydrocarbons (ozone precursors) in ambient (TD–GC dual FID)

DB-1 column showing

retention of C₇ to C₉



PLOT column showing
retention of C₂ to C₆
compounds



- **Target compounds:**
 - C₂ to C₁₀ hydrocarbons (ozone precursors)
 - ‘Ozone precursors’ focusing trap at –30°C and flow path at 120°C
- **Performance in field operation:**
 - Detection limits: <0.05 ppb
 - Retention time stability: <0.2% RSD across all compounds
 - Standard reproducibility: 0.2–5% RSD
 - Excellent peak shape for splitless injection

Can canisters do everything?



- ✓ Great for C₂ to C₁₂ compounds
- ✓ Suitable for rapid transfer (not storage) of ultra-volatile reactive compounds such as H₂S
- ✓ Ideal for simple grab-sampling

- ✗ NOT suitable for compounds with volatility less than C_{10/12}
- ✗ NOT suitable for high-concentration samples
- ✗ Time-weighted average sampling is NOT as easy with a canister

How thermal desorption works

1. Tube desorption and inlet split



Sample tubes and traps can contain multiple sorbents, for analysis of a wide range of analytes.

During desorption of the tube the flow of analytes can be split and re-collected onto a clean sorbent tube.

Sample tube heated in flow of carrier gas and analytes swept onto an electrically cooled focusing trap, typically held between ambient and -30°C .

How thermal desorption works

2. Trap desorption and outlet split



Focusing trap rapidly heated (up to 100° C/s) in a reverse flow of carrier gas ('backflush' operation), to transfer the analytes to the GC column.

Extended re-collection

Quantitative re-collection allows:

- Valuable samples to be re-analysed, overcoming the historic 'one-shot' limitation of thermal desorption.
- Allows complete analyte transfer to be validated, ensuring compliance with standard methods.
- Aids method development and troubleshooting.

Sorbent tubes: Sampling techniques


- Active sampling
- Passive sampling

Which method to use?


Pumped sampling

Sample collection

- Sample (e.g. air) is collected



- Compounds of interest are adsorbed on the sorbent surface



- Lighter gases such as nitrogen, argon and carbon dioxide pass through



Flow rate = 20–100 mL/min

Volume = 500 mL to 100 L

Important – do not exceed breakthrough volume for a compound on a given sorbent

Air monitoring - pumped

Sorbent selection for both tubes and focusing trap are very important


Semi-volatile compounds – **Weak sorbent**

Helps prevent retention of unwanted compounds



Very volatile compounds – **Strong sorbent**

Prevents **breakthrough** of light compounds




A complex example (US EPA TO-17)





1	Propene	30	Carbon tetrachloride
2	Dichlorodifluoromethane	31	Benzene
3	1,2-Dichloro-1,1,2,2-tetrafluoroethane (Freon® 114)	32	n-Heptane
4	Chloromethane	33	Trichloroethene
5	Chloroethane	34	1,2-Dichloropropane
6	Buta-1,3-diene	35	1,4-Dioxane
7	Vinyl chloride	36	Bromodichloromethane
8	Bromomethane	37	cis-1,3-Dichloropropene
9	1,2-Dichloroethane	38	4-Methylpentan-2-one
10	Trichlorofluoromethane	39	Toluene
11	Ethanol	40	trans-1,3-Dichloropropene
12	1,1-Dichloroethene	41	1,1,2-Trichloroethane
13	1,1,2-Trichloro-1,2,2-trifluoroethane (Freon® 113)	42	Tetrachloroethene
14	Acetone	43	Hexan-2-one
15	Carbon disulfide	44	Dibromochloromethane
16	Isopropanol	45	1,2-Dibromoethane
17	Dichloromethane	46	Chlorobenzene
18	Methyl tert-butyl ether	47	Ethylbenzene
19	cis-1,2-Dichloroethene	48	m-, p-Xylene
20	n-Hexane	49	
21	1,1-Dichloroethane	50	o-Xylene
22	Vinyl acetate	51	Styrene
23	trans-1,2-Dichloroethene	52	Bromoform
24	Butan-2-one	53	1,1,2,2-Tetrachloroethane
25	Ethyl acetate	54	Trimethylbenzene
26	Tetrahydrofuran	55	Trimethylbenzene
27	Chloroform	56	1-Ethyl-4-methylbenzene
28	1,1,1-Trichloroethane	57	Dichlorobenzene
29	Cyclohexane	58	Dichlorobenzene
		59	α-Chlorotoluene
		60	Dichlorobenzene
		61	1,2,4-Trichlorobenzene
		62	Hexachlorobuta-1,3-diene

Splitless desorption of 'Air toxics' tube loaded with 1 L of 1 ppb std
GC/MS

Source: Markes Application Note TDTs 86





Not just volatiles...



Passive (diffusive) sampling

- Diffusive sampling = a simple and cost-effective method of collecting the large number of samples required in many air monitoring programmes
- Vapours migrate across the air gap at a constant 'uptake rate'
- Diffusive sampling is a slow process – typically sample for days



When should I use diffusive sampling?

- ✓ You know the compound that you are looking for
- ✓ There is a validated uptake rate available for that compound
- ✓ The test atmosphere is not heavily contaminated with a wide range of other organic compounds at much higher concentrations
- ✓ The expected concentration of analyte in the atmosphere is such that the desired sampling time (usually between 4–8 hours (occupational) and 1–4 weeks (environmental)) will result in a mass on the tube which is above the limit of detection of the TD–GC(MS) method
- ✓ You are looking for several compounds of the same volatility

- ✗ You are using a multi-bed sorbent tube
- ✗ You are sampling a completely unknown atmosphere
- ✗ You want to sample two (or more) compounds of widely differing volatilities (e.g. acetone and toluene)
- ✗ There are no uptake rates available for the compounds of interest

EPA 325 – Refinery perimeter monitoring

New federal regulation (CFR 40) to be implemented mid 2015, compliance within 3 years

- Requires continuous monitoring of vapor-phase organics around the boundary of oil refineries. Target analytes:
 - Benzene
 - Hazardous air pollutants (HAPs)
 - VOCs in refinery air (light/middle distillate fuels)
- **US EPA Methods:**
 - 325 A (Sampling): 2-week passive sampling using sorbent tubes.
 - 325 B (Analysis): TD–GC(MS) determination (MS recommended)



~~Which sampling method is best?~~

How do I get the right results?

- What are the compounds of interest?
 - What is the expected concentration?
 - How long is the monitoring period?
 - Is the matrix compatible?
 - Etc.....
-
- You can not determine which method is best until you know the sampling situation.
 - **Define the problem, before finding the solution,** or as Einstein is quoted as having said

“ If I had one hour to save the world, I would spend fifty-five minutes defining the problem and only five minutes finding the solution.”



Questions?

nwatson@markes.com

Acknowledge the work
of the Markes
Application Team

Contact us/Learn More



1-866-483-5684



enquiries@markes.com



www.markes.com



@MarkesInt



Markes International



Google plus



TheMarkes1

