



... efficient, environmentally friendly and cost effective approach

***Single Analysis of Volatile and Semi-Volatile Organics in Air.***  
***The Cost Effective***  
***Green Solution***

*Roberta Provost, Air Method Development Chemist  
Pace Analytical Services  
Lee Marotta, Sr Field Application Scientist  
PerkinElmer Instruments*

- ▶ Introduction
- ▶ Industry Problem/Solution
  - Combining a two analysis solution into one analysis
- ▶ EPA Method TO-17 and Thermal Desorption
- ▶ Analytical Data and Performance
- ▶ Results from Site Sampling (MGP site)
- ▶ Air Monitoring Study Summary

- ▶ 2010: Soil Vapor Intrusion (SVI) Tube (patented)
  - ( $C_3$  to  $C_{26}$ )
  - Combines VOC & SVOC from the seven VOA gases to pyrene
  - Thank you to CARO Analytical Services for their help
- ▶ 2011: XRO-640 (patent pending)
  - ( $C_6$  to  $C_{44}$ )
  - Residue in Liquefied Petroleum Gas (LPG)
  - Combines VOC & SVOC from BTEX to benzo(g,h,i)perylene
  - Thank you to Alberta Innovates Tech Futures for their help
- ▶ 2013: XRO-440 (patent pending)
  - ( $C_4$  to  $C_{44}$ )
  - Combines VOC & SVOC from 1,3-butadiene to benzo(g,h,i)perylene
  - Thank you to Pace Analytical Services for their help





... efficient, environmentally friendly and cost effective approach

***Success***

***Move from a Two Analysis***

***to One Analysis***

***The cost effective,***

***environmentally friendly***

***Approach***

PM10 (particulates)



TO-13 (PAHs)



TO-15 (VOCs)



# Why the XRO-440 Tube was Designed?

- Eliminate liquid extractions
  - Save time and \$\$\$
  - Improve productivity and efficiency
  - Enhance recoveries
  - A Greener analysis

**Two Analysis now ONE**


**Improve Productivity  
Cost Effective  
Environmentally Friendly**

TO - 13



TO - 15

TO - 17



# Only Four Canisters to a Box .... YIKES



... can fit over 500 Tubes!



# The Analytical Solution for Air Monitoring



... efficient, environmentally friendly and cost effective approach


## *The TO-17 Solution for the Analysis of VOCs and SVOCs in air*





# Research required *NEW* TD Tube Design

Volatiles	Semi-volatiles
1,3-Butadiene	Napthalene
Benzene	2-Methylnapthalene
Toluene	1-Methylnapthalene
Ethyl Benzene	Acenaphthylene
Xylenes	Acenaphthene
	Fluorene
	Phenathrene
	Acenaphthene
	Fluoranthene
	Pyrene
	Benzo(a)anthracene
	Chrysene
	Benzo(b)fluoranthene
	Benzo(k)fluoranthene
	Benzo(e)pyrene
	Benzo(a)pyrene
	Indeno(1,2,3-c,d)pyrene
	Dibenz(a,h)anthracene
	Benzo(g,h,i)perylene



... started with these targets for MGP sites

# The Analytical Solution for Air Monitoring



... efficient, environmentally friendly and cost effective approach

*Operation*



# Step 1: Sample Tube (Primary) Desorption

IS and/or surrogate spike (optional)  
Impedance check (optional)  
Ambient purge (at least 1 min)

Desorption via  
✓ Temperature  
✓ Flow  
✓ Time

Inert gas


Inlet split (Optional)

Desorb flow

Peltier cooled to -35°C  
No need for liquid cryogen



## Step 2: Transfer Sample to Analytical Column







... efficient, environmentally friendly and cost effective approach

## ***Analytical Performance***

Data was acquired in  
Simultaneous Full Scan/ SIM  
Mode. Only results from Full  
Scan are reported




# Total Ion Chromatogram (TIC)




Target Compound	Range 0.2 to 50ng	Reporting Limit (ug/m <sup>3</sup> ) 45L sample volume	Precision (%RSD) n=6
1,3-Butadiene	0.9961	0.0111*	1.89
Benzene	0.9971	0.0044	0.90
Toluene	0.9991	0.0044	0.94
Ethyl Benzene	0.9989	0.0044	0.77
m & p - Xylenes	15.54%	0.0044	0.95
o - Xylene	0.9994	0.0044	1.57
Naphthalene	25.07%	0.0044	0.92
2-Methylnaphthalene	11.79%	0.0044	1.69
1-Methylnaphthalene	19.05%	0.0044	0.65
Acenaphthylene	11.32%	0.0044	1.87
Acenaphthene	14.40%	0.0044	1.48
Fluorene	20.96%	0.0044	2.27
Phenanthrene	8.13%	0.0044	1.67
Anthracene	15.54%	0.0044	2.27
Fluroanthene	7.23%	0.0044	1.41
Pyrene	22.44%	0.0044	1.24
Benzo[a]anthracene	18.93%	0.0044	2.04
Chrysene	19.21%	0.0044	1.92
Benzo[b&k]fluoranthene	16.21%	0.0044	5.96
Benzo[e]pyrene	16.61%	0.0044	0.80
Benzo[a]pyrene	10.86%	0.0044	0.99
Indeno[1,2,3-c,d]pyrene	20.28%	0.0044	1.78
Dibenz[a,h]anthracene	0.9951	0.0044	1.21
Benzo[g,h,i]perylene	0.9952	0.0044	1.97



Compound 12 name: o - Xylene  
Coefficient of Determination: 0.999930  
Calibration curve: 0.292679 \* x + 0.00222540  
Response type: Internal Std ( Ref 1 ), Area \* ( IS Conc. / IS Area )  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



Compound 22 name: Pyrene  
Coefficient of Determination: 0.999243  
Calibration curve:  $1.90286 \cdot x + -0.239747$   
Response type: Internal Std ( Ref 4 ), Area \* ( IS Conc. / IS Area )  
Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None



# The Analytical Solution for Air Monitoring




... efficient, environmentally friendly and cost effective approach




***Breakthrough and Recovery  
Experiments and Results***

- ▶ Occurs when target compounds are not adsorbed by adsorbents
- ▶ EPA TO-17 definition: **“The volume sampled when the amount of analyte collected in a back-up sorbent tube reaches a certain percentage (typically 5%) of the total amount collected by both sorbent tubes”**



- 8000ng on tube
- 85% humidified N<sub>2</sub> for 100 min at a flow of 100 mL/min
- Total volume 10 L sampling
- All studies performed in triplicate

- ▶ A primary TD tube was attached to a gaseous standard to continuously deliver target compounds (mimics a real-world sampling event)
- ▶ A BT tube was attached and monitored on a regular basis
- ▶ Ultimately, the primary tube was loaded with >200mg analyte with no detectable breakthrough




non-detectable breakthrough on all compounds first target 1,3-butadiene

## ► Recovery/Carryover Experiments

- Analyzed spiked tube (50ng)
- Analyzed trap
- Analyzed valve
- Re-analyzed spiked tube

Target Analyte	Trap Test	Tube Test	Valve Test
Benzene	nd	nd	nd
Toluene	nd	nd	nd
Ethylbenzene	nd	nd	nd
m&p-Xylene	nd	nd	nd
o-Xylene	nd	nd	nd
Naphthalene	nd	nd	nd
2-Methylnaphthalene	nd	nd	nd
Acenaphthylene	nd	nd	nd
Acenaphthene	nd	nd	nd
Fluorene	nd	nd	nd
Phenanthrene	nd	nd	nd
Fluoranthene	nd	nd	nd
Chrysene	nd	nd	nd
Benzo[a]pyrene	nd	nd	nd
Indeno[1,2,3-cd]pyrene	nd	nd	nd
Benzo[g,h,i]perylene	nd	nd	nd





# The Analytical Solution for Air Monitoring




... efficient, environmentally friendly and cost effective approach

## ***Site Study, 2012-2013***


TO-17 data was collected in  
simultaneous Full Scan/SIM mode  
only Full Scan data is presented



- ▶ Compare TO-13 / TO-15 to TO-17 from an active MGP remediation site
- ▶ 72-hour sample collection
- ▶ Continuous sampling for six weeks
- ▶ One sample location selected (AMS-03)
- ▶ Three 72-hour samples from each site were selected for comparison
- ▶ 14 sampling sets







- ▶ Two types of tubes investigated (XRO-440 and XRO-640)
- ▶ Each type of tube was sampled in duplicate
- ▶ One of the duplicates had filter attached and analyzed
- ▶ A breakthrough tube was attached to every tube sampled



## 72 Hour Sampling Duration

- ▶ TO-13 = ~1,000,000 Liters (1000 m<sup>3</sup>)
- ▶ TO-15 = 6 Liters
- ▶ TO-17 = ~45 Liters



# Moisture: Hydrophobic adsorbents

Tube	Time for Dry Purge
Sample Tube Type 1	nd added moisture on tube
Sample Tube Type 2	2 minute Dry Purge

Only slight water retention with  
45L sample volume!





- ▶ There was non-detectable targets on the breakthrough tubes from the site studies with an average of 45 liter sample volume


Target Analyte	ug/m3 (first tube)	ug/m3 (second tube)	equates to ng on tube	% relative dif
1,3-Butadiene	0.0193	0.0200	0.87	3.6
Benzene	0.3827	0.4090	17.2	6.6
Toluene	1.3329	9.7246	60.0	152
Ethyl Benzene	0.1543	0.2136	6.94	32.2
m & p - Xylenes	0.4413	0.6081	19.9	31.8
o - Xylene	0.1434	0.1586	6.45	10.1
Naphthalene	3.1182	3.4084	140	8.9
2-Methylnaphthalene	0.6185	0.6083	27.8	1.7
1-Methylnaphthalene	0.2647	0.3138	11.9	17.0
Acenaphthylene	0.0656	0.0492	2.95	28.6
Acenaphthene	0.3022	0.2660	13.6	12.8
Fluorene	0.1238	0.1768	5.57	35.3
Phenanthrene	0.0547	0.0931	2.46	52.0
Anthracene	0.0803	0.0915	3.61	13.1
Fluroanthene	0.0040	0.0050	0.18	21.7
Pyrene	0.0032	0.0050	0.14	45.4
Benzo[a]anthracene	0.0067	0.0054	0.30	21.0
Chrysene	0.0046	0.0033	0.21	32.7
Benzo[b&k]fluoranthene	0.0044	nd	0.20	
Benzo[e]pyrene	0.0044	nd	0.20	
Benzo[a]pyrene	0.0074	nd	0.33	
Indeno[1,2,3-c,d]pyrene	nd	nd		
Dibenz[a,h]anthracene	nd	nd		
Benzo[g,h,i]perylene	nd	nd		

Sample ID Analyte	AMS-03-081313		AMS-03-081613		AMS-03-082213	
	TO13 & TO15	TO17	TO13 & 15	TO17	TO13 & 15	TO17
Benzene	0.52	0.54	1.20	0.56	1.2	0.48
Ethyl Benzene	1.4	0.29	1.4	0.29	3.2	0.78
Toluene	2	1.9	1.6	1.1	4.4	0.86
m,p-Xylene	2.8	0.87	2.8	0.83	6.5	0.91
o-Xylene	1.4	0.22	1.4	0.24	3.2	0.26
Naphthalene (TO15)	3.2	6.9	3.3	1.1	3.9	2.4
Naphthalene (TO13)	1.2	6.9	0.95	1.3	0.18	2.4
2-Methylnaphthalene	0.48	0.74	0.51	0.41	0.10	0.67
1-Methylnaphthalene	0.25	0.36	0.26	0.25	0.088	1.2
Acenaphthylene	0.023	0.12	0.0094	0.047	0.17	0.40
Acenaphthene	0.27	0.25	0.30	0.21	0.17	1.7
Fluorene	0.14	0.10	0.17	0.081	0.13	0.45
Phenanthrene	0.16	0.077	0.17	0.019	0.24	0.24
Anthracene	0.023	0.11	0.0061	0.0258	0.016	0.36
Fluoranthene	0.023	0.0092	0.017	0.0047	0.13	0.013
Pyrene	0.023	0.0043	0.0088	0.0047	0.027	0.0018
Benzo(a)anthracene	0.023	0.0051	0.0011	0.0047	0.00074	0.00029
Chrysene	0.023	0.0051	0.0016	0.0047	0.0014	0.0046
Benzo(b+k)fluoranthene	0.046	0.010	0.0027	0.0094	0.00091	0.0092
Benzo(e)pyrene	0.023	0.0051	0.0011	0.0047	0.00050	0.0046
Benzo(a)pyrene	0.023	0.0051	0.0010	0.0066	0.00050	0.0046
Indeno(1,2,3-cd)pyrene	0.023	0.0051	0.0009	0.0065	0.00050	0.0046
Dibenz(a,h)anthracene	0.023	0.0051	0.00050	0.0140	0.00050	0.0046
Benzo(g,h,i)perylene	0.023	0.0051	0.0013	0.0220	0.00050	0.015


Note: yellow cells are non-detect with the reporting limit value for that target in the cell

AMS-03-081313


Non-detect TO15/TO13 TO-17



Note: dotted, unfilled bars are non-detects. Value represents reporting limit.



Note: dotted, unfilled bars are non-detects. Value represents reporting limit.



Note: dotted, unfilled bars are non-detects. Value represents reporting limit.




... efficient, environmentally friendly and cost effective approach

## ***Other Projects and Site Studies***



# Site Experiments

Site	Investigating (Research)	Sample Volume
Indoor Air	Comparing PAHs in Sub-Slab to Indoor Air	1 liter
Industrial Sub-Slab (MI)	2-Methyl Naphthalene was the compound of interest comparing recovery using TO-15 to TO-17 (XRO-440 tube)	1 liter
Industrial Sub-Slab (IN)	Same site sampling onto SVI tube and XRO-440 tube	1 liter



Soil Vapor Intrusion™ Tubes

## Indoor Air Site – Soil Vapor Intrusion Concern (conc $\mu\text{g}/\text{m}^3$ )

Location	Stock Room		Sewing Room		Maint Room
	Sub-slab	Indoor Air	Sub-slab	Indoor Air	Drain
<b>Target Analyte (SVOC)</b>					
Naphthalene	0.251	3.96	0.286	1.15	1.46
2-Methylnaphthalene	0.345	0.802	0.414	0.252	0.900
1-Methylnaphthalene	0.284	0.408	0.234	0.138	0.529
Acenaphthylene	0.0773	nd	nd	nd	0.180
Acenaphthene	0.157	nd	nd	nd	nd
Fluorene	0.170	0.124	0.166	0.210	0.119
Phenanthrene	0.158	0.111	nd	1.40	0.0794
Anthracene	nd	0.0335	nd	0.240	nd
Fluoroanthene	0.0208	nd	0.0336	0.145	nd
Pyrene	0.0286	0.0192	nd	0.238	nd
Benzo[a]anthracene	nd	nd	nd	0.0334	nd
Chrysene	nd	nd	nd	0.0315	nd
Benzo[b]fluoranthene	nd	nd	nd	nd	nd
Benzo[k]fluoranthene	nd	nd	nd	nd	nd
Benzo[e]pyrene	nd	nd	nd	nd	nd
Benzo[a]pyrene	nd	nd	nd	nd	nd
Indeno[1,2,3-c,d]pyrene	nd	nd	nd	nd	nd
Dibenz[a,h]anthracene	nd	nd	nd	nd	nd
Benzo[g,h,i]perylene	nd	nd	nd	nd	nd

## Results from Site in Michigan: Soil Gas (XRO-440 Tubes)

Component	Reporting Limit (RL)	001B	004B	022B	002B	009B	016B
Naphthalene	0.30	1.39	1.52	0.81	0.453	0.414	1.1
2-Methylnaphthalene	0.30	0.818	0.253	0.189	0.247	0.129	1.2
1-Methylnaphthalene	0.30	0.415	0.146	0.116	0.14	0.117	0.631
Acenaphthylene	0.30	0.0682	0.0338	0.0458	0.043	0.0232	0.114
Acenaphthene	0.30	0.0527	0.0178	0.0246	0.015	0.00897	0
Fluorene	0.30	0.0542	0.0162	0.0455	0.0427	0.0312	0.051
Phenanthrene	0.30	0.212	0.174	0.214	0.211	0.247	0.546
Anthracene	0.30	0.063	0.0416	0.049	0.0455	0.0619	0.0933
Fluoranthene	0.30	0.079	0.0873	0.0892	0.116	0.136	0.128
Pyrene	0.30	0.0976	0.05	0.0558	0.0638	0.0576	0.057
Benzo(a)anthracene	0.30	0.165	0.0465	0.0729	0.0935	0.109	0.103
Chrysene	0.30	0.108	0.0623	0.0768	0.0833	0.0682	0.0576
Benzo(b)fluoranthene	0.30	0.568	0.547	0.179	0.45	0.435	0.356
Benzo(k)fluoranthene	0.30	0.426	0.413	0.164	0.371	0.25	0.13
Benzo(e)pyrene	0.30	0.193	0.148	nd	0.23	0.229	0.195
Perylene	0.30	nd	0.4	nd	nd	nd	nd
Benzo(a)pyrene	0.30	0.247	nd	nd	nd	nd	nd
Ideno(1,2,3-cd)pyrene	0.30	nd	nd	nd	nd	nd	nd
Dibenz(a,h)anthracene	0.30	nd	nd	nd	nd	nd	nd
Benzo(g,h,i)perylene	0.30	nd	nd	nd	nd	nd	nd

... five (5) samples of 28 are shown. All had similar results

# Indiana Site: TO-17 (SVI vs XRO-440 Tubes)

Sample			VP - 105		VP - 106		VP - 111		VP - 114		VP - 115	
Compound	SVI RL (ng)	XRO-440 RL (ng)	SVI	XRO - 440	SVI	XRO - 440	SVI	XRO - 440	SVI	XRO - 440	SVI	XRO - 440
Naphthalene	5.0	0.30	2.04	0.999	2.16	1.87	2.25	2.1	2.15	1.1	8.46	6.55
2-Methylnaphthalene	2.0	0.30	1.72	0.439	1.37	1.32	0.612	0.71	0.474	0.742	20.9	20
1-Methylnaphthalene	2.0	0.30	1.11	0.262	1.08	0.78	0.604	0.49	0.38	0.714	12.8	13
Acenaphthylene	2.0	0.30	0.0278	0.0381	0.0327	0.0828	0.0238	0.0514	0.0776	0.0932	1.86	0.704
Acenaphthene	2.0	0.30	0	0.051	0.186	0.118	0.136	0.0382	0.101	0.0465	1.2	0.947
Fluorene	2.0	0.30	0.28	0.0678	0.281	0.106	0.213	0.0676	0.172	0.043	0.795	0.416
Phenanthrene	2.0	0.30	0.249	0.196	0.242	0.347	0.255	0.237	0.179	0.147	0.42	0.358
Anthracene	nd	0.30	nd	nd	nd	0.0622	nd	nd	nd	nd	nd	0.136
Fluoranthene	nd	0.30	nd	0.0576	nd	0.0591	nd	0.0579	nd	0.0539	nd	0.0363
Pyrene	nd	0.30	nd	0.162	nd	0.161	nd	0.15	nd	0.101	nd	0.173
Benzo(a)anthracene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Chrysene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	0.0386	nd	nd
Benzo(b)fluoranthene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Benzo(k)fluoranthene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Benzo(e)pyrene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Perylene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Benzo(a)pyrene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Ideno(1,2,3-cd)pyrene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Dibenz(a,h)anthracene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd
Benzo(g,h,i)perylene	nd	0.30	nd	nd	nd	nd	nd	nd	nd	nd	nd	nd

... five (5) samples of 18 are shown. All had similar results

- ▶ EPA Method TO-17 is performance-based, guidance method
  - **Section 2.5 states:** “...This method provides performance criteria to demonstrate acceptable performance of the method (or modifications of the method) for monitoring a compound or set of compounds.”
  
- ▶ EPA has seen this data and has given verbal acceptance stating that TO-17 is performance based so targets may be included as long as criteria is met
  - U.S. Environmental Protection Agency  
Office of Air Quality Planning and Standards  
Ambient Air Monitoring Group C304-06  
Research Triangle Park, NC 27711

- ▶ Analytical performance proves concept
- ▶ Site data suggests this is a better alternative
- ▶ One analysis instead of two:
  - Reduce sampling and analytical costs and disposal
  - Save on shipping and labor costs
  - Enhance productivity and efficiency
  - Increase profits
  - **Better for our environment ... A Greener analysis**
- ▶ More data is available

- ▶ Pace Analytical Services
  - Amy Jacobson, manager
  - Nathan Eklund, Program Manager
  - Kesler Krieg, Sales Manager for Air
- ▶ PerkinElmer Instruments
  - James Day, Sr Service Engineer

Thank you



# The Analytical Solution for Air Monitoring



... efficient, environmentally friendly and cost effective approach

***Questions please ????***

Contact Info: [lee.marotta@perkinelmer.com](mailto:lee.marotta@perkinelmer.com)

