

## Quantification of anti-mTNF-alpha ADCC activity using *iLite*<sup>®</sup> ADCC Assay Ready Cells

*For research and professional use only. Not for use in diagnostic procedures.*

*This application note contains a suggested protocol and performance data. Each individual laboratory must set up their own method and perform relevant validations.*

### Background

Antibody-dependent cell-mediated cytotoxicity (ADCC) is a mechanism whereby pathogenic cells are lysed by lymphocytes, most often Natural Killer (NK) cells. The mechanism involves binding of antibodies to surface antigens on the pathogen. Crosslinking of these antibodies to NK cells through the binding of the Fc-portion to Fc receptors on the NK cells leads to activation of the NK cell and formation of an immune synapse with the pathogenic cell. The NK cell releases cytotoxic granules containing granzymes and perforin into the synapse, leading to apoptosis of the targeted cell (1).

Monoclonal anti-TNF-alpha inhibitors are one of the most common types of antibody drugs on the market. The mechanism by which these drugs act is most often blocking of soluble TNF-alpha or the TNF receptor. Since TNF-alpha is produced through cleavage of a membrane bound form of the protein, there is a risk that anti-TNF-alpha antibodies bind the membrane bound TNF-alpha and trigger effector functions, such as ADCC. By using the *iLite* ADCC mTNF-alpha Activity Assay, the risk of an anti-TNF-alpha antibody drug inducing ADCC can be quantified.

### Principle of the assay


The *iLite*<sup>®</sup> ADCC Assay Ready Cells are engineered cells that enable antibody-dependent cell-mediated cytotoxicity (ADCC) to be examined through the specific expression of Firefly luciferase. When the antibodies of interest bind to the antigens on the surface of the target cell, the target-bound antibodies will be presented to the Fc receptors (FcγRIIIa) on the effector cell. When the Fc-portion of the target-bound antibodies binds to the receptor, multiple cross-linking of the two cell types occurs. This will initiate a signaling cascade which triggers the expression of Firefly luciferase (FL) in the effector cell. In this application note, we describe the use of an effector cell line (*iLite*<sup>®</sup> ADCC Effector (V) Assay Ready Cells) that over-express FcγRIIIa and contain the FL reporter gene that responds to the principal transcription factors that mediate signaling from the FcγRIIIa receptor, together with a positive target cell line which over-expresses membrane bound TNF-alpha (*iLite*<sup>®</sup> mTNF-alpha (+) Target Assay Ready Cells). *iLite*<sup>®</sup> ADCC Effector (V) Assay Ready Cells also contain the Renilla Luciferase (RL) reporter gene, under the control of a constitutive promoter, that allows anti-TNF drug-induced FL activity to be normalized with respect to the constitutive expression of RL. This renders assay results independent of variations in cell number, serum matrix effects, or lysis of the effector cells by the target cells. In addition, we also describe the use of a negative control in the form of a target cell line without expression of membrane bound TNF-alpha (*iLite*<sup>®</sup> mTNF-alpha (-) Target Assay Ready Cells). The Firefly luciferase signal can be measured in a luminometer following addition and incubation of luciferase substrate. The Firefly luciferase signal is proportional to the functional activity of Infliximab in the sample (Fig.1).

## Material and equipment needed

| Material and equipment  | Suggested supplier  | Reference   |
|---|---|---|
| <i>iLite</i> <sup>®</sup> ADCC Effector (V) Assay Ready Cells                       | Svar Life Science   | BM5001  |
| <i>iLite</i> <sup>®</sup> mTNF-alpha (+) Target Assay Ready Cells                   | Svar Life Science   | BM5013  |
| <i>iLite</i> <sup>®</sup> mTNF-alpha (-) Target Assay Ready Cells                   | Svar Life Science   | BM5014  |
| Diluent (RPMI 1640 + 9% heat inactivated FBS + 1% Penicillin Streptomycin)          | Gibco   | 61870 (RPMI)<br>26140-079 (FBS)<br>15140-122<br>(Penicillin-Streptomycin) |
| Infliximab or analogues   | Janssen Biotech   | NA  |
| Firefly/Renilla luciferase substrate  | Promega   | E2940, Dual-Glo <sup>®</sup> Luciferase Assay System                      |
| Plate; White walled micro well plate suitable for luminescence                      | PerkinElmer   | 6005680   |
| Microplate Luminometer with appropriate reading software – no filter on luminometer | Contact Svar Life Science for list of recommended suppliers | NA  |
| Incubator, 37°C with 5% CO <sub>2</sub>   | NA  | NA  |
| Water bath, 37°C  | NA  | NA  |
| Single-channel and multi-channel pipettes with polypropylene disposable tips        | NA  | NA  |
| Polypropylene tubes or plate for dilution   | NA  | NA  |
| Single-use polypropylene reservoir  | NA  | NA  |
| Plate shaker  | NA  | NA  |
| Timer   | NA  | NA  |

## Preparation of calibrators (example given with infliximab)

The ADCC effect of the infliximab antibody from Janssen Biotech has successfully been measured in combination with a mix of ADCC Effector (V) Assay Ready Cells and mTNF-alpha (+) Target Assay Ready Cells. As a negative control, a combination of ADCC Effector (V) Assay Ready Cells and mTNF-alpha (-) Target Assay Ready Cells was used. In the present assay an Effector: Target ratio of 4:1 has been used. The optimal ratio is dependent on the antibody and target cells used and should be determined each time a new assay is set up. The table below shows recommended dilutions of Infliximab when making an 8-point calibration curve.



**Figure 1.** Example of Infiximab calibration curve using Firefly Luciferase substrate Dual-Glo® Luciferase Reagent. Values are shown as mean of triplicate  $\pm$  SD and values on x-axis are given as **final concentration** in the wells before addition of Dual-Glo® Luciferase Reagent.

| Calibrator | Infiximab  |                        |
|------------|------------|------------------------|
|            | Calibrator | solution conc. (ng/mL) |
| 1          |            | 1 000                  |
| 2          |            | 333                    |
| 3          |            | 111                    |
| 4          |            | 37                     |
| 5          |            | 12                     |
| 6          |            | 4.1                    |
| 7          |            | 1.4                    |
| 8          |            | 0                      |

**Table 1.** Suggested calibrator **solution concentrations** for Infiximab.

## Protocol

### Assay preparation and incubation

1. Design a plate layout.
2. Dilute calibrators, controls and samples to fall within the expected **in assay values (= final concentration)** of 0-500 ng/mL.
3. Add 40  $\mu$ L calibrators, controls and samples in duplicate to assigned wells.
4. Thaw the vial of ADCC Effector (V) Assay Ready Cells and the vials of mTNF-alpha (+) Target Assay Ready Cells and mTNF-alpha (-) Target Assay Ready Cells in a 37°C water bath with gentle agitation<sup>1</sup>.
5. Mix the cell suspensions very carefully **at least 10 times with a pipette** in order to ensure a homogeneous distribution of cells.
6. Dilute 200  $\mu$ L of the ADCC Effector (V) Assay Ready Cells and 200  $\mu$ L the mTNF-alpha (+) Target Assay Ready Cells with 3.44 mL Diluent. The total volume of the diluted ADCC Effector (V) /Target mTNF-alpha (+) Assay Ready Cells mixture is 3.84 mL.
7. In a separate tube, dilute 50  $\mu$ L of the ADCC Effector (V) Assay Ready Cells with 50  $\mu$ L of the mTNF-alpha (-) Target Assay Ready Cells with 860  $\mu$ L Diluent. The total volume of the diluted ADCC Effector (V) /Target mTNF-alpha (-) Assay Ready Cells mixture is 960  $\mu$ L.
8. Add 40  $\mu$ L of the diluted cells to each well to be tested.
9. Place the lid on the plate and mix on a plate shaker at **minimum of 750 rpm** for 10 sec. Alternatively, mix the cell suspensions very carefully in the wells by pipette. Insufficient mixing can cause reduced assay sensitivity.
10. Incubate for 4 hours at 37°C with 5% CO<sub>2</sub>.

### Adding substrate solutions


11. Equilibrate the plate and the substrate solutions to room temperature.
12. Prepare the **Firefly luciferase substrate** in accordance with the supplier's instructions and add 80  $\mu$ L per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.

<sup>1</sup> The color of the cell suspension might vary from salmon pink to a more yellowish red. The difference in color is due to small variations in pH of the solution and do not affect function of the cells.

13. If appropriate, prepare the **Renilla luciferase substrate** in accordance with the supplier's instructions and add 80  $\mu\text{L}$  per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.

## Normalization

The reporter gene used for normalization, Renilla luciferase, is under the control of a tyrosine kinase promoter, and is thus constitutively expressed. Unspecific effects such as serum matrix effects or differences in cell number can be obviated by relating the specific Firefly signal with the Renilla normalization signal through simple division.



## Precautions

- This application note is intended for professional laboratory research use only. The data and results originating from following the Application Note should not be used either in diagnostic procedures or in human therapeutic applications.
- Use and handle the material and instruments referenced according to the suppliers'/manufacturers' instructions or product specifications accompanying the individual material and instruments.
- Dispose of all sample specimens, infected or potentially infected material in accordance with good microbiological practice. All such materials should be handled and disposed as though potentially infectious.
- Residues of chemicals and preparations are generally considered as biohazardous waste and should be inactivated prior to disposal by autoclaving or using bleach. All such materials should be disposed of in accordance with established safety procedures.


## Proprietary Information

In accepting delivery of *iLite*<sup>®</sup> Assay Ready Cells the recipient agrees not to sub-culture these cells, attempt to sub-culture them or to give them to a third-party recipient, and only to use them directly in assays. *iLite*<sup>®</sup> cell-based products are covered by patents which are the property of Svar Life Science

AB and any attempt to reproduce the delivered *iLite*<sup>®</sup> Assay Ready Cells is an infringement of these patents.

## QUICK GUIDE

### Quantification of anti-mTNF-alpha ADCC activity using *iLite*<sup>®</sup> ADCC Assay Ready Cells



## Troubleshooting and FAQ

Please consult the Svar Life Science website [www.svarlifescience.com](http://www.svarlifescience.com)

## References

1. **Weiner GJ.** *Building better monoclonal antibody-based therapeutics.* Nat Rev Cancer 15(6): 361-70 (2015).