## EDIA™ Anti-CCP

For professional use only  
Usage réservé aux professionnels  
Sólo para uso profesional  
Nur für den Fachgebrauch  
Solo per uso professionale

Document No. LABEL-DOC-0040 2.0

### EDIA™ Anti-CCP

English:	page.....	2
Français:	page.....	16
Español:	página.....	30
Deutsch:	Seite .....	44
Italiano:	pagina.....	59

FCCP 100

## **INTENDED USE**

The EDIA™ Anti-CCP test kit is an enzyme-linked immunosorbent assay (ELISA) for detection and semi-quantitation of IgG antibodies to Cyclic Citrullinated Peptides (CCP) in human serum or plasma. The assay is used to detect antibodies in a single specimen. The results of the assay are to be used as an aid to the diagnosis of Rheumatoid Arthritis (RA), in conjunction with other laboratory and clinical findings. The analysis should be performed by trained laboratory professionals. For in vitro diagnostic use.

## **SUMMARY AND EXPLANATION OF THE TEST**

Rheumatoid Arthritis (RA) is one of the most common systemic autoimmune diseases. The aetiology of the disease, which affects up to 1-2% of the world population, is unknown. The diagnosis of RA depends primarily on clinical manifestation of the disease. The only serological test routinely used is the determination of the presence of rheumatoid factors (RF) in the serum. RF are antibodies directed to the constant region of immunoglobulins of the IgG class. However, these antibodies are also present in relatively high percentages in other autoimmune diseases, infections and in up to 15% of healthy individuals.

Antibodies of a more specific nature have also been found in sera of RA patients (see (1) for an overview). Anti-perinuclear factor (APF) antibodies are reported to be present in around 50% of RA patients with a specificity of over 70% (2). A number of cyclic synthetic peptides not related to filaggrin or other known proteins have been described which are specifically recognized by autoantibodies in sera from RA patients (3). These peptides were subsequently used in an ELISA for the detection of RA-specific autoantibodies (3). Clinical evaluation studies showed that the ELISA was positive in a significant number of well-defined RA patient sera with a high specificity against disease controls (3-8). A diagnostic and prognostic value for the measurement of the anti Cyclic Citrullinated Peptides (anti-CCP) antibodies was found in relation to joint involvement and radiological damage in early RA (7, 9-14). The EDIA™ anti-CCP assay offered by Svar Life Science is based on highly purified synthetic peptides containing citrulline residues and is an additional tool in the diagnosis of RA. This anti-CCP kit contains improved synthetic peptides selected on the basis of performance in the detection of RA autoantibodies (8-14).

## **PRINCIPLE OF THE ASSAY**

The EDIA™ anti-CCP antibody kit is based on an ELISA method. The test utilizes microtitre wells coated with citrullinated synthetic peptides (antigen). Diluted patient serum or plasma is applied to the wells and incubated. If specific antibodies are present, they will bind to the antigen in the wells. Unbound material is washed away and any bound antibody is detected by adding alkaline phosphatase labelled antibodies to human IgG, followed by a second washing step and incubation with a chromogenic substrate.

The presence of reacting antibodies will result in the development of colour, which is proportional to the quantity of bound antibody, and this is determined photometrically

## KIT COMPONENTS

<b>A</b>	IgG Conjugate	1 x 15 mL	Alkaline phosphatase-labelled goat polyclonal antibody to human IgG, Tris buffer, protein stabiliser, <0.1% (w/v) sodium azide. <b>Ready-to-use.</b>
<b>B</b>	Substrate	1 x 15 mL	Mg <sup>2+</sup> , phenolphthalein monophosphate (PMP), buffer solution. <b>Ready-to-use.</b> Do not expose to light during storage. 
<b>C</b>	Stop solution	1 x 15 mL	Sodium hydroxide, EDTA, carbonate buffer (pH>10). <b>Ready-to-use.</b> 
<b>D</b>	Wash Buffer Concentrate (16x)	2 x 25 mL	Borate buffer, 0.4% (w/v) sodium azide. <b>Dilute before use.</b> 
<b>E</b>	CCP-Coated Wells in a Frame	96 wells	Coated with CCP-antigen, in a resealable foil pack with desiccant.
<b>F</b>	Sample diluent concentrate (5x)	1 x 25 mL	Phosphate buffer, protein stabiliser, 0.5% (w/v) sodium azide. <b>Dilute before use.</b> 
<b>1-6</b>	Anti-CCP Calibrators	6 x 1.0 mL	Human plasma, buffer, <0.1% (w/v) sodium azide. 0, 2, 8, 30, 100, 300 U/mL. <b>Ready-to-use.</b>
<b>7</b>	Anti-CCP Reference Control	1 x 1.5 mL	Human plasma, buffer, <0.1% (w/v) sodium azide. <b>Ready-to-use.</b>
<b>+/-</b>	Positive Control Negative Control	1 x 1.5 mL 1 x 1.5 mL	Human plasma, buffer, <0.1% (w/v) sodium azide. <b>Ready-to-use.</b>
	Instruction for use (IFU)		

## STORAGE OF REAGENTS

### Handling and Procedural Notes

1. Store kit components at 2-8° C and use until the expiry date on the labels. Do not use expired reagents.
2. Do not mix different lot numbers.
3. Do not freeze kits.
4. Wash Buffer Concentrate and Sample Diluent Concentrate must be diluted before use. All other reagents are ready-to-use.
5. Diluted Wash Buffer and diluted Sample Diluent are stable at 2-8° C for up to 6 months if microbial contamination is avoided.
6. Replace surplus microtitre wells in the foil pack and store within the resealable aluminium foil bag with the desiccant at 2-8° C, until required.
7. Do not expose Substrate to light during storage.
8. Avoid contamination of reagents. Use a new disposable pipette tip for each reagent or sample manipulation.
9. Optimal results will be obtained by strict adherence to this protocol. Careful pipetting and washing throughout the procedure are necessary to maintain precision and accuracy.

### Indications of Deterioration

The Substrate should be pale yellow in colour. Pink colouring indicates contamination and the reagent must be discarded. Turbidity or precipitation in any component indicates deterioration and the component should be discarded.

**Sample Collection and Storage**

The assay is recommended for serum/plasma (EDTA and Heparin) samples; do not use grossly haemolysed lipaemic or turbid samples. Thoroughly mix thawed samples before assay and avoid repeated freeze/thawing. Store samples for a maximum of 48 hours at 4-8° C. For prolonged storage freeze samples at -20° C.

**WARNINGS AND PRECAUTIONS**

**For in vitro diagnostic use only.**

**Safety Precautions**

1. Adhere strictly to the instructions in this booklet, particularly for handling and storage conditions.
2. Controls and calibrators contain serum of human origin. Although tested against and confirmed negative for HIV 1+2, HCV, HbsAg and HIV-1 Ag, this material must be treated as potentially infectious. - The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
3. Do not pipette by mouth.
4. Do not smoke, eat, drink or apply cosmetics in areas where kits and samples are handled.
5. Any skin complaints, cuts, abrasions and other skin lesions should be suitably protected.
6. The Calibrators, Controls, Conjugate, Sample Diluent Concentrate and Wash Buffer Concentrate contain sodium azide which can react with lead and copper plumbing to form highly explosive metal azides. On disposal, drain with large quantities of water to prevent azide build-up.
7. The Stop Solution contains sodium hydroxide. Avoid contact with skin, eyes and mucous membranes. Spillage should be mopped up with copious amounts of water. If contact with skin or eyes occurs, irrigate with water and seek medical attention immediately.
8. The substrate contains PMP, Bronidox L and Diethanolamine. Avoid contact with skin, eyes and respiratory system. If contact with skin, eyes or respiratory system occurs, rinse with water and seek medical advice.
9. Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

**Warning**

**B.**

SUBS

Contains: Diethanolamine

H319:	Causes serious eye irritation.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313:	If eye irritation persists: Get medical advice/attention.

C.

SOLN	STOP
------	------

**Warning**

Contains: Sodium hydroxide

H315:	Causes skin irritation.
H319:	Causes serious eye irritation.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352:	IF ON SKIN: Wash with plenty of soap and water.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P332+P313:	If skin irritation occurs: Get medical advice/attention.
P337+P313:	If eye irritation persists: Get medical advice/attention.

D. and F.

BUF	WASH	16X
DIL	SPE	5X

**Warning**

Contains: Sodium azide

H302:	Harmful if swallowed.
EUH032:	Contact with acids liberates very toxic gas.
H412:	Harmful to aquatic life with long lasting effects.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P301+P312:	IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
P273:	Avoid release to the environment.

## PREPARATION

### **Materials/Equipment Required but not Provided**

1. 96 well plate/strip reader with 550nm filter (540-565nm is acceptable).
2. Precision pipettes to dispense 10  $\mu$ L, 100  $\mu$ L, 1 mL. Automatic pipette to dispense 100  $\mu$ L. Automatic pipette to dispense 200  $\mu$ L for manual washing, automatic plate washer optional.
3. Glass/plastic measuring cylinders: 1 $\times$ 100 mL, 1 $\times$ 400 mL.
4. 1 mL volume containers.
5. Distilled/deionised water.
6. Paper towels.
7. Timer for 30 and 60 minutes intervals.

**Preparation for the Assay**

Before use, bring all kit components, including the microtitre wells, to room temperature (18-25° C) Mix reagents by gentle inversion.

**Do not dilute Reference, Positive and Negative Controls.**

Dilute the following reagents and samples and mix thoroughly.

Reagent	Volume	Add
Wash Buffer Concentrate	1 vial	375 mL distilled/deionised water
Sample Diluent Concentrate	1 vial	100 mL distilled/deionised water

Sample	Volume	Add
Sample	10 µL	1 mL diluted Sample Diluent

Calculate the number of microtitre wells required for the current assay, and retain these in the microtitre frame. Return surplus wells to the foil pack with desiccant, and store at 2-8° C until required. Ensure that all wells are securely held within the microtitre frame, with the assay identification tab along the bottom edge below row H.

**ASSAY PROTOCOL**

**Qualitative protocol:** run Reference Control, Positive and Negative Controls, and samples.

**Semi-quantitative protocol:** run Calibrators (1-6), Positive and Negative Controls, and samples.

- Reference wells for identification.
- Pipette 100 µL Reference, Positive and Negative Controls/Calibrators in duplicate, and pre-diluted (1:101) patient samples into appropriate wells. Remember to change pipette tips between additions. This step should not exceed 15 minutes for any one set of Calibrators/Controls/samples.
- Incubate 60±10 minutes at 18-25° C.
- Decant well contents by quick inversion over a sink suitable for the disposal of biological materials, bearing in mind the potential infective hazard of the samples. Blot inverted wells with paper towels.
- Wash wells **three times** with a minimum of 200 µL diluted Wash buffer. **Decant and blot after each wash step.**
- Add 100 µL Conjugate to each well.
- Incubate 30±5 minutes at 18-25° C.
- Repeat steps 4 and 5.
- Add 100 µL Substrate to each well.
- Incubate 30±5 minutes at 18-25° C. **Do not decant.**
- Add 100 µL Stop Solution to each well, in the same order and rate as the Substrate. Tap wells gently to mix.
- Read wells within 24 hours at 550nm (540-565nm)

## CALCULATION AND INTERPRETATION OF RESULTS

Consider each assay separately when calculating and interpreting results.

### Qualitative protocol

Calculate the absorbance value (optical density) ratio for the Positive and Negative Controls, and for each sample.

$$\text{Absorbance Ratio} = \frac{\text{Sample or Control Absorbance Value}}{\text{mean Reference Control Absorbance Value}}$$

Users should calculate a cut-off between positive and negative samples that is specific to their patient populations. Results from the patient populations used in the Svar Life Science clinical trial suggest the following cut-off:

<u>Absorbance Ratio</u>	<u>Results Interpretation</u>
<0.95	Negative
≥0.95 to ≤1.0	Borderline - recommend repeat testing
>1.0	Positive

### Semi-quantitative protocol

Plot the mean absorbance value of each Calibrator against log<sub>10</sub> concentration (see following table) on suitable graph paper. Concentrations of Controls and samples can then be read from the calibration curve; a typical plot employing 4-parameter logistic (4PL) is shown below for reference purposes, it must not be used for interpreting results. Samples with absorbances above Calibrator 6 (300 U/mL) are outside the range of the assay, and should be reported as >300U/mL or diluted and re-assayed, correcting for this further dilution factor.


Note: As in any assay measuring antibodies, this assay determines the activity of the antibody present in the sample, rather than the concentration. Activity can be affected by a number of parameters, such as antibody avidity.

Samples with results ≤5 U/mL are defined as negative. Samples >5 U/mL are defined as positive.

Calibrator Concentrations

Calibrator Number	Concentration U/mL
1	0
2	2
3	8
4	30
5	100
6	300

Typical Standard Curve

## QUALITY CONTROL

Ensure that adequate maintenance and calibration of the plate-reader is performed according to the manufacturer's instructions, and that the correct wavelength is employed.

Users should ensure they are fully acquainted with the instructions for the assay, particularly the Warnings and Precautions section, and the Handling and Procedural Notes.

Users should demonstrate that they can obtain performance specifications for precision and reportable range of test results comparable to those established by the manufacturer before reporting patient test results. It is recommended that the pre-diluted Positive and Negative Controls are run in duplicate in all assays to monitor the quality of the test procedure. Run the ready-to-use Reference Control in duplicate in all qualitative assays.

Assuming the precision specifications described by the manufacturer are met, failure of any Control to meet the Control ratio specifications below renders the assay invalid and patient results should not be reported. The operator may repeat the assay, having reviewed the procedure, or contact the distributor/manufacturer. If repeating the assay, prepare a fresh dilution of the sample. Laboratories may wish to include in-house controls in each assay run. Store such control material at or below - 20° C and avoid repeat freeze/thaw cycles. Preservatives such as sodium azide at 0.1% (w/v) will not affect sample results.

Levels of analytes identified in particular diseases are those established by the manufacturer for specific populations, and may not necessarily mirror the literature. Incidence levels, their relationship for specific diseases, reference ranges, and appropriate cut-off points should all be calculated for the specific populations serviced by users.

### Control Ratio Specifications

Protocol	Specifications
Qualitative (ratios)	$\frac{\text{Positive Control Absorbance}}{\text{Reference Control Absorbance}}$ See Positive Control label
	$\frac{\text{Negative Control Absorbance}}{\text{Reference Control Absorbance}}$ <0.95
Semi-quantitative	See Positive Control label for acceptable expected range (U/mL)
	Negative Control concentration < 3 U/mL

## EXPECTED VALUES

The EDIA™ anti-CCP ELISA measures antibodies against synthetic peptides with citrulline residues (anti-CCP). The EDIA™ anti-CCP ELISA is calibrated in the semi-quantitative assay in arbitrary units using a positive patient serum pool as calibrator. The calibration curve ranges from 0-300 U/mL. These values have been chosen arbitrarily by Svar Life Science since no generally recognised (inter)national standard exists for expressing the titre of anti-CCP antibodies. The specificity and sensitivity were evaluated with 416 RA patients, 531 diseased non-RA patients (including other autoimmune, a range of inflammatory and infectious diseases) and 262 healthy controls, samples from adult subjects. The sensitivity was 76%. The specificity was 98% with diseased non-RA patients and 99% with apparently healthy individuals.

To determine the cut-off level, a Receiver Operator Characteristic (ROC) analysis was performed. The analysis revealed a point with the lowest possible sum of false-positives and false-negative samples and the cut-off was accordingly set at 5 U/mL.

*Reference Range*  
 $\leq 5\text{U/mL}$  = Negative  
 $> 5\text{U/mL}$  = Positive

## PERFORMANCE CHARACTERISTICS

**Table 1. Percent agreement** of the EDIA™ anti-CCP compared to an alternative anti-CCP ELISA. A total of 678 frozen retrospective sera were assayed. 416 samples from RA patients and 262 samples were from apparently healthy blood donors.

EDIA™ anti-CCP	Alternative ELISA		
	Positive	Negative	Total
Positive	317	2	319
Negative	5	354	359
Total	322	356	678

Positive Percent Agreement:  $317/322 = 98.4\%$  95% CI = 96.4 - 99.5%  
 Negative Percent Agreement:  $354/356 = 99.4\%$  95% CI = 98.0 - 99.9%  
 Overall Percent Agreement:  $671/678 = 99.0\%$  95% CI = 97.9 - 99.6%

The 95% confidence interval (CI) was calculated using the exact method.

**Table 2. Clinical sensitivity and specificity.** A total of 1209 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results.

	n	negative	positive
<b>Blood donors</b>	262	260	2
<b>RA</b>	416	99	317
<b>Infectious diseases</b>	86	85	1
<b>Inflammatory diseases (nonRA)</b>	445	437	8
<b>TOTAL</b>	1209	881	328

(Data on file)

### Clinical sensitivity

RA  $317/416 = 76.2\%$  95% CI = 72.1 - 80.3%

### Clinical specificity

Blood donors  $260/262 = 99.2\%$  95% CI = 97.3 - 99.9%  
 Infectious diseases  $85/86 = 98.8\%$  95% CI = 93.7 - 100%  
 Inflammatory diseases (nonRA)  $437/445 = 98.2\%$  95% CI = 96.5 - 99.2%  
 Combined non-RA groups  $782/793 = 98.6\%$  95% CI = 97.5 - 99.3%

The 95% confidence interval (CI) was calculated using the exact method.

**Table 3. Intra-assay precision** was determined by testing six different samples eight times each.

	<b>High</b> (U/mL)	<b>Medium</b> (U/mL)	<b>Low</b> (U/mL)
Mean	173.9	34.0	9.9
S.D.	13.8	0.6	0.2
% C.V.	7.9	1.9	2.1
	<b>Low</b> (U/mL)	<b>Low</b> (U/mL)	<b>Low</b> (U/mL)
Mean	11.8	7.8	9.7
S.D.	0.5	0.1	0.4
%C.V.	4.0	1.9	4.4

**Table 4. Inter-assay precision** was determined by testing six different samples eight times each. Results were obtained for three different runs.

	<b>High</b> (U/mL)	<b>Medium</b> (U/mL)	<b>Low</b> (U/mL)
Mean.	183.8	36.6	9.3
S.D.	19.5	3.0	0.9
% C.V.	10.6	8.2	9.8
	<b>Low</b> (U/mL)	<b>Low</b> (U/mL)	<b>Low</b> (U/mL)
Mean	11.9	7.8	10.6
S.D.	0.8	0.7	0.9
%C.V.	6.3	9.5	8.9

**Table 5. Batch to batch variation** was determined by testing six different samples eight times each. Results were obtained for three different batches.

	<b>High</b> (U/mL)	<b>Medium</b> (U/mL)	<b>Low</b> (U/mL)
Mean	232.5	41.6	11.5
S.D.	30.9	4.5	1.3
% C.V.	13.3	10.8	11.2
	<b>Low</b> (U/mL)	<b>Low</b> (U/mL)	<b>Low</b> (U/mL)
Mean	14.1	9.8	13.0
S.D.	1.2	1.0	1.2
%C.V.	8.3	10.4	9.3

**Table 6. Dilution recovery** was determined by testing five serial dilutions of three different patient samples.

Sample	Dilution	Mean Measured Concentration (U/mL)	Calculated Concentration (U/mL)	Dilution Corrected % Recovery
1	1/100	205.0	205.0	100
	1/200	110.5	102.5	108
	1/400	47.3	51.3	92
	1/800	24.8	25.6	97
	1/1600	10.8	12.8	84
Sample	Dilution	Mean Measured Concentration (U/mL)	Calculated Concentration (U/mL)	Dilution corrected % recovery
2	1/100	138.9	138.9	100
	1/200	70.3	69.5	101
	1/400	40.4	34.7	116
	1/800	18.3	17.4	105
	1/1600	8.7	8.7	100
Sample	Dilution	Mean Measured Concentration (U/mL)	Calculated Concentration (U/mL)	Dilution corrected % recovery
3	1/100	47.3	47.3	100
	1/200	26.7	23.6	113
	1/400	13.0	11.8	110
	1/800	6.3	5.9	107
	1/1600	3.0	3.0	103

**Detection Limit**

The detection limit of the assay was determined by running the zero calibrator 12 times on three different lots. The detection limit of 0.5 U/mL was calculated by finding the mean plus two standard deviations.

**Interference study**

Three low positive samples were spiked to the following concentrations in diluted serum samples; Bilirubin F at 0.188 mg/dL, Bilirubin C at 0.2 mg/dL, Haemoglobin at 453 mg/dL, Chyle at 0.24 U/dL and Rheumatoid Factor at 200 IU/mL. The data indicates that the assayed concentrations do not interfere with the anti-CCP results.

To assess potential cross-reactivity of CCP IgG antigen with other autoantibodies, a total of 405 samples of different aetiology were assayed. Samples from patients diagnosed with Crohn's disease, Colitis ulcerosa, SLE, Sjögren's syndrome, Osteoarthritis, Scleroderma, Multiple sclerosis, MCTD, Inflammatory bowel disease, Polymyositis/Dermatomyositis, nonRA autoimmune patients and samples reacting with MPO-ANCA, PR3-ANCA and ds-DNA. Data indicates that the assayed autoantibodies show no significant cross-reactivity.

## LIMITATIONS OF USE

1. A positive result must be used in conjunction with clinical evaluation and other diagnostic procedures. The values obtained from this assay are intended to be an aid to diagnosis only. Each physician must interpret the results in conjunction with the patient's history, physical findings and other diagnostic procedures.
2. Elevated anti-CCP antibodies may be seen in individuals with no evidence of clinical disease. Also, some individuals with RA may have undetectable antibodies. Anti-CCP antibody levels do not necessarily correlate to disease state.
3. Because anti-CCP antibody levels do not necessarily correlate to disease state treatment should not be initiated or changed based on a positive result. Clinical findings should be taken into account for all treatment decisions.
4. Monitoring anti-CCP antibody levels for progression and or remission of RA has not been established.
5. The performance characteristics for this assay have not been established for paediatric specimens. The diagnostic value of anti-CCP antibodies has not been determined for juvenile arthritis.

## REFERENCES

1. Van Boekel M, et al. Autoantibody systems in Rheumatoid Arthritis: specificity, sensitivity and diagnostic value. *Arthritis Res.* 4, 87-93, 2002.
2. Nienhuis R, et al. A new serum factor in patients with Rheumatoid Arthritis. The anti perinuclear factor. *Ann. Rheum. Dis.* 23, 302-305, 1964.
3. Schellekens G, et al. Citrulline is an essential constituent of antigenic determinants recognized by Rheumatoid Arthritis-specific autoantibodies. *J. Clin. Invest.* 101, 273-281, 1998.
4. Van Jaarsveld C, et al. The prognostic value of the antiperinuclear factor, anti-citrullinated peptide antibodies and rheumatoid factor in early Rheumatoid Arthritis. *Clin. Exp. Rheumatol.* 17, 689-697, 1999.
5. Schellekens, G, et al. The diagnostic properties of Rheumatoid Arthritis antibodies recognizing a cyclic citrullinated peptide. *Arthritis Rheum.* 43, 155-163, 2000.
6. Bizzaro N, et al. Diagnostic accuracy of the anti-citrulline antibody assay for Rheumatoid Arthritis. *Clinical Chemistry* 47, 1089-1093, 2001.
7. Visser H, et al. How to diagnose Rheumatoid Arthritis early? A predication model for persistent (erosive) arthritis. *Arthritis Rheum.* 46, 357-365, 2002.
8. Van Venrooij W, et al. Anti-citrullinated protein/peptide antibody and its role in the diagnosis and prognosis of early Rheumatoid Arthritis. *Neth. J. Med.* 60, 383-388, 2002.
9. Vossenaar, E, et al. Anti-CCP antibodies, a highly specific marker for (early) Rheumatoid Arthritis. *Clin. Applied Imm. Rev.* 4, 239-262, 2004.
10. Meyer O, et al. Anticitrullinated protein/peptide antibody assay in early Rheumatoid Arthritis for predicting five year radiographic damage. *Ann. Rheum. Dis* 62, 120-126, 2003.
11. Rantapää-Dahlqvist S, et al. Antibodies against citrullinated peptide and IgA rheumatoid factor predict the development of Rheumatoid Arthritis. *Arthritis Rheum.* 48, 2741-2749, 2003.
12. Forslind K, et al. Prediction of radiological outcome in early RA in clinical practice: role of antibodies to citrullinated peptides (Anti-CCP) *Ann. Rheum. Dis.* 63, 1090-1095, 2004.
13. Kastbom A, et al. Anti-CCP antibody test predicts the disease course during three years in early Rheumatoid Arthritis (the TIRA project). *Ann. Rheum. Dis.* 63, 1085-1089, 2004.
14. van Gaalen F, et al. Autoantibodies to cyclic citrullinated peptides predict progression to Rheumatoid Arthritis in patients with undifferentiated arthritis. *Arthritis Rheum.* 50, 709-715, 2004.

## **SUMMARY OF PROTOCOL**

1. Dilute samples 1:101. Do not dilute Calibrators, Reference, Positive and Negative Controls.
2. Add 100  $\mu$ L of Reference, Positive and Negative Controls/Calibrators in duplicate and pre-diluted samples into referenced wells.
3. Incubate  $60\pm 10$  minutes at 18-25° C.
4. Wash wells 3 times.
5. Add 100  $\mu$ L of Conjugate to each well.
6. Incubate  $30\pm 5$  minutes at 18-25° C.
7. Wash wells 3 times.
8. Add 100  $\mu$ L of Substrate to each well.
9. Incubate  $30\pm 5$  minutes at 18-25° C.
10. Add 100  $\mu$ L of Stop Solution to each well.
11. Read the absorbance at 550nm.

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## BIBLIOGRAPHIE

1. Van Boekel M, et al. Autoantibody systems in Rheumatoid Arthritis: specificity, sensitivity and diagnostic value. *Arthritis Res.* 4, 87-93, 2002.
2. Nienhuis R, et al. A new serum factor in patients with Rheumatoid Arthritis. The anti perinuclear factor. *Ann. Rheum. Dis.* 23, 302-305, 1964.
3. Schellekens G, et al. Citrulline is an essential constituent of antigenic determinants recognized by Rheumatoid Arthritis-specific autoantibodies. *J. Clin. Invest.* 101, 273-281, 1998.
4. Van Jaarsveld C, et al. The prognostic value of the antiperinuclear factor, anti-citrullinated peptide antibodies and rheumatoid factor in early Rheumatoid Arthritis. *Clin. Exp. Rheumatol.* 17, 689-697, 1999.
5. Schellekens, G, et al. The diagnostic properties of Rheumatoid Arthritis antibodies recognizing a cyclic citrullinated peptide. *Arthritis Rheum.* 43, 155-163, 2000.
6. Bizzaro N, et al. Diagnostic accuracy of the anti-citrulline antibody assay for Rheumatoid Arthritis. *Clinical Chemistry* 47, 1089-1093, 2001.
7. Visser H, et al. How to diagnose Rheumatoid Arthritis early? A predication model for persistent (erosive) arthritis. *Arthritis Rheum.* 46, 357-365, 2002.
8. Van Venrooij W, et al. Anti-citrullinated protein/peptide antibody and its role in the diagnosis and prognosis of early Rheumatoid Arthritis. *Neth. J. Med.* 60, 383-388, 2002.
9. Vossenaar, E, et al. Anti-CCP antibodies, a highly specific marker for (early) Rheumatoid Arthritis. *Clin. Applied Imm. Rev.* 4, 239-262, 2004.
10. Meyer O, et al. Anticitrullinated protein/peptide antibody assay in early Rheumatoid Arthritis for predicting five year radiographic damage. *Ann. Rheum. Dis* 62, 120-126, 2003.
11. Rantapää-Dahlqvist S, et al. Antibodies against citrullinated peptide and IgA rheumatoid factor predict the development of Rheumatoid Arthritis. *Arthritis Rheum.* 48, 2741-2749, 2003.
12. Forslind K, et al. Prediction of radiological outcome in early RA in clinical practice: role of antibodies to citrullinated peptides (Anti-CCP). *Ann. Rheum. Dis.* 63, 1090-1095, 2004.
13. Kastbom A, et al. Anti-CCP antibody test predicts the disease course during three years in early Rheumatoid Arthritis (the TIRA project). *Ann. Rheum. Dis.* 63, 1085-1089, 2004.
14. van Gaalen F, et al. Autoantibodies to cyclic citrullinated peptides predict progression to Rheumatoid Arthritis in patients with undifferentiated arthritis. *Arthritis Rheum.* 50, 709-715, 2004.

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**SVAR LIFE SCIENCE AB**