Instruction

DIASTAT[®] Anti-dsDNA

For professional use only
Usage reserve aux professionnels
Sólo para uso profesional
Nur für den fachgebrauch
Solo per uso professionale
Endast för professionell användning

Document No. LABEL-DOC-0024 2.0

DIASTAT[®] Anti-dsDNA

English: page 2
Français: page 14
Español: página 27
Deutsch: Seite 40
Italiano: pagina 52
Svenska: sida..... 64

INTENDED USE

The DIASTAT® anti-dsDNA test is a quantitative/qualitative enzyme-linked immunosorbent assay (ELISA) for the detection of IgG and IgM autoantibodies specific for double-stranded DNA (dsDNA) in human serum or plasma. It is intended to aid in the diagnosis of systemic lupus erythematosus, and to monitor dsDNA antibody levels during treatment and remission and is not definitive in isolation. Autoantibody levels represent one parameter in a multicriterion diagnostic process.

INTRODUCTION

Systemic rheumatic diseases are autoimmune disorders such as systemic lupus erythematosus (SLE), rheumatoid arthritis and mixed connective tissue disease. A general feature of systemic rheumatic diseases is the presence of circulating antibodies to a variety of cellular antigens. The detection and serological characterisation of specific autoantibodies plays an important role in the differential diagnosis of these diseases^{1,2}.

Antibodies against dsDNA are highly specific for SLE and are detected at a high frequency in untreated patients with active disease³. The presence of dsDNA antibodies in SLE is one of the criteria for disease classification by the American Rheumatism Association⁴. Several investigators have reported a correlation between disease activity and fluctuations in antibody levels^{5,6}, consequently the monitoring of dsDNA activity is considered of use in the management of SLE patients⁷.

A number of techniques are available for detecting dsDNA autoantibodies, including Farr radioimmunoassay and the *Crithidia luciliae* immunofluorescence assay. These tests have inherent drawbacks in complexity, lack of sensitivity and reproducibility. ELISAs offer simplicity, sensitivity, reproducibility and objectivity over other methods⁸ and can be standardised by reference to the World Health Organisation (WHO) anti-dsDNA standard, Wo/80⁹.

PRINCIPLE OF THE ASSAY

The wells of the microtitre strips are coated with a preparation of calf thymus selected for its high dsDNA content. During the first incubation, specific autoantibodies in diluted serum or plasma bind to the antigen-coated surface. The wells are then washed to remove unbound components. In the second incubation, the Conjugate, enzyme-labelled antibodies to human IgG and IgM, binds any surface-bound autoantibodies. After further washing, specific autoantibodies are traced by incubation with the Substrate. Addition of Stop Solution terminates the reaction, resulting in a coloured end-product. The amount of Conjugate bound is measured in absorbance units. In the qualitative protocol, the amount of Conjugate bound by the sample is compared with that bound by the Reference Control. In the quantitative protocol, the concentration of anti-dsDNA autoantibody can be estimated by interpolation from a dose-response curve based on Standards. The Standards are referenced against WHO Reference Preparation Wo/80 and reported in International units per mL (IU/mL).

KIT COMPONENTS

A	IgG/IgM Conjugate	1 × 15mL	Alkaline phosphatase-labelled antibodies to human IgG and IgM, Tris buffer, protein stabiliser, <0.1% (w/v) sodium azide. Ready-to-use.	
B	Substrate	1 × 15mL	Mg ²⁺ , phenolphthalein monophosphate (PMP), buffer solution. Ready-to-use. Do not expose to light during storage.	
C	Stop Solution	1 × 15mL	Sodium hydroxide, EDTA, carbonate buffer (pH >10). Ready-to-use.	
D	Wash Buffer Concentrate (16X)	2 × 25mL	Borate buffer, 0.4% (w/v) sodium azide. Dilute before use.	
E	dsDNA-Coated Wells and Strip Holder	12 × 8 well microtitre strips	Coated with dsDNA antigen, in a resealable foil pack with desiccant. Colour-coded LIME GREEN . Individual wells can be broken off from each microtitre strip.	
F	Sample Diluent 2 Concentrate (5X)	1 × 25mL	Phosphate buffer, protein stabiliser, 5% (w/v) Triton X-100, 0.5% (w/v) sodium azide. Dilute before use.	
1-5	Anti-dsDNA Standards	5 × 1.0mL	Human plasma, buffer, <0.1% (w/v) sodium azide. 0, 10, 50, 150, 300 IU/mL. Ready-to-use.	
6	Anti-dsDNA Reference Control	1 × 1.5mL	Human plasma, buffer, <0.1% (w/v) sodium azide. Ready-to-use.	
+/-	Positive Control Negative Control	1 × 0.2mL 1 × 0.1mL	Human plasma, <0.1% (w/v) sodium azide. Dilute 1:101 with diluted Sample Diluent 2 before use, as for samples.	
	Pack Leaflet			

STORAGE OF REAGENTS

Opened Kit Stability

A kit was opened and reused on three occasions over a three months period with no adverse effect on kit performance.

Handling and Procedural Notes

1. Store kit components at 2-8° C and use until the expiry date on the labels. Do not use expired reagents.
2. Do not mix different lot numbers.
3. Do not freeze kits.
4. Wash Buffer Concentrate, Sample Diluent 2 Concentrate and Positive and Negative Controls must be diluted before use. All other reagents are ready-to-use.
5. Diluted Wash Buffer and diluted Sample Diluent 2 are stable at 2-8° C for up to 6 months if microbial contamination is avoided.
6. Replace surplus microtitre strips in the foil pack and store with the desiccant at 2-8° C, until required.
7. The plate holder is adapted for use with snappable wells only.
8. Do not expose Substrate to light during storage.
9. Avoid contamination of reagents. Use a new disposable pipette tip for each reagent or sample manipulation.

Indications of Deterioration

The Substrate should be pale yellow in colour. Pink colouring indicates contamination and the reagent must be discarded. Turbidity or precipitation in any component indicates deterioration and the component should be discarded.

Sample Collection and Storage

The assay is recommended for serum/plasma samples; do not use lipaemic, haemolysed or turbid samples. Thoroughly mix thawed samples before assay and avoid repeated freeze/thawing. Do not heat-inactivate samples, this may yield false positive results.

Samples must be assayed within 12 hours of collection, or may be stored at or below -20° C.

WARNINGS AND PRECAUTIONS

For in vitro diagnostic use only.

Safety Precautions

1. Adhere strictly to the instructions in this booklet, particularly for handling and storage conditions.
2. Standards and Controls contain human plasma tested by FDA-cleared assays for hepatitis B surface antigen, HCV, HIV antigen and HIV antibodies and found to be non-reactive/negative. As no known test offers complete assurance that infectious agents are absent, Standards and Controls should be considered potentially infectious and handled with the same precautions as any other potentially biohazardous material. The CDC/NIH Health Manual "Biosafety in Microbiological and Biomedical Laboratories", 3rd edition, 1993, describes how these materials should be handled in accordance with Good Laboratory Practice. This is applicable in the USA.
3. Do not pipette by mouth.
4. Do not smoke, eat, drink or apply cosmetics in areas where kits and samples are handled.
5. Any skin complaints, cuts, abrasions and other skin lesions should be suitably protected.
6. The Standards, Controls, Conjugate, Sample Diluent 2 Concentrate and Wash Buffer Concentrate contain sodium azide which can react with lead and copper plumbing to form highly explosive metal azides. On disposal, drain with large quantities of water to prevent azide build-up.
7. The Stop Solution contains sodium hydroxide. Avoid contact with skin, eyes and mucous membranes. Spillage should be mopped up with copious amounts of water. If contact with skin or eyes occurs, irrigate with water and seek medical attention immediately.
8. The substrate contains PMP, Bronidox L, and Diethanolamine. Avoid contact with skin, eyes and respiratory system. If contact with skin, eyes or respiratory system occurs, rinse with water and seek medical advice.
9. Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

Warning

B.

SUBS

Contains: Diethanolamine

H319:	Causes serious eye irritation.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313:	If eye irritation persists: Get medical advice/attention.

**Warning****C.**

SOLN	STOP
------	------

Contains: Sodium hydroxide

H315:	Causes skin irritation.
H319:	Causes serious eye irritation.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352:	IF ON SKIN: Wash with plenty of soap and water.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P332+P313:	If skin irritation occurs: Get medical advice/attention.
P337+P313:	If eye irritation persists: Get medical advice/attention.

**Warning****D. and F.**

BUF	WASH	16X
DIL	SPE	5X

Contains: Sodium azide

H302:	Harmful if swallowed.
EUH032:	Contact with acids liberates very toxic gas.
H412:	Harmful to aquatic life with long lasting effects.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P301+P312:	IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
P273:	Avoid release to the environment.

PREPARATION

Materials/Equipment Required but not Provided

1. 96 well plate/strip reader with 550nm filter (540-565nm is acceptable).
2. Precision pipettes to dispense 10 μ L, 100 μ L, 1mL. Automatic pipette to dispense 100 μ L. Automatic pipette to dispense 200 μ L for manual washing, automatic plate washer optional.
3. Glass/plastic measuring cylinders: 1 \times 100mL, 1 \times 400mL.
4. 1mL volume containers.
5. Distilled/deionised water.
6. Paper towels.
7. Timer for 30 and 60 minutes intervals.

Preparation for the Assay

Allow all kit components, including the microtitre strips, to warm up to 18-25° C for 30-60 minutes before use. Mix reagents by gentle inversion.

Do not dilute the Reference Control.

Dilute the following reagents and mix thoroughly.

Reagent	Volume	Add
Wash Buffer Concentrate	1 vial	375mL distilled/deionised water
Sample Diluent 2 Concentrate	1 vial	100mL distilled/deionised water
Positive and Negative Controls/samples	10µL	1mL diluted Sample Diluent 2

Microtitre wells are supplied in strips of eight. If other than a multiple of eight wells are required, proceed as follows.

1. Remove strip from holder by pushing underside of wells.
2. Snap off required number of wells.
3. Hinge rectangular hole into bottom edge (to H) of the holder groove.
4. Ensure the square hole, with nick on left, is firmly held along the top edge (row A).

ASSAY PROTOCOL

Qualitative protocol: run Reference Control, Positive and Negative Controls, and samples.

Quantitative protocol: run Standards (1-5), Positive and Negative Controls, and samples.

1. Reference wells for identification.
2. Pipette 100µL Reference Control/Standards in duplicate, pre-diluted Positive and Negative Controls, and pre-diluted patient samples into appropriate wells. Remember to change pipette tips between additions. This step should not exceed 15 minutes for any one set of Standards/Controls/samples.
3. Incubate 60±10 minutes at 18-25° C.
4. Decant strip contents by quick inversion over a sink suitable for the disposal of biological materials, bearing in mind the potential infective hazard of the samples. Blot inverted strips well with paper towels.
5. Wash wells **three times** with a minimum of 200µL diluted Wash Buffer. **Decant and blot after each wash step.**
6. Add 100µL IgG/IgM Conjugate to each well.
7. Incubate 30±5 minutes at 18-25° C.
8. Repeat steps 4 and 5.
9. Add 100µL Substrate to each well.
10. Incubate 30±5 minutes at 18-25° C. **Do not decant.**
11. Add 100µL Stop Solution to each well, in the same order and rate as the Substrate. Tap wells gently to mix.
12. Read strips within 24 hours at 550nm (540-565nm).

CALCULATION AND INTERPRETATION OF RESULTS

Consider each assay separately when calculating and interpreting results.

Qualitative Protocol

Calculate the absorbance value (optical density) ratio for the Positive and Negative Controls, and for each sample.

$$\text{Absorbance Ratio} = \frac{\text{Sample or Control Absorbance Value}}{\text{mean Reference Control Absorbance Value}}$$

Users should calculate a cut-off between positive and negative samples that is specific to their patient populations. Results from the patient populations used in the Svar Life Science clinical trial suggest the following cut-off:

Absorbance Ratio

<0.95

≥0.95 to ≤1.0

>1.0

Result Interpretation

Negative

Borderline - recommend repeat testing

Positive

Quantitative Protocol

Plot the mean absorbance value of each Standard against log₁₀ Standard concentration (see following table) on suitable graph paper. Concentrations of Controls and samples can then be read from the standard curve; a typical plot is shown below for reference purposes, it must not be used for interpreting results. 4-parameter logistic (4PL), log/logit or spline curve fits are also satisfactory.


Samples with absorbances above Standard 5 (300 IU/mL) are outside the range of the assay, and should be reported as >300 IU/mL, diluted and re-assayed, correcting for this further dilution factor.

NB: As in any assay measuring antibodies, this assay determines the activity of the antibody present in the sample, rather than the concentration. Activity can be affected by a number of parameters, such as antibody avidity.

Standard Concentrations

Standard Number	Concentration IU/mL
1	0
2	10
3	50
4	150
5	300

Typical Standard Curve

QUALITY CONTROL

Ensure that adequate maintenance and calibration of the plate-reader is performed according to the manufacturer's instructions, and that the correct wavelength is employed.

Users should ensure they are fully acquainted with the instructions for the assay, particularly the Warnings and Precautions section, and the Handling and Procedural Notes. Users should demonstrate that they can obtain performance specifications for precision and reportable range of test results comparable to those established by the manufacturer before reporting patient test results. It is recommended that the pre-diluted Positive and Negative Controls are run in duplicate in all assays to monitor the quality of the test procedure. Run the ready-to-use Reference Control in duplicate in all qualitative assays.

Assuming the precision specifications described by the manufacturer are met, failure of any Control to meet the Control ratio specifications below renders the assay invalid and patient results should not be reported. The operator may repeat the assay, having reviewed their procedure, or contact the distributor/manufacturer. If repeating the assay, prepare a fresh dilution of each Control and sample. Laboratories may wish to include in-house controls in each assay run. Store such control material at or below -20° C and avoid repeat freeze/thaw cycles. Preservatives such as sodium azide at 0.1% (w/v) will not affect sample results.

Levels of analytes identified in particular diseases are those established by the manufacturer for specific populations, and may not necessarily mirror the literature. Incidence levels, their relationship to specific diseases, reference ranges, and appropriate cut-off points should all be calculated for the specific populations serviced by users.

Control Ratio Specifications

Protocol	Specifications
Qualitative (ratios)	$\frac{\text{Positive Control Absorbance}}{\text{Reference Control Absorbance}}$ see Positive Control label
	$\frac{\text{Negative Control Absorbance}}{\text{Reference Control Absorbance}} < 0.95$
Quantitative	See Positive Control label for acceptable expected range (IU/mL)
	Negative Control concentration <30 IU/mL

EXPECTED VALUES

The following is a guide to interpreting results. It is recommended that users establish reference ranges for the populations served by their own laboratories.

Anti-dsDNA Range IU/mL	Suggested Interpretation
<30	Negative for antibodies to dsDNA.
≥30 to 50	Borderline for antibodies to dsDNA - possible SLE or other connective tissue disease. Consider other testing for differential diagnosis, e.g. DIASTAT® anti-Sm assay. A repeat dsDNA test should be carried out on a subsequent sample.
>50 to 300	Positive for antibodies to dsDNA - probably SLE.
>300	Strong positive for antibodies to dsDNA - highly diagnostic for SLE.


Asymptomatic Samples

192 serum samples from asymptomatic apparently healthy donors were assayed for anti-dsDNA antibody. The negative cut-off was defined as the mean value plus two standard deviations and was calculated to be 30 IU/mL. 190/192 samples (99%) were negative for antibodies to dsDNA using this value.

The asymptomatic reference range was established as <30 IU/mL.

SLE Samples

257 serum samples were obtained from patients with clinically proven SLE. Of these, 142 (55%) were positive by the DIASTAT anti-dsDNA test. The overall frequency distribution of anti-dsDNA antibodies in these SLE patients is illustrated opposite.

**Non-SLE Samples**

A further 119 samples from non-SLE disease groups were assayed. The distribution of results is given in the following table.

Disease	Negative 0-30 IU/mL	Weak Positive 30-50 IU/mL	Positive 50-300 IU/mL
Rheumatoid Arthritis	32	2	1
Polymyositis	22	1	-
Sjögren's Syndrome	25	1	-
Viral Infections	27	1	-
Hypergammaglobulinaemia	6	1	-

PERFORMANCE DATA**Correlation**

338 serum samples (100 asymptomatic, 238 confirmed SLE patients) were assayed in parallel using the DIASTAT® anti-dsDNA kit with a cut-off value of 50 IU/mL and another ELISA with a value of 80 IU/mL. Weak positives as defined by both assays were excluded. The results are shown below.

		DIASTAT		
		+ve	-ve	Total
ELISA	+ve	93	12	105
	-ve	7	226	233
	Total	100	228	338

Overall agreement = 94%.

Good correlation was obtained with DIASTAT® anti-dsDNA and both the Farr radioimmunoassay and Crithidia luciliae immunofluorescence assay.

Dilution Characteristics

Four dilutions of four patient samples were assayed in three assays using three kit batches. The following table shows the mean values obtained and the dilution-corrected recovery.

Sample	Dilution	Mean Value IU/mL	Dilution Corrected % Recovery
1	A	171	100
	A/2	79.5	93
	A/4	39.8	93
	A/8	19.7	92
2	A	487	100
	A/2	211	87
	A/4	116	96
	A/8	57.6	95
3	A	146	100
	A/2	67.1	92
	A/4	33.3	91
	A/8	18.1	99
4	A	154	100
	A/2	89.1	116
	A/4	44.6	116
	A/8	26.5	137

Imprecision

- Intra-assay imprecision** determined by testing three controls, with replication of 12.

Control	Mean Value IU/mL	%CV
1	22.4	6.0
2	66.9	4.6
3	105	4.5

- Inter-assay imprecision** determined by testing three controls in 50 assays, using five operators and three kit batches.

Control	Mean Value IU/mL	%CV
1	22.2	8.7
2	65.3	9.1
3	105	8.1

LIMITATIONS OF USE



1. Although the presence of high titres of antibodies to dsDNA is indicative of SLE, the data must be considered in light of other clinical and laboratory findings.
2. Some individuals may have high levels of anti-dsDNA antibodies with little or no evidence of clinical disease. By contrast, some patients with SLE may have undetectable levels of these antibodies.
3. While IgG anti-dsDNA antibodies predominate in SLE, IgM anti-dsDNA antibodies can also be found in conjunction with IgG in some samples. Samples that contain high levels of IgM and low levels of IgG may not produce positive results in tests that only measure IgG antibodies¹⁰. This assay is configured to detect both IgG **and** IgM anti-dsDNA antibodies.
4. For repeat patient sampling, e.g. for monitoring, the same type of sample (serum or plasma) should be used throughout the study period.

REFERENCES

1. Tan EM. Autoantibodies to Nuclear Antigens (ANA): Their Immunobiology and Medicine. *Adv Immunol*, **33**, 167-240, 1982,
2. Morrow J, Isenberg D. *Autoimmune Rheumatic Diseases*, Blackwell Scientific Publications, Oxford 1987.
3. Tan EM. Antinuclear Antibodies: Diagnostic Markers for Autoimmune Diseases and Probes for Cell Biology. *Adv Immunol*, **44**, 93-151, 1989.
4. Tan EM, Cohen AS, et al. The 1982 Revised Criteria for the Classification of Systemic Lupus Erythematosus. *Arth Rheum*, **25**, 1271-1277, 1982.
5. Swaak AJG, Grownwold J, et al. Prognostic Value of anti-dsDNA in SLE. *Ann Rheum Dis*, **41** 388-395, 1982.
6. Smeenk R, Brinkman K, et al. Antibodies to DNA in Patients with Systemic Lupus Erythematosus: Their Role in the Diagnosis, the Follow-up and the Pathogenesis of the Disease. *Clin Rheum*, **9** (1), 100-110, 1990.
7. Borg EJ, Horst G, et al. Measurement of Increases in anti-dsDNA Antibody Levels as a Predictor of Disease Exacerbation in SLE. *Arth Rheum*, **33**, 634-643, 1990.
8. Halbert SP, Karsh J, Anken M. Studies on Autoantibodies to Deoxyribonucleic Acid and Deoxyribonucleoprotein with Enzyme-Immunoassay (ELISA). *J Lab Clin Med*, **97**, 97-111, 1981.
9. Feltkamp TEW, Kirkwood TBL, et al. The First International Standard for Antibodies to Double Stranded DNA. *Ann Rheum Dis*, **47**, 740-746, 1988.
10. Isenberg DA, Schoenfeld Y, Schwartz RS. Multiple Serologic Reactions and their Relationship to Clinical Activity in Systemic Lupus Erythematosus, Sjögren's Syndrome and Rheumatoid Arthritis. *Arthritis Rheum*, **27**, 132-138, 1984.

SUMMARY OF PROTOCOL

1. Dilute samples and Positive and Negative Controls 1:101. Do not dilute Standards or Reference Control.
2. Add 100µL of Reference Control/Standards in duplicate, pre-diluted Positive and Negative Controls and samples into referenced wells of the microtitre strip.
3. Incubate 60±10 minutes at 18-25° C.
4. Wash strips 3 times.
5. Add 100µL of IgG/IgM Conjugate to each well.
6. Incubate 30±5 minutes at 18-25° C.
7. Wash strips 3 times.
8. Add 100µL of Substrate to each well.
9. Incubate 30±5 minutes at 18-25° C.
10. Add 100µL of Stop Solution to each well.
11. Read absorbance at 550nm.

	Batch code / Code du lot / Código de lote / Chargenbezeichnung / Codice del lotto / Satsnummer
	Catalogue number/ Référence du catalogue / Número de catálogo / Bestellnummer / Numero di catalogo / Katalognummer
	Use by/ Utiliser jusque / Fecha de caducidad / Verwendbar bis / Utilizzare entro / Hållbar till
	Temperature limitation/ Limites de temperature / Limite de temperatura / Zulässiger Temperatur-bereich / Limite di temperatura / Temperaturgränser
	Biological risks/ Risques biologiques / Riesgo biológico/Biogefährdung /Rischio biologico / Biologisk risk
	Consult instructions for use/ Consulter les instructions d'utilisation / Consulte las instrucciones de uso/ Gebrauchsanweisung beachten / Consultare le istruzioni per l'uso / Se bruksanvisning
	In Vitro Diagnostic Medical Device/ Dispositif medical de diagnostic in vitro / Producto sanitario para diagnóstico in vitro / In Vitro Diagnostikum / Dispositivo medico-diagnostico in vitro / In vitro diagnostika
	Manufacturer/ Fabricant / Fabricante / Hersteller / Fabbicante / Tillverkare
	Contains sufficient for <n> tests/ Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Ausreichend für "n" Ansätze / Conteúdo suficiente para <n> ensaios / Innehållet tillräckligt för n tester
	Warning / Attention / Atención / Achtung / Attenzione / Varning
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive / Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro / Konform mit Richtlinie 98/79/EG zu In-vitro- Diagnostika / Conformità alla direttiva 98/79/CE relativa ai dispositivi medico- diagnostici in vitro / Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.
	Conjugate/ Conjugué / Conjugado / Konjugat / Conjugato / Konjugat
	Substrate/ Substrat / Substrato / Substrat / Substrato / Substrat
	Stop solution/ Solution d'Arrêt / Solución de Parada / Stopplösung/ Soluzione bloccante / Stopplösning

BUF	WASH	16 x	Wash buffer concentrate (16 X)/ Concentré tampon de (16X lavage)/ Concentrado de Búfer de lavado (16 X) / Washpuffer-Konzentrat (16 X) / Tampone di lavaggio concentrato (16 X) / Tvättbuffert koncentrat (16 X)
Ag			dsDNA-coated wells and strip holder / Cupules enduites de ADNdb et Portebandes / Soporte para Bandas y Vasos Recubiertos con dsDNA / dsDNA-beschichtete Vertiefungen und Streifenrahmen / Pozzetti rivestiti di dsDNA e supporto per strip / dsDNA-klädda brunnar och striphållare
DIL	SPE	5 X	Sample Diluent Concentrate (5 X) / Concentré diluent pour échantillons (5 X)/ Concentrado de Diluente de Muestra (5 X) / Probediluens Konzentrat / Diluente per campioni concentrato (5 X) / Provspädningsbuffert koncentrat (5 X)
CAL	X		Anti-dsDNA Standards 1-5/ Etalons anti-ADNdb 1-5 / Estandares Anti-dsDNA 1-5 / Anti-dsDNA Standards 1-5 / Standard anti-dsDNA 1-5 / Anti-dsDNA-standarder
CONTROL	REF		Anti-dsDNA Reference Control/ Témoïn de référence anti-ADNdb / Control de Referencia Anti-dsDNA / Anti-dsDNA Referenzkontrolle / Controllo di riferimento anti-dsDNA / Anti-dsDNA referenskontroll
CONTROL	+		Positive Control/ Témoïn positive / Control Positivo / Positiv-Kontrolle / Controllo Positivo / Positiva kontroll
CONTROL	-		Negative Control/ Témoïn negatif / Control Negativo / Negativ-Kontrolle / Controllo negativo / Negativ kontroll

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REFERENSER

1. Tan, E.M. 1982. Autoantibodies to nuclear antigens (ANA): Their immunobiology and medicine. *Adv. Immunol.* 33: 167-240.
2. Morrow, J. and Isenberg, D. 1987. In "Autoimmune Rheumatic Diseases," Blackwell.
3. Tan, E.M. 1989. Antinuclear antibodies: Diagnostic markers for autoimmune diseases and probes for cell biology. *Ad. Immunol.* 44: 93-151.
4. Tan, E.M., *et al.* 1982. The 1982 revised criteria for the classification of systemic lupus erythematosus. *Arth. Rheum* 25:1271-1277.
5. Swaak, A.J.G., *et al.* 1982. Prognostic value of anti-dsDNA in SLE. *Ann. Rheum. Dis.* 41: 388-395.
6. Smeenk, R., *et al.* 1990. Antibodies to DNA in patients with systemic lupus erythematosus. Their role in the diagnosis, the followup and the pathogenesis of the disease. *Clin. Rheum.* 9: 100-110.
7. Borg, E.J., *et al.* 1990. *Measurement of increases in anti-dsDNA antibody levels as a predictor of disease exacerbation in SLE.* *Arth. Rheum.* 33: 634-643.
8. Halbert, S.P., Karsh, J., Anken, M. 1981. Studies on autoantibodies to deoxyribonucleic acid and deoxyribonucleoprotein with enzymeimmunoassay (ELISA). *J. Lab. Clin. Med.* 97: 97-111.
9. Feltkamp, T.E.W., Kirkwood, T.B.L., Maini, R.N., Aarden, L.A. 1988. The first international standard for antibodies to double stranded DNA. *Ann. Rheum. Dis.* 47: 740-746.
10. Data on file and available on request.
11. Isenberg, D.A., Schoenfeld, Y., Schwartz, R.S. 1984. Multiple serologic reactions and their relationship to clinical activity in systemic lupus erythematosus, Sjögren's Syndrome and rheumatoid arthritis. *Arthritis Rheum.* 27: 132-138.
12. HHS Publication No. 93-8395, 3rd ed. May 1993. *Biosafety in Microbiological and Biomedical Laboratories.* Washington, DC: US Department of Health and Human Services.
13. *Manual Guide – Safety Management No. CDC-22 Decontamination of Laboratory Sink Drains to Remove Azide Salts.* April 30, 1976. Atlanta, GA: Centers for Disease Control.

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SVAR LIFE SCIENCE AB