Instruction

WIESLAB® Vasculitis screen

Enzyme immunoassay for detection of
autoantibodies against GBM, Proteinase 3 (PR3-ANCA)
and Myeloperoxidase (MPO-ANCA)

Microtitration 96 wells
Store the kit at +2-8° C
For in vitro diagnostic use only

Document No. LABEL-DOC-0049 2.0

WIESLAB® Vasculitis screen

English:	page	2
Français:	page	14
Español:	página	21
Deutsch:	Seite	28
Italiano:	pagina	35
Português:	página	42
Dansk:	side.....	50
Norsk:	side.....	57
Svenska:	sida.....	64

INTENDED USE

The Wieslab® Vasculitis Screening Test Kit is an enzyme-linked immunosorbent assay (ELISA) for the qualitative detection of antibodies to glomerular basement membrane (GBM), Proteinase-3 (PR3) and Myeloperoxidase (MPO) in human sera. The assay is used to detect antibodies in a single serum specimen. The results of the assay are to be used as an aid to the diagnosis of reno-pulmonary syndromes and rapidly progressive glomerulonephritis, especially Goodpasture syndrome (GP), Wegener's granulomatosis (WG) and microscopic polyangiitis (MP). The assay is intended for use in patients with signs and symptoms consistent with GP, WG, and MP. It is not intended for screening a healthy population.

FOR IN VITRO DIAGNOSTIC USE.

A positive result should always be confirmed by a quantitative assay.

Summary and explanation

Goodpasture syndrome is characterised by lung haemorrhage, renal failure and the presence of anti-GBM antibodies. The diagnosis is based on clinical signs of lung haemorrhage and rapidly progressive glomerulonephritis, but there are several other autoimmune diseases which may present with similar symptoms and less than one third of patients with reno-pulmonary syndromes have antibodies against the Goodpasture antigen, the majority having either proteinase 3 (PR3)-ANCA or myeloperoxidase (MPO)-ANCA (1, 2).

Thus patients with reno-pulmonary syndromes can be identified by three autoantibodies: i.e. anti-GBM, PR3-ANCA or MPO-ANCA (2).

The molecular nature of the Goodpasture antigen is described in (3) and is the NC1 domain of the alpha 3 chain of type IV collagen. The antigen is characterised by a restricted tissue distribution, occurring mainly in kidney and lungs.

ANCAs (anti-neutrophil cytoplasmic antibodies) are a family of autoantibodies related to vasculitis and inflammatory disorders. The first method to detect ANCA was indirect immunofluorescence (IIF) performed on ethanol fixed granulocytes (4). This method yields two patterns, a cytoplasmic staining of the granulocyte denoting the presence of c-ANCAs, and a perinuclear staining denoting the presence of p-ANCAs. IIF was followed by ELISAs using the purified proteins (5). The most important antigens are proteinase 3 (PR3) and myeloperoxidase (MPO). PR3 is a serine protease with a molecular weight of 29kD, and MPO is a dimer with a molecular weight of 140kD. Thus antibodies to proteinase 3 are termed PR3-ANCA, and antibodies to myeloperoxidase are termed MPO-ANCA.

Antibodies to PR3 and MPO can be detected either by a direct ELISA method or capture technique where monoclonal antibodies capture the antigens. The capture technique is therefore an alternative method to detect ANCA. It has been described that the capture technique for PR3-ANCA has a higher sensitivity and correlates to severity of disease (6, 7). It is uncommon with relapse of disease if the capture ELISA is negative (8, 9, 10, 11).

Approximately 80-90% of WG patients manifest PR3-ANCA and 5-15% MPO-ANCA. In microscopic polyangiitis (MP) most patients with active MP are characterised by positive ANCA test results, MPO-ANCA being more frequent than PR3-ANCA.

Principle of the Wieslab® Vasculitis screen

The wells of a microtiterplate are coated with purified GBM antigen (Bovine) (row A and B), anti-PR3 monoclonal antibody and purified proteinase 3 (Human Neutrophil) (row C and D), anti-MPO monoclonal antibody and purified myeloperoxidase (Human Neutrophil) (row E and F), monoclonal mouse antibody (row G) and HSA (Human Serum Albumin) (row H). G and H are control wells.

During the first incubation, specific antibodies in diluted serum will bind to the antigen coating. The wells are then washed to remove unbound antibodies and other components. A conjugate of alkaline phosphatase labelled (goat) antibodies to human IgG binds to the antibodies in the wells in the second incubation. After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. Addition of stop solution terminates the reaction, resulting in a coloured endproduct. The amount of bound antibodies correlates to the colour intensity and is measured as absorbance (optical density (OD)).

Warnings and precautions

- For in vitro diagnostic use.
- Optimal results will be obtained by strict adherence to this protocol.
- Do not use components past the expiration date and do not mix components from different lots.
- The human serum components used in the preparation of the controls in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The stop solution contains 10 mM EDTA.
- The concentrations of anti GBM and ANCA in a given specimen determined with assays from different manufactures can vary due to differences in assay methods and reagent specificity.
- Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

BUF	WASH	30X
DIL		
CONJ		

CONTROL	+
CONTROL	-
SOLN	STOP
SUBS	pNPP

Warning

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

H317:	May cause an allergic skin reaction.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P302+352:	IF ON SKIN: Wash with plenty of soap and water.
P333+313:	If skin irritation or rash occurs: Get medical advice/attention.

Specimen collection

The Wieslab® Vasculitis screen assay is for use with serum. Handle as if capable of transmitting infections agents.

Avoid using sera which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8°C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20°C or colder.

Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

Kit components and storage of reagents

- One frame with 96 wells coated with purified GBM antigen (Bovine), monoclonal anti-

proteinase 3/proteinase 3 (Human Neutrophil) and monoclonal anti-myeloperoxidase/myeloperoxidase (Human Neutrophil) and control wells with monoclonal mouse antibody or human serum albumin: one lid sealed in a foil pack with a dry pack.

- 3 mL negative control (NC) containing human serum in diluent (green colour).
- 3 mL positive control (PC) containing human serum in diluent (red colour).
- 13 mL conjugate containing alkaline phosphatase-labelled (Goat) antibodies to human IgG (blue colour).
- 32 mL Diluent (Dil) containing PBS (red colour).
- 13 mL Substrate pNPP.
- 13 mL Stop Solution.
- 30 mL Wash solution 30x concentrated.

All reagents in the kit are ready for use except wash solution and should be stored at + 2-8° C.

Materials or equipment required but not provided

Microplate reader with filter 405 nm. Distilled water for preparing more washing solution. Washer for strips, absorbent tissue, plastic tubes to dilute the patient sample and a timer. Precision pipettes with disposable tips.

Do not use reference wavelength (620 nm or similar) when reading the plate. Do not use reference wavelength subtracted OD values for the calculations.

PROCEDURE

All solutions should be used at room temperature. Do not open the foil packaged plate until it has been room tempered as damp from condensation may have a negative affect on the antigen. Incubate all steps at room temperature (20-28°C) with lid. Remove only the number of wells needed for testing, resealing the aluminium package carefully. In the event that the wells become separated they can be reassembled in the correct order. Each well in a column is individually marked with different numbers of notches, one notch in well A, two notches in well B and so on.

Preparation of washing solution

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37°C water bath until the crystals have dissolved before dilution of wash solution. Dilute 10 mL of the 30x concentrated wash solution in 290 mL distilled water. When stored at 2-8°C, the diluted wash solution is stable until the date of expiration of the kit.

Dilution of serum and incubation

Dilute the patient's serum 1/20 with diluent (950 µL diluent + 50 µL serum).

Pipet 100 µL/well of negative control (NC), positive control (PC) according to the diagram below. Diluted patient serum (P) is added to each well in the Sample Column. Incubate for 10 minutes.

	Reference Column	Sample Column	Antigen coat per row
A	NC	P	GBM
B	PC	P	GBM
C	NC	P	PR3
D	PC	P	PR3
E	NC	P	MPO
F	PC	P	MPO
G	NC	P	Mouse monoclonal
H		P	HSA

After serum incubation

Wash 4 times with 300 µL washing solution/well, filling and emptying the wells each time; after the last wash, empty the wells by tapping the wells on an absorbent tissue.

Adding conjugate

Add 100 µL conjugate to each well. Incubate for 10 minutes.

After conjugate incubation

Wash as before.

Adding substrate

Add 100 µL substrate pNPP to each well, incubate for 20 minutes.

Adding stop solution

Add 100 µL stop solution to each well and read the absorbance at 405 nm on a microplate reader within 2 hours.

Do not use reference wavelength (620 nm or similar) when reading the plate.

Calculations

Do not use reference wavelength subtracted OD values for the calculations.

An optical density (OD) ratio for each patient sample is calculated as follows:

$$\text{OD ratio} = \frac{\text{OD patient sample}}{\text{OD of negative control}}$$

The patient sample is negative if the OD ratio is < 3.0 and positive if the OD ratio is ≥ 3.0.

Quality Control

All calculations should only be made on the 405 nm data without reference wavelength OD subtraction.

The OD for the NC for each antigen should be < 0.4.

See Certificate for acceptable expected range (Ratio) for the PC.

The range is different for each antigen.

If the PC or NC values are not within their respective ranges, the test should be considered invalid and the test should be repeated.

Control wells: The control well with mouse monoclonal antibody (row G) should have a negative OD ratio i.e. < 3.0. The OD value for the patient sample in the HSA well (row H) should be < 0.4. If the OD value for the patient sample is ≥0.4 in the HSA or positive on more than one antigen the sample may contain unspecific reactivity. In some samples unspecific reactivity will disappear with a higher dilution as used in a quantitative test.

Dual PR3/MPO positives: Occasionally patient samples are positive for both PR3 and MPO. This is extremely rare. Careful evaluation of such results are then necessary. A simultaneous rise in OD compared to NC or even positive ratio result in the mab-control well (row G) strongly indicates a false positive reaction due to the presence of antibodies to the mouse monoclonals used to capture the antigens. A further indication of that false reaction is that the sample is negative in the GP-coated wells. In these situations the test results can not be used. It is recommended to retest the sample by the use of Wieslab Capture PR3/MPO-ANCA quantitative methods where a higher sample dilution is used.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off.

Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to NCCLS C24-A for guidance on appropriate QC practices.

Interpretation of results

A sample with an OD ratio of:

< 3.0 = Negative

≥ 3.0 = Positive

Owing to the high sensitivity of the screening kit, in a few percent of cases samples with an OD ratio ≥ 3.0 may be found to be negative in a quantitative assay.

A positive test result should always be confirmed by a quantitative assay.

Limitations

The individual patient's OD ratio can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardization of results. The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used in combination with clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive, with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies.

Immunosuppressive therapy should not be started on basis of a positive anti-GBM, ANCA result. Initiation or changes in treatment should not be based on changes in anti-GBM, ANCA concentration alone, but rather on careful clinical observation.

Expected results

ANCA and anti-GBM are rarely found in normal healthy individuals. The anti-GBM, ANCA screening kit was tested with 120 normal sera. 120 were found to be negative. The annual incidence of Goodpasture syndrome is one per million and one new patient with PR3-ANCA is expected per 100,000 inhabitants per year. Around 10% of patients with WG are negative in both IIF and ELISA.

Anti-GBM was tested with 40 GP sera, 40 were found to be positive.

The MPO-ANCA was tested with 40 WG sera, 3 were found to be positive. The MPO-ANCA was tested with 40 MP sera, 28 were found to be positive.

The PR3-ANCA was tested with 40 WG sera, 38 were found to be positive. The PR3-ANCA was tested with 40 MP sera, 9 were found to be positive. (Table 1)

Performance characteristics

Table 1. Clinical sensitivity and specificity. A total of 288 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total number	Negative <3			Positive ≥ 3.0		
		Anti-GBM	Capture-PR3	Capture-MPO	Anti-GBM	Capture-PR3	Capture-MPO
GP	40	0	36	29	40	4	11
WG	40	40	2	37	0	38	3
MPA	40	40	31	12	0	9	28
Blood Donors	120	120	120	120	0	0	0
RA	20	20	20	20	0	0	0
SLE	8	8	8	8	0	0	0
UICo	20	20	20	20	0	0	0
Totalt	288	248	237	246	40	51	42

GP = Goodpasture syndrome

WG = Wegener's granulomatosis

MPA = microscopic polyangiitis

RA (anti-CCP positive samples) = rheumatoid arthritis

SLE (anti-Sm positive samples) = systemic lupus erythematosus

UICo = ulcerative colitis

GBM = anti-glomerular basement membrane

Clinical sensitivity

Anti-GBM

GP 40/40 = 100 % 95% CI = 91.2 – 100%

Capture-PR3 ANCA

WG 38/40 = 95.0 % 95% CI = 83.1 – 99.4%

MP 9/40 = 22.5 % 95% CI = 10.8 – 38.4%

Capture-MPO ANCA

WG 3/40 = 7.5 % 95% CI = 1.8 – 20.4%

MP 28/40 = 70.0 % 95% CI = 53.5 – 83.4 %

Clinical specificity

Anti-GBM

SLE 8/8 = 100 % 95% CI = 63.1 – 100%

RA 20/20 = 100 % 95% CI = 83.2 – 100%

UICo 20/20 = 100 % 95% CI = 83.2 – 100%

Blood Donors 120/120 = 100 % 95% CI = 97.0 – 100%

Capture-PR3 ANCA

SLE 8/8 = 100 % 95% CI = 63.1 – 100%

RA 20/20 = 100 % 95% CI = 83.2 – 100%

UICo 20/20 = 100 % 95% CI = 83.2 – 100%

Blood Donors 120/120 = 100 % 95% CI = 97.0 – 100%

Capture-MPO ANCA

SLE 8/8 = 100 % 95% CI = 63.1 – 100%

RA 20/20 = 100 % 95% CI = 83.2 – 100%

UICo 20/20 = 100 % 95% CI = 83.2 – 100%

Blood Donors 120/120 = 100 % 95% CI = 97.0 – 100%

Table 2. Relative sensitivity and specificity of the Wieslab® Vasculitis screen kit compared to an alternative qualitative ELISA. A total of 288 frozen retrospective sera were assayed on the Wieslab® Vasculitis screen kit and Wieslab® anti-GBM, ANCA screen.

Anti-GBM		Wieslab® anti-GBM, ANCA screen	
		Negative	Positive
Wieslab® Vasculitis screen kit	Negative	248	0
	Positive	0	40

Total 288 serum samples.

Capture-PR3		Wieslab® anti-GBM, ANCA screen	
		Negative	Positive
Wieslab® Vasculitis screen kit	Negative	232	4
	Positive	3	46

Total 285 serum samples; 3 borderline samples in reference kit were excluded from calculations.

Capture-MPO		Wieslab® anti-GBM, ANCA screen	
Wieslab® Vasculitis screen kit		Negative	Positive
	Negative	238	5
	Positive	0	42

Total 285 serum samples; 3 borderline samples in reference kit were excluded from calculations.

Positive Percent Agreement

Anti-GBM	= 40/40 = 100 %	95% CI = 91.2 – 100%
Capture-PR3 ANCA	= 46/50 = 92.0%	95% CI = 80.8 – 97.8%
Capture-MPO ANCA	= 42/47 = 89.4%	95% CI = 76.9 – 96.5%

Negative Percent Agreement

Anti-GBM	= 248/248 = 100 %	95% CI = 98.5 – 100%
Capture-PR3 ANCA	= 232/235 = 98.7 %	95% CI = 96.3 – 99.7%
Capture-MPO ANCA	= 238/238 = 100 %	95% CI = 98.5 – 100%

Overall Percent Agreement

Anti-GBM	= 288/288 = 100%	95% CI = 98.7 – 100%
Capture-PR3 ANCA	= 278/285 = 97.5%	95% CI = 95 – 99%
Capture-MPO ANCA	= 280/285 = 98.2%	95% CI = 96 – 99.4%

The 95% confidence interval (CI) was calculated using the exact method.

Table 3. Batch to batch variation was determined by testing six different samples in eight replicates on three different batches.

	1	2	3	4	5	6
Anti-GBM						
Mean OD ratio	3.9	5.9	4.4	6.5	8.3	6.5
SD	0.50	0.44	0.51	0.89	0.58	0.75
% CV	13	7	12	14	7	12
Capture-PR3						
Mean OD ratio	2.9	3.0	3.3	4.4	6.8	6.5
SD	0.71	0.51	0.43	0.50	0.98	1.48
% CV	24	17	13	11	14	23
Capture-MPO						
Mean OD ratio	3.1	3.2	3.1	6.8	5.6	6.5
SD	0.33	0.25	0.19	0.77	0.38	0.48
% CV	11	8	6	11	7	7

Table 4. Inter-assay precision was determined by testing six different samples in eight replicates at three separate occasions.

	1	2	3	4	5	6
Anti-GBM						
Mean OD ratio	3.4	5.8	4.3	6.3	8.9	6.3
SD	0.31	0.20	0.15	0.21	0.22	0.25
% CV	9	3	4	3	2	4
Capture- PR3						
Mean OD ratio	4.0	3.3	3.4	4.4	6.9	6.9
SD	0.76	0.25	0.21	0.52	0.30	0.54
% CV	19	8	6	12	4	8
Capture-MPO						
Mean OD ratio	3.5	3.6	3.1	7.3	5.6	6.9
SD	0.19	0.41	0.20	0.36	0.55	0.39
% CV	5	11	6	5	10	6

Table 5. Intra-assay precision was determined by testing six different samples in eight replicates at one occasion.

	1	2	3	4	5	6
Anti-GBM						
Mean OD ratio	4.5	6.9	4.8	7.7	4.8	7.7
SD	0.21	0.15	0.28	0.37	0.28	0.38
% CV	5	2	6	5	6	5
Capture-PR3						
Mean OD ratio	5.1	4.0	4.7	7.5	4.7	7.5
SD	0.27	0.11	0.13	0.34	0.13	0.33
% CV	5	3	3	5	3	4
Capture-MPO						
Mean OD ratio	4.2	4.2	4.5	9.3	4.5	9.3
SD	0.10	0.18	0.15	0.21	0.15	0.21
% CV	2	4	3	2	3	2

Troubleshooting

Problem:	Possible causes:	Solution:
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagents. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water used to prepare solution. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto-dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

References

1. Saxena R, Bygren P, Arvastson B, Wieslander J. Circulating autoantibodies as serological markers in the differential diagnosis of pulmonary renal syndrome. *J Intern Med.* 1995 Aug;238(2):143-52.
2. Westman KW, Bygren PG, Eilert I, Wiik A, Wieslander J. Rapid screening assay for anti-GBM antibody and ANCA; an important tool for the differential diagnosis of pulmonary renal syndromes. *Nephrol Dial Transplant.* 1997 Sep;12(9):1863-8.
3. Gunnarsson A, Hellmark T, Wieslander J. Molecular properties of the Goodpasture epitope. *J Biol Chem.* 2000 Oct 6;275(40):30844-8.
4. van der Woude FJ, Rasmussen N, Lobatto S, Wiik A, Permin H, van Es LA, van der Giessen M, van der Hem GK, The TH. Autoantibodies against neutrophils and monocytes: tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet.* 1985 Feb 23;1(8426):425-9.
5. Segelmark M, Westman K, Wieslander J. How and why should we detect ANCA? *Clin Exp Rheumatol.* 2000 Sep-Oct;18(5):629-35.
6. Westman KW, Selga D, Bygren P, Segelmark M, Baslund B, Wiik A, Wieslander J. Clinical evaluation of a capture ELISA for detection of proteinase-3 antineutrophil cytoplasmic antibody. *Kidney Int.* 1998 May;53(5):1230-6.
7. Csernok E, Holle J, Hellmich B, Willem J, Tervaert C, Kallenberg CG, Limburg PC, Niles J, Pan G, Specks U, Westman K, Wieslander J, De Groot K, Gross WL. Evaluation of capture ELISA for detection of antineutrophil cytoplasmic antibodies directed against proteinase 3 in Wegener's granulomatosis: first results from a multicentre study. *Rheumatology (Oxford).* 2004 Feb;43(2):174-80. Epub 2003 Oct 29.
8. Westman KW, Bygren PG, Olsson H, Ranstam J, Wieslander J. Relapse rate, renal survival, and cancer morbidity in patients with Wegener's granulomatosis or microscopic polyangiitis with renal involvement. *J Am Soc Nephrol.* 1998 May;9(5):842-52.
9. Gisslén K, Wieslander J, Westberg G, Herlitz H. Relationship between anti-neutrophil cytoplasmic antibody determined with conventional binding and the capture assay, and long-term clinical course in vasculitis. *J Intern Med.* 2002 Feb;251(2):129-35.
10. Segelmark M, Phillips BD, Hogan SL, Falk RJ, Jennette JC. Monitoring proteinase 3 antineutrophil cytoplasmic antibodies for detection of relapses in small vessel vasculitis. *Clin Diagn Lab Immunol.* 2003 Sep;10(5):769-74.
11. Westman KW, Selga D, Isberg PE, Bladström A, Olsson H. High proteinase 3-anti-neutrophil cytoplasmic antibody (ANCA) level measured by the capture enzyme-linked immunosorbent assay method is associated with decreased patient survival in ANCA-associated vasculitis with renal involvement. *J Am Soc Nephrol.* 2003 Nov;14(11):2926-33.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Explicação dos símbolos Forklaring til symboler. Symboler som brukes på etiketter. Förklaringar till symboler.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Inneholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medisinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigénio. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Fortynning. Diluent. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Solução de lavagem concentrada 30x. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentrert 30 ganger. Tvättbuffert 30x konc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
SOLN STOP	Stop solution. Solution d'Arrêt. Solución de Parada. Stopplösung. Soluzione bloccante. Solução de paragem. Stoppløsning. Stoppløsning. Stoppløsning
CONTROL X	Control. Contrôle. Controllo. Kontrolle. Controllo. Controllo. Kontroll. Kontroll. Kontroll.

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SVAR LIFE SCIENCE AB