Instruction

WIESLAB® Anti-GBM, ANCA screen kit

Enzyme immunoassay for detection of
autoantibodies against GBM, Proteinase 3 (PR3-ANCA)
and Myeloperoxidase (MPO-ANCA)

Double microtitration strips (6x16) 96 wells
Store the kit at +2-8° C
For in vitro diagnostic use only

Document No. LABEL-DOC-0027, 2.0
December, 2018


WIESLAB® Anti-GBM, ANCA screen kit

English: page	2
Français: page	13
Español: página	25
Deutsch: Seite	32
Italiano: pagina	39
Dansk: side.....	46
Norsk: side.....	53
Svenska: sida.....	60

REF

GCP 100

IVD

INTENDED USE

The Wieslab® Anti-GBM, ANCA Screening Test Kit is an enzyme-linked immunosorbent assay (ELISA) for the qualitative detection of antibodies to glomerular basement membrane (GBM), Proteinase-3 (PR3) and Myeloperoxidase (MPO) in human sera. The assay is used to detect antibodies in a single serum specimen. The results of the assay are to be used as an aid to the diagnosis of reno-pulmonary syndromes and rapidly progressive glomerulonephritis, especially Goodpasture syndrome (GP), Wegener's granulomatosis (WG) and microscopic polyangiitis (MP). The assay is intended for use in patients with signs and symptoms consistent with GP, WG, and MP. It is not intended for screening a healthy population.

FOR IN VITRO DIAGNOSTIC USE.

A positive result should always be confirmed by a semi-quantitative assay.

Summary and explanation

Goodpasture syndrome is characterised by lung haemorrhage, renal failure and the presence of anti-GBM antibodies. The diagnosis is based on clinical signs of lung haemorrhage and rapidly progressive glomerulonephritis combined with findings of anti-GBM antibodies. As there are several other autoimmune diseases which may present with similar symptoms, this kit is a useful aid in differentiation. Less than one third of patients with reno-pulmonary syndromes have antibodies against the Goodpasture antigen, the majority having either proteinase 3-ANCA or myeloperoxidase-ANCA (1). Indirect immunofluorescence was formerly used to detect anti-GBM antibodies. When the first ELISA based on a collagenase digest was published in 1981, assays using crude extracts were the only alternative. In 1984, the specific antigen was shown to derive from the C-terminal domain of type IV collagen (α 3 chain) (2, 3) and sensitive and specific assays were subsequently developed (4). The molecular nature of the Goodpasture antigen is reviewed in (5). The antigen is characterised by a restricted tissue distribution, occurring mainly in kidney and lungs.

ANCAs (anti-neutrophil cytoplasmic antibodies) are a family of autoantibodies related to vasculitis and inflammatory disorders. Since 1985, when c-ANCA was shown to be related to Wegener's granulomatosis (6), interest in ANCAs has steadily increased, and today these antibodies are considered to be major diagnostic tools for the investigation of systemic vasculitis.

The first method to detect ANCA was indirect immunofluorescence (IIF) performed on ethanol fixed granulocytes (7). This method yields two patterns, a cytoplasmic staining of the granulocyte denoting the presence of c-ANCAs, and a perinuclear staining denoting the presence of p-ANCAs. IIF was followed by ELISAs using the purified proteins (8).

The granulocyte is full of granules each with many different proteins. It was early shown that antibodies from systemic vasculitis patients bind to the alpha fraction containing the azurophil granules. The most important proteins were proteinase 3 (PR3)(9,10) and myeloperoxidase (MPO)(11,12). PR3 is a serine protease with a molecular weight of 29kD, and MPO is a dimer with a molecular weight of 140kD. Thus antibodies to proteinase 3 are termed PR3-ANCA, and antibodies to myeloperoxidase are termed MPO-ANCA.

Approximately 80-90% of WG patients manifest PR3-ANCA and 5-15% MPO-ANCA. One category of vasculitis is microscopic polyangiitis (MP). Most patients with active MP are characterised by positive ANCA test results, MPO-ANCA being more frequent than PR3-ANCA (13).

The Wieslab® anti-GBM, ANCA screening kit has been clinically evaluated (14).

Principle of the Wieslab® ANCA screen kit

The wells of a double microtiterstrip are coated with purified (Bovine) GP antigen (A1, A2, B1 and B2), purified (Human Neutrophils) proteinase 3 (C1, C2, D1 and D2), purified (Human Neutrophils) myeloperoxidase (E1, E2, F1 and F2), HSA (Human Serum Albumin) (G1, G2) and empty (H1, H2). G and H are blank wells. During the first incubation, specific antibodies in diluted serum will bind to the antigen coating. The wells are then washed to remove unbound antibodies and other components. A conjugate of alkaline phosphatase labelled (goat) antibodies to human IgG binds to the antibodies in the wells in the second incubation. After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour intensity and is measured as absorbance (optical density (OD)).

Warnings and precautions

- For in vitro diagnostic use.
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of anti GBM and ANCA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity
- Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

BUF	WASH	30X
DIL		
CONJ		

CONTROL	+
CONTROL	-
SUBS	pNPP

Warning

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

- H317: May cause an allergic skin reaction.
 P264: Wash hands thoroughly after handling.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P302+352: IF ON SKIN: Wash with plenty of soap and water.
 P333+313: If skin irritation or rash occurs: Get medical advice/attention.

Specimen collection

The anti-GBM, ANCA screen assay is for use with serum. Handle as if capable of transmitting infections agents.

Avoid using sera which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8 °C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20 °C or colder.

Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

Kit components and storage of reagents

- One frame with doublestrips (6x (2x8) wells) coated with Goodpasture antigen (Bovine source), purified proteinase 3 (Human Neutrophil source) and myeloperoxidase (Human Neutrophil source), one lid sealed in a foil pack with a dry pack.
- 0.75 mL negative control (NC) containing human serum in diluent (green colour).
- 0.75 mL positive control (PC) containing human serum in diluent (red colour).
- 13 mL conjugate containing alkaline phosphatase-labelled (Goat) antibodies to human IgG (blue colour).
- 32 mL Diluent (Dil) containing PBS (reddish colour).
- 13 Substrate pNPP.
- 30 mL Wash solution 30x concentrated.

All reagents in the kit are ready for use except wash solution and should be stored at + 2-8 ° C. Remove only the number of wells needed for testing, resealing the aluminium package carefully.

Materials or equipment required but not provided

Microplate reader with filter 405 nm. NaCl, Tween 20 and distilled water for preparing more washing solution.

Washer for strips, absorbent tissue, tubes and a timer. Precision pipettes with disposable tips.

Do not use reference wavelength (620 nm or similar) when reading the plate. Do not use reference wavelength subtracted OD values for the calculations.

PROCEDURE

All solutions should be used at room temperature.

Remove only the number of wells needed for testing, resealing the aluminium package carefully.

If the assay is performed without a shaker, incubation times should be adjusted as follows: serum incubation 10 minutes, no change; conjugate incubation 20-30 minutes; substrate incubation 20-30 minutes.

Preparation of washing solution

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37°C water bath until the crystals have dissolved before dilution of wash solution. Dilute 10 mL of the 30x concentrated wash solution in 290 mL distilled water. When stored at 2-8 ° C, the diluted wash solution is stable until the date of expiration of the kit.

Dilution of serum and incubation

1. Pipet 75 µL diluent (Dil) into each well of a double microtitre strip (not G2+H2).
2. Pipet 25 µL negative control (NC) into each of wells A2, C2 and E2.
3. Pipet 25 µL positive control (PC) into each of wells B2, D2 and F2.
4. Pipet 25 µL patient serum into each well of column 1 according to the schema below. Put on the lid and incubate for 10 minutes on a shaker.

	Column 1	Column 2	Antigen coat
A	75 µL Dil + 25 µL serum	75 µL Dil + 25 µL NC	GP antigen
B	75 µL Dil + 25 µL serum	75 µL Dil + 25 µL PC	GP antigen
C	75 µL Dil + 25 µL serum	75 µL Dil + 25 µL NC	PR3 antigen
D	75 µL Dil + 25 µL serum	75 µL Dil + 25 µL PC	PR3 antigen
E	75 µL Dil + 25 µL serum	75 µL Dil + 25 µL NC	MPO antigen
F	75 µL Dil + 25 µL serum	75 µL Dil + 25 µL PC	MPO antigen
G	75 µL Dil + 25 µL serum		HSA as blank
H	75 µL Dil + 25 µL serum		Empty as blank

After serum incubation

Wash 4 times with washing solution, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding conjugate

Add 100 µL conjugate to each well. Put on the lid and incubate for 10 minutes on a shaker (without a shaker, incubate 20-30 minutes).

After conjugate incubation

Wash 4 times with washing solution, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding substrate

Add 100 µL substrate pNPP to each well, incubate for 10-20 minutes (without a shaker, incubate 20-30 minutes). Read the absorbance at 405 nm on a microplate reader at 10-20 minutes after substrate addition and 20-30 minutes after substrate addition without a shaker.

Do not use reference wavelength (620 nm or similar) when reading the plate.

Calculations

Do not use reference wavelength subtracted OD values for the calculations.

An optical density (OD) ratio for each patient sample is calculated as follows:

$$\text{OD ratio} = \frac{\text{OD patient sample}}{\text{OD of negative control}}$$

The patient sample is negative if the OD ratio is < 3.0. Equivocal if the OD ratio is 3.0-4.0 and positive if the OD ratio is > 4.0.

Quality Control

All calculations should only be made on the 405 nm data, without reference wavelength OD subtraction.

The OD for the NC for each antigen should be < 0.3.

The OD for the PC for each antigen should be > 0.9

If the PC or NC values are not within their respective ranges, the test should be considered invalid and the test should be repeated.

The OD value for the patient sample in the HSA and in the uncoated wells should be < 0.3.

If the OD value for the patient sample is >0.3 in some of the HSA or uncoated wells or positive on more than one antigen the sample may contain unspecific reactivities.

Sometimes unspecific reactivities will disappear with a higher dilution.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off. It is recommended to assay an additional control at the assay cut-off.

Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to NCCLS C24-A for guidance on appropriate QC practices.

Interpretation of results

A sample with a OD ratio of:

< 3.0 = **Negative** 3.0 – 4.0 = **Equivocal**; Retest, if still equivocal retest by an alternative method
> 4.0 = **Positive**

Owing to the high sensitivity of the screening kit, in a few percent of cases samples with an OD ratio > 3 may be found to be negative in a quantitative assay.

An equivocal or a positive test result should always be confirmed by a quantitative assay.

Limitations

The individual patient's OD ratio can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardization of results. The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used in combination with clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive, with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies.

Immunosuppressive therapy should not be started on basis of a positive anti-GBM, ANCA result.

Initiation or changes in treatment should not be based on changes in anti-GBM, ANCA concentration alone, but rather on careful clinical observation.

Expected results

ANCA and anti-GBM are rarely found in normal healthy individuals. The anti-GBM, ANCA screening kit was tested with 133 normal sera. 126 were found to be negative, 7 sera were in the equivocal range. The annual incidence of Goodpasture syndrome is one per million and one new patient with PR3-ANCA

is expected per 100,000 inhabitants per year. Around 10% of patients with WG are negative in both IIF and ELISA.

Anti-GBM was tested with 38 GP sera, 38 were found to be positive.

The MPO-ANCA was tested with 42 WG sera, 4 were found to be positive. The MPO-ANCA was tested with 43 MP sera, 20 were found to be positive. The PR3-ANCA was tested with 42 WG sera, 39 were found to be positive. The PR3-ANCA was tested with 43 MP sera, 21 were found to be positive. (Table 1)

Performance characteristics

Table 1. Clinical sensitivity and specificity. A total of 326 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total	Negative < 3			Equivocal 3-4			Positive >4		
		GBM	PR3	MPO	GBM	PR3	MPO	GBM	PR3	MPO
Blood donors: (NS)	131	128	131	127	3	0	4	0	0	0
WG:	42	-	3	37	-	0	1	-	39	4
MP:	43	-	20	23	-	2	0	-	21	20
SLE:	31	31	31	24	0	0	2	0	0	*5
RA:	41	41	41	40	0	0	1	0	0	0
GP:	38	0	-	-	0	-	-	38	-	-

WG = Wegener's granulomatosis,
SLE = systemic lupus erythematosus

MP = microscopic polyangiitis
GBM = glomerular basement membrane

RA = rheumatoid arthritis

* All 5 sera was positive in quantitative MPO-ANCA test.

Clinical sensitivity (Equivocal samples are not included in the calculation)

PR3-ANCA

WG = 39/42 = 92.9 % 95% CI = 84.9 – 100%. MP = 21/41 = 51.2 % 95% CI = 35.6 – 66.8%

MPO-ANCA

WG = 4/41 = 9.8 % 95% CI = 4.9 – 19.0%. MP = 20/43 = 46.5 % 95% CI = 31.3 – 61.7%

Anti-GBM

GP = 38/38 = 100 % 95% CI = 92.2 – 100%

Clinical specificity (Equivocal samples are not included in the calculation)

PR3-ANCA

SLE = 31/31 = 100 % 95% CI = 90.4 – 100%. NS = 131/13 = 100 % 95% CI = 97.7 – 100%
RA = 41/41 = 100 % 95% CI = 92.7 – 100%.

MPO-ANCA

SLE = 24/29 = 82.8 % 95% CI = 68.7 – 96.8% RA = 40/40 = 100 % 95% CI = 92.6 – 100%
NS = 127/127 = 100 % 95% CI = 97.6 – 100%

Anti-GBM

SLE = 31/3 = 100 % 95% CI = 90.4 – 100%. RA = 41/41 = 100 % 95% CI = 92.7 – 100%
NS = 128/128 = 100 % 95% CI = 97.6 – 100%

Table 2. Relative sensitivity and specificity of the Wieslab® anti-GBM, ANCA screen kit compared to an alternative semi-quantitative ELISA. A total of 216 frozen retrospective sera were assayed on the Wieslab® anti-GBM, ANCA screen kit and a semi-quantitative ELISA for PR-3 and

MPO. Also, a total of 169 frozen retrospective sera were assayed on the Wieslab® anti-GBM, ANCA screen kit and an semi-quantitative ELISA for anti-GBM. The following table summarises the results.

Wieslab® anti-GBM, ANCA screen										
Semi-quantitative ELISA		Negative < 3			Equivocal			Positive >4		
		GBM	PR3	MPO	GBM	PR3	MPO	GBM	PR3	MPO
MPO-ANCA	Positive	0	0	1	0	0	0	0	0	23
PR3-ANCA	Positive	0	1	0	0	0	0	0	59	0
Anti-GBM	Positive	0	-	-	0	-	-	37	-	-
	Negative	128	152	182	3	2	4	0	0	0
	Equivocal	0	1	4	0	0	1	1	1	1
	Total	128	154	187	3	2	5	38	60	24

Relative sensitivity (Equivocal samples are not included in the calculation)

Relative sensitivity PR3-ANCA = 59/60 = 98.3% 95% CI = 95.0 – 100%
 Relative sensitivity MPO-ANCA = 23/24 = 95.8% 95% CI = 87.7 – 100%
 Relative sensitivity anti-GBM = 37/37 = 100 % 95% CI = 92.0 – 100%

Relative specificity (Equivocal samples are not included in the calculation)

Relative specificity PR3-ANCA = 152/152 = 100 % 95% CI = 98.0 – 100%
 Relative specificity MPO-ANCA = 182/182 = 100 % 95% CI = 98.4 – 100%
 Relative specificity anti-GBM = 128/128 = 100 % 95% CI = 97.7 – 100%
 The 95% confidence interval (CI) was calculated using the exact method.

Table 3. Batch to batch variation was determined by testing three different samples. Results were obtained for 4 different batches.

Sample	Mean	SD	CV%	Sample	Mean	SD	CV%	Sample	Mean	SD	CV%
PR3	OD ratio			MPO	OD ratio			GBM	OD ratio		
2	37.5	3.8	10	3	16.8	2.2	13	1	21.8	2.1	9
5	24.3	3.6	15	6	31.5	1.9	6	4	29.0	2.2	4
8	35.5	2.4	7	9	28.3	1.5	5	7	33.8	2.4	7

Table 4. Inter-assay precision was determined by testing one sample. Results were obtained for six different runs.

Sample	Mean	SD	CV %	Sample	Mean	SD	CV %	Sample	Mean	SD	CV%
PR3	OD ratio			MPO	OD ratio			GBM	OD ratio		
PK	27.3	2.9	11	PK	17.3	3.8	21	PK	16.5	1.1	6
K5	12.1	1.2	10	K6	16.5	0.84	5	K4	3.6	0.25	7

Table 5. Intra-assay precision was determined by testing one sample in 22 wells.

Sample	Mean	SD	CV %	Sample	Mean	SD	CV %	Sample	Mean	SD	CV%
PR3	OD			MPO	OD			GBM	OD		
PK	1,3	0,07	6	PK	1,8	0,06	3	PK	1,2	0,19	16
K5	1,3	0,06	5	K6	1,46	0,07	5	K4	0,6	0,03	5

Troubleshooting

Problem:	Possible causes:	Solution:
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagents. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water used to prepare solution. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto-dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

REFERENCES. RÉFÉRENCES. REFERENCIAS. LITERATUR. BIBLIOGRAFIA. LITTERATUR. REFERANSER. REFERENSER

- 1. van der Woude FJ, Rasmussen N, Lobatto S, Wiik A, Permin H, Van Es L et al. Autoantibodies against neutrophils and monocytes: Tool for diagnosis and marker of disease activity in Wegener's granulomatosis. Lancet 1985;i:425-429.**
- 2. Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): Manual of biological markers of disease. Dordrecht, Kluwer Academic Publishers; 1993:A9,1-14.
- 3. Rasmussen N, Sjolín C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). J Immunol Methods 1990;127:139-145.
- 4. Wieslander J.** How are anti neutrophil cytoplasmic autoantibodies detected? Am J Kidney Dis 1991;18:154-158.
- 5. Jenne DE, Tschopp J, Ludemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. Nature 1990;346:520.
- 6. Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
- 7. Falk RJ, Jennette JC.** Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. N Engl J Med 1988;318:1651-1657.
- 8. Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): Manual of Biological Markers of Disease. Dordrecht, Kluwer; 1994:1-9.
- 9. Jennette C, Falk R.** Small vessels vasculitis. N Engl J Med 1997, 337, 1512-1525.
- 10. Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. Kidney Int 1995, 48, 844-850.
- 11. Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The T, van der Hem GK et al.** Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. Lancet 1990; 336:709-711.
- 12. Savige J, Davis D, Falk DJ, Jennette JC, Wiik A.** ANCA and associated diseases; A review of clinical and laboratory features. Kidney Int 2000, 57, 846-862.
- 13. Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. Am J Clin. Pathol 1999, 111, 507-513.
- 14. Westman K, Bygren P, Eilert I, Wiik A, Wieslander J.** Rapid screening assay for anti-GBM antibody and ANCAs; an important tool for the differential diagnosis of pulmonary renal syndroms. Nephrol Dial transplant 1997, 12, 1863-1868.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklrning der Symbole. La spiegazione di simboli. Forklaring til symboler. Symbolforklaring. Frklaringar till symboler.

	Batch code. Numro de lot. Nmero de lote. Chargen-Nummer. Numero di lotto. Cdigo do lote. Partinummer. Lot nummer. Satsnummer.
	Catalogue number. Rfrence catalogue. Nmero de catlogo. Katalog-Nummer. Numero di catalogo. Nmero catalogo. Katalognummer.
	Use-by date. Date de premption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udlbsdato. Utlpsdato. Anvnd fre.
	Temperature limit. Seuls de tempratures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Frvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biolgico. Biologische Gefhrdung. Rishio biologico. Risco biolgico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instrues de utilizao. Se brugsanvisning. Se bruksanvisningen. Ls instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif mdical de diagnostic in vitro. Producto sanitario para diagnstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos mdicos para diagnstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atensin. Achtung. Attenzione. Ateno. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend fr 96 Tests. Contenuto sufficiente per 96 test. Inholder tilstrkkelig for 96 test. Inholder tilstrekkelig for 96 test. Innehller tillrckligt fr 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformment  la directive europenne 98/79/CE relative aux dispositifs mdicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformit alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rdets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. verensstmmer med direktiv 98/79/EG fr medicintekniske produkter.

Ag	Antigen. Antigène. Antigeno, Antigene. L'antigene. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Diluent. Fortynning. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Conjugato. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x conc. Soluzione di lavaggio 30x conc. Vaskebuffert 30x conc. Vaskeløsning 30x kons. Tvättbuffert 30x conc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
CONTROL X	Control. Kontrolle. Contrôle. Controllo. Kontrol. Kontroll. Kontroll. Kontroll.

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- 1. van der Woude FJ, Rasmussen N, Lobatto S, Wiik A, Permin H, Van Es L et al.** Autoantibodies against neutrophils and monocytes: Tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet* 1985;i:425-429.
- 2. Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): *Manual of biological markers of disease*. Dordrecht, Kluwer Academic Publishers; 1993:A9,1-14.
- 3. Rasmussen N, Sjolín C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990;127:139-145.
- 4. Wieslander J.** How are anti neutrophil cytoplasmic autoantibodies detected? *Am J Kidney Dis* 1991;18:154-158.
- 5. Jenne DE, Tschopp J, Ludemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990;346:520.
- 6. Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Dordrecht, Kluwer; 1994:1-9.
- 7. Falk RJ, Jennette JC.** Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 1988;318:1651-1657.
- 8. Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Dordrecht, Kluwer; 1994:1-9.
- 9. Jennette C, Falk R.** Small vessels vasculitis. *N Engl J Med* 1997, 337, 1512-1525.
- 10. Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48, 844-850.
- 11. Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The T, van der Hem GK et al.** Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. *Lancet* 1990; 336:709-711.
- 12. Savige J, Davis D, Falk DJ, Jennette JC, Wiik A.** ANCA and associated diseases; A review of clinical and laboratory features. *Kidney Int* 2000, 57, 846-862.
- 13. Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin. Pathol* 1999, 111, 507-513.
- 14. Westman K, Bygren P, Eilert I, Wiik A, Wieslander J.** Rapid screening assay for anti-GBM antibody and ANCAs; an important tool for the differential diagnosis of pulmonary renal syndroms. *Nephrol Dial transplant* 1997, 12, 1863-1868.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklrning der Symbole. La spiegazione di simboli. Forklaring til symboler. Symbolforklaring. Frklaringar till symboler.

	Batch code. Numro de lot. Nmero de lote. Chargen-Nummer. Numero di lotto. Cdigo do lote. Partinummer. Lot nummer.Satsnummer.
	Catalogue number. Rfrence catalogue. Nmero de catlogo. Katalog-Nummer. Numero di catalogo. Nmero catalogo. Katalognummer.
	Use-by date. Date de premption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udlbsdato. Utlpsdato. Anvnd fre.
	Temperature limit. Seuls de tempratures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Frvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biolgico. Biologische Gefhrdung. Rishio biologico. Risco biolgico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instrues de utilizao. Se brugsanvisning. Se bruksanvisningen. Ls instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif mdical de diagnostic in vitro. Producto sanitario para diagnstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos mdicos para diagnstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atensin. Achtung. Attenzione. Ateno. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend fr 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrkkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehller tillrckligt fr 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformment  la directive europenne 98/79/CE relative aux dispositifs mdicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformit alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rdets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. verensstmmes med direktiv 98/79/EG fr medicintekniske produkter.

Ag	Antigen. Antigène. Antigeno, Antigene. L'antigene. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Diluent. Fortynning. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Conjugato. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x conc. Soluzione di lavaggio 30x conc. Vaskebuffert 30x conc. Vaskeløsning 30x kons. Tvättbuffert 30x conc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
CONTROL X	Control. Kontrolle. Contrôle. Controllo. Kontrol. Kontroll. Kontroll. Kontroll.

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1. **van der Woude FJ, Rasmussen N, Lobatto S, Wiik A, Permin H, Van Es L et al.** Autoantibodies against neutrophils and monocytes: Tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet* 1985;i:425-429.
2. **Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): *Manual of biological markers of disease*. Dordrecht, Kluwer Academic Publishers; 1993:A9,1-14.
3. **Rasmussen N, Sjolín C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990;127:139-145.
4. **Wieslander J.** How are anti neutrophil cytoplasmic autoantibodies detected? *Am J Kidney Dis* 1991;18:154-158.
5. **Jenne DE, Tschopp J, Ludemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990;346:520.
6. **Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Dordrecht, Kluwer; 1994:1-9.
7. **Falk RJ, Jennette JC.** Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 1988;318:1651-1657.
8. **Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Dordrecht, Kluwer; 1994:1-9.
9. **Jennette C, Falk R.** Small vessels vasculitis. *N Engl J Med* 1997, 337, 1512-1525.
10. **Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48, 844-850.
11. **Cohen Tervaert JW, Huitema MG, Hene RJ, Sluiter WJ, The T, van der Hem GK et al.** Prevention of relapses in Wegener's granulomatosis by treatment based on antineutrophil cytoplasmic antibody titre. *Lancet* 1990; 336:709-711.
12. **Savigne J, Davis D, Falk DJ, Jennette JC, Wiik A.** ANCA and associated diseases; A review of clinical and laboratory features. *Kidney Int* 2000, 57, 846-862.
13. **Savigne J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin. Pathol* 1999, 111, 507-513.
14. **Westman K, Bygren P, Eilert I, Wiik A, Wieslander J.** Rapid screening assay for anti-GBM antibody and ANCAs; an important tool for the differential diagnosis of pulmonary renal syndroms. *Nephrol Dial transplant* 1997, 12, 1863-1868.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklrning der Symbole. La spiegazione di simboli. Forklaring til symboler. Symbolforklaring. Frklaringar till symboler.

	Batch code. Numro de lot. Nmero de lote. Chargen-Nummer. Numero di lotto. Cdigo do lote. Partinummer. Lot nummer.Satsnummer.
	Catalogue number. Rfrence catalogue. Nmero de catlogo. Katalog-Nummer. Numero di catalogo. Nmero catalogo. Katalognummer.
	Use-by date. Date de premption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udlbsdato. Utlpsdato. Anvnd fre.
	Temperature limit. Seuls de tempratures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Frvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biolgico. Biologische Gefhrdung. Rishio biologico. Risco biolgico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instrues de utilizao. Se brugsanvisning. Se bruksanvisningen. Ls instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif mdical de diagnostic in vitro. Producto sanitario para diagnstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos mdicos para diagnstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atensin. Achtung. Attenzione. Ateno. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend fr 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrkkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehller tillrckligt fr 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformment  la directive europenne 98/79/CE relative aux dispositifs mdicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformit alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rdets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. verensstmmer med direktiv 98/79/EG fr medicintekniske produkter.

Ag	Antigen. Antigène. Antigeno, Antigene. L'antigene. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Diluent. Fortynning. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Conjugato. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffert 30x konc. Vaskeløsning 30x kons. Tvättbuffert 30x konc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
CONTROL X	Control. Kontrolle. Contrôle. Controllo. Kontrol. Kontroll. Kontroll. Kontroll.

SVAR LIFE SCIENCE AB
Lundavägen 151, SE-212 24 Malmö, Sweden
Phone: +46 40 53 76 00
E-mail: info@svarlifescience.com
www.svarlifescience.com