Instruction

**WIESLAB® Anti-GBM semi quantitative kit
(Goodpasture syndrome)**

Enzyme immunoassay for detection of autoantibodies
against glomerular basement membrane

Microtitration strips (12x8) 96 wells
Store the kit at +2-8° C
For in vitro diagnostic use only

Document No. LABEL-DOC-0036
December 2018

**WIESLAB® Anti-GBM semi quantitative kit
(Goodpasture syndrome)**

English:	page	2
Français:	page	12
Español:	page	23
Deutsch:	seite.....	30
Italiano:	pagina	37
Português	pagina	44
Dansk:	side.....	51
Norsk:	side.....	57
Svenska:	sida.....	64

REF

GP 104X

IVD

INTENDED USE

The Wieslab® anti GBM test kit is an enzyme-linked immunosorbent assay (ELISA) for detection and semi-quantitation of IgG antibodies to glomerular basement membrane (GBM) in human sera. The assay is used to detect antibodies in a single serum specimen. The results of the assay are to be used as an aid to the diagnosis of Goodpasture syndrome. The analysis should be performed by trained laboratory professionals.

FOR IN VITRO DIAGNOSTIC USE.

Summary and explanation

Goodpasture syndrome is characterised by lung haemorrhage, renal failure and the presence of anti-GBM antibodies. Diagnosis is based on clinical signs of lung hemorrhage and rapidly progressive glomerulonephritis combined with findings of anti-GBM antibodies.

As there are several other autoimmune diseases which may present with similar symptoms, this kit is a useful aid in differentiation between these diseases.

Less than one third of patients with reno-pulmonary syndromes have antibodies against the Goodpasture antigen, the majority having either Proteinase 3-ANCA or Myeloperoxidase-ANCA (1). Historically, indirect immunofluorescence was used to detect anti-GBM antibodies. When the first ELISA based on a collagenase digest was published in 1981 (2), assays using crude extracts were the only alternative. In 1984 the specific antigen was shown to derive from the C-terminal domain of type IV collagen (alpha 3 chain) (3, 4, 5, 6, 7, 11, 12, 13) and sensitive and specific assays were subsequently developed (8). The molecular nature of the Goodpasture antigen is reviewed in (9, 14). The antigen is characterised by a restricted tissue distribution occurring mainly in kidney and lungs 10).

The sensitivity and specificity of ELISA to detect anti-GBM antibodies is high (15). False positive reactions occur, mainly in SLE and other diseases with polyclonal activation, most of which yield low titers of 10-20 U/ml.

Principle of the Wieslab® anti-GBM kit

The wells of the microtitre strips are coated with purified GBM antigen (alfa 3 chain). During the first incubation, specific antibodies in diluted serum, will bind to the antigen coating.

The wells are then washed to remove unbound antibodies and other components.

A conjugate of alkaline phosphatase-labelled antibodies to human IgG binds to the antibodies in the wells in this second incubation.

After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour intensity and is measured in terms of absorbance (optical density (OD)). The absorbance is then calculated against a calibrator curve and the results are given in arbitrary U/ml.

Warnings and precautions

- For in vitro diagnostic use.
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of anti-GBM in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.

- Material safety data sheets for all hazardous components contained in this kit are available on request from Svar Life Science.

BUF	WASH	30X
DIL		
CONJ		
CAL		

CONTROL	+
CONTROL	-
SUBS	pNPP

Warning

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

H317:	May cause an allergic skin reaction.
P264:	Wash hands thoroughly after handling.
P280:	Wear protective gloves/protective clothing/eye protection/face protection.
P302+352:	IF ON SKIN: Wash with plenty of soap and water.
P333+313:	If skin irritation or rash occurs: Get medical advice/attention.

Specimen collection

The anti-GBM assay is for use with serum. Handle as if capable of transmitting infectious agents. Avoid using sera which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8° C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20° C or colder. Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

Kit components and storage of reagents

- One frame with strips (12x8) coated with GBM antigen (alfa 3 chain) one lid sealed in a foil pack with a dry pack.
- 1,5 ml negative control (NC) containing human serum in diluent.
- 1,5 ml positive control (PC) containing human serum in diluent.
- 13 ml conjugate containing alkaline phosphatase-labelled antibodies to human IgG (blue colour).
- 2 x 32 ml Diluent (Dil) containing PBS (red colour).
- 13 ml Substrate pNPP.
- 30 ml wash solution 30x concentrated.
- Five calibrators containing human serum in diluent. 1.5 ml Cal 1 = 320 U/ml, 1.5 ml Cal 2 = 160 U/ml, 1.5 ml Cal 3 = 80 U/ml, 1.5 ml Cal 4 = 40 U/ml, 1.5 ml Cal 5 = 10 U/ml.

All reagents in the kit are ready for use except wash solution and should be stored at 2-8° C.

Remove only the number of strips needed for testing, resealing the aluminium package carefully.

Materials or equipment required but not provided

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.

PROCEDURE

All solutions should be used at room temperature. Incubate in all steps at room temperature (20-25° C) with lid. Incubate serum for 30 minutes, conjugate for 30 minutes and substrate for 60 minutes (\pm 10 minutes).

Preparation of washing solution

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37°C water bath until the crystals have dissolved before dilution of wash solution. Dilute 10 ml of the 30x concentrated wash solution in 290 ml distilled water. When stored at 2-8° C, the diluted wash solution is stable until the date of expiration of the kit.

Dilution of serum and incubation

Dilute the patient's serum 1/80 with diluent (395 μ l diluent +5 μ l serum).

Pipet 100 μ l/well in duplicate of diluent (as a blank), Calibrator 1, 2, 3, 4, 5, NC, PC and diluted patient's serum (P) according to the diagram below. Incubate for 30 minutes.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Dil	Cal 4	P1									
B	Dil	Cal 4	P1									
C	Cal 1	Cal 5	P2									
D	Cal 1	Cal 5	P2									
E	Cal 2	NC	etc									
F	Cal 2	NC										
G	Cal 3	PC										
H	Cal 3	PC										

After serum incubation

Wash 3 times with 300 μ l washing solution / well, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding conjugate

Add 100 μ l conjugate to each well. Incubate for 30 minutes.

After conjugate incubation

Wash as before.

Adding substrate solution

Add 100 μ l substrate pNPP to each well, incubate for 60 minutes (\pm 10 minutes).

Read the absorbance at 405 nm on a microplate reader.


Calculations

Subtract the OD value for the blank from the other OD values.

Construct a calibrator curve by plotting the OD against the U/ml values of the 5 calibrators. The five calibrators provided have been assigned arbitrary values of 320 U/ml for calibrator 1, 160 U/ml for calibrator 2, 80 U/ml for calibrator 3, 40 U/ml for calibrator 4 and 10 U/ml for calibrator 5. Read the U/ml value of the patient from the constructed curve. Values greater than 320 should be reported as >320, or reassay them with a higher dilution.

Arbitrary U/ml have been adopted by Wieslab®, as no generally recognised international standard exists for expressing GBM titres.

Example:	Calibrator	U/ml	Absorbance
	1	320	2.3
	2	160	1.8
	3	80	1.2
	4	40	0.7
	5	10	0.18

A sample with an absorbance value of 0.7 will read on the X-axis as having 40 U/ml of anti-GBM.

Important: The curve is an example and should not be used for actual patient interpretation.

Quality Control

The OD for the negative control should be less than that of calibrator 5.

The OD for calibrator 1 should be > 1,0

The value for the positive and negative controls see lot certificate.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off. It is recommended to assay an additional control at the assay cut-off.

If any of the values are not within their respective ranges, the test should be considered invalid and the test should be repeated. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to NCCLS C24-A for guidance on appropriate QC practices.

Interpretation of results

< 10 U/ml = **Negative**

10-20 U/ml = **Equivocal**; Retest, if still equivocal retest by an alternative method or test a new sample

>20 U/ml = **Positive**

Limitations

The individual patient's antibody titre can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardisation of results.

The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used as an adjunct to clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive for anti-GBM antibodies with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies.

Immunosuppressive therapy should not be started on basis of a positive anti-GBM result. Initiation or changes in treatment should not be based on changes in anti-GBM concentration alone, but rather on careful clinical observation.

Expected results

Anti-GBM are rarely found in normal healthy individuals. The anti-GBM ELISA was tested with 120 normal sera. 120 were found to be negative. The annual incidence of Goodpasture syndrome is one per million inhabitants per year.

Anti-GBM antibodies are found in the sera of most Goodpasture syndrome patients. The anti-GBM ELISA was tested with 40 GP sera, 38 (95%) were found to be positive.

Performance characteristics

Table 1. Clinical sensitivity and specificity. A total of 199 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total number	Negative <10 U/ml	Equivocal 10-20 U/ml	Positive >20 U/ml
Blood donors:	80	80	0	0
WG / MP	27	27	0	0
SLE:	19	19	0	0
Others:	33	33	0	0
Anti-GBM:	40	0	2	38

WG = Wegener's granulomatosis, MP = microscopic polyangiitis

SLE = systemic lupus erythematosus Others = RA, UC, etc.

GBM = glomerular basement membrane

Clinical sensitivity (Equivocal samples excluded from calculations)

GBM = $38/38 = 100\%$ 95% CI = 90.7-100 %

Clinical specificity (Equivocal samples excluded from calculations)

WG/MP = $27/27 = 100\%$ 95% CI = 87.2-100 %

SLE = $19/19 = 100\%$ 95% CI = 82.4-100 %

Others = $33/33 = 100\%$ 95% CI = 89.4-100 %

Donors = $80/80 = 100\%$ 95% CI = 95.5-100 %

The 95% confidence interval (CI) was calculated using the exact method.

Table 2. Relative sensitivity and specificity of the Wieslab® anti-GBM kit compared to GBM IFA.

A total of 68 frozen retrospective sera were assayed. The following table summarises the results.

Wieslab® anti-GBM				
GBM IFA	Positive	Equivocal	Negative	Total
Positive	55	3	0	58
Negative	1	1	8	10
Total	56	4	8	68

Sera falling in the equivocal range were excluded from the following calculations.

Relative sensitivity = $55/55 = 100\%$ 95% CI = 93.5-100 %

Relative specificity = $8/9 = 88.9\%$ 95% CI = 51.8-99.7 %

Relative accuracy = $63/64 = 98.4\%$ 95% CI = 91.6-100 %

Table 3. Relative sensitivity and specificity of the Wieslab® anti-GBM kit compared to an alternative ELISA.

(Wieslab® 60x60x60) A total of 159 frozen retrospective sera were assayed. The following table summarises the results.

Wieslab® anti-GBM					
		Positive	Equivocal	Negative	Total
Alternative ELISA	Positive	35	0	0	35
	Equivocal	3	1	0	4
	Negative	0	1	119	120
	Total	38	2	119	159

Sera falling in the equivocal range were excluded from the following calculations.

Relative sensitivity = $35/35 = 100\%$ 95% CI = 90.0-100 %

Relative specificity = $119/120 = 99.2\%$ 95% CI = 95.4-100 %

Relative accuracy = $154/155 = 99.4\%$ 95% CI = 96.5-100 %

The 95% confidence interval (CI) was calculated using the exact method.

Table 4. Batch to batch variation was determined by testing four different samples in duplicate. Results were obtained for four different batches.

Sample	Mean value	SD	CV %
1	58 U/ml	5	9
2	45 U/ml	7	15
3	24 U/ml	2	9
PK	40 U/ml	5	11

Table 5. Inter-assay precision was determined by testing two different samples in duplicate. Results were obtained for seven different runs.

Sample	Mean value	SD	CV %
1	40 U/ml	2.6	7
2	18 U/ml	1.9	11

Table 6. Intra-assay precision was determined by testing one samples in 70 wells.

Sample	Mean value	SD	CV %
1	67 U/ml	5.4	8

Table 7. Linearity. The values were determined for serial two-fold dilutions of four positive sera. The values were compared to log 2 of dilution by standard linear regression. The data indicates that the assay has a linear relationship with serum dilution.

Serum	neat	1:2	1:4	1:8	1:16	r
1	60	36	18	9		0,977
2	217	126	82	71	26	0,988
3	98	60	34	23	6	0,998
4	69	42	19	12		0.999

Troubleshooting

Problem:	Possible causes:	Solution:
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagents. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water used to prepare solution. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5% 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

REFERENCES. RÉFÉRENCES. REFERENCIAS. LITERATUR. BIBLIOGRAFIA. LITTERATUR. REFERANSER. REFERENSER

1. **Saxena R, Bygren P, Arvastson B, Wieslander J.** Circulating autoantibodies as serological markers in the differential diagnosis of pulmonary renal syndrome. *J Intern Med.* 1995. 238:143-152.
2. **Wieslander J, Bygren P, Heinegård D.** Anti-basement membrane antibody: Immunoenzymatic assay and specificity of antibodies. *Scand. J. Clin. Lab. Invest.* 1981. 41:763-772.
3. **Wieslander J, Bygren P, Heinegård D.** Isolation of the specific glomerular basement membrane antigen involved in Goodpasture syndrome. *Proc Natl Acad Sci U. S. A.* 1984a. 81:1544-1548.
4. **Wieslander J, et al** Goodpasture antigen of the glomerular basement membrane: localization to noncollagenous regions of type IV collagen. *Proc Natl Acad Sci U. S. A.* 1984b. 81:3838-3842.
5. **Segelmark M, Butkowski R, Wieslander J.** Antigen restriction and IgG subclasses among anti-GBM autoantibodies. *Nephrol Dial Transplant* 1990.5: 991-996.
6. **Hellmark T, Johansson C, Wieslander J.** Characterization of anti-GBM antibodies involved in Goodpastures syndrome. *Kidney Int.* 1994. 46: 823-829.
7. **Hellmark T, Brunmark C, Trojner J, Wieslander J.** Epitope mapping of anti-GBM antibodies with synthetic peptides. *Clin. Exp. Immunol.* 1996, 105 504-510.
8. **Saxena R, Isaksson B, Bygren P, Wieslander J.** A rapid assay for circulating anti-glomerular basement membrane antibodies in Goodpasture syndrome. *J. Immunol. Methods* 1989 118:73-78.
9. **Hellmark T, Segelmark M, Bygren P, Wieslander J.** Glomerular basement membrane antibodies. In "Autoantibodies" Eds Peter JB, Shoenfeld Y, Elsevier Science 1996, 291-298.
10. **Kleppel MM, et al.** Human tissue distribution of novel basement membrane collagen. *Am J Path.* 1989. 134 4: 813-825.
11. **Hellmark T, et al.** Identification of a clinically relevant immunodominant region in Goodpasture disease. *Kidney Int* 1999, 55, 936-944.
12. **Segelmark M, Burkhardt H, Wieslander J.** Goodpasture disease: Characterization of a single conformational epitope as the target of pathogenic autoantibodies. *J. Biol. Chem.* 1999, 274, 25862-25868.
13. **Gunnarsson A, et al.** Molecular properties of the Goodpasture epitope. *J Biol Chem* 2000, 275, 30844-30848.
14. **Hudson B, et al..** Alport's Syndrome Goodpasture's Syndrome and type IV collagen *N Engl J Med* 2003, 348, 2543-2556.
15. **Segelmark M, Hellmark T, Wieslander J..** The prognostic significance in Goodpasture's disease of specificity, titre and affinity of anti-glomerular basement membrane antibodies. *Nephron Clin Pract* 2003, 94, 59-68.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklrning der Symbole. La spiegazione di simboli. Explicao dos smbolos Forklaring til symboler. Symboler som brukes p etiketter. Frklaringar till symboler.

	Batch code. Numero de lot. Nmero de lote. Chargen-Nummer. Numero di lotto. Cdigo do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Rfrence catalogue. Nmero de catlogo. Katalog-Nummer. Numero di catalogo. Nmero catalogo. Katalognummer.
	Use-by date. Date de premption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udlbsdato. Utlpsdato. Anvnd fre.
	Temperature limit. Seuls de tempratures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Frvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biolgico. Biologische Gefhrdung. Rishio biologico. Risco biolgico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instrues de utilizao. Se brugsanvisning. Se bruksanvisningen. Ls instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif mdical de diagnostic in vitro. Producto sanitario para diagnstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos mdicos para diagnstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atensin. Achtung. Attenzione. Ateno. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend fr 96 Tests. Contenuto sufficiente per 96 test. Inneholder tilstrkkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehller tillrckligt fr 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformment  la directive europenne 98/79/CE relative aux dispositifs mdicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformit alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rdets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. verensstmmer med direktiv 98/79/EG fr medicintekniske produkter.

Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigénio. Antigen. Antigen. Antigen.
DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Diluente. Fortynning. Diluent. Spädningsbuffert.
CONJ	Conjugate. Conjugué. Conjugado. Konjugat. Coniugato. Conjugado. Konjugat. Konjugat. Konjugat.
BUF WASH 30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Solução de lavagem concentrada 30x. Vaskebuffer 30x koncentreret.. Vaskebuffer, konsentret 30 ganger. Tvättbuffert 30x konc.
SUBS pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrat pNPP. Substrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP. Substrat pNPP.
CAL X	Calibrator. Etalon. Calibrador. Calibratore. Calibrador. Kalibrator. Kalibrator. Kalibrator. Kalibrator.
CONTROL X	Control. Contrôle. Controllo. Kontrolle. Controllo. Controllo. Kontroll. Kontroll. Kontroll.

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SVAR LIFE SCIENCE AB