

Quantification of IL-12 inhibitor activity using *iLite*® IL-12 Assay Ready Cells

For research and professional use only. Not for use in diagnostic procedures.

*This application note contains a suggested protocol and performance data.
Each individual laboratory must set up their own method and perform relevant validations.*

Background

Interleukin 12 (IL-12) is a 70 kDa cytokine primarily released upon antigenic stimulation of macrophages, neutrophils and dendritic cells. As T cell-stimulating factor, IL-12 plays a critical role in the regulation of T-helper 1 cell responses. IL-12 is composed of two subunits, p35 and p40, covalently linked by a single disulfide bond. The p40 subunit, which binds to the receptor chain IL-12R β 1, is shared with another heterodimeric cytokine, IL-23. However, the two cytokines exert distinct non-redundant biological functions (1). Therapeutic agents targeting both IL-12 and IL-23 cytokines are currently used to treat psoriasis and psoriatic arthritis, and related agents are in clinical testing for a variety of inflammatory disorders (2).

Principle of the assay

The *iLite*® IL-12 Assay Ready Cells are engineered cells optimized to express Firefly luciferase under the control of an IL-12 responsive promoter. When IL-12 binds to the hetero-dimeric cell surface receptor composed of IL-12R β 1 and IL-12R β 2 it activates the expression of the IL-12 regulated Firefly luciferase reporter gene. The Firefly luciferase signal can be measured in a luminometer following addition and incubation of luciferase substrate. The Firefly luciferase signal is proportional to the functional activity of IL-12 in the sample. In the presence of inhibitory activity against IL-12, the amount of free IL-12 is reduced, resulting in a decreased stimulation of Firefly luciferase expression. Accordingly, the Firefly luciferase signal is inversely proportional to the amount of inhibitory activity in a sample.

The *iLite*® IL-12 Assay Ready Cells can be utilized as a highly sensitive assay for quantification of IL-12 inhibitor activity in test samples, including human serum.

Material and equipment needed

Material and equipment	Suggested supplier	Reference
<i>iLite</i> ® IL-12 Assay Ready Cells	Svar Life Science	BM4012
Diluent (RPMI containing 9% heat inactivated FBS + 1% Penicillin-Streptomycin).	Gibco	61870-044 (RPMI 1640) 26140-079 (FBS) 15140-122 (Penicillin-Streptomycin)
Ustekinumab or analogues	NA	NA
IL-12 or analogues	R&D	219-IL
Firefly/Renilla luciferase substrate	Promega	E2920, Dual-Glo Luciferase Assay System
Plate; White walled micro well plate suitable for luminescence	Revvity	6055680
Microplate Luminometer with appropriate reading software – no filter on luminometer	Contact Svar Life Science for list of recommended suppliers	NA
Incubator, 37 °C with 5% CO ₂	NA	NA
Water bath, 37 °C	NA	NA

Single-channel and multi-channel pipettes with polypropylene disposable tips	NA	NA
Polypropylene tubes or plate for dilution	NA	NA
Single-use polypropylene reservoir	NA	NA
Plate shaker	NA	NA
Timer	NA	NA

Protocol

Preparation of IL-12 inhibitor

Ustekinumab has successfully been used to inhibit IL-12 and the IL-12 regulated Firefly luciferase expression in *iLite*® IL-12 Assay Ready Cells (refer to the table and graph below).




Figure 1. Example of IL-12 inhibitory curve

Final IL-12 6 ng/mL	Ustekinumab
	Suggested calibrator solution concentrations, ng/mL
A	200 000
B	33 332
C	5 556
D	924
E	156
F	25.6
G	4.4
H	0

Table 1. Suggested calibrator solution concentrations for IL-12 inhibitor

Assay preparation and incubation

1. Design a plate layout. It is recommended to perform the test at least in duplicates.
2. Perform a serial dilution of the reference IL-12 inhibitor. Ensure matrix consistency between reference inhibitor solutions, control solutions, and sample solutions.
3. Add 20 µL of the reference IL-12 inhibitor dilutions, controls and samples to assigned wells (final concentration will be a quarter of solution concentration).
4. Add 20 µL of 24 ng/mL IL-12 to all wells (final concentration will be 6 ng/mL IL-12).
5. Place the lid on the plate, mix and incubate the plate for 30 minutes at 37 °C with 5% CO₂.
6. Thaw the vial of *iLite*® IL-12 Assay Ready Cells in a 37°C water bath with gentle agitation. The cell suspension is mixed very carefully ten times with pipette to ensure a homogeneous distribution of cells.
7. Dilute 250 µL cells with 5.75 mL Diluent.
8. Add 40 µL diluted cell suspension to each well.
9. Place the lid on the plate, mix and incubate for 5 hours at 37 °C with 5% CO₂.

Adding substrate solutions

10. Equilibrate the plate and the substrate solutions to room temperature.
11. Prepare the **Firefly luciferase** substrate according to the supplier's instructions and add 80 µL per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.
12. If appropriate, prepare the **Renilla luciferase** substrate according to the supplier's instructions and add 80 µL per well. Mix and protect the plate from light. After 10 minutes incubation at room temperature read in a luminometer.

Precautions


- This application note is intended for professional laboratory research use only. The data and results originating from following the Application Note should not be used either in diagnostic procedures or in human therapeutic applications.
- Use and handle the material and instruments referenced according to the suppliers'/manufacturers' instructions or product specifications accompanying the individual material and instruments.
- Dispose of all sample specimens, infected or potentially infected material in accordance with good microbiological practice. All such materials should be handled and disposed as though potentially infectious.
- Residues of chemicals and preparations are generally considered as biohazardous waste and should be inactivated prior to disposal by autoclaving or using bleach. All such materials should be disposed of in accordance with established safety procedures.

Proprietary Information

In accepting delivery of *iLite*[®] Assay Ready Cells the recipient agrees not to sub-culture these cells, attempt to sub-culture them or to give them to a third-party recipient, and only to use them directly in assays. *iLite*[®] cell-based products are covered by patents which are the property of Svar Life Science AB and any attempt to reproduce the delivered *iLite*[®] Assay Ready Cells is an infringement of these patents

QUICK GUIDE

Quantification of IL-12 inhibitor activity using *iLite*[®] IL-12 Assay Ready Cells



Troubleshooting and FAQ

Please consult the Svar Life Science website www.svarlifescience.com

References

1. **Hsieh CS, Macatonia SE, Tripp CS, Wolf SF, O'Garra A, Murphy KM.** (April 1993). *Development of TH1 CD4+ T cells through IL-12 produced by Listeria-induced macrophages.* Science 260: 547–5499 (1993).
2. **Teng MW, Bowman EP, McElwee JJ, Smyth MJ, Casanova JL, Cooper AM, Cua DJ.** *IL-12 and IL-23 cytokines: from discovery to targeted therapies for immune-mediated inflammatory diseases.* Nature Medicine 21: 719–729 (2015).