trophon® CPP

Manufacturer:
Nanosonics Limited
14 Mars Road, Lane Cove
NSW 2066, Australia
+61 2 8063 1600
info@nanosonics.com.au

REF

EC REP


Nanosonics Europe GmbH
Poppenbueteler Bogen 66
22399 Hamburg - Germany
+49 40 46856885

nanosonics

www.nanosonics.com.au

CE

L00032 4.0 01/2017

ENGLISH

Product Introduction

trophon® Curved Probe Positioner (CPP) is an accessory to the trophon EPR to improve the positioning of approved curved probes in the device chamber.

trophon CPP is to be used only with trophon EPR and approved curved probes. Refer to trophon EPR Validated Probes List on the Nanosonics website for more information on approved curved probes for use.

Do not use the trophon CPP if discoloured, cracked or broken.

This product includes:

- Two (2) trophon Curved Probe Positioners
- One (1) Instructions For Use (IFU)

Set Up & Use

WARNING: Take care not to touch the inside chamber when inserting and removing trophon CPP as the chamber could be hot to touch.

1. Remove the trophon CPP from packaging.
2. Wearing clean gloves, open the chamber door and refer to figure 1 and 2 for correct orientation and location in the device chamber.
3. Gently insert the trophon CPP into the device as shown in figure 3 and slide until it can't be pushed in any further.
4. Once positioned correctly, place the pre-cleaned and dry approved curved probe into the device as shown in figure 4 ensuring the probe tip is positioned towards the right side of the chamber.
5. Ensure the probe is not touching the chamber walls and is above the embossed line as per figure 5.
6. Figure 6 highlights incorrect placement of a curved probe in the device chamber. Refer to the trophon EPR User Manual for additional information.

NOTE: When the trophon CPP is not in use, remove and store in a safe and dry location.

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