

Application Note


One-step immunostaining with Nano-Secondaries

Introduction: What is a Nano-Secondary?

ChromoTek's Nano-Secondaries are a novel class of secondary antibodies for higher resolution, cleaner images, and faster immunostaining. Nano-Secondaries consist of Nanobodies/ V_HHs that bind to primary antibodies with high affinity in a species and subclass specific manner. Nano-Secondaries are conjugated to Alexa Fluor® dyes.

Introduction: One-step immunostaining

Because the Nano-Secondaries are monovalent and bind with high specificity and affinity to their target IgGs, they can be simultaneously incubated with the primary antibody. This results in a one-step immunostaining. It saves incubation time and reduces washing steps, and hand-on time. Simultaneous incubation also supports multiplexing, live-cell immunostaining, and improves cell viability for flow cytometric analysis.



One-step staining protocol for immunofluorescence detection


The protocol below provides guidelines for one-step immunostaining of cultured adherent mammalian cells. For immunostaining of other cell types, tissues or whole organs please adjust accordingly.

Note, the efficiency of one-step immunostaining depends on the primary antibody. Carefully follow the supplier's recommendations for the specific primary antibody you are using, especially regarding cell fixation (PFA or acetone/methanol), permeabilization reagent, blocking buffer composition, incubation time (1 h or overnight) and temperature (room temperature or 4°C). It helps a lot to pre-test your primary antibody with respect to their optimal dilution, isotype, and correct staining pattern using conventional two-step immunostaining before moving to a faster one-step protocol.

| Step | Basic condition | Alternative condition |
|----------------------------------|------------------------|-------------------------------------|
| 1. Fix cells | 4% PFA-PBS – 10 min RT | (acetone/) methanol - 10 min on ice |
| 2. Wash cells 3-5 times with PBS | | |


| | | |
|---|---|--|
| 3. Permeabilize cells | 0.1% - 0.5% Triton-X100 in PBS - 5 min RT | saponin, digitonin, NP40 in PBS or no permeabilization |
| 4. Wash cells 3-5 times with PBS | | |
| 5. Block cells | 4% BSA-PBS - 20 min RT | 5% normal goat serum, dry milk, cold fish skin gelatin |
| 6. Incubate cells with primary antibody and Nano-Secondary in blocking solution | 1 h RT | overnight at 4°C; 15-30 min at 4°C if live-cell |
| 7. Wash cells 3-5 times with PBS | | |
| 8. Counterstain with DAPI, mount, image | | |

Simultaneous one-step immunostaining vs. sequential immunostaining. HeLa cells were immunostained with different primary rabbit antibodies (left anti-Actin, right anti-Lamin) and secondary alpaca anti-rabbit IgG, recombinant V_HH, Alexa Fluor[®] 647 (1:1,000). Scale bar, 20 μm.



What is a Nanobody or V_HH?

In addition to conventional IgG antibodies, alpacas also possess heavy chain only IgGs. These antibodies lack the C_H1 domain of the heavy chain and are devoid of light chains. Their antigen binding domain is built up solely by their heavy chain and is called V_HH or Nanobody. Nano-Secondaries are Nanobodies against rabbit IgG or mouse IgG subclasses that are chemically conjugated to Alexa Fluor dyes.



Products

| Nano-Secondaries | Product Size | Product Code |
|---|--------------|--------------------------|
| Alpaca anti-rabbit IgG, recombinant VHH, Alexa Fluor® 488, 568, 647 | 10 µl | srb1AF488/568/647-1-10 |
| | 100 µl | srb1AF488/568/647-1-100 |
| Alpaca anti-mouse IgG1, recombinant VHH, Alexa Fluor® 488, 568, 647 | 10 µl | sms1AF488/568/647-1-10 |
| | 100 µl | sms1AF488/568/647-1-100 |
| Alpaca anti-mouse IgG2b, recombinant VHH, Alexa Fluor® 488, 568, 647 | 10 µl | sms2bAF488/568/647-1-10 |
| | 100 µl | sms2bAF488/568/647-1-100 |
| Alpaca anti-mouse IgG3, recombinant VHH, Alexa Fluor® 647 | 10 µl | sms3AF647-1-10 |
| | 100 µl | sms3AF647-1-100 |

| Related Products | Product Size | Product Code |
|---|--------------|--------------|
| Mouse monoclonal IgG1 antibody [28A] to Halo-tag | 20 µl | 28a8-20 |
| | 100 µl | 28a8-100 |
| GFP antibody rabbit polyclonal [PABG1] | 20 µl | PABG1-20 |
| | 100 µl | PABG1-100 |

For product details, information, and ordering visit www.chromotek.com.

Contact: support@chromotek.com

For research use only. Not for diagnostic or therapeutic use.

ChromoTek is a registered trademark of ChromoTek GmbH. Nanobody is a registered trademark of Ablynx, a Sanofi company. Alexa Fluor® is a registered trademark of Life Technologies Corporation, a part of Thermo Fisher Scientific Inc. Other suppliers' products may be trademarks or registered trademarks of the corresponding supplier each. Statements on other suppliers' products are given according to our best knowledge.