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## **Instruction Manual**

*Please read before use!*

# **opus® dispenser**

*Motorisches Dosieren*

*Motorised dispensing*

D9580002 / 0109

Hirschmann Laborgeräte  
GmbH & Co. KG  
Postfach 11 53 - Hauptstr. 7 - 15  
D- 74246 Eberstadt  
Germany  
Tel. +49(0)71 34/511-0  
Fax +49(0)71 34/511-90  
[www.hirschmannlab.de](http://www.hirschmannlab.de)  
[info@hirschmannlab.de](mailto:info@hirschmannlab.de)

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## 1. Warning and safety instructions

### **Warning!**

Read the instruction manual prior to using the unit and observe all safety instructions. Any person working in the laboratory should be familiar with these safety instructions prior to using the unit, and the safety instructions should be accessible at all times. However, the instruction manual cannot illustrate every safety problem which may occur. Every user is personally responsible for observing the safety and health protection regulations and determining restrictions prior to using the unit. These restrictions should then be observed.

- Observe all general safety instructions for the laboratory (e.g. wearing of protective clothing, safety glasses and protective gloves when using certain liquids.
- Observe the instructions and regulations of reagent manufacturers.
- The unit should only be used for dispensing liquids, and it is imperative that specified usage prohibitions and

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*restrictions be observed (consult the manufacturer if necessary).*

- Prior to use, check all connections for leaks and ensure that the ejection unit and suction hose are securely fitted. A loosely-fitted ejection unit could cause splashing when the piston unit is lowered.*
- It is imperative that the closure cap be removed from the ejection hose prior to dispensing. Never activate the unit while the closure cap is fitted.*
- Ensure that neither you yourself nor other persons are at risk. Never work facing in the direction of another person. Avoid splashing and use suitable collecting vessels. Always dispense onto the interior wall of the collecting vessel rather than directly onto its base or the contents contained therein.*
- Never use force!*
- Seal the ejection hose with the closure cap after use. Caution! Reagent can drip out!*
- Use only original accessories and original spare parts.*

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- Allow components to cool to room temperature after sterilisation (see chapter 13).
- In case of malfunction, e.g.
  - piston unit stiffness or seizure
  - Sticky or leaking suction/ejection or recirculation valve
  - Leaks
  - Component breakage Stop working immediately. Observe the instructions on cleaning (see chapter 12) and problem solutions (see chapter 17) before continuing to use the unit. Send to manufacturer for repair if necessary. Observe chapter 15 in this respect!
- Do not deposit the unit or power supply unit in water or other liquids to avoid the risk of electric shock!
- Ensure that the mains power voltage corresponds to the details on the power supply unit.
- The power supply unit should not be used any further if the cable or power supply unit are damaged in any way.
- Never touch the plug on the power supply unit with wet hands!

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- The power supply unit should only be used indoors.  
Protect against damp!

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**Warning!**

Liability is excluded for the unit or any subsequent damage if the unit is opened. Avoid incorrect use!

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**Warning!**


Always use an original power supply unit and original accessories and spare parts to avoid damage to the unit.

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## 2. Individual components

1. Ejection unit (complete)
2. Recirculation valve
3. Suction valve
4. Ejection valve
5. Basic unit
6. Control unit
7. Power supply unit

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### 3. Function overview

- The unit is a motorised liquid dispensing system. The basic unit (A 45 thread) and adapter variants (see accessories) facilitate the use of commercially-available reagent bottles.
- The system employs the piston stroke principle and functions via ball valves. The piston unit is advanced by a DC motor with gear drive and increments. Actuation is realised with a microprocessor control. Power is supplied by a power supply unit. The unit can be optionally powered with a battery pack.

#### **Note:**

The upper part of the basic unit is not sealed against gas and damp. It is therefore not permissible to operate the unit in caustic fumes.

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- The suction valve is opened when the piston unit is drawn upwards. Liquid is suctioned into the cylinder from a storage bottle through the suction hose. The suction valve closes when suctioning is completed.
- The ejection valve is opened when the piston unit is lowered and the liquid in the cylinder is ejected.

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## 4. Prohibited use

### **The unit should not be used with the following:**

- Liquids which attack the materials Halar (ECTFE), FEP, PTFE, borosilicate glass 3.3 or Hastelloy
- Suspensions (e.g. active carbon suspension), as solid particles can block the valves
- Hydrofluoric acid, as it attacks glass
- Easily-combustible media (note the flashpoint)
- Use in potentially-explosive environments

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## 5. Restrictions on use

- The working temperature range for the unit and liquid lies between +4 °C and +40 °C.
- The unit should only be used and optional batteries loaded in areas where it can be assured that no hazardous potentially-explosive atmospheres occur.
- Do not deposit in potentially-explosive areas.

### **Warning – important information!**

Never use force where the piston unit has seized. Terminate the process immediately and realise basic cleaning (see chapter 12).

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
## 6. Preparing the unit

### **Warning, please note before assembly!**

- Check all parts for evidence of damage!
- Flawless and safe operation of the unit is only assured if the undamaged original accessories included in the delivery are used.
- Never use damaged accessories or accessories from other manufacturers.

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
### 6.1. Securing the suction hose

- Hold the suction hose beside the storage bottle. The hose end with the union nut should face towards the base of the bottle.
- Cut the hose at a slant level to the lower edge of the locking ring on the storage bottle.

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
- Push the suction hose on up to the stop.

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- Then tighten the union nut onto the suction valve.

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

## 6.2. Securing the recirculation hose

- Insert the recirculation hose firmly into the recirculation opening in the valve block

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
## 6.3. Securing the ejection unit

- Ejection unit with:
  - hose guide
  - union nut
  - ejection hose
  - closure cap
- Push the ejection unit on up to the stop.

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- Tighten the union nut firmly and ensure that it is securely fitted.

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
**Warning!**

- Ensure that the ejection unit is securely fitted
- Use only original ejection unit  
(see accessories – spare parts)
- Do not connect a deformed or damaged ejection unit.

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#### **6.4. Screwing the unit onto the bottle**

- Screw the unit onto the storage bottle (use the threaded adapter where appropriate).

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#### **6.5. Turning the ejection unit to the correct position**

The unit and ejection unit can be aligned in practically any direction relative to each other, thanks to two rotation axes integrated in the unit.

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### **Warning!**

- *Observe all safety instructions, usage prohibitions and restrictions on use.*
- *Wear protective clothing, safety glasses and protective gloves*

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### **6.6. Bottle label aligning unit**

- *Turn the unit into the desired position relative to the bottle label*
- *The unit should not be positioned directly at the right-hand stop of the rotation axis (see chapter 12 for dismantling position).*

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### **6.7. Battery installation (optional)**

The control unit should be opened to facilitate installation of the battery pack. All cables should be removed from the control unit prior to commencing battery installation. All 4 screws at the rear of the control unit should then be loosened (Phillips screw PZ1). The control unit should then be rested on its feet again. The electronic unit should be tilted away over the upper edge of the housing after removing the cover. Attention should be paid to the wiring during this. The battery pack is then fixed in the centre of the housing floor with the assistance of the double-sided adhesive tapes provided. Ensure that the connection cables for the batteries are in the upper right-hand corner during this. The battery connection plug should be connected to the respective connection on the PCB. The PCB can then be positioned on its base, the cover replaced and screws hand-tightened again. Attention should once more be paid to the wiring during this.

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## **6.8. Battery charging (optional)**

The unit is, in principle, ready for immediate use when connected to the mains power supply. However, the battery should be charged prior to initial commissioning if operation with the battery option is chosen.

The battery is charged when the power supply unit is connected, regardless of whether the unit is activated or deactivated. Charging is indicated by the battery indicator on the display and by the flashing red LED on the display module if the unit is switched off. Charging is terminated automatically when the battery is fully charged (battery overloading is impossible!) and a green LED flashes on the display module. The unit is also ready for operation during battery charging.

A warning message appears on the display when the battery capacity is no longer adequate to ensure correct operation. Work may only be continued in this case after the power supply unit is connected.

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
**Warning:**

The unit should only be operated with the original power supply unit. The mains voltage must correspond to the voltage specified on the power supply unit! (Max. deviation  $\pm 10\%$ ).

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## **7. Commencing work**

*Connect the basic unit and control unit with the cable provided (blue 5-pole connector and socket). Establish the power supply by connecting the power supply unit (or, optionally, the battery pack).*

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
Stopp  
Stop

Not all buttons can be used when certain operating conditions prevail, particularly if a process (dispensing, venting, etc.) is currently active. These buttons are then temporarily deactivated. Inactive buttons are indicated with a grey symbol frame.

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## 7.2. Activating/Deactivating the unit

The unit is activated by pressing the touch screen for approx. 2 seconds.

The following screen then appears:

To switch off the unit, activate the “Menu”  button and then select “OFF”. Acknowledge with .

- If the unit is only supplied with power via the optional-battery, it deactivates automatically after 5 minutes in the absence of entries by the user or if no process is currently active.
- All user configurations are retained if the unit is deactivated automatically or with the “Off” button, so that work can be continued immediately by pressing the touch screen again.

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### **Warning!**

- Safety instructions must be observed (see chapter 1.)!
- Usage prohibitions and restrictions of use must be observed (see chapter 6.)
- The unit should only be utilised with the ejection unit turned away from the user and other persons. The ejection unit should and valve block should be positioned turned away to the left or right of the body for this purpose.

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### 7.3. LED status

Basic unit: LED flashes green: Unit ready for operation  
LED flashes red: Malfunction  
(see: chapter 17)

Control unit: LED flashes green: Power supply unit  
plugged in or power supply unit plugged  
in and battery fully charged (only when the  
unit is deactivated)  
LED flashes red: Battery charging (only  
when the unit deactivated)

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
### 7.4. Closure cap removal

- Position a collecting vessel under the ejection unit
- Pull of the closure cap as shown in the illustration. **Warning: reagent can drip out during this operation!**

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- Push the closure cap to the rear as shown in the illustration.
- The closure cap can remain in the retention position while work is being realised

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### **7.5. Venting the unit**

- The recirculation lever should be positioned at  
Recirculation (the tip of the recirculation lever  
points downwards) if air bubbles visibly form in  
the cylinder chamber (e.g. after screwing onto  
a new storage vessel and subsequent filling).

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The unit is then automatically vented using the

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“Venting” symbol  in the main window.

1 mL is suctioned and the cylinder drained when the symbol is pressed. 40 % of the nominal volume is then suctioned and dispensed again. The symbol  is converted into a “Stop” symbol  during venting. The venting process can be interrupted by pressing the Stop symbol.

Venting is not compulsory if the unit was correctly deactivated after its last use and the control unit not detached from the basic unit. However, venting is recommended in this case. It is imperative in some cases that venting be realised first when the system is activated. This is indicated in that only the “Venting”  and “Menu”  button are active.

## 8. Dispensing

A previously-configured volume is emitted during

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dispensing. This volume should not exceed ten times the value of the nominal piston volume (i.e. maximum 200 mL in the case of a 20 mL piston and maximum 500 mL for a 50 mL piston). Multiple dispensing operations are possible in this respect.

### 8.1. Parameter configuration

The following parameters can be configured: the volume, suctioning and ejection speed, the number of repetitions and the cycle time. All menus are laid out with the same schematic:

### 8.2. Volume configuration

Parameters are displayed in the main window.  
- An input field in which the parameter (in this case the

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volume) is displayed at the top opens when the volume display is pressed. The current value is shown immediately beneath. Details of the minimum and maximum value are displayed adjacent to this, and these values should not be exceeded or undershot.

- Inputting of the new value is realised via the numerical sequence, and incorrect entries can be corrected with the “DEL” field.

The new value is imported by pressing the “OK”  symbol, provided it is valid. Please note the min. and max. values in the case of an invalid value.

The new value is rejected by pressing the “Abort”  symbol and the main window opens.

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### 8.3. Configuring the suctioning and ejection speed

It is possible to configure the suctioning and ejection speed independently of each other. The current values are illustrated in [mL/s] on the right of the main screen beside the respective symbols. The speeds configured can vary slightly, depending on the ambient and operating conditions.

- An input field is opened when the symbol for suctioning or ejection speed  is pressed  in which the parameter is displayed in the upper part. The current value is shown immediately beneath. Details of the minimum and maximum value are displayed adjacent to this, and these values should not be exceeded or undershot.
- Inputting of the new value is realised via the numerical sequence, and incorrect entries can be corrected with the “DEL” field.

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- The new value is imported by pressing the “OK”  symbol, provided it is valid. Please note the min. and max. values in the case of an invalid value.
- The new value is rejected by pressing the “Abort”  symbol and the main window opens.

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#### 8.4. Serial dispensing

A volume which has already been configured once can be consecutively dispensed automatically several times if required. Dispensing is realised within a previously-configured cycle time. The number of cycles  and cycle time settings are utilised to configure a process of this kind.

 The cycle time describes the time interval from the start of the dispensing procedure until the start of the next procedure. It therefore encompasses suctioning, ejection and, depending on the application, a subsequent waiting time. The minimum possible cycle time is therefore dynamically

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dependent on speed settings and dispensing volume. A serial dispensing procedure only occurs if the number of cycles entered is greater than 1. As soon as this is the case, a check is conducted for every configuration change to ascertain whether the configured cycle time  is adequate to realise a dispensing procedure. The minimum cycle time is automatically adapted to the new conditions if this is not the case, and a message informs the user of this adaptation.

### **8.5. Configuring the number of cycles**

the frequency with which the procedure should be realised can be configured here. The current value is shown on the right of the main screen beside the respective symbol. The value always indicates the number of procedures remaining if a procedure is started.

- An input field is opened when the symbol  for the number of repetitions is pressed, and the parameter is displayed in the upper part of this window. The current value

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is shown immediately beneath. Details of the minimum and maximum value are displayed adjacent to this, and these values should not be exceeded or undershot. 1 to 9999 repetitions are possible.

- Inputting of the new value is realised via the numerical sequence, and incorrect entries can be corrected with the “DEL” field.
- The new value is imported by pressing the “OK”  symbol and provided it is valid. Please note the min. and max. values in the case of an invalid value.
- The new value is rejected by pressing the “Abort”  symbol and the main window opens.

## 8.6. Configuring the cycle duration

The current value is shown on the right of the main screen in [s] beside the respective symbol. The display commences counting backwards to the start of the next procedure if a process is started.

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The cycle time for individual dispensing procedures is input by pressing the  symbol for the cycle time. An input field opens containing the current value. Details of the minimum and maximum value are displayed adjacent to this, and these values should not be exceeded or undershot. The cycle time can be up to 9999 seconds. The min. value depends on the current configurations (see chapter 8.4).

- Inputting of the new value is realised via the numerical sequence, and incorrect entries can be corrected with the “DEL”field.
- The new value is imported by pressing the “OK”  symbol, provided it is valid. Please note the min. and max. values in the case of an invalid value.

### **8.7. Starting/Stopping dispensing**

After all parameters are defined, the procedure can be started by pressing “Start” .

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This is converted immediately after into a “Stop” button . The dispensed volume is displayed during dispensing.

Press “Stop” once to interrupt the procedure. The residual volume in the piston is displayed and marked with the “RV” symbol to the left of the volume indicator on the display. The residual volume in the unit is dispensed by pressing “Start”  once again. If serial dispensing is interrupted with the “Stop” button, it can be continued at this point by pressing “Start” or pressing the “Number of cycles”  button which resets the start settings.

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## **9. Stepper**

*The Stepper is a special type of opus dispenser. Up to 10 individual partial volumes can be dispense from a cylinder charge in stepper mode. Partial volumes can be dispensed manually or time-controlled.*

*Time-controlled dispensing is defined with the delta time and cycle time parameters. The delta time describes the interval between the dispensing start of a partial volume and the automatic dispensing start of the next partial volume. The cycle time encompasses the complete interval (cylinder filling, automatic dispensing of partial volumes and, where applicable, waiting period) that passes between the start of the complete volume dispensing and the start of the next dispensing operation.*

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Press the “Menu” symbol  to switch to the stepper function. It is possible in the “Function” submenu to change between the “Stepper” and “Dispenser” functions at the press of a button.

### **9.1. Stepper parameters**

The following parameters can be configured for the stepper:

- the different partial volumes (a maximum of ten)
- The suctioning and ejection speed
- The number of repetitions
- The cycle time
- The delta time

Only the overall volume is shown in the main window. The next or active volumes to be dispensed are always displayed during the procedure. The V1-V10 display to the right

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of the volume indicator signals in the main window which is the current dispensing volume. All parameters can be configured as described in 8 (with the exception of the partial volumes and delta time).

## 9.2. Configuring partial volumes

An input field opens when the volume display is pressed, and the individual partial volumes 1 to 10 can be configured here. Pressing one of the partial volumes opens a menu which can be operated as described in 8.1. Where a partial volume is set to 0, this partial volume and all those subsequent to it are deactivated. It is therefore also possible to generate a stepper program with less than 10 partial volumes. To reuse these partial volumes, a volume can be entered again by pressing the last empty box. The settings are imported by pressing the “OK”  symbol and you return to the main window.

### 9.3. Configuring the delta and cycle time

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### 9.3. Configuring the delta and cycle time

As with dispensing, the delta and cycle time is dynamically dependent on the speed settings and dispensing volume. Press the  symbol to enter the delta and cycle time for individual dispensing procedures. The delta time is configured first. An input field opens in which the current value is displayed. Details of the minimum and maximum value are displayed adjacent to this, and these values should not be exceeded or undershot. The delta time can be up to 9999 seconds. The min. value depends on the current configurations (see chapter 8.4). The largest of the configured partial volumes serves as a basis here. Where different partial volumes are configured, the smaller volumes follow after a suitable waiting period.

- Inputting of the new value is realised via the numerical sequence, and incorrect entries can be corrected with the “DEL” field.
- If “MAN” is pressed, the individual

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partial volumes are dispensed consecutively by pressing the Start button .

- The new value is imported by pressing the “OK”  symbol, provided it is valid. Please note the min. and max. values in the case of an invalid value.

You return immediately to the main screen if the delta time is manually configured. It is otherwise possible to configure the cycle time .

The cycle time can be up to 9999 seconds. The min. value depends on the current configurations (see chapter 8.4).

- Inputting of the new value is realised via the numerical sequence, and incorrect entries can be corrected with the “DEL” field.
- The cycle time is set to manual by pressing “MAN”. Each new cycle must be started by pressing “Start”  in this case.
- The new value is imported by pressing the “OK”, provided it is valid. Please note the min. and max. values in the case of an invalid value.

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Provided the delta time setting is not set to manual, a check is conducted for every configuration change to ascertain whether the times configured are adequate to realise a dispensing procedure. The minimum delta time or cycle time is automatically adapted to the new conditions if this is not the case, and a message informs the user of this adaptation.

#### **9.4. Starting/Stopping dispensing**

After all parameters are defined, the procedure can be started by pressing “Start” . This is converted immediately after into a “Stop” button. . The volume dispensed is displayed during dispensing. Press “Stop” once to interrupt the procedure. The residual volume in the piston is displayed and marked with the “RV” symbol to the left of the volume indicator on the display. The residual volume in the unit is dispensed by pressing “Start” once again.

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
## 10. Special functions

These can be activated with the  "Menu" button.

### 10.1. Main screen

The unit is set to default mode by pressing the "Default settings" button. All user-specified settings are deleted. Resetting can be acknowledged or interrupted in the subsequent message. You can return to the main screen via the "Main" menu item.

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## 10.2. Language selection

The user can choose from four different display languages:

- German
- English
- French
- Spanish

The display language is configured in the submenu "Language"

The unit imports the new language after one of the language fields is pressed.

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### 10.3. Quick-Cal

The unit is configured by the manufacturer for a reference temperature of 20 °C, taking the cubic thermal expansion coefficient for borosilicate glass (DIN ISO 4787) into consideration (Ex 20 °C). Quality 3 water is used for this adjustment, conforming to DIN ISO 3696.

Customer calibration enables the user to calibrate his unit exactly to suit a special application case. 20 mL or 50 mL (depending on the unit) of the medium are dispensed for this purpose and then weighed on an adequately-accurate

weighing scale, with the actual volume being calculated taking the density into consideration.

The actual volume can now be entered in the “CAL” submenu. It is important to note here that the value calibrated by the customer should not deviate more than 10 % from the nominal value (an error message is generated if the deviation is too high).

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After successful activation of the value calibrated by the customer, the right-hand display changes to “ON” when “CAL” is pressed, and the “CAL” indication also appears in the main window adjacent to the volume. The unit now automatically corrects each volume to be dispensed on the basis of the calibration realised.

The value calibrated by the customer can also be deactivated again by pressing the “CAL” button again, and the unit now continues working with the default setting.

The calibration value entered once is retained after the unit is switched off and can be activated or deactivated at any time.

#### **10.4. Function**

It is possible in the “Function” submenu to change between the “Stepper” and “Dispenser” functions at the press of a button (see chapter 8 and 9).

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## 10.5. Method saving

It is possible to save all settings relating to volumes (or partial volumes), numbers of cycles , cycle time = delta time , current customer value calibration, active or inactive customer value calibration, suctioning and ejection speed as a method.

Nine memory locations are available for this purpose.

Press “Save method” in the “Method” submenu to save a method. A numerical field appears in which you can select a memory location (1-9) in which the method should be saved. Memory locations which are already occupied are indicated with red numbers, but these can be written over. The memory location selected appears above the numerical field. After pressing “OK”  the method is saved to the selected location. The indication M1-9 (depending on the saved method) now appears in the main window to the left of the volume indication.

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## 10.6. Activating the method

In order to activate a saved method, press “Active method” in the “Method” submenu. A numerical field appears in which the method to be activated can be selected. Unoccupied memory locations have a grey background and cannot be selected. Pressing “OK”  imports the new method and you change to the main window. The indication M1-9 (depending on the activated method) now appears in the main window to the left of the volume indication.

### Note

The M1-9 indication disappears on the display again as soon as a value has been changed after saving or activation of a method, as work is no longer realised according to this method. Methods for the stepper function and dispenser function are saved to a separate memory location. This means that each function has its own 9 memory locations. For example, a method saved in dispenser mode cannot be activated in stepper mode. The same applies if work is realised alternately with a 20 and 50 mL basic unit using the same control unit.

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## **10.7. Switches**

*The connection of additional devices enables the external actuation of certain buttons on the main screen. The following devices are currently available for this purpose:*

- Pedal switch*
- Hand switch*
- Mouse with two buttons*
- Control cable*

*The auxiliary device is connected to the opus control unit. Only one auxiliary device can be connected to each control unit. The pedal and hand switch each have one button function, while the mouse and control cable can be configured with two button functions.*

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- Start
- Entlüften

In order to configure this function, the auxiliary device should be first connected to the black interface provided on the control unit. Then press "One button" or "Two buttons" in the "Buttons" submenu. A list in these buttons indicates which auxiliary devices can be configure under "One button" or "Two buttons". A further submenu opens in which the function can be assigned to the button by pressing one of the symbols.

The selection now appears with a red margin. The procedure is then acknowledged with "OK"  .

The following functions can be assigned to the dispenser/stepper (see 7.2):

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*The assigned button function can be actuated from the external unit and the touch screen as soon as an external device is connected.*

### **10.8. Info**

*Serial and revision numbers for both the control unit and basis unit can be viewed in the “Info” submenu.*

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## **11. Interface**

*The unit can be operated from a PC via a serial or USB interface. A separate cable is available for this purpose which should be connected to the red interface at the rear of the control unit. The respective interface protocol can be obtained from the manufacturer.*

## **12. Cleaning**

*Your unit should be cleaned regularly to ensure a long service life.*

*Cleaning of the unit is imperative in the following cases:*

- If the piston unit has seized*
- After dispensing of reagents whose use is restricted (see 4 and 5)*
- Prior to changing reagent*
- Prior to sterilisation*

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- Prior to replacing suction / ejection or recirculation valves
- Prior to longer storage periods

### **Warning!**

- Put on protective clothing, safety glasses and protective gloves prior to commencing cleaning.
- Safety instructions (see 1) must be observed.
- The suction hose, ejection unit, glass cylinder and valves are filled with liquid! All openings should be turned away from the body / from persons.

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
## **12.1. Service/Cleaning**


- Terminate dispensing, fit the closure cap and deactivate the unit.
- Place the unit mounted on a storage bottle into a cleaning container.

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
- Using suitable gloves, unscrew the unit from the storage bottle and remove it (the suction hose no longer sticks into the liquid)

**Warning: Liquid may emerge from the suction hose!**

- Hold the ejection unit over the storage bottle opening.


- Pull off the closure cap and push to the rear into the retention position.

- Activate the unit, press the “Venting”  button repeatedly until no more liquid emerges.

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
## 12.2. Basic cleaning with the venting function

- Immerse the suction hose in a suitable cleaning solution.

Press the "Venting"  button to purge the unit completely.

### **Caution:**


- Residual liquid may cause splashing
- Hold away from body at all times

- Pull the suction hose out of the cleaning agent and purge the hose and unit with the "Venting" function until they are completely empty (until no more liquid emerges).
- Remove the following components from the unit:
- Unscrew the ejection unit from the unit.
- Pull off the recirculation hose.
- Release the union nut on the suction hose and pull off the suction hose.

Immerse these components in a suitable cleaning solution.

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### 12.3. Intensive cleaning

- Clean first as described in 12.1. and 12.2.
- Remove the valve block, cylinder and piston:
- Turn the upper part of the unit up to the right-hand stop (the point is aligned with the right-hand edge of the inscription on the chrome ring), then draw 6 mm upwards and turn to the left up to the left-hand stop. Now pull the valve block and cylinder downwards and away from the upper part of the unit.

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
#### **Note!**

The piston can be removed from the spindle if necessary. The spindle is held in position with a finger for this purpose, and the piston can be pushed to the side by hand to disconnect it from the spindle. For greater ease of access, the inspection window can also be removed if necessary.

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
If the basic unit is in the cleaning position, position mounting can be facilitated by moving the spindle slowly to the lower stop by pressing the "Move" button in the message on the control unit for a longer period of time. If the spindle is already in the lowermost position, it can be raised easily by pressing the button again for a longer period of time.

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- Check the piston for damage to the seal lips (replace if necessary) and deposit in the cleaning solution with the components listed in 12.2.
- Clean the cylinder with a brush and cleaning solution and rinse with a rinsing solution.
- Remount the piston, valve block and cylinder:
- Push the piston onto the spindle
- Fit the upper part of the unit onto the valve block and cylinder. Align the marking on the upper part of the unit with the left-hand edge of the marking on the chrome ring.

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- Push the upper part of the unit 6 mm onto the valve block and cylinder up to the first stop and turn to the right-hand stop, pushing a further 6 mm together subsequent to this.
- Now turn the ejection unit to the desired position.
- Check the function of the unit.

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**Warning!**

Insert the piston carefully into the glass cylinder, as the seal lip can otherwise be damaged. This will cause the system to leak!

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### 13. Sterilisation

The components in the unit conveying media are resistant to steam sterilisation up to 121 °C and an exposure time conforming to DIN te = 20 min.

**Warning:**

Only the following components should be steam sterilised:

- Valve block with cylinder
- Suction hose
- Ejection unit
- Recirculation hose
- Piston

The upper part of the unit, control unit and power supply unit should never be sterilised with steam!  
(See 2 for overview of individual components)

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### **13.1. Preparation**

- Cleaning should be realised prior to sterilisation (see 12).

### **13.2. Dismantling**


**Note:**

Dismantling can only be realised in the home position (the piston unit is in the lowest position). The unit should only be dismantled when deactivated!


- In addition to the components listed in 12.1., the other unit components conveying media described in 12.2. are dismantled and cleaned.

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

### 13.3. Loosening the suction/ejection and recirculation valves

- Disconnect the suction hose.
- Disconnect the ejection unit.
- Loosen the suction valve with the assembly wrench by turning it counterclockwise a half a revolution.

- Loosen the ejection valve with the assembly wrench by turning it counterclockwise a half a revolution.

- Loosen the recirculation valve with the assembly wrench by turning it counterclockwise half a revolution.

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
- Grasp the cover laterally and pull off.

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- Open the nut in a counterclockwise direction by turning it a half a revolution and, loosen the tap by applying pressure to the nut and dismantle the piston.

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
- Disconnect the piston from the spindle by pushing it to the side (see 12.3).

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**Important:**

The piston should be cleaned in a suitable solution prior to sterilisation. Then push the piston into the glass cylinder to a point where the piston protrudes slightly out of the glass cylinder.

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
### 13.4. Sterilisation procedure

Realise steam sterilisation of all components conveying media (up to 121 °C, exposure time conforming to DIN  $t_e = 20$  min).

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
#### **Warning!**

- Allow the unit to cool to room temperature and only then commence reassembly.
- Check all components for damage and replace if necessary.

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### 13.5. Assembly

- Remount the piston, valve block and cylinder:  
Fit the upper part of the unit onto the valve block and cylinder. Align the marking on the upper part of the unit with the left-hand edge of the marking on the chrome ring. Push the upper part of the unit 6 mm onto the valve block and cylinder up to the first stop and turn to the right-hand stop, pushing a further 6 mm together subsequent to this. Now turn the ejection unit to the desired position.

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



### **13.6. Tightening the suction/ejection and recirculation valve**

- Tighten the suction valve firmly with the assembly wrench by turning it in a clockwise direction.

- Tighten the ejection valve firmly with the assembly wrench by turning it in a clockwise direction.


- Tighten the recirculation valve firmly with the assembly wrench by turning it in a clockwise direction.

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- Tighten the nut on the valve tap by hand by turning it in a clockwise direction.

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### 13.7. Function inspection

- Activate the unit.
- Function inspection with the "Venting" symbol  . The venting function should be realised (can even be realised without liquid).
- Please check "Problem solutions" in case of malfunctions (see chapter 17) or contact the manufacturer.

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
## 14. Valve replacement

Leaking or damaged suction, ejection and recirculation valves should be replaced immediately (see accessories and spare parts in 19).



### 14.1. Suction valve replacement

- Realise cleaning (see 12).
- Unscrew the suction hose.
- Stand the unit on its head to prevent the valve ball falling out during replacement.
- Unscrew the suction valve with the assembly wrench and remove the valve ball.
- Insert the ball again and screw on the new suction valve, tightening with the assembly wrench.
- Push on the suction hose up to the stop and tighten the union nut firmly (see 6.1.).

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### **14.2. Ejection valve replacement**

- Realise cleaning (see 12).
- Unscrew the ejection unit.
- Screw out the ejection valve with the assembly wrench and remove.
- Screw in the new ejection valve and tighten firmly with the assembly wrench.
- Push on the ejection unit up to the stop and tighten the union nut firmly (see 13.6.). Check to ensure a secure fitting!

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### **14.3. Recirculation valve replacement**

- Unscrew the recirculation valve with the assembly wrench and remove.
- Screw in the new recirculation valve and tighten firmly with the assembly wrench.

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## **15. Repair by manufacturer**

- A defective unit should be cleaned carefully prior to sending for repair (see chapter 12).
- Units will only be repaired if they are cleaned by the user and are no longer contaminated with any substances which could pose a risk to health.
- Uncleaned units will be returned unrepaired to the sender.
- Describe the malfunction accurately.
- Specify the liquid used (important!)
- The guarantee is rendered null and void if the unit has been subject to interventions other than those described (e.g. application of force).

## **16. Volumetric inspection**

The parameters specified in chapter 18 for the regularity (R%) and coefficient of variation (CV%) of a unit can be checked as follows:

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### **1. Weighing procedure (see for example DIN EN ISO 8655-6)**

- Dispense the test volume with the unit (deionised water).
- Weigh the quantity dispensed on the analysis scale.
- Convert this quantity into volume.

**Note: consideration should be given to temperature and air pressure.**

### **2. Repetitive weighing**

- The weighing procedure described in 1 should be repeated at least five times to obtain reliable details on the regularity and coefficient of variation of a unit.

### **3. Calculating the regularity and coefficient of variation**

- The regularity (R%) and coefficient of variation (CV%) of the unit should be calculated on the basis of statistical quality control formulas:

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$$V = m \cdot Z$$

**Conversion of the mass determined to volume**

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$$R\% = \frac{100 \cdot \bar{V} - V_s}{V_0}$$

**Systematic measuring deviation (regularity R%)**

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$$s = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}}$$

**Random measuring deviation  
(coefficient of variation CV%)**

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### **Legend**

|     |   |  |
|-----|---|--|
| $m$ | = | weighed value determined in g  |
| $V$ | = | volume determined in mL  |
| $Z$ | = | correction factor<br>(e.g.: 1.0029 $\mu\text{l}/\text{mg}$ at 20 °C and 101.3 kPa 1,0040 $\mu\text{l}/\text{mg}$ at 25 °C and 101.3 kPa) |

|           |   |                       |
|-----------|---|-----------------------|
| $\bar{V}$ | = | mean volume value     |
| $n$       | = | number of measurement |
| $V_s$     | = | test volume           |
| $V_0$     | = | nominal volume        |
| $s$       | = | standard deviation    |

This brief guide merely illustrates the principal procedure and calculation for calibrating the unit. The precise procedure is described in detail in DIN EN ISO 8655.

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## **17. Problem solutions**

### **17.1. Messages**

The red flashing LED in the basic unit indicates that the unit is in a non-operational condition. The control unit also generates a suitable warning message:

| <i>Messages</i>                          | <i>Reason</i>                                    | <i>Solution</i>   |
|--|--|---|
| <i>Piston blockage, check basic unit</i> | <i>- Recirculation lever not applied to stop</i> | <i>- Turn the recirculation lever completely to the front or rear</i>     |
|  | <i>- Valves are blocked</i>                      | <i>- Clean or replace valves</i>  |
|  | <i>-- Liquid viscosity is too high</i>           | <i>- Select slower speed, use dispensing tip instead of titration tip</i> |
|  | <i>- Closure cap is not removed</i>              | <i>- Remove closure cap</i>   |

|  |   |   |
|--|---|---|
| <i>Unit is in cleaning position, valve block is not fitted</i> | <ul style="list-style-type: none"> <li>- Upper part of basic unit is turned to the right-hand stop (dismantling position)</li> <li>- Upper part of basic unit has beendismantled</li> </ul> | <ul style="list-style-type: none"> <li>- Turn the upper part slightly to the left</li> <li>- Assembly basic unit again</li> </ul> |
| <i>No basic unit identified</i>                                | <ul style="list-style-type: none"> <li>- Basic unit not connected to control unit</li> </ul>  | <ul style="list-style-type: none"> <li>- Check cable connection</li> </ul>  |
| <i>No piston fitted</i>  | <ul style="list-style-type: none"> <li>- Basic unit was assembled without piston</li> </ul>   | <ul style="list-style-type: none"> <li>- Dismantle basic unit and fit piston again</li> </ul>                                     |

*Special cleaning position case:*

*If the basic unit is in the cleaning position, pressing and holding the "Move" button in the message on the control unit for a longer period will move the spindle slowly down to the lower stop to facilitate piston fitting. If the spindle is already at the lowest point, pressing the button for a longer period will move it slightly upwards again.*

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---

## 17.2. Other problem solutions

If the unit can no longer be operated, it can be deactivated by pressing and holding (approx. 10 seconds) in the upper right-hand corner of the touch screen (battery symbol position) for a longer period. All changes realised since the last activation are lost as a result, so this function should only be used in emergencies.

| <i>Problem</i>  | <i>Reason</i>   | <i>Solution</i>   |
|---|---|---|
| <i>Liquid emerges from the gap between the suction / ejection valve / recirculation valve and valve block</i> | <i>- Valves loose or damaged</i>  | <i>- Tighten the valves firmly with the assembly wrench or replace if necessary. Warning: the suction / ejection / and recirculation valve should be tightened firmly with the assembly wrench after sterilisation.</i> |
| <i>Liquid emerges between the ejection valve and ejection unit</i>  | <i>- Ejection unit loose or damaged</i><br><i>- Incorrect ejection unit</i> | <i>- Mount original ejection unit correctly</i><br><i>- Use only original ejection unit.</i>  |
| <i>Unit is suctioning in air</i>  | <i>- Unit not vented</i><br><i>- Suction hose loose</i>                     | <i>- Vent unit</i><br><i>- Tighten suction hose firmly</i>  |

- Suction hose and/or union nut  
damaged

- Replace suction hose

---

Liquid emerges from recirculation valve  
(rear of unit)

- Ejection unit loose or damaged

- Tighten recirculation valve firmly with  
the assembly wrench or replace

---

Liquid emerges between the valve shaft  
and valve block (only visible when  
cover is removed)

- Valve shaft loose

- Tighten nut on valve shaft firmly

---

Erroneous dispensing volume

- Ejection unit loose or damaged

- Mount original ejection unit correctly.  
Use only original ejection unit.

- Suction hose loose or damaged

- Tighten suction hose firmly

- Suction / ejection / and recirculation  
valves loose or damaged

- Tighten the valves firmly with the  
assembly wrench or replace if  
necessary. Realise cleaning

---

|  |  |   |
|--|--|---|
| <i>Liquid over piston unit</i>               | - Damaged seal lip on piston unit                                  | - Replace piston unit   |
| <i>Visible air bubbles in glass cylinder</i> | - Unit not vented<br>- Suction hose or union nut loose, deformed   | - Vent unit<br>- Tighten suction hose firmly.<br>Change original suction hose |
| <i>Inaccurate dispensing volume</i>          | - Piston unit seal lip damaged<br>- Ejection unit loose or damaged | - Replace complete piston unit<br>- Mount ejection unit correctly             |
| <i>Erroneous values on display</i>           | - Unit error   | - Reset unit to default settings using the button in the "Main" menu item     |

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### **17.3. Contacting the manufacturer**

*Should the measures above fail to completely remedy the malfunction, the dealer or manufacturer should be contacted. Please specify the serial number and revision number of the basic unit and control unit. Serial numbers are embossed on unit components and can be viewed in the Info menu. Revision numbers can only be viewed in the Info menu.*

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### **18.1. Basic unit**

Volume range:

0.1 mL – 200 mL (maximum piston lift 20 mL, automatic addition), resolution: < 1  $\mu$ l.

1 mL – 500 mL (maximum piston lift 50 mL, automatic addition), resolution: < 2  $\mu$ l. Volume accuracy:

dispenser: 50 mL: R = 0.6 %, CV = 0.2 %

20 mL: R = 0.6 %, CV = 0.2 %

stepper: 50 mL: R = 0.5 %, CV = 0.25 %

20 mL: R = 0.5 %, CV = 0.3 %

\* Values specified for regularity (R) and coefficient of variation (CV) relative to maximum volume

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## **18.2. Control unit**

*Individual Dispensing:*

*volume of 0.1 mL – 200 mL (20 mL unit) and  
1 mL – 500 mL (50 mL unit).*

*Time-controlled serial dispensing:*

*volume of 0.1 mL – 200 mL (20 mL unit) and  
1 mL – 500 mL (50 mL unit).*

*Number of repetitions (1 to 9999)*

*Cycle time (1 s to 9999 s)*

*Stepper function: Dispensing of up to 10 partial volumes  
(can vary) from a single cylinder charge manually or via  
serial dispensing functions*

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*Conformity:*

*Conformity certificate conforming to German calibration regulations (Appendix 12).*

*CE:*

*This symbol confirms that the unit meets the stipulated requirements of pertinent EU directives and is tested in conformance with these.*

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## **20. Accessories - Spare parts**

|   |           |
|---|-----------|
| <i>System platform</i>                                | 958 60 01 |
| <i>Data power cable</i>                               | 958 70 01 |
| <i>Battery pack</i>                                   | 958 71 00 |
| <i>Touch screen film</i>                              | 958 72 00 |
| <i>Valve block</i>                                    | 958 60 20 |
| <i>Valve block</i>                                    | 958 60 50 |
| <i>Pedal switch</i>                                   | 956 40 02 |
| <i>Ejection unit, flexible, spiral shaped</i>         | 931 30 26 |
| <i>Connection cable for serial interface (9-pole)</i> | 956 40 05 |
| <i>Connection cable for USB</i>                       | 956 40 06 |

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## 21. Clearance certificate form

for opus® basic unit .....(SNR) verwendetes Medium: .....

for opus® operation unit "dispenser" .....(SNR)

|  | <b>Yes</b>               | <b>No</b>                | <b>Please define in case of YES</b> |
|--|--------------------------|--------------------------|-------------------------------------|
| Blood, bodily fluids, patholog. samples:                 | <input type="checkbox"/> | <input type="checkbox"/> | _____                               |
| Other hazardous biologic substances:                     | <input type="checkbox"/> | <input type="checkbox"/> | _____                               |
| Biologically-degradable materials which can pose a risk: | <input type="checkbox"/> | <input type="checkbox"/> | _____                               |
| Chemicals containing substances harmful to health:       | <input type="checkbox"/> | <input type="checkbox"/> | _____                               |
| Radioactive materials:                                   | <input type="checkbox"/> | <input type="checkbox"/> | _____                               |
| Other hazardous and toxic substances:                    | <input type="checkbox"/> | <input type="checkbox"/> | _____                               |

### **We hereby declare that**

- the unit has been carefully cleaned and decontaminated prior to shipment
- no risks resulting from contamination exist
- we are aware that liability claims can be made against us in the event of damage being caused by contaminated units.

### **Address of sender:**

Name: \_\_\_\_\_

Date: \_\_\_\_\_

Signature: \_\_\_\_\_