

## Comprehensive Genomic Profile Report

Powered by pierianDx

PierianDx  
77 Maryland Plaza  
St. Louis, MO 63108



PATIENT	DOB	DISEASE	MEDICAL RECORD #	REPORT DATE	REPORT STATUS
John Doe	02/04/1981	Non-small cell Lung Cancer	6563465346	02/18/2019	Final

### Report Summary

2 IA	0 IB	1 IIC	0 IID	High	Stable	13
GENOMIC FINDINGS BY TIER + LEVEL				TMB	MSI	TRIALS

#### GENOMIC FINDINGS

##### Tier I - Strong Clinical Significance

VARIANT	LEVEL	CLINICAL IMPACT
<b>NCOA4-RET</b> fusion	A	May benefit from - Cabozantinib, Vandatinib in <i>non-small cell lung cancer</i>
<b>KRAS</b> p.G12D c.35G>A	A	Not likely to benefit from - Erlotinib, Gefitinib, Afatinib, Osimertinib, Dacomitinib in <i>non-small cell lung cancer</i>  Unfavorable prognosis in - in non-small cell lung cancer

##### Tier II - Potential Clinical Significance

VARIANT	LEVEL	CLINICAL IMPACT
<b>PDGFRA</b> p.D842V c.2525A>T	C	May benefit from - Dasatinib in <i>gastrointestinal stromal tumor</i>  Not likely to benefit from - Sunitinib, Imatinib in <i>gastrointestinal stromal tumor</i>

##### Other Biomarkers

VARIANT	LEVEL	VALUE	CLINICAL IMPACT
<b>TMB</b>	high	24 mut/Mb	May benefit from - Nivolumab, Nivolumab + Ipilimumab in <i>non-small cell lung cancer</i>
<b>MSI</b>	stable	5% Unstable Sites	



# Report Optimization: Maximizing the Clinical Utility of Next Generation Sequencing

Eric Vail, MD

Director, Molecular Pathology  
Cedars-Sinai Medical Center



Rakesh Nagarajan, MD, PhD  
Executive Chairman and Founder  
PierianDx



Moderated by: OnLive

# Agenda

## 1 **The Evolving Landscape of Precision Medicine and Oncology Testing**

Rakesh Nagarajan, MD, PhD - PierianDx

## 2 **The Utility and Associated Challenges of Genomic Profiling**

Eric Vail, MD - Cedars-Sinai Medical Center

## 3 **Optimizing Clinical Utility & Value of Genomic Reporting**

Eric Vail, MD - Cedars-Sinai Medical Center

## 4 **Q&A**

Eric Vail, MD and Rakesh Nagarajan, MD, PhD

# Rapid Rise of Precision Medicine

## Numbers and Milestones

**284** Total # of FDA approved pharmacogenomic biomarkers; 108 in oncology<sup>1</sup>

**25** Personalized medicine approvals in 2018 (42% of NMEs)<sup>2</sup>

**10** Cancer-related personalized medicine approvals in 2018<sup>2</sup>


**2nd** Approval of a cancer drug (Vitrakvi) indication based on biomarker, not tumor type<sup>2</sup>

**3/16/18** CMS announces national coverage policy for NGS diagnostics

## Hospitals billing for comprehensive genomic profiling (CGP)<sup>3</sup>



## Most of the CGP volume is going to a few independent reference labs.<sup>3</sup>



Sources:

1. FDA - Table of Pharmacogenomic Biomarkers in Drug Labeling. Data as of June 2018

2. Personalized Medicine at FDA: A Progress & Outlook Report, Personalized Medicine Coalition, 2018.

3. Boston Healthcare Associates analysis of claims data 2016.

# Biomarker and NGS Testing

## Biomarker Testing

Progressive increase in **oncology** due to tumor agnostic biomarkers developed to inform targeted and immune therapies

**Standardized testing algorithms** drive biomarker testing for common tumors (NSCLC, CRC, breast) at diagnosis at many healthcare settings

## NGS Testing


More laboratories performing NGS for common tumor types with approved therapies

**Tumor agnostic** markers, such as MSI, TMB, and NTRK fusions are significant driver of more NGS testing

**Improving payor coverage** with recent Medicare coverage for FDA-approved NGS tests

## Approved

### Microsatellite Instability, NTRK Fusions




## Emerging

### Tumor Mutational Burden





# Oncology Testing is Evolving



# Overcoming the Challenges



Scarcity of genomic analysts



Growing number of assays from multiple vendors




NGS Test Reimbursement




Amount of data to curate and interpret from publications



Rapidly changing therapies, guidelines, clinical trials



# Wisdom in Every Report



**Today:** 45+ health systems and laboratories in partner sharing network

**2019:** PierianDx signs multi-year deal with Illumina to support cancer research and diagnostics

**2018:** Moffitt 1st to launch TruSight™ Tumor 170 clinically

**2014:** PierianDx established; Moffitt Cancer Center 1st to go live

**2011:** PierianDx technology developed at Washington University in St. Louis

# Technology Enabled Services

## Clinical Genomics Workspace

All-in-one informatics and reporting software




## Validation & Interpretation Services

Experienced team to **fast-track growth**



## Medically Powered Knowledgebase

The largest opt-in content **sharing network**



## Laboratory Services

**Turnkey**, validated assays and informatics

# Utility and Challenges Associated with Genomic Profiling

# Biomarker Detection



Traditionally for tumor characterization, targets were identified by limited assays

- IHC
- FISH
- Allele specific PCR


Highly specific but provide limited information

Requires multiple tests for complete analysis


- Increase in markers increases costs and tissue requirements




No Mutation Detected



V600E Mutation Detected



# Predictive Molecular Biomarkers in Oncology




- Allows for large amounts of sequencing relatively cheaply and quickly

## Commercially available for over a decade

- Significantly improved in that time
- Currently a very reliable technology

## Technical cost now less than analytical

- Depending on panel size can have a very large amount of variants to report



Hotspot Genes (87)				Full-Length Genes (48)		Copy Number Genes (43)		Gene Fusions (Inter- and Intra-genic)	
AKT1	FGFR8	MDM2	SAC	ARID1A	PALB2	AKT1	MDM2	AKT2	NF1
AKT2	FLT3	MEK	STAT3	ATM	PIK3R1	AKT2	MDM4	AKK	NOTCH1
AKT3	FGFR2	MTOR	TEK	ATR	PMS2	AKT3	MEK1	AK	NOTCH4
ALK	GATA2	MYC	TOP2	ATRX	POU5	ALK	MYC	AKL	NRG1
AR	BAP1	MYCN	U2AF1	BAP1	PTCH1	AR	MYCL	BRAP	NRK1
ARAF	GNAQ	MYD88	XPO1	BRCA1	PTEN	ARL	MYCN	BRCA1	NTN2
ATL	GNA3	NF2Y2		BRCA2	RAD50	BRAP	NTN1	BRCA2	NTN3
BLAP	HSP3A	NRAS		CDK12	CDK5	CDK12	NTN2	CDKN2A	NUM1
BTX	HIST1H18	NTN1		CDKN18	NDS18	CDK2D	NTN3	EGFR	PDGFR
CB1	HNF1A	NTN2		CDKN1A	NDS18C	CDK2D	PDGFR	ERBB2	PDGFRB
CDK6	HNR4S	NTN3		CDKN2B	NDS19	CDK2E1	PDGFRB	ERBB4	PIK3CA
CDK8	IDH1	PDGFR		CHUK1	NR1	CDK2	PIK3CA	ERG	PRKG
CDK9	IDH2	PDGFRB		CHUKBP	NR3A3	CDK4	PIK3CB	ESR1	PRKAC
CHUK2	JAK2	PIK3CB		CHUK2	STT2	CDK6	PRKAC	ETV3	PRKACB
CSF1R	JAK2	PIK3CB		PANCO2	SLK	CDK9	NFICTOR	ETV4	PTEN
CTNNB1	JAK3	PPP2R1A		PANCI	SMARCA4	ERBB2	TEK	ETV5	NDS18
DCK2	KDR	PTEN1		PRKX1	SMARCB1	ESR1		FGFR1	RAF1
EGFR	KIT	RAC1		MLH1	STK11	FGFR18		FGFR2	NR1
ERBB2	KIT2RN	RAF1		NF1B1A	TP53	FGFR3		FGFR3	NR1A
ERBB3	KRAS	RET		TP53	TSC1	FGFR		FGFR	RET
ERBB4	MAP2K1	RHE		MSH6	TSC2	FGFR2		FLT3	ROS1
ERCC2	MAP2K1	RHOA		NRN		FGFR3		JAK2	RSP02
ESR1	MAP2K2	ROS1		NP1		FGFR4		KRAS	RSP03
EGH	MAP2K4	SPB1		NP2		FLT3		MDM4	TEK
FGFR1	MAPK1	SHAD4		NOTCH1		FGFR1			
FGFR2	MAP2	SHAD		NOTCH2		KIT		MYB	
FGFR3	MDM4	SPOR		NOTCH3				MYBL1	




# Therapy Selection

15 genes and MSI are FDA approved predictive biomarkers

- *ALK, BRAF, BRCA1/2, EGFR, ERBB2, FGFR2/3, KRAS, NRAS, NTRK1/2/3, PIK3CA, ROS1*
- MSI and *NTRK* pan-tumor, remainder tissue specific

Multiple others have NCCN emerging therapy designation and/or high levels of clinical evidence to support off-label use

- *ATM, MET, PALB2, RET, TMB*



# Additional Utility of Genomic Profiling

## Clinical Trial Recruitment

- Currently ~5% of patients are added to a clinical trial

## Germline variant detection

- Impact for entire families

## Diagnostic precision

- Many tumors defined by their molecular event
- Sub-classification of tumors also increasingly driven by genomic profile
- Heme leading but solid tumors will follow

## Tissue of origin (primary vs. metastasis)



National Comprehensive  
Cancer Network®

**“NCCN believes that best management of any patient with cancer is in a clinical trial. Participation in clinical trials is especially encouraged.”**

*--NCCN Guidelines, 2014*

# Challenges with Genomic Reporting



Pathologists



Clinicians



Patients



Regulators

## Competing audiences

- Different educational backgrounds even within groups
- Different “wants” (sometimes opposing)

## Lack of standardization

- Genomic nomenclature can be tricky especially with older “pre-NGS” genes/variants
- Reporting of genomic changes, VAF, CN, fusion transcripts differ lab to lab
- What level of evidence should be used for clinical recommendations?

## Report findings vs. Analyze findings

- Most labs at least report some analysis of variants
- Matched therapies, clinical trials, guidelines and/or global analysis all can be added on

# NTRK Fusion Reporting

*NTRK* fusions can be reported many different ways

- Note hyphen between fusion partner and gene
- In this example, fusion is referred to as “rearrangement”
- Every laboratory reports *NTRK* fusions slightly differently

## GENOMIC VARIANTS

### Somatic - Potentially Actionable

### Variant Allele Fraction

**TP53** p.Y220C Missense variant - LOF

90,3% 

**FGFR1** Copy number gain

**EIF4B - NTRK3** Chromosomal rearrangement

### Somatic - Biologically Relevant

**KAT6A** Copy number gain

### Germline - Pathogenic / Likely Pathogenic

No pathogenic germline variants were found in the limited set of genes we report on.

## IMMUNOTHERAPY MARKERS

### Tumor Mutational Burden

### Microsatellite Instability Status

**5.0 m/MB** 81st percentile

**Stable** Equivocal High

## TREATMENT IMPLICATIONS

**Larotrectinib (TRK Inhibitor)**  
 ✓ FDA approved, current diagnosis

EIF4B - NTRK3 Chromosomal rearrangement  
 Consensus, solid tumors: [FDA](#)

# NTRK Fusion Reporting

## High Impact Results

BIOMARKER	METHOD	RESULT	THERAPY ASSOCIATION		BIOMARKER LEVEL*
KRAS	NGS	Mutation Not Detected	BENEFIT	cetuximab, panitumumab	Level 1
NRAS	NGS	Mutation Not Detected			
NTRK1	RNA-Seq	Fusion Detected	BENEFIT	larotrectinib	Level 1
BRAF	NGS	Mutation Not Detected			

\* Biomarker reporting classification: Level 1 - highest level of clinical evidence and/or is endorsed by standard clinical guidelines; Level 3 - potential clinical significance


- Larotrectinib is included in the therapy association
- Prominent use of the word “fusion”
- Note hyphen between fusion partner and gene

## Important Note

A TPM3-NTRK1 gene fusion was detected in this tumor. This fusion has been previously reported in several different tumors (Ardini 2014 Mol Oncol 8:1495; Chiang 2018 Am J Surg Pathol 42:791).

# NTRK Fusion Reporting

- Does not include the fusion partner
- Larotrectinib is listed as therapy
- Prominent use of the word "fusion"

THERAPY CONSIDERATIONS FOR COLORECTAL CANCER			
	Markers Identified	Therapies in Colorectal Cancer	Therapies in Other Tumor Types
Level 1	 NTRK1 Fusion	larotrectinib <sup>1</sup>	No therapies in other tumor types
Level 2	No FDA Level 2 markers with immunotherapy or targeted therapy considerations for the tumor type tested or other tumor types were identified		
Level 3	No FDA Level 2 markers with immunotherapy or targeted therapy considerations for the tumor type tested or other tumor types were identified		
FDA Evidence Levels: 1) Companion diagnostic; 2) Practice guidelines, clinically validated; 3) Clinically significant, analytically validated with clinical or mechanistic rationale (clinical trials, off-label therapies, or peer reviewed evidence)			

Clinical Trial Markers Identified	
Immunotherapy	Targeted Therapy
Level 3 No immunotherapy clinical trials identified	NTRK1 Fusion
See clinical trials page 2	

Therapeutic Markers NOT Identified	
Immunotherapy	Targeted Therapy
CD8 (IHC) Minimally Infiltrating CD8 (RNA-Seq) 18% (Noninflamed) MSI Stable PD-L1 (IHC 22C3) 0% TPS PD-L1 (RNA-Seq) 42% Rank TMB 1.8/Mb (Low)	BRAF V600E KRAS mutation NRAS mutation

PATHOLOGIST SUMMARY INTERPRETATION
<p><b>MOLECULAR SUMMARY:</b> This NTRK fusion positive, microsatellite stable colorectal cancer with a very low mutational burden is non-inflamed with a moderately low number of CD8+ T-cells in a minimally infiltrating pattern, a high number of FOXP3+CD4+ T-cells, and negative expression of PD-L1 by IHC (TPS=0%; 22C3 clone). No dominant immunosuppressive mechanism is identified in this tumor and the overall assessment is a somewhat non-immunogenic tumor. The APC c.3958_3962dup (S1321RfsX2) mutation identified in this tumor is an uncommon variant in this gene.</p> <p><b>LIKELIHOOD OF RESPONSE BASED ON EVIDENCE IN CURRENT LITERATURE:</b> From a targeted therapy perspective, the NTRK1 fusion identified in this case is level 1 clinical evidence for treatment with a TRK inhibitor, such as larotrectinib or entrectinib. TRK inhibitor response rates for NTRK fusion positive tumors are &gt;75% regardless of tumor histology. For a review of this topic, see the recent publication by Cocco et al. (PMID: 30333516)</p> <p>From an immunotherapeutic perspective, response to PD-1 axis checkpoint blockade in this case is not favorable. Recommendation is augmentation of PD-1 axis checkpoint blockade with chemotherapy or radiation (see clinical trials section of the report).</p>

# NTRK Fusion Reporting

## Companion Diagnostic (CDx) Associated Findings

GENOMIC FINDINGS DETECTED	FDA-APPROVED THERAPEUTIC OPTIONS
<b>KRAS</b> wildtype (codons 12 & 13)	Erbitux® (Cetuximab)
<b>KRAS/NRAS</b> wildtype (codons 12, 13, 59, 61, 117, & 146 in exons 2, 3, & 4)	Vectibix® (Panitumumab)

- Note hyphen between fusion partner and gene
- Larotrectinib is not listed as therapy on this page
- Prominent use of the word “fusion”

### OTHER ALTERATIONS & BIOMARKERS IDENTIFIED

Results reported in this section are not prescriptive or conclusive for labeled use of any specific therapeutic product. See *professional services* section for additional information.

**Microsatellite status** MSI-High §  
**Tumor Mutational Burden** 34 Muts/Mb §  
**ASXL1** G645fs\*58  
**CIC** G136fs\*8  
**CREBBP** H2384fs\*12  
**FLCN** H429fs\*39

**NTRK1** TPM3(NM\_152263)-NTRK1(NM\_002529) fusion (T10\*; N10) §  
**RNF43** R286W  
**RNF43** G659fs\*41  
**SOX9** Q357\*  
**TSC2** G654fs\*2

§ Refer to appendix for limitation statements related to detection of any copy number alterations, gene rearrangements, MSI or TMB result in this section.

Please refer to appendix for Explanation of Clinical Significance Classification and for variants of unknown significance (VUS).



# NTRK Fusion Reporting

- Note hyphen between fusion partner and gene
- Larotrectinib is listed as therapy with clinical benefit
- Prominent use of the word “fusion”

BIOMARKER FINDINGS		THERAPIES WITH CLINICAL BENEFIT (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL BENEFIT (IN OTHER TUMOR TYPE)
<b>Microsatellite status - MSI-High</b>		Nivolumab	none
10 Trials see p. 15		Pembrolizumab	
<b>Tumor Mutational Burden - TMB-High (34 Muts/Mb)</b>		none	Atezolizumab
			Avelumab
			Durvalumab
			Nivolumab
			Pembrolizumab
10 Trials see p. 17			
GENOMIC FINDINGS		THERAPIES WITH CLINICAL BENEFIT (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL BENEFIT (IN OTHER TUMOR TYPE)
<b>NTRK1 - TPM3-NTRK1 fusion</b>		Larotrectinib	Crizotinib
8 Trials see p. 19			
<b>RNF43 - G659fs*41, R286W</b>		none	none
2 Trials see p. 21			
<b>TSC2 - G654fs*2</b>		none	none
10 Trials see p. 22			



# NTRK Fusion Reporting


## NTRK1 Point Mutation vs. Fusion

- Do not confuse NTRK1 point mutation with "fusion."
- Point mutation is not part of the indication for larotrectinib.
- Anecdotal findings suggest that larotrectinib does not work in these patients.


## Variants of Unknown Significance

Gene	Variant
<i>NOTCH3</i>	P.R544C chr19:g.15298126G>A NM_000435.2:c.1630C>T
<i>NTRK1</i>	p.Q80R chr1:g.156834172A>G NM_002529.3:c.239A>G
<i>TSC2</i>	p.P330L chr16:g.2110684C>T NM_000548.3:c.989C>T
<i>ARID1A</i>	p.G1254S chr1:g.27099881G>A NM_006015.4:c.3760G>A


# Why is this Important?



**B** NCCN-directed therapy



**C** EGFR inhibitor receipt status



# How can Pathology Help?

1

## Comprehensive analysis and reporting

- Utilize EMR in crafting reports for true personalization
- Provide updates on novel therapeutics
- Hand match clinical trials in and outside network
- Explain discordant results

2

## Outreach and education

- Both clinicians and general pathologists
- Builds personal relationships and trust in the laboratory
- Generally enhances patient care

# Structure of Cedars-Sinai Reports

First page has everything really important for patient care

Summary of results with therapies and guidelines

## Molecular pathology interpretation

- Provides personalized, in-depth analysis
- Lists relevant therapies and guidelines
- Suggests relevant clinical trials

Remaining (3-5) pages has everything else

- Clinical trial details
- Variant details
- VUS
- Regulatory requirements (disclaimers, methodology etc.)

CS - Comprehensive Cancer Panel					
SPECIMEN INFORMATION					
Specimen Type:	Tissue specimen from bronchus			Source:	RIGHT MAIN STEM BRONCHUS
Accession Number:	[REDACTED]			Block:	A1
Percent Tumor Nuclei:	35				
Clinical Indication:	Bronchial tissue diffusely invaded by poorly differentiated, mucinous adenocarcinoma (lung primary)				
RESULT SUMMARY					
Variant Detected	FDA Approved Therapy Within Indication	FDA Approved Therapy Outside Indication	Resistance to Therapies	Guidelines	Clinical Trial Opportunity
<i>EML4, ALK</i> Fusion	Brigatinib, Crizotinib, Alectinib, Lorlatinib, Ceritinib	No	No	Yes - See Variant Details	Yes - see below
<i>TP53</i> p.G245S 16.2%	No	No	No	No	Yes - see below
IMMUNOTHERAPY					
TMB Low: 2 muts/Mb					
Tumor mutational burden (TMB) is a measurement of the amount of nonsynonymous somatic mutations present within a tumor sample. Tumors that have low TMB rarely respond to immunotherapy (PMID 28835386, 29657128, 30643254).					
MOLECULAR PATHOLOGY INTERPRETATION					
This [REDACTED] with bronchial tissue diffusely invaded by poorly differentiated, mucinous adenocarcinoma (lung primary). FISH testing was negative for rearrangements in <i>ALK</i> and <i>ROS</i> and amplification of <i>MET</i> . The information provided in this report relates to the currently known clinical significance and possible therapeutic implications of mutations detected in the above specimen from the right main stem bronchus.					
An <i>EML4/ALK</i> fusion and an inactivating <i>TP53</i> mutation were detected. This fusion is known to be responsive to TKI therapy. In the treatment naive setting the NCCN recommends alectinib (preferred), crizotinib, brigatinib, and ceritinib ( <i>ALK</i> inhibitors, FDA approved for <i>ALK</i> -positive metastatic NSCLC) for patients with <i>ALK</i> -positive disease (category 1) (NCCN NSCLC v3.2019). Additionally, RNA based detection of fusions is considered to be significantly more sensitive than FISH and would explain the discrepancy in this case (PMID: 27769042). Of note, in <i>ALK</i> -rearranged/ <i>TP53</i> co-mutated NSCLC patients who received systemic therapy, median progression-free survival and overall survival were significantly lower compared with <i>TP53</i> wild-type patients (PMID: 30165392). Clinical trials are available for both of these mutations.					
No other clinically significant mutations were detected (see gene list below).					

# Explanation of Rare Variants

## Sequencing produces many rare variants

- Clinicians are often unprepared to evaluate them
- Need for descriptive reporting and direct consultation
- Often have molecular studies that seem discrepant
- EMR is crucial in these cases

## Results Summary

Variant Detected	FDA Approved Therapy Within Indication	FDA Approved Therapy Outside Indication	Resistance to Therapies	Guidelines	Clinical Trial Opportunity
<i>EGFR</i> Exon 18 to 25 Kinase Domain Duplication	Osimertinib, Erlotinib, Afatinib, Gefitinib	No	No	Yes - see variant details below	Yes - see below

## 49 y/o m with metastatic (bone & brain) NSCLC

- *EGFR* (Sanger), *BRAF* (PCR) and *PD-L1* (IHC) negative at OSH
- Remaining biomarkers QNS

## Came to Cedars-Sinai for treatment

- Repeated biopsy to complete molecular profile

## Found to have an *EGFR* KDD (exon 18-25)

- Rare activating *EGFR* variant (~1% of *EGFR* mut)
- Sparse clinical evidence suggests sensitivity to TKI therapy (particularly afatinib)

## Discussed findings with clinician prior to and during tumor board

- Explained discrepancy with prior *EGFR* study
- Osimertinib ultimately chosen (Brain penetration)

**Patient started therapy and has had complete symptom relief and radiographic shrinkage of the brain and bone mets**

# Updating Targeted Therapies

## Molecular therapeutics is advancing rapidly

- Novel therapies are continually being approved
- Often therapies that work on a target in one tumor type work in another
- No patient should miss out if available

73 y/o female with GBM diagnosed August of last year

## Found to have an FGFR3-TACC3 fusion

- At the time no targeted therapies, just trials

## Received standard of care (TMZ, radiation)

- Developed aplastic anemia and was no longer eligible for clinical trial
- Considering hospice

## Erdafitinib (FGFR TKI) approved recently for FGFR2/3 mutated urothelial carcinoma

- Queried database for patients with those mutations (pan-cancer)
- Reached out to clinician and amended report to include this as off label therapy
- Patient was able to recently obtain and begin treatment

Variant Detected	FDA Approved Therapy Within Indication	FDA Approved Therapy Outside Indication	Resistance to Therapies	Guidelines	Clinical Trial Opportunity
<b>FGFR3, TACC3</b> Fusion	No	No	No	No	Yes - see below
<b>TP53</b> p.G245C 49.4%	No	No	No	No	Yes - see below
<b>CDK4</b> Amplification	No	No	No	No	Yes - see below
<b>FGFR3</b> Amplification	No	No	No	No	Yes - see below
<b>CCND2</b> Amplification	No	No	No	No	Yes - see below

# Clinical Trial Matching

Most patients have a variant with an open related clinical trial

Including clinical trials on the report is an area of disagreement between clinicians

- Some utilize it, some skip right over
- How can we make higher yield?

Curating them takes a lot of work

- Easily searchable/automatable databases help
- Still large manual portion

Will include trials in the first page interpretation when especially clinically relevant

- In system or local area a big plus
- Good published clinical data
- Example: EGFR/ERBB2 exon 20 insertions referred to trial treating with poziotinib at an outside institution.

## POTENTIAL CLINICAL TRIALS

Clinical Trials associated with this patient's genomic profile and tumor type are displayed below.

TITLE	TRIAL IDENTIFIER	PHASE	VARIANT
Randomized Blinded Phase III Assessment of Second or Third-Line Chemotherapy With Docetaxel + Plinabulin Compared to Docetaxel + Placebo in Patients With Advanced Non-Small Cell Lung Cancer and With at Least One Measurable Lung Lesion	<a href="https://clinicaltrials.gov/show/NCT02504489">NCT02504489</a> <a href="https://clinicaltrials.gov/show/NCT02504489">https://clinicaltrials.gov/show/NCT02504489</a>	III	<b>ALK, TRMT61B</b> ALK-TRMT61B fusion transcript
Randomized Blinded Phase III Assessment of Second or Third-Line Chemotherapy With Docetaxel + Plinabulin Compared to Docetaxel + Placebo in Patients With Advanced Non-Small Cell Lung Cancer and With at Least One Measurable Lung Lesion	<a href="https://clinicaltrials.gov/show/NCT02504489">NCT02504489</a> <a href="https://clinicaltrials.gov/show/NCT02504489">https://clinicaltrials.gov/show/NCT02504489</a>	III	<b>EML4, ALK</b> EML4-ALK fusion transcript
Targeted Therapy in Children and Adolescents With Recurrent, Progressive and Unresectable Inflammatory Myofibroblastic Tumor With the Inhibitor of Tyrosine Kinase -Crizotinib	<a href="https://clinicaltrials.gov/show/NCT03874273">NCT03874273</a> <a href="https://clinicaltrials.gov/show/NCT03874273">https://clinicaltrials.gov/show/NCT03874273</a>	II/III	<b>ALK, TRMT61B</b> ALK-TRMT61B fusion transcript
Targeted Therapy in Children and Adolescents With Recurrent, Progressive and Unresectable Inflammatory Myofibroblastic Tumor With the Inhibitor of Tyrosine Kinase -Crizotinib	<a href="https://clinicaltrials.gov/show/NCT03874273">NCT03874273</a> <a href="https://clinicaltrials.gov/show/NCT03874273">https://clinicaltrials.gov/show/NCT03874273</a>	II/III	<b>EML4, ALK</b> EML4-ALK fusion transcript
PHASE I SAFETY, PHARMACOKINETIC AND PHARMACODYNAMIC STUDY OF PF-02341066, A C-MET/ HGFR SELECTIVE TYROSINE KINASE INHIBITOR, ADMINISTERED ORALLY TO PATIENTS WITH ADVANCED CANCER	<a href="https://clinicaltrials.gov/show/NCT00585195">NCT00585195</a> <a href="https://clinicaltrials.gov/show/NCT00585195">https://clinicaltrials.gov/show/NCT00585195</a>	I	<b>EML4, ALK</b> EML4-ALK fusion transcript

# Tumor Mutation Burden (TMB)

Gold standard is T/N exome


- Mostly done clinically with tumor only panels

Currently reported as a category and a number

- High -  $\geq 10$  mutations/Mb
- Intermediate - 5-9 mutations/Mb
- Low -  $\leq 4$  mutations/Mb

“One size fits all approach”

Does this really translate from tumor to tumor?






## IMMUNOTHERAPY

TMB Low: 2 muts/Mb

Tumor mutational burden (TMB) is a measurement of the amount of nonsynonymous somatic mutations present within a tumor sample. Tumors that have low TMB rarely respond to immunotherapy (PMID 28835386, 29657128, 30643254)



# TMB Reporting




# Visibility of Reporting

## Reports are often PDFs

- Not discrete (unsearchable)
- Disconnected from the histopathology
- Easily lost in the shuffle of the dreaded media tab

May lead to situations where molecular results are never reviewed!

## Multiple possible solutions



Direct discrete data  
integration into EMR



Molecular  
consult notes



EMR  
alerts

# Summary


- Genomic profiling provides **high value for patients and laboratories**
- Molecular oncology reporting is **highly complex**
- Differences in reporting have real-world clinical impact. **Standardization and integration with EMR can help mitigate challenges.**
- When in doubt, **pick up the phone!**



# Questions and Answers

## Eric Vail, MD

Director, Molecular Pathology  
Cedars-Sinai Medical Center



## Rakesh Nagarajan, MD, PhD

Executive Chairman and Founder  
PierianDx

