

# 3pp® Wrist P.O.P.™ (Point of Pressure) Splint

## Application Sheet



1.

Apply base sleeve



2.

Apply pads to pressure strap. Use either 1/8" pads, 1/4" pads or one of each



3.

Apply strap with pad pressing down on the ulnar styloid



4.

Pull the strap around so the second pad pushes up under the radius

Watch the Wrist P.O.P. video on:  
[www.3pointproducts.com](http://www.3pointproducts.com)  
[www.ohmyarthritis.com](http://www.ohmyarthritis.com)



Wear with the Design Line™ cover

©3-Point Products, Inc.