

## On-Resin Disulfide Bridge Formation

### INTRODUCTION

Many naturally occurring peptides contain intra-disulfide bridges, which play an important role in biological activities. There are many ways to form a disulfide bridge in the solution phase and solid phase. Formation of a disulfide bridge in the solution phase is well known and widely used in the peptide community. However, a method for on-resin disulfide bridge formation has also been developed and would be a great tool for library and screening purposes<sup>1,2</sup>. In this application note, on-resin disulfide bridge formation is demonstrated. The commercially available peptides human amylin (1-13) and oxytocin were synthesized by oxidative cyclization with Thallium(III)<sup>3</sup> or Iodine<sup>4</sup> prior to cleavage.

Amylin (1-13), Human: H-KCNTATCATQRLA-OH, Disulfide Bridge, Cys<sup>2</sup>-Cys<sup>7</sup>

Oxytocin: H-CYIQNCPLG-NH<sub>2</sub>, Disulfide Bridge, Cys<sup>1</sup>-Cys<sup>6</sup>

For on-resin disulfide bridge formation, the Ac<sub>m</sub> protecting group was used to protect the cysteine side chain during the synthesis. The Ac<sub>m</sub> protecting group is stable to TFA, but is removed oxidatively with TI(III) or I<sub>2</sub> during disulfide formation. Linear and cyclic amylin and oxytocin were prepared and analyzed by reverse phase HPLC and mass spectrometry.

### METHOD

**Ac<sub>m</sub>-protected Linear Peptide Synthesis:** linear human amylin (1-13), (H-KC\*NTATC\*ATQRLA-OH)<sup>5</sup> and oxytocin (H-C\*YIQNC\*PLG-NH<sub>2</sub>)<sup>5</sup> were synthesized on a *Prelude*<sup>TM</sup> peptide synthesizer under the following conditions:

<sup>1</sup> S. Gazal, G. Gellerman, E. Glukhov, C. Gilon. *J. Peptide Res.*, **58**, 527 (2001).

<sup>2</sup> M. C. Munson and G. Barany. *J. Am. Chem. Soc.*, **115**, 10203 (1993).

<sup>3</sup> N. Fujii, A. Otaka, S. Funakoshi, K. Bessho, T. Watanabe, K. Akaji, and H. Yajima. *Chem. Pharm. Bull.*, **35**, 2339 (1987).

<sup>4</sup> B. Kamber, A. Hartmann, K. Eisler, B. Riniker, H. Rink, P. Sieber, and W. Rittel. *Helv. Chim. Acta*, **63**, 899 (1980).

<sup>5</sup> C\* represent Ac<sub>m</sub> protected Cys.

Scale: 40 μmol; Resin: Fmoc-Rink-MBHA (0.47 mmol/g); Deprotection: 20% piperidine in DMF, 3 min then 20 min; Coupling: 1:1:2 AA/HCTU/NMM in DMF, 2 x 45 min. Cleavage: 92.5:2.5:2.5:2.5 TFA/EDT/H<sub>2</sub>O/TIS, 2 hours.

**Cyclic Peptide Synthesis:** Peptides were synthesized the same as the linear peptides, except the disulfide bridge was formed prior to cleavage.

- **Disulfide Bridge Formation with Iodine:** Treat resin with I<sub>2</sub> (10 eq.) in DMF/H<sub>2</sub>O (4:1) for 40 min. Wash resin with DMF x 2, 2% Ascorbic acid in DMF x 2, DMF x 5 and DCM x 2.
- **Disulfide Bridge Formation with TI(CF<sub>3</sub>CO<sub>2</sub>)<sub>3</sub>:** Treat resin with TI(CF<sub>3</sub>CO<sub>2</sub>)<sub>3</sub> (1.2 eq.) in DMF for 40 min x 2. Wash resin with DMF x 6 and DCM x 6.


**Analysis:** Peptides were analyzed on a Varian Microsorb C-18 column (4.5 x 250 mm) on a Varian Pro-Star HPLC using an aqueous acetonitrile, 0.1% TFA buffer system with an increasing gradient of 5-60% acetonitrile over 55 minutes. Detection was at 214 nm. Mass analysis was performed using a Perseptive Biosystems MALDI-TOF mass spectrometer.

### RESULTS/DISCUSSION

Mass spectrometry confirmed the successful synthesis of both linear and cyclic peptides. Ac<sub>m</sub>-protected linear human amylin (1-13) and oxytocin results are shown in Table 1, while cyclic peptide results are shown in Table 2.

**Table 1:** Mass spectrometry analysis of Ac<sub>m</sub>-protected linear peptides.

| Peptide  | Expected m/z | Observed m/z            |
|----------|--------------|-------------------------|
| Amylin   | 1522         | 1547 [Na <sup>+</sup> ] |
| Oxytocin | 1154         | 1176 [Na <sup>+</sup> ] |



**Figure 1:** Linear peptide: (a) amylin, (b) oxytocin. Peptide cyclized with iodine: (c) amylin, (d) oxytocin. Peptide cyclized with thallium: (e) amylin, (f) oxytocin.

**Table 2:** Mass spectrometry analysis of cyclic peptides.

| Peptide  | Expected m/z | Observed m/z              |
|----------|--------------|---------------------------|
| Amylin   | 1378         | 1378.9                    |
| Oxytocin | 1007         | 1029.6 [Na <sup>+</sup> ] |

HPLC results are shown in Figure 1. From these results, it is clear that cyclization with thallium produced a purer product than iodine for both peptides.

**CONCLUSION**

This application note demonstrates that on-resin disulfide bridge formation using Thallium or Iodine would be good tool for synthesizing a disulfide bridge-containing library.