

Multiplex Peptide Synthesizer



www.gyrosproteintechnologies.com

GYROS PRCTEIN Technologies

Symphony X[®] Peptide Synthesizer

Ultra PurePep™ Pathway "inside": The proprietary ultra PurePep Pathway is enhanced with a dedicated microfluidic matrix valve block pathway for each reaction vessel resulting in pure peptides.

- Maximum Throughput and Production
- Fast, Efficient Couplings
- High Quality Crude Peptides
- Extremely Flexible
- Automated Cleavage

8 SOLVENT POSITIONS

- Maximize flexibility
- 20 L solvent container (other sizes available upon request)
- Automated cleaning routines
- No calibration needed for viscosity
- Sensor driven constant volume deliveries

12 INDEPENDENT REACTION VESSEL PAIRS

- Up to 24 RVs with in situ activation
- Optional IR heating and/or UV monitoring in RV position 1
- Up to 12 RVs with pre-activation
- Adjustable nitrogen mixing
- Independent operation
- Top-down washing

IR HEATING

- Oscillation mixing ensures no hot spots
- Test the same protocol with and without heat
- Optimize process variables
- Precise temperature control



28 and 12 Additional Amino Acid Positions

- 10, 120 and 400 mL sizes available
- Pre-packed amino acid bottles available
- Single-Shot[™] deliveries
- Precise optical sensor delivery measurements
- Automated cleaning routines

REAL TIME UV MONITORING

- Monitors reaction in vessel and not the waste stream
- Extends deprotection time when needed
- Synthesize the right peptide the first time
- Increases yield and purity
- Reduces solvent consumption

High Quality Performance



Twenty-four 42-mers in 14 hc	ours
Deprotection Time:	2 min
Coupling time:	5 min
Total Synthesis Time, 24 Peptides:	14 hours

Twenty-four 68-mers in 25 hours					
Deprotection Time:	4 min				
Coupling time:	5 min				
Total Synthesis Time, 24 Peptides:	25 hours				

Innovative Features

ULTRA PUREPEP PATHWAY



This proprietary matrix valve block system was specifically designed to minimize dead volumes, eliminating cross contamination and delivering peptides with high crude purity.

AUTOMATED CLEAVAGE

- Customizable cleavage programs
- TFA resistant vessels
- Automated cleaning post cleavage
- 50 mL collection tube capacity
- Adjustable start times

	Step 1. Select the Amins Acid file and the Solvent file for th Amine Acid file: 94				PIC-RIA					alvent the		-	2.0					
Aniesi Acie Mer			44.44						and the local diversion of the local diversio		44.55	*						
Step 2: Select program(s). Irre-Synthesis Program				ogum gaswell.prg							Sele	ect file						
Delault Program					ga landem.prg							Select Ne						
Post dynthesis Program				m	gaend prg							Select Ne						
itep 3.	Select a	sequence.																
					Seques					hth.seq				Sele	ect file			
							ter leng	th seque			the N	program step						
a a a a a a a a a a a a a a a a a a a	ndem M		8		Legues				ect	1.500				Sele	ext Ne			
redness															Length	Clemina A		
PV	7 5	1 1 3 2	¥ 3 4	6											18	COOM	94.12	
n/ 6	6 H W														34		1140.3	
	4	K	<				c	pole: 1				>		>>		CONPU	1140.	
tep 5	Optional	Change t		m for the		nce cyde.												
	min	- 0			×						٠							
															Select Re			
							Propert	·							Select Ne			
	2		•		•					60,00	andem p	-						
8	New File		Op	en file		De	iete file					Save Ar	file			140		
	Print			Prop	ear.			Seq Edit			RV Op	erations		Retu		Mair Mer	N	
• *			***	815	P11	817	***		81.33		₩11	Al hearth			March 1	man 2		
		Rv Automa	ted Open		Re1	Re 1		8410	Re Cycl	e Progres	RV 13				RV Manual Op	venitions		
11		RV Automa	ted Open Synthe	con	Per t	Cycle		Pro-			11	Time	Rep		Rir Manual Op Action	entions Oterve	pe .	
	Pet		Synthe	100m	~	Cycle			Rir Cycl Step	e Progres	11	Tre			Action	Otero		
	745			100m	~ ~	Cycle	Prog	part .	Rir Cycl SDep	e Progres	11	Tere	141		Action	Cleave	wage	
	2.ert 2.ert		Synthe	atom sis Isyn		Cycle	20 1 3	pan	Rir Cycl Step	e Progres	11	Time			Action	Otero	wage	
2 9			Synthe Act JN	itions isis isin isin	~ ~ ~	Cycle 1 or 1 11 or 21 11 or 21 11 or 21	Pag.	gan Kara Kara Kara	Ref Cycl 50ep 2 of 1 2 of 2 2 of 2 2 of 2	e Progres	11	Time	1.41 1.41 1.41 1.41		Action cone cone cone cone	Cleave	wage	
- 9 - 9	745		System Act, (N 64, 1412	atons Isis Isin Isin Isin	* * *	Cycle 141 1143 1143	F 1 3 3 3	parti Kury Kury Kury	Re Cycl 50ep 2 of 1 2 of 1 2 of 1	e Progres	11	Time	1.41 1.41 1.41		Action cone cone cone	Clean Load Clean Load Clean Load Clean	wage wage	
- 9 - 9	tet		System Actija 64,144	atons Isis Isin Isin Isin		Cycle 1 41 11 41 11 41 11 41 11 41 11 41 11 41 11 41	F	ram darg darg darg darg darg	Ref Cycl 50ep 2 of 2 3 of 2 3 of 2 2 of 2 3 of 2 2 of 2 2 of 2 2 of 2	e Progres	11	Time	141 141 141 141 141		Action done done done done done	Cleans Load Clea	wage wage	
1 51 2 51 3 51 4 51	745		System Act, (N 64, 1412	ations isis isin isin isin isin	* * * * *	Cycle 1 of 1 11 of 11 11 of 11 11 of 11 11 of 11 11 of 11 11 of 11 11 of 11	F 1 3 1 3 3	ram darg darg darg darg darg darg	Rv Cycl Slop 2 47 2 2 47 2 2 47 2 2 47 2 2 47 2 2 47 2	e Progres	11	Time	1.41 1.41 1.41 1.41 1.41 1.41 1.41		Action 4006 8006 8006 8006 8006	Clean Load Clean Load Clean Load Clean	wage wage wage	
1 51 2 51 3 51 4 51 5 51	145 145 145 145		System Active (4,144) (4,144) (4,144) (4,144)	attens Isis Isin Isin Isin Isin		Cycle 1 41 11 41 11 41 11 41 11 41 11 41 11 41 11 41		ram dang dang dang dang dang dang dang dang	Rir Cycl 50rg 2 or 1 2 or 2 3 or 2 2 or 2 2 or 2 2 or 2 2 or 2 2 or 2 2 or 2	e Progres	11	Time	Lark Lark Lark Lark Lark Lark Lark Lark		Action core core core core core core core	Clear Lear Des Lear Des Lear Des Lear Des	wage wage wage	
	1443 2441 2443 2443 2443		Synthe AC2_04 64,1442 64,1442 64,1442 64,1442 64,1442	atom 199 1990 - 1990 - 1990 - 1990 - 199		Cyth 1 41 11 411 11 41 11 41 1		rury dary dary dary dary dary dary	Ref Cycl 5010 2 47 1 2 47 2 3 47 1 2 47 2 3 47 2 2 47 2 2 47 2 2 47 2 2 47 2	e Progres	11	Time	Lark Lark Lark Lark Lark Lark Lark		Action itone itone itone itone itone itone itone	Lead Oran Lead Oran Lead Oran Lead Oran Lead Oran Lead Oran	wage wage wage	
	145 145 145 145		System Active (4,144) (4,144) (4,144) (4,144)	atom 199 1990 - 1990 - 1990 - 1990 - 199		Cycle 1 41 11 41 11 41 11 41 11 41 11 41 11 41 11 41	N 1 3 3 3 3 3 3	ram dang dang dang dang dang dang dang dang	Rir Cycl 50rg 2 or 1 2 or 2 3 or 2 2 or 2 2 or 2 2 or 2 2 or 2 2 or 2 2 or 2	e Progres	11	Time	Lark Lark Lark Lark Lark Lark Lark Lark		Action core core core core core core core	Clear Lear Des Lear Des Lear Des Lear Des	wage wage wage	
1 54 2 54 3 54 4 54 5 54 5 54 5 54 5 54 5 54 5	1443 2441 2443 2443 2443		Synthe AC2_04 64,1442 64,1442 64,1442 64,1442 64,1442	atom 198 1991 - 1994 - 1994 - 1994 - 199		Cytle 1 or 1 31 or 31 31 or 33 31 or 33	Per 0 00 0 00 0 00 0 00 0 00 0 00 0 00 0	9981	Per Cycl 5000 2 of 2 2 of 2 0 0 0 0 0 0 0 0 0 0 0 0 0	e Progres	11	Tree	1.011 1.011 1.011 1.011 1.011 1.011 1.011 1.011		Actor 0005 0005 0005 0005 0006 0006 0006	Lead Oran Lead Oran Lead Oran Lead Oran Lead Oran Lead Oran	waje waje waje	
1 9 2 9 3 9 4 9 4 9 2 9 2 9 2 9 2 9	tert tert tert tert tert tert		5ynthe Ac2,34 (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344)	ston 199 - 199 - 199 - 199 - 199 - 199 - 199 - 199 - 199 -		Cycle 1 or 1 31 or 33 10 or 33	Per 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	2201 1023 1023 1023 1023 1023 1023 1023	Riv Cycl 5000 2 47 2 3 47 2 4 5 47 4 4 5 47 4 5 47 4 4 5 47 4 5 47 4 5 47 4 5 47 4 5 47 4 5 47 4 5	e Progres	11	Tree	Lark Lark Lark Lark Lark Lark Lark		Action 0.006 0.006 0.006 0.006 0.006 0.006	Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear	waji waji waji waji	
1 9 2 9 3 9 4 9 4 9 2 9 2 9 2 9 2 9	1443 1745 2441 2443 2443		Synthe Ac7,34 (4,345) (4,35) (4,35	ston 199 - 199 - 199 - 199 - 199 - 199 - 199 - 199 - 199 -		Cytle 1 or 1 31 or 31 31 or 33 31 or 33	Per 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	9981	Per Cycl 5000 2 of 2 2 of 2 0 0 0 0 0 0 0 0 0 0 0 0 0	e Progres	11	Time	1.011 1.011 1.011 1.011 1.011 1.011 1.011 1.011		Actor 0005 0005 0005 0005 0006 0006 0006	Clean	waji waji waji waji	
2 9 3 9 4 9 5 9 6 9 7 9 8 9 9 9 9 9 9 9 9 9 9 9 9 9 9	tert tert tert tert tert tert		5ynthe Ac2,34 (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344) (94,344)	atons Isla Islan - San -		Cyth 141 147 147 147 147 147 147 147 147 147		para dang dang dang dang dang dang dang dan	Religid 50ep 2472 2472 2472 2472 2472 2472 2472 2472 2472 2472	e Progres	11	Time	Lori Lori Lori Lori Lori Lori Lori Lori		Actor 0006 0006 0006 0006 0006 0006 0006	Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear	wage wage wage wage wage	
1 9 2 9 3 9 4 9 5 9 5 9 4 9 7 9 5	tart tart tart tart tart tart tart		Synthe Ac2-JN (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2) (e,1er2) (e,1er2 (e,1er2) (e,1er2	ations 198 1991 - 5991 - 5991 - 1991 - 991 - 991 - 991 - 991 - 991 -		Cycle 1441 11443 11443 11443 11443 11443 11443 11443 11443 11443 11443 11443 11443			Reflect Step 2 41 2 42 2 42 2 42 2 42 2 42 2 42 2 42	e Progres	11		Lori Lori Lori Lori Lori Lori Lori Lori		Actem 0.04 0.04 0.04 0.04 0.04 0.04 0.04 0.0	Cleara Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear	waya waya waya waya waya waya waya	
b 99 2 99 3 99 4 99 5 99 5 99 6 99 6 99 6 99 6 99 6 99 6	145 145 145 145 145 145 145 145		System AC3_3A (44,144)(44,144) (44,144)(44,144) (44,144)(44,144) (44,144)(44,144)(44,144)(44,144)(44,144)(44,144	ations ispin - ispin - ispin - ispin - yin -		Cycle 1 441 1143 1143 1143 1143 1143 1144 1144		2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	Re Cycl Step 2 er 2 2 er 2 2 2 er 2 2 er 2 2 er 2 2 er 2 2 er 2 2 er 2 2 2 er 2 2 2 er 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	e Progres	11	Tere	Lori Lori Lori Lori Lori Lori Lori Lori		Actem 0.016 0.016 0.016 0.016 0.016 0.016 0.016 0.016 0.016 0.016	Cleara Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear	waje waje waje waje waje	
b 99 2 99 3 99 4 99 5 99 5 99 6 99 6 99 6 99 6 99 6 99 6	145 145 145 145 145 145 145 145 145		Synthe Ac2-JN (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2 (e,1er2) (e,1er2) (e,1er2 (e,1er2) (e,1er2	ations ispin - ispin - ispin - ispin - yin -		Cycle 1441 11443 11443 11443 11443 11443 11443 11443 11443 11443 11443 11443 11443			Reflect Step 2 41 2 42 2 42 2 42 2 42 2 42 2 42 2 42	e Progres	11	Tee Internet	Lori Lori Lori Lori Lori Lori Lori Lori		Actem 0.04 0.04 0.04 0.04 0.04 0.04 0.04 0.0	Cleara Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear	waje waje waje waje waje	
b 99 2 99 3 99 4 99 5 99 5 99 6 99 6 99 6 99 6 99 6 99 6	144 144 144 144 144 144 144 144 144		System AC3_3A (44,144)(44,144) (44,144)(44,144) (44,144)(44,144) (44,144)(44,144)(44,144)(44,144)(44,144)(44,144	ations ispin - ispin - ispin - ispin - yin -		Cycle 1441 1142 1142 1142 1142 1142 1142 114			Recycl Soep 2 41 2 42 2 42 2 42 2 42 2 42 2 42 2 42	Reg		See .	Lock Lock Lock Lock Lock Lock Lock Lock		Actim 0.006 0.006 0.006 0.006 0.006 0.006 0.006 0.006 0.006	Cleara Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear Land Clear	waje waje waje waje waje	
b 99 c 90 c 90 c 90	144 144 144 144 144 144 144 144 144		5 yethe AC3_(A) (4,140	ations ispin - ispin - ispin - ispin - yin -		Cycle 1441 1148 114800 11480 114800 114800 114800 11480000000000			Ref Cycl Step 2 443 2 442 2 443 2 443 2 443 2 443 2 443 2 443 2 443 2 443 2 443	Reg		See	Lot		Actim 0.005 0.005 0.005 0.005 0.005 0.005 0.005 0.005 0.005	Cleare Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear Lead Clear	waje waje waje waje waje	

SOFTWARE FEATURES

- Select sequence(s) and assign a program for each cycle
- Monitor progress of all reactions
- Customizable amino acid and solvent tables
- Good Laboratory Practice data recording, diagnostics and error reporting
- Flexible restart options for interrupted syntheses

SINGLE-SHOT

All positions can be used in Single-Shot mode. 12 positions may only be used as Single-Shots. Add expensive special reagents to any reactor without priming or wasting a drop. Save money with no manual interuption required.

SAFE-RESPONSE FEATURE

Keep your reactions safe during unexpected interruptions:

- Drain active reaction vessels if needed
- Wash resin with solvent of your choice
- Prevent side reactions and undesired by products

High Quality Chemicals, Amino Acids & Accessories

MAXIMIZE YOUR PERFORMANCE

Reliable Chemistries:

- Amino Acids, resins, activators, solvents and more:
 - Guaranteed weight, purity and quantity
 - All of our chemicals are 100% QC tested
 - 100% fully synthetic amino acids
 - Available in custom pack sizes

Total Solution Delivering Uncompromising Purity:

- Instrument Performance
 - Symphony X peptide synthesizers are optimized using our high quality chemistries
 - Minimize down time caused by clogging resulting from impure and unstable reagents
- Better Synthesis Results
 - Synthesize your peptide right the first time
 - Consistent quality reduces synthesis variability
 - Better yield and purity with high quality reagents
- Rapid and Reliable Setup
 - Start syntheses quickly with convenient, affordable, pre-weighed chemicals

STRENGTHEN YOUR SYMPHONY X

- Symphony X Start-Up Kits
 - Out-of-the box peptide synthesis reagents
 - 24 Plastic reaction vessels for simultaneous use
 - Pre-packed amino acids for test peptides
- Symphony X Pre-Packed Bottles
 - Ready-to-use amino acid bottles measured and QC tested
 - Sizes available in 5 mmol, 10 mmol and 20 mmol quantities
- Bulk Packaging
 - 25 gram, 100 gram, kilogram and custom lots available
 - Pseudoproline dipeptides available in bulk quantities
 - Deprotect and base solutions, pre-mixed and ready to run
 - TFA for cleaving and premium wash solvents available in various quantities

CUSTOMIZE YOUR SYMPHONY X

 We supply reaction vessels and amino acid bottles to provide added functionality and flexibility

Visit www.gyrosproteintechnologies.com and our webshop for all your peptide synthesis needs - amino acids, resins, coupling agents, solvents and more.

Symphony X Specifications

Number of reaction vessels	24 (12 with pre-activation)
Synthesis scale range	0.005* mmol - 24 mmol (up to ~2 g of resin per RV) *minimum deliveries 1 mL or 0.5 mL with Single-Shot
Reaction vessel volume	Plastic, disposable - 10 mL, 45 mL sizes Borosilicate glass - 10 mL, 40 mL
Number of solvent positions/sizes	1 primary solvent, 20 L capacity. Other sizes available
	7 user-defined solvent positions, 0.5 $$ - 8 L capacities, larger capacity vessels available
Number of amino acid positions	Up to 40 amino acid positions 10 mL, 120 mL and 400 mL bottles
Chemistries supported	Fmoc, t-Boc, organic, peptoid, combinatorial, branched, PNA
Activation	DIC/HOBt, DIC/OxymaPure, HBTU, HATU, HCTU, TBTU, PyBOP
Throughput	User defined
Fluid measurement technique	Fixed loop measurement with optical sensors
Fluid transfer method	Positive pressure with nitrogen
Agitation method	Nitrogen bubbling and/or oscillation mixing can be adjustable and programmable
Optional heating method	Induction heating with IR pyrometer temperature sensing
Cleavage	Automatic, programmable (for Fmoc chemistry only)
Waste container	(2) 20 L D.O.T rated containers with over-flow sensor in cap
Control interface	Proprietary Symphony X control software Internal CPU and touchscreen monitor included
Reporting	Real-time log file updated for every instrument function, print to external file, pull up on screen or printer
Power	115 V/60 Hz or 230 V/50 Hz
Dimensions	45"W, 29"D, 66"H (114 cm W, 74 cm D, 168 cm H)
Weight	550 lbs. (250 kg)
Warranty	One year, parts and labor
Optional	Infrared (IR) Heating
	UV monitoring

ORDERING INFORMATION

Catalog #	Description
SYMX-110	Symphony X Peptide Synthesizer, 115 V
SYMX-230	Symphony X Peptide Synthesizer, 230 V

Gyros Protein Technologies logo, Intellisynth, Prelude, PS3, PurePep, Single-Shot, Sonata, Symphony and Tribute are trademarks of Gyros Protein Technologies Group.

COMU, OxymaPure, PyClock, K-Oxyma and Py-Oxim are trademarks owned by Luxembourg Bio Technologies Ltd. Learn more at http://www.luxembourg-bio.com/

All other trademarks are the property of their respective owners. $\ensuremath{\mathbb{O}}$ Gyros Protein Technologies 2017

D0029505/082021/B

For more information: peptides@gyrosproteintech.com www.gyrosproteintechnologies.com

Office:

4675 South Coach Drive Tucson, AZ 85714, U.S.A. Phone: 520-629-9626 Toll Free: 800-477-6834