

# THE IMPACT OF OCCLUSIONS ON DOSE ACCURACY

**INSTECH**

ITO Conference, Cologne, September 2015

# 3 Syringe Pumps

INSTECH




<b>Model</b>	OrchesTA Model 100	Baxter AS50	Harvard Apparatus PHD2000
<b>Stated Accuracy</b>	$\pm 3\%^*$	$\pm 3\%^*$	$\pm 0.35\%^*$

\* All excluding syringe variability

DISCLAIMER:  
Instech sells [OrchesTA](#) and [Harvard Apparatus](#) pumps.  
The Baxter AS50 is discontinued.

# Occlusions in Infusion Studies

Forgot to  
open pinch  
clamp



Rat infusion set  
dead volume:  $\sim 350\mu\text{l}$


Kink in  
external tube

Kink in catheter  
or tip occlusion


# What Can Happen?

## PLANNED DOSE

3ml BD plastic syringe




# What Can Happen?



# What Can Happen?

Catheter occlusion  
15 minutes into study



Scenario	ml into animal	ml left in syringe
Alarm, back off, fix, restart	2.0	0.0
Alarm, fix, bolus, restart	2.0	0.0
No alarm, pressure build up	0.5	1.5
No alarm, <b>leak</b>	0.5	0.0

# Results: OrchesTA Model 100


INSTECH



Scenario	Alarm after 18sec Auto back off Fix occlusion Restart
Vol. dispensed (error)	1.97ml (-1.5%)
Vol. delivered (error)	1.97ml (-1.5%)

# Results: Baxter AS50

INSTECH




Scenario	Alarm after <b>2.0min</b> Manual back off Fix occlusion Restart
Vol. dispensed (error)	1.59ml (-20%)
Vol. delivered (error)	1.59ml (-20%)



# Results: Harvard PHD2000

INSTECH



Leak

Scenario	No alarm <b>Leak</b> after 50min
Vol. dispensed (error)	1.95ml (-2.3%)
Vol. delivered (error)	0.5ml (-75%)

# Summary



Model	OrchesTA Model 100	Baxter AS50	Harvard Apparatus PHD2000
Occlusion Sensor	In-line force transducer	Motor current monitor	None
Time to alarm	18sec	120sec	n/a
Dose Error	-1.5%	-20%	-75%
Reporting Error*	0%	0%	-75%

\* If syringes are weighed to confirm every dose. Does not include weighing error.

## BACKUP DATA FOR "IMPACT OF OCCLUSIONS ON DOSE ACCURACY"

Presentation by Jeffrey Sommers (Instech Laboratories) at ITO Conference in Cologne Germany, 7 September 2015

### Experiment Parameters

Intended dose: 2.0ml infusion over 1hr  
 Equipment used: Pump, plastic syringe (3 or 20ml), Instech KVAH95T rat tether kit, VAH95AB rat harness, 15cm 3Fr PU catheter  
 Occlusion simulation: Catheter tip occluded at 15min (after 0.5ml infused) with pinch clamp. On occlusion alarm, follow pump operating manual to relieve pressure; undo pinch clamp; restart.

### Results

PUMP	SYRINGE	OCCLUSION SETTING	TIME TO SENSE OCCLUSION (MIN)	SYRINGE WEIGHT (g)			DISPENSE ERROR	LEAK?	DOSE ERROR	
				START	END	CHANGE				
OrchesTA model 100	B-D 3ml	normal	Trial 1:	0.3	5.93	3.95	1.98	-1.0%	no	-1.0%
			Trial 2:	0.3	5.86	3.90	1.96	-2.0%	no	-2.0%
			Trial 3:	0.3	5.77	3.80	1.97	-1.5%	no	-1.5%
			Average:	0.3			1.97	-1.5%		<b>-1.5%</b>
Medfusion 2001	B-D 3ml	medium	Trial 1:	1.2	6.25	4.43	1.82	-9.0%	no	-9.0%
			Trial 2:	1.2	6.40	4.60	1.80	-10.0%	no	-10.0%
			Trial 3:	1.2	6.32	4.53	1.79	-10.5%	no	-10.5%
			Average:	1.2			1.80	-9.8%		<b>-9.8%</b>
SAI 3D	B-D 3ml	medium	Trial 1:	2.1	6.40	4.60	1.80	-10.0%	no	-10.0%
			Trial 2:	2.1	6.05	4.27	1.78	-11.0%	no	-11.0%
			Trial 3:	2.1	6.31	4.52	1.79	-10.5%	no	-10.5%
			Average:	2.1			1.79	-10.5%		<b>-10.5%</b>
Baxter AS50	B-D 3ml	medium	Trial 1:	2.1	5.75	4.18	1.57	-21.5%	no	-21.5%
			Trial 2:	1.8	6.05	4.42	1.63	-18.5%	no	-18.5%
			Trial 3:	2.1	5.91	4.33	1.58	-21.0%	no	-21.0%
			Average:	2.0			1.59	-20.3%		<b>-20.3%</b>
Harvard PHD2000	B-D 3ml	n/a	Trial 1:	n/a	6.03	4.09	1.94	-3.0%	yes @ 48min	-75%
			Trial 2:	n/a	6.11	4.15	1.96	-2.0%	yes @ 54min	-75%
			Trial 3:	n/a	5.85	3.89	1.96	-2.0%	yes @ 48min	-75%
			Average:				1.95	-2.3%		<b>-75%</b>
Harvard 11 Elite	B-D 3ml	1% (force limit)	Trial 1:	did not alarm	6.22	4.26	1.96	-2.0%	yes @ 48min	-75%
			Trial 2:	did not alarm	5.99	4.02	1.97	-1.5%	yes @ 54min	-75%
			Trial 3:	did not alarm	5.88	3.92	1.96	-2.0%	yes @ 50min	-75%
			Average:				1.96	-1.8%		<b>-75%</b>
OrchesTA model 100	B-D 20ml	normal	Trial 1:	11.0	28.57	26.78	1.79	-10.5%	no	-10.5%
			Trial 2:	11.0	26.78	24.97	1.81	-9.5%	no	-9.5%
			Trial 3:	11.0	24.97	23.17	1.80	-10.0%	no	-10.0%
			Average:	11.0			1.80	-10.0%		<b>-10.0%</b>
Medfusion 2001	B-D 20ml	medium	Trial 1:	13.2	23.15	21.57	1.58	-21.0%	no	-21.0%
			Trial 2:	14.4	21.57	20.01	1.56	-22.0%	no	-22.0%
			Trial 3:	14.7	20.01	18.42	1.59	-20.5%	no	-20.5%
			Average:	14.1			1.58	-21.2%		<b>-21.2%</b>
Baxter AS50	B-D 20ml	medium	Trial 1:	15.0	26.97	25.96	1.01	-49.5%	no	-49.5%
			Trial 2:	15.0	25.96	24.99	0.97	-51.5%	no	-51.5%
			Trial 3:	15.0	24.99	23.89	1.10	-45.0%	no	-45.0%
			Average:	15.0			1.03	-48.7%		<b>-48.7%</b>
SAI 3D	B-D 20ml	medium	Trial 1:	23.0	29.54	28.79	0.75	-62.5%	no	-62.5%
			Trial 2:	22.5	28.79	27.90	0.89	-55.5%	no	-55.5%
			Trial 3:	21.6	27.90	27.28	0.62	-69.0%	no	-69.0%
			Average:	22.4			0.75	-62.3%		<b>-62.3%</b>

### DOSE ERROR, 3ML SYRINGE

