

DATE

25-March-2020

TO

Barco Healthcare Display System Customers

FROM

Albert Xthona & Inge Haesaerts, Product Managers

Addendum to Barco Display Cleaning Instructions

Dear Customers,

Many people would like to use stronger concentrations of alcohol-based products to clean their Coronis and Nio display systems. In response, we have tested higher concentrations of alcohol and you should now feel free to use up to 70% solutions of ethanol or isopropyl alcohol as shown below in the annotated section from the User Guide.

Cleaning your display

5.1 Cleaning instructions

To clean the display

Clean the display using a sponge, cleaning cloth or soft tissue, lightly moistened with a recognized cleaning product for medical equipment. Read and follow all label instructions on the cleaning product. In case of doubt about a certain cleaning product, use plain water.

Do not use following products:

- Alcohol/solvents at higher concentration > ~~5%~~ 70%
- Strong alkalis lye, strong solvents
- Acid
- Detergents with fluoride
- Detergents with ammonia
- Detergents with abrasives
- Steel wool
- Sponge with abrasives
- Steel blades
- Cloth with steel thread

 **CAUTION:** Take care not to damage or scratch the front glass or LCD. Be careful with rings or other jewelry and do not apply excessive pressure on the front glass or LCD.

 **CAUTION:** Do not apply or spray liquid directly to the display as excess liquid may cause damage to internal electronics. Instead, apply the liquid to a cleaning cloth.

This is applicable to models with {MDCC, MDCG, MDMC, MDMG, MDNC, MDNG, MDRC}, e.g. MDMC-12133, MDCC-6530, manufactured after 1-January 2008.

Regards,
Albert & Inge

FDR Go PLUS & FDR Visionary Suite

OEM Recommended Cleaning Instructions

To avoid secondary infection

Outline

This document is Shimadzu's voluntary cleaning instructions for the equipment proposed to minimize the risk of secondary infection. Every country must deal with regulations separately. Standards required for specific clinical use differs from place to place.

General Cleaning

Before using disinfectant agents, it is necessary to clean the surface of the equipment in advance. Please remove the dirt on the surface using a cloth without lint, or a soft toothbrush that does not damage the surface. If dirt is difficult to remove, please dilute dishwasher detergent and soap with lukewarm water, and use them.

*Please make sure that the dishwasher detergent and soap used for cleaning is completely wiped off before disinfectant treatment.

Disinfectant Agents (Disinfecting)

The following agents can be used. But there are some restrictions of use according to the target parts as stated in table below. Please make sure that disinfectant agent is wiped off after disinfecting is completed.

(A) Chlorine disinfectants

- Sodium dichloroisocyanurate solution (1% maximum)
- Sodium hypochlorite solution (1% maximum)

(B) Alcohol disinfectants

- Commercially available isopropyl alcohol solution (Up to 99 wt% can be used)
- Rubbing alcohol (76.9-81.4 vol% Ethanol, Isopropyl alcohol as an additive)

Table 1: Restriction of usable disinfectant according to the target parts

Parts	Usable disinfectant	Risk
Plastic parts	Chlorine disinfectant only.	Possibility of deterioration if alcohol disinfectant is used. Especially acrylic resin parts <ul style="list-style-type: none"> ● General x-ray system Acrylic table top (transparent resin part) ● Collimator (front transparent resin part) ● Touch screen for mobile equipment ● dosimeter
Metal parts without paint	Alcohol disinfectant only.	Possibility of deterioration if chlorine disinfectant is used.

 Caution

Please keep use of disinfectant agent to a minimum. Due to long-term disinfection, fading and cracking may occur in the exterior or resins. For example, rubber and plastics may deteriorate. If the product changes due to disinfecting, immediately stop using the product and request repair to your service representative.

 Caution

Please do not spray disinfectant directly on the equipment. If disinfectant penetrates the device, it may cause malfunction or other problems. Please turn off the power during disinfecting. When turning on the power after disinfecting, please thoroughly ventilate the inside of the examination room. If flammable gas remains in the room, there is a potential for danger of ignition, smoke, explosion, or electric shock when turning on the system.

Examples of agents that can be used

List of disinfectant agents are extracted from the US Environmental Protection Agency (EPA) list of approved disinfectants for use against a new coronavirus infection (COVID -19).

This list does not guarantee the efficacy of disinfection.

Product name	company	Active ingredient	EPA registration number
Sani-Cloth Bleach Germicidal Disposable Wipe	DPI	Sodium hypochlorite	9480-8
PURELL Professional Surface Disinfectant Wipes	GOJO Industries	Ethanol	84150-1
Micro-kill Bleach Germicidal Bleach Wipes	Medline Industries	Sodium hypochlorite	37549-1
Clorox Healthcare® Bleach Germicidal Wipes	Clorox Professional Products Company	Sodium hypochlorite	67619-12

This instruction is applicable to the products manufactured by Medical Systems Div., Shimadzu Corporation only.