

# INTERPRETING SINGLE ANTIBIOTIC RESULTS

**Veterinary BioFilm Susceptibility Report**

Lab Number

Patient

Order Date

Strain ID **95.00% Escherichia coli**


Gram **Negative**

Breed **Golden Retriever**

Site Swab


Clinic Name

Veterinarian



**BecScreen**

**Golden Retriever**



Antibiotics	SIR Breakpoints	Minimum Inhibitory Concentration (MIC)		Minimum Biofilm Eradication Concentration (MBEC)	
Amikacin	S: ≤16, I: 32, R: ≥64	≤8	<b>S</b>	16	<b>S</b>
Amoxicillin + Clavulanate	S: ≤4/2, I: 8/4-16/8, R: ≥32/16	>32/16	R	>32/16	R
Ampicillin	S: ≤8, I: 16, R: ≥32	>32	R	>32	R
Azithromycin	S: ≤1, R: ≥2	8	R	>8	R
Cefazolin	S: ≤8, R: ≥32	8	<b>S</b>	>8	R
Cefovecin	S: ≤2, R: ≥4	8	R	>8	R
Cefpodoxime	S: ≤2, I: 4, R: ≥8	>8	R	>8	R
Ceftiofur	S: ≤2, I: 4, R: ≥8	>8	R	>8	R
Chloramphenicol	S: ≤8, I: 16, R: ≥32	≤4	<b>S</b>	>32	R
Ciprofloxacin	S: ≤1, I: 2, R: ≥4	≤0.5	<b>S</b>	≤0.5	<b>S</b>
Doxycycline	S: ≤4, I: 8, R: ≥16	≤2	<b>S</b>	≤2	<b>S</b>
Enrofloxacin	S: ≤0.5, I: 1-2, R: ≥4	≤0.5	<b>S</b>	≤0.5	<b>S</b>
Gentamicin Sulfate	S: ≤2, I: 4, R: ≥8	≤2	<b>S</b>	4	I
Imipenem	S: ≤1, I: 2-8, R: ≥16	16	R	>16	R
Marbofloxacin	S: ≤1, I: 2, R: ≥4	≤0.5	<b>S</b>	≤0.5	<b>S</b>
Minocycline	S: ≤4, I: 8, R: ≥16	4	<b>S</b>	>16	R
Orbifloxacin	S: ≤1, I: 2-4, R: ≥8	≤1	<b>S</b>	≤1	<b>S</b>
Piperacillin	S: ≤16, I: 32-64, R: ≥128	>128	R	>128	R
Tetracycline	S: ≤4, I: 8, R: ≥16	≤2	<b>S</b>	>16	R
Ticarcillin / Clavulanate	S: ≤16/2, I: 32/2-64/2, R: ≥128/2	≤16/2	<b>S</b>	128/2	R
Tobramycin	S: ≤4, I: 8, R: ≥16	≤2	<b>S</b>	4	<b>S</b>
Trimethoprim / Sulfamethoxazole	S: ≤2/28, R: ≥4/76	≤0.5/9.5	<b>S</b>	4/76	R

FIGURE 1: These results indicate that there are 6 antibiotic options that can provide possible\* treatment for the above infection in a biofilm state (**SEE THOSE HIGHLIGHTED IN GREEN UNDER MBEC HEADER**)

The normal therapeutic dosage for a patient of this weight/size would be utilized, unless directed otherwise by your pharmacist or reference guide.

Note the vast increase in suggested antibiotics using standard culture and sensitivity (MIC), that would have suggested ineffective antibiotics.

\*unless contraindicated



S: Susceptible    I: Intermediate    R: Resistant

BecScreen is a qualitative, in-vitro diagnostic test designed for use in determining antimicrobial susceptibility of both planktonic and biofilm isolates. Minimum Inhibitory Concentration (MIC) results are for planktonic cells susceptibility and Minimum Biofilm Eradication Concentration (MBEC) results are for biofilms susceptibility. Results for each individual antibiotic may be applied to other antibiotics in the same class.



# INTERPRETING ANTIBIOTIC COMBINATIONS RESULTS

Veterinary BioFilm Susceptibility Report					
Lab Number		Breed		Cocker Spaniel	
Patient	Spot	Site Swab	Skin		
Order Date	Strain ID		Clinic Name		
Staphylococcus aureus					
Gram Positive					
Antibiotics	SIR Breakpoints	Minimum Inhibitory Concentration (MIC)	Minimum Biofilm Eradication Concentration (MBEC)		
Amoxicillin:Clavulanate + Clindamycin	1:0.5/0.5 - 2:1/1	≤1:0.5/0.5	S	>2:1/1	R
Amoxicillin:Clavulanate + Imipenem	1:0.5/1 - 2:1/2	≤1:0.5/1	S	>2:1/2	R
Amoxicillin:Clavulanate + Orbifloxacin	1:0.5/1 - 2:1/2	≤1:0.5/1	S	2:1/2	I
Cefalexin + Clindamycin	1/0.5 - 2/1	≤1/0.5	S	>2/1	R
Cefalexin + Doxycycline	1/4 - 2/8	≤1/4	S	>2/8	R
Cefalexin + Gentamicin	1/2 - 2/4	≤1/2	S	≤1/2	S
Cefalexin + Imipenem	1/1 - 2/2	≤1/1	S	>2/2	R
Cefalexin + Orbifloxacin	1/1 - 2/2	≤1/1	S	2/2	I
Cefalexin + Trimethoprim/Sulfamethoxazole	1/2:38 - 2/2:38	≤1/2:38	S	>2/2:38	R
Chloramphenicol + Clindamycin	8/0.5 - 16/1	≤8/0.5	S	>16/1	R
Chloramphenicol + Erythromycin	8/0.5 - 16/1	≤8/0.5	S	>16/1	R
Chloramphenicol + Gentamicin	8/2 - 16/4	≤8/2	S	≤8/2	S
Chloramphenicol + Orbifloxacin	8/1 - 16/2	≤8/1	S	>16/2	R
Chloramphenicol + Trimethoprim/Sulfamethoxazole	8/2:38 - 16/2:38	≤8/2:38	S	>16/2:38	R
Chloramphenicol + Vancomycin	8/2 - 16/4	≤8/2	S	≤8/2	S
Clindamycin + Ampicillin	0.5/1 - 1/2	≤0.5/1	S	>1/2	R
Clindamycin + Imipenem	0.5/1 - 1/2	≤0.5/1	S	1/2	I
Clindamycin + Orbifloxacin	0.5/1 - 1/2	≤0.5/1	S	1/2	I
Clindamycin + Vancomycin	0.5/2 - 1/4	≤0.5/2	S	≤0.5/2	S
Doxycycline + Chloramphenicol	4/8 - 8/16	>8/16	R	>8/16	R
Doxycycline + Clindamycin	4/0.5 - 8/1	≤4/0.5	S	>8/1	R
Doxycycline + Erythromycin	4/0.5 - 8/1	≤4/0.5	S	>8/1	R
Doxycycline + Imipenem	4/1 - 8/2	≤4/1	S	>8/2	R
Erythromycin + Cefalexin	0.5/1 - 1/2	≤0.5/1	S	>1/2	R
Gentamicin + Amoxicillin/Clavulanate	2/1:0.5 - 4/2:1	≤2/1:0.5	S	≤2/1:0.5	S
Gentamicin + Cefalexin	2/1 - 4/2	≤2/1	S	≤2/1	S
Gentamicin + Ampicillin	2/1 - 4/2	≤2/1	S	4/2	I
Gentamicin + Clindamycin	2/0.5 - 4/1	≤2/0.5	S	≤2/0.5	S
Gentamicin + Doxycycline	2/4 - 4/8	≤2/4	S	4/8	I
Gentamicin + Minocycline	2/4 - 4/8	≤2/4	S	4/8	I
Gentamicin + Orbifloxacin	2/1 - 4/2	≤2/1	S	≤2/1	S
Minocycline + Cefalexin	4/1 - 8/2	≤4/1	S	>8/2	R
Minocycline + Orbifloxacin	4/1 - 8/2	≤4/1	S	8/2	I
Orbifloxacin + Ampicillin	1/1 - 2/2	≤1/1	S	2/2	I
	1/4 - 2/8	≤1/4	S	2/8	I
	1/1 - 2/2	≤1/1	S	>2/2	R
Amoxicillin/Clavulanate	2:38/1:0.5 - 2:38/2:1	≤2:38/1:0.5	S	>2:38/2:1	R
Amoxicillin	2:38/1 - 2:38/2	≤2:38/1	S	>2:38/2	R
Clindamycin	2:38/0.5 - 2:38/1	>2:38/1	R	>2:38/1	R
Doxycycline	2:38/4 - 2:38/8	≤2:38/4	S	≤2:38/4	S
Imipenem	2:38/2 - 2:38/4	≤2:38/2	S	>2:38/4	R
Orbifloxacin	2:38/1 - 2:38/2	>2:38/2	R	>2:38/2	R
Vancomycin	2:38/1 - 2:38/2	≤2:38/1	S	≤2:38/1	S
	2:38/2 - 2:38/4	≤2:38/2	S	≤2:38/2	S
	2/1:0.5 - 4/2:1	≤2/1:0.5	S	≤2/1:0.5	S
	2/1 - 4/2	≤2/1	S	≤2/1	S
	2/1 - 4/2	≤2/1	S	≤2/1	S

FIGURE 2: These results indicate that there are 14 antibiotic options that can provide possible\* treatment for the above infection in a biofilm state (SEE THOSE HIGHLIGHTED IN GREEN UNDER MBEC HEADER)

The normal therapeutic dosage for a patient of this weight/size would be utilized, unless directed otherwise by your pharmacist or reference guide.

Note the vast increase in suggested antibiotics using standard culture and sensitivity (MIC), that would have suggested ineffective antibiotics.

\*unless contraindicated

S: Susceptible I: Intermediate R: Resistant  
 BecScreen is a qualitative, *in-vitro* diagnostic test designed for use in determining antimicrobial susceptibility of both planktonic and biofilm states. Minimum Inhibitory Concentration (MIC) results are for planktonic cells susceptibility and Minimum Biofilm Eradication Concentration (MBEC) results are for biofilm susceptibility. Results for each individual antibiotic may be applied to other antibiotics in the same class.