





Care@Home™

CMS

Reference Manual



August 2021



Table of Contents

- 1. Overview.....7
 - 1.1. The CP Configuration File.....7
 - 1.2. System Requirements8
- 2. CMS Deployment Process.....9
 - 2.1. Installing CMS9
 - 2.2. Identifying the Communications Port.....10
 - 2.3. Configuring Application Settings.....11
 - 2.3.1 Understanding the *config.ini* File.....12
 - 2.3.2 Understanding the *setup.ini* File14
 - 2.3.3 Understanding the *amigo.ini* File17
 - 2.3.4 Understanding the *DTMFConfig.ini* File.....20
- 3. Understanding the CMS Connection.....22
 - 3.1. Establishing a Local Connection23
 - 3.1.1 Local Connection to the CP.....23
 - 3.2. Establishing a Remote Connection24
 - 3.2.1 Establishing Remote Connection with API26
 - 3.2.2 Establishing Remote Connection with SMS.....27
 - 3.2.3 Establishing Remote Connection with PSTN Modem.....28
 - 3.3. Setting up CP Connectivity32
 - 3.3.1 Uploading CP Settings Procedure32
 - 3.3.2 Reconnecting to the CP35
- 4. Getting Started37
 - 4.1. Main Action Icons37
 - 4.2. File Menu Functions39
 - 4.2.1 Open Function40
 - 4.2.2 Quit Function40
 - 4.2.3 Print... Function.....40
 - 4.2.4 Print Log... Function.....40
 - 4.3. The Main Navigator Menu.....41
- 5. The Control Panel Module43
 - 5.1. Account Number43
 - 5.2. ARC Account Number44
 - 5.3. Service Type.....44
 - 5.4. More Options.....45
 - 5.4.1 Photo Config47
 - 5.4.2 Out of Home49

- 5.4.3 Analytics51
- 5.4.4 Periods 52
- 5.4.5 Rules.....54
- 5.4.6 PERS.....58
- 5.4.7 DTMF Config 60
- 5.4.8 Log63
- 5.4.9 VPD Configuration Parameters 64
- 5.4.10 EPA Configuration Parameters68
- 5.4.11 RF Settings.....70
- 5.5. Periodic Test70
- 5.6. Mains Failure70
- 5.7. Temperature71
- 5.8. Emergency Pendant..... 72
- 5.9. Medical Alarm 72
- 5.10. Password 73
- 5.11. Speakerphone..... 73
- 5.12. Remote Call-In74
- 5.13. Sound..... 75
- 5.14. Emergency..... 75
- 5.15. Inactivity.....76
- 5.16. Activity Timer76
- 5.17. Smoke Detector..... 77
- 5.18. The Control Panel Module Buttons 77
- 6. The Communications Module79
 - 6.1. General Communications 80
 - 6.1.1 Link Lost or Restored 80
 - 6.1.2 IP ACK Timeout Configuration 81
 - 6.1.3 Communication 81
 - 6.1.4 DNIS82
 - 6.1.5 Phone Number82
 - 6.1.6 Authentication Password83
 - 6.2. Cellular Communications.....83
 - 6.2.1 Module..... 84
 - 6.2.2 Provider.....85
 - 6.2.3 APN85
 - 6.2.4 Audio Codec86
 - 6.2.5 SIM86
 - 6.2.6 Dialer Configurations 87

- 6.2.7 Call Configurations..... 87
- 6.3. PSTN Communications.....88
- 6.4. Ethernet Communications..... 90
- 7. The Dialer Module 92
 - 7.1. Destination Definitions by CP Type..... 92
 - 7.2. Destination Definitions for an ES6502HC CP..... 94
 - 7.2.1 Defining Communications Types and Channels..... 94
 - 7.2.2 Selecting Communication Channels 97
 - 7.2.3 Defining Receiver Address Attributes 99
 - 7.3. Destination Definitions for a C7000 CP- Pre Version 6.1101
 - 7.3.1 Defining Communications Types and Channels..... 103
 - 7.3.2 Selecting Communication Channels104
 - 7.3.3 Defining Receiver Address Attributes 106
 - 7.4. Destination Definitions for a C7000 CP - Version 6.1 and Later 108
 - 7.4.1 Defining Communications Types and Interfaces.....110
 - 7.4.2 Selecting Interfaces by CP Type111
 - 7.4.3 Defining Receiver Address Attributes 113
 - 7.5. Managing the Destination List 115
 - 7.5.1 Adding a Dialer Destination 115
 - 7.5.2 Updating a Dialer Destination 116
 - 7.5.3 Deleting a Dialer Destination 116
 - 7.6. The Dialing Sequence 116
 - 7.6.1 Dialing Sequence Dual Mode 120
- 8. The Activity Sensors Module 122
 - 8.1. Defining Sensor Device Attributes..... 123
 - 8.2. Managing Activity Sensors..... 125
 - 8.2.1 Adding an Activity Sensor126
 - 8.2.2 Updating an Activity Sensor 127
 - 8.2.3 Deleting an Activity Sensor 127
 - 8.3. Zone Status128
- 9. The Safety & SOS Module.....129
 - 9.1. Defining Attributes for Safety & SOS Devices..... 130
 - 9.2. Specific Safety and Emergency Attributes 131
 - 9.3. Managing Safety & SOS Devices..... 135
 - 9.3.1 Adding a Safety or SOS Device136
 - 9.3.2 Updating a Safety or SOS Device.....136
 - 9.3.3 Deleting a Safety or SOS Device 137
 - 9.4. RF Input Status 137

10.	Voice & LED Module	139
10.1.	General Configuration Parameters	141
10.2.	Reminder Schedule Types.....	144
10.2.1	One-Time Schedule Type.....	145
10.2.2	Weekly Schedule Type	145
10.2.3	Day of the Month Schedule Type	145
10.2.4	Day of Week Schedule Type	145
10.3.	Medication Reminders.....	146
10.4.	Appointment Reminders.....	147
10.5.	Ride Reminders	148
10.6.	Alarm Reminders.....	148
10.7.	Pendant Test Reminder.....	149
11.	Mobile Operators	152
11.1.	APN Information	152
11.2.	Managing the Mobile Operator List.....	153
11.2.1	Adding a Mobile Operator	153
11.2.2	Updating a Mobile Operator	154
11.2.3	Deleting a Mobile Operator	154
11.3.	Load APN Information.....	155
12.	The Custom Labels Module	157
12.1.	Defining Custom Labels	157
12.2.	Managing Custom Labels.....	158
13.	The Log Events Module	159
13.1.	Log Events Module Actions	160
13.2.	Number of Events	160
13.3.	Read Log	160
13.4.	Reset Log	160
13.5.	Save Log	161
13.6.	Open Log	161
13.7.	Print Log.....	162
13.8.	Read Log File.....	163
14.	Using CMS for Care@Home™ Implementation	164
14.1.	Loading Factory Settings.....	164
14.2.	Verifying Factory Settings.....	165
14.2.1	Account Number	165
14.2.2	ARC Account Number	165
14.2.3	APN Settings.....	166
14.2.4	Customizing the CP Configuration.....	167

- 14.3. Configuring the Care@Home™ CP..... 168
 - 14.3.1 Dialer Module..... 168
 - 14.3.2 Redefining Dial Sequence Retries..... 170
 - 14.3.3 Activity Sensors Module..... 171
 - 14.3.4 Safety and SOS Module 175
 - 14.3.5 Log Events Module.....176
- 15. Building a Kit-It Configuration177
 - 15.1. Introduction to Kitting.....177
 - 15.2. Understanding a Kit-It Configuration178
 - 15.3. Creating a Predefined Kit Configuration Using CMS in Kit-It Mode179
- 16. Remote Boot Software Updates.....182
 - 16.1. Launching the Remote Boot Feature182
 - 16.2. Remote Boot Setup183
 - 16.3. The Software/Firmware Update Procedure 184
- 17. Using the Remote Multiple Device Manager Tool 189
 - 17.1. Prerequisites 189
 - 17.2. Upgrading the CPs190
 - 17.3. Managing RMDM Processing 194
 - 17.4. Upgrade Progress Reporting195
 - 17.4.1 Panel List Details Box..... 196
 - 17.4.2 Current CMS task 196
 - 17.4.3 Session Statistics197
 - 17.4.4 Export feature 198
- Appendix A USB to RS232 Adapters 199
 - The USB Gear Serial Adapter 199
 - The VScom USB-COM Mini.....200
- Appendix B Voice Announcements.....201
- Appendix C Alerts by Type203

1. Overview

The Care@Home™ CMS application provides a tool to configure the control panel (CP) parameters that define the operating environment for the Care@Home™ functionality. The Care@Home™ CMS tool allows you to customize Care@Home™ as required for a resident's needs, before and after its installation on the premises.

The Care@Home™ CMS application supports the following panels:

- ES6502HC
- C7000

All references to the CP in this guide refer to both the ES6502HC and the C7000 unless otherwise noted.

When availability criteria are noted for a feature, if a panel type is not mentioned, the feature is fully supported by both panel types.

For more information about the C7000, refer to the ESUGSC077 Care@Home™ C7000 User Guide.

You can modify the panel parameter values as well as update the CP software and the firmware for its peripheral devices. These modifications can be performed locally on the resident's premises or remotely by the monitoring station.



NOTE: This user guide refers to Care@Home™ CMS version 6.3.20.

1.1. The CP Configuration File

Each CP contains a set of parameters defining the CP configuration and other information for the firmware of the CP and its peripheral devices. This data defines and controls the operation of Care@Home™.



CMS provides you with the tools to save this set of configured parameters in an external CP configuration file, with the extension **.cpf**.

The purpose of the **.cpf** file is two-fold:

- The **.cpf** file may be used as an offline backup to store a given set of configuration parameters and peripheral firmware information. The **.cpf** file can be kept indefinitely.
- The **.cpf** file allows you to distribute a given set of configuration parameters and peripheral firmware by updating and managing multiple CPs, using the CMS warehouse locally and the RMDM remotely. Refer to 17 Using the Remote Multiple Device Manager Tool on page 189.

For information on creating predefined CP configuration files to update multiple CPs, refer to 15 Building a Kit-It Configuration on page 177.

The tools to manage a CP configuration are:

- : Transfer data from the CP to CMS. Refer to 4.1 Main Action Icons on page 37.
- : Transfer data from CMS to the CP. Refer to 4.1 Main Action Icons on page 37.
- **Save** (in **File** menu): Saves the configuration parameter settings from CMS to a CP configuration file (**.cpf**). Refer to 4.2 File Menu Functions on page 39.
- **Open** (in **File** menu): Opens a CP configuration file (**.cpf**) and uploads the configuration parameter settings to CMS, and at the end of the session, can download changes to the CP. Refer to 4.2 File Menu Functions on page 39.

The data downloaded from the CMS includes:

- The CP parameter values
- The peripheral device parameter values
- The operation-related configuration values

1.2. System Requirements

The following are the PC minimum requirements for the CMS application:

- 1 GHz processor
- 1 GB of RAM
- 1 GB free disk space
- Windows 8, 8.1, and 10

2. CMS Deployment Process

Implementation of the Care@Home™ CMS application includes:

- Installing CMS
- Setting administrator privileges
- Identifying COM port
- Configuring application settings




2.1. Installing CMS

Contact Essence professional services for the latest version of the CMS installation file.

To install Care@Home™ CMS:


1. Copy the CMS installation execution file, **Care@Home CMS <release number> setup.exe**, to the PC.
2. Run the CMS installation execution file and follow the on-screen instructions.

The Care@Home™ CMS applications are installed, and their icons appear on your desktop:

-  - Care@Home™ CMS
-  - Care@Home™ CMS Warehouse
-  - Care@Home™ Remote Multiple Device Manager (RMDM)


2.2. Identifying the Communications Port

Associate the communications (COM) port used to connect to the CP with the COM port selected in the **Connect** window of the Care@Home™ CMS application.




To identify the COM port number used to connect to the CP:

1. Click  .

2. Click **Control Panel** on the Windows  menu.

3. Click **Device Manager** on the **Control Panel**  menu. The **Device Manager** window appears.

- Under **Ports (Com & LPT)**, locate the communications port to which the CP cable is to be connected.



2.3. Configuring Application Settings

To configure application settings, update your configuration **.ini** files, located at the following default path:


[c:\Program Files \(x86\)\Essence\Care@Home CMS\](c:\Program Files (x86)\Essence\Care@Home CMS\)

Use a text editor in administrator mode for editing any of the following files:

- [config.ini](#)
- [setup.ini](#)
- [amigo.ini](#)
- [DTMFConfig.ini](#)

2.3.1 Understanding the *config.ini* File


The **config.ini** file defines the CP communication connections from modems and ports to the CP.




The parameter values for the local COM port connection between the PC on which the Care@Home™ CMS is installed, and the CP are defined in the **config.ini** file.

If the COM port is defined incorrectly, the configuration file must be edited, using a text editor, to correct the assignment statements accordingly.

For example, the user chooses to use serial port COM5 to connect the CP to the local PC.




In the configuration file, **config.ini**, the SERIAL1 port is defined as COM1.




The user must edit the configuration file to redefine the SERIAL1 port as COM5.

To edit the *config.ini* file:

1. Open the *config.ini* file in a text editor.
2. For [Config], modify: SERIAL1=COM5->Cable
3. For [SERIAL1], modify: ComNumber=5
4. Save the changes.



When the user opens the **Connect** window, the first COM port listed in the *config.ini* file appears as the default COM port.



The table below defines the parameters included in the *config.ini* file.

Sections	Parameter	Description or Value
[ALL]	TRACELEVEL	Not in use
[Config]	SERIAL<n>=Com<m>	Section name=port name
	MODEM<n>=Com<m>	Modem port
[SERIAL<n>]	Type=0	Serial connection section
	ComNumber=<n>	COM port number
	Baud=<nnnnn>	Baud rate overridden by CMS
	DataBits=8	Standard

Sections	Parameter	Description or Value
	Parity=0	Standard
	StopBits=1	Standard
	Note: Amigo only Channel=<n>	4 - Amigo
[MODEM<n>]	Type=1	Modem connection
	ComNumber=<n>	COM port number
	Baud=<nnnnnn>	Baud rate overridden by CMS
	DataBits=8	Standard
	Parity=0	Standard
	StopBits=1	Standard
	Channel=<n>	0 – GSM; 1 – PSTN
	ATCommands=\Modem_PSTN. ini	Modem .ini file

2.3.2 Understanding the *setup.ini* File

The **setup.ini** file assigns the parameter values which determine how the CMS and the CMS warehouse work with the CP.

The following table defines the parameters included in the **setup.ini** file:

Category	Parameter	Values	Description	Data UM
[Setup]	EntriesStartAdd	&H11C000	Not in use	
	EntriesEndAdd	&H11FFFF	Not in use	
	TimeOutACKFrm	15	Not in use	seconds
	TimeOutCMDFrm	20	Not in use	seconds
	UserUpdate	5	Not in use	minutes
	LogColor	● 0 - Black	Background color of the Log file: <ul style="list-style-type: none"> ● 0 = Black ● FF = Red ● FF00 = Green ● FFFF = Yellow ● FF0000 = Blue ● FF00FF = Magenta ● FFFF00 = Cyan ● FFFFFFFF = White 	hexadecimal RGB format
	IgnorePanelDifferences	0 - Default 1-9999	Default: Download of another CP type .cpf file is not allowed Any non-zero numeric value: Download of another CP type .cpf file is allowed	
	ApiSleep	100	Not in use	msecs
	SyncSleepTime	3000	Standard	msecs
	HideDeleteConfirmation	1	● 1 – True (Hidden) ● 0 – False (Visible)	0/1

Category	Parameter	Values	Description	Data UM
[MODEM]	KeepAliveInterval	5	For remote connections, CMS sends "Keep Alive" messages to the CP to keep the connection open.	seconds
	UserActionTimeout Note: <ul style="list-style-type: none"> ● Warning issued: "No user actions were made in the last x-1 minute (sec). Link will disconnect within 1 minute. ● Before connection is closed: "Link will close automatically." 	4	When no data has been sent between CMS and the CP for n minutes, other than keep-alive messages, CMS closes the remote connection.	minutes
[RemoteBoot]	SizeOfBootBlock	512	Standard – maximum packets Number of packets CMS sends during firmware file update before requesting the packets' status	packets
	BootDelayGSM	10	Packet transmission delay time – Not in use	msecs
	BootDelayPSTN	10	Packet transmission delay when updating a file using PSTN modem	msecs
	BootDelaySERIAL	35	Packet transmission delay when updating a file using local serial modem	msecs
	BootDelayIP	50	Packet transmission delay time – Not in use	msecs
	BootDelayGPRS	50	Packet transmission delay time – Not in use	msecs
	BootDelayAmigo	50	Packet transmission delay when updating a file over a CMS Amigo connection	msecs
	GSMDelayType	1	Not in use	1/2
[PSTN_Double_Ring]	DelayBeforeRing Note: Used to connect CMS to panels sharing a landline with non-CP devices – Prevents CP from answering first to calls for non-CP devices	13000	Required delay for the CP to receive two rings from CMS before disconnecting, and to respond automatically on next CMS call.	msecs
	WaitWindow	25	Delay for CMS second call	seconds
	DisconnectPolicy	1	<ul style="list-style-type: none"> ● 0 = Close COM port ● 1 = Hang up (ATH) command 	0/1
[GSMDData]	DelayBeforeFirstBuffer	7000	Not in use	msecs
	DelayBetweenBuffers	20	Not in use	msecs
[IP]	DelayBeforeFirstBuffer	7000	Not in use	msecs
	DelayBetweenBuffers	150	Not in use	msecs
[Warehouse]	DefaultChangeAccount	1	Not in use	
	DefaultPort	COM1	Default COM port to use in CMS warehouse	
	MruMax	5	Not in use	

Category	Parameter	Values	Description	Data UM
	ShowPanel	1 - Update	Whether CMS warehouse displays the kit management screen after successfully updating system configurations and files for current CP: <ul style="list-style-type: none"> ● 0 – CP update ends, button label Next Panel ● 1 – Manage CP devices before ending CP update, button label Update 	1/0
	showCMSFormAndIDOrigin	0	Not in use	
	rsu	1	Not in use	
	LogPath	Default=null	Path to CMS warehouse log files and output CPF files Default path: (%appdata%\Essence\Care@Home\Warehouse\Logs)	
	SelectVoice	0	Not in use	
	SelectFirmware	0	Not in use	
	ShowSimNumber	0	Whether the phone (SIM or landline) number is available for configuration	1/0?
	AutoIncrementAccount	0	Whether CMS warehouse increments (plus 1) the CP account number for subsequent successful CP updates	1/0?
	AllowEmptyAccount	1	Whether CMS warehouse assigns an account number to a CP without an instruction to do so: <ul style="list-style-type: none"> ● 0 – assign account number ● 1 – only with instruction 	1/0?
	ShowCombineDnisAccountWarning	1	Whether to check if CP is configured to combine DNIS with the CP account number <ul style="list-style-type: none"> ● 0 – No ● 1 – Yes 	1/0
	EnableChangeSIM	0	Whether the default should be to change the SIM (marked) <ul style="list-style-type: none"> ● 0 – Clear ● 1 - Marked 	1/0
	EnableChangeAccount	0	Whether the default should be to change the CP account number (marked) <ul style="list-style-type: none"> ● 0 – Clear ● 1 - Marked 	1/0
	EnableConfiguration	1	Not in use	
	EnableVerify Note: Verifying the burned FW doubles the FW update time	0	Whether the default should be to verify that the burned FW matches the FW file content (marked) <ul style="list-style-type: none"> ● 0 – Clear ● 1 - Marked 	1/0
	EnableVoice	0	Not in use	
	EnableFirmware	1	Not in use	
[PSTN]	DelayBeforeFirstBuffer	7000	CMS Wait period between establishing a PSTN modem connection and sending the first command to the CP.	msecs
	DelayBetweenBuffers	1	CMS wait period before each command is sent to a CP on a PSTN connection.	msecs

Category	Parameter	Values	Description	Data UM
[barcode]	defaultUseInPERS	1	Whether the barcode reader should be used by default when adding a device on a PERS panel ● 0 – No ● 1 – Yes	1/0
	defaultUseInAdlFamily	0	Whether the barcode reader should be used by default when adding a device on a Family/PERS-E panel ● 0 – No ● 1 - Yes	1/0
	defaultUseInAdlPro	0	Whether the barcode reader should be used by default when adding a device on a Pro panel ● 0 – No ● 1 - Yes	1/0

For example, to configure the CMS and CMS warehouse to use a barcode reader as the default for adding devices, edit your **setup.ini** file, using a text editor in administrator mode.

To edit the configuration file:

1. Open the **setup.ini** file in a text editor in administrator mode.
2. Change the barcode parameter: **defaultUseInPERS/AdlFamily/AdlPro** = 1, depending on the service package.
3. Save the changes.

2.3.3 Understanding the *amigo.ini* File

The **amigo.ini** file assigns the parameter values which configure:

- Amigo communication connections
- Remote connections that work either with SMS or with API.


The following table defines the parameters included in the **amigo.ini** file:


Parameter	Communications	Description or Value	Data UM
waitForPgToFinish	FTP by SMS	Note: When sending commands directly by SMS: Wait for the CP to download a file by FTP and place an acknowledgment file on the FTP server (PG command).	Time (msecs)
waitForPfToFinish	FTP by API	Note: When sending commands via API Wait for CP to download a file by FTP and send an event to the Care@Home™ server with the result (PF command).	Time (msecs)
panelTimeout	N/A	After CP successfully connects to CMS Amigo and no conversation has occurred within a pre-defined period, timeout takes effect. Default time period: 6 minutes	minutes
waitForPanelToCall	Amigo	Wait time:30 seconds - how long the CMS waits for the CP to establish connection	seconds

Parameter	Communications	Description or Value	Data UM
amigoPort	Amigo	Port # 3000 Note: This TCP port must be open on your firewall	
amigoAddress	Amigo	The external IP address or domain name of the computer where the CMS is located	
amigoPath	Amigo	64 bit: c:\Program Files (x86)\Essence Security\CMS Amigo\ 32 bit: c:\Program Files\Essence Security\CMS Amigo\ 	
waitBeforeStartingAmigo	Amigo	1000 milliseconds - default	
source	SMS	One of the source codes assigned to the Essence aerialink account: "Please enter source number for SMS provider"	
smsPassword	SMS	The password provided for the aerialink service	
smsUsername	SMS	The username provided for the aerialink service	
smsSPLink	SMS	URL address to the aerialink service	
noCapi=1	Connect API	0 for using cApi, 1 for using SMS direct	0/1
backofficePassword	API	"Please enter backoffice password"	
BackofficeEmail	API	"Please enter backoffice username"	
capiTimeout	API	50	seconds
cApiUrl	API	http://localhost/Connect.API/Commands/PanelSessionStart.json	

For example, to establish a remote cellular connection with SMS, your **amigo.ini** file must include the SMS-related parameters.

To establish a remote cellular connection with SMS:

1. Click . The Care@Home™ **Connect** window appears on the CMS home page.
2. Select **COM16** for the amigo remote cellular connection protocol.



3. Enter the CP international SIM phone number, without leading zeroes or a "+".

4. Enter the four digit CP DTMF code. The default code is "1234".
5. Click . The **Panel Info** window appears.
6. Click to upload the CP data to CMS.

```

amigo.ini

[all]
cApiUrl = http://localhost/Connect.API/Commands/PanelSessionStart.json
capiTimeout = 50 ; seconds
backofficeEmail = CmsOperatorUser
backofficePassword = 11111111111111111111
smsSpLink = https://apix.aerialink.net/v4/messages/
smsUsername = <aerialink user>
smsPassword = <aerialink password>
source = 1111111111
amigoPath = C:\Program Files (x86)\Essence Security\CMS Amigo
amigoAddress = 222.99.222.111 ; Domain|IP
amigoPort = 3000 ; port
waitBeforeStartingAmigo = 1000 ; milliseconds
waitForPanelToCall = 30 ; seconds
panelTimeout = 6 ; optional minutes
noCapi = 1 ; 0 for using cApi, 1 for using SMS direct
waitForPftoFinish = 5000 ; Time in milliseconds for PF command to be sent and finished
waitForPgtoFinish = 5000 ; Time in milliseconds for PG command to be sent and finished
    
```

The configuration file *config.ini* includes definitions for the amigo connection:

```

[Config]
SERIAL1=COM6->Cable
SERIAL2=COM2->Cable
SERIAL3=COM1->Cable
SERIAL4=COM5->Cable
SERIAL5=COM99->TACTICAL
SERIAL6=COM16->Amigo
[SERIAL6]
Type=0 ;*** 0 = SERIAL , 1 = MODEM ***
ComNumber=16
Baud=115200
DataBits=8
Parity=0
StopBits=1
Channel=4 ;*** 4 - Amigo ***
    
```



NOTE: All functionality is available using the remote connection. However, the response time for this type of connection is slower than direct connection types.

2.3.4 Understanding the *DTMFConfig.ini* File

The DTMF configuration file (*DTMFConfig.ini*) defines the default DTMF codes for each of the DTMF templates in the USA and Europe.

Parameter Values

Digits	DTMF Key code representation
1-9	DTMF key codes 1-9
10	DTMF key code 0
11	DTMF key code *
12	DTMF key code #
0	Default DTMF Key

The following table defines the parameters in the *DTMFConfig.ini* file.

Section	Parameter	Description/Action
[GENERAL]	default = USA	Default is DTMF codes for the USA
[USA/EU/UMO]	Half Duplex Speaker: ● Operator speaks ON ● Microphone OFF	Talk
	Half Duplex Mic: ● Operator listens ON ● Speaker OFF	Listen
	Full Duplex Inc Speaker	Increase speaker volume: 1 - 5 levels
	Full Duplex Dec Speaker	Decrease speaker volume: 1 - 5 levels
	End Call	Disconnect the call
	Alternate End Call	Alternate disconnect (UMO only)
	Reset Panel	Note: PSTN CP only Perform reset to the panel after 5 presses
	Switch to Next Voice Device	Next voice extender device
	Switch to CP	Voice extender switch to CP
	Go to Full Duplex	Switch to full duplex: talk and listen

```
*DTMFConfig.ini - Notepad
File Edit Format View Help
[GENERAL]
default=USA
[USA]
;Possible values
HalfDuplexSpeaker=1
HalfDuplexMic=3
FullDuplexIncSpeaker=2
FullDuplexDecSpeaker=10
EndCall=9
AlternateEndCall=0
ResetPanel=8
SwitchToNextVoiceDevice=6
SwitchToCP=5
GotoFullDuplex=4
[EU]
HalfDuplexSpeaker=8
HalfDuplexMic=7
FullDuplexIncSpeaker=3
FullDuplexDecSpeaker=1
EndCall=10
AlternateEndCall=0
ResetPanel=2
SwitchToNextVoiceDevice=6
SwitchToCP=5
GotoFullDuplex=4
[UMO]
HalfDuplexSpeaker=8
HalfDuplexMic=7
FullDuplexIncSpeaker=3
FullDuplexDecSpeaker=1
EndCall=11
AlternateEndCall=10
ResetPanel=2
SwitchToNextVoiceDevice=6
SwitchToCP=5
GotoFullDuplex=4
```

USA DTMF

EU DTMF

UMO DTMF


Refer to 5.4.7 DTMF Config on page 60.

3. Understanding the CMS Connection

The **Connect** window enables the CMS connection as follows:

- To allow the technical support personnel to perform maintenance and error recovery
- To enable support-level access to configure the CP parameters
- To upload configuration parameter values and monitoring information from the CP to Care@Home™ CMS
- To activate the Remote Boot feature: Refer to 16 Remote Boot Software Updates on page 182.

The **Connect** window appears automatically when you open Care@Home™ CMS.



The following table describes the **Connect** window parameters with instructions how to apply them.

Name	Definition and Instructions
User	Enter your support-level username to enable configuration of support-allowed CP parameters. Note: This is for Essence technical support use.
Password	Enter your support-level password to enable configuration of support-allowed CP parameters. Note: This is for Essence technical support use.
Protocol	Upload Select this protocol to establish a connection to the CP to maintain the CP settings.
	Remote Boot Select this protocol to activate the Remote Boot feature to update the software and firmware of the CP and the peripheral devices. NOTE: The Remote Boot window appears when this option is selected.
	Remote Software Update using FTP NOTE: This option is only valid for Care@Home™ PERS.
Connect	Select the connection method. If you are using a cable connection, select a COM port.

Name	Definition and Instructions
Phone	Enter the phone number only for a remote connection.



NOTE: Consult Essence professional services to receive your support-level user credentials.

The following table describes the functionality of the action buttons at the bottom of the **Connect** window.

Action	Instructions
Connect	Click to establish a connection with the CP according to the protocol selected.
Abort	Click to abort the connection process. The button is disabled until the connection process begins.
Clear	Click to clear the message field of all displayed information, such as system messages and error messages.
Exit	Click to exit the Connect window. Care@Home™ CMS enters “offline” mode. Offline mode allows you use the features in the File Menu, such as open and view the .cpf file. Refer to 4.2 File Menu Functions on page 39.

3.1. Establishing a Local Connection

The Care@Home™ CMS can connect to the CP using a local COM port.

3.1.1 Local Connection to the CP

To connect the CP to the PC, use an Essence ES7000BACMS cable with a serial COM port (RS232). Contact Essence professional services personnel for information about this cable.



When there is no RS232 COM port available on the PC, use an USB-RS232 adapter, such as:

- ESCL05019 mini-USB adapter

- USB Gear Serial adapter

Refer to Appendix A USB to RS232 Adapters on page 199.



NOTE: The CP must be turned on **before** connecting the cable to the RS232 COM port on the local computer.

For the RS232 Com port connection, there is a set of beeps that are sounded for different connect or disconnect actions as follows:

- A beep emitting a long rising pitch indicates a successful plug-in connection.
- A beep emitting a double low octave pitch indicates an unsuccessful plug in.
- A beep emitting a long falling pitch indicates a successful disconnect.

When the CP beeps indicate a “disconnect” or an unsuccessful plug in, you must re-connect to the CP.

Refer to 3.3.2 Reconnecting to the CP on page 35.

3.2. Establishing a Remote Connection

For a remote a connection, the CP can use the following communication channels:

- Cellular and Ethernet channels use Amigo communications application
- PSTN with modem connection



NOTE: All functionality is available using the remote connection. However, the response time for this type of connection is slower than direct connection types.


The user should use the Windows **Remote Desktop** utility to connect remotely to the CMS.

The Windows **Remote Desktop** connection requires the user to connect to a server and to enter user credentials for accessing the server via Amigo communications.




NOTE: Contact Essence professional services to receive the server IP address and user credentials.

After establishing the remote connection with the server, the **Remote Desktop** screen displays the Care@Home™ CMS shortcut. Double-click to open CMS.



The **Connect** window appears.




The following table provides instructions for configuring the parameters on the **Connect** window for a remote connection.

Name		Definition and Instructions	Data Values	Required/Optional	Default Value
User		For service provider support user or Essence technical support use only		Optional	Blank
Password		For service provider support user or Essence technical support use only		Optional	Blank
Protocol: Upload		Select this option to establish a connection with the CP to maintain the panel settings.	Radio group	Optional	Selected
Protocol: Remote Boot		Select this option to activate the Remote Boot feature to update the CP and peripheral firmware. NOTE: The Remote Boot window appears when this option is selected.	Radio group	Optional	Not Selected
Connect		Select the remote connection method via either: <ul style="list-style-type: none"> ● Amigo communications application ● PSTN with a modem 	Dropdown List	Required	Amigo: COM16 PSTN: COM9
Remote connection via	Panel ID	Enter the account number assigned to the CP	Up to 8 digits	Required	None
	Phone number	For SMS	19 chars		


3.2.1 Establishing Remote Connection with API

Perform the following procedure to establish a remote connection from the CP to the CMS, after connecting and logging onto the remote server via Remote Desktop:


1. Select the COM16 port (default) for the Amigo remote connection.
2. Enter the CP account number.
3. Click .

System messages are displayed in the message window, describing the status of the connection process.

The figure below illustrates the messages displayed that explain that the system opens the Amigo virtual port and calls the CP using the Connect API Web Service. When the CP responds, the system displays messages that the “COM is OK” and that the system is verifying the communications connection.



When the connection is established, the **Panel Info** window appears displaying the resident’s current status in the window title of the **Panel Info** window.

In the **Panel Info** window, click , to upload the CP data to CMS.

Refer to 3.3.1 Uploading CP Settings Procedure on page 32.



NOTE: All functionality is available using the remote connection. However, the response time for this type of connection is slower than direct connection types.

3.2.2 Establishing Remote Connection with SMS


A remote connection with SMS is used when the CP is not yet registered on the Care@Home™ server. This remote connection is also used with registered CPs.

Remote connection with SMS is used depending on the service package assigned to the CP. CPs with the PERS service package connect to the CMS via remote connection with SMS. CPs with either PERS-E, Family, or Pro service package connect to the CMS via remote connection with either SMS or API.


Before establishing a remote cellular connection with SMS, ensure that the following files are located on your computer, depending on your computer operating system:

- The amigo configuration file (**amigo.ini**) is located at the following path:
 - 64-bit OS: **C:\Program Files (x86)\Essence\Care@Home CMS**

1. Plug the modem into an available USB port on the local PC. The Windows **Device Manager** identifies the device as a modem.




2. Double-click on the U.S. Robotics line. The U.S. Robotics **Modem Properties** window appears.
3. Click the **Modem** tab. The COM port to which the modem is assigned is displayed. In the example, the **Device Manager** assigned the modem to **Port "COM9"**.




4. To change the COM port:
 - a. Click the Advanced tab in the U.S. Robotics Modem Properties window. The Advanced window appears.

- b. Click Advanced Port Settings. A window appears displaying the Advanced Settings for COM9.



- c. Select the COM port you want to assign to the modem from the **COM Port Number** dropdown list.



- d. Click **OK**.

3.2.3.2 Configuring Landline Modem COM Port

Continuing with the example of establishing the remote PSTN connection to the CP, COM9 is the COM port used by the PSTN landline modem.

To edit the configuration file, **config.ini**, to assign COM9 to the modem:

1. Open the configuration file, config.ini, in a text editor.
2. For [Config], modify: MODEM4=COM9->MODEM PSTN
3. For [MODEM4], modify: ComNumber=9
4. Save the file changes.

The figure below shows the configuration file, **config.ini**, revised to assign the COM port, **COM9**, to the landline PSTN Modem.

```


*Config Robots.ini - Notepad
File Edit Format View Help
[ALL]
TRACELEVEL=1
useXML=0
ShowStatusBar=0
useBarcode=1
[Config]
SERIAL1=COM1->Cable
SERIAL2=COM2->Cable
MODEM1=COM1->MODEM EICON GSM
MODEM2=COM2->MODEM EICON PSTN
MODEM3=COM1->MODEM SIEMENS GSM
MODEM4=COM9->MODEM PSTN ; USRobotics ←
[SERIAL1]
Type=0 ;*** 0 = SERIAL , 1 = MODEM ***
ComNumber=1
Baud=9600
DataBits=8
Parity=0
StopBits=1
.
.
.
[MODEM4]
Type=1 ;*** 0 = SERIAL , 1 = MODEM ***
ComNumber=9
Baud=115200
DataBits=8
Parity=0
StopBits=1
Channel=1 ;*** 0 = GSM , 1 = PSTN ***
ATCommands=\\MODEM_PSTN_USRobotics.ini ←
    
```

3.2.3.3 Connecting the Remote Landline Modem

Continuing with the example of establishing the remote PSTN connection to the CP, connect the CP to the remote PSTN Landline modem.



To remotely connect the Care@Home™ CMS to the CP:



1. Double click . Care@Home™ CMS opens and the **Connect** window appears.
2. Click the dropdown list of the **Connect** field.
3. Select the COM port that connects to the PSTN modem. In the example, the COM port for the PSTN Modem is COM9.
4. Enter the telephone (landline) number of the remote PSTN landline CP. Enter the phone number in international telephone number format:

Country code||Area code||Phone number

5. Click Connect...



3.3. Setting up CP Connectivity

Connect the Care@Home™ CMS to the CP, to use CMS to manage your Care@Home™ CP and peripheral devices functionality.


3.3.1 Uploading CP Settings Procedure

To establish a connection to the CP:




1. Double click . The CMS home page opens and the **Connect** window appears.
2. Select **Upload**.
3. Select the COM port according to the connection method required:
 - **Local connection:** refer to 3.1 Establishing a Local Connection on page 23.
 - **Remote Cellular connection:** refer to 3.2 Establishing a Remote Connection on page 24.
4. Enter a phone number, if required by the connection method chosen.
5. Click . When the connection is established, the **Panel Info** window appears.
6. Click to upload the CP data to CMS.

Following the successful completion of the connection procedure, the **Panel Info** window appears.



The following table describes the information and functionality of the **Panel Info** window.

Name	Data/Buttons	Definition or Instructions	
Firmware	Serial	Serial number of the CP	
	Main CPU	Firmware version of the CP	
Actual Panel Time	Date	The date the CP connected to CMS, in DD/MM/YYYY format	
	Time	The time the CP connected to CMS, in HH24:MI:SS format	
Commands	Radio button commands	Test Mode	For service provider system administrator or Essence technical support use only
		Operational (default)	Panel is in normal operation mode.
	Push button commands	Reset Panel	Click to restart the CP. The restart occurs at the end of the current session.
		Refresh	Click to load new status information of the CP to the Panel Info window.
		Click to upload the CP configuration parameter values to the CMS from the CP. The Remote Panel window appears displaying a progress bar of the data transfer process.	

Name	Data/Buttons	Definition or Instructions
GSM Coverage		Current strength of the GSM reception signal, relevant to 3G, 2G, and LTE models


When the connection is established, the connection information is displayed on the left-most part of the CMS home page.

Reset Panel

Account Number
25250001

Panel Firmware
Serial: 703F56C3
CPU: 2.07.1.230.06.04

Panel Connect Time
Date: 22/08/2018
Time: 13:44:00

Gsm Coverage
 100%

Panel hardware
Version: 2.A

<= Panel CPU Version =>

Panel Type

<= ES6502HC


C7000 =>

Reset Panel

Account Number
00008156

Panel Firmware
Serial: N/A
CPU: 6.01.1.231.05.08

Panel Connect Time
Date: 26/07/2018
Time: 13:01:00

Gsm Coverage
 100%

Panel hardware
Version: C7000-v1.1, 2.B


The **Panel Connection** data allows you to identify which type of panel you have, whether ES6502HC or C7000 and which panel CPU version is installed.



NOTE: Remember: When availability criteria are noted for a feature:

- If a panel type is **not** mentioned, the feature is fully supported by the panel type.
- If a panel type is mentioned, compare the panel version listed to the **Panel CPU Version** displayed with the **Panel Connection** data.

Along with the service package, these criteria can help you understand if a specific feature or functionality is available on your Care@Home™.

Name	Definition and Instructions
Reset Panel	Click to restart the CP. The restart occurs at the end of the current session.
Account Number	Identification number currently assigned to the CP
Panel Firmware	Serial ES6502HC/ C7000 panel serial number
	CPU ES6502HC / C7000 panel firmware version (Panel CPU Ver) NOTE: Where the panel version criteria is listed, to determine if a feature or function is included in your Care@Home™, compare the criteria to this parameter (Panel CPU Version or Panel CPU Ver)
Panel Connect Time	Date The date the CP connected to CMS, in DD/MM/YYYY format
	Time The time the CP connected to CMS, in HH24:MI:SS format
GSM Coverage	 Current strength of the GSM reception signal, relevant to 3G, 2G, and LTE models
Panel Hardware	Hardware revision: If the panel is a C7000, the “C7000” should appear as the prefix to the hardware version

You can update the configuration parameter values of the CP by:

- Clicking **Reset Panel** on the **Panel Info** window
- Clicking **Reset Panel** on the left side of the CMS home page
- Selecting the **Open** function in the **File** menu on the CMS home page

The **Reset Panel** function requests restarting the CP when the CMS session ends.

At the end of the CMS session, the CP integrates the updated configuration parameter values into the CP functionality.

3.3.2 Reconnecting to the CP

When the local connection to the CP fails, you must re-establish the connection to the CP.

There are several instances when it is necessary to re-connect to the CP:

- The CP is not connected to a working USB cable.
 The messages in the **System Messages** window show the communication checking process. This process terminates with an **Authentication Time-Out** message.
 CMS issues the **Authentication Time-Out** message.
- The user does not wait the acceptable time-period before attempting to reconnect to the CP.

As in the previous instance, the connection process terminates with an **Authentication Time-Out**. CMS issues an error message.

- The COM port specified is being used by another device.

The COM port chosen in the **Connect** window, from the **Connect** dropdown list, is unavailable.


The messages in the **System Messages** window show that the attempt to access the COM port failed.

- The COM port chosen in the **Connect** window is invalid.

The COM port chosen in the **Connect** window, from the **Connect** dropdown list, is incorrectly defined in the configuration file.


As in the previous instance, the messages in the **System Messages** window show that the attempt to access the COM port failed.

To reconnect the CP to the local PC:

1. Click .
2. Remove the USB cable from the USB connection on the CP and **wait at least 40 seconds**. This allows the CP to reset itself.



NOTE: It is crucial to wait at least 40 seconds.


3. Reattach the USB cable to the USB connection on the CP. A successful plug-in beep is sounded. If the beep is not sounded, try again.
4. For information about the CP beeps, refer to 3.1.1 Local Connection to the CP on page 23.
5. Click .

4. Getting Started




Double click . The Care@Home™ CMS home page appears.








The Care@Home™ CMS home page provides a **File** menu, action icons, and the **Main** navigator menu to operate within the application.




4.1. Main Action Icons

The following is a list of the action icons on the CMS toolbar, along with the functionality description of each icon.


Action Icon	Action Definition and Instructions
	Exit the Care@Home™ CMS application.

Action Icon	Action Definition and Instructions
	CMS properties: Revising COM port properties is reserved for Essence technical support personnel use only.
	Disconnect from the CP.
	Connect to the CP: Opens the Connect window. Refer to 3 Understanding the CMS Connection on page 22.
	Download the configuration parameter settings from CMS to the CP.
	Upload the configuration parameter settings from the local configuration file (.cpf file) to CMS.
	Refresh: Load new status information to CMS.
	Monitor activity sensors, safety devices, and emergency devices connected to the CP.




Click  or the **Monitor** sub-menu on the menu bar to access the monitoring windows. The monitoring windows display the status information for the devices connected to the CP:

- The **Activity Sensors - Monitor** window is a consolidated view of the information about the sensor devices as displayed in the **Activity Sensors** module.



ID	Serial Number	Device Type	Label	BYPASS	OPEN	BATTERY	TAMPER	SUPERVISO...	RF Level (%)
1	AB928457	IPD	Living Room	ACTIVE	OK	OK	OK	OK	80
2	18493751	IPD	Dining Room	ACTIVE	OK	OK	OK	OK	80
3	A298347C	PIR	Childrens Room	ACTIVE	OK	OK	OK	OK	80
4	38947538	MGLS	Main Entry	ACTIVE	OK	OK	OK	OK	80
5	34975394	IPD	Courtyard	ACTIVE	OK	OK	OK	OK	80
6	23896419	MGLS		ACTIVE	OK	OK	OK	OK	80
7	0745B488	MGLS		ACTIVE	OK	OK	OK	OK	80

- The **Safety & SOS - Monitor** window is a consolidated view of the information about the safety and emergency devices as displayed in the **Safety and SOS** module.




NOTE: The statuses are either OK, KO, ACTIVE, or BYPASSED except the battery status. Values for the battery status are:



- OK
- LOW
- EMPTY

4.2. File Menu Functions



The following table explains the functions in the **File** menu which allow management of the CMS parameters file (**.cpf** file) for the Care@Home™ CMS.

Name	Definition and Instructions
Open	Locate and “open” the local .cpf file on the PC to be uploaded to CMS, using the Windows CMS File Dialog box.
Save	Save the current configuration parameter settings to a local .cpf file using the same file name. If the .cpf file is a new file, use the Save As... function.
Save as...	Save the current configuration parameter settings to a new local .cpf file with a new file name.
Print...	Produce an external file or hardcopy of the configuration parameter settings. Refer to 4.2.3 Print... Function on page 40.

Name	Definition and Instructions
Print Log...	Refer to 13.7 Print Log on page 162.
Print Setup	Set the printer settings to produce a hardcopy of the configuration parameter settings.
Quit	Select this function to exit Care@Home™ CMS. A window appears to confirm the exit request. Refer to 4.2.2. Quit Function on page 40.

4.2.1 Open Function

The **Open** function not only locates and opens the **.cpf** file but also automatically uploads the configuration parameter settings to the CP, at the end of the session.

4.2.2 Quit Function

The **Quit** function verifies that the changes made during the CMS session were updated and/or saved, before allowing you to exit CMS.

To exit CMS:

1. In the **File** menu, click **Quit**.
2. If your CMS is connected to a CP, a message appears prompting you to confirm the exit request.
3. Click **Yes** to exit CMS.
4. If open issues remain, the application displays several caution messages to prompt you to confirm that you are aware of these issues.
5. Click **OK** to close these messages to exit CMS.

4.2.3 Print... Function

The **Print...** function opens the Print Panel Configuration window. This function provides the tools to choose the media format to print the CP data. You can choose whether to produce an Excel file, a **.csv** (comma separated values) file, or a hardcopy of the configuration parameter settings.


4.2.4 Print Log.. Function

The **Print Log...** function has same functionality as in the **Log Events** module. Refer to 13.7 Print Log on page 162.

4.3. The Main Navigator Menu

The **Main** navigator menu provides a guide to access the modules in the CMS application.

The following is a list of the CMS application modules, including a description of each module's purpose:



- **Control Panel** – This module provides the tools to configure the Care@Home™ CP parameters.
Refer to 5 The Control Panel Module on page 43.
- **Communications** – This module allows for configuration management of multiple communication channels, allowing for remote activation or redundancy in case of channel failure or other connectivity issue.
Refer to 6 The Communications Module on page 79.
- **Dialer** – This module provides the tools to define the communication methods, such as phone numbers, IP addresses, and ports, and to manage the flow required to communicate with a monitoring station.
Refer to 7 The Dialer Module on page 92.
- **Activity Sensors** – This module provides the tools to manage and control the activity sensors installed on the resident's premises.
Refer to 8 The Activity Sensors Module on page 122.
- **Safety and SOS** – This module provides the tools to manage and control the safety and emergency devices installed on the resident's premises.
Refer to 9 The Safety & SOS Module on page 129.
- **Voice & LED** – This module provides the operator with a tool to define which type(s) of announcements and reminders are to be used, to manage the vocal announcements and reminders for the resident, and to control the CP LED indicators.

The reminders announce:

- When the residents should take their medication
- The schedule time of an appointment
- The time transportation should arrive to take a resident to an appointment or meeting

- Wake-up calls like an alarm clock
- Random requests from the operator for the residents to test their emergency buttons

Refer to 10 Voice & LED Module on page 139.

- **Mobile Operators** – This module provides the tools to provide APN information for the SIM card in the Care@Home™ CP.

Refer to 11 Mobile Operators on page 152.


- **Custom Labels** - This module provides the tools to define and manage the personalized labels, defined by the user, to label areas on the resident's premises not included in the list of possible areas provided by the application.

Refer to 12 The Custom Labels Module on page 157.


- **Log Events** – This module provides the tools to view the log of events and activities recorded by the CP for a given time-period.

Refer to 13 The Log Events Module on page 159.

5. The Control Panel Module



Click **Control Panel** on the **Main** navigator menu. The **Control Panel** module is the tool that allows the configuration of the settings for the Care@Home™ CP installed on the resident’s premises.



5.1. Account Number

The **Account Number** is a unique 8-digit integer, used as the required CP identification number, assigned to a CP by the monitoring station.



NOTE: When the GSM configuration parameter, **Combine DNIS with Account No**, is marked, the account number is a concatenation of the DNIS with up to the last five digits of the assigned number.

5.2. ARC Account Number

The ARC account number is a unique 16-digit integer. It is an optional customer-defined identifier. The ARC account number is sent to the monitoring station using the following protocols:

- SCAIP – ADL version 2.5.14 and later
- CPC – PERS version 2.6 and later

5.3. Service Type

A Care@Home™ service type is the product offering of Care@Home™ that provides a level of security and monitoring for the resident.

The supported Care@Home™ service types by panel and software version are:

Service Type	ES6502HC		C7000	
	Earlier than Panel CPU 2.05	Panel CPU 2.05 or later	Prior to Care@Home 6.02	Care@Home 6.02 and later
PERS	√	√	√	√
PERS-E				√
Pro	√	√		√
Family		√		√

Select the service type assigned to the CP. CMS activates the functionality specific for the service type selected.



NOTE: For ES6502HC, select **Family** for the **PERS-E** service type and disable Report Interval.

When uploading CP parameter values from a local **.cpf** file or from a CP connected to the CMS and the parameter values are defined for a specific service type, the **Service Type** parameter is automatically displayed and the relevant Care@Home™ functionality activated. The CMS interface changes accordingly.


Click  to display additional CP configuration parameters.

5.4. More Options

The CP extended configuration includes standard and support-level parameter categories. **More Options** displays the additional developer-level categories if you connected to the CP with developer credentials.



NOTE: For ES6502HC, access to **More Options** is available for CP software version 2.03 or later.



The following table lists the standard and developer-level parameter categories.

Tab name	Access	Description	Refer to	Tab Visible			
				ES6502HC		C7000	
				PERS	Pro& Family	Not PERS or PERS-E	All
Photo Config	Standard	Definitions for the camera sensors configuration	5.4.1 Photo Config on page 47		√	√	
Out of Home	Support	Definitions of Out of Home monitoring and alert notification	5.4.2 Out of Home on page 49		√	√	
Analytics	Standard	Other additional configuration parameters	5.4.3 Analytics on page 51		√	√	
Periods	Support	Definitions of periods for rules' definitions	5.4.4 Periods on page 52		√	Not PERS or PERS-E	

Tab name	Access	Description	Refer to	Tab Visible			
				ES6502HC		C7000	
				PERS	Pro& Family	Not PERS or PERS-E	All
Rules	Support	Rules' criteria settings	5.4.5 Rules on page 54		√	Not PERS or PERS-E	
PERS	Standard	Additional communications protocol, CP battery management, and debug option for technical support	5.4.6 PERS on page 58	Panel CPU Ver 2.03 and later	Panel CPU Ver 2.04 and later		√
DTMF Config	Standard	Definitions for DTMF keys configuration	5.4.7 DTMF Config on page 60	Panel CPU Ver 2.03 and later	Panel CPU Ver 2.04 and later		√
Log	Standard	Log transmission and transfer parameters	5.4.8 Log on page 63		Panel CPU Ver 2.05 and later & GSM-enabled		√
VPD Config	Standard	VPD configuration parameters at CP level, and a migration tool for full VPD functionality.	5.4.9 VPD Configuration Parameters on page 64	Panel CPU Ver 2.03 and later	Panel CPU Ver 2.04 and later		√
EPA Config	Standard	Configuration for EPA at CP level, for all EPAs in your Care@Home™	5.4.10 EPA Configuration Parameters on page 68	Panel CPU Ver 2.03.4 and later	Panel CPU Ver 2.05.6 and later		√
RF	Standard	For Essence technical support only	5.4.11 RF Settings on page 70		Panel CPU Ver 2.05.6 and later		√

To revise the CP configuration:

1. Click . The Control Panel Configuration Extension appears.



2. If you connected to the CP without support credentials, the **Control Panel Configuration Extension** appears with the standard tabs.


3. If you connected to the CP with support credentials, the **Control Panel Configuration Extension** appears with both the standard tabs and the developer-level tabs.
4. Refer to 3 Understanding the CMS Connection on page 22.
5. Click the tab you want to edit and enter the values for the parameters on the tab.
6. On each tab, click **Update** to add the revised configuration parameter values to the CP configuration.
7. Click **Close** or click X on window to exit the **Control Panel Configuration Extension**.

5.4.1 Photo Config



NOTE: The **Photo Config** parameters are available for Family and Pro only

The **Photo Config** tab allows you to configure the camera detectors (IPDs) which include infrared sensors.



When triggered, each of the IPDs photographs five consecutive photos of the triggering event. If both cameras are enabled, the photos are taken in the order defined by the configuration parameters. The purpose of the **Photo Config** tab is to configure this camera-related functionality.



NOTE: IPDs are optional.

If IPDs are installed, mark the relevant checkboxes to enable the following capabilities:

- **Alarm** – If enabled, the CP instructs the IPDs to take photos when an alarm is triggered. The default is disabled (clear).
- **Photo Request** - If enabled, the CP instructs the IPDs to take photos when an event is detected by any of the configured sensors. The default is disabled (clear).
- **SOS** - If enabled, the CP instructs the IPDs to take photos when an emergency event is triggered. The default is disabled (clear).

If IPDs are installed, define the **default order** in which the IPDs photograph events when the camera option is triggered by the CP or pendant.

- Camera 1 – first priority
- Camera 2 - second priority

The numeric value selected for the camera parameter is the device identification number assigned to an IPD in your Care@Home™ configuration.

If the camera parameter is set to **None**, the camera option is disabled. The default configuration is **None**.

You can override the default order per device in the Activity Sensors and the Safety and SOS modules.


Refer to 8 The Activity Sensors Module on page 122.

Refer to 9 The Safety & SOS Module on page 129.

5.4.2 Out of Home



NOTE: The **Out of Home** parameters are available only for service types Family and Pro.



The **Out of Home** tab allows you to configure the monitoring and alert notifications for the **Out of Home** event. The parameters that define an **Out of Home** event are used when applying the **Out of Home** non-standard rule type. This tab requires support-level access.



Caution: Set these parameters carefully to avoid false alerts.


The following table describes the **Out of Home** configuration parameters.

Parameter		Description	Values
<Out of Home feature activation>	Enabled	Mark to activate Out of Home .	Marked=Enabled Clear= Disabled
	Enabled Alert	Mark to activate the Out of Home alert.	Marked=Enabled Clear= Disabled
<Out of Home> Entry	Delay	Select the amount of time to wait, from the time the Entry event occurs, before activating detection.	Dropdown list: 0 - 50 minutes
	Duration	Select the amount of time to detect movement before deciding whether the resident is at home. Note: If the duration has expired and no detections occur, the resident is still Out of Home .	Dropdown list: 5 – 720 minutes

Parameter		Description	Values
	Threshold	Select the number of detections that must occur to trigger a Back at Home alert. Note: A Back at Home alert is triggered if the number of detections occurs, even if the duration time has not expired.	Dropdown list: 1 - 16 detections
<Out of Home> Exit	Delay	Select the amount of time to wait, from the time the Exit event occurs, before activating detection.	Dropdown list: <ul style="list-style-type: none"> ● ES6502HC ADL firmware, CP version 2.09.6 and later (except 2.14): 0-254 minutes in increments of 1 minute. ● In all other releases (including C7000): 0-50 in increments of 5 minutes
	Duration	Select the amount of time to detect movement before deciding whether the resident is not at home. Note: If there are no detections, the resident is Out of Home .	Dropdown list: 5 - 720 minutes
	Threshold	Select the number of detections that must occur to not trigger an Out of Home Exit alert.	Dropdown list: 1 - 16 detections

5.4.3 Analytics

The **Analytics** tab allows you to configure service level functionality, as needed.



The following table describes the **Analytics** configuration parameters.


Name	Definition and Instructions	Data Values	ES6502HC		C7000		
			PERS	Pro & Family	Not PERS	All	
Care@Home:	Report Interval	To configure the interval for reporting activities, select the reporting frequency in minutes Note: For PERS-E, this parameter can be disabled. On the C7000 the parameter is disabled by default.	Dropdown List "Default", 10 - 180 with increments of 10. For C7000, and ES6502HC with CPU version between 2.09.1 and 2.14, list includes "Disable"	N/A	√	√	
	Sleeping interval	Select the amount of time of the PIR "blocking window". This is the time the PIR "sleeps" between detections. Note: 10 seconds is recommended Note: Parameter appears for support-level user only	Dropdown List "Default", 10 - 300 with increments of 10. Default: 160-seconds	N/A	√	√	

Name	Definition and Instructions		Data Values	ES6502HC		C7000	
				PERS	Pro & Family	Not PERS	All
STUN Cellular	Enable	Mark to enable the STUN cellular communications type	Marked=Enabled Clear= Disabled	N/A	√	√	
	Interval	Select the interval in seconds	Integer (5): 25 - 43200 sec	N/A	√	√	
STUN Ethernet	Enable	Mark to enable the STUN Ethernet communications type	Marked=Enabled Clear= Disabled	N/A	N/A	√	
	Interval	Select the interval in seconds	Integer (5): 25 - 43200 sec	N/A	N/A	√	
Multiple Detectors:	Note: ES6502HC with panel software version 2.09.6 through 2.14 – Pro service package Mark to limit one detector per activity type		Marked=Enabled Clear= Disabled	N/A	√		

5.4.4 Periods





NOTE: The **Periods** parameters are available only for service types Family and Pro.



The **Periods** tab allows you to define the periods to use when defining rules. This tab requires support-level access.

You can mark a period of 24 hours or define the beginning and end of a period, up to 24 hours.

The following table describes all the parameters and functionality of the **Periods** tab.

Parameter	Description	Values	Notes
ID	The identification code for the period: The code assigned is the number of the row where the period is defined.	Display Only You can define up to 16 periods.	The Rule definitions refer to this identification code.
24 Hours	Mark to define a 24-hour period.	Marked = 24 hours	Default: Clear- no period
Period Start	Select the Period begin time.	Dropdown list: 00:00 – 23:45 in 15-minute intervals	
Period End	Select the Period finish time.	Dropdown list: 24-hours of 15-minute intervals beginning with the first 15-minute interval following the begin time	

Click **Clear** to delete the selected period definition. A prompt appears to confirm your request.


Click **Clear All** to delete all the period definitions. A prompt appears to confirm your request.

Click **Update** to save all the period definitions.

5.4.5 Rules



NOTE: The **Rules** parameters are available only for service types Family and Pro.



The **Rules** tab allows you to define the rules to trigger activity alerts. This tab requires support-level access.

When you define a rule, you assign the rule to a rule type. Each rule type has attributes that determine the rule’s behavior. The rule can trigger alerts. For a list of rules types see page 55.

For example, if a rule is assigned the “Wandering” type, the rule behavior definition includes:


- The **period** in which the activity of “opening of the entrance door” is tracked
- The monitoring of the activity during the **period** on each of the marked **days of the week**
- The **duration** after which an alert is broadcast if no activity is detected

There are rule types that are triggered by the number of PoSAs identified within a given period.

The **Pro** and **Family** service types interpret movement detections differently. **Pro** tracks resident “activities” whereas **Family** tracks resident “presence”. The algorithm determines “presence” by recognizing “sustained activity” detections.

A **Period of Sustained Activity** (PoSA) is a period during which an activity sensor detects multiple consecutive movements. For example, if the resident is in the dining room, a PoSA is recognized when the activity sensor in the dining room detects two or more movements of the resident at or in the vicinity of the dining room table.

The following table describes the parameters of the **Rules** tab.

Parameter	Description	Values	Notes
Rule ID	The identification code for the rule: The code assigned is the number of the row where the rule is defined.	Display Only	Maximum: 128 rules.
Period ID	Select the identification code of the time frame in which the rule should monitor for specific activity or no activity.	Dropdown list of Periods defined in the Periods .	Refer to 5.4.4 Periods on page 52
Type	Note: Option Possible fall is not available for Family service type Select the rule type associated with the rule being defined.	Dropdown list of the rule types.	For a list of rules types see page 55.
Group / Device	Choose either an activity group of devices or a specific device whose detections are used by the rule to trigger alert notifications	Options: ● Group ● Device Options are enabled based on the Type selected.	Group functionality allows you to manually define a group of devices for rule usage only.
Group ID / Device ID	Select activity group identification code Select device identification code	Dropdown list is determined by value of previous parameter.	
Activity	When Group is chosen, select the activity type assigned to the activity group.	Dropdown list of activities according to the default device configuration	Enabled only for Group
Value	The field changes based on the Type selected.	Dropdown list: ● Detections ● Durations (min) ● Durations (min) with days of the week (DoWs) ● Durations (sec) ● PoSAs ● Expiration (min)	Depends on the selected rule
<Days of the week>	Mark the days of the week to activate the rule.	Sunday through Saturday 	Shown only for rules Wandering and Unexpected Entry/Exit

Click **Clear All** to delete all the rule definitions. A prompt appears to confirm your request.

Click **Update** to save all the rule definitions.

The following describes the rule types and the parameters required for each rule type.

Name	Description	Group (Gr)/ Device (Dv)	Value
Excessive Activity	Note: For Family or for Pro with Panel CPU Version 2.4.5 and later, only with rule Wake & Well An alert is triggered when more than the expected number of activity detections is made, during the chosen period.	Gr or Dv If Gr, select Activity type	Detections (1-255)

Name	Description	Group (Gr)/ Device (Dv)	Value
Extreme Low Activities	Note: Family only An alert is triggered when the resident is on the premises and no activity detections are made for longer than the expected duration, during the chosen period.	N/A	Duration (min)
High Number of PoSA	Note: Family only An alert is triggered when more than the expected number of PoSAs is identified, during the chosen period.	Gr or Dv If Gr, select Activity type	PoSAs (1-255)
High Sum of PoSA Durations	Note: Family only An alert is triggered when the total duration of the PoSAs identified during the chosen period exceeds the expected duration.	Gr or Dv If Gr, select Activity type	Duration (min)
Low Activity	An alert is triggered when the resident is on the premises and the number of activity detections made, during the chosen period, is lower than the expected number of detections.	Gr or Dv If Gr, select Activity type	Detections (1-255)
Low Number of PoSA	Note: Family only An alert is triggered when less than the expected number of PoSAs is identified, during the chosen period.	Gr or Dv If Gr, select Activity type	PoSAs (1-255)
Low Sum of PoSA Durations	Note: Family only An alert is triggered when the total duration of the PoSAs identified during the chosen period is below the expected duration.	Gr or Dv If Gr, select Activity type	Duration (min)
Magnet Open Too Long	An alert is triggered when a door on the premises is open for longer than the expected duration, during the chosen period.	Dv	Duration (min)
Possible Fall	Note: Pro only An alert is triggered when the resident is on the premises and no activity detections have been made in the expected duration following the last activity detection, during the chosen period.	Dv	Duration (min)
Short Duration Between Extremes	Note: Family only An alert is triggered when the resident is on the premises and during the chosen period, either: <ul style="list-style-type: none"> ● No activity detections have been made by a specific device -or- ● The time between detections, following the last activity is shorter than the expected duration. 	Dv	Duration (min)
Sustained Activity	Note: Family only During the chosen period, PoSAs are identified when two or more activities are detected by the same activity sensor, prior to the configured expiration of a PoSA timer.	Gr or Dv If Gr, select Activity type	Expiration (min)
Unexpected Entry/Exit Timer	Note: For Family or for Pro with Panel CPU Version 2.4.5 and later An alert is triggered when a door on the premises is used, during the chosen period, on given weekdays. Detections made during the expected duration are ignored.	Gr or Dv If Gr, select Activity type	Duration (sec) and DoWs


Name	Description	Group (Gr)/ Device (Dv)	Value
Wandering	<p>Note: Family only</p> <p>An alert is triggered when a door on the premises is used, during the chosen period, on given weekdays and no detections are made for longer than the expected duration.</p>	N/A	Duration (min) and DoWs

5.4.6 PERS



NOTE: For ES6502HC, **PERS** parameters are available:

- Family or Pro with **Panel CPU Version 2.04** or later
- PERS with **Panel CPU Version 2.03** or later



The following table describes the PERS parameters that define PERS-related functionality at CP level.

Name	Definition and Instructions		Data Values	ES6502HC		C7000
				PERS Firmware	ADL Firmware	
Battery Management	Set Maximum Charging Level	Mark to allow setting the maximum battery charging level. Clear to allow the default maximum battery charging level.	Set level = Marked Use default level = Cleared	Panel CPU Ver 2.03.3 or later	N/A	√
	Charge Battery Up To	Select the maximum battery charging level.	Dropdown list			
Contact ID Acct No. with Group	Note: For PSTN models only Mark to enable the usage of the Contact ID group or partition as an extension of the account field in the CP configuration.		Enable = Marked Disable = Cleared	PSTN enabled & Panel CPU Ver 2.14.1 or later	N/A	✗

Name	Definition and Instructions	Data Values	ES6502HC		C7000
			PERS Firmware	ADL Firmware	
Carrier Redundancy	<p>Note: For CP models with GPRS SIA IP and ethernet communications channels only</p> <p>Mark to enable the CP to send messages to the same third-party receiver, using two different IP addresses.</p>	<p>Enable = Marked</p> <p>Disable = Cleared</p>	GSM enabled & Panel CPU Ver 2.03.3 or later	N/A	√
PERS Signaling Mode	<p>Note: For PSTN only</p> <p>Mark to force the CP to use only the Voice Signaling dialer. This option is used when only the Voice Signaling dialer can reach the monitoring station.</p> <p>Note: The CP sends only medical alarm (MA) messages.</p>	<p>Enable = Marked</p> <p>Disable = Cleared</p>	PSTN enabled & Panel CPU Ver later than 2.05 or Panel CPU Ver from 2.03.1.44 until 2.04	N/A	X
Send Device ID in Medical Alarm Restore Message	<p>Mark to add the ID of the device, that triggered the medical alarm, to the message, when one of the following occurs:</p> <ul style="list-style-type: none"> ● DTMF code 9 entered. ● Call disconnected without DTMF code. ● Panel RESET button pressed. 	<p>Enable = Marked</p> <p>Disable = Clear</p> <p>Default: Clear</p>	Disabled and unmarked		Enabled for CP versions: 6.01.9 through 6.02 and 6.03
Send Medical Alarms During Medical Alarm Window	<p>Within a single medical alarm window:</p> <p>Mark to send an alarm for each medical alarm triggered.</p> <p>Clear to send only one alarm, ignoring additional medical alarms.</p>	<p>Enable = Marked</p> <p>Disable = Clear</p> <p>Default: Marked</p>			
Debug For Essence Professional Services personnel only	<p>Disable</p> <p>Select this option to exit debug mode and resume regular operations.</p>	Dropdown list	Panel CPU Ver 2.03 or later	Panel CPU Ver 2.04 or later	√
	<p>Continuous Debug Logging:</p> <p>Select this option to enable the CP to log events continuously, sending the log files to the debug server whenever the buffer is full.</p>				
	<p>One-time Debug Logging:</p> <p>Select this option to enable the CP to log events only until the buffer is full. The CP sends the log file to the debug server and stops logging events.</p>				
	<p>Click to trigger the CP to:</p> <ul style="list-style-type: none"> ● Send the current log file to the debug server ● Stop logging events. 				


Name	Definition and Instructions		Data Values	ES6502HC		C7000
				PERS Firmware	ADL Firmware	
Emergency Call Retry	Enable	Mark to enable automatic redial when an emergency call is terminated prematurely.	Enable = Marked Disable = Cleared Default: Marked Default for panels not supporting the feature: Unmarked and disabled	N/A	N/A	CP version 6.04 and later.
	Number of Retries	Number of times to retry the emergency call Active when Enable is marked.	1-3			
	Minimum Call Duration	Minimum length of call, in seconds, for call to be considered completed. Active when Enable is marked.	5-30 seconds (5 second intervals)			
Control Panel Emergency Button	Enable	Mark to enable the CP to generate an emergency alarm when the Emergency Button is pressed.	Enable = Marked Disable = Cleared Default: Marked Default for panels not supporting the feature: Marked and disabled	N/A	N/A	CP version 6.05.2 and later.

5.4.7 DTMF Config

When you:

- Upload information from a CP or open a **.cpf** file, you see the DTMF keys as configured in the CP.

- Use a legacy **.cpf** file and there are no DTMF codes in the file, the downloaded DTMF codes will be the default codes from the **DTMFConfig.ini** file.



The purpose of the **DTMF Configurations** tab is to define the telephone keys to use for conversation flow, through the CP, between the service provider operators and the resident.

Refer to 2.3.4 Understanding the *DTMFConfig.ini* File on page 20.

The following table describes the configuration parameters on the **DTMF Configurations** tab.

Name	Names, Descriptions, and/or Instructions Note: Names are subject to panel software versions	ES6502HC		C7000
		PERS	Family or Pro / Family, Pro, & PERS from Panel CPU Ver 2.09.6	
DTMF Templates	Choose the template you need: ● USA ● EU ● UMO (Verklizan)	√	√	√
DTMF Keys	Switch to Full Duplex	√	√	√
	Half duplex: Operator speaks	√	√	√
	Half duplex: Operator listens	√	√	√
	Increase speaker volume / Full duplex: Increase speaker volume/ Full duplex: Increase speaker volume (GSM only)	√	√	√
	Decrease speaker volume / Full duplex: Decrease speaker volume/ Full duplex: Decrease speaker volume (GSM only)	√	√	√
	Disconnect the call	√	√	√

Name	Names, Descriptions, and/or Instructions Note: Names are subject to panel software versions	ES6502HC		C7000
		PERS	Family or Pro / Family, Pro, & PERS from Panel CPU Ver 2.09.6	
	Note: For Verklizan UMO only Alternate disconnect is "0"	N/A	N/A	C7000 CPU Ver 6.02.2 or later
	Note: For VPD voice extender feature Switch to CP	N/A	Panel CPU Ver 2.08 or later	√
	Note: For VPD voice extender feature Switch to next voice device	N/A	Panel CPU Ver 2.08 or later	√
	Perform reset to the panel after 5 presses	PSTN only	N/A	PSTN only

The following table lists the default values for each of the DTMF templates:

Action	USA	Europe	UMO	Comments
Talk	1	8	8	
Listen	3	7	7	
Increase speaker volume: 1 - 5 levels	2	3	3	
Decrease speaker volume: 1 - 5 levels	0	1	1	
Disconnect the call	9	0	*	
Alternate disconnect (UMO only)	No action	No action	0	Verklizan UMO only
Perform reset to the panel after 5 presses	8	2	2	PSTN CPs only
Next voice extender device	6	6	6	
Voice extender switch to CP	5	5	5	
Full duplex: Talk and listen	4	4	4	




NOTE: Click **Update** to add the revised DTMF parameter values to the CP configuration.

5.4.8 Log



NOTE: For ES6502HC, when GSM is enabled, **LOG** parameters are available only:

- For service types Family or Pro
- With **Panel CPU Version 2.05** or later




The **Log** tab allows you to send event log files automatically to a specified FTP Server for archiving purposes. The archived remote log files are assigned the extension **.rlg**.

The remote log files (**.rlg**) are retrieved from the FTP server upon request. Once retrieved onto the local PC, the **Log Events** module can locate and display the remote log files.

Refer to 13 The Log Events Module on page 159.

The following table describes the **Log** configuration parameters as they appear on the **Log** tab.

Parameter	Description	Values	
FTP	Domain Name	Mark if a URL is required rather than an IP address.	
	<FTP address>	IP: Enter a valid IP address.	Marked = URL Clear= IP (Default)
		URL: Enter a URL that references an Internet resource.	Four 3-digit strings with values from 0 - 255
	Port	Enter the port number used with the IP address entered in the previous field.	Character Integer (5) Values: 1-65536

Parameter	Description		Values								
	Path	Enter the path to the FTP Server.	Character								
	Username	Enter your username for your FTP login credentials.	Character								
	FTP Password	Enter your FTP password for your FTP login credentials.	Character								
		Click to view your encoded FTP password.	N/A								
Shipment	Enable Automatic Shipment	Mark to activate the automatic shipment of your event log file to the FTP server. The log file is automatically transmitted when the event log file reaches 90-100% capacity.	Marked = Enabled Clear= Disabled (Default)								
	Only over Ethernet	Mark to enable automatic shipment using only Ethernet communications. Note: If Ethernet communications is enabled but cellular communications is not enabled, this option default is enabled and display only.	Marked = Enabled Clear= Disabled (Default)								
	Enable Periodic Shipment	Mark to activate the Periodic transmission of your event log file to the FTP Server. The log file is automatically transmitted according to the periodic interval defined. Note: This field is activated only if the previous field is marked.	Marked = Enabled Clear= Disabled (Default)								
	Periodic Interval Note: These fields are displayed only if the previous field is marked.	<table border="1"> <tr> <td>Days</td> <td>Enter a time frame of days up to 90 days.</td> <td>Integer: 0-90</td> </tr> <tr> <td>Hours</td> <td>Enter a time frame in hours up to 24 hours.</td> <td>Integer:0-23</td> </tr> <tr> <td>Minutes</td> <td>Enter a time frame in minutes up to 60 minutes.</td> <td>Integer:0-59</td> </tr> </table>	Days	Enter a time frame of days up to 90 days.	Integer: 0-90	Hours	Enter a time frame in hours up to 24 hours.	Integer:0-23	Minutes	Enter a time frame in minutes up to 60 minutes.	Integer:0-59
Days	Enter a time frame of days up to 90 days.	Integer: 0-90									
Hours	Enter a time frame in hours up to 24 hours.	Integer:0-23									
Minutes	Enter a time frame in minutes up to 60 minutes.	Integer:0-59									


5.4.9 VPD Configuration Parameters



NOTE: For ES6502HC, **VPD** configuration parameters are available for:

- Family or Pro with **Panel CPU Version** 2.04 or later
- PERS with **Panel CPU Version** 2.03 or later

The purpose of the **VPD Configuration** tab is to configure VPD-related functionality at CP level. For example, the functionality defined by the following parameters applies to all the VPDs installed on the premises.



The following table defines the VPD configuration parameters at panel level:

Name	Definition and Instructions	Data Values	ES6502HC		C7000
			PERS	Family or Pro	
Detect Phrase Once	Do not recognize the trigger phrase again after the trigger phrase has already been recognized.	Enable = Marked Disable = Clear	√	√	√
Max Wait Between Phrases	Enter the maximum time (in seconds) required to wait from the first time the trigger phrase is recognized until the next time the trigger phrase is recognized.	Number in seconds Values: 10 – 600 Default: 10	√	√	√
Default Phrase Score	Score for at least one phrase For the VPD to trigger an emergency alarm, the emergency phrase must be: <ul style="list-style-type: none"> ● Detected twice, and ● One occurrence of the phrase must at least match this setting Refer to 5.4.9.1 VPD Quality Configuration on page 66	1 to 9	N/A	Panel CPU Ver 2.08 or later	√
	Score for both phrases For the VPD to trigger an emergency alarm, the emergency phrase must be: <ul style="list-style-type: none"> ● Detected twice, and ● Both occurrences of the phrase must at least match this setting Refer to 5.4.9.1 VPD Quality Configuration on page 66				
Phrase Detection Sound	Mark if you want the VPD to beep when an emergency alarm is triggered by voice detection.	Enable = Marked Disable = Clear	N/A	Panel CPU Ver 2.08 or later	√

Name	Definition and Instructions	Data Values	ES6502HC		C7000
			PER S	Family or Pro	
Voice Extender	Enable Voice Extender Feature Mark to enable the voice extender feature capability.	Enable = Marked Disable = Clear	N/A	Panel CPU Ver 2.08 or later	√
	Enable Switching Between Voice Devices Mark to enable switching between devices. Only available if the voice extender feature capability is enabled.	Enable = Marked Disable = Clear			
	Route to Alarm Triggering Device Mark to have communication from the monitoring station automatically routed to the VPD which triggers an emergency call. Only available if the voice extender feature capability is enabled.	Enable = Marked Disable = Clear			
	Enable LED Mark if you want the LED to light on the VPD when the operator is either listening or speaking	Enable = Marked Disable = Clear			

5.4.9.1 VPD Quality Configuration

The VPD can be configured to require the emergency phrases said by the resident to match a specified quality. The VPD triggers an emergency call when it detects two occurrences of the emergency phrase, and:

- Both occurrences of the phrase must at least match the **Score for both phrases** quality setting
- One occurrence of the phrase must at least match the **Score for at least one phrase** quality setting
- Detect Phrase Once is disabled

For example, assume the following:

- **Score for both phrases** is set to 3
- **Score for at least one phrase** is set to 7

The following table summarizes possible outcomes given different qualities of phrases said by the resident:

Resident Phrases	Both Phrases: 3	One Phrase: 7	Outcome
<ul style="list-style-type: none"> ● Quality 2 ● Quality 2 	✗	✗	No emergency call
<ul style="list-style-type: none"> ● Quality 3 ● Quality 4 	✓	✗	No emergency call
<ul style="list-style-type: none"> ● Quality 3 ● Quality 7 	✓	✓	Emergency call
<ul style="list-style-type: none"> ● Quality 7 ● Quality 7 	✓	✓	Emergency call




NOTE: If the current VPD firmware version does not support phrase scores (the **Score for both phrases** and **Score for at least one phrase** settings), the phrase scores configuration goes into effect automatically when the VPD firmware version is upgraded.

5.4.10 EPA Configuration Parameters



NOTE: For ES6502HC, **EPA** configuration parameters are available:

- Family or Pro with **Panel CPU Version** 2.05.6 or later
- PERS with **Panel CPU Version** 2.03.4 or later



The purpose of the **EPA Configuration** tab is to configure EPA-related functionality at CP level. For example, the functionality defined by the following parameters applies to all the EPAs installed on the premises.

The following table defines the EPA configuration parameters at panel level:

Name	Definition and Instructions	Data Values	ES6502HC		C7000
			PERS	Family or Pro	
Usage Note: Usage data includes: <ul style="list-style-type: none"> The step count: reported for each 8-hour interval within the last 24 hours. Confirmation that the EPA is being worn: reported hourly. 	Enable Mark to allow the EPA to send usage data to the CP at the frequency specified. Default: Clear	Enable = Marked Disable = Clear	√	√	√
	Usage Alert Frequency Select the frequency at which the usage data is sent to the CP. Default: 8 hours	Dropdown list values(hours): 8, 16, 24	√	√	√
BLE	Enable Mark to enable Bluetooth Low Energy (BLE) communications with mobile devices. Default: Clear	Enable = Marked Disable = Clear	N/A	Panel CPU Ver 2.07 or later	Not PERS
	Enable Periodic Test Mark to enable periodic connectivity test for an EPA connected to a mobile device via BLE. Default: Clear Note: Parameter disabled when BLE-Enable is not marked.	Enable = Marked Disable = Clear	N/A	Panel CPU Ver 2.07 or later	Not PERS
Fall Detection	Wake Up Sensitivity Select the sensitivity level for triggering a detection The higher the level – the more sensitive Note: It is highly recommended to use the default level of 5	Dropdown list: Default, 1-7	√	√	√
	Enable for New Device Mark to enable accelerometer features for EPA devices recently added to Care@Home™ . Default: Marked	Enable = Marked Disable = Clear	N/A	Panel CPU Ver from 2.06.6 until 2.07 or Panel CPU Ver 2.09 or later	√
	Enable Free Fall Mark to enable the fall test of an EPA's fall detection. Fall tests generate a medical alarm (MA). Default: Clear	Enable = Marked Disable = Clear	N/A	Panel CPU Ver 2.09.5 and later	C7000 CPU Ver 6.01.5 and later

5.4.11 RF Settings



NOTE: ES6502HC supports **RF** settings only for service types Family or Pro with **Panel CPU Version 2.05.4** or later

This **RF Settings** parameter is for Essence internal use only.

5.5. Periodic Test

The following parameters combined define the reporting frequency for forwarding a periodic system test report to the monitoring station.

Name	Definition and Instructions	Data Values	Factory Settings
Enable Night Period	NOTE: Only for an ES6502HC CP with PSTN interface and PERS service type with SW version >= 2.03.2, or C7000 with PSTN. Mark to allow scheduling a period test at night.	Enable = Marked Disable = Clear	Disable = Clear
Days	The part of the reporting frequency that defines number of days	Number (2) Values: 1-31 Unit of measure: days	7 days
Hours	The part of the reporting frequency that defines number of hours	Number (2) Values: 0-23 Unit of measure: hours	0 hours
Minutes	The part of the reporting frequency that defines number of minutes	Number (2) Values: 0-59 Unit of measure: minutes	0 minutes

5.6. Mains Failure

The following table describes the CP configuration parameters for broadcasting power-outage events.

Name	Definition and Instructions	Data Values	Required/Optional	Factory Settings
Message Enable	Mark to allow the broadcast of power-outage events to the monitoring station.	Enable = Marked Disable = Clear	Optional	Enable = Marked
Message Delay	A fixed time-delay before broadcasting a power-outage event.	N/A	Display Only	4 minutes

Name	Definition and Instructions	Data Values	Required/Optional	Factory Settings
Random Range Limit	A 2-minute time-interval during which the power-outage event is broadcast randomly. This is used to avoid overloading the monitoring station.	N/A	Display Only	2 minutes

5.7. Temperature



NOTE: For ES6502HC, the **Measurement Interval** parameter is available with **Panel CPU Version 2.05** or later.

The following table describes the CP configuration parameters for extreme temperature.

Name	Definition and Instructions	Data Values	Factory Settings
Extreme Temperature Enable	Mark to activate the capability to trigger the safety alarm when a device reports a temperature that passes the given thresholds.	Enable = Marked Disable = Clear	Disable = Clear
Freeze Threshold	Set the required low threshold for cold temperatures.	Dropdown list Values: 30-75 at 5° increments. Unit of measure: Fahrenheit	40°F
Heat Threshold	Set the required high threshold for hot temperatures.	Dropdown list Values: 75-120 at 5° increments Unit of measure: Fahrenheit	100°F
Measurement Interval	<p>Note: For:</p> <ul style="list-style-type: none"> ● C7000 CP ● ES6502HC CP with Family or Pro service types with CP SW version >= 2.05 ● Only for MDS, SPB, and VPD only Set the room temperature measurement interval for the device.	Values: 0-255 minutes Unit of measure: minutes	

5.8. Emergency Pendant

The following table describes the CP configuration parameters for using any of the following emergency pendant devices: EP, EPP and EPA.

Name	Definition and Instructions	Data Values	Factory Settings
Min Press Duration	Set the minimum press time to trigger an emergency event using an emergency pendant	Dropdown List Values: 0.1-12 with increments of 0.5 Unit of measure: seconds	0.3 seconds
Pairing Min Press	Set the minimum time to press the button to pair the emergency pendant with the CP. Note: To reset a pendant that was pre-registered with the panel, press for a duration between the Pairing Min Press and the Pairing Max Press.	Dropdown List Values: 4-15 Unit of measure: seconds	4 seconds
Pairing Max Press	Set the maximum time to press the button to pair the emergency pendant with the CP. Note: To reset a pendant that was pre-registered with the panel, press for a duration between the Pairing Min Press and the Pairing Max Press.	Dropdown List Values: 4-23 with increments of 0.5 Unit of measure: seconds	23 seconds
Supervision Report See note below	The maximum number of hours within which the CP must receive at least one supervisory report from the pendant. Note: If no report is received, the CP sends an event message to the monitoring station that the emergency pendant is lost.	A dropdown list appears for some configurations, otherwise the value is set to 335 hours. Dropdown values: 5, 40, 50, 60, 70, 80, 90, 100, 200, 335, 400	335 hours



NOTE: For ES6502HC, the **Supervision Report** parameter is hardcoded as 335 hours with **Panel CPU Version 2.03** or later.

5.9. Medical Alarm

The following table describes the CP configuration **Medical Alarm** parameters for detection of:


- Extreme temperature
- Medical alarm
- Fire, smoke or gas
- Entrance door left open

Name	Definition and Instructions	Data Values	Factory Settings
Delay	Set the duration before the medical alarm is triggered, to allow you to cancel the medical alarm, if triggered by mistake. Cancel the medical alarm by pressing the CP RESET button.	Number (3) Values: 0-180 Unit of measure: seconds	10 seconds
Voice Window	Set the time frame for the CP to automatically answer an incoming call. The window allows the operator to call the resident after a call was terminated accidentally, (not by DTMF) after the medical alarm was triggered.	Number (2) Values: 0-60, 120, 180, 240. "Disable" appears when value is 0. "Unlimited" option available for some configurations. Unit of measure: minutes	10 minutes

5.10. Password

Enter the required 4-digit DTMF code used to allow the CP to receive remote calls. Default: 1234. Notify Essence technical support if you use a code other than the default.

5.11. Speakerphone

Mark to activate the option to answer incoming calls using the  which triggers a medical alarm. This parameter activates the speaker in the CP. The default value is **Enabled**.

5.12. Remote Call-In

The following table describes the CP configuration parameters for incoming calls on the CP.

Name	Definition and Instructions	Data Values	ES6502HC			C7000
			PERS	Pro	Family	
Enable	Mark to allow the CP to automatically answer incoming calls.	Enable = Marked Disable = Cleared Default: Cleared	√	√	√	√
Rings	Set the number of rings for the CP to wait before automatically answering a remote call. The user must enter the DTMF code to start a half-duplex conversation.	Number (2) Values: 4-12	√	√	√	√
Mute Rings	Mark to silence the CP ringer for remote calls.	Enable = Marked Disable = Cleared Default: Cleared	Panel CPU Ver is <ul style="list-style-type: none"> ● Later than or equal to 1.03.45 but earlier than 2 ● Later than 2.03 	Panel CPU Ver is 2.04 or later	Panel CPU Ver is 2.04 or later	√

5.13. Sound

The following table describes the CP configuration parameters for the buttons, sirens, and voice announcements.

Name	Definition and Instructions	Data Values	ES6502HC		C7000	
			PERS	Not PERS	PERS	Not PERS
Announcement Volume	Set the volume of the button tones and the vocal announcements. Default: 7	Dropdown List Values: 1-7 where 1 is Low.	✓	✓	✓	✓
Medical Siren Duration	Set the duration for sounding the medical alarm. Default: 5 seconds	Values: 0-180 or 0-30 seconds	0-30 seconds for Panel CPU Ver <ul style="list-style-type: none"> ● Later than or equal to 1.03.45 but earlier than 2 ● Equal to 2.03 or later Else 180 seconds	0-180 seconds	0-30 seconds	0-180 seconds
Safety Siren Duration	Set the duration for sounding the safety siren.	Values: Default =180 seconds, 1-250 seconds, or Disable	N/A	Supported from Panel CPU Ver 2.07 or later	✓	✓
Announcement While Connecting	Mark to enable repeating a "Call in progress" announcement, to assure residents that medical alarms are being handled.	Enable = Marked Disable = Cleared Default: Cleared	Disabled and clear		Enabled for CP versions: 6.01.9 through 6.02 and 6.03	

5.14. Emergency

The following table describes the CP configuration parameters for receiving an emergency call.

Name	Definition and Instructions	Data Values	ES6502HC			C7000
			PERS	Pro	Family	
Conversation Mode	Select the starting mode of the conversation flow between the resident and the monitoring station during an emergency call.	Full Duplex	✓	✓	✓	✓
		Half Duplex	Panel CPU Ver earlier than 2.03.5	✓	✓	
		Half Duplex – Talk	Panel CPU Ver 2.03.5 or later			✓
		Half Duplex - listen	Panel CPU Ver 2.03.5 or later			✓
Speaker Volume	Set the CP speaker volume level for the beginning of the conversation, for a half-duplex emergency call.	Dropdown list Values: 1-7	Panel CPU Ver is: <ul style="list-style-type: none"> ● Later or equal to 1.03.45 but earlier than 2 ● Later than 2.03 	Panel CPU Ver 2.04 or later	Panel CPU Ver 2.04 or later	✓

5.15. Inactivity



NOTE: For ES6502HC, the **Inactivity** feature is available and enabled only if your service type is PERS.

The **Inactivity** feature uses a global parameter that allows you to redefine the functionality of the devices in the **Activity Sensors** module, such as the motion detector (PIR) and the camera detector (IPD), to detect inactivity.

When the feature is enabled, the devices act as inactivity sensors. If no activity is detected during the defined **Time Window**, the CP sends an inactivity event to the monitoring station.

The following table describes the CP configuration **Inactivity** parameters.

Name	Definition and Instructions	Data Values	Required /Optional	Default
Enable	Mark to activate inactivity monitoring.	Enable = Marked Disable = Cleared	Optional	Disable = Cleared
Time Window	Define the duration to monitor inactivity.	Dropdown list Values 8-48 Unit of measure: Hours	Required	8 hours

5.16. Activity Timer



NOTE: For ES6502HC, the **Activity** feature is available and enabled only if your service type is PERS.

The purpose of the activity timer parameters is to monitor whether the resident is active.

The resident is expected to press the CP **RESET** button during defined time intervals. The CP reminds the resident to press **RESET**, 15 minutes prior to the end of the time interval. The CP issues an alarm if the resident does not press **RESET** during the defined time interval.



NOTE: If activity is not verified, an inactivity event is sent to the monitoring station.

The following table describes the CP configuration **Activity Timer** parameters.

Name	Definition and Instructions	Data Values	Required /Optional
Enable	Mark to activate the ability to monitor whether the resident is active within the premises. NOTE: If marked, the CP emergency call button LED is lit Blue : <ul style="list-style-type: none"> ● For ES6502HC, for the day if both Timer 1 or Timer 2 are clear. ● For C7000, only during the periods defined by enabled Timer 1 or Timer 2 parameters 	Enable = Marked Disable = Cleared Default: Cleared	Optional
Timer 1	Mark to activate the first time-interval when resident activity is monitored.	Enable = Marked Disable = Cleared Default: Cleared	Optional
	< Start / End Interval > Define the first time-interval, defined by a start time and an end time, when resident activity is monitored. NOTE: Overlapping time periods are not acceptable.	Time-Picker field Format: HH:MM AM/PM	Required if Enabled
Timer 2	Mark to activate the second time-interval when user activity is monitored.	Enable = Marked Disable = Cleared Default: Cleared	Optional
	< Start / End Interval > Define the second time-interval, defined by a start time and an end time, when resident activity is monitored. NOTE: Do not enter overlapping time periods.	Time-Picker field Format: HH:MM AM/PM	Required if Enabled

5.17. Smoke Detector



NOTE: For ES6502HC, the **Smoke Detector** parameter is supported.

For service type Family and Pro with **Panel CPU Version** 2.09.9 or later, the parameter is enabled. By default, it is cleared.

Mark to activate the option to initiate a call to the monitoring station when smoke is detected. This parameter activates the speaker in the CP. The default value is **Enabled**.

5.18. The Control Panel Module Buttons


The **Control Panel** module buttons are located at the bottom of the screen. The functionality of the module buttons is:

- **Clear** – Restores the parameter values to the values of the last update or the original values from the CP upload.


- **Update** - Updates all the new and revised configuration parameter values to the **Control Panel** module.

You must click **Update** to ensure that the modifications are integrated into the CP configuration, before downloading or saving.

6. The Communications Module



Click **Communications** on the **Main** navigator menu. The **Communications** module allows for configuration management of multiple communication channels, enabling remote activation or redundancy, in case of channel failure or other connectivity issue. The multiple channels support two-way data communications between the Care@Home™ CP, installed on the resident’s premises, and the monitoring station.



The following are the available CP communication channels:

- Public Switched Telephone Network (PSTN)
- Cellular
- Ethernet

The **.cpf** file received from Essence for each CP is pre-configured with the matching communications channel parameters.

6.1. General Communications



NOTE: To save changes, remember to click **Update** before exiting the module or the application.

General Communications defines the communications-related parameters for the C7000 configured for multiple communications interfaces and ES6502HC panels with cellular capability.

6.1.1 Link Lost or Restored

When the CP is configured with multiple communication interfaces, if there is no link for a specific communication interface, this configuration defines whether warning events should be sent over other available link(s).

Name	Definition and Instructions	Data Values
Event	Whether a link event is meant to be sent, if there is no link on the communication interface: If sent, whether for an active or enabled link.	<ul style="list-style-type: none"> ● Sent for each active link ● Sent for each enabled link ● Never sent
Event Message Delay	Length of message delay, if the delay frame is enabled. Note: If the warning event is never sent, the delay frame is disabled.	0-60 minutes
Event Random Range Limit	Length of random delay before sending link events, if there are regional communications issues.	0-60 minutes
Voice & LED	Whether CP should warn user with vocal announcement and/or LEDs, if there is no link on the communication interface.	<ul style="list-style-type: none"> ● Warn when both links are down ● Warn when enabled link is down ● Warn when active link is down

6.1.2 IP ACK Timeout Configuration

The following table describes the acknowledgement timeout parameter.

Name	Definition and Instructions	Data Values
ACK Timeout	Set the wait time to receive an acknowledgement message.	Number (3) Values: 15-300 seconds

6.1.3 Communication

The following table describes the communication parameter.

Name	Definition and Instructions	Data Values	Required/Optional
Enable PSTN-Cellular Redundancy	Note: Not configurable	Disabled= Clear	N/A

6.1.4 DNIS

The following describes the general communication parameters for the DNIS.

Name	Definition and Instructions	Data Values	Factory Settings
DNIS (Dialed Number Identification Service)	<p>Note: Only when cellular or Ethernet communications is enabled</p> <p>Enter the DNIS number, included in the SIA IP and SIA DC-09 messages sent from the CP to the monitoring station.</p> <p>The number is used by the monitoring station to differentiate between multiple service providers.</p>	Hexadecimal	00000
Combine DNIS with Account No.	<p>Note: Relevant for:</p> <ul style="list-style-type: none"> ● C7000 (cellular / Ethernet onboard) ● ES6502HC (cellular enabled): <ul style="list-style-type: none"> ■ Panel CPU version 2.03.7 through 2.04 ■ C Panel CPU version 2.07 and later <p>Mark to combine the last five digits of the account number with the DNIS, SIA IP communication channel only. Can also be used with DC-09, using the entire account number</p> <p>Format: <DNIS> <Last 5 digits (Account #)></p> <p>Note: Combining the DNIS with the account number increases the number of available account numbers.</p>	Combine = Marked DNIS only = Cleared	Clear

6.1.5 Phone Number

The following describes the phone number of the SIM or landline.

Name	Definition and Instructions	Data Values	ES6502HC			C7000
			PERS	Pro	Famil y	
Phone Number	<p>Enter the cellular number of the SIM card installed in the CP. For the C7000 PSTN/Ethernet model, enter the landline number.</p> <p>Note: If the CP configuration, uploaded to the CMS warehouse, includes a SCAIP dialer with callback enabled, the callback SIM number entered during CMS warehouse processing overrides this SIM number.</p>	International cellular phone number or landline for C7000 PSTN/Ethernet model	Panel CPU Ver 2.03.2 or later	Panel CPU Ver 2.05.4 or later	Panel CPU Ver 2.05.4 or later	√

6.1.6 Authentication Password




NOTE: Applies to

- ES6502HC:
 - Cellular enabled
 - Panel CPU version 2.05.6.15 through 2.14
- C7000

The following table describes the authentication password parameter.

Name	Definition and Instructions	Data Values	Required/Optional
Password	SCAIP password	Enable = Marked Disable = Clear	Optional

6.2. Cellular Communications



Cellular Communications contains cellular configurations, defining the communications-related parameters for both:

- C7000 CP
- ES6502HC CP

The communications-related parameters are relevant to the ES6502HC CP and to C7000 CPs configured for a cellular communication interface.



NOTE: All cellular communications parameters are only required when cellular communication is enabled.

6.2.1 Module

The following table describes the module parameters.

Name	Definition and Instructions	Data Values
Module	Module identification for cellular communications. Display only.	Includes CMS module code plus the module version: major, minor, subversion and build.
IMEI Number	International Mobile Equipment Identity (IMEI) number. Display only.	IMEI number

The module for the ES6502HC panel is always 3G/2G

The following table lists the modules by GSM band for the C7000 panel:

CP Cellular Modules	Band	Band Data Allowed <i>Note:</i> Only if cellular communication is enabled
C1NA	3G	4G, 3G (N.America LTE)
C1EUR	3G	3G, 2G (Europe 3G)
C1AP	4G	4G, 3G (APAC LTE)
C1EU	4G	4G, 3G, 2G (Europe LTE)
C1LA	4G	4G, 3G, 2G (Latin America LTE)

6.2.2 Provider

The following table describes the GSM configuration parameters for the cellular provider.

Name	Definition and Instructions	Data Values	Required/Optional	Factory Settings
Lock	Mark to allow the selection of a specific cellular provider. Note: This parameter disables SIM card roaming.	Enable roaming = Marked Disable roaming = Clear	Optional	Disable = Clear
Name	Choose the cellular provider for the CP SIM. Note: The available options can be configured using <i>Mobile Operators.ini</i>	Dropdown list per resident's region or N/A	Required	N/A
Switch Automatically	Mark to switch automatically between providers if Lock is not marked. C7000 version 6.03.2 and later only.	Enable automatic switching = Marked Disable automatic switching = Clear	Optional	Disable = Clear
Blacklisted Provider Code	Indicates blocked cellular provider.	Cellular provider code	Read-only	N/A

6.2.3 APN

The following table describes the cellular configuration parameters for APN credentials.

Name	Definition and Instructions	Data Values
APN	The Access Point Name (APN) for your SIM card service provider.	The CP loads the actual APN information from the Mobile Operators module, according to the APN name from the SIM-card inserted in the CP.
User Name	The APN user name assigned by your service provider	
Password	The APN password assigned by your service provider.	



NOTE: If the APN information does not allow the CP to connect, the CP should request other APN information from the network. When the CP finally connects successfully, the CP saves the APN information retrieved from the mobile operator module in these parameters.

6.2.4 Audio Codec



NOTE: The Audio Codec parameter is available parameter is available for:

- C7000 panels
- ES6502HC ADL firmware version between
 - 2.06.6 and 2.07 (excluding 2.07).
 - 2.09 and 2.14 (excluding 2.14).

The following table describes the GSM configuration parameter for the GSM/UMTS audio codecs settings.

Name	Definition and Instructions	Data Values	Factory Settings
Codecs	Configure the GSM/UMTS audio codecs required setting.	Values: 0-255	On initialization, Codecs is set automatically: If the audio language is either Slovenian or Hebrew, Codecs is set to 95. Otherwise, Codecs is set to 0.



NOTE: Consult with Essence technical support personnel before changing the GSM/UMTS audio codecs setting.

6.2.5 SIM

The following describes the GSM configuration parameters for the SIM.

Name	Definition and Instructions	Data Values
ICCID	Unique SIM card identifier. Display only.	ICCID number
PIN Code Enable	Mark PIN Code check box to require a PIN code entry.	Enabled=Marked Disabled=Clear
PIN Code	PIN for the CP SIM card.	Number (4)

6.2.6 Dialer Configurations

The following describes the GSM configuration parameters for the dialer configurations.

Name	Definition and Instructions	Data Values	Factory Settings
Min Reception Level	Set the minimum GSM reception level.	Number (2) Values: 1-31	10
Number of Retries	Enter the number of times the dialer redials a phone number.	Number (3) Values: 0-255	2
Retry Timeout	Enter the time interval, in seconds, between redials.	Number (3) Values: 0-255 Unit of measure: seconds	1


6.2.7 Call Configurations

The following describes the GSM configuration parameters for the call configurations.

Name	Definition and Instructions	Data Values	Factory Settings
Ring Tones	Select the ring tones for when the CP receives a call.	Dropdown list Values, "Default", 1-18	"Default"
DTMF blocking period	NOTE: For ES6502HC, same for PSTN and GSM calls. Select the period duration (in seconds) that the control panel ignores DTMF from the beginning of the emergency audio call	Default or 1-60 seconds	"Default"

Refer to 2.3.4 Understanding the *DTMFConfig.ini* on page 20.

6.3. PSTN Communications




PSTN Communications contains PSTN configurations, defining the communications-related parameters for both:

- C7000 CP
- ES6502HC CP

The communications-related parameters are relevant to the ES6502HC CP and to C7000 CPs that are configured for a PSTN communication interface.

Name	Definition and Instructions		Data Values	ES6502HC			C7000
				PERS	Pro	Family	
Enable	<p>Note: For CP with PSTN interface Mark to activate control panel landline communications</p>		Enable = Marked Disable = Clear	√			√
Configurations	Tone/Pulse	Select the dialing method of the resident's telephone	Dropdown List Tone or Pulse	When PSTN communication is enabled			√
	Ack Timeout	Set the wait time to receive an acknowledgement message	Number (3) Values: 1-180 in seconds	When PSTN communication is enabled			√
	DTMF Configuration	Select the external dual-tone multi-frequency (DTMF) configuration file name. Note: Essence recommends using the Profile-1 settings. To define your own settings, use the Profile-1 settings as a guiding template.	Dropdown list of available DTMF configuration files: <ul style="list-style-type: none"> ● Default from Panel ● Profile-1 ● <Additional DTMF profiles as defined in the <i>PstnDtmf.ini</i> file> 	When PSTN communication is enabled and SW version >= 1.03			√
	DTMF Blocking Period	This feature is used during an emergency call, for a given time-period, to intercept the DTMF codes sent to the CP. The purpose is to prevent the CP from acting according to these DTMF codes. Note: The DTMF codes are sent when the monitoring station is transferring the call to the operator assigned to answer the emergency call.	Dropdown List Values: 1-60 seconds and "Default" Note: Default – The blocking period is disabled. All DTMF codes are received by the CP.	When PSTN communication is enabled: same configuration is used for PSTN and GSM, since only one interface can be enabled at the same time			Separate configurations for cellular and PSTN interfaces

6.4. Ethernet Communications



Ethernet Communications contains Ethernet configurations, defining the communications-related parameters only for the C7000 CP configured for an Ethernet communication interface.




NOTE: Ethernet configurations and related parameters are only relevant to C7000 CPs.

The following parameters are mandatory if the C7000 CP is configured for an Ethernet interface.

Name	Definition and Instructions		Data Values	Comment
Ethernet	Enable	Note: For CP with Ethernet interface Mark to activate control panel full Ethernet communications	Enable = Marked ● Disable = Clear	●
	MAC Address	The Ethernet mac address assigned to the C7000 CP.	Hexadecimal Format: FF-FF-FF-FF-FF-FF	Display only
IP Setup Configurations	Automatic IP	Mark to activate automatic IP setup configuration by DHCP.	Enable = Marked Disable = Clear	Default: clear
	IP Address	Enter the CP's IP address on the local Ethernet network	4 3-digit strings	Mandatory

Name	Definition and Instructions		Data Values	Comment
	Subnet Mask	Enter the subnet mask for your router to work, according to your router network settings	Values from 0- 255	
	Default Gateway	Enter your router's IP address		
	Primary DNS	Enter the IP address for the domain name system server.		
	Secondary DNS	Enter the IP address for the backup domain name system server.	4 3-digit strings ● Values from 0-255	● Optional

7. The Dialer Module



Click **Dialer** on the **Main** navigator menu. The **Dialer** module is the tool that enables the definition of communication methods and the management of the flow required to communicate with a monitoring station or Care@Home™ servers. The module allows you to manage the destination list, comprised of domain names, IP addresses, and in some cases, phone numbers, along with the connect sequences for communication between the Care@Home™ CP and the monitoring station or Care@Home™ servers.



NOTE: If Care@Home™ is configured to work **without** the PERS receiver, do not use the PERS Dialer functionality such as:

- The Comm Type: Message
- The **Sequence Protocols:** PERS Supervisory and PERS Alert

7.1. Destination Definitions by CP Type

A destination is defined as one of the following:

- Phone number
- IP address and a port number
- Domain name (URL) and a port number

Destination attributes include, among others, the communications type, communication channels (interfaces), and protocols, defined for each destination, depending on the CP model type.

The CP model types are:

- **ES6502HC** – This legacy CP model allows for either PSTN or GSM communications
- **C7000 Pre Version 6.1** – This upgraded CP model allows you to select one of the following specific communication channels, alongside the traditional communication channels like PSTN.
 - LTE (4G)
 - 3G

- 2G
- **C7000 Version 6.1 and later** – From this version, the CP is upgraded to support the following multiple communication channel types, known as interfaces:
 - PSTN
 - Cellular
 - Ethernet

This version allows for multiple communication interfaces to be enabled simultaneously.

This upgraded CP model supports the Ethernet interface in addition to the LTE (4G), 2G and 3G interfaces, alongside the traditional communication interfaces.




NOTE: The **Message** communication type in this control panel is redefined as multiple interfaces with multiple protocols for each interface.


7.2. Destination Definitions for an ES6502HC CP

This legacy control panel model supports either PSTN or GSM communications.

The following shows ES6502HC dialers:



7.2.1 Defining Communications Types and Channels



The following table describes the basic attributes required to add a destination to the Destination list.

Name	Definition and Instructions	Data Values	ES6502HC	
			PERS	Pro & Family
Destination ID	Internal identification of a dialer	Values: 1-16	✓	✓
Phone Num/ IP :Port	The phone number or IP/URL plus port number entered for the destination ID receiver address	Receiver address	✓	✓
Comm Type	Select the communications type for the Destination. NOTE: Use the "Message" Comm Type only if the Care@Home™ system is configured to work with a PERS receiver.	Analytics Message	N/A	✓
		Message	✓	✓
		Data Stream	N/A	✓
		Video	N/A	✓
		Voice	N/A	GSM enabled & Panel CPU Ver 2.05.2 or later
		STUN	N/A	✓
Comm Channels	Select the broadcast method to use with the communications type entered in the previous field.	Debug Printing	Panel CPU Ver 2.03 or later	Panel CPU Ver 2.04 or later
		Dropdown list changes according to the selected Comm Type.	✓	✓

The following table lists the communication channels available for each communication type.

Comm Type	Description	Comm Channels		Emergency Number
		ES6502HC		
		PERS	Family or Pro	
Analytics Message	For sending event messages to the Care@Home™ Server	N/A	<ul style="list-style-type: none"> ● GPRS ● GPRS 2G with Panel CPU Ver 2.04 or later 	Emergency number allowed only for first 10 rows with IP address/domain name Note: Number erased only using delete key

Comm Type	Description	Comm Channels		Emergency Number
		ES6502HC		
		PERS	Family or Pro	
Message	For sending event messages to the PERS receiver. Note: Use this Comm Type only if Care@Home™ is configured to work with a PERS receiver.	<ul style="list-style-type: none"> ● GPRS with Panel CPU Ver earlier than 2.05.4 ● GPRS SIA IP GSM-enabled & Panel CPU Ver 2.03.3 or later ● GPRS SCAIP GSM-enabled & Panel CPU Ver 2.03.3 or later ● SIA DC-09 GSM-enabled & Panel CPU Ver 2.03.3 or later ● SIA-PSTN-enabled ● Contact ID - PSTN-enabled & Panel CPU Ver 2.03 or later ● Voice Signaling – PSTN-enabled & Panel CPU Ver 2.03 or later ● CPC – PSTN-enabled & Panel CPU Ver 2.03.2 or later ● BS 8521 – PSTN-enabled and panel CPU version 2.16 or later ● 4x2 – Panel version 2.14 and later ● T2i SMS - Panel CPU version 2.11 and later (Essence technical support only) 	<ul style="list-style-type: none"> ● GPRS with Panel CPU Ver earlier than 2.05.4 ● GPRS SIA IP enabled & Panel CPU Ver 2.05.4 or later ● GPRS SCAIP - enabled & Panel CPU Ver 2.05.4 or later ● GPRS SIA DC-09 enabled & Panel CPU Ver 2.05.4 or later ● GPRS BS 8521 - Ver 2.09.6 or later ● T2i SMS - Panel CPU version 2.11 and later (Essence technical support only) 	Emergency number allowed only for first 10 rows Note: Number erased only using delete key
Data Stream	For sending files to the Care@Home™ server that contain: <ul style="list-style-type: none"> ● Detections from Family and Pro activity sensors ● Panel and peripheral information 	N/A	<ul style="list-style-type: none"> ● GPRS ● 2G with Panel CPU Ver 2.04 or later 	
Video	For sending photos taken by the camera detectors (IPD)	N/A	<ul style="list-style-type: none"> ● GPRS ● 2G with Panel CPU Ver 2.04 or later 	
Voice	For calling the monitoring station Fail-Safe number	Care@Home 2.5.14 or later and panel software version 2.11 or later	GSM	Phone number required Note: Number erased only using delete key

Comm Type	Description	Comm Channels		Emergency Number
		ES6502HC		
		PERS	Family or Pro	
STUN	For sending messages, via UDP, to the Care@Home™ server that enables the CP to receive commands via UDP	N/A	● GPRS	
Debug Printing	For Essence technical support only	<ul style="list-style-type: none"> ● Panel CPU Ver 2.03 or later ● GPRS – GSM enabled ● PSTN - PSTN enabled 	<ul style="list-style-type: none"> ● Panel CPU Ver 2.04 or later ● GPRS 	

7.2.2 Selecting Communication Channels

The available **communication channels** are the same for either communication type method:

- **PSTN** –Phone number.
- **GPRS** – Options:
 - IP address
 - Domain address (URL) with a port number, such as 3061
 - Phone number for GPRS BS8521 and T2i SMS

To select your **communication channels**, choose the broadcast method that matches your CP communications method.

7.2.2.1 Configuring PSTN Communications Method

For PSTN communications method, the available broadcast methods are:

- SIA DC-03
- Contact ID
- Voice Signaling
- CPC default/NEAT
- GPRD BS 8521
- 4x2

When choosing any of these methods, the **Dialer** module prompts you to enter the landline phone number of the monitoring station receiver.

When choosing **Voice Signaling**, the **Dialer** module prompts you to enter the **direct** phone number of the monitoring station. The Voice Signaling mode is used to ensure that event messages reach the monitoring station, even if the receiver cannot interpret messages sent by other protocols. For


ES6502HC panels with SW version \geq 2.16 (PERS 2.4) – the phone number in the Voice Signaling destination is for emergency alarms – but additional dialing scenarios can also be entered.

7.2.2.2 Configuring GSM Communications Method

For GSM communications method, the following are the available cellular communication protocol methods, between the CP and the monitoring station:


- GPRS SIA IP
- GPRS SCAIP
- GPRS SIA DC-09
- GPRS BS 8521
- T2i SMS

When choosing the **GPRS** communication channel, the **Dialer** module prompts you to enter receiver address information, except when using GPRS BS 8521 or T2i SMS.




7.2.3 Defining Receiver Address Attributes

Most of the communication channels are based on receiver address attributes.



For the **Message** communication type, the communication channels require the same attributes except for **SCAIP**.



The following table lists and defines the receiver address attributes.

Name	Definition and Instructions	Data Values	Required/Optional
Domain Name	Mark if a URL is required rather than an IP address.	Enable = Marked Disable = Clear	Optional
URL	When Domain Name is marked, enter the web address that refers to the IP address of the monitoring station receiver. NOTE: Domain name can be up to 32 ASCII characters only	Character “<domain name>” For example: “www.essence-grp.com”	Required if Domain Name marked
IP address	When Domain Name is not marked, enter the IPv4 receiver address.	A set of 4 3-digit strings Values from 0- 255	Required
Port*	Enter the port number used by the monitoring station receiver. The default ports are listed in the table below.	Integer (5) Values: 1-65535	Required

Name	Definition and Instructions	Data Values	Required/Optional
Protocol**	<p>NOTE: Not applicable for Debug Printing</p> <p>Select the IP communications protocols:</p> <ul style="list-style-type: none"> ● TCP ● UDP 	Dropdown list depends on communication channel	Required
Emergency Number	<p>NOTE: Comm Type Analytics Message and Message with GSM enabled. For C7000, GSM enabled is not required</p> <p>Enter a phone number to call in case of an emergency, after sending the emergency event to the receiver and/or the Care@Home™ backend.</p>	<p>19-character string. Accepts numeric digits, "+" and "#" only.</p> <p>NOTE: Available only for the first 10 rows</p> <p>Note: Number erased only using delete key</p>	Optional

* - The following table lists the receiver address ports defined as defaults for the communication channels:

Interfaces	Default Port
SIA-IP	3061
SCAIP	5060 – default SIP port
SIA DC-09	Clear – no default
Analytics Alert	10060
STUN	3478
Video	5666
Data stream	
Debug printing	

** - The following receiver address protocols are allowed for the communication type and/or channels:

Communication Type	Communication Channel	Receiver Address Protocols
		ES6502HC
Analytics Message Alert		TCP/UDP
Message	SIA IP	UDP
	SCAIP	UDP
	SIA DC-09	UDP/TCP
Data stream		TCP/UDP
Video		TCP/UDP
Voice		
STUN		UDP
Debug printing		N/A (TCP)

If the chosen protocol is SCAIP, click **More** to enter additional SCAIP receiver information.

Name	Definition and Instructions	Data Values
Callback	<p>Mark if you want the monitoring station to call you back, in case of an emergency.</p> <p>Note: The cellular phone number must be entered. Refer to 6.1.5 Phone Number on page 82.</p>	Enable = Marked, Disable = Clear Optional
Request Delay	Select time frame in seconds for request delay.	Dropdown list: 5 - 19 seconds Optional
Request Duration	Select time frame in minutes for request duration.	Dropdown list: Values: 3 - 9 minutes Optional


Name	Definition and Instructions	Data Values
Server User Name	Enter your user name on your server.	Required If SCAIP, authentication password is entered. Refer to 6.1.6 Authentication Password on page 83.
SCAIP Flavor	Note: For ES6502HC, this parameter is enabled and required for service types Family or Pro with Panel CPU Version 2.05.6.15 or later	Drop down list values: <ul style="list-style-type: none"> ● Default ● ENAI - Panel CPU version 2.05.6 and later ● VITAL - Panel CPU version 2.09.1 ● Answerlink - Panel CPU version 2.09.5 and later ● NEAT - Panel CPU version 2.11 and later ● VICTRIX - Panel CPU version 2.11 and later

Click **Update** to add the IP addresses and more to the destination list of the **Dialer** module along with the attributes entered for each destination.

7.3. Destination Definitions for a C7000 CP- Pre Version 6.1

Functionality of the C7000 CP is defined in relationship to Care@Home™ version 6.0. For example, prior to Care@Home™ version 6.0, only GPRS communication channels were available for **Message** dialers.


The following shows C7000 dialers prior to Care@Home™ version 6.1:



NOTE: Prior to Care@Home™ version 6.1, the C7000 supported only the PERS service type.

The following sections discuss functionality in relation to Care@Home™ pre-version 6.1.

7.3.1 Defining Communications Types and Channels



The following table describes the basic attributes required to add a destination to the Destination list.

Name	Definition and Instructions	Data Values	C7000 Before ver 6.1
ID	Internal identification of a dialer	Values: 1-16	✓
Phone Num / IP: Port	The phone number or IP/URL plus port number entered for the destination ID receiver address	Receiver address	✓
Comm Type	Select the communications type for the Destination. NOTE: Use the "Message" Comm Type only if the Care@Home™ system is configured to work with a PERS receiver.	Message	✓
		Voice	✓
		Debug Printing	✓
Comm Channels	Select the broadcast method to use with the communications type entered in the previous field.	Dropdown list changes according to the selected Comm Type .	✓

The following table lists the communication channels available for each communication type.

Comm Type	Description	C7000 Before ver 6.1 Comm Channels	Emergency Number
Analytics Message	For sending event messages to the Care@Home™ Server	<ul style="list-style-type: none"> ● GPRS ● GPRS 2G 	Emergency number allowed only for first 10 rows Note: Number erased only using delete key
Message	For sending event messages to the PERS receiver. Note: Use this Comm Type only if Care@Home™ is configured to work with a PERS receiver.	<ul style="list-style-type: none"> ● SIA IP ● SIA DC-09 ● SCAIP ● SIA ● Contact ID ● CPC ● Voice Signaling ● Debug Printing ● 4x2 ● BS 8521:2009 	Emergency number allowed only for first 10 rows Note: Number erased only using delete key
Data Stream	For sending files, that contain the detections from activity sensors, to the Care@Home™ server	<ul style="list-style-type: none"> ● GPRS ● GPRS 2G 	
Video	For sending photos taken by the camera detectors (IPD)	<ul style="list-style-type: none"> ● GPRS ● GPRS 2G 	
Voice	For calling the monitoring station Fail-Safe number	<ul style="list-style-type: none"> ● GSM ● PSTN 	Non-emergency phone number required Note: Number erased only using delete key
STUN	For sending messages, via UDP, to the Care@Home™ server that enables the CP to receive commands via UDP	<ul style="list-style-type: none"> ● GPRS ● GPRS 2G 	
Debug Printing	For Essence technical support only	<ul style="list-style-type: none"> ● GPRS ● GPRS 2G ● PSTN 	

7.3.2 Selecting Communication Channels

The available **communication channels** are:

- **PSTN** – Enter a phone number.
- **LTE**– Enter either an IP address or a domain address (URL) with a port number, such as 3061.
- **GSM** – Enter a phone number

To select your **communication channel**, choose the broadcast method that matches your CP communications method.

7.3.2.1 Configuring PSTN Communications Method

For PSTN communications method, the available broadcast methods are:

- SIA
- Contact ID
- Voice Signaling
- CPC

When choosing any of these methods, the **Dialer** module prompts you to enter the landline phone number of the monitoring station receiver.

When choosing **Voice Signaling**, the **Dialer** module prompts you to enter the **direct** phone number of the monitoring station. The **Voice Signaling** mode is used to ensure that the medical alarm (MA) event messages reach the monitoring station. **Voice Signaling** dialer works only in half-duplex mode.

7.3.2.2 Configuring Cellular Communications Method


For cellular communications method, the following are the available cellular communication protocol methods, between the CP and the monitoring station:

- SIA IP
- SCAIP
- SIA DC-09
- BS 8521


When choosing the cellular interfaces, the **Dialer** module prompts you to enter receiver address information, except for BS 8521.

7.3.3 Defining Receiver Address Attributes

Most of the communication channels are based on receiver address attributes.



For the **Message** communication type, the interfaces require the same attributes except for **SCAIP**.



The following table lists and defines the receiver address attributes.

Name	Definition and Instructions	Data Values	Required/Optional
Domain Name	Mark if a URL is required rather than an IP address.	Enable = Marked Disable = Clear	Optional
URL	When Domain Name is marked, enter the DNS that refers to the IP address of the monitoring station receiver. NOTE: Domain name can be up to 32 ASCII characters only	Character “<domain name>” For example: “www.essence-grp.com”	Required if Domain Name marked
IP address	When Domain Name is not marked, enter the IPv4 receiver address.	A set of 4 3-digit strings Values from 0- 255	Required
Port*	Enter the port number used by the monitoring station receiver. The default ports are listed in the table below.	Integer (5) Values: 1-65536	Required
Protocol**	NOTE: Not applicable for Debug Printing Select the IP communications protocols: ● TCP ● UDP	Dropdown list depends on communication channel	Required

Name	Definition and Instructions	Data Values	Required/Optional
Emergency Number	<p>NOTE: Comm Type Analytics Message and Message when using SIA IP, SIA DC-09, or SCAIP.</p> <p>Enter a phone number to call in case of an emergency, after sending the emergency event to the receiver and/or the Care@Home™ backend.</p>	<p>19-character string. Accepts numeric digits, “+” and “#” only.</p> <p>NOTE: Available only for the first 10 rows</p> <p>Note: Number erased only using delete key</p>	Optional

* - The following receiver address ports are defined as defaults for the communication channels:

Communication Channel	Default Port
SIA-IP	3061
SCAIP	5060 – default SIP port
SIA DC-09	Clear – no default
Analytics Alert	10060
STUN	3478
Video	5666
Data stream	
Debug printing	

** - The following receiver address protocols are allowed for the communication type and/or channels:

Communication Type	Communication Channel	Receiver Address Protocols
		C7000
Analytics Message Alert		TCP/UDP
Message	SIA IP	UDP
	SCAIP	UDP
	SIA DC-09	TCP/UDP
Data stream		TCP/UDP
Video		TCP/UDP
Voice		N/A
STUN		UDP
Debug printing		N/A (TCP)

If the chosen protocol is SCAIP, click **More** to enter additional SCAIP receiver information.


Name	Definition and Instructions	Data Values
Callback	<p>Mark if you want the monitoring station to call you back, in case of an emergency.</p> <p>Note: The GSM SIM number must be entered. Refer to 6.1.5 Phone Number on page 82.</p>	Enable = Marked, Disable = Clear Optional
Request Delay	Select time frame in seconds for request delay.	Dropdown list: 5 - 19 seconds Optional
Request Duration	Select time frame in minutes for request duration.	Dropdown list: Values: 3 - 9 minutes Optional
Server User Name	Enter your user name on your server.	Required If SCAIP, authentication password is entered. Refer to 6.1.6 Authentication Password on page 83

Name	Definition and Instructions	Data Values
SCAIP Flavor	Select the advanced classification element (flavor) that classifies how the Social Care Alarm Internet Protocol (SCA-IP) sends event messages between the CP and the receiver over the communication network.	Drop down list values: Default ENAI VITAL - CP SW version >= 6.1.4 Answerlink - CP SW version >= 6.1.4

Click **Update** to add the IP addresses and more to the destination list of the **Dialer** module along with the attributes entered for each destination.

7.4. Destination Definitions for a C7000 CP - Version 6.1 and Later

Functionality of the C7000 CP is defined in relationship to Care@Home™ version 6.1. The following shows C7000 dialers for Care@Home™ version 6.1 and later:



Version 6.1 supports the following communication channel types, known as interfaces:

- PSTN
- LTE (4G)
- Ethernet


Multiple communication interfaces can be enabled simultaneously.



NOTE: The **Message** communication type in this control panel is redefined as multiple interfaces with multiple protocols for each interface.

The following sections discuss functionality in relation to Care@Home™ version 6.1.

7.4.1 Defining Communications Types and Interfaces



The following table describes the basic attributes required to add a destination to the Destination list.

Name	Definition and Instructions	Data Values	C7000 Ver 6.1 and later	
			PERS-E & PERS	Family & Pro
ID	Internal identification of a dialer	Values: 1-16	✓	
Phone Num / IP: Port	The phone number or IP/URL plus port number entered for the destination ID receiver address	Receiver address	✓	
Comm Type	Select the communications type for the Destination. NOTE: Use the "Message" Comm Type only if the Care@Home™ system is configured to work with a PERS receiver.	Analytics Message	✓	
		Message		✓
		Data Stream	✓	
		Video	✓	
		Voice		✓
		STUN	✓	
Interfaces	Select the broadcast method to use with the communications type entered in the previous field.	Debug Printing		✓
		Dropdown list changes according to the selected Comm Type .		✓

The following table lists the interfaces available for each communication type.

Comm Type	Description	C7000 Interfaces Ver 6.1 and later	Emergency Number
		Family, Pro, PERS & PERS-E	
Analytics Message	For sending event messages to the Care@Home™ Server	<ul style="list-style-type: none"> ● LTE ● Ethernet 	Emergency number allowed only for first 10 rows Note: Number erased only using delete key

Comm Type	Description	C7000 Interfaces Ver 6.1 and later		Emergency Number
		Family, Pro, PERS & PERS-E		
Message	For sending event messages to the PERS receiver. Note: Use this Comm Type only if Care@Home™ is configured to work with a PERS receiver.	● Interfaces	● Protocols	Emergency number allowed only for first 10 rows Note: Number erased only using delete key
		● LTE (4G)	● SIA DC09 ● SIA IP ● SCAIP ● BS 8521	
		● Ethernet	● SIA DC09 ● SIA IP ● SCAIP	
		● PSTN	● SIA DC03 ● Contact ID ● Voice Signaling ● PSTN CPC ● BS 8521	
Data Stream	For sending files, that contain the detections from activity sensors, to the Care@Home™ server	● LTE ● Ethernet		
Video	For sending photos taken by the camera detectors (IPD)	● LTE ● Ethernet Note: Only for Family & Pro		
Voice	For calling the monitoring station Fail-Safe number	● GSM ● PSTN Note: Only for Family, Pro, & PERS-E		Non-emergency phone number required Note: Number erased only using delete key
STUN	For sending messages, via UDP, to the Care@Home™ server that enables the CP to receive commands via UDP	● LTE ● Ethernet Note: Supports all interfaces using IP receiver address		
Debug Printing	For Essence technical support only	● LTE ● Ethernet ● PSTN		

7.4.2 Selecting Interfaces by CP Type

The available **interfaces** are the same for either communication type method:

- **PSTN** – Enter a phone number.
- **LTE**– Enter either an IP address or a domain address (URL) with a port number, such as 3061.
- **Ethernet** - Enter either an IP address or a domain address (URL) with a port number, such as 5060.

To select your **Interface**, choose the broadcast method that matches your CP communications method.

7.4.2.1 Configuring PSTN Communications Method

For PSTN communications method, the available broadcast methods are:

- SIA

- Contact ID
- Voice Signaling
- CPC
- BS8521

When choosing any of these methods, the **Dialer** module prompts you to enter the landline phone number of the monitoring station receiver.


When choosing **Voice Signaling**, the **Dialer** module prompts you to enter the **direct** phone number of the monitoring station. The **Voice Signaling** mode is used to ensure that the medical alarm (MA) event messages reach the monitoring station. **Voice Signaling** dialer works only in half-duplex mode.

7.4.2.2 Configuring Cellular Communications Method

For cellular communications method, the following are the available cellular communication protocol methods, between the CP and the monitoring station:

- SIA IP
- SCAIP
- SIA DC-09
- BS8521

When choosing the cellular interfaces, the **Dialer** module prompts you to enter receiver address information.



7.4.2.3 Configuring Ethernet Communications Method


For the ethernet communications methods the following protocols are available:

- SIA IP


- SCAIP
- SIA DC-09

7.4.3 Defining Receiver Address Attributes

Most of the interfaces are based on receiver address attributes.



For the **Message** communication type, the interfaces require the same attributes except for **SCAIP**.



The following table lists and defines the receiver address attributes.

Name	Definition and Instructions	Data Values	Required/Optional
Domain Name	Mark if a URL is required rather than an IP address.	Enable = Marked Disable = Clear	Optional
URL	When Domain Name is marked, enter the web address that refers to the IP address of the monitoring station receiver. NOTE: Domain name can be up to 32 ASCII characters only	Character "<domain name>" For example: "www.essence-grp.com"	Required if Domain Name marked

Name	Definition and Instructions	Data Values	Required/Optional
IP address	When Domain Name is not marked, enter the IPv4 receiver address.	A set of 4 3-digit strings Values from 0- 255	Required
Port*	Enter the port number used by the monitoring station receiver. The default ports are listed in the table below.	Integer (5) Values: 1-65536	Required
Protocol**	NOTE: Not applicable for Debug Printing Select the IP communications protocol. <ul style="list-style-type: none"> <input type="radio"/> TCP <input type="radio"/> UDP 	Dropdown list depends on communication channel	Required
Emergency Number	NOTE: Comm Type Analytics Message and Message with GSM enabled. For C7000, GSM enabled is not required Enter a phone number to call in case of an emergency, after sending the emergency event to the receiver and/or the Care@Home™ backend.	19-character string. Accepts numeric digits, "+" and "#" only. NOTE: Available only for the first 10 rows Note: Number erased only using delete key	Optional

* - The following receiver address ports are defined as defaults for the interfaces:

Communication Types and Protocols	Default Port
SIA-IP	3061
SCAIP	5060 – default SIP port
SIA DC-09	Clear – no default
Analytics Alert	10060
STUN	3478
Video	5666
Data stream	
Debug printing	

** - The following receiver address protocols are allowed for the communication type and/or channels:

Communication Type	Interfaces/Protocols	Receiver Address Protocols
		C7000 Version 6.1 and later
Analytics Message Alert		TCP/UDP
Message	SIA IP	UDP
	SCAIP	UDP
	SIA DC-09	TCP/UDP
Data stream		TCP/UDP
Video		TCP/UDP
Voice		N/A
STUN		UDP
Debug printing		N/A (TCP)

If the chosen protocol is SCAIP, click **More** to enter additional SCAIP receiver information.

Name	Definition and Instructions	Data Values
Callback	Mark if you want the monitoring station to call you back, in case of an emergency. Note: The GSM SIM number must be entered. Refer to 6.1.5 Phone Number on page 82.	Enable = Marked, Disable = Clear Optional
Request Delay	Select time frame in seconds for request delay.	Dropdown list: 5 - 19 seconds Optional
Request Duration	Select time frame in minutes for request duration.	Dropdown list: Values: 3 - 9 minutes Optional
Server User Name	Enter your user name on your server.	Required If SCAIP, authentication password is entered. Refer to 6.1.6 Authentication Password on page 83.
SCAIP Flavor	For service types Family or Pro with panel CPU version 2.05.6.15 or later	Drop down list values: Default ENAI VITAL -CP SW version >= 6.1.4 Answerlink - CP SW version >= 6.1.4

Click **Update** to add the IP addresses and more to the destination list of the **Dialer** module along with the attributes entered for each destination.

7.5. Managing the Destination List


The **Dialer** module allows you to manage the **Destination** list:

- Add a new destination in a blank row in the list of destinations.
- Update the attributes of an existing destination.
- Delete an existing destination from the **Destination** list.

7.5.1 Adding a Dialer Destination

To add a new destination to the **Destination** list:

1. Select an empty row in the **Destination** list table.
2. Enter the destination information in the fields to the right of the **Destination** list table:
 - Comm Type
 - Comm Channel
 - Receiver Address Information
3. Right-click on the row selected in Step 1. A window appears displaying the functions available to manage the **Destination** list.



4. Click **Update** to add the new destination to the **Destination** list. The selected row is populated with the information entered in Step 2.

7.5.2 Updating a Dialer Destination

To modify an existing destination in the **Destination** list.

1. Select a row in the **Destination** list with an existing destination.
2. Revise the destination information in the fields to the right of the **Destination** list.
3. Right-click on the row selected in Step 1. A window appears displaying the functions available to manage the **Destination** list.
4. Click **Update**. The selected row is populated with the information revised in Step 2.

7.5.3 Deleting a Dialer Destination

To delete an existing destination from the **Destination** list:


1. Right-click on the row, in the **Destinations** window, that you want to delete from the Destination list. A window appears displaying the functions available to manage the destination list.
2. Click **Delete**. A window appears prompting you to confirm.
3. Click **Yes**.

7.6. The Dialing Sequence

The table below lists the communications type assigned to which dialing sequences.

Dialing Sequence	Communications Type
Message	<ul style="list-style-type: none"> ● Analytics Alert (not PERS) ● PERS Supervisory/Alert ● For ADL firmware – Panel CPU version 2.11 and later or PERS firmware – Panel CPU version 2.17 and later <ul style="list-style-type: none"> ■ PERS User Alert ■ PERS External Alert ■ PERS Technical Alert
Video	Video (Family and Pro only)
STUN	STUN (not PERS)
Data	Data Stream for Activity Files (not PERS)
	Debug Printing for Debug Printing (not PERS)
Fail-Safe Number	Voice – <ul style="list-style-type: none"> ● C7000 – not PERS ● ES6502HC: <ul style="list-style-type: none"> ■ Panel CPU version 2.11 through 2.14 ■ Family/Pro and panel CPU version 2.05.2 and later


CMS verifies that the destinations chosen for the dialing sequence are of the same communications type as the dialing sequence protocol.




ID	Phone Num / IP:PORT	Comm Type	Comm Channel	Emergency Number
1	192.168.005.105:1348	Analytics Message	GPRS	9725367840
2	212.056.122.048:1648	Analytics Message	GPRS	9725476387
3	www.amazon.co.uk:7845	Analytics Message	GPRS	00140188224
4	204.111.065.035:2315	Data Stream	GPRS	
5	www.pashbar.co.il:2584	Video	GPRS	
6	212.029.208.141:5666	Data Stream	GPRS	0
7	www.bingo.com:4718	STUN	GPRS	
8	212.029.208.141:3478	STUN	GPRS	0
9	111.185.222.010:6225	Video	GPRS	

For example, the destination list includes destinations of different communications types. Destinations 1, 2 and 3 are of communications types valid for the **Analytics Message** dialing sequence.

When choosing a destination with a communications type invalid for the selected dialing sequence, an error message appears.



Each communications type has its own protocols.



The following table lists the protocols assigned to each dialing sequence type:

Dialing Sequence	Protocol	ES6502HC		C7000	
		PERS	Family or Pro	Not PERS	All
Message	Analytics Alert	N/A		√	
	PERS Supervisory		PERS prefix added		PERS prefix added
	PERS Alert		PERS prefix added		PERS prefix added
	User Alert	Panel CPU version 2.11 through 2.14 or 2.17 and later	Versions 2.11-2.13.xx	N/A	N/A
	Technical Alert	Panel CPU version 2.11 through 2.14 or 2.17 and later	Versions 2.11-2.13.xx	N/A	N/A
	External Alert	Panel CPU version 2.11 through 2.14 or 2.17 and later	Versions 2.11-2.13.xx	N/A	N/A
Video	Comfort		√	√	
STUN	STUN		√	√	
Data	Activity Files	N/A	√	√	
	Debug Printing	Panel CPU Ver 2.03 or later	√		√
Fail-Safe Number	Failsafe	N/A	GSM enabled with Panel CPU Version 2.05.2 or later	√	Panel CPU version 2.11 through 2.14

View the dialing sequences and their protocols:

Analytics Message with Dual Mode

Dial Sequence

Message | Video | STUN | DATA | Fail-Safe Number

ID	Protocol	Sequence	Retries
1	PERS Supervisory	, 1, 2	1
2	PERS User Alert	, 1, 1, 1, 2, 2, 2	1
3	PERS External Alert	, 1, 1, 1, 2, 2, 2	1
4	PERS Technical Alert	, 1, 1, 1, 2, 2, 2	1
5	Analytics Alert	, 3	4

Protocol: PERS Technical Alert

Sequence:

No. of Retries:

Video Dialing Sequence Protocol

Dial Sequence

Message | Video | STUN | DATA | Fail-Safe Number

ID	Protocol	Sequence	Retries
1	Comfort	, 14, 14	1

Protocol: Comfort

Sequence:

No. of Retries:

STUN Dialing Sequence Protocol

Dial Sequence

Message | Video | STUN | DATA | Fail-Safe Number

ID	Protocol	Sequence	Retries
1	STUN	9	1

Protocol: STUN

Address:

No. of Retries:

Data Dialing Sequence

Dial Sequence

Message | Video | STUN | DATA | Fail-Safe Number


ID	Protocol	Sequence	Retries
1	Activity Files	, 6	3
2	Debug Printing		

Protocol: Activity Files

Sequence:

No. of Retries:

Fail Safe Dialing Sequence




Select a protocol and choose the destinations listed in the **Destination** list window for the relevant dialing sequence, as defined below.

Name	Definition and Instructions	Data Values	Required or Optional	Default Value
Sequence (Address for Fail Safe and STUN)	Enter the IDs for the relevant destinations, in the order they must be contacted. Refer to the list of Destination ID numbers. NOTE: Only choose an ID number for which an IP address or a URL is listed.	Number Values: 1-16	Required	None
No. of Retries	Enter the number of times to retry the sequences until at least one destination responds. NOTE: Recommended: Configure STUN only for a single retry. Not available for Fail Safe.	Dropdown List Values: 1-8	Required	1

Click **Update** to add the dialing sequence to the selected protocol.

7.6.1 Dialing Sequence Dual Mode

Dual mode is when you configure your system to use both a third-party PERS receiver and the Care@Home™ server. In this configuration, the events supported by the PERS receiver, such as emergency alerts, are sent first to the PERS receiver and then to the Care@Home™ server.




For example, the CP sends events to the monitoring station using Dual Mode, as follows:


1. The CP sends the events, via the **PERS Alert** protocol, to the monitoring station PERS receiver according to the details of Destination ID #2. This may be an IP address and port number, domain name or telephone number.

2. The CP attempts to communicate with the receiver twice, according to the number of retries in the example.
3. If an acknowledgement is received and the event is an emergency alert, call the destination monitoring station emergency number. The **PERS Alert** succeeded.
4. Otherwise, if no acknowledgement is received, the **PERS Alert** failed.
5. The CP sends the events, via the **Analytics Alert** protocol, to the monitoring station Care@Home™ server, according to the details of Destination ID #3
6. The CP attempts to communicate with the Care@Home™ server for Destination ID #3.
7. If an acknowledgement is received and the event is an emergency alert:
 - If the **PERS Alert** failed, the CP calls the destination monitoring station emergency number.
 - If the **PERS Alert** succeeded, the CP does not call the destination monitoring station emergency number.
8. The dialing sequence in dual mode is completed.
9. If no acknowledgement is received, the CP sends the events, via the **Analytics Alert** protocol, to the monitoring station Care@Home™ server, according to the details of Destination ID #4.
10. The CP attempts to communicate with the Care@Home™ server for Destination ID #4.
11. If an acknowledgement is received and the event is an emergency alert:
 - If the **PERS Alert** failed, the CP calls the destination monitoring station emergency number.
 - If the **PERS Alert** succeeded, the CP does not call the destination monitoring station emergency number.
12. The dialing sequence in dual mode is completed.
13. If no acknowledgement is received, repeat steps 4 through 9, according to the number of retries, in this case, one additional retry.
14. If an acknowledgement is received and the event is an emergency alert:
 - If the **PERS Alert** failed, the CP calls the destination monitoring station emergency number.
 - If the **PERS Alert** succeeded, the CP does not call the destination monitoring station emergency number.
15. The dialing sequence in dual mode is completed.
16. If no acknowledgement is received using both **PERS Alert** and **Analytics Alert** protocols:
 - a. The CP calls the Fail-safe number.
 - b. The CP issues the following vocal announcement:
 - c. "Attention, communication with the monitoring station is lost."

8. The Activity Sensors Module



Click **Activity Sensors** on the **Main** navigator menu. The **Activity Sensors** module is the tool that allows management of devices for the Care@Home™ CP installed on the resident’s premises.



All the sensor devices connected to the CP are listed by identification number in this module. For each device, the following information is displayed:

- **Identification number**- The internal identification number assigned to the sensor device synchronized with the Care@Home™ web application
- **Serial Number**
- **Device Type**
- **Mode** – The characteristics of the device type functions, such as “Chime”, “Bypass”, and “Inactivity”
- **Status** – The status of the sensor device, such as “Tamper”, “Open” and “Low Bat”

8.1. Defining Sensor Device Attributes

The following table describes the attributes for the sensor devices.

Name	Definition and Instructions	Data Values	Required/Optional	Default Value
Device ID	The identification number assigned by the CP for each sensor device.	<ul style="list-style-type: none"> ● ES6502HC: <ul style="list-style-type: none"> ■ Panel CPU version 2.08 and earlier or 2.14 and later – 1-55 ■ Panel CPU version 2.08 and later – 1-20. ● C7000 - Integer from 1-20 	N/A	Display Only
RF	The radio frequency (RF) level of the selected sensor device. It is the received signal strength indicator (RSSI) level of the selected device.	Percent	N/A	Display Only
Device Type	The device type of the selected sensor device. The sensor devices are: <ul style="list-style-type: none"> ● MGLS – Door or bed sensor ● PIR – Motion detector ● IPD – Camera detector (optional) 	String	N/A	Display Only
Version	The version of the firmware of the device Note: Parameter enabled for ES6502HC – CP SW version >=2.03	Number	N/A	Display Only
Zone Label	Note: Parameter enabled for ES6502HC - CP SW version >= 2.0. Choose the label that describes the location where the sensor was installed. Note: The dropdown list includes both system-defined location labels and user-defined custom labels.	Dropdown list	Optional	None
First priority IPD	Note: Parameter enabled for Family and Pro only. Choose the first camera to photograph the events. Choose NONE for no camera coverage.	Dropdown list	Optional	None
Second priority IPD	Note: Parameter enabled for Family and Pro only. Choose the second camera to photograph the events. Choose NONE for no camera coverage.	Dropdown list	Required	None
Enable LEDs	Note: MGLS device only Mark to enable the device LEDs to light up when door is open or closed.	Enable = Marked Disable = Clear	Optional	Disable = Clear
Resets possible fall	Note: Pro only To prevent false fall detections, clear for sensors in a combined bathroom/restroom. Note: Parameter appears for support-level user only			

Name	Definition and Instructions	Data Values	Required/Optional	Default Value
Resets wandering alert	<p>Note: Parameter enabled for ES6502HC, for Family or Pro with Panel CPU Version 2.05 or later and for C7000, only for Family (not PERS-E) or Pro.</p> <p>When marked, prevents a detection from the device from triggering a wandering alert.</p> <p>Example: When there is more than one resident, we may want a wandering alert to be sent if the only movement is by the resident in the bedroom. In this case, clear the checkbox for the bedroom.</p> <p>Note: Parameter appears for support-level user only</p>			
Activity Type	<p>Choose the type of activity the sensor is configured to detect.</p>	<p>For PERS, only MGLS for Bed sensor or Unspecified</p> <p>For Family, Pro or PERS-E:</p> <ul style="list-style-type: none"> ● Refrigerator door ● Front door ● Bedroom ● Restroom ● Bathroom ● Living room ● Dining room ● Bed (ES6502HC CPU version 2.09.6 through 2.14 only) ● Unspecified 	Required	None
Activity Group	<p>Note: Parameter enabled for ES6502HC, for Family or Pro and for C7000, only for Family (not PERS-E) or Pro.</p> <p>Allows devices with the same activity type to be grouped.</p> <p>Note: Parameter appears for support-level user only</p>	<p>Enable = Marked</p> <p>Disable = Clear</p>	Optional	Disable = Clear

Name	Definition and Instructions	Data Values	Required/Optional	Default Value
Send Detections	<p>Note: Parameter enabled for ES6502HC, for Family or Pro and for C7000, only for Family (not PERS-E) or Pro. Mark to enable the sensor device to send all its event detections to the CP.</p>	<p>Enable = Marked Disable = Clear</p>	Optional	Disable = Clear
External Magnet	<p>Note: Parameter enabled for ES6502HC, for Family or Pro and for C7000, only for Family, PERS-E or Pro. NOTE: MGLS device only Mark to allow third-party magnetic devices on metallic surfaces to be installed in your Care@Home™, without effecting RF transmission.</p>	<p>Enable = Marked Disable = Clear</p>	Optional	Disable = Clear
(Enable modes)	<p>Activate or deactivate modes that alter the functionality of the activity sensor devices.</p> <p>Chime- Mark to sound the chime, if the sensor device detects activity.</p> <p>Bypass – Mark to ignore the detections sent by the sensor device, including the ability to sense extreme temperatures. The sensor device remains connected to the CP.</p> <p>Note: Bypass Extreme Temperature can be marked even if Bypass is clear.</p> <p>Inactivity – Mark to enable the sensor to send an Inactivity event to the monitoring station, when no activity occurs during a given time-frame.</p>	<p>Enable = Marked Disable = Clear</p> <p>For the following set of checkboxes:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Chime <input type="checkbox"/> Bypass <input type="checkbox"/> Bypass Extreme Temperature <input type="checkbox"/> Inactivity 	Optional	None

Click **Update** to add the sensor device attributes to the CP.

8.2. Managing Activity Sensors



NOTE: For ES6502HC, you can manage activity sensors from CP SW version 2.03.

The **Activity Sensors** module allows you to manage the Activity Sensor list manually:

- Add a device to the list.
- Update the attributes of a device in the list.
- Delete a device from the list.


To delete a device, you require privilege access usually reserved for a system administrator or Essence technical support personnel.

8.2.1 Adding an Activity Sensor

Below is an illustration of adding an activity sensor to the Activity Sensors list.

To add an activity sensor to the **Activity Sensors** list:

1. Right-click on the **Activity Sensors** list. A window appears displaying the functions available to manage the list.



2. Click **Add**. The **Device Details** window appears.
3. If you have a barcode reader:
 - a. Mark Use Barcode.
 - b. Using the barcode reader, scan the barcode on the device packaging. The scanned data appears in the **Barcode** field. The serial number and device type from the scanned data are loaded to their respective fields.
4. Refer to 2.3.2 Understanding the *setup.ini* File on page 17.
5. If you do not have a barcode reader:
 - a. Clear Use Barcode.
 - b. Enter the last eight characters on the barcode label.
 - c. Select the device type from the dropdown list.
6. Click Add Device.


8.2.2 Updating an Activity Sensor

To update a row in the **Activity Sensor** list with an existing activity sensor:

1. Select the row in the **Activity Sensor** list with the activity sensor requiring update.

ID	Serial Number	Device Type	Mode	Status	Label
1	70091E64	MGLS		NL	Kitchen
2	700A037A	MGLS		NL	Main Entry
3	002134BC	PIR		NL	Bedroom
4	002133D5	PIR		NL	Restroom
5	00213476	PIR		NL	Bathroom
6	00582B41	IPD	I	NL	Living Room
7	002133DC	PIR	I	NL	Dining Room
8	09AB9D32	IPD		NL	Hall
9	124ABC68	MGLS		NL	Choose Name
10	34374BC8	IPD		NL	Choose Name

2. Modify the values of the activity sensor attributes requiring revision.
3. Click **Update** in the function window or at the bottom of the screen.



8.2.3 Deleting an Activity Sensor


To delete an existing activity sensor, from the **Activity Sensor** list, includes removing all the data that has been accumulated for that sensor. This involves “erasing” all the event history not only for the activity sensor requiring deletion but also for all the activity sensors installed on the premises.

To delete an activity sensor from the **Activity Sensor** list:

1. Right-click on the row in the Activity Sensor list with the activity sensor requiring deletion. A window appears asking to confirm.
2. Click **Delete**. A window appears prompting you to confirm.
3. Click **Yes**.

8.3. Zone Status


The **Zone Status** box displays the status (normal = **Green**; problem = **Red**) of the selected sensor device.




The following is the list of the available statuses:

- **Tamper** – The selected sensor device was touched, moved, or installed incorrectly.
- **Low Battery** – The battery of the selected sensor device requires replacement.
- **Bypass** – The selected sensor device appears to be deactivated.
- **Open** –MGLS devices, “open/closed”
- **Supervisory** – Whether a problem with the selected sensor device has occurred. If communication between the CP and the device has been lost.
- **Learned** – Whether the selected sensor device has been paired with the control panel.

9. The Safety & SOS Module



Click **Safety & SOS** on the **Main** navigator menu. The **Safety & SOS** module is the tool that allows management of the safety and emergency devices for the Care@Home™ CP installed on the resident’s premises.



All the safety and emergency devices connected to the CP are listed by identification number in the Safety & SOS devices list. For each device, the following information is displayed:

- **ID** - The internal identification number assigned to the emergency device synchronized with the Care@Home™ web application

- **Serial Number**
- **Device Type**
- **Mode** – Whether to disregard notifications broadcast by the device

9.1. Defining Attributes for Safety & SOS Devices

For each safety and emergency device, there are attributes that define:

- General, common device functionality
- Functionality specific and unique for that device type

The following table describes the general attributes of the safety and emergency devices.


Name	Definition and Instructions
Device ID	The identification number assigned by the CP for each safety or emergency device. <ul style="list-style-type: none"> ● C7000 – 1-10 ● ES6502HC – up to version 2.03 – 1-8 ● ES6502HC – version 2.03 to 2.03.2 – 1-16 ● ES6502HC – version >= 2.03.2 – 1-10
RF	The RF level of the selected safety and emergency device relative to the received signal strength indicator (RSSI) scale. <p>Note: An RF level of 80% and above is recommended. If the RF-level is lower than 50%, install the device in a better reception area, or move the control panel.</p>
Hardware Type	The code for the hardware type of the selected device: <ul style="list-style-type: none"> ● MDS – Multidimensional fall detector ● SPB – Stationary panic button ● EP and EPP – Emergency pendant devices (panic button) ● SK2 - Smoke detector ● VPD – Voice panic detector <p>Note: For VPD, the hardware revision number is displayed</p> <ul style="list-style-type: none"> ● EPA – Emergency pendant with fall detection ● MGLS – Magnetic transmitter input device that connects non-RF or other devices to Care@Home™, such as a water detector
Device Type	<p>Note: For ES6502HC, attribute enabled only if Hardware Type is MGLS and CP version >= 2.03</p> The description of the selected safety or emergency device. When Hardware Type is MGLS, available device type options are: <ul style="list-style-type: none"> ● Gas Detector ● Fire Detector ● Water Detector ● Panic Device

Name	Definition and Instructions
Version	<p>Note: For ES6502HC, attribute displayed only if Hardware Type is MGLS and when CP version >= 2.03</p> <p>The firmware version installed on the device.</p> <p>Note: For the VPD and MDS, two firmware versions are displayed:</p> <ul style="list-style-type: none"> ● Host ● VPD – Voice detection ● MDsense – fall detection

9.2. Specific Safety and Emergency Attributes


Each safety and emergency device type attributes define how the device alerts the resident, caregiver, and/or monitoring station, when a situation occurs that might threaten the safety or welfare of the resident.

MDS Attributes



Device ID: 1 RF: N/A
 Hardware Type: MDS
 Device Type: Fall Detector
 Version: T
 Device Label: Choose Name
 Bypass Mode
 Bypass Extreme Temperature
 First priority IPD: NONE
 Second priority IPD: NONE
 Enabled LED
 Update

VPD Attributes



Device ID: 2 RF: 98 %
 Hardware Type: VPD < Rev. 3 >
 Device Type: VPD
 Version: T 03.03
 04.03
 Location: Other Room
 Device Label: Choose Name
 Bypass Mode
 Enable Answer Call
 Bypass Extreme Temperature
 Enable VPD LED
 Score for at least one phrase: 6
 Score for both phrases: 6
 Voice Triggering: Enable
 Update

EP and EPP Attributes

Device ID: 1	RF: 80 %
Hardware Type: EP	
Device Type: Panic Device	
Version: 1.01	
Device Label: Choose Name	
<input type="checkbox"/> Bypass Mode	
<input checked="" type="checkbox"/> Enable Answer Call	
<input type="checkbox"/> Bypass Extreme Temperature	
First priority IPD: NONE	
Second priority IPD: NONE	
Update	

Device ID: 1	RF: 80 %
Hardware Type: EPP	
Device Type: Panic Device	
Version: 1.06	
Device Label: Choose Name	
<input type="checkbox"/> Bypass Mode	
<input checked="" type="checkbox"/> Enable Answer Call	
<input type="checkbox"/> Bypass Extreme Temperature	
First priority IPD: NONE	
Second priority IPD: NONE	
Update	

EPA Attributes

Device ID: 2	RF: N/A
Hardware Type: EPA	
Device Type: EPA	
Version:	
Device Label: Choose Name	
<input type="checkbox"/> Bypass Mode	
<input checked="" type="checkbox"/> Enable Answer Call	
<input type="checkbox"/> Bypass Extreme Temperature	
Update	


SPB Attributes

Device ID: 7	RF: N/A
Hardware Type: SPB	
Device Type: SPB	
Version: 1.01	
Device Label: Cash Register	
<input type="checkbox"/> Bypass Mode	
<input checked="" type="checkbox"/> Enable Answer Call	
<input type="checkbox"/> Bypass Extreme Temperature	
First priority IPD: 16	
Second priority IPD: 9	
<input checked="" type="checkbox"/> Enable SPB LED	
Update	

Smoke Detector Attributes

Device ID: 3	RF: 80 %
Hardware Type: SK2	
Device Type: Smoke Detector	
Version: 1.01	
Device Label: Kitchen	
<input type="checkbox"/> Bypass Mode	
<input type="checkbox"/> Bypass Extreme Temperature	
First priority IPD: 1	
Second priority IPD: 2	
Safety Standard: EN	
Update	

Water Detector Attributes



The following table describes the attributes specific to the hardware types of the safety and emergency devices.

Name	Definition and Instructions	Data Values	Required/Optional	Factory Settings
Location	<p>Note: VPD only</p> <p>Choose the label that describes the location where the VPD is installed. Used to inform emergency response operator.</p> <p>Note: For ES6502HC, enabled for panel CPU version 2.08 through 2.14.</p>	Dropdown list	Required	Unspecified
Device Label	<p>Note: Disabled for the EP, EPP, or the EPA devices</p> <p>Choose the label that describes the location where the safety or emergency device is installed.</p> <p>Note: The dropdown list is the same as the dropdown list for the Zone Label attribute, in the Activity Sensors module. The list includes both system-defined location labels and user-defined custom labels.</p> <p>Refer to 12 The Custom Labels Module on page 157.</p> <p>Note: Enabled for:</p> <ul style="list-style-type: none"> ● C7000 ● ES6502HC – CP SW version >= 2 	Dropdown list	Required	None
<Bypass options>	<p>Mark to activate or deactivate modes that alter the functionality of the safety and emergency devices.</p>			
	<p>Bypass mode</p> <p>Mark to deactivate the sensor device. The sensor device remains connected to the CP.</p>	Disable= Clear	Optional	Disable = Clear
	<p>Enable Answer Call via Buttons</p> <p>NOTE: Not for MDS, SK2 or the MGLS with a non-Panic device</p> <p>Mark to answer incoming calls using the panic button.</p>			Enable = Marked
	<p>Bypass Extreme Temperature</p> <p>Mark to enable the CP to send extreme temperature reports to the monitoring station, if temperature readings exceed the defined thresholds.</p>			Disable = Clear

Name	Definition and Instructions	Data Values	Required/Optional	Factory Settings
First Priority IPD	<p>NOTE: Not for EP, EPP and EPA</p> <p>Choose the first camera to photograph the events. Choose NONE to disable the camera option.</p> <p>Note: Parameter enabled for ES6502HC, for Family or Pro and for C7000, not PERS</p>	Dropdown list	Optional	None
Second Priority IPD	<p>NOTE: Not for EP, EPP, and EPA</p> <p>Choose the second camera to photograph the events. Choose NONE to disable the camera option.</p> <p>Note: Parameter enabled for ES6502HC, for Family or Pro and for C7000, not PERS</p>	Dropdown list	Optional	None
Enable LEDs Enable VPD LED Enable SPB LED	<p>NOTE: MDS, SPB, VPD, and MGLS only</p> <p>Mark to activate the LED display for the device. When the device is triggered, the LED is lit. If disabled, the device continues to function. However, the LED is not lit.</p> <p>For the VPD and SPB, the factory setting is enabled (Marked).</p>	Enable = Marked Disable= Clear	Optional	Disable = Clear
Score for at least one phrase	<p>NOTE: VPD only</p> <p>For the VPD to trigger an emergency alarm, the emergency phrase must be:</p> <ul style="list-style-type: none"> ● Detected twice, and ● At least one occurrence of the phrase must match this setting, and ● Detect Phrase Once is cleared <p>Refer to 5.4.9.1 VPD Quality Configuration on page 66</p>	1-9		
Score for both phrases	<p>NOTE: VPD only</p> <p>For the VPD to trigger an emergency alarm, the emergency phrase must be:</p> <ul style="list-style-type: none"> ● Detected twice, and ● Both occurrences of the phrase must at least match this setting <p>Refer to 5.4.9.1 VPD Quality Configuration on page 66</p> <p>Note: Applies to C7000 and ES6502HC with panel CPU version 2.08 through 2.14.</p> <p>For earlier versions or for PERS firmware, panel CPU version 2.08 through 2.14, CMS assigns the parameter the same value as Score for at least one phrase</p>	1-9		
Safety Standard	<p>NOTE: SK2 Only:</p> <p>Select the safety standard to apply to the smoke detector. The standards relate to the event reporting frequency. The safety standards are:</p> <p>Default –out-of-the box factory settings designed to preserve battery life.</p> <p>EN – European standards</p> <p>UL – USA standards</p> <p>ULC – Canadian standards</p>	Dropdown list: <ul style="list-style-type: none"> ● Default ● EN ● UL ● ULC 	Required	Default

Name	Definition and Instructions	Data Values	Required/Optional	Factory Settings
Enable Fall Detection	<p>Note: EPA and MDS only. For Essence technical support</p> <p>Mark to allow the device to detect falls and send fall alerts and usage data to the CP.</p> <p>When adding an EPA, the default value is determined by the Enable for New Device parameter in the EPA Config tab in the Control Panel Configuration Extension window.</p> <p>Refer to 5.4.10 EPA Configuration Parameters on page 68.</p> <p>The default for MDS is enabled.</p>	<p>Enable = Marked</p> <p>Disable= Clear</p>	Optional	
Enable debug mode	<p>NOTE: EPA and MDS only: For Essence technical support.</p> <p>Mark to allow the device to run in debug mode.</p>	<p>Enable = Marked</p> <p>Disable= Clear</p>	Optional	Disable = Clear
Voice Triggering	<p>NOTE: For VPD and MDS only. For Essence technical support</p> <p>Select ENABLE to activate the Intelligent Voice Activation™ (IVA) technology, to recognize the trigger phrase.</p>	<p>Dropdown list</p> <p>Enable</p> <p>Disable</p>	Required	Enable
Detector Delay	<p>NOTE: MDS Only. For use by Essence Technical Support.</p> <p>Select the number of seconds between fall and alarm.</p>	5-60 seconds	Required	

Click **Update** to save the safety and emergency device attributes to the CP.

9.3. Managing Safety & SOS Devices



NOTE: For ES6502HC, you can manage safety and SOS devices for PERS with **Panel CPU Version 2.03** or later.


The **Safety & SOS** module allows you to manage the **Safety & SOS** device list manually:

- Add a device to the list.
- Update the attributes of a device in the list.
- Delete a device from the list.

9.3.1 Adding a Safety or SOS Device

To add a safety or SOS device:

1. Right-click on the safety & SOS device list. A window appears displaying the functions available to manage the list.



2. Click **Add**. The **Device Details** window appears.
3. If you have a barcode reader:
4. Mark Use Barcode.
 - a. Using the barcode reader, scan the barcode on the device packaging. The scanned data appears in the **Barcode** field. The serial number and device type from the scanned data are loaded to their respective fields.
5. Refer to 2.3.2 Understanding the *setup.ini* File on page 17.
6. If you do not have a barcode reader:
 - a. Clear Use Barcode.
 - b. Enter the last eight characters on the barcode label.
 - c. Select the device type from the dropdown list.
7. Click Add Device.

9.3.2 Updating a Safety or SOS Device


To update a row in the Safety & SOS device list with an existing safety or emergency device:

1. Select the row in the Safety & SOS device list with the safety or emergency device requiring update.

2. Modify the values of the safety or emergency device attributes requiring revision.

ID	Serial Number	Device Type	Mode	Label
1	AB754CD6	Panic Device	Enable	
2	000001C2	VPD	Enable	Bathroom
3	BC485A34	Panic Device	Enable	
4	E363C167	Water Detector	Enable	Laundry
5	A3837D37	Smoke Detector	Enable	Kitchen
6	BC45D134	VPD	Enable	
7				

3. Click **Update** in the function window or at the bottom of the screen.



9.3.3 Deleting a Safety or SOS Device

To delete a safety or emergency device from the **Safety & SOS** device list:

1. Right-click on the row in the **Safety & SOS** device list with the safety or emergency device you want to delete. A window appears prompting you to confirm the delete request.
2. Click **Delete**. A window appears prompting you to confirm.
3. Click **Yes**.

9.4. RF Input Status

The RF input status reports the status of a selected safety or emergency device (normal = **Green**; problem = **Red**).




The following is the list of the available statuses:

- **Tamper** – The selected device was touched, moved, or installed incorrectly.
- **Low battery** – The battery of the selected device requires replacement.
- **Dead battery** – For EPA, EPP, and MDS devices only: The battery requires immediate replacement.
- **Active/Bypass** – Alert reports from the selected device are sent or ignored.
- **Supervisory** – Whether communication between the device and the control panel has been lost.
- **Learned** - Whether the selected device has been paired with the control panel.
- **External Power** – MDS only. Connected to external power. Lights red if MDS was connected to power and then disconnected.

10. Voice & LED Module



NOTE: For ES6502HC, the **Voice and LED module** parameters are available for PERS service type with **Panel CPU Version 2.03** or later.



Click **Voice & LED** on the **Main** navigator menu. The **Voice & LED** module allows you to configure:

- The receipt of vocal announcements for technical, medical, and safety issues
- The CP to “vocally” remind the resident about scheduled activities.
- The CP LED indicators
- The pendant test operational parameters

The vocal reminder is characterized by the following attributes:

- The purpose of the reminder announcement
- Additional information related to the purpose of the reminder
- The scheduled frequency of the reminder announcement (schedule type)
- The period of the reminder announcement
- The time of the reminder announcement


Vocal reminders use the following schedule types:

- **One Time**
- **Weekly**
- **Monthly (Day of the month)**
- **Monthly (Day of the week)**

For each schedule type, you can choose the “date” and time for the reminder to be announced. The time configuration for a reminder announcement is restricted to the hour with 15-minute increments.


Refer to 10.2 Reminder Schedule Types on page 144.

The structure for defining a reminder includes the parameters that define the attributes of a vocal reminder. For example, below illustrates the structure to define a weekly vocal reminder.



The following table describes the parameters that define the attributes of a vocal reminder.

Parameter	Definition/Instructions	Attributes
(Reminder) Enabled	Mark to activate the vocal reminder.	Enable = Marked Disable = Clear
(Schedule Type)	Select the type of schedule you require for the vocal reminder. Refer to 10.2 Reminder Schedule Types on page 144.	Dropdown list: <ul style="list-style-type: none"> ● One-Time ● Weekly ● Monthly (Day of Month) ● Monthly (Day of Week)
Report	Mark to send a Missed Reminder report to the monitoring station when a reminder is not acknowledged, or an EP test is incomplete.	Enable = Marked Disable = Clear
Period: Day and Time	Day: Enter either the day of the week, the day of the month, or the date when you require the reminder announcement.	The format of the date is dependent on the schedule type. Refer to 10.2 Reminder Schedule Types on page 144.
	Time: Enter the time when you require the reminder announcement	Time in <i>hh:mm</i> format
Days/Weeks	Mark one or more days/weeks when you require the reminder announcement.	Per day/week: Enable = Marked Disable = Clear
Set Default	Click to reset to default settings. The Reminder structure is reset to factory settings, disabled, and displayed in grey.	Push Button
Time to Announce	Enter the time to inform the resident when an event is scheduled.	Time in <i>hh24:mi</i> format





The following are the available types of vocal reminders:

- **Medication:** Refer to 10.3 Medication Reminders on page 146.
- **Appointments:** Refer to 10.4 Appointment Reminders on page 147.
- **Rides:** Refer to 10.5 Ride Reminders on page 148.
- **Alarms:** Refer to 10.6 Alarm Reminders on page 148.

- **Pendant Test:** Refer to 10.7 Pendant Test Reminder on page 149.

10.1. General Configuration Parameters

The **Voice & LED** module includes the global parameters that configure the scope of the vocal announcements and reminder definitions and allows you to manage the CP **Trouble**  LED and other LED indicators.



Voice & LED General Configuration

The following table defines the **Voice & LED** configuration parameters.

Field	Definition/Instructions	Attributes	Required/Optional
Silence Voice Announcements	Mark to silence the vocal announcements. Note: If enabled, all vocal announcement groups are silenced.	Enable = Marked Disable = Clear	Optional
Language	The language of the vocal announcements is determined by the voice file downloaded to the panel. Contact Essence for the voice file. Supported languages: <ul style="list-style-type: none"> ● English ● Spanish ● French ● Portuguese ● German ● Hebrew ● Chinese ● Slovenian ● Danish ● Swedish ● Finnish ● Norwegian ● Welsh ● Dutch – supported for C7000 and for ES6502HC with panel CPU version 2.05.4 through 2.14. 	N/A	Display Only
Enable Groups	Clear the checkboxes of the group of announcements you want to silence. Refer to Appendix B Voice Announcements on page 201.	Enable = Marked Disable = Clear	Optional
	Technical – Clear to silence vocal messages about technical problems such as tamper and low battery power.		
	Safety Alarm – Clear to silence vocal messages for safety alarms such as fire alarm and flood alarm.		
	Medical Alarm – Clear to silence vocal messages about panic and emergency alarms.		
	Vocal Reminders – Clear to silence vocal reminder messages.		
	Functional – Clear to silence vocal messages for operations performed using the CP and the peripherals.		
Vocal Reminders: Enable	Mark to allow the use of the Vocal Reminders feature.	Enable = Marked Disable = Clear	Optional
Send Missed Reminder report	Mark to allow the broadcast of Missed Reminder reports. Each reminder type has a Report parameter that, if marked, allows sending a Missed Reminder report to the monitoring station if a reminder is not acknowledged or the pendant test is incomplete. Note: If this parameter is not marked, none of the reminders can send Missed Reminder reports.	Allow = Marked Not allow = Clear	Optional
Timeout	Select the timeframe for the continued repetition of the reminder announcements. At the end of the timeframe, if the resident has not pressed RESET, a Missed Reminder report is sent to the monitoring station, if Report is marked. Note: For the pendant test this is the test interval within which the resident must press their emergency pendant button.	15 minutes or 30 minutes	Required

Field	Definition/Instructions	Attributes	Required/Optional
Snooze Time	Select the wait interval for rebroadcasting reminder announcements. Note: If you choose "No Repetition", the reminder announcement is not broadcast again.	Dropdown list values: No Repetition, 1 - 30 minutes, cannot be larger than Timeout.	Required

The following table defines the parameters for configuring the CP's **Trouble LED**.









NOTE: For ES6502HC, the **Panel's Trouble LED** parameters are available for **Panel CPU Version 2.03** or later.

Parameter	Definition/Instructions	
Disable	Mark to deactivate the CP Trouble LED . When an alert is issued for the alert categories, the LED remains unlit. Values: Marked = LED parameters disabled; Clear = LED parameters enabled.	
Enable Groups	Mark the group of trouble states for which you want to light the CP trouble LED . Values: Enable = Marked; Disable = Cleared	
	Group	Description
	Communication loss	When the CP loses communications with the monitoring station
	Supervisory loss	When the CP loses communications with a peripheral
	Panel is in low/missing battery	When the CP battery is low or missing
	Peripheral tamper is on	When a peripheral is tampered
	Panel's power supply in trouble	When there is trouble with your power supply
	Low peripheral battery	When a peripheral reports low battery
	Jamming detection	When RF communications detects jamming
Connectivity test failed	Requires panel with cellular interface. When the cellular connectivity test, that verifies communication between the CP and the Care@Home™ server, failed. For both panel types, for the PERS service type, this parameter is clear and disabled.	

The following table defines when to mark the CP LED indicators.



NOTE: For ES6502HC, the **Panel's Trouble LED** parameters are available for ES6502HC panel CPU version 2.05.2 through 2.14

LED	Definition	Attributes	Required / Optional
GSM Communication LED 	The LED lights up Green whenever the CP sends information to the Care@Home™ server.	Enable = Marked Disable = Clear	Optional
Power LED 	The LED lights up Green when the CP is connected to a power outlet. The LED blinks when working off the backup battery.		
Emergency Duration 	The LED lights up Red during a voice window for an emergency call. If disabled, the LED is lit Red for the duration of the voice window.		
	Enter a value greater than 0 for the time frame for the LED to be lit Red during the voice window.	Seconds	Required only if Emergency Duration is marked
Connectivity Test 	Note: For ES6502HC, enabled for 2.05.3 <= CP SW version < 2.14... Not enabled for PERS on any control panel. Requires panel with cellular capability. The Communication LED is Orange for the cellular connectivity test - Refer to ESUGSC003 Care@Home™ CP User Guide.	Enable = Marked Disable = Clear	Optional
Fixed Blue LED 	Note: Enabled for: <ul style="list-style-type: none"> ● ES6502HC: panel CPU version 2.05.3 through 2.14 ● C7000 The LED surrounding the button is lit Blue all the time.		
Blinking LED 	Note: Enabled for: <ul style="list-style-type: none"> ● ES6502HC: Panel CPU version 2.05.4 through 2.14 ● C7000 The LED surrounding the button blinks Blue during CP voice announcements.		

10.2. Reminder Schedule Types

The following topics describe the different schedule types of the voice announcements and how to apply the schedule types in setting up the voice announcements.

- **One Time:** Refer to 10.2.1 One-Time Schedule Type on page 145.

- **Weekly:** Refer to 10.2.2 Weekly Schedule Type on page 145.
- **Monthly (Day of Month):** Refer to 10.2.3 Day of the Month Schedule Type on page 145.
- **Monthly (Day of Week):** Refer to 10.2.4 Day of Week Schedule Type on page 145.

10.2.1 One-Time Schedule Type

For the **One-Time Schedule Type**:

1. Choose a specific date for the reminder to be announced.
2. Choose a specific time for the reminder to be announced.

On the chosen date, at the time specified, the reminder is announced.

10.2.2 Weekly Schedule Type

For the **Weekly Schedule Type**:

1. Choose a specific time for the reminder to be announced.
2. Choose one or more days of the week when the reminder is to be announced.

Every week, on the chosen day or days of the week, at the time specified, the reminder is announced.

10.2.3 Day of the Month Schedule Type

For the **Monthly (Day of Month) Schedule Type**:

1. Choose a specific day of a month for the reminder to be announced.
2. Choose a specific time for the reminder to be announced.

Every month, on the chosen day, at the time specified, the reminder is announced.

10.2.4 Day of Week Schedule Type

For the **Monthly (Day of Week) Schedule Type**:

1. Choose a day of the week for the reminder to be announced.
2. Choose a specific time for the reminder to be announced.
3. Choose the week or weeks of the month for the reminder to be announced.

Every month, during the chosen week or weeks, on the chosen day of the week, at the time specified, the reminder is announced.

10.3. Medication Reminders

The **Voice & LED** module allows you to define up to four scheduled **medication** reminders.


The reminders are identified by the labels: MI1 – MI4

The purpose of the **medication** reminders is to broadcast vocal announcements through the CP to remind the residents to take their medication.

The text of the message is: "**Attention!** It is time to take your medication. Please take them and press the **Reset** button"

To acknowledge the reminder, the resident must press the blinking CP **RESET** button, within the configured **Timeout** timeframe. If the resident does not press **RESET** before the **Timeout** timeframe ends, a **Missed Reminder** report is sent to the monitoring station, if **Report** is enabled.

For information about configuring the **Timeout** parameter, refer to 10.1 General Configuration Parameters on page 141.



Refer to 10.2 Reminder Schedule Types on page 144.

10.4. Appointment Reminders

The **Voice & LED** module allows you to define up to four scheduled **Appointment** reminders.

The reminders are identified by the labels: MI5 – MI8

The purpose of the **appointment** reminders is to broadcast vocal announcements through the CP to remind the residents to prepare for their appointments.


Each reminder is defined with an additional parameter: **Time to Announce**.

Enter the scheduled time of the appointment as the **Time to Announce**, which is added to the text of the vocal announcement.

The text of the message is: "**Attention!** Your appointment is at <*time to announce*>. Please press the **Reset** button."

To acknowledge the reminder, the resident must press the blinking CP **RESET** button, within the configured **Timeout** timeframe. If the resident does not press **RESET** before the **Timeout** timeframe ends, a **Missed Reminder** report is sent to the monitoring station, if **Report** is enabled.

For information about configuring the **Timeout** parameter, 10.1 General Configuration Parameters on page 141.



ID	Enabled	Schedule Type	Time to announce	Report
5	<input checked="" type="checkbox"/>	One time	13:00	<input checked="" type="checkbox"/> 15-01-2020 14:00
6	<input checked="" type="checkbox"/>	Weekly	09:30	<input checked="" type="checkbox"/> 10:00 Days <input type="checkbox"/> Sun <input checked="" type="checkbox"/> Mon <input checked="" type="checkbox"/> Tue <input checked="" type="checkbox"/> Wed <input checked="" type="checkbox"/> Thu <input checked="" type="checkbox"/> Fri <input type="checkbox"/> Sat
7	<input checked="" type="checkbox"/>	Weekly	12:30	<input checked="" type="checkbox"/> 14:00 Days <input type="checkbox"/> Sun <input checked="" type="checkbox"/> Mon <input type="checkbox"/> Tue <input checked="" type="checkbox"/> Wed <input type="checkbox"/> Thu <input type="checkbox"/> Fri <input type="checkbox"/> Sat
8	<input checked="" type="checkbox"/>	Monthly (Day of Week)	15:30	<input checked="" type="checkbox"/> Wed 16:30 Weeks <input type="checkbox"/> 1 <input checked="" type="checkbox"/> 2 <input type="checkbox"/> 3 <input checked="" type="checkbox"/> 4 <input type="checkbox"/> Last

Refer to 10.2 Reminder Schedule Types on page 144.

10.5. Ride Reminders

The **Voice & LED** module allows you to define up to four scheduled **ride** reminders.

The reminders are identified by the labels: MI9 – MI12

The purpose of the **ride** reminders is to broadcast vocal announcements through the CP to remind the residents to prepare for their rides.


Each reminder is defined with an additional parameter: **Time to Announce**.

Enter the scheduled time of the ride as the **Time to Announce**, which is added to the text of the vocal announcement.

The text of the announcement is: "**Attention!** Your ride will arrive at <*time to announce*>. Please press the **Reset** button."

To acknowledge the reminder, the resident must press the blinking CP **RESET** button, within the configured **Timeout** time frame. If the resident does not press **RESET** before the **Timeout** timeframe ends, a **Missed Reminder** report is sent to the monitoring station, if **Report** is enabled.

For information about configuring the **Timeout** parameter, 10.1 General Configuration Parameters on page 141.



ID	Enabled	Schedule Type	Time to announce	Report
9	<input checked="" type="checkbox"/>	One time	14:30	<input checked="" type="checkbox"/> 15-01-2020 14:45
10	<input checked="" type="checkbox"/>	Weekly	14:45	<input checked="" type="checkbox"/> 15:00 Days <input type="checkbox"/> Sun <input type="checkbox"/> Mon <input type="checkbox"/> Tue <input checked="" type="checkbox"/> Wed <input type="checkbox"/> Thu <input type="checkbox"/> Fri <input type="checkbox"/> Sat
11	<input checked="" type="checkbox"/>	Monthly (Day of Month)	09:30	<input checked="" type="checkbox"/> 27 10:00
12	<input checked="" type="checkbox"/>	Monthly (Day of Week)	10:15	<input checked="" type="checkbox"/> Wed 11:00 Weeks <input type="checkbox"/> 1 <input checked="" type="checkbox"/> 2 <input type="checkbox"/> 3 <input checked="" type="checkbox"/> 4 <input type="checkbox"/> Last

Refer to 10.2 Reminder Schedule Types on page 144.

10.6. Alarm Reminders

The **Voice & LED** module allows you to define up to four scheduled **alarm** reminders.


The reminders are identified by the labels: MI13 – MI16

The purpose of the **alarm** reminders is to use the CP as an alarm clock to remind the residents of scheduled times, such as when to wake up in the morning or when to wake up from their other rest periods.

The text of the announcement is: "**Attention!** Alarm reminder. Please press the **Reset** button."

To acknowledge the reminder, the resident must press the blinking CP **RESET** button, within the configured **Timeout** time frame. If the resident does not press **RESET** before the **Timeout** time frame ends, a **Missed Reminder Report** is sent to the monitoring station, if the **Report** feature is enabled.

For information about configuring the **Timeout** parameter, 10.1 General Configuration Parameters on page 141.



Refer to 10.2 Reminder Schedule Types on page 144.

10.7. Pendant Test Reminder

The **Voice & LED** module allows the operator to define one scheduled **Pendant Test** reminder.

The reminder is identified by the label: MI17

The purpose of the **Pendant Test** reminder is to broadcast vocal announcements through the CP to remind the resident to perform a test of their emergency pendant.

The pendant test supports the EP, the EPA, the EPP, and the EPA-BC.



NOTE: For ES6502HC, the **Pendant Test Duration** parameter is available for:

- C7000
- ES6502HC panel CPU version 2.03.2 through 2.04, or 2.14 and later


To configure pendant test parameters:

On the **General** tab under **Vocal Reminders**:

1. Ensure that **Report Missed Reminders** is marked.
2. Set the **Timeout** pendant test interval.
3. If you want the pendant test to repeat if the pendant test fails to complete, set **Snooze Time** for the wait interval between pendant tests.
4. If you do not want the pendant test to repeat if the pendant test fails to complete, set **Snooze Time** to "No Repetition".

On the **Pendant Test** tab:

1. Mark **Report** to allow sending a **Missed Reminder** report when the pendant test fails.
2. Set the **Pendant Test Duration** to the wait interval, from the schedule date and time to begin the **Pendant Test** interval. The wait interval can be either 1, 2, or 3 hours.



When the duration time finishes:

- The CP **RESET** button starts blinking
- The **Timeout** test interval begins.
- The pendant test reminder announcement is broadcast once

The text of the reminder is:

"Attention! It is time to test your pendant. Please press your pendant button."

To complete the pendant test, the resident must press their emergency pendant button.

If the resident does not press the pendant button before the test interval ends, a **Missed Reminder** report is sent to the monitoring station, if **Report** is enabled.


Refer to 10.2 Reminder Schedule Types on page 144.

11. Mobile Operators




NOTE: The **Mobile Operators** module parameters are available for panels with cellular capabilities:

- C7000
- ES6502HC panel CPU version 2.04 through 2.14



Click **Mobile Operators** on the **Main** navigator menu. The **Mobile Operators** module is the tool that allows management of the mobile operators APN information for Care@Home™ control panels installed on the resident’s premises.



11.1. APN Information

Access Point Name (APN) information provides access to the cellular communications channel. This information is provided by your SIM-card’s Cellular Operator or your Service Provider.

- APN Name
- User name
- Password

The following table provides instructions for entering the APN information for each Mobile Operator.

Field	Definition and Instructions	Attributes
ID	Automatic allocation of an internal identification code for each Mobile Operator.	N/A
Mobile Operator Code	Enter the Mobile Operator code.	Integer (6)
APN	Enter the Access Point Name (APN) name from your SIM card or cellular provider.	Mandatory. Up to 32 ASCII characters
User	Enter the APN user name.	Optional. Up to 28 ASCII characters
Password	Enter the APN password.	Optional. Up to 12 ASCII characters

11.2.Managing the Mobile Operator List

The **Mobile Operators** module allows you to manage the APN information for a list of mobile operators. The **Mobile Operators** module supports APN information for up to 10 mobile operators.

The module provides you with the capability to perform the following functions:


- Add a new mobile operator in a blank row in the Mobile Operators list.
- Update a row with an existing mobile operator.
- Delete a row with an existing mobile operator.

11.2.1 Adding a Mobile Operator

The **Update** function allows you to **Add** a new mobile operator by updating an empty row.

1. Select an empty row from 1-10 in the Mobile Operator list window.
2. Enter the destination information in the fields to the right of the Mobile Operator list window:
 - Mobile Operator Code
 - APN
 - User
 - Password

- Right-click on the row in the Mobile Operator list selected in Step 1. A function window appears displaying the functions available to manage the Mobile Operator list.
- Click **Update** in the function window or **Update** at the bottom of the screen, to **add** the new Mobile Operator to the Mobile Operator list. The selected row is populated with the information entered in Step 2.



11.2.2 Updating a Mobile Operator

The **Update** function allows you to revise an existing destination in the Mobile Operator list.


- Select a row in the Mobile Operator list window with an existing destination.
- Revise the APN information in the fields to the right of the Mobile Operator list window:
 - Mobile Operator Code
 - APN
 - User
 - Password
- Right-click on the row, in the Mobile Operator list window, selected in Step 1. A function window appears.
- Click **Update** in the function window or click **Update** at the bottom of the screen. The selected row is populated with the information entered in Step 2.

11.2.3 Deleting a Mobile Operator

The **Delete** function allows the user to remove an existing mobile operator from the Mobile Operator list.

- Right-click on the row, in the Mobile Operator list, that is to be deleted. A function window appears.

- Click **Delete** in the function window or **Delete** at the bottom of the screen. A window appears prompting you to confirm the delete request.




- Click **Yes**.

11.3. Load APN Information

APN information is found in a **.cpf** file.

You can upload the APN information from another **.cpf** file by using **Load APNs From CPF**. The APN information is loaded directly into the **Mobile Operators** module.



To load APN information from an external file:

- Enter the CMS **Mobile Operators** module.
- Click **Load APNs From CPF**. An **Open CMS File** browser window appears.
- Click the **.cpf** file that includes the APN information.
- Click **Open** to load the APN information from the **.cpf** file into the **Mobile Operators** module. The Mobile Operators list displays the updated APN information.




NOTE: Only APN information is uploaded from the selected external **.cpf** file into the CMS application.

5. To save the APN information:


- Download the CMS information to the CP.
- Save the CMS information to a CMS parameter file (**.cpf**).

Refer to 4.1 Main Action Icons on page 37.


12. The Custom Labels Module



Click **Custom Labels** on the Main navigator menu. The **Custom Labels** module is the tool that allows management of the personalized labels for Care@Home™ CMS installed on the resident’s premises.



12.1. Defining Custom Labels



The following table describes the parameters needed to define the personalized labels to the CMS application.

Name	Definition	Data Values	Factory Settings
Label ID	The internal identification of a label in the Custom Label list.	<ul style="list-style-type: none"> ● For ES6502HC, <ul style="list-style-type: none"> ■ ES6502HC PERS firmware – Panel CPU version 2.03.2 through 2.04, or 2.14 and later – 1-20 labels. ■ ES6502HC ADL firmware – panel CPU version 2.04 through 2.14 – 1-20 labels. ■ Otherwise, 1 - 50 labels ● For C7000, 1 - 20 labels 	
Encoding	The encoding character set with which to enter the personalized text for the custom label.	Dropdown list – one option	Latin-1
Label Text	Free form text up to 31 characters.	Aphanumeric	

12.2. Managing Custom Labels

To add a custom label:

1. Select an empty row in the **Custom Label** list.
2. Enter the name of the label in the **Label Text** field.
3. Click **Update**.

The update process adds the custom labels to the dropdown list for the following parameters in other modules:

- The **Zone Label** parameter in the **Activity Sensors** module
- The **Device Label** parameter in the **Safety and SOS** module


To update an existing label:

1. Select the row in the **Custom Label** list with the label you want to edit.
2. Modify the text in the **Label Text** field.
3. Click **Update**.


To delete an existing label:

1. Select the row in the **Custom Label** list with the label you want to delete.
2. Remove the text in the **Label Text** field.
3. Click **Update**.

13. The Log Events Module



Click **Log Events** on the **Main** navigator menu. The **Log Events** module is the tool that allows management of the event log files for the Care@Home™ CP installed on the resident’s premises.



Care@Home™ provides the ability to capture and read a CP event log file (.lgf file). Essence Professional Support personnel use the log file to analyze system behaviors and malfunctions.

13.1. Log Events Module Actions

The log files are stored in the CP. The **Log Events** module action buttons allow the user to view, store, print, or reinitialize the events log file in the CP.



13.2. Number of Events

The **Number of Events** parameter allows the user to enter the number of records to retrieve from the CP cyclical event log file. The records retrieved are the last events recorded by the log file.

The log file can store up to 945 event records. Since the log file is cyclical, these records are of the last 945 events recorded by the CP and its peripheral devices.


13.3. Read Log

Click **Read Log** to retrieve the number of records entered in the **Number of Events** parameter. The records are retrieved from the CP event log file in “Last In First Out” (LIFO) order. The records are displayed in the **Log Events** main page in “descending order” by the date and time the event was reported.

Note: Time gaps may appear in the log when the panel’s clock is updated.

13.4. Reset Log

Click **Reset Log** to empty the current CP event log file. You are prompted to confirm the reset log request.




Reset Log Confirm Prompt

The contents in the CP event log file are deleted. The event log file is ready to log new CP events, issued by the CP, the sensors, and other peripheral devices connected to the CP.

13.5. Save Log

Click **Save Log** to create an external copy of the retrieved records from the event log file in the CP. The records displayed on the **Log Events** screen are saved to an external file.

A Windows **Save Log** dialog box appears, to allow the user to specify the log file name with the extension **.lgf** and the location to store the log file on the local PC.



Save Log Window




NOTE: Use the **Save Log** feature to create the external **.lgf** file, when requested by the Essence Professional Service personnel to send them the CP events log file.


13.6. Open Log

Click **Open Log** to open the Windows **Open Log** dialog box, to locate an event log file (**.lgf**) on the local PC, “open” it, and display its contents on the **Log Events** page.

For example, the records in the event log file are retrieved and displayed on the **Log Events** page.



Open Log Window



Log File Opened on the Log Events Screen

13.7. Print Log

Click **Print Log** to open the **Print Log** window, to allow the user to choose whether to produce an Excel file, or a .csv file.

The log file is formatted in Windows EXCEL to properly view and understand the event log records.

For example, an event log file is formatted and printed as an Excel file.

	A	B	C	D	E	F	G	H
1	ACCOUNT: 00000259	CPU A: 2.03.230.05.12	CPU C: 246.00	DATE: 01-01-15	HOUR: 14:44:00			
2								
3	ID	Event Description	Date	Device Type	Device ID	Additional Information		
4	1	Provider time request	updating time	01-01-15 09:58			Day=1	Month=1 Hour=9
5	2	Status Panel		01-01-15 11:27	RKP	0	Duration=0	
6	3	Local Maintenance Mode		01-01-15 11:27	Control Panel	1	Partition FULL	RF device Normal A
7	4	Local Serial Interface		01-01-15 11:27			Start Maintenance Machine	Normal Activation
8	5	GPRS data modem		01-01-15 11:07			Failed to start connection	Message_group_number=0 RSSI=22
9	6	GPRS data modem		01-01-15 11:06			Failed to start connection	Message_group_number=0 RSSI=22
10	7	GPRS data modem		01-01-15 11:06			Failed to start connection	Message_group_number=0 RSSI=21
11	8	Send SIA Message System Automatic Test		01-01-15 11:05	Control Panel		Message_group_number=1	

Print Log Function Results

13.8. Read Log File

Click **Read Log File** to open the Windows **Open** dialog box to locate the remote log file (*.rlg*) on the debug server, “open” it, and display its contents on the **Log Events** page.

The *.rlg* file stores dialer debug data written to the remote log file. Refer to 5.4.8 Log on page 63.

14. Using CMS for Care@Home™ Implementation

Care@Home™ CMS provides the user interface to verify that Care@Home™ has been predefined according to the resident's needs. CMS has the tools to modify and redefine the configuration parameter values of the CP and its peripherals, as needed. These tools allow the user to define additional peripheral devices to Care@Home™.

To maintain Care@Home™ CP configuration parameters, first the user must upload the CP's default factory settings (configuration parameter values) to CMS. The user can use the tools provided by CMS to maintain the configuration parameters.

Follow the following example of using CMS to implement Care@Home™ :

A Care@Home™ kit is assembled according to a predefined default configuration.

For example:

- One CP
- Five motion detectors (PIRs)
- Two magnetic sensors (MGLSs)

The Care@Home™ CP is initialized with the default factory settings for the CP configuration parameter values. For example, the account number should be set for the account number assigned to the resident for whom the Care@Home™ is allocated.

The default peripheral devices are defined in the Care@Home™ CP. Optional peripheral devices can be added to the configuration to provide additional capabilities to Care@Home™ , such as

- A camera detector (IPD)
- Safety devices, such as a smoke alarm (SK2) or a water detector (MGLS)
- Emergency devices, such as a multidimensional fall detector (MDS) or a voice panic detector (VPD)

The Care@Home™ default configuration is flexible and allows for expansion. Care@Home™ can accept additional peripheral devices as needed.

14.1>Loading Factory Settings


The factory settings are defined in the Care@Home™ CP, provided by the manufacturer. The CP configuration parameter values must be uploaded to CMS for the user to verify that the factory settings are correct.

To upload the CP configuration parameter values:

1. Connect the CP to power and switch it on.

2. Connect a serial cable to the COM port on the PC and the CP mini-USB. Wait for the rising tone to confirm connection. Refer to 3.1.1 Local Connection to the CP on page 23.



3. Double click  .
4. On the **Connect** screen, select the COM port from **Com Port** list. Refer to 2.2 Identifying the Communications Port on page 10.



5. Click  .
6. Click  to upload the configuration parameter values from the CP to CMS.

14.2. Verifying Factory Settings

Following the upload of the CP configuration parameter values, CMS is prepared for the user to verify that the factory settings are defined correctly in the CP.

The user can view the factory settings of the CP in the CMS **Control Panel**. Refer to 5 The Control Panel Module on page 43.

Enter changes if needed and click **Update**. The following are configuration parameter values to verify.

14.2.1 Account Number

The CP account number should be set to the account number assigned to the resident for whom Care@Home™ is allocated.

Verify that value displayed in the **Account Number** parameter is the account number assigned to the resident. If not, modify the value in the parameter field accordingly.

Refer to 5.1 Account Number on page 43.


14.2.2 ARC Account Number

The ARC account number is an optional identifier. If using the ARC number, verify that it is correct. Refer to 5.2 ARC Account Number on page 44.

14.2.3 APN Settings

Control panels with cellular capabilities must have at least one mobile operator defined. This mobile operator should be the cellular service for the SIM card used by the CP. Refer to 6.1.5 Phone Number on page 82.

The cellular service for the SIM card inserted in the CP determines the APN information required. To access the APN settings, refer to 6.2.3 APN on page 85.



APN Information in the Communications Module

The purpose of the **Mobile Operators** module is to provide the APN information for the chosen cellular service. For example, the APN information in the **Mobile Operators** module is needed for the

cellular service used by the SIM card. The CP automatically retrieves the APN information associated with the cellular service used by the SIM card.

ID	Operator Code	APN Name	User	Password
1	42501	uinternet	orange	orange
2	42502	internet	Comfort	internet
3	42503	internet.pelephone.net.il	pcl@3g	pcl
4	42502	internet.golantelecom.net.il		
5	42502	net.hotm	hot502	2468
6	42507	net.hotm	hot507	1234
7	42508	internet.golantelecom.net.il		
8	42503	internet.ri	ri@3g	ri
9	42501	uwap.orange.co.il		
10	42515	uwap.orange.co.il	orange515	1234

Mobile Operators APN Information


Refer to 11 Mobile Operators on page 152.

14.2.4 Customizing the CP Configuration

The CP configuration covers all aspects of the CP functionality for itself and for its peripheral devices. The CP module allows the user to customize the capabilities of the CP functionality according to the user’s needs.

For example, the user may want the medical siren on the premises to sound for a specific duration before the CP contacts the emergency call center.

The **Sound** fields manage the medical siren on the premises. Modify the number of seconds for the siren to sound the alarm.




Sound	
Announcement Volume	1
Medical Siren Duration	15 (sec)
Safety Siren Duration	Default (sec)
Announcement While Connecting	<input type="checkbox"/>

The Medical Siren Duration

Another example of customizing the CP configuration is redefining the number of rings the CP waits before automatically answering a call from the emergency call center.

The **Remote Call-In** fields manage the receipt of external phone calls from the emergency call center.



Remote Call-In Rings

Choose the number of rings, from the dropdown list, for the CP to wait until automatically answering external calls to the CP.



NOTE: Remember to click **Update**, at the bottom of the **Control Panel** module page, upon completion of any modifications made in this module.

14.3. Configuring the Care@Home™ CP

The user can use CMS to configure the CP. Options include:

- Dialing
- Working with peripherals
- Event logging

14.3.1 Dialer Module


As part of the factory settings, the manufacturer assigns the ports and the default server’s IP address to each of the communication channels. This is to ensure that the CP communicates with the default server.

The CMS **Dialer** Module provides the user with the tools to maintain these settings for communications from the CP. The ports and IP addresses along with other dialer destination information can be viewed and verified in the **Dialer** module.

For example, shows the factory settings for the communication channels included in a Care@Home™ configuration. Each communication channel is assigned dialer destination information, such as:

- Destination ID

- IP Address for the Default Server
- Port
- Comm Type
- Comm Channel



Factory Settings in the Dialer Module

Refer to 7 The Dialer Module on page 92.

The user can maintain the destination information for each destination in the **Destinations** list using the tools provided in the **Dialer** module.



NOTE: Click **Update** on the bottom right of the **Destinations** list, to save the changes made to each line in the **Destinations** list.

The following are examples for maintaining dialer destination information.

14.3.1.1 Altering IP Addresses

Although the port assignments usually remain unchanged, the IP addresses associated with the ports can be altered. Altering the IP addresses is only necessary if the default server is replaced by another server and network.



NOTE: Do not change the port assignments without first consulting Essence technical support personnel.

14.3.1.2 Add Dialer Destinations


Additional dialer destinations can be added to the **Dialer** module.

For example, an additional server can be configured if the user wanted to provide a backup server in the event the default server fails to respond. The communication channels associated with an additional server must be added to the **Dialer Destinations** list.

In this case, the user can use the tools provided in the **Dialer** module to add new destinations to the **Dialer Destinations** list.

14.3.1.3 Edit Dialer Information

Emergency numbers or IP addresses can be assigned to **Analytics Message Comm Type** destinations

1-10. This is generally the destination the CP calls when  is pressed. When using SCAIP, the receiver may instruct the control panel to use a different number.

For example, the user should verify that a valid emergency call number is defined for the **Analytics Message Comm Type**.

If no emergency call number is defined for the **Analytics Message Comm Type** or the emergency call number defined is invalid, the user should redefine the emergency call number.

In this case, the user can use the tools provided in the **Dialer** module to edit existing destinations in the **Dialer Destinations** list.


Refer to 7.5 Managing the Destination List on page 115 and

A Fail-Safe number can also be assigned to ensure proper communication. Refer to 7.6 The Dialing Sequence on page 116

14.3.2 Redefining Dial Sequence Retries

The factory settings include defining the dialing sequences for each **Comm Type**. The **Dialer** module provides the tools to redefine the **Dial Sequence** parameters for each **Comm Type**. For example, the

Dial Sequence parameters include the number of retries when attempting to communicate with a **Comm Type**.



Redefine Number of Retries

Choose a number to redefine the number of retries allowed to communicate through a **Comm Type**. Click **Update** to save the modification.

Refer to 7.6 The Dialing Sequence on page 116.

14.3.3 Activity Sensors Module


As part of the factory settings, the manufacturer defines each of the peripheral devices in the CP configuration as activity sensors. Each peripheral device is defined according to the device Care@Home™ activity.

The CMS **Activity Sensors** module provides the user with the tools to maintain the definitions and characteristics of each of the CP’s activity related peripheral devices, like motion, door, and bed sensors.. The activity sensor definitions and characteristics information can be viewed and verified in the **Activity Sensors** module.

The factory settings for the peripheral devices are included in the CP configuration. Each peripheral device is defined as an activity sensor. The definitions and characteristics information defined for each peripheral device determines the device’s purpose or role in the Care@Home™ configuration.

The definitions and characteristics are displayed in the parameters, such as:

- Device ID
- Serial Number
- Device Type
- Activity Type



Activity Sensors Module

Refer to 8 The Activity Sensors Module on page 122.

The user can maintain the definitions and characteristics information for each activity sensor in the activity sensor list using the tools provided in the **Activity Sensors** module.

The following are examples for maintaining activity sensor information.

14.3.3.1 Adding Activity Sensors

The user can use the tools provided in the **Activity Sensors** module to add new activity sensors.

For example:

- The user may require additional peripheral devices in other rooms on the premises.
- An optional camera detector (IPD) is installed as part of the CP configuration.

All these peripheral devices must be defined as activity sensors in the **Activity Sensors** module

14.3.3.2 Altering the Activity Type

The default Pro CP configuration is comprised of two MGLS devices and five PIR devices. Each device is predefined for its activity type, which determines the purpose of the device and where the device is installed.



NOTE: Consult Essence technical support personnel before changing activity sensors' activity type definitions.

If the addition of devices to the CP configuration changes a device's role or purpose, the activity type for that device must be redefined.

If a camera detector(s) is added to the CP configuration, the characteristics of the PIR devices can be revised to enable camera coverage. When camera coverage is enabled, the IPD photographs events if a rule, like possible fall, is detected.

The user can choose the device identification number for the IPD that is to be triggered:

- First Priority IPD
- Second Priority IPD

Choosing **NONE** for no camera coverage for the device.


Refer to 8.1 Defining Sensor Device Attributes on page 123.

The user can use the tools provided in the **Activity Sensors** module to edit the existing activity sensors.

14.3.3.3 Resetting the Send Detections Capability

Characteristics of a device can be altered. Characteristics are defined by the parameter fields to the right of the **Activity Sensor** list.

For example, it is important to set the characteristic parameter that determines whether to send detections from a device to the server. Mark **Send Detections** and click **Update**, to enable sending detections from the device to the server.



Device Characteristic Parameters



NOTE: Click **Update**, at the bottom of the **Activity Sensors** module page, to save the changes made to each line in the **Activity Sensors** list.

Refer to 8.2 Managing Activity Sensors on page 125.

14.3.3.4 Configuring a Camera Detector

If a camera detector (IPD) is added to the Care@Home™ configuration, the user must not only define the IPD in **Activity Sensors** module but also configure the IPD to be triggered automatically by specific events or other triggers that occur in Care@Home™.

The tool for configuring an IPD is in the **Photo Configuration** tab of the **More Options** window of the **Control Panel** module.





Photo Configuration

To configure IPD camera coverage:


1. Mark **Alarm** to automatically trigger the IPD when the alarm is sounded.
2. Mark **Photo Request** to enable the control center to request a photograph when handling an alarm.
3. Mark **SOS** to automatically trigger the IPD when  is pressed.
4. If **SOS** is marked, for camera coverage priority:
 - Select the first camera detector for **Camera 1** or **NONE** for no camera coverage.
 - Select the second camera detector for **Camera 2** or **NONE** for no camera coverage.
5. Click **Update**.

14.3.4 Safety and SOS Module

If the user decides to add safety devices or emergency devices or both, as part of the resident’s CP configuration, the manufacturer will include the definition of these devices in the factory settings.

The devices are defined as safety or emergency devices in the CMS **Safety and SOS** module. The Safety and SOS definitions and characteristics information can be viewed and verified in the **Safety and SOS** module.

The safety and SOS devices are defined in the **Safety and SOS** module. To the right of the list of safety and SOS devices, the user can find the definitions and characteristics for each of the safety and SOS devices.



Safety and SOS Module

Refer to 9 The Safety & SOS Module on page 129.

The **Safety and SOS** module provides the tools for user to maintain the safety and emergency devices.

Refer to 9.3 Managing Safety & SOS Devices on page 135.

14.3.5 Log Events Module

The purpose of the **Log Events** module is to allow the user to view the log of events that occurred during Care@Home™ operations.

If an error occurs in the CP, the user can retrieve information from the events log file that may be able to explain the error and aid in recovery and solving the problem.

The user can:

- Download the log file from the CP (Read Log)
- Display the log file (Read Log).
- Save the log file to an external file in a readable format (Save Log)

Refer to 13 The Log Events Module on page 159.

15. Building a Kit-It Configuration

You can use the Kit-It tool to make configurations for the kitting process.

The following explains:

- Kitting
- Kit-It configurations
- How to create a predefined kit configuration using CMS in Kit-It mode

15.1. Introduction to Kitting

Kitting is when warehouse personnel assemble and customize Care@Home™ kits for your customers' service level requirements. The kits include a CP and its peripherals. The purpose of kitting is to update the CP with the updated system parameter data and a revised device inventory, sometimes with additional devices.

Kitting is needed for the following scenarios:

- A new CP – changes to software and the device inventory needed for a service level
- A CP needs a software upgrade – no changes to the device inventory
- An CP needs customization - changes to the device inventory needed
- A CP is inactive or needs a service level upgrade - changes to software and the device inventory needed for a service level



NOTE: Consult with your Essence account manager when creating a pre-defined configuration file for one or more active CPs.

The kitting process involves running updates for CP firmware and vocal messages, and uploading a predefined kit configuration, prepared by the account manager, to update the CP, followed by changing the CP's device inventory by adding, removing, and switching devices, by serial number and activity assignment. The process per CP ends by downloading the customized configuration to the CP and saving the customized configuration to a **.cpf** file.

Connect the CP to CMS on your PC and use the CMS to make changes to the CP software and device inventory. Refer to 14 Using CMS for Care@Home™ Implementation on page 164. To add **optional** devices to the CP's device inventory, use CMS in **Kit-It** mode.

Refer to 15.3 Creating a Predefined Kit Configuration Using CMS in Kit-It Mode on page 179.

Each CP update requires a predefined kit configuration. The kitting process for the various scenarios is generally the same. The difference is whether the CP's device inventory is being changed or fully redefined.

When the CP software update and the CP's device inventory changes are complete, download the customized configuration to the CP. Any optional devices added to the configuration become basic devices in the CP device inventory.

For information about warehouse kitting process, refer to ESUGSC030 Care@Home CMS Warehouse Guide.

15.2. Understanding a Kit-It Configuration

A predefined kit configuration includes the CP's system parameter updates, the basic device inventory required for a service level, and sometimes additional devices to allow for expanding the device inventory.

The following table lists the kit configuration content.

Name	Definition	Refer to:
CP parameters	System parameters, vocal reminders definitions, and mobile operator information that define the CP and peripheral functions of the residents' Care@Home™: Updates the CP configuration and downloads to the kit configuration file. Account number, ARC account number and control panel phone number are not overwritten.	5 The Control Panel Module on page 43 10 Voice & LED Module on page 139 11 Mobile Operators on page 152 6 The Communications Module on page 79
Dialers	Communication methods used to manage communication between a CP and its monitoring station: If the configuration includes GPRS SCAIP communications channels with callback enabled, a callback SIM number is required for the warehouse CP update. Note: Dialer functionality is dependent on your CP model	7.1 Destination Definitions by CP Type on page 92 7.6 The Dialing Sequence on page 116
Rules and Periods	Used to monitor a resident and analyze the resident's activity, specifically to identify and report behavior abnormalities:	5.4.4 Periods on page 52 5.4.5 Rules on page 54
Basic devices	Activity sensors and safety and SOS devices included in a CP's device inventory, as required by a service level: Updates the CP device inventory and downloads to the kit configuration file.	8.2.1 Adding an Activity Sensor on page 126 9.3.1 Adding a Safety or SOS Device on page 136
Optional devices	During kitting, if there are no optional devices in the predefined kit configuration, devices cannot be added to the kit configuration.	

The device categories are:

- Activity sensors
- Safety and SOS devices

Add camera detectors used to capture comfort videos as optional devices.



NOTE: The predefined kit configuration file (**.cpf**) settings must be compatible with the CP's type and firmware version and communication capabilities.

Predefined kit configuration files (**.cpf**) are located in the following default path:

[C:\Program Files \(x86\)\Essence\Care@Home CMS\Input](C:\Program Files (x86)\Essence\Care@Home CMS\Input)

There are several ways to create a predefined kit configuration:

- Using CMS, save your current or updated CP configuration to an external (**.cpf**) file.
Refer to 4.2 File Menu Functions on page 39.
- Using CMS in **Kit-It** mode, create a predefined kit configuration for a service level, including basic and optional devices.
Refer to 15.3 Creating a Predefined Kit Configuration Using CMS in Kit-It Mode on page 179.
- Using CMS warehouse **Kit-It**, use the **Kit-It** functions to create a new kit configuration by reinitializing and/or customizing a predefined kit configuration.

For information about CMS warehouse **Kit-It**, refer to ESUGSC030 Care@Home™ CMS Warehouse Guide.

15.3. Creating a Predefined Kit Configuration Using CMS in Kit-It Mode

You can create a predefined kit configuration using CMS in **Kit-It** mode. Contact Essence technical support for your **Kit-It** user credentials.

1. Before you begin, prepare the following:
 - Serial cable
 - The service level required
 - Optional device types required by the service provider for customization





NOTE: For service level kit configurations, use fictitious serial numbers for basic and optional devices.

2. To create a predefined kit configuration using CMS in **Kit-It** mode:

3. Connect the CP to power and switch it on.
4. Connect a serial cable to the COM port on the PC and the CP mini-USB. Wait for the rising tone to confirm connection.



5. Double click  .
6. On the **Connect** screen, enter your **Kit-It** user credentials.
7. From **Connect** drop down list, select the COM port used to connect to the serial cable.
8. Click **Connect**. When the connection is established, you are in **Kit-It** mode.



Indicates **Kit-It** mode

Working in Kit-It Mode

9. If the configuration is for an existing customer, enter the CP account number. Refer to 5.1 Account Number on page 43.
10. If you want to use the dialer callback feature, enable **Callback** on a GPRS SCAIP dialer and click **Update**. Depending on your CP model, refer to 7.1 Destination Definitions by CP Type on page 92.



NOTE: When running the CMS warehouse with **Configuration File** marked, you must enter the callback SIM number. When downloading the CMS warehouse kit configuration to the CP, the callback SIM number overrides the phone number. Refer to 6.1.5 Phone Number on page 82.

11. In the **Activity Sensors** module and the **Safety and SOS** module, add basic devices for the service level and optional devices if required by the service provider:
 - a. For Mode, click Basic or Optional.
 - b. Add device.
 - ◆ Refer to 8.2.1 Adding an Activity Sensor on page 126.

- ◆ Refer to 9.3.1 Adding a Safety or SOS Device on page 136.
 - c. For each optional device, enter a **Label**.
 - d. Click **Update**.
12. Save the updated CP configuration as a predefined kit configuration (**.cpf**) file, with a filename that includes the service level and the service provider identification. Both the basic and the optional devices are included in the configuration file.



16. Remote Boot Software Updates

Software and firmware updates for the CP and its peripherals are performed by the **Remote Boot Feature**. The **Remote Boot** feature allows the management of the CP, the voice announcements, and the CP software and peripheral firmware updates.

16.1. Launching the Remote Boot Feature

To launch the **Remote Boot** feature:




1. Double click . Care@Home™ CMS opens the home page and the **Connect** window appears.
2. Select the **Remote Boot** protocol.
3. Select the COM port according the connection method.
4. Click .
5. Following the successful completion of the connection procedure, the **Remote Boot** window appears.

16.2. Remote Boot Setup

The Remote Boot setup is required to prepare for the software and firmware updates.

The Remote Boot Setup includes:

- Setting the **Boot Options**
- Selecting the update (.esi) file



The Remote Boot Window

The following are the **Boot Options**:

- **Verify** – If marked, the CP automatically verifies that the software/firmware (.esi) file transfer completed successfully.
- **Automatic Firmware Update** – If marked, the CP automatically performs the software/firmware update. When not marked, use the **Reprogram** feature to perform the software/firmware update.
- **Target** – Allows you to specify the target device for which the **Remote Boot** feature is to update the firmware.

The following table lists the target devices and their usage.

Target Device	Upgrade
Main CPU	CP software
IPD CPU	The camera detector
SK2 CPU	The smoke detector

Target Device	Upgrade
EP CPU	Portable emergency button
EPP CPU	Portable emergency button
EPA CPU	Portable emergency button with fall detection
SPB CPU	Stationary panic button
Voice File	Recordings of voice announcements
VPD CPU	<ul style="list-style-type: none"> ● The voice panic detector host ● Voice processing unit
MDS CPU	<ul style="list-style-type: none"> ● Multidimensional fall detector host ● Multidimensional sensor


Following the definition of the **Boot Options**, use  to locate and select the update file.

16.3. The Software/Firmware Update Procedure

To update the target devices' software/firmware:


1. Define the **Boot Options**, located on the right-side of the **Remote Boot** window.

2. Click . The **Open ESI File** window appears.





Open ESI File Window

- Using the browser, select the update (.esi) file, received from the PSO team, for the target device chosen in the **Boot Option** setup. The file details appear in the **Details** box at the bottom of the **Remote Boot** window.



ESI File Details

- Click . The application verifies that the .esi file is the file for the target device. If the file is not for the target device, the download process aborts and an error message appears.



- If the file is for the target device, the update download process begins.


Information about the progress of the download process appears in the **View System Log** box and in the **Details** box. Note that there is a progress bar at the bottom of the window.



The information displayed includes:


- **File Name** - The name of update file
- **File Size** – The amount of data to be transferred
- **Sending** – The amount of data transferred
- **Time Left** – The estimated time remaining
- **Speed (bps)** – The data transfer rate

If you need to abort the update process, click **Stop**. The update process aborts and an error message appears.




To restart the update process:







1. Click .
2. Disconnect from the CP.
3. Re-connect to the CP.
4. Launch the **Remote Boot** feature
5. Click .
6. When the update process ends, **Success** confirmation appears along with the information displayed in both the **View System Log** box and in the **Details** box.



7. The CMS automatically disconnects from the CP and issues a message.



The following table summarizes the functionality of the **Remote Boot** action buttons and boxes.

Name	Definition and Instructions	Required or Optional	Default Value
 Select File	Select the file for the target device. A browser dialog window appears to allow location and selection ('open') of the relevant update file (.esi file).	Optional	Enabled
 Start Download	Click to perform the update file transfer and update the software/firmware, according to the Boot Options .	Optional	Disabled
 Reprogram	Click to manually initiate the software/firmware update.	Optional	Disabled
 Verify Firmware	Click to verify that the existing software for the Main CPU updated successfully. NOTE: This button is enabled only when the Main CPU is selected.	Optional	Disabled
 Stop	Click to stop the Remote Boot process.	Optional	Disabled
 Exit	Click to exit the Remote Boot window.	Optional	Enabled
View System Log	Displays the system messages and the error messages issued during the Remote Boot processing.	N/A	Display Only
Details	When a file is selected, the details of the file selected appear in this field.	N/A	Display Only

17. Using the Remote Multiple Device Manager Tool

The Care@Home™ CMS Remote Multiple Device Manager (RMDM) utility can for update software and firmware, and account numbers for deployed PSTN CPs:

- Firmware for any of the following target devices:
 - CP CPU
 - VPD
 - EPP
 - EPA
 - Smoke Detector (SK2)
- Voice files for CP
- CP parameters from a predefined set of values you configured in the Care@Home™ CMS application.

The Care@Home™ CMS RMDM utility uses a PSTN landline modem to connect remotely to the CPs in the panel list. Essence recommends that you use the USRobotics Dial-up External Modem (product code: USR5637). Refer to 3.2.3.1 Installing the Landline Modem on page 28.




NOTE: This reference manual refers to version 1.2.0 of the Care@Home™ CMS RMDM utility.

17.1. Prerequisites

Before you begin the update of the CPs, do the following:

1. Establish the remote PSTN connection using the PSTN modem.
2. Refer to 3.2.3 Establishing Remote Connection with PSTN Modem on page 28.
3. Create a **Panel List** file (**.csv**), using a text editor, with the phone numbers and account numbers for the CPs to be updated.
4. The file structure includes a single column of landline phone numbers. Each phone number must be followed by four commas. Make a second column for account numbers.
5. Enter the phone numbers in the format you would use for dialing the CP.

6. For example, CP phone numbers are in the US format with a leading “1” for calling phones outside your area code.



7. Locate one of the following files as the input file:
 - Firmware file (.esi) with which to upgrade the CPs or peripherals
 - Configuration file (.cpf) with which to update the CPs configuration parameters
8. Identify the target device type to be upgraded.
9. Determine the time periods for running the upgrade process.



NOTE: Since the process calls the resident landline phone number, it is important to define time periods that will not disturb the resident.


10. Define the optimum number of retries and time to wait between retry attempts.

17.2. Upgrading the CPs

The RMDM utility is designed for prolonged run times. If a CP is busy or offline, the RMDM allows you to rerun the update process just for the CPs that were not updated.


To update the CPs, the peripherals, or the CP configuration parameters:

1. Ensure that you have performed the prerequisite tasks. Refer to.




2. Double-click . The main window of the Care@Home™ RMDM utility appears.


- For the **Panel List** field, click . The **Open CSV File** window appears.




- Select the **Panel List** file (.csv). The path and name of the selected file appear in the **Panel List** field. The **Panel List** file details appear in the box to the right of the **Panel List** field.




- For the **Input File**, click . The **Open input file** window appears.




6. If you want to upgrade the CP or a peripheral, select a firmware upgrade file (.esi) for the target device.
7. If you want to update the CP configuration parameters, select a configuration file (.cpf).
8. If you want to upgrade the CP or a peripheral, select the device type to upgrade from the **Target** dropdown list.




9. If you want to update the CP configuration parameters, select **Configuration** from the **Target** dropdown list.
10. Mark the **Upgrade** period(s) that you defined when preparing the prerequisites. Refer to 17.1 Prerequisites on page 189.
 - If no **Period** checkboxes are marked, you can run the upgrade process nonstop until all the target devices associated with the CPs in the **Panel List** are upgraded.
 - If the **Period** checkboxes are marked such that one interval overlaps the other, the following error message appears:



- If the end time of a period is earlier than the start time of that same period, the following error message appears:




11. Enter the number of **Retries** and **Wait Time between Retries** as you defined when preparing the prerequisites. Refer to 17.1 Prerequisites on page 189.
12. Ensure that the values for all the parameters described in Steps 1-10 are entered and completed correctly.
13. Click **Execute**. If you entered all the parameters correctly, the utility begins the upgrade process.




If you did not enter the file parameters correctly, error messages are displayed:


- If you did not enter the **Panel List** file, the following error message appears:




- If the utility cannot locate the **Panel List** file you entered, the following error message appears:



- If you did not enter the **Input file**, the following error message appears:




- If the utility cannot locate the selected **Input file**, the following error message appears:







17.3. Managing RMDM Processing

The RMDM utility provides you with the tools to manage and control the upgrade process.

The following table provides the operating instructions for using the **Session Control Actions** to manage the upgrade process.

Session Control Actions	Instructions	Confirmation Windows
Execute	Click to initiate RMDM upgrade processing.	
Pause	Click to suspend processing. The confirmation window appears, prompting you to confirm the Pause action. Click Yes to terminate the current CMS task.	 <p>A confirmation dialog box titled "Remote Multiple Device Manager" with the text: "Are you sure you want to pause the session? This will terminate the current CMS task." It has "Yes" and "No" buttons.</p>

Session Control Actions	Instructions	Confirmation Windows
Exit	<p>Click to stop processing and close the utility main window. The confirmation window appears, prompting you to confirm the Exit action.</p> <p>Click Yes to terminate the current session.</p>	
New Session	<p>Following a Pause action, if you wish to dismiss the previous upgrade processing and restart the upgrade process of the Panel List file from the beginning.</p> <p>The confirmation window appears, prompting you to confirm the New Session action.</p> <ul style="list-style-type: none"> ● Click Yes to start a new session. ● Click No to retain the previous processing. 	
Continue Last Session	<p>Following a Pause action, the message appears in Red, prompting you to decide which action to choose.</p> <p>Mark option if you wish to resume processing of the last session.</p>	
Upgrade Failed Panels	<p>This option is enabled following the completion of the upgrade process for the Panel List file.</p> <p>If the results of the upgrade process include upgrade failure to one or more CPs, the option to restart the upgrade process is enabled.</p> <p>Mark if you wish to run the upgrade process again for the CPs that failed the upgrade process.</p>	

17.4. Upgrade Progress Reporting

The RMDM main window reports the progress of the process and the upgrade results for each CP as it is upgraded.

There are four areas of progress reporting:

- [Panel List Details box](#)
- [Current CMS task](#)
- [Session Statistics](#)
- [Export feature](#)

17.4.1 Panel List Details Box


The contents of the **Panel List** file are displayed in the box to the right of the **Panel List** field. The utility updates the status of the upgrade of each CP in the **Panel List** file.

The following table describes the information reported in the **Panel List** details box.

Column Names	Description	Notes
Phone	The landline phone number of the CP to be upgraded with the chosen firmware	
Status	The status of the CP Upgrade processing Caution: If residents answer the ringing of their landlines when the utility connects with the CPs, the upgrade process for each of the CPs fails.	<ul style="list-style-type: none"> ● Successful ● Fail ● In Progress ● Pending
Description	Additional information about the Upgrade process such as: <ul style="list-style-type: none"> ● The number of the retry currently running ● The reason for a failed Upgrade process 	
Target	Upgrade target: CPU or voice file	
File Name	File used for upgrade	
Account	Updated account number	

17.4.2 Current CMS task

The **Current CMS Task** information provides you with the real-time status of the **Upgrade** processing of the current CP.



For example, the **Current CMS Task** information shows:


- The phone number currently being processed
- The statistics about the current processing of the firmware file
- The version the CPU of the current CP
- Progress of the current upgrade process

The following table describes the information reported for the **Current CMS Task**.

Task Information	Description
Panel phone number	The phone number of current CP in process
Remaining time	The time remaining to transfer the firmware file
Sent	The current number of bytes transferred
File size	The total number of bytes to transfer
CMS status	Additional information about the CP, such as the version of the CP CPU software
CMS	The progress bar that represents the status of the current Upgrade process.

17.4.3 Session Statistics

The **Session Statistics** reports the overall progress of the **Upgrade** process within the session.



For example, the **Session Statistics** show:

- How many CPs have been processed?
- How many CPs are left to process?
- The overall progress of the upgrade process

The following table describes the information reported for the **Session Statistics**.

Statistics	Description
Success	The number of successful CP upgrades
Fail	The number of failed CP upgrades
Pending	The number of CP upgrades remaining
Current	The position in of the current CP in the Panel List file
CMS	The progress bar that represents the overall status of the Panel List Upgrade process.

17.4.4 Export feature

The **Export** feature allows you to generate a report of the CPs **Upgrade** process in the **Panel List**.

To generate the list, click **Export**. The utility generates a **.csv** file which includes the information concerning the current **Panel List** file upgrade process.

The following shows an example of an exported CSV file of the upgrade process report.

	A	B	C	D	E	F
1	Phone	Status	Description	Target	Filename	Account
2	15168682211	Success		Main CPU	c:\upgrades\PERS_2_3_1_230_1_31.esi	2256358
3	17186652214	Fail	Busy			2264897
4	12126695525	InProgress				2236559
5						

The following table describes the information reported in the upgrade process report.

Column Name	Description
Phone	The landline phone number of the CP processed
Status	The status of the Upgrade processing of the CP
Description	Additional information about the Upgrade process such as: <ul style="list-style-type: none"> ● The number of the retry currently running ● The reason for the failed Upgrade process
Target	The CP or peripheral device for which the upgrade process is running
ESI Path	The path and name of the firmware file (.esi)

Appendix A USB to RS232 Adapters

If the PC on which Care@Home™ CMS is installed has no available RS232 COM ports, the user is instructed to use a RS232-USB adapter to connect the CP to the PC.

Essence technical support recommends the following two alternatives:

- USB Gear Serial Adapter
- VScom USB-COM Mini

These are the descriptions and special features of the adapters.

The USB Gear Serial Adapter



The following is a list of some of the features of this USB to serial adapter:

- Over 500 kbps of data transfer capability
- Works with Modems, PDA's, Cell phones, Digital Cameras and more serial legacy devices
- LED visuals for transmit and receive between devices
- Compatible with: Windows 98, 2000, 2003, ME, XP, Vista, Windows 7 (32/64Bit) Microsoft Certified Drivers Provided/ MAC 8.6, 9.0, 9.1/Linux
- The USB Gear's 12" USB Serial Converter, with its powerful Hi-Speed DB9 Serial adapter, comes with software drivers and a manual. The customer can download the software drivers.

The VScom USB-COM Mini



The VScom USB-COM Mini adapter provides a high-speed RS232 serial port with 128 byte FIFO to allow for high speed communications, even for heavy loaded systems. The RS232 serial port on this adapter also enables high speed settings such as 500,000 bps.

The following is a list of the key features of this USB to serial adapter:

- 1 x RS232 DB9 male port

- Speed up to 921.6 kps

- No external power supply

The customer can access VS Com and more information about this adapter by following this link:
<http://www.vscom.de/619.htm>

Appendix B Voice Announcements

The following table includes examples of the texts for the vocal announcements listed by group.

Announcement Group	Announcement Text
Vocal Reminders	"Attention, it is time to take your medication. Please take them and press the Reset button"
	"Attention, your appointment is at [hour, minute AM/PM]. (Please press the Reset button)"
	"Attention, your ride will arrive at [hour, minute AM/PM]. Please press the Reset button)"
	"Attention, alarm reminder. (Please press the Reset button)"
	"Attention, it is time to test your pendant. Please press pendant button"
Medical	"Emergency alarm"
	"Activity alarm"
	"Emergency alarm – fall detected"
	"Call in progress"
Safety	"Flood alarm"
	"Fire alarm"
	"Temperature alarm"
	"Gas alarm"
Functional	"Cellular signal strength is x percent"
	" Control panel on"
	"Control panel shutting down"
	"Device added successfully"
	"Device already paired with control panel"
	"Memory full, additional device cannot be added"
	"Could not add device, please try again"
	"Pendant recognized"
	"Activity timer on"
	"Activity timer off"
	"Activity timer will expire soon. Please press the reset button."
	"Signaling mode on"
	"Signaling mode off"
	"To erase all devices press the pairing button, to exit press the RESET button"
	"Successful operation - All devices were erased from control panel"
	"Device paired with Control Panel"
	"Attention, a device is tampered – control panel"
	"Control panel entered test mode"
	"Control panel exited test mode"
	"Test successful"

Announcement Group	Announcement Text	
	"Bed sensor added successfully" – ES6502HC ADL only	
Technical	"Attention, communication with monitoring station is lost"	
	"Battery missing in control panel"	
	"Attention, low battery in console"	
	"Attention, a device is tampered"	
	"Attention, phone line is disrupted"	
	"Could not add device, please try again"	
	"Attention, control panel switched to backup battery"	
	"Attention, low battery in a device was detected"	
	"Attention, control panel backup battery needs to be replaced"	
	"Attention, supervision loss of a device was detected"	
	"Attention, RF jamming Was detected"	
	"RF jamming restored"	
	"Control panel performed software reset"	
	"Malfunction alarm"	
	"Installation test starts"	
	"Installation test ends"	
	"Attention, supervision loss of"	"motion detector"
	"input magnet detector "	
	"magnet sensor"	
	"fixed emergency button"	
	"portable emergency button"	
	"smoke Detector"	

Appendix C Alerts by Type

The specific alerts included in each alert type are found in the following table:

Alert Type	Alert Description
External	Sensor registered an extremely hot temperature.
External	High temperature alarm cancelled by CP RESET button
External	No movement detected for longer than expected, or no button pressed during the expected time frame
External	Sensor registered an extremely cold temperature
External	Cold temperature alarm cancelled by CP RESET button
External	Fire or smoke detected, and alarm triggered
External	Fire or smoke alarm cancelled by CP RESET button
External	Fire or smoke sensors restored, and alarms reset
External	Temperature sensor restored, and extreme heat alarm reset
External	Temperature sensor restored, and extreme cold alarm reset
External	Activity detected after a previous no activity alarm
External	Missed reminder: the CP RESET button wasn't pressed after a vocal reminder
External	Flood detected, and alarm triggered
External	Flood sensors restored, and alarms reset
External	Gas leakage detected, and alarm triggered
External	Gas sensors restored, and alarms reset
External	The resident remains in a location for longer than expected, according to their daily pattern of activity
External	The resident is away from the premises for longer than expected
External	The resident returned after being away for longer than expected
External	The entrance door is opened at an unexpected time, suggesting, for example, that the resident is wandering outside
External	A door is left open for longer than expected
Supervisory	The CP sends periodic communication test
Technical	CP main power restored.
Technical	CP lost main power supply; switched to backup battery
Technical	Medical emergency alarm device communication with CP restored
Technical	Medical emergency alarm device lost communication with CP
Technical	CP powered up/CP battery power is restored after battery was empty or missing
Technical	Device was tampered with
Technical	Device restored, and tamper reset
Technical	Batteries replaced in peripheral device
Technical	Low or empty battery in peripheral device
Technical	Empty or missing battery in CP
Technical	CP battery power recharged
Technical	CP communication with receiver lost.
Technical	Low battery in CP
Technical	Communication restored with fire or smoke sensors
Technical	Communication lost with fire or smoke sensors

Alert Type	Alert Description
Technical	Communication lost with temperature sensor
Technical	Local programming started
Technical	Local programming ended
Technical	Remote programming started
Technical	Remote reset
Technical	Remote programming ended successfully
Technical	Remote programming failed
Technical	A problem was detected with the control panel's power supply
Technical	CP main power supply restored
Technical	CP time was manually changed
Technical	Communication lost with flood sensors
Technical	Communication restored with flood sensors
Technical	Communication lost with gas sensors
Technical	Communication restored with gas sensors
Technical	Radio interference
Technical	Radio interference restored
Technical	CP powered up
Technical	PSTN link down
Technical	PSTN link restored
Technical	Cellular link down
Technical	Cellular link restored
Technical	Ethernet link down
Technical	Ethernet link restored
Technical	Communication lost with pendant
Technical	Communication restored with pendant
Technical	Low or empty battery in pendant
User Alert	Emergency call ended
User Alert	Emergency call started
User Alert	Either emergency alarm or fall detection alert triggered
User Alert	Medical alarm situation cleared by ARC
User Alert	Medical alarm cancelled by CP RESET button
User Alert	Emergency alarm from pendant
User Alert	Fall detected
User Alert	Emergency alarm triggered by voice
User Alert	Cord pulled (emergency alarm)
User Alert	Emergency alarm from fixed trigger

Legal Notice

Usage of this document, and all information (including product information) provided within, are subject to the following terms and conditions, and all applicable laws. If you do not agree with these terms, please do not access or use the remainder of this document.

This document contains highly confidential information, which is proprietary to Essence SmartCare Ltd. and/or its affiliates (hereafter, "Essence"). No part of this document's contents may be used, copied, disclosed or conveyed to any third party in any manner whatsoever without prior written permission from Essence.

The information included in this document is intended for your knowledge and for negotiation purposes only. Essence makes no implicit representations or warranties with respect to such information.

The information included in this document is subject to change without notice. Any decision to rely on the information contained herein shall be at your sole responsibility, and Essence will not accept any liability for your decision to use any information or for any damages resulting therefrom.

Certain laws do not allow limitations on implied warranties or the exclusion or limitation of certain damages. If these laws apply to you, some or all of the above disclaimers, exclusions, or limitations may not apply to you.

By using the information contained herein, you agree that the laws of the State of Israel, without regard to principles of conflict of laws, will govern any dispute of any sort that might arise between you and Essence regarding the information contained herein, and any such dispute shall be settled exclusively in the competent courts of Tel Aviv-Jaffa, Israel.

All registered or unregistered trademarks, product names, logos and other service marks mentioned within this document are the property of Essence, or their respective owners. Nothing contained herein shall be construed as conferring by implication, estoppels, or otherwise any license or right, either express or implied, under any patent or trademark of Essence or any third party. No use of any trademark may be made without the prior written authorization of Essence.

This document and all of its contents are protected intellectual property of Essence. Any copying, reprinting, reuse, reproduction, adaptation, distribution or translation without the prior written permission of Essence is prohibited.

Please check your End User License Agreement (EULA) for terms and conditions.

© 2021 All rights reserved to Essence SmartCare Ltd.

For more information, please contact:

Essence SmartCare Ltd.
12 Abba Eban Avenue,
Ackerstein Towers Bldg. D
Herzliya Pituach, 4612001 Israel
www.essence-grp.com
Tel: +972-73-2447777
Fax: +972-9-7729962