



Pillar® ONCO/Reveal™ Multi-Cancer with CNV & RNA Fusion Panel

The ONCO/Reveal Multi-Cancer with CNV and RNA Fusion Panel is a combined DNA/RNA multi-cancer panel. The assay combines the DNA-based ONCO/Reveal Multi-Cancer v4 with CNV Panel with the ONCO/Reveal Multi-Cancer RNA Fusion v2 Panel allowing for joint sequencing. The assay uses proprietary Stem-Loop Inhibition-Mediated amplification (SLIMamp®) technology, a tiled amplicon-based library prep chemistry for efficient single-tube target enrichment.

Variants and CNVs detected from DNA

ABL1	CDKN2A	FBXW7	GNAS	KIT●	NPM1	PTPN11	SRC
AKT1	CSF1R	FGFR1●	HNF1A	KRAS●	NRAS	RAC1	STK11
ALK	CTNNB1	FGFR2●	HRAS	MAP2K1	NTRK1	RB1	TP53
APC	DDR2	FGFR3●	IDH1	MET▲	NTRK2	RET	VHL
ATM	EGFR▲	FLT3●	IDH2	MLH1	NTRK3	ROS1	
BRAF	ERBB2▲	FOXL2	JAK2	MPL	PDGFRA●	SMAD4	
CCNE1●	ERBB4	GNA11	JAK3	MYC▲	PIK3CA●	SMARCB1	
CDH1	EZH2	GNAQ	KDR●	NOTCH1	PTEN	SMO	

CNVs detected and verified by NIST reference standard are indicated by ▲
CNVs can also be detected in genes indicated by ●

Fusions and expression insights detected from RNA

Driver gene fusions (fusion partners not listed)

ALK	EGFR	FGFR2	MET	NTRK1	NTRK3	PPARG	RAF1	ROS1
BRAF	ERG	FGFR3	NRG1	NTRK2	PBX1	PRKACA	RET	TFE3

Expression imbalance

ALK	FGFR3	NRG1	NTRK1	NTRK2	NTRK3	PBX1	RET	ROS1
-----	-------	------	-------	-------	-------	------	-----	------

Expression control

HMBS	TBP
------	-----

Simple NGS library prep workflow

Maintain control of samples and results with single-tube, tiled amplification that can be performed in-house by any NGS lab

Sensitive and robust chemistry

Achieve variant detection as low as 1% VAF[†] without UIDs[‡] even with limited DNA input or poor sample quality


Reduced fully-loaded lab costs

Improve lab efficiency and reduce “no calls”, repeat testing, and difficult interpretation decisions

*Content is based on data from ongoing clinical trials along with sources that include the College of American Pathologists (CAP), the Association for Molecular Pathology (AMP), the National Comprehensive Cancer Network (NCCN), and the Catalog of Somatic Mutations in Cancer (COSMIC) database.
[†]VAF, variant allele frequency; [‡]UID, unique ID; also known as unique molecular ID (UMI)

For Research Use Only. Not for use in diagnostic procedures.

Simple, one-day workflow



Panel specifications

	ONCO/Reveal Multi-Cancer v4 with CNV	ONCO/Reveal Multi-Cancer RNA Fusion v2
Enrichment chemistry	Multiplex PCR using tiled amplicons	
Number of pools	2 pools	
Number of genes/amplicons	60/341	18/>80 plus MET 14 exon skipping
Number of targets	Hotspots in 60 genes; CNVs for 14 genes	Fusions in 18 driver genes; expression for 11 genes
Variant types	SNVs, small and medium indels, and CNVs	Fusion RNA transcripts
Average amplicon size	125bp (range 86bp – 185bp)	120bp
Recommended input range	5ng to 80ng DNA	10ng to 50ng RNA
Sample types	DNA from tissue or blood; FFPE	RNA from FFPE or tissue
Mapping rate	99.3% ± 0.3%	n/a
% on-target aligned reads	99.5% ± 0.1 %	n/a
Coverage uniformity (% targets with >0.2X mean coverage)	98.2% ± 0.7%	n/a
Total assay time (from DNA to sequencer)	<8 hours	<9-10 hours
Sequencing platforms	Illumina®, MGISEQ™ and Ion Torrent™ platforms	

Order or learn more at
pillar-biosciences.com



For Research Use Only. Not for use in diagnostic procedures.
 ©2019 Pillar Biosciences. An ISO 13485:2016 certified company. Pillar®, SLIMamp®, PiVAT® and ONCO/Reveal™ are trademarks of Pillar Biosciences, Inc. Illumina® is a trademark of Illumina, Inc. Ion Torrent™ is a trademark of Thermo Fisher Scientific. MGISEQ™ is a trademark of MGI Tech Co, Ltd.
 MK-0024-2

Pillar Biosciences, Inc.
 12 Michigan Dr
 Natick, MA 01760
 (800) 514-9307
info@pillar-biosciences.com