


Polymer studies by QCM-D


Polymers are versatile building blocks within applications ranging from plastics used in everyday life to sophisticated biopolymers for medical devices. Quartz Crystal Microbalance with Dissipation monitoring, QCM-D, has been used for the characterization of a wide variety of polymers.

Introduction

Control of polymer build-up, conformation and degradation are key features in the study and development of polymer based materials. QCM-D provides this control through measuring changes in both mass and softness (related to energy dissipation) of polymer layers in real-time.



[Figure 1]: Illustration of polymers adsorbing in a brush (top) or pancake (bottom) conformation.



[Figure 2]: Thickness of multilayer after each polyelectrolyte deposition, measured with QCM-D.


Grafted polymer conformation and phase transition [1]

Polymers grafted to a surface can adopt several different conformations. The ability to control these transitions increases the applicability of the material. In this study QCM-D was chosen specifically for the sensitivity of the technique to study changes in polymer conformation. It was possible to control the conformation of the polymers in situ and in real-time by altering temperature and grafting density and also to identify different phases in the polymer adsorption process.

Build-up and control of polyelectrolyte multilayers [2]

Polyelectrolyte multilayers are used to modify surfaces in a vast number of applications thanks to their flexibility and robustness. Even though they are simple to build, the mechanisms of how different conditions affect the build-up is not well understood. QCM-D is an ideal tool to get better understanding of this as it can study the build-up of hydrated films in real-time. In this study polyelectrolyte multilayers were evaluated in terms of thickness per layer and the influence of salt concentration on build-up rate and total thickness.

Biopolymer grafting and cross-linking through specific interaction [3]



[Figure 3]: Illustration showing the frequency and dissipation shifts as hyaluronan (blue in illustration) is immobilized and later cross-linked by tumor necrosis factor-stimulated gene-6 (TSG-6) (green ovals).

References:

- [1] Zhang, G., Wu, C., Quartz crystal microbalance studies on conformational change of polymer chains at interface. *Macromolecular Rapid Communications*, 30 (4-5), 328-335, 2009
- [2] Zan, X., Peng, B., Hoagland, D.A., Su, Z., Polyelectrolyte uptake by PEMs: Impact of salt concentration. *Polym. Chem.*, 2, 2581-2589, 2011
- [3] Baranova, N. S., Nilebäck, E., Haller, F. M., Briggs, D. C., Svedhem, S., Day, A. J., Richter, R. P., The inflammation-associated protein TSG-6 cross-links hyaluronan via hyaluronan-induced TSG-6 oligomers. *J. Biol. Chem.*, 286 (29), 25675-25686, 2011