




*The most significant
advancements in TMD therapy in
over 35 years.*

It's Time for

MAPA[®]

In Your Practice

Dr. Danny K. Crout, DMD, MS has partnered with Pittman Dental Laboratory exclusively to produce the MAPA[®], or Maxillary Anterior Passive Appliance. The MAPA[®] is a multi-use appliance with multiple functionality for patients. It can be used as a deprogrammer, TMD/TMJ treatment device, night guard, and a device to protect newly placed anterior restorations.



The MAPA[®] is an every-day appliance, meaning that the patient doesn't wear it only at night. This allows the patient to become familiar with wearing the MAPA[®] as an every-day appliance, not a night only appliance which may hinder sleep comfort. Most patients become used to wearing the MAPA[®] within 10 days. The material is clear, making the MAPA[®] almost impossible to see.

For more information about incorporating the MAPA[®] into your practice, please visit www.MAPAbiteguard.com or www.pittmandental.com

**Act now and get 50% off your first MAPA[®]
with Pittman Dental Laboratory
(code:MAPA50)**



Impression Instructions and Techniques

Place tip of the tongue to the soft palate.
Practice biting on the second molars
(pretend if not present). Goal is to get
condyles in the TMJ fossae to rotate in their
correct orthopedic position.



With tip of tongue on soft palate, bite on
cotton rolls, 4 folded gauzes (2x2's) or most
commonly on 2 folded gauzes placed on
molars (keep off 2nd bicuspids). Goal is to
record bite registration just inside freeway
space anteriorly from second bicuspid to
second bicuspid.



Using a stable bite registration material,
record from 2nd bicuspid around to
opposite 2nd bicuspid.



For more information regarding the MAPA®, visit
www.MAPAbiteguard.com. Send your impressions and/or models
with bite registration to us and we will take care of the rest!