Rapid and sensitive detection of EGFR C797S mutations using a blood-based droplet digital PCR assay



Background

The adoption of third-generation tyrosine kinase inhibitors (TKIs) to treat non-small cell lung cancer (NSCLC) for EGFR positive patient cases has led to the emergence of acquired resistance pathways. Approximately 40% of EGFR^{T790M}-positive NSCLC cases display an acquired mutation of amino acid 797 (C797S) that mediates this resistance, which can develop after a median of 10 months following treatment. The effectiveness of treatment for patients with EGFR^{C797S} may depend on whether the mutation is present in *cis* or *trans* with EGFR^{T790M}, the type of sensitizing mutation (EGFR^{del19} or EGFR^{L858R}), and previous treatment(s). Patients with EGFR^{del19} or EGFR^{L858R} in conjunction with EGFR^{T790M} and EGFR^{C797S} remain resistant to osimertinib. We have developed a blood-based test that can detect two of the most common C797S nucleotide mutations (T>A and G>C) using droplet digitalTM PCR (ddPCRTM) technology and assays. Analytic sensitivity and specificity were assessed using synthetic DNA designed to mimic the EGFR^{C797S} variants and their locus detected in the assay. Normal healthy donor samples as well as reference positive and negative NSCLC donor samples were assessed for clinical specificity and sensitivity. Finally, the precision of the assay was evaluated with both clinical and analytical samples. Specifically, we evaluated the assay at high, medium, and low mutation frequencies over three consecutive days, including repeat runs on one day, and with multiple operators. This assay is capable of accurately and precisely detecting multiple EGFR^{C797S} variants from blood specimens in the clinical laboratory. Consistent with other ddPCR blood-based EGFR variant assays we have developed, the limit of detection was 0.02% for C797S variants. The EGFR C797S assay described here is approved by NYS-CLEP and may be of utility as part of the molecular diagnoses for physicians treating patients with osimertinib.



Figure 1. Overview of the GeneStrat[®] Test Sample Processing for cfDNA in the Biodesix Laboratory. Patient sample testing is initiated when whole blood is drawn and shipped to the Biodesix Laboratory. Patient samples for testing are accessioned into the Laboratory Information Management System (LIMS). Following nucleic acid extraction, samples are processed using the Bio-Rad QX200[™] Droplet Digital PCR (ddPCR) system. Numbers of mutant and wild-type copies are assessed using QuantaSoft[™] Software, and Test Result Reports (TRRs) are generated from the secure LIMS.



Figure 2. GeneStrat Test Turn Around Time (TAT) from sample receipt to result generation. TAT was compiled over a 34 month period (n=9,820). Data excludes weekends, holidays, and samples held for >24 hours due to incomplete clinical information on the Test Request Form.

Table 1. assay.

Nucleotide Sequence	Amino Acio
c.2390G>C	p.C7
c.2389T>A	p.C7

Table 2. Validated variants within the EGFR C797S cis/trans assay.

Nucleotide Sequence	Amino Aci
c.2390G>C	p.C
c.2389T>A	p.C
c.2369C>T	p.T7

Biodesix Inc., 2970 Wilderness Place, Boulder CO 80301



		R
Sample Number		
Concordance*		
Mean	Mutant Copies	
Weatt	Wild-type Copies	
	Mutant Copies	
Range	Wild-type Copies	





Data Analysis Test results to using Bio-Rad physician QuantaSoft within 72 hours Software GG Accession No: EDKA30170116-002
Physician: Dr. Physicia a 1000 7000 W00 4000 Y000 W00 Wild-type Gary Pestano, Ph.D., New York Laboratory Director Validated variants within EGFR C797S **COSMIC ID** Seauence COSM5945664 COSM6493937 COSMIC ID Sequence 797S COSM5945664 797S COSM6493937 790M COSM6240

Victoria Edwards, Leisa Jackson, Jordan M Reese PhD, Hestia Mellert PhD, and Gary A Pestano PhD

© 2019 Biodesix, Inc. All rights reserved.

Abstract # 411

EGFR C797S Multiplex Assay GFR C797S Variant Copies EGFR WT Copies Result						
0	482	Negative				
0	728	Negative				
0	390	Negative				
1.8	832	Negative				
0	736	Negative				
0	674	Negative				
0	654	Negative				
0	688	Negative				
0	450	Negative				
88	528	Positive				
Ο	0	Negative				