

Real-World Performance of Blood-Based Host Immune Profiling in First-line Immunotherapy Treatment in Advanced Stage Non-Small Cell Lung Cancer

Abstract # 9545
Poster # 311

R. Brian Mitchell¹, Patricia L. Rich², Paul R. Walker³, Eric S. Schaefer⁴, Jiahuai Tan⁵, Nagaprasad Nagajothi⁶, John W. Dubay⁷, Jason K. Boyd⁸, David N. Oubre⁹, Nadeem Ikhtlaque¹⁰, Ray D. Page¹¹, Mazen Khalil¹², Suman K. Sinha¹³, Hafez Halawani¹⁴, S. Scott Boniol¹⁵, Edgardo S. Santos¹⁶, James M. Orsini¹⁷, Emily K. Pauli¹⁸, Wallace L. Akerley¹⁹, INSIGHT investigators

¹Virginia Cancer Institute, Midlothian, VA; ²Cancer Treatment Centers of America, Atlanta, GA; ³East Carolina University, Greenville, NC; ⁴Highlands Oncology Group, Fayetteville, AR; ⁵Karmanos Cancer Institute, Detroit, MI; ⁶Mercy Medical Cancer Center, Canton, OH; ⁷Lewis & Faye Manderson Cancer Center, Tuscaloosa, AL; ⁸Southeastern Medical Oncology Center, Goldsboro, NC; ⁹Pontchartrain Cancer Center, Mandeville, LA; ¹⁰Franciscan St. Francis, Indianapolis, IN; ¹¹Center for Cancer and Blood Disorders, Fort Worth, TX; ¹²Clopton Clinic, Jonesboro, AR; ¹³Christus Health, Tyler, TX; ¹⁴Christus Cabrini Cancer Center, Alexandria, LA; ¹⁵Christus Cancer Treatment Center, Shreveport, LA; ¹⁶Lynn Cancer Institute, Boca Raton, FL; ¹⁷Essex Oncology of North Jersey, Montclair, NJ; ¹⁸Clearview Cancer Institute, Huntsville, AL; ¹⁹Huntsman Cancer Institute, Salt Lake City, UT.

Introduction

- Immune checkpoint inhibition (ICI) has improved outcomes for many patients with treatment-naïve advanced non-small cell lung cancer (NSCLC).
- However, better biomarkers are needed to predict patient response and guide treatment decisions considering added toxicity and higher cost of combination treatments.
- A prospectively designed, observational registry study (NCT03289780) assesses the ability of a clinically validated, blood-based, host immune classifier (HIC) to predict ICI therapy outcomes.
- The proteomic classifier was based on MALDI-TOF mass spectrometry and a machine learning algorithm for classification (K-Nearest Neighbor).

Methods


INSIGHT study overview:

- All enrolled subjects are tested and designated HIC-Hot (HIC-H) or HIC-Cold (HIC-C) prior to therapy initiation.
- Subject characteristics, therapeutic decisions, staging, disease monitoring metrics and available biomarker data have been collected.
- Subject follow-up occurs for up to 18 months.
- First subject first visit for was May 11, 2016
- Last subject first visit for the first 2000 subject cohort was January 15, 2019
- Database lock for interim analysis was March 2020

Statistical Methods:

A pre-specified interim analysis was performed on the first 2000 subjects enrolled in the study with at least 1 year follow up. Overall survival (OS) in months (mo) is summarized as median and 95% confidence interval (CI) and as Kaplan-Meier plots and compared between HIC-defined subgroups or between therapies by Cox Proportional hazard ratios and P values. Data from the subgroup (N=877) comprising subjects with advanced stage (IIIB and higher) NSCLC treated with first-line regimens are presented here.

Figure 1. Study subject population summary.



Results




Figure 2. Kaplan-Meier plot of OS for first-line advanced stage NSCLC subjects receiving platinum-doublet chemo (PTCT) or ICI-containing therapy (ICI). Und: Undefined.

Table 1. Subject demographics by HIC subgroup.

	HIC-Hot (N=622)	HIC-Cold (N=255)
Age (years)	P=0.075	
Mean (SD)	68.4 (10.3)	67.0 (10.0)
Median (Range)	68 (35 – 92)	66 (36 – 95)
Gender, N (%)	P=0.004	
Female	290 (46.6%)	92 (36.1%)
Male	332 (53.4%)	163 (63.9%)
Histology, N (%)	P=0.006	
Adenocarcinoma	430 (69.1%)	143 (56.1%)
Squamous	116 (18.6%)	72 (28.2%)
NSCLC Other	76 (12.2%)	40 (15.7%)
Disease Stage at Study Entry, N (%)	P=0.749	
Stage IIIB	78 (12.5%)	34 (13.3%)
Stage IV	544 (87.5%)	221 (86.7%)
ECOG Status At Study Entry, N (%)	P=0.075	
0	184 (29.6%)	38 (14.9%)
1	299 (48.1%)	133 (52.2%)
2	117 (18.8%)	66 (25.9%)
3	22 (3.5%)	13 (5.1%)
4	0	5 (2.0%)
Smoking Status, N (%)	P=0.062	
Current Smoker	230 (37.0%)	102 (40.0%)
Former Smoker	308 (49.5%)	133 (52.2%)
Never Smoker	84 (13.5%)	20 (7.8%)
PD-L1 Expression, N (%)	P=0.809	
Negative (< 1%)	144 (41.9%)	51 (39.2%)
Low (≥ 1% and < 50%)	81 (23.5%)	30 (23.1%)
High (≥ 50%)	119 (34.6%)	49 (37.7%)
N/A	278	125
Treatment Regimen, N (%)	P=0.022	
ICI Monotherapy	80 (12.9%)	37 (14.5%)
ICI Combination	116 (18.6%)	51 (20.0%)
Pt-based Chemotherapy	284 (45.7%)	108 (42.4%)
Other Chemotherapy	13 (2.1%)	6 (2.4%)
TKI	35 (5.6%)	5 (2.0%)
Radiation	25 (4.0%)	7 (2.7%)
Other Active Treatment	14 (2.3%)	2 (0.8%)
No Active Treatment	55 (8.8%)	39 (15.3%)

N/A: PD-L1 expression levels not reported in the study (excluded from percentage and Chi²), Pt: platinum, TKI: tyrosine kinase inhibitor

Results




Figure 3. Kaplan-Meier plot of OS for subjects with advanced stage NSCLC without active treatment by HIC result. HR: Hazard ratio by Cox proportional hazards.




Figure 4. Kaplan-Meier plot of OS by HIC for subjects with advanced stage NSCLC receiving first-line (A) platinum (pt) chemotherapy, (B) any ICI-containing regimen, (C) ICI monotherapy or (D) ICI and chemotherapy combinations. NR: Not reached, Und: Undefined.

Table 2. PD-L1 expression of ICI-receiving subjects by HIC subgroup.

PD-L1 Expression, N (%)	HIC-Hot (N=126)	HIC-Cold (N=64)	P Value by Chi ²
Negative (< 1%)	13 (10%)	9 (14%)	0.694
Low (≥ 1% and < 50%)	32 (25%)	14 (22%)	
High (≥ 50%)	81 (65%)	41 (64%)	

Results

Table 3. Multivariate analysis of OS of subjects receiving an ICI-containing regimen in first-line, including HIC classification and therapy type.

Covariate	vs.	HR (95% CI)	CPH p-value
HIC-Cold	HIC-Hot	0.37 (0.25-0.53)	<0.0001
ECOG PS 0	PS 1	1.48 (0.90-2.43)	0.125
	PS 2+	2.62 (1.52-4.50)	0.0005
Adenocarcinoma	Squamous Cell	1.42 (0.86-2.34)	0.168
	NSCLC Other	0.54 (0.30-1.00)	0.047
ICI+Chemotherapy	ICI Monotherapy	1.67 (1.11-2.51)	0.013
	With Other	1.94 (0.66-5.74)	0.232
PD-L1 high (≥ 50%)	PD-L1 low (≥ 1% and < 50%)	1.54 (0.91-2.60)	0.104
	PD-L1 negative (< 1%)	1.88 (1.07-3.32)	0.028
	N/A	0.99 (0.62-1.61)	0.987
Age < 65	Age ≥ 65	1.51 (1.07-2.15)	0.020
Never Smoker	Ever Smoker	1.81 (0.89-3.91)	0.112
Stage IIIB	Stage IV	1.50 (0.72-3.14)	0.283
Female	Male	1.02 (0.71-1.46)	0.933




Figure 5. Kaplan-Meier plot of OS for subjects with advanced stage NSCLC receiving first-line (A) ICI monotherapy or (C) ICI and chemotherapy combinations by PD-L1 expression and (B,D) PD-L1 high subgroup by HIC. Neg: Negative, NR: Not reached, Und: Undefined, Tx: Therapy.

Conclusions

- In a real-world clinical setting, overall survival of patients with treatment-naïve advanced stage NSCLC receiving platinum-based chemotherapy did not differ significantly from patients receiving any type of ICI containing regimen.
- Profiling the immune response to lung cancer using the proteomic host immune classifier identifies a chronic inflammatory disease state, HIC-Cold, in almost 30% of subjects which is associated with poor prognosis:
 - in the absence of treatment
 - independent of treatment type
 - independent of PD-L1 and other clinical factors
- Subjects receiving ICI and chemotherapy in combination experienced longer overall survival than subjects receiving ICI alone, with and without adjusting for PD-L1 expression and other clinical factors.
- In subjects with high PD-L1 expression:
 - HIC-Hot receiving ICI alone experienced similar overall survival to that of triplet combination.
 - HIC-Cold had shorter overall survival on ICI alone than on ICI combination indicating that HIC might identify a disease state less responsive to ICI alone.



Blood based HIC is independent of PD-L1 expression and may provide clinically meaningful information for selection of ICI treatment type or best supportive care in patients with advanced stage NSCLC.