




JSR Life Sciences

Lifetime study for JSR Amsphere A3 – Amsphere™ A3 protein A resin

Many factors are required to assess and implement a protein A resin into a mAb purification process. Dynamic binding capacity, purity, operating flow rate and caustic compatibility are a few among many evaluation parameters. Resin lifetime is a key parameter that needs to be assessed as it is directly related to both resin volume and process costs. The resin should maintain performance parameters over the intended number of process cycles. The following describes the effects of process cycling and caustic CIP on the purification performance of Amsphere A3 over 100 cycles.

Materials and Methods

- Biosimilar mAb, 1.1 g/l titer, pH 8.45, from CHO cell line
- Initial HCP level 350,000 ppm
- Clarified by centrifugation and 0.22 micro filtration
- HCP determination by Cygnus Kit (F550)
- DNA determination by Picogreen® assay (Invitrogen)
- Leached protein A determination by Cygnus Kit (F740)



Amsphere A3 is a new protein A resin designed with a surface modified base bead and alkali-resistant optimized ligand.

Protein A ligand

- High DBC via controlled conformation and orientation
- High alkaline stability from protein engineering

Surface modification

- Low HCP levels by surface hydrophilization


Base bead formulation

- High DBC at high flow rate
- Excellent pressure and flow properties via rigid crosslinking

TABLE 1: CIP CYCLE TEST CONDITIONS

CIP CYCLE TEST CONDITIONS	
Number of cycles	100 cycles
Column	0.5 x 5 cm (1 ml)
Sample	Biosimilar mAb (Trastuzumab)
Loading amount	>10% BT (Run 1, 11, 21...) 80% load of 10% BT DBC (Run 2, 12, 22...)
Loading frequency	Added to every 10 cycles
Residence time	5 min

FIGURE 1: TEST CONDITIONS OF CHROMATOGRAPHY PROCESS



Results

1. DYNAMIC BINDING CAPACITY




FIGURE 2: DBC IS MAINTAINED AT INITIAL LEVELS FOR 100 CYCLES

2. mAb RECOVERY




FIGURE 3: YIELD IS ABOVE 80% AND CONSTANT FOR 100 CYCLES

3. HCP REDUCTION




FIGURE 4: HCP REDUCTION IS APPROXIMATELY 2.0 LRV AND IS CONSTANT FOR 100 CYCLES

4. DNA REDUCTION






FIGURE 5: DNA LRV IS CONSISTENT AT APPROXIMATELY 2.9 THROUGHOUT THE 100 CYCLES


5. LEACHED PROTEIN A INTO TARGET POOL




6. ELUTION POOL



7/8. CHROMATOGRAMS COMPARISON



9. ELUTION PROFILE



Summary

The CIP test results indicate stable performance over 100 operating cycles, including NaOH exposure of 0.1 M NaOH every cycle and 0.5 M NaOH every tenth cycle. The mAb recovery, HCP reduction, DNA reduction and elution pool pH are consistent across the 100 cycles. The chromatograms from the first to the 101st cycle also show nearly identical response both in terms of elution profile and column volume.

This data set indicates stable resin performance can be expected after multiple cycles.



JSR Life Sciences

Amsphere™ is a global trademark of JSR Corporation
©2017 JSR Corporation – All rights reserved.

For more information visit www.jsrlifesciences.com

NORTH AMERICA

JSR Life Sciences
JSR Micro, Inc.
1280 North Mathilda Avenue
Sunnyvale, CA 94089
408-543-8800
bioprocess.us@jsrlifesciences.com

ASIA

JSR Corporation
1-9-2 Higashi-Shimbashi
Minatoku, Tokyo
105-8640 Japan
+81-3-6218-3557
bp@jls.jsr.co.jp

EUROPE

JSR Life Sciences
JSR Micro NV
Technologielaan 8
3001 Leuven
Belgium
+32-16-668-721
bioprocess.eu@jsrlifesciences.com