



Purification of antibody fragments and single domain antibodies with Amsphere™ A3 Protein A resin

Platteau, Gerald^{1*}; Ströhlein, Guido¹; Gaspariunaite, Vaiva²; Vincke, Cécile²; Sterckx, Yann²; Muyldermans, Serge²

1 JSR Life Sciences – JSR Micro NV, Technologielaan 8, 3001 Leuven, Belgium

2 Laboratory of Cellular and Molecular Immunology, Vrije Universiteit Brussel (VUB) – Pleinlaan 2, 1050 Elsene, Belgium

* Contact: Gerald Platteau gerald.platteau@jsrlifesciences.com

Abstract

The wide range of antibody fragment variants for which Amsphere A3 can be used as a capture resin is described, together with insight to the related binding mechanisms obtained via x-ray crystallography.

Introduction

The standard capture purification of full-size, classical antibodies is typically performed with Protein A affinity chromatography. Binding of the Protein A affinity ligand to the “Fc” region of antibodies (Abs) (abbreviations are defined at end of paper) takes place at the juncture of the constant domains 2 and 3 of the Ab heavy chains (Lewis *et al.* 2008). This high affinity binding primarily involves hydrophobic interactions. Antibodies that belong to the same subclass have greater than 95% homologous Fc-regions, allowing Protein A to be a capture platform for a wide range of antibodies (Ghose *et al.* 2005).

Following the wave of successful commercial monoclonal antibody products, various forms of antibody fragments are now becoming an important class of next-generation therapeutic proteins. This includes Fabs and fusion proteins of the Fab variable domains. From the variable domain of the heavy-chain antibodies of camelids, the VHH sdAbs have been derived. These VHHs represent some of the smallest antigen binding antibody-derived proteins. As such they are more stable than full size mAbs, can be produced in microbial organisms, and offer higher target binding events per gram of product. Due to the lack of an Fc region, these antibody fragments cannot be captured with most engineered Protein A affinity ligands. However, Amsphere™ A3 Protein A ligand exhibits a high affinity for VHH single domain antibodies (Figure 1).




Figure 1: DBC values at 3min residence time for a bivalent VHH (VHH sdAb X from Ablynx™) of Amsphere A3 and four other commercially available Protein A affinity resins (noted in regard to base matrix type). Feed material: HCCF from *Pichia Pastoris* culture. Column: 0.5 cm diameter; 5.0 cm bed height.

Compared to 4 other commercially available PrA resins, Amsphere A3 has the highest DBC. The resin that comes closest has a 40% lower DBC, and given its glass matrix, is not as caustic stable. Results in Figure 1 were obtained with a default protocol for Protein A affinity purification, which included VHH elution at pH 3 (Table 1).

Step	Solution - buffer	Flow (cm/hr)	Res. time (min)	Length (CV)
Equilibration	20mM sodium phosphate pH 7.5	300	1.0	5.0
Load	clarified feed Ablynx from <i>Pichia Pastoris</i> ; adjusted	100	3.0	Until >10% breakthrough
Wash 1	20mM sodium phosphate pH 7.5			5.0
Wash 2	20mM sodium phosphate; 0.5M NaCl; pH 7.5			5.0
Wash 3	20mM sodium phosphate pH 7.5			3.0
Elution	50mM sodium acetate pH 3.0			5.0
CIP part 1	0.1M NaOH (elution buffer replacement)	55	2.0	3.0
CIP part 2	0.1M NaOH (15 min contact time)	300	5.5	3.0
Equilibration	20mM sodium phosphate pH 7.5	300	1.0	5.0

Table 1: Protocol used for determining the DBC values in Figure 1.

Figure 2 shows the chromatographic performance of a batch purification run of the VHH with Amsphere A3, demonstrating a low elution volume with good protein recovery and HCP clearance.




Figure 2: Chromatogram of the purification of a VHH sdAb with Amsphere A3. Blue: UV280 absorbance. Brown: conductivity (mS/cm). Purple: pH. Experiments were performed on an ÄKTA™ avant 25 from GE Healthcare. A Pichia pastoris 2nd Generation HCP ELISA kit (Cygnus Technologies, cat no F640) was used for HCP measurement.

It is of interest to delve further into the interaction between the Amsphere A3 ligand and the VHH sdAbs which make the results in Figure 1 possible. Since binding to the Fab region of the mAb Herceptin® was also demonstrated (data not shown), the mechanism of binding to the VH domain was also studied. Conditions for binding of VHH to Amsphere A3 were evaluated. DBC data for VHHs were compared to those for mAbs. The use of Amsphere A3 for antibody fragment capture was compared against other currently available resin types. Finally, the range of antibody fragment variants for which Amsphere A3 can be used as capture resin is described.

Binding sites in the Amsphere A3 ligand and VH(H) molecules

In collaboration with the Laboratory of Cellular and Molecular Immunology at the Vrije Universiteit Brussel (VUB), the crystal structure of a VHH-PrA complex was obtained. Figure 3 provides a schematic overview of the different steps in the process to obtain the crystal structure.




Figure 3: Overview of the different steps in the process of obtaining the crystal structure of a complex formed between a VHH sdAb and the monomeric Amsphere Protein A ligand.

The process (Figure 3) resulted in the crystal structure shown in Figure 4. The structure of both proteins was resolved to < 2 Ångstroms and is well defined in the complex. On the left, we see the 3 helices of the Protein A monomer. On the right, we have the typical structure of a VHH, with 2 sheets of 4 and 5 β -strands respectively. The sheet with 5 β -strands corresponds to the side in a VH that normally interacts with the VL domain in a classical antibody. The 4 β -stranded sheet is solvent exposed and does not participate in antigen recognition. The 3 CDRs are the loops on top of the VHH molecule.




Figure 4: Crystal structure of Amsphere A3 Protein A ligand monomer in complex with VHH. The three helices of Protein A ligand are on the left and the β -stranded structure of the VHH is on the right. The CDRs are the loops on the top of the VHH. The amino acids involved in Protein A ligand and VHH association are shown in stick representation and colored by element (N in blue, O in red, C in grey).



Binding sites for the VHH are located in helix 2 and 3 of the Protein A ligand. 7 amino acids in the Protein A are identified as most important for the interaction with the VHH. Where the binding sites for the Fc region are in helix 1 and 2 and a glutamine in helix 2 is the only amino acid involved in both Fc and VHH binding (Graille *et al.* 2000).

Interaction sites for Protein A are located in framework regions 1 and 3 of the VHH. This is the side of the sheet with 4 β -strands. 8 amino acids of the VHH face the Protein A domain and these are mainly polar and/or charged residues: Ser17, Arg19, Lys65, Thr69, Ser71, Gln82, Asn84 and Ser85

In order to verify whether there is a similar mechanism for Fab binding, an overlay was made of the VHH used in this project, with the published crystal structure of the VH domain of a human Fab fragment (Figure 5). It is clear that these 2 domains have an excellent overlap. All amino acids interacting with the Protein A in the VHH are shared between the VHH and this VH domain. All these amino acids were also conserved in the VH domain of Herceptin (sequence alignment not shown). This strongly suggests a similar binding mechanism for both VHH and VH domains. In the VH domain of a classical antibody, the Protein A-binding region does not interact with the VL domain and is not involved in antigen binding.




Figure 5: Protein A's-eye view of the overlay of the VHH used in our project with the published crystal structure of the VH domain of the Fab fragment of a human IgM (Graille *et al.* 2000). The VHH of our complex is shown in orange, the VH domain in cyan. The CDRs are located on top and are purple colored. Side chains of amino acids that interact with the Protein A are shown in stick representation and colored by element (N in blue, O in red, C in grey).

Operating conditions

For a VHH sdAb, a range of binding conditions was tested. Related chromatography conditions are the same as in Table 1. Only the conditions of the feed were modified, as shown in Table 2. It should be noted that these conditions, representing a range of feed pH (6.9 - 7.5), conductivity (8.3 - 43.5 mS/cm) and target concentration (1.7 - 4.2 g/L), all resulted in significant DBC values.



FEED					
Type	pH	Cond. (mS/cm)	Conc. (g/L)	DBC at 10% breakthrough (mg/ml)	Capacity based on recovery in eluate (mg/ml)
HCCF (*)	6.9	8.8	4.08	17.0	16.5
HCCF higher pH	7.5	8.8	4.23	16.0	14.1
HCCF diluted	6.9	9.1	1.71	14.9	14.0
HCCF diluted high cond.	6.9	43.5	2.16	13.3	14.3
HCCF with 100mM thioglycerol	6.9	8.3	4.15	15.1	13.4

Table 2: DBC results of Amsphere A3 for VHH sdAb X from Ablynx, for a range of feed conditions. (*): HCCF = clarified feed VHH sdAb X from Ablynx (*Pichia Pastoris* supernatant). Column: 0.5 cm diameter, 5.0cm bed height (pre-packed column from Repligen).

Since the interaction with Protein A is based on charged residues, the interaction was expected to be salt dependent. It was therefore double checked if the VHH would elute at high salt concentration (Table 3 and Figure 6).

Step	Solution - buffer	Flow (cm/hr)	Res. time (min)	Length (CV)
Equilibration	20mM sodium phosphate pH 7.5	300	1.0	5.0
Load	VHH HCCF (*)	100	3.0	3.2
Wash 1	20mM sodium phosphate pH 7.5			7.0
Wash 2	20mM sodium phosphate; 1.0M NaCl; pH 7.5			5.0
Wash 3	20mM sodium phosphate pH 7.5			3.0
Elution	50mM sodium acetate pH 3.0			12.0
CIP part 1	0.1M NaOH (elution buffer replacement)	150	2.0	3.0
CIP part 2	0.1M NaOH (15 min contact time)	55	5.5	3.0
Equilibration	20mM sodium phosphate pH 7.5	300	1.0	5.0

Table 3: Protocol used to investigate conductivity dependency of binding of VHH sdAb X to Amsphere A3. (*): HCCF = clarified feed VHH sdAb X from Ablynx (*Pichia Pastoris* supernatant) with conc. 4.08 g/l – pH 6.9 and 8.8 mS/cm.




Figure 6: Chromatogram of the purification of a VHH sdAb with Amsphere A3. The sdAb is still bound to the column after applying a wash step with 1M NaCl. Blue: UV280 absorbance. Brown: conductivity (mS/cm). Purple: pH.

The sdAb is still bound to the column after applying a wash step with 1M NaCl. The binding of a Fab fragment, created by papain digest of Herceptin, was found to be salt dependent (data not shown). It is therefore recommended to use low conductivity conditions (< 5.0 mS/cm) for initial binding studies of VH domain-related targets, but possibly scout higher conductivity conditions.

Evaluating and interpreting DBC values

Compared to full size mAbs, the g/L capacity values shown for monovalent and bivalent VHHs are about 2-3 times lower (Figure 7). However, the MW of a full size IgG (approx. 150 kDa) is about 10 times larger than that of a monovalent VHH (approx. 15 kDa). The data in Figure 7, reflects a 5 minute residence time for the mAbs and only 1 minute for the VHH targets. No significant difference was found in DBC at 1 versus 5 minutes for the smaller VHH targets.




Figure 7: Blue bars: DBC values of Amsphere A3 for 3 monoclonal antibodies and 2 VHH sdAbs, expressed in gram molecules bound per liter resin at the point of 10% breakthrough. Green points: molecular weight in kDa.

The Figure 7 equivalent molar DBC data, representing the number of molecules that are bound to the resin, expressed in $\mu\text{mol/liter}$ are given in Figure 8. About 3 to 4 times more monovalent VHHs bind per multimeric Protein A ligand than for mAbs. Such results may reflect steric hindrance spatial limitations regarding binding of multiple targets per Protein A ligand.




Figure 8: Blue bars: DBC values of Amsphere A3 for 3 monoclonal antibodies and 2 VHH sdAbs, expressed in μmol molecules bound per liter resin at the point of 10% breakthrough. Green points: DBC values of Amsphere A3 for 3 monoclonal antibodies and 2 VHH sdAbs, expressed in Figure 7 in gram molecules bound per liter resin at the point of 10% breakthrough.



Comparison to other resins for antibody fragment capture

Two large resin categories are currently used for antibody fragment capture. Ion exchange and mixed mode resins separate the target molecules from impurities based on their surface charge at a given pH and salt concentration. With these resins it is possible to obtain higher binding capacities. However, the advantage of Protein A compared to these resin types is that the process development time is much shorter. Protein A chromatography requires shorter process development work. Also, since it is an affinity step, purity is already very high after the capture step, combined with product recoveries over 95%. Linked to the limited process development is the fact that affinity chromatography has the potential for use as a platform capture step.

The second resin category are affinity resins produced specifically to bind a certain antibody domain, like Protein G or Protein L (Rodrigo *et al.*, 2015). Amsphere A3 has the advantage over these resins of good caustic stability, allowing a high number of runs and thus lowering the overall cost per gram purified product. It may also allow target elution at higher pH (e. g. 3 versus 2). Table 4 gives an overview of chemical stability results related to exposure times of Amsphere A3 in different sodium hydroxide concentrations. Depending on the desired CIP and sanitization regime, the resin lifetime in number of cycles can be computed from the data in the table. The 80% DBC retention is given as it is often used as an endpoint for determining useful resin lifetime.

NaOH concentration (mol/l)	Exposure time at 22 °C (hrs)	
	At 90% DBC retention	At 80% DBC retention
0.1	106	162
0.2	37	68
0.5	18	27

Table 4: Caustic stability data for Amsphere A3. Exposure times for different sodium hydroxide concentrations are shown after which the DBC for polyclonal IgG at 4 min residence time is still 90 and 80% of the starting value.

Figure 9 summarizes the above noted advantages of Amsphere A3 over other resin types in regard to primary capture of antibody fragments.




Figure 9: Scheme showing the advantages of using Amsphere A3 over other resin types currently used for capture of antibody fragments.

Platform scope for antibody fragment capture

The studies noted above provided some insight into the molecular basis of binding of Amsphere A3 to VH(H) antibody formats. Figure 10 summarizes the 2 variable domain binding targets for Amsphere A3.




Figure 10: Overview of the variable domain binding to Amsphere A3

The literature reflects binding of Protein A to VH is restricted to molecules of the human VH3 subfamily, with no examples from any of the other gene families. In almost every case, reduced or eliminated binding can be linked to mutations in the described residues. Specifically for VHH sdAbs, 110 molecules were screened in collaboration with Ablynx NV, of which 99% showed binding to the Amsphere A3 Protein A ligand. The fact that the interacting residues are in the framework regions, (so not involved in antigen binding), and that there is no interaction with VL in a Fab, opens the possibility of mutagenizing non-Protein A interacting Ab fragments into binders (Fridy *et al.* 2015; Henry *et al.* 2016).

Figure 11 shows the various forms of antibody fragment formats for which Amsphere A3 can be a suitable capture resin. These can be divided into 3 major categories: VHH single domain antibodies, and constructs containing a VH3 domain like Fabs and scFvs.




Figure 11: Overview of the different types of antibodies and antibody fragments. The green boxes indicate the various forms of antibody fragment formats for which Amsphere A3 can be a suitable capture resin: Molecules containing a VH3 domain and VHH single domain antibodies. (original figure: Holliger & Hudson. 2005)



Conclusions

Besides the use for Fc containing antibodies and constructs, the Protein A resin Amsphere A3 can be used for purification of antibody fragment formats, like the VHH single domain antibodies, and constructs containing a VH3 domain like Fabs and different fusion proteins such as single chain variable fragments. An overview of these binding targets is given in Figure 12. The well-known robustness, high selectivity, and scalability and regulatory acceptance of Protein A based resins make Amsphere A3 an ideal platform technique, minimizing DSP development time.

Abbreviation list

Ab:	Antibody	mAb:	Monoclonal Antibody
CDRs:	Complementarity-Determining Regions	MW:	Molecular Weight
COGs:	Cost Of Goods	PD:	Process Development
DBC:	Dynamic Binding Capacity	PrA:	Protein A
DSP:	Downstream Processing	scFv:	single-chain variable fragment
Fab:	Fragment Antigen-Binding	sdAb:	Single Domain Antibody
Fc:	Fragment crystallizable	VH:	Variable domain of the Heavy chain
HCCF:	Harvested Cell Culture Fluid	VH3:	Variable domain of the Heavy chain of the gene family 3
HCP:	Host Cell Protein	VHH:	Variable domain of the Heavy chain of a (camelid) Heavy chain antibody
IgG:	Immunoglobulin G	VL:	Variable domain of the Light chain
kDa:	kilodalton		
LRV:	Log Reduction Value		

References

- Fridy *et al.* 2015. Engineered high-affinity nanobodies recognizing staphylococcal Protein A and suitable for native isolation of protein complexes. *Analytical Biochemistry* 477, 92–94.
- Ghose *et al.* 2005. Antibody Variable region Interactions with Protein A: Implications for the Development of Generic Purification Processes. *Biotechnology and Bioengineering* vol. 92, nr. 6.
- Graille *et al.* 2000. Crystal structure of a *Staphylococcus aureus* protein A domain complexed with the Fab fragment of a human IgM antibody: Structural basis for recognition of B-cell receptors and superantigen activity. *Proc. Natl. Acad. Sci.* 97, 5399-5404.
- Henry *et al.* 2016. A Rational Engineering Strategy for Designing Protein A-Binding Camelid Single-Domain Antibodies. *PLOS ONE* | DOI:10.1371/journal.pone.0163113.
- Holliger & Hudson. 2005. Engineered antibody fragments and the rise of single domains. *Nature biotechnology* volume 23, number 9.
- Lewis *et al.* 2008. A Common Theme in Interaction of Bacterial Immunoglobulin-binding Proteins with Immunoglobulins Illustrated in the Equine System. *The Journal of Biological Chemistry* vol. 283, NO. 25, pp. 17615–17623, June 20.
- Rodrigo *et al.* 2015. Antibody fragments and their purification by Protein L affinity chromatography. *Antibodies*, vol. 4, pp. 259-277.



Contact info

NORTH AMERICA

JSR Life Sciences
1280 North Mathilda Avenue
Sunnyvale, CA 94089
408-543-8800
bioprocess.us@jsrlifesciences.com

ASIA

JSR Corporation
1-9-2 Higashi-Shimbashi
Minatoku, Tokyo
105-8640 Japan
+81-3-6218-3557
bp@jls.jsr.co.jp

EUROPE

JSR Life Sciences
JSR Micro NV, Technologielaan 8
3001 Leuven
Belgium
+32-16-668-721
bioprocess.eu@jsrlifesciences.com

Amsphere™ is a registered trademark of JSR Corporation.

©2017 JSR Life Sciences All rights reserved.

For more technical and ordering information on Amsphere A3, please visit www.jsrlifesciences.com.