


HemaSpot-SE™

Serum collection at your fingertips!



Stick finger with lancet.

Add 3 drops of blood.


Wait 2 minutes.

Close lid and ship.


Shaping the Future of Biosampling

HemaSpot-SE is an easy to use patent-pending device that simplifies blood sample collection and provides robust storage/shipping of blood samples at ambient temperatures.

Sample Analysis



Samples punched (positions A, B, C, D) from blood separated and dried on HemaSpot-SE show high DNA levels only at position A. RNA is found in both cell (A,B) and non-cell (C,D) positions. The small molecule Nifedipine shows similar levels for all positions.



For Research Use Only. Not For Use in Diagnostic Procedures


Call us at **512.827.9627** for a sample or email us at info@SpotOnSciences.com. For more information on any of our products or services please visit www.SpotOnSciences.com

HemaSpot-SE™

Serum collection at your fingertips!

Innovative Sampling

HemaSpot-SE's patent-pending blood separation device simplifies sample collection and improves separation of blood.



Applications

- Clinical Trials
- Medical Research
- Military field medicine
- Disaster & public health crises
- Disease surveillance
- Epidemiology & population studies
- Biobanking
- Biomarker studies
- Space medicine

And many more uses!

Features & Benefits

HemaSpot-SE's solution improves collection techniques and results by providing unique features:

- **Robust**
Protects the sample from contamination
- **Safe**
Reduces biohazard exposure
- **Fingerstick Simplicity**
No Special training required
- **Easy**
Ambient temperature storage and shipping
- **Fast**
Takes only 5 minutes from sampling to shipping
- **Stable**
Sample stays stable at ambient temperatures

Contact Us

Call us at 415-347-3412 for samples or email us at info@SpotOnSciences.com For more information on any of our products or services please visit www.SpotOnSciences.com.