# SISU, PRODUCT DROUGTALOG

MADE IN THE USA



**MORE PROTECTION.** LESS MOUTHGUARD.<sup>™</sup>

sisuguard.com



SMALL

Up to 5' (150 cm)

# UP TO \$50,000 DENTAL WARRANTY ON ALL OF OUR MOUTHGUARDS

CUSTOM FIT AND REMOLDABLE UP TO 20 TIMES

## AERO FOR TEAM SPORTS

At a lean 1.6mm, the ultralight design of the Aero is 50% thinner than conventional sports mouth guards, providing maximum comfort and convenience. Combined with its quick, custom and remoldable fit, you'll entirely forget it's there.

For smaller mouths and children

Available in 13 colors:















### MAX FOR HIGH-IMPACT SPORTS

At 2.4mm thin, the SISU Max Mouthguard offers the most protection of any SISU guard, but is still 30% thinner than conventional sports mouth guards. It's customfitted so you'll be able to talk, breathe and drink better than ever before. It's really the best of both worlds.

> **ONE SIZE** Over 5' (150 cm) For average and larger mouths







Made in the USA



#### GO FOR ESSENTIAL PROTECTION

Smaller in size and price, the all new SISU GO was made to be both discreet and comfortable while still providing the necessary canine to canine coverage to protect where most dental injuries happen. The SISU protection athletes expect, at a more affordable price point.

**ONE SIZE** Over 4' (120 cm)

One size fits most



Available in 5 colors:



Made in the USA

### SENSE ELECTRONIC HEAD IMPACT TRACKER

With SISU Sense<sup>™</sup>, we've taken our mission to protect the athlete to the next level. By infusing our HitSense<sup>™</sup> Technology into our industry-leading mouthguard, SISU Sense measures every hit an athlete absorbs, then transmits the data to our proprietary app. The result: protection on the field today, and tomorrow, that makes sense.





Made in the USA



1-time use

# ACCESSORIES



#### Made in the USA

## FLAT2FIT<sup>TM</sup> TECHNOLOGY

SISU Mouthguards use Flat2Fit technology to create a completely custom fit to your mouth. The material comes flat in the packaging, and within minutes, transforms into a perfectly fit mouth guard ready to protect you from whatever comes your way.



Remoldable up to 20 times



# POWERED BY DIFFUSIX **TECHNOLOGY**

## THE SISU ADVANTAGE:

SISU Mouthguards are made from a non-compressible, perforated material that distributes impact force over the entire surface of each mouthguard. With their custom, remoldable fit and unmatched strength, SISU Mouthguards are able to better protect your teeth while providing you the freedom to easily talk, breathe and drink.



Conventional mouthguards are made form a material that is too soft to fully protect your teeth from harm.



With our one-of-a-kind Diffusix<sup>™</sup> Technology, SISU Mouthguards are made from a material that absorbs and diffuses impact throughout the entire mouthguard, meaning that much less impact force is transferred to your teeth.



## **QUICK-FIT GUIDE**

#### Notice:

If you have braces, **DO NOT** attempt to fit your SISU Mouthguard at home. Instead, contact your orthodontist (see below for more info). If you don't have braces, let's get started.

- 1 Wash your hands and the SISU Mouthguard with soap and water. Heat water in a small bowl until hot to the touch (160°F / 71°C). Place your SISU Mouthguard into the hot water until soft and pliable, approx. 30-60 seconds (use utensil to determine softness).
- 2 Once soft, use the utensil to carefully lift the SISU Mouthguard out of the water. DO NOT let the softened mouthguard fold on itself the material will glue together, making it difficult to reshape.
- **3** Balancing the pliable guard with your thumbs, lean in toward the mirror and open your mouth just wide enough to insert the guard. Insert the SISU Mouth-guard and bite down gently on the bite pad in the middle of the guard.
- 4 Hold the mouthguard in place using your teeth and tongue. With your fingers, fold the guard straight up. Be sure to mold the guard up to the gum line on the front and sides of your teeth to create a solid shield.
- **5** Once the SISU Mouthguard is in the proper position, close your lips over the guard and gently create suction until the material has hardened to a snug, perfect fit. Leave the guard in your mouth for 3-5 minutes to allow time to set.









## **GOT BRACES?**

If you have braces or are undergoing other orthodontic treatment, **DO NOT** attempt to fit your SISU Mouthguard at home. Instead have your dentist or orthodontist fit your SISU Mouthguard for you! He/she will use the appropriate barriers between the guard and braces so that the material, while pliable, doesn't bond to your brackets.

We also recommend that no matter your child's age, all athletes use an adult-sized SISU Mouthguard whenever wearing braces in order to ensure a proper fit.

Visit sisuguard.com/fitting-with-braces for more information.



#### THE OFFICIAL MOUTHGUARD OF

# **MOUTHGUARD COMPARISON**



	AERO			MAX	GO	SENSE
Sizes	Small	Medium	Large	One Size	One Size	One Size
Athlete Height	Up to 5' (150 cm)	5′ - 6′ (150 cm - 180 cm)	Over 6′ (180 cm)	Over 5′ (150 cm)	Over 4' (120 cm)	Over 5′ (150 cm)
Works with braces	<b>No</b> (Use Medium or Large)	Yes	Yes	Yes	Yes	Yes
Thickness	1.6 mm			2.4 mm	1.6 mm	2.4 mm
Dental warranty	Up to \$50,000			Up to \$50,000	Up to \$50,000	Up to \$50,000
Uses	All sports			High-impact sports	Essential Protection	Head Impact Tracking
Best for	Lacrosse Field Hockey Ice Hockey Football Basketball Roller Derby Team Sports And more!			Rugby MMA Boxing Karate/BJJ Hockey Football Impact Sports And more!	Basketball Soccer Cycling Water Polo Baseball Volleyball Action Sports And more!	Rugby MMA Lacrosse Roller Derby Hockey Football Impact Sports And more!
MSRP	\$24.99			\$34.99	\$11.99	\$99.99

Research suggest the taller you are, the larger your dental arch. And thus, a bigger mouthguard is required as height increases. Our new and improved sizing system uses the findings of extensive research and testing to suggest the proper guard for you based on how you measure up.

\*SISU Mouthguards are BPA, latex, PVC, and phthalate free, and FDA compliant.

## CONTACT INFORMATION

Questions? Please contact our team at 800.444.0570 or send an email to customerservice@sisuguard.com

sisuguard.com