



ADULT TOOTHBRUSHES:

Not all toothbrushes are the same. We offer a variety of clinically proven products to meet the needs of those with healthy gingiva, restorations, orthodontics, periodontal disease or post surgery.



SUNSTAR

G·U·M
HEALTHY GUMS. HEALTHY LIFE.®

Technique[®] PRO

45°

0.01



GUNZ
dental

P 1800 025 300
F 1800 025 202
www.gunz.com.au

HEAL

SUNSTAR




Extremely tapered bristle design

Unique bristles tapered to a 0.01mm tip to gently remove biofilm more effectively¹


7x along the sulcus




7x beneath the sulcus



1.6x between the teeth



0.01



45°

45° angle Quad-Grip handle for optimal Bass/Modified Bass brushing technique



Technique Pro Manual:

ORDER CODES:

- GUM525 Technique Pro Soft Compact Head Pkt 1
- GUM526 Technique Pro Medium Compact Head Pkt 1
- GUM526CY Technique Pro Soft Compact Head Bulk (Box 96)
- GUM528CY Technique Pro Medium Compact Head Bulk (Box 96)

The GUM® Technique PRO manual toothbrush is designed for optimal gum health control. Allows patients to master the proper brushing technique while benefiting from superior efficacy and gentle cleaning :

- 1) Below the gum line
- 2) Along the gum line
- 3) Between teeth

Feature	Benefit
Extremely tapered bristles	Deeper penetration below the gumline (subgingival), along the gumline (gingival margin site) and between teeth (interproximal) to remove plaque in areas not effectively reached by ordinary bristles
Bi-level bristles	Multi-level, effective cleaning: <ul style="list-style-type: none"> • Below the gumline, along the gumline and between teeth (thanks to extremely tapered bristles) • Exposed tooth surfaced (thanks to regular bristles)
Angled bristles	<ul style="list-style-type: none"> • Enhance between teeth penetration and cleaning by striking the surface in different directions
Patented, no slip Quad-Grip handle	<ul style="list-style-type: none"> • Promotes the modified Bass technique (recommended by dentists) - Helps brushing at 45° angle for optimal below the gumline cleaning.

¹ than conventional bristles.
1. In-Vitro Test, YRC Inc., September 2008



P 1800 025 300

F 1800 025 202

www.gunz.com.au



S-Class Soft Assorted Colour Range



4 Row Assorted Colour Range



CareDent

Classic
design, great
quality
brush range

Caredent Adult Toothbrushes:

ORDER CODES:

CAR1015 Caredent Toothbrush S-Class Soft x48 Assorted
CAR3050 Caredent Toothbrush 4 Row Adult x72 Assorted

Classic designed adult toothbrush featuring a compact head to reach all those areas in the mouth that other brushes cannot reach. Made with premium Dupont Tynex bristles, this brush has a timeless design and highly recommended by dentists, hygienists and dental professionals

Pink Ribbon:

ORDER CODES:

CAR1040 Caredent Toothbrush Pink Ribbon Retail x12
CAR1045 Caredent Toothbrush Pink Ribbon Cello Wp x48

Caredent is an official partner of the Cancer Council, and is proud to introduce the Pink Ribbon Toothbrush. Each brush purchased is a contribution to Cancer Council. The Pink Ribbon brush is available in Pink for her and Blue for him. An elegant crystal clear handle with a TPE ribbon design, a small compact head, an easy grip handle together with quality end rounded bristles.



P 1800 025 300

F 1800 025 202

www.gunz.com.au



SEE THE GUM SENSIVITAL IN ACTION:

<http://bit.ly/1C4gmbT>



SensiVital Toothbrush:

ORDER CODE: GUM509 SensiVital Toothbrush x1

The GUM® SensiVital Toothbrush gently removes plaque for an absolute clean and smooth sensation; particularly recommended for sensitive gums. Features include:

- Dual-level bristle trim for 2 levels of cleaning
- Longer micro-thin bristles reach deep into interdental spaces and along the gum line without harming gums
- Shorter bristles offer maximum stability to clean the entire tooth surface
- Ultra Soft bristles (0.15mm) with tapered tip are designed to be ultra-gentle on sensitive teeth and gums



Post-Surgical Toothbrush:

ORDER CODE: GUM317 Delicate, Post-Surgical Toothbrush x1

After oral surgery it is important to keep the mouth clean, especially to help the wound heal uneventfully. The GUM® Post-Surgical toothbrush has been designed with those instructions in mind with very gentle but highly effective removal of plaque and food debris from around the wound site without dislodging sutures.

Care Sensitive Toothbrush:

ORDER CODE: CAR6593 Care Sensitive U/Soft Toothbrush x12

Caredent ultrasoft sensitive toothbrush with micro bristles. A large TPE thumb grooved handle for ease of use and a comfortable grip. Great for those patients with special needs.



P 1800 025 300

F 1800 025 202

www.gunz.com.au

CareDent

Brushes for those on the go!



HastyPasty Travel Brush Colour Range

Travel Brushes:

ORDER CODE:

CAR6400 Caredent HastyPasty Travel Brush. Pkt 100.

Simple but effective pre-pasted fresh mint toothbrushes from CareDent. All individually wrapped for added hygiene benefits. Great for those travelling or simply on the go.

ORDER CODE:

CAR7146 Caredent TraveLine Toothbrush. Pkt 48 Assorted.

A timeless travel brush with a plastic cap which extends to a handle. Quality bristles and great for those travelling on business, holidays, or weekends away.



TraveLine Brush Colour Range

Brushing with a new v old manual toothbrush: The general rule is to replace manual toothbrushes every 2 or 3 months. Some professionals suggest otherwise.

READ MORE:

<http://bit.ly/1tiPxI7>



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Battery Powered dual-action



Spinbrush™ ProClean

ORDER CODE: ARM2304 Spinbrush ProClean Medium x1
ORDER CODE: ARM2305 Spinbrush ProClean Refill Med, Pkt 2

SpinBrush™ ProClean™ battery powered toothbrushes feature dual-action technology. A circular oscillating head spins to clean and polish even between teeth. The secondary brush bristles move up and down to scrub away plaque and massage along the gumline. Together both heads deliver up to 3,700 brush strokes per minute for a deeper more effective clean.

- Removes up to 70% more plaque in hard-to-reach places
- Features dual-action technology for a combination of spinning and side-to-side brushing head movements
- Ergonomic handle for superior comfort and control
- Replacement heads available
- Battery-powered, includes 2 replaceable AA batteries



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Battery Powered whitening



Spinbrush™ ProWhitening

ORDER CODE: ARM2306 Spinbrush ProWhitening Med x1
ORDER CODE: ARM2307 Spinbrush ProWhitening
Refill Med, Pkt 2

SpinBrush™ ProWhitening™ is powered by patented Dual Action Technology - simultaneous side-to-side spinning upper head and up and down motions on the secondary head - works with specialty whitening bristles to gently polish away surface stains as it cleans.

- Whiter teeth in just 7 days
- Innovative whitening bristle system
- Ergonomic handle for superior comfort and control
- Replacement heads available
- Includes 2 replaceable AA batteries




P 1800 025 300

F 1800 025 202

www.gunz.com.au



Battery Powered sensitive



Spinbrush™ Pro Sensitive

ORDER CODE: ARM2333 Spinbrush ProSensitive Brush x1
ORDER CODE: ARM2334 Spinbrush ProSensitive
Refill Med, Pkt 2

SpinBrush® ProSensitive is a battery powered toothbrush that is gentle on the gums. Specially designed for sensitive teeth and gums, the brush has been made with soft-flex bristles that are longer and more flexible for a gentler clean. Spinbrush ProSensitive has also been equipped with a gum guard motor ideal for sensitive teeth and gums and a precision brush head for control around sensitive areas.

The highest levels of plaque in the mouth: Studies have reported that the area in the mouth most likely to have been missed with tooth brushing is the lower right lingual.

READ MORE:

<http://bit.ly/1BtfVsF>



P 1800 025 300

F 1800 025 202

www.gunz.com.au



SPECIALTY TOOTHBRUSHES:

Clinically proven products to meet the needs of patients undergoing more specialised care – such as orthodontics, implants and periodontal surgery.




GUNZ
dental

P 1800 025 300
F 1800 025 202
www.gunz.com.au



Orthodontic Toothbrush:

ORDER CODE: GUM124 Orthodontic V Trim, Pkt 6

The GUM® Orthodontic toothbrush has been developed for safe and effective brushing of teeth fitted with orthodontic appliances including braces, brackets, tubes, and wires.

The unique V-cut bristle design perfectly adapts to the anatomy of the teeth with braces and gently cleans in and around the orthodontic appliances. This special configuration of bristles reaches into the crevices between the gums and the brackets, as well as below the brackets, under the wires, and between the brackets themselves. This makes it easy for your patient to care for those hard-to-reach areas that you as a dental professional are most concerned about. The GUM® Orthodontic toothbrush keeps braces as well as teeth and gums clean, helping to prevent dental decay and gum diseases.

The sleek handle offers excellent adaptation to the size of the hand and adds comfort and control while brushing. It helps to position the GUM® Orthodontic toothbrush at the correct 45-degree angle for optimal and easy cleaning.

Brushing twice daily with the GUM® Orthodontic toothbrush will keep your patient's mouth clean and healthy while they acquire the best results from orthodontic treatment!



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Post-Surgical Toothbrush:

ORDER CODE: GUM317 Post Surgical Brush x1

After oral surgery it is important that patients keep their mouth clean, especially to help the wound heal uneventfully. As an oral surgeon, dentist or dental hygienist you will have given your patient individualised advice about how to look after the surgical wound. The GUM® Post-Surgical toothbrush has been designed with those instructions in mind – to help keep the healing wound clean.

The main purpose of the GUM® Post-surgical toothbrush is to be very gentle but highly effective in removing dental plaque and food debris from near the healing wound and any sutures keeping the wound closed. The brush is designed to be used until the surgical site is fully healed. It is used in addition to the patient's normal oral hygiene routine of toothbrushing and between-teeth cleaning for all teeth and gums not involved in the surgery. This keeps the rest of the mouth clean and healthy and helps reduce the buildup of potentially infectious bacteria around the surgical site.

Being so soft and gentle, the GUM® Post-Surgical toothbrush can be used around mouth irritations and ulcers and to clean exposed roots that may have sensitivity.

Hygienists Shouldn't Fear Patients with Implants: Hygienists must treat patients with titanium implants almost every day. Don't risk scratching the titanium abutment.

READ MORE:

<http://bit.ly/1uVg5Kj>



P 1800 025 300

F 1800 025 202

www.gunz.com.au



End-Tuft Toothbrush:

ORDER CODE: GUM308 End-Tuft Toothbrush x1

The GUM® End Tuft toothbrush is clinically designed for hard-to-reach or simply difficult-to-clean areas where the dental plaque buildup will become harmful to teeth or gums.

The small round brush head comprises seven tufts of tightly packed soft nylon bristles, trimmed so the bristles in the centre can reach deeper into small spaces. The brush handle is ergonomically designed for a firm grip, giving the control and precision necessary to clean where most other cleaning aids cannot reach.

The GUM® End Tuft toothbrush is recommended for patients who have the following:

- Exposed spaces between the roots of molar teeth due to periodontitis and gum recession, keeping furcations clean
- Complex fixed dentures, bridges or dental implants that are inherently hard to clean
- Orthodontic appliances
- Difficulty in cleaning the back surfaces of the last teeth at the end of the dental arch at the back of the mouth
- Other hard-to-reach areas that they only your dentist or dental hygienist can identify

Use the GUM® End Tuft toothbrush in addition to normal toothbrush and between-teeth cleaning aids to clean those hard-to-reach and difficult-to-clean areas.



P 1800 025 300

F 1800 025 202

www.gunz.com.au

SUNSTAR

G·U·M
 HEALTHY GUMS. HEALTHY LIFE.™

premier®

“Perfect
Partners”
for denture
care



Denture Toothbrush:

ORDER CODE: GUM201 Denture Toothbrush x1

The GUM® Denture brush is recommended for the daily care of removable dentures and acrylic retainers. The GUM® Denture brush consists of two differently configured brush heads: a flat bristled head for smooth surfaces and a single-tufted head for hard-to-reach areas. The strong flat bristles are extremely resilient and efficiently remove dental plaque, food particles and other debris from the inner and outer surfaces. The single-tufted head consists of tapered bristles, which provide superior cleaning in the hard-to-reach areas regular toothbrushes cannot reach.

The GUM® Denture brush has a comfortable ergonomic handle with lever-grip superbly adapted for people with limited dexterity. The textured handle provides wide thumb support and sturdy handgrip, which assists in easy and precise brushing.



Premier ProClean:

ORDER CODE: PRE11103 Pro Clean Tartar & Stain Remover Pkt 25 sachets

The ProClean™ offers dental professionals a highly effective tartar and stain remover for dental appliances. Use in-office and dispense for patient home use, ProClean is ideal for patients who have retainers and night guards. Since these appliances are worn when sleeping, the patient usually does not bring these to the dental office to be cleaned. A brief 2 hour soak once a week with Pro Clean is recommended. One box of 25 sachets is ideal between 6 monthly hygiene visits.



GUNZ
dental

P 1800 025 300

F 1800 025 202

www.gunz.com.au



SUNSTAR



Dual action
the key to a
clean, fresh
mouth

Tongue Cleaner:

ORDER CODE: GUM760 Tongue Cleaner x1

The GUM® dual action tongue cleaner efficiently removes dental plaque, food debris, and dead cells from the surface of the tongue. One side of the tongue cleaner consists of two rows of bristles that reach into the tiny folds and crevices on the tongue surface to dislodge the dental plaque and debris. The other side of the cleaner consists of two rows of scrapers that remove the loosened plaque and debris for thorough cleaning.

The GUM® dual action tongue cleaner has a narrow head that is comfortable and minimises sensitivity and the gag reflex. The long handle provides easy access to the far-back areas of the tongue. The comfort grip handle design enables firm control and precision guidance for both brushing and scraping actions.

Dual action for thorough, deep, and gentle tongue cleaning! Include tongue cleaning as part of your twice daily toothbrushing and once daily between-teeth cleaning regimen.



P 1800 025 300

F 1800 025 202

www.gunz.com.au

CareDent

Caredent specialty toothbrush range



Orthodontic V Toothbrush Denture Brush

Caredent Orthodontic V Brush:

ORDER CODE: CAR1016 Caredent Orthodontic V x12

Designed for orthodontics. This unique toothbrush comes with a compact head together with special manufactured super soft.

Caredent Denture Brush:

ORDER CODE: CAR5370 Caredent Denture Brush x6

Clever ergonomically designed denture brush with extra control in mind during the design process. Featuring two brush heads and a rubber polishing tip helps polish the denture surface. Long multi tuft bristles for effective cleaning. Small angle trimmed bristles help clean hard to reach areas and grooves. Large easygrip handle provides the extra control. Great value.

Clinical Abstract: The effect of topical fluorides on decalcification in patients with fixed orthodontic appliances.

READ MORE:



<http://bit.ly/1EMUZxY>



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Gum Disclosing Tablets:

ORDER CODE: GUM824 Red-Cote Disclosing Tablets x12/blister

ORDER CODE: GUM800-BULK Red-Cote Disclosing Tablets x250

The GUM® Red Cote disclosing tablets highlight harmful bacterial plaque on teeth and show where more brushing and flossing are needed. Cherry flavour.

Caredent Disclosing Tablets:

ORDER CODE: CAR6630 Caredent Disclosing Tablets x100

See the plaque you need to brush away. An old favourite. CareDent disclosing tablets are an excellent aid to help identify plaque. Safe and simple. Erythrosine based. Great for Adults and Kids!

Tepe Disclosing Tablets:

ORDER CODE: ODP84676 Tepe Plaque Disclosing Tabs x10

TePe PlaqSearch™ – efficient plaque control. To keep gums and teeth healthy it is important to regularly remove plaque through toothbrushing and interdental cleaning.

TePe PlaqSearch™ colours plaque to show where you need to brush better.

Quintessence International Clinical Study: Closing the gap between oral hygiene and minimally invasive dentistry; a review on the resin infiltration technique.

READ MORE:

<http://bit.ly/1EMWRXj>



P 1800 025 300

F 1800 025 202

www.gunz.com.au



CHILDREN'S TOOTHBRUSHES:

Make brushing fun! Help promote a lifetime of happy healthy oral care with this great range of manual and power toothbrushes designed especially for children.




GUNZ
dental

P 1800 025 300
F 1800 025 202
www.gunz.com.au



SUNSTAR



0-6: Good oral hygiene habits start young

0-2
years

Baby Toothbrush:

ORDER CODE: GUM213 Gum Baby Toothbrush 0-2 yrs. Pkt 12

- Extra small head for better fit in baby's mouth
- Longer handle to facilitate parent's hold of the brush
- Gentle feathered bristles thoroughly remove plaque from tooth surfaces, interdental spaces and gum line
- Split tips reach into pits and fissures, where cavities can easily start

3-6
years

Kid's Toothbrush:

ORDER CODE: GUM214 Gum Kids Toothbrush 3-6 yrs. Pkt 12

- Small head adapted to child mouth size
- Small, rounded handle to allow for all types of grip by the young child
- Gentle feathered bristles thoroughly remove plaque from tooth surfaces, interdental spaces and gum line
- Split tips reach into pits & fissures, where cavities can easily start
- Comes with hygiene cap



P 1800 025 300

F 1800 025 202

www.gunz.com.au



SUNSTAR

G·U·M
HEALTHY GUMS. HEALTHY LIFE.®

7-10+: Good oral hygiene habits start young

7-12 years

Junior Toothbrush:

ORDER CODE: GUM215 Junior Toothbrush 7-12 yrs. Pkt 12

- Longer and thin handle neck for better access to back teeth
- Handle with thumb pad for easy grip and better control of brush

10+ years

Technique Junior:

ORDER CODE: GUM221 Technique Jnr Brush 10+ yrs. Pkt 12

- Designed to help develop a proper brushing technique
- Soft bristles with Dome Trim® design
- Quad-Grip handle guides the child's hand to correctly position the toothbrush at a 45° angle for the most effective below-the-gum line cleaning



45°

45° angle Quad-Grip handle for optimal Bass/Modified Bass brushing technique




GUNZ
dental

P 1800 025 300

F 1800 025 202

www.gunz.com.au



CareDent

Brushes
designed
with fun in
mind!



Sparkle Toothbrush Assorted Colour Range



Sparkle Toothbrush with Timer



Octopus Toothbrush Colour Range

Sparkle Toothbrush:

ORDER CODE: CAR2048 CareDent junior Sparkle x72 Asst
ORDER CODE: CAR5360 CareDent Sparkle with Timer x6

A colourful classic. A well designed toothbrush kids will love. Available in 6 vibrant colours. Perfect for boys and girls. Compact head with soft bristles for thorough and gentle clean. 2 minute timer packs option also available.

Octopus Toothbrush:

ORDER CODE: CAR5350 CareDent Octopus x48

Cool and quirky toothbrush for kids. Designed with fun in mind, encouraging kids to brush. Lights flash for one minute to indicate the time required for correct upper and lower brushing.




GUNZ
dental

P 1800 025 300

F 1800 025 202

www.gunz.com.au

CURADEN  SWISS
LEADING ORAL CARE PEOPLE

For perfect
and painless
cleaning
everytime!



CuraKid Super Soft Colour Range

CuraKid Super Soft:

ORDER CODE: CURCK4260 Super Soft Toothbrush x1

Ultra gentle Curen® bristles guarantee perfect and painless cleaning of the child's small teeth. No risk of injury due to hard or sharp surfaces, thanks to the small, tapered head - tufted with ultra gentle Curen filaments and fully protected with soft rubber coating.

The anti-slip handle fits perfectly in the hands of both the small child and the parents. The practical flat stand for hygienic storage for the toothbrush.

Mum's Oral Health Predicts Infants

Oral Health: Mothers share many health outcomes with their children including their oral bacteria.

READ MORE:

<http://bit.ly/1y1CCAaw>



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Fun power brushes for boys & girls!



Spinbrush™ MyWay!™

ORDER CODE: ARM2320 Spinbrush My Way Girls, 1 pce
ORDER CODE: ARM2321 Spinbrush My Way Boys, 1 pce

The Arm & Hammer™ Spinbrush™ My Way!™ Crazy Wrapz™ for Girls & Boys is the battery-powered toothbrush that makes brushing fun, while letting kids express themselves. The brush comes with 141 stickers for kid-friendly personalisation. When kids have fun, the brushing gets done!

- Includes 27 fun, water-safe stickers
- Kids brush 38% longer when using a Spinbrush™ battery-powered toothbrush vs. using an ordinary manual brush*
- Fixed lower bristles with oscillating upper bristles
- Smaller toothbrush head for smaller mouths
- Ergonomically designed with a comfortable grip
- Battery-powered: includes 2 replaceable AAA batteries

* Children under 3 should use under adult supervision



P 1800 025 300

F 1800 025 202

www.gunz.com.au



FLOSS RANGE:

A comprehensive range of flossing options with something for all oral care needs, including tight contacts, below the gum line, implant and orthodontic care.



P 1800 025 300
F 1800 025 202
www.gunz.com.au

SUNSTAR



Perfect for
tight
interdental
contacts



Gum Expanding Floss:

ORDER CODES:

GUM2030 Expanding Floss (30m) x1

GUM2030-BULK Expanding Floss (10m) x72

Thin and easy to insert even between the tightest contacts. The floss expands during use allowing the floss to adapt to varying sizes of interdental spaces for more surface contact. Gentle on gingival tissue, perfect for sensitive teeth & gums.



Gum Butler Weave Floss:

ORDER CODES:

GUM1155 Butlerweave Floss waxed (50m) x1

GUM1855 Butlerweave Floss waxed, Mint (50m) x1

Unique woven floss offers great durability. The patented "woven" design allows the floss to spread out on tooth surface to clean and remove plaque effectively. The wide, flat profile of the flossing fibres allow it to glide between tight contacts.



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Special fibres for orthodontic & implant



Gum Fine Floss:

ORDER CODE: GUM1555 Gum Fine Floss (50m) x1

Thin, shred-resistant floss for cleaning between tightly spaced teeth and below the gum line. Special 'S' twist fibres for both durability and spreadability to effectively remove and carry away plaque.



Ortho+Floss

Gum Ortho+Floss:

ORDER CODES: GUM3220 Ortho+ Floss. Pkt 6

Ideal for orthodontic appliances. Built in threader allows easy insertion of the floss under and around the orthodontic appliances. Floss expands to adapt to different sizes of interdental spaces for more tooth surface contact.



Implant+Floss

Gum Implant+Floss:

ORDER CODES: GUM3200 Implant+ Floss. Pkt 6

Ideal for implants, bridges, crowns and wide interdental spaces. With built in threader.



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Easy Flossers



250
USES

Flosbrush



Gum Easy Flossers:

ORDER CODE: GUM890 Easy-Flossers, Pkt 6

Convenient and disposable, special floss slides easily between teeth and resists shredding, breaking or sagging. Multiple filaments provide a textured surface for more tooth contact for effective plaque removal. Angled-neck to access posterior teeth.

Gum Flosbrush Automatic:

ORDER CODE: GUM874M Flosbrush Automatic x1 (250 uses)

Easy to use. Just click the lever to advance a clean piece of floss. Angled design reaches back teeth. Unique and innovative design with extra shred-resistant floss fibre. Interior chambers separate new and used floss. With hygiene cap.




Gum Eez-Thru Floss Threaders

ORDER CODE:

GUM840 Eez-Thru Floss Threaders, Pkt 25

GUM840-BULK Eez-Thru Floss Threaders, (100 envelopes of 5)

Versatile nylon loop helps patients thread floss through spaces where most other products cannot reach; under bridges between connected crowns and around orthodontic and implant appliances.



CareDent

For your
everyday
flossing
needs!



PerioTape Minis



PerioTape Regular



PerioTape Clam

Caredent PerioTape:

ORDER CODES:

- CAR3004 Caredent PerioTape 4m Minis x48
CAR3100 Caredent PerioTape 100m Clam x6
CAR3030 Caredent PerioTape 30m x6pk

Premium PTFE tape. Available in a great looking easy grip pack featuring an easy to use flip top lid or colourful Mini's. PerioTape is a unique shed resistant fibre that slides easily between your teeth without shredding.

Caredent EeziFlosser:

ORDER CODE: CAR5515 Caredent EeziFlosser 24x6

Caredent EeziFlossers are manufactured using a quality polyethylene fibre. The ergonomically designed handle incorporates a large thumb sized pad handle for a comfortable grip. Each flosser incorporates a toothpick at the tip for an additional benefit.



Award Winning Caries Detection, Right Inside Your Surgery: Oral disease continues to present a major public health issue. Early detection is the key to reversing this trend.

READ MORE:

<http://bit.ly/1z4hRak>


GUNZ
dental

P 1800 025 300

F 1800 025 202

www.gunz.com.au

CareDent

For your
everyday
flossing
needs!



Denta Floss Nylon Regular & Mint

Caredent Denta Floss:

ORDER CODES:

CAR4010 Caredent Denta Floss Nylon Wax 50m x6

CAR5888 Caredent Denta Floss Nylon Mint 50m x6

Every day nylon floss coated with patented ptfе and wax, helps prevent decay and maintains healthy gums. Vibrant blue coloured pack, easy to carry and easy to use. Great value makes this a great starting point for those seeking an alternative to the more expensive brands of floss and tape.

Gunz Dental also stock the following flossing products from the Johnson & Johnson preventative care range:



DPA31036 REACH Cleanpaste Flosser Starter pack x 1

DPA31037 REACH Cleanpaste Daily Flosser Refill pack x 28

DPA35701 Dentotape 20 metres

DPA35700 Floss, Waxed- 200m x 6

DPA35702 Floss Waxed Mint 6 x 50 metres

DPA35704 Johnson & Johnson Floss Waxed 6 x 50 metres



P 1800 025 300

F 1800 025 202

www.gunz.com.au

ORTHODONTIC CARE:

A range of products to help prevent tooth decay and gum disease by keeping orthodontic appliances as well as teeth and gums clean.



P 1800 025 300
F 1800 025 202
www.gunz.com.au



SUNSTAR



Features unique V-cut bristle design



Gum Orthodontic V Trim Toothbrush:

ORDER CODE: GUM124 Orthodontic V Trim toothbrush. Pkt 6

The GUM® Orthodontic toothbrush has been developed for safe and effective brushing of teeth fitted with orthodontic appliances including braces, brackets, tubes, and wires.

The unique V-cut bristle design perfectly adapts to the anatomy of the teeth with braces and gently cleans in and around the orthodontic appliances. This special configuration of bristles reaches into the crevices between the gums and the brackets, as well as below the brackets, under the wires, and between the brackets themselves—the places that your dental professional is most concerned about. The GUM® Orthodontic toothbrush keeps braces as well as teeth and gum clean, helping to prevent dental decay and gum diseases. The sleek handle offers excellent adaptation to the size of the hand and adds comfort and control while brushing. It helps to position the GUM® Orthodontic toothbrush at the correct 45-degree angle for optimal and easy cleaning.



V-Trim cut

- ADAPTS AROUND ORTHO BRACKETS AND WIRES - The unique V-cut head helps bristles adapt around brackets and wires.
- CLEANS TEETH AND GUMS - Designed to clean the orthodontic appliance as well as the teeth and gums, removing dental plaque from the hard to reach areas where dental decay and gum disease may occur
- HANDLE DESIGN PROVIDES CONTROL and POSITIONING - Handle enables the correct 45-degree angle and thumb position while giving comfort and control while brushing.
- LONG THIN NECK - Enables thorough cleaning even in the more difficult to reach areas at the back of the mouth



P 1800 025 300

F 1800 025 202

www.gunz.com.au

SUNSTAR




Floss for orthodontic comfort



Gum ORTHO+Floss:

ORDER CODE: GUM3220 ORTHO+Floss PK 6

3 in 1 floss. Built in stiff plastic threader with a spongy brush to sweep out food caught under archwires and around brackets. Soft nylon floss section for regular interdental plaque removal. Easy to dispense and cut with new thumb control wheel incorporated into the dispenser.



**SEE THE NEW ADA PATIENT INFORMATION
VIDEO - GET IT FOR YOUR WAITING ROOM:**

<http://bit.ly/1BtsF2w>



P 1800 025 300

F 1800 025 202

www.gunz.com.au

SUNSTAR



Wax for orthodontic comfort



Gum Orthodontic Wax:

GUM723 Orthodontic Wax, Unflavoured x 1

GUM724 Orthodontic Wax, Mint Flavoured x 1

GUM® Orthodontic Wax has been formulated to relieve discomfort by adhering to orthodontic appliances such as braces and wires, and to the auxiliary hooks, tubes, and ligatures. GUM® Orthodontic Wax is also suitable for patients wearing partial dentures where wire clasps may be causing discomfort, sore spots, or abrasions. The wax is easily applied.



Caredent Orthodontic Wax:

CAR6620 Caredent Ortho Wax Mint x50

CAR6625 Caredent Ortho Wax Reg x50

Braces can cause sores and irritation of the soft tissues of the mouth. Caredent Ortho Wax helps relieve wire and bracket irritation. Wax comes in a convenient carry case. A natural wax product. Easy to use. Cut off the required amount of wax, roll it into a ball and apply it on the bracket or brace which is irritating the soft tissue.



P 1800 025 300

F 1800 025 202

www.gunz.com.au

INTERDENTAL CLEANING:

A great range of products designed to fit all kinds of interdental spaces. Featuring unique technologies that are flexible for easy access to ensure reduced bacterial contamination of the gingiva.



P 1800 025 300
F 1800 025 202
www.gunz.com.au



SEE GUM SOFT-PICKS IN ACTION:

<http://bit.ly/1LdwNsW>



NEW! Gum Soft-Picks:


ORDER CODES:

GUM632M40 Soft-Picks Regular with Travel Case. Pkt 40

GUM634M40 Soft-Picks Large with Travel Case. Pkt 40

Flexible and soft rubber bristles with fluoride slide between teeth to safely remove plaque and food particles. Wire free technology ensures that it does not damage the gingiva.


CLINICALLY PROVEN EFFICACY
Reduces gingivitis and removes plaque between teeth as effectively as floss. 7



P 1800 025 300

F 1800 025 202

www.gunz.com.au



SEE GUM TRAV-LER IN ACTION:

<http://bit.ly/1z4q8Lf>



Gum Trav-ler Interdental:

ORDER CODES:

GUM1312 Trav-ler 0.6 mm (6) (6brushes-6caps)

GUM1314 Trav-ler 0.8 mm (6) (6brushes-6caps)

GUM1412 Trav-ler 0.9 mm (6) (6brushes-6caps)

GUM1414 Trav-ler 1.1 mm (6) (6brushes-6caps)

GUM1512 Trav-ler 1.2 mm (4) (4brushes-4caps)

GUM1514 Trav-ler 1.3 mm (4) (4brushes-4caps)

GUM1612 Trav-ler 1.4 mm (6) (6brushes-6caps)

GUM1614 Trav-ler 1.6 mm (6) (6brushes-6caps)

GUM1618 Trav-ler 2.0 mm (4) (4brushes-4caps)

UCHX-treated bristles (unique patented technology)

inhibit bacterial growth. Different sized brushes to fit all kinds of interdental spaces. Flexible neck for easy access and added protection due to nylon coated wire. Cap can be used as handle for easier posterior reach.



ISO standardization for manual interdental brushes states that it should be measured by the minimum diameter of a hole which the interdental brush can pass through without stem deformation



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Gum Proxabrush Classic:

ORDER CODES:


- GUM605 Proxabrush Classic Plastic Handle w 612 & 614
- GUM412 Proxabrush Refill Brushes Ultra Fine, Cylindrical, Pkt 8
- GUM414 Proxabrush Refill Brushes Ultra Fine, Tapered, Pkt 8
- GUM512 Proxabrush Refill Brushes Extra Fine, Cylindrical, Pkt 8
- GUM514 Proxabrush Refill Brushes Extra Fine, Tapered, Pkt 8
- GUM612 Proxabrush Refill Brushes Fine, Cylindrical, Pkt 8
- GUM614 Proxabrush Refill Brushes Fine, Tapered, Pkt 8

Traditional handle and refill system for all kinds of interdental spaces.

Gum Stimulator:

ORDER CODE: GUM600 Stimulator Handle and Tip x1

The first and original stimulator for gingival massage and interproximal cleaning.



SEE THE GUM STIMULATOR TIP IN ACTION: <http://bit.ly/1LdvOsR>



P 1800 025 300

F 1800 025 202

www.gunz.com.au

CareDent

Easy grip
for easy
plaque
removal



Caredent ProximAll colour range

ProximAll Pro Brush:

ORDER CODES:

CAR6704	ProximAll Brush 1XXXXF Pnk 4pkx12
CAR6706	ProximAll Brush 2 XXXF Ora 4pkx12
CAR6708	ProximAll Brush 3 XXF Red 4pkx12
CAR6710	ProximAll Brush 4 XF Blue 4pkx12
CAR6712	ProximAll Brush 5 Fine Yel 4pkx12
CAR6714	ProximAll Brush 6 Med Grn 4pkx12
CAR6724	ProximAll Pro Brush 1 XXXXF Pink
CAR6726	ProximAll Pro Brush 2 XXXF Oran50
CAR6728	ProximAll Pro Brush 3 XXF Red50
CAR6730	ProximAll Pro Brush 4 XF Bluex50
CAR6732	ProximAll Pro Brush 5 Fine Yellx50
CAR6734	ProximAll Pro Brush 6 Med Greenx50
CAR6736	ProximAll Pro Brush 7 Large Purpx50

Interdental brush made with quality German plastic coated wire and Tynex Dupont bristles. A rubber thumb grooved handle for a comfortable grip. Cap extends the handle for additional grip, offering additional benefits for those with dexterity issues.

Caredent Interbrush:

ORDER CODES:

CAR7115	InterBrush Sm Profx10 Box36 Blue
CAR7120	InterBrush Med Profx10 Box36 Orange
CAR7125	InterBrush Lg Profx10 Box36 Purple
CAR7130MIX	InterBrush Mix Prof 36pks(20)
CAR7200	InterBrush Sm #30 x 6 Blue
CAR7210	Interbrush Med #30 Box6 Orange
CAR7220	InterBrush Lg #30 Box6 Purple

A wire free patented interdental brush from Caredent.




GUNZ
dental

P 1800 025 300

F 1800 025 202

www.gunz.com.au



Good access
for implants
and braces

TePe Interdental Brushes:

ORDER CODES:

ODP84271	Interdent Brush x25 T2271 Pink	ODP84656	Interdent Brush 8x10 0.4 Pink
ODP84273	Interdent Brush x25 T2273 Red	ODP84658	Interdent Brush 8x10 0.5 Red
ODP84274	Interdent Brush x25 T2274 Blue	ODP84659	Interdent Brush 8x10 0.6 Blue
ODP84275	Interdent Brush x25 T2275 Yel	ODP84700	Angle Interdent Brush x6 Pink
ODP84276	Interdent Brush x25 T2276 Grn	ODP84701	Angle Interdent Brush x6 Orange
ODP84277	Interdent Brush x25 T2277 Prp	ODP84702	Angle Interdent Brush x6 Red
ODP84643	Interdent Brush 8x10 0.7 Yellow	ODP84703	Angle Interdent Brush x6 Blue
ODP84644	Interdent Brush 8x10 0.8 Green	ODP84704	Angle Interdent Brush x6 Yellow
ODP84645	Interdent Brush 8x10 1.1 Purple	ODP84705	Angle Interdent Brush x6 Green
ODP84646	Interdent Brush 8x10 1.3 Grey	ODP84709	Angle Interdent Brush x6 Assort
ODP84647	Interdent Brush 8x10 1.5 Blk		

The interdental brushes can also be used for cleaning around implants or braces. TePe Interdental Brushes Original and Extra soft have a straight neck, while TePe Angle has a pre-angled neck for extra good access between the back teeth and from the inside. Plastic coating on the wire for safe and gentle cleaning.



Curaprox Prime Series:

ORDER CODES:

CURCPS011	Curaprox 1/Br Prime 1.1mm 4
CURCPS06	Curaprox 1/Br Prime 0.6mm 5
CURCPS07	Curaprox 1/Br Prime 0.7mm 5
CURCPS08	Curaprox 1/Br Prime 0.8mm 5
CURCPS09	Curaprox 1/Br Prime 0.9mm 5

Made with patented CURAL surgical wire offering extreme durability with the thinnest available wire cores. Curaprox Prime respects the papilla (which floss cannot) and is equally suited to young and elderly patients.



P 1800 025 300

F 1800 025 202

www.gunz.com.au



ORAL THERAPEUTICS:

Mouth rinses, gels, toothpastes and dissolvables – for the prevention and remedy of diseases of the teeth and oral cavity.



P 1800 025 300
F 1800 025 202
www.gunz.com.au

Alcohol free daily mouth rinses



Curasept ADS205



Curasept ADS920



Curasept ADS220



Curasept ADS212

Curasept Mouth Rinse Range:

ORDER CODES:

- CURADS205 Curasept Chlorhex+ M/Rinse Fl 200ml
- CURADS920 Curasept Chlorhex M/Rinse.20% 900ml
- CURADS220 Curasept Chlorhex M/Rinse.20% 200ml
- CURADS212 Curasept Chlorhex M/Rinse.12% 200ml

All Curasept chlorhexidine rinses are alcohol free, non-staining and they will not affect taste perception. They offer excellent protection and improved efficacy through better patient compliance.

6 Steps to making Dental Hygiene Dept. a Profit Centre: No longer a loss leader, dental hygiene treatments can add to your bottom line without working extra hours.

READ MORE:

<http://bit.ly/15FNA6U>

Bleaching and Caries Control in Elderly Patients: The use of 10% carbamide peroxide applied nightly can help elderly patients with oral hygiene care.

READ MORE:
<http://bit.ly/15Lqzj5>



Curasept ADS705



Curasept ADS712



Curaprox Enzykal

CURADEN + SWISS
 LEADING ORAL CARE PEOPLE

Excellent everyday toothpastes/ gels

Curaden Toothpaste/Gel:

ORDER CODE:

CURADS705 Curasept TPaste Gel .05% Ch/F 75ml

Curasept Toothpaste with 0.05% sodium fluoride and 0.05% chlorhexidine for daily use. Especially suited to patients with orthodontic appliances or bands but recommended to all as an excellent every day toothpaste. An ideal gel/toothpaste containing NO sodium lauryl sulphate ("SLS"). Recommended where the patient is using chlorhexidine mouthwashes.

ORDER CODE:

CURADS712 Curasept TPaste Gel .12%ChI 75ml

Curasept gel/toothpaste with 0.12% Chlorhexidine Diguconate. An ideal gel/toothpaste containing NO sodium lauryl sulphate ("SLS"). Recommended where the patient is using chlorhexidine mouthwashes.

ORDER CODE:

CURENZY Curaprox Enzykal + FIT/Paste 75ml

Enzykal is the new toothpaste from Curaprox. It contains no sodium lauryl sulphate ("SLS") but it is high in fluoride at 950ppm. It does not contain chlorhexidine, but it is completely compatible with the Curasept chlorhexidine mouth rinses.



P 1800 025 300

F 1800 025 202

www.gunz.com.au



Why not consider giving you patients a protective fluoride treatment at the next Hygiene appointment?

SUGGESTED ORAL CARE TREATMENT

REGIME IN ACTION: <http://bit.ly/1z4vSVn>



Curaden Gingival Gel:

ORDER CODE:

CURADS350

Curasept Chlorhex Gin Gel .50% 30ml

Curasept 0.50% chlorhexidine gel in a 30ml tube. The gel is smooth flowing and easy to use. It is, of course, non-staining and alcohol free. Popular with Periodontists and general dentists alike. Chlorhexidine is the gold standard in chemical plaque control. For application in the spaces between your teeth, apply the gel by applying to interdental brush and insert through the space. No need to rinse afterwards.

Germiphene Gel:

ORDER CODE: WN999035

Germiphene Gel 7 Neut Fluor 450ml

2% sodium fluoride gel ideal for use with any patient. Gluten Free/Dye Free Fluoride uptake within 4 minutes. No sour or bitter aftertaste due to the neutral pH. Excellent post whitening desensitiser. Ideal for root caries, dry mouth, decay, recurrent decay and orthodontic decalcification.



P 1800 025 300

F 1800 025 202

www.gunz.com.au

CareDent

Fast acting relief for dry mouth & ulcers



Caredent Rapids:

ORDER CODES:

CAR6845 Caredent Rapids for Dry Mouth and Ulcers.

Contains 12 foil packs of each.

CareDent's new mouth strip technology offers fast acting relief to those that suffer from the common problem of a dry mouth and mouth ulcer. Rapids use the latest dissolvable technology and provide instant relief. The dispenser kit is ideal for the surgery counter or waiting room allowing easy storage as well as encouraging patients requesting more details on the products. The side windows will be of a great help to all practice managers when working on the stock count.



P 1800 025 300

F 1800 025 202

www.gunz.com.au