


Safest and highest quality
safety system designed
and constructed to
eliminate leaks!

Neonatal / Pediatric Enteral Safety Feeding System



FDA approved

NeoDevices Enteral-only Syringe

- BD Syringe
- Includes syringe label

ND-S3EL	3cc	100/box	ND-S20EL	20cc	100/box
ND-S5EL	5 cc	100/box	ND-S30EL	30cc	50/box
ND-S10EL	10cc	100/box	ND-S60EL	60cc	40/box

NeoDevices Polyurethane Feeding Tubes Enteral Only - Safety System

- Designed for long-term OG/NG use
- Lipid Resistant
- Includes a variety of hub connections

FTS5.0P-EO	5.0 Fr	40cm	FTL6.5P-EO	6.5 Fr	90cm
FTM5.0P-EO	5.0 Fr	60cm	FTS8.0P-EO	8.0 Fr	50cm
FTL6.5P-EO	5.0 Fr	90cm	FTM8.0P-EO	8.0 Fr	60cm
FTS6.5P-EO	6.5 Fr	40cm	FTL8.0P-EO	8.0 Fr	90cm
FTM6.5P-EO	6.5 Fr	60cm	FTL10.0P-EO	10.0 Fr	42cm

NeoDevices Enteral Extension Sets with Orange Stripe Tubing

- DEHP Free - All of our products contain non- DEHP plasticizers.
- Labeling - Each set has a visual reminder label "For Enteral Feeding Only"
- Lipid Resistant
- Priming Volume Approx 2.0ml.

60" Orange Stripe Tubing -
Differentiates Enteral Tubing
from IV Tubing

ND60EO	Set with Non-IV Enteral Stepped Adapter and Enteral Only Connector, 50/box
ND60EOY	Bifurcated Set with Non-IV Stepped Adapter and Enteral Only Connector with 2" Orange Stripe Tubing with Enteral Only Connector, 50/box

NeoDevices Silicone Feeding Tubes Enteral Only - Safety System

- Designed for long-term OG/NG use
- Lipid Resistant
- Includes a variety of hub connections

FTS5.0S-EO	5.0 Fr	40cm	FTL6.5S-EO	6.5 Fr	90cm
FTM5.0S-EO	5.0 Fr	60cm	FTS8.0S-EO	8.0 Fr	40cm
FTL5.0S-EO	5.0 Fr	90cm	FTM8.0S-EO	8.0 Fr	60cm
FTS6.5S-EO	6.5 Fr	40cm	FTL8.0S-EO	8.0 Fr	90cm
FTM6.5S-EO	6.5 Fr	60cm			

FDA approved

