

Salter Labs®
Quality Care Begins With Quality Products
Salter Family of Oxygen Cannulas



In our continuing commitment to manufacture the finest products the marketplace has to offer, Salter Labs® provides the most complete range of cannulas.

Our cannula line includes an extensive selection of oxygen cannulas ranging in sizes from Premature to Adult, covering flow rates from 3 LPM up to 15 LPM. Also available, is a series of demand cannulas for use with dual port oxygen conserving devices, available in sizes from Pediatric to Adult (not pictured.)

Salter Labs® offers four focus Salter-Style® cannulas:

- **Salter-Style® Micro Cannulas*** for flows up to 3 LPM
- **Salter-Style® Original Cannulas*** for flows up to 6 LPM
- **Salter-Style® Quiet Cannulas*** for flows up to 8 LPM
A quieter version for higher flows
- **Salter-Style® High Flow Cannulas*** for flows up to 15 LPM

In addition, our cannulas are available with a wide range of tubing lengths and an assortment of end fittings.

Many of the Salter-Style® cannulas are also available in the TLCannula™ version, which incorporates pre-installed soft foam ear cushions on the headset tubing to reduce over the ear irritation and discomfort.

1610 - Premature
Up to 3 LPM

1611 - Neonate
Up to 3 LPM

1601 - Infant
Up to 3 LPM

1615 - Intermediate Infant
Up to 3 LPM

1602 - Pediatric/Small Adult
Up to 6 LPM









1616 - Micro, Adult
Up to 3 LPM

1600 - Original, Adult
Up to 6 LPM

1600Q - Quiet, Adult
Up to 8 LPM

1600HF - High Flow, Adult
Up to 15 LPM

1056 - Conventional Adult
Up to 6 LPM



100 W. Sycamore Road - Arvin, CA 93203
 Phone: 661-854-3166
 1-800-421-0024
 1-800-235-4203
 Fax: 661-854-3850
 Toll Free Fax: 1-800-628-4690
 Fax: 661-854-3850
 E-mail: salterlabs@us.salterlabs.com
 www.salterlabs.com



ISO 13485

VITAL CARE PRODUCTS, INC.
 800 837 5511 + sales@vitalcareproducts.com



slc-05