




Premature infants have special skin care needs  
 Casmed products are designed to accommodate those unique needs

## Neonatal Intensive Care CASMED Monitoring Supplies



### Klear-Trace® Electrodes

- Proprietary hydrogel adhesive provides excellent skin contact for a high-quality signal
- Easy application & removal without irritating babies' skin
- Radiolucency lets electrodes stay in place during x-ray
- Typically maintains adhesive quality for 3-7 days
- 100% tested
- 7/8" x 1 1/2"



4800	Circle	300/case	(3/pouch)
4810	Rectangle	300/case	(3/pouch)
2000	Rectangle	300/case	(3/pouch)
2500	Rectangle	250/case	(5/pouch)
2500-22	Rectangle	250/case	(5/pouch)

### NeoGuard® Temperature Probes

- Single-patient-use probe designed for length of stay improves cost recovery
- Low profile design and thin, pliable wire enhance flexibility and patient comfort
- Each probe packaged with a NeoGuard Reflector
- Packaged 40/case



4501	Replaces Air Shields	Probe #6820970
4503	Replaces Ohmeda	Probe #LA003

### NeoGuard® Reflectors

- Holds NeoGuard Temp Probe securely to skin while shielding it from radiant heat
- Hydrogel adhesive allows probe to be removed and repositioned without irritating babies' skin

4200	1.5" Heart	200/case
4220	1" Heart	100/case

### NeoGuard® Limboard® Arm Boards


- Soft foam construction for comfortable, effective IV stabilization
- Proprietary hydrogel strips hold Imboard in place without harming the skin
- Can be reshaped multiple times as IV sties are repositioned



4400	1" x 4" with gel strips	48/case
4401	1.5" x 4.5" with gel strips	48/case
4400B	1" x 4" w/o gel strips	48/case
4401B	1.5" x 4.5" w/o gel strips	48/case

### NeoGuard® Hydrogel Adhesive Tape

- 25/box (2/pouch)
- |       |               |
|-------|---------------|
| 4400S | 0.5" x 5.125" |
| 4401S | 0.635" x 6.5" |



VITAL CARE PRODUCTS, INC.  
 3046 Brecksville Road + Suite 2 + Richfield, Ohio 44286  
 800 837 5511 + www.vitalcareproducts.com

Critical Care Solutions Delivered Full Circle

