



KaVo QUATTROcare Plus

Protect Your Investment

- Easy to use
- Safe controlled cleaning
- Automatic maintenance for up to 4 handpieces
- Chuck maintenance adapter



4 Easy Steps

- 

1
Position the instruments
- 

2
Close the door
- 

3
Press Start and everything is DONE.
- 

4
Weekly chuck maintenance

QUATTROcare plus Spray

- For cleaning and care of instruments / air motors in the KaVo QUATTROcare
- Designed specifically to withstand sterilization
- Relubrication after sterilization is unnecessary, eliminating the risk of recontaminating sterilized handpieces

Go to WWW.TRYKAVO.COM to start your free and easy trial today!



KaVo. Dental Excellence.

888-ASK-KAVO · www.kavousa.com

KaVo QUATTROcare Plus

The user friendly QUATTROcare Plus system is designed to clean, lubricate and purge your instruments with the simple push of a button. Loading and operation is easy as well as the automatic care and maintenance procedure. Using the QUATTROcare Plus system means less wear, less repairs, less staff down time, smoother operation, and increased longevity of your handpieces. The QUATTROcare Plus can be used for virtually all manufacturer's high-speeds, low-speed contra-angles, straight handpieces, air driven scalers, and air motors. Complete handpiece maintenance is available with the new chuck care feature. Once per week, cleaning and care of handpiece chucks, maintains their gripping force, reliability and long service life.

QUATTROcare Plus & Accessories

| Product | Description | Part No. |
|--------------------------------------|---|-----------|
| QUATTROcare Plus | Automatic handpiece maintenance and purging system | 1008.3805 |
| Multiflex adapter | Adapter for KaVo Cellular Optic High-Speeds, SONICflex Air Scaler, INTRAFlex Low-Speeds | 0411.7991 |
| E-Type Adapter | Adapter for INTRAMatic Low-Speed and surgical attachments | 0411.7981 |
| Adapter for Comfortdrive | Attaches Comfortdrive to QUATTROcare Plus | 1005.1707 |
| INTRA Head adapter | Attaches INTRA Head to QUATTROcare Plus | 0411.7941 |
| Bumper sheath for INTRA head adapter | | 0411.5532 |
| Fuse for QUATTROcare Plus | | 1008.6112 |
| QUATTROcare Plus Spray | For cleaning and care of instruments / air motors in the QUATTROcare (6 cans) | 1005.4524 |



QUATTROcare Plus Spray

- Dissolves inorganic residues
- Propellout removes dissolved inorganic residues out of the instruments
- Lubricates inner parts of the instruments
- **Must be used exclusively in the QUATTROcare Plus**



QUATTROcare Plus Spray expands within your handpiece to both clean and lubricate.



When you purchase the KaVo QUATTROcare unit, you will receive up to 2.5 year Warranty on all your handpieces

Go to WWW.TRYKAVO.COM to start your free and easy trial today!

KaVo Dental
11727 Fruehauf Drive
Charlotte, NC 28273
888.ASK.KAVO (888.275.5286)
Fax 1.704.927.0816
www.kavousa.com

©2012 KaVo Dental Corporation. 906.8988 - rev 1 - 1.12



KaVo. Dental Excellence.
888-ASK-KAVO · www.kavousa.com